



TiLite ZR ZRA - Order Form



Dealer Name*:

Account #*:

City, State, Zip:

Phone:

Email:

ATP:

Rep Name:

P.O. Number:

Quote Number:

Client Reference*:

Client Weight:

*Required field

Special considerations that need to be addressed here (e.g., diagnosis):

Frame Style

Z2FS2	ZRA Titanium Frame Full adjustability. 265lb weight limit. 150lb weight limit for seat sizes 10"x10" to 12"x12". For rear seat width and seat depth $B \geq 15"$ use K0005. For rear seat width and seat depth $< 15"$ use E1235.	K0005, E1235	\$ 4,500
Z2FS1	ZR Titanium Frame Ultimate TiFIT Minimal adjustability. 265lb weight limit. 150lb weight limit for seat sizes 10"x10" to 12"x12". For rear seat width and seat depth $B \geq 15"$ use K0005. For rear seat width and seat depth $< 15"$ use E1235.		\$ 750

Frame Upgrades

ZHDR1	Heavy Duty Frame 300lb weight limit. No charge when 20" rear seat width and 20" seat depth are selected. Includes billet aluminum front wheels. Requires tension adjustable by straps back upholstery, Shadow, Spinergy SPOX or LX rear wheels. Only available with 0° or 2° camber. Not available with Slipstream or Glide forks, carbon fiber camber tube , velcro adjustable upholstery, swing away tubular armrests or swing away flip back tubular armrests.		\$ 485
-------	--	--	---------------

Transit Options

ZTTDO3	Transit Tie-Downs Not available with ZR. Tested to WC-19. Consists of four tie-down brackets, two pelvic belt mounting brackets and reinforced frame. Please refer to the ZRA Transit Tie-Down Option Brochure for important information regarding transit tie-down options.	E1023	\$ 363.00
ZWTB	Transport Brackets Consists of four tie-down brackets. Unoccupied use only. Please refer to the ZR Transport Bracket Supplement for important information regarding transit options.	E1022	\$ 379.00

SmartDrive Power Assist System

ZMAXSD1	Permobil SmartDrive MX2+ Not available with rear seat heights less than 13.5" or 20" rear wheels. Maximum user weight capacity of 331 lbs.	E0986	\$ 7,529.00
ZMAXSD18	Switch Control Buttons		No charge
ZMAXSD8	SwitchControl Button with Mono Jack Option		No charge
ZMAXSD19	Buddy Button For use with SwitchControl button with mono jack option.		No charge
ZMAXSD4	Carrying Bag		\$ 84.00
ZMAXSD10	Universal Padded Handle		\$ 22.00

Frame Finishes

ZFF1	Satin Finish Available with titanium frames only.	No charge
	ZSCBT1 Buffing Pads for Satin Frame	\$ 22.00
ZFF	Powder Coated Frame Select below.	No charge
Frame Color:	ZFF8 High Gloss Black ZFF15 Tangerine Metallic ZFF19 Black Lava ZFF57 Indigo Purple ZFF60 Silver Streak ZFF64 Daffodil Yellow ZFF10 Red Pepper ZFF30 Acid Green Metallic ZFF27 Titanium Grey ZFF40 Hyper Blue ZFF63 Lavendar Dusk ZFF29 Neon Pink ZFF17 Midnight Blue Pearl ZFF59 Pomegranate Red ZFF41 Flame Orange ZFF58 Black Galaxy ZFF14 Ocean Blue Metallic ZFF18 White Pearl ZFF62 Blue Mist ZFF42 OD Green ZFF61 Seafoam Green	
ZFF3	Tattoos Tattoos on front frame tubes. Titanium frame finished in satin. Tattoos will be trimmed and adjusted to fit your frame.	\$ 375
	ZFF65 Liquid ZFF66 Bloom ZFF67 Rebel ZFF68 Vector	



Liquid



Bloom



Rebel



Vector

Color Anodize Package

Standard:	Black	Standard
Anodized Package:	Gold Green Purple Orange Blue Red	\$ 71.00

Includes caster mount, caster fork, caster caps and ZR camber tube bracket. If no anodize colour is selected, default will be black.

TiFit: Tailored Individualized Frame

Measurements shown below are our most popular selections. Not seeing what you want? Contact customer service or view our Creative Mobility Division (CMD) catalog [here](#) to view additional options. Please refer to Reference Diagrams at end of form for more information.

Rear Seat Width 'A' | Measure from outside of seat tube at back post to the outside of the opposite seat tube at back post.

ZSW1:	10"	11"	12"	13"	14"	15"	16"	17"	18"	19"		
ZSW2:	20"										E2201	\$ 536.00

HCPCS code only applies to 20" rear seat width. >20" sizes available through CMD.

Front Seat Width (Seat Taper) 'T' | Measure from inside of front frame tube to inside of opposite front frame tube. Required with V-front end, not with standard front end.

ZFSW1	No Taper 2.5" narrower than rear seat width.	Standard
ZFSW3	1.5" Taper 4" narrower than rear seat width. Includes tension adjustable by straps seat upholstery. Not with 10" and 11" seat widths. Bolt-on adjustable sling available through our Creative Mobility Division #20.	\$ 304.00

Seat Depth 'B' | Measure from front of back post to the front edge of seat sling. Frame bend will start 0.5" in front of seat sling unless Custom Frame Depth is specified or the combination of specifications affect the rolling dynamics of the chair.

ZSD1:	10"	11"	12"	13"	14"	15"	16"	17"	18"	19"		
ZSD2:	20"										E2203	\$ 568.00
ZSD3 Custom Frame Depth:	Specify:											No charge

Measure from front of seat sling to beginning of frame bend. Ranges from 0.5"-2" and can be given as specific as 0.25" increments. A custom frame depth outside the specified range may be available through Creative Mobility Division. A custom frame depth of 20" will have an additional charge.

Ergonomic Seat 'U' | Measure from front of back post to beginning of upward bend of seat tube.

ZERG1	ZR Ergonomic Seat Comes standard with tension-adjustable by strap seat upholstery. Bolt-on adjustable sling available through our Creative Mobility Division #20. Solid seat available through our Creative Mobility Division CMD #16.	\$ 369.00
	5" 6" 7" 8" 9"	
ZERG2	ZRA Ergonomic Seat Center of gravity adjustment may be affected by this option. Not with V-front or seat taper, seat back that locks when folded, desk armrests, depth adjustable backrests, fold down rigid sideguards or fendered adjustable sideguards. Bolt-on adjustable sling available through our Creative Mobility Division #20. Solid seat available through our Creative Mobility Division CMD #16. Minimum seat depth 10" for 6". Minimum seat depth 12" for 7". Minimum seat depth 13" for 8".	\$ 433.00
	6" 7" 8"	

Front Seat Height 'C' | Measure from floor to top of seat tube at front edge of seat sling.

ZFSH1:	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	21"
---------------	-----	-------	-----	-------	-----	-------	-----	-------	-----	-------	-----

ZRSH1:	12.5"	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"
	18.5"	19"	19.5"	20"								

Rear Wheel	26"	15.5"-20"
	25"	14.5"-20"
	24"	14.5"-20"
	22"	14"-19"
	20"	12.5"-18.5"

Front Angle or Overall Frame Length | *Select one. If Overall Frame Length is selected, front angle will be bent to achieve requested length.*

ZFA1	Frame Angle 'L' <i>Measure from frame front to floor. 85° not available with 6" front wheels. 90° not available with 5" or 6" front wheels.</i>					Standard
	70°	75°	80°	85°	90°	

Footrest

Seat to Footrest 'E':	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"
-----------------------	-------	-----	-------	-----	-------	-----	-------	-----	-------	-----	-------	-----

Measure from front edge of seat sling to top rear of footrest. Measurement must be at least 2" shorter than Front Seat Height 'C'.

ZFTR3	Angle Adjustable Footrest		Standard
ZFTR1	Titanium Footrest w/Flat ABS Cover		\$ 217.00
ZFTR2	Titanium Open Loop Footrest		\$ 277.00
ZFTR12	Angle Adjustable Flip-Back Footrest	Not available with 10" or 11" rear seat width or heavy duty frame.	K0040 \$ 244.00

Seat to Footrest 'E': 7"-9" 9"-11" 11"-13" Set at

Measure from front edge of seat sling to top rear of footrest.

ZFTR4	Angle Adjustable High Mount Footrest	K0040	\$ 109.00
-------	--------------------------------------	-------	-----------

Footrest Width | Select standard front end 'I' if no seat taper 'T' has been selected. Select V-front end 'V' if a seat taper 'T' has been selected.

ZFE1 **No V-Front Taper 'I'** Find your rear seat width 'A' in the chart below and select one of the options to the right. Measure from inside of front frame tube to inside of opposite front frame tube. Standard

Rear Seat Width 'A'	10"	7.5"		
	11"	8.5"		
	12"	9.5"	8"	
	13"	10.5"	9"	
	14"	11.5"	10"	
	15"	12.5"	11"	9"
	16"	13.5"	12"	10"
	17"	14.5"	13"	11"
	18"	15.5"	14"	12"
	19"	16.5"	15"	13"
20"	17.5"	16"	14"	

No Taper

1.5" Taper

3.5" Taper

Footrest Width 'I'

ZFE2	V-Front Taper 'V' <i>Measure from inside of front frame tube to inside of opposite front frame tube 2.5" above footrest. Titanium footrests have 1" of adjustment up and down. Not with ergonomic seat frames.</i>	\$ 315.00
------	---	-----------

9" 10" 11" 12" 13" 14" 15"

Back Supports

ROHO Back Supports

SBRS14	Agility Active Aluminum, Min Contour with Foam	<i>Includes fixed hardware. Select shell height.</i>	E2620	\$ 738.00
SBRS15	Agility Active Aluminum, Mid Contour with Foam	<i>Includes fixed hardware. Select shell height.</i>	E2620	\$ 738.00
	8" 10"			
SBRS7	Agility Carbon Min Contour Back with Foam/Air Insert	<i>Includes fixed hardware. Select shell height.</i>	E2620	\$ 738.00
ACB1	Agility Carbon Min Upgrade	<i>Required with Agility Carbon Backs. Not available with bolt-on push handles.</i>		\$ 464.00
SBRS9	Agility Carbon Mid Contour Back with Foam/Air Insert	<i>Includes fixed hardware. Select shell height.</i>	E2620	\$ 738.00
ACB2	Agility Carbon Mid Upgrade	<i>Required with Agility Carbon Backs. Not available with bolt-on push handles.</i>		\$ 464.00
	10" 13" 16"			

Comfort Back Supports

SBRS22	Comfort Embrace Back	<i>Replaces Comfort Acta-Embrace Back. Includes 2-point quick release Compass hardware. Select shell height.</i>	E2613	\$ 700
	10" 12" 14" 16"			
	Comfort-Tek Cover	Diamond Spacer Cover		
	Privacy Shield			\$ 52
SBRS23	Comfort Relief Back	<i>Replaces Comfort Acta-Relief Back. Includes 2-point quick release Compass hardware. Select shell height.</i>	E2615	\$ 886
	10" 12" 14" 16"			
	Comfort-Tek Cover	Diamond Spacer Cover		
	Privacy Shield			\$ 52

Back Upholstery

ZBUP2	Velcro Adjustable	<i>Not available with fold-down push handles.</i>		Standard
ZBUP3	NEW Tension Adjustable By Straps, Padded Nylon	<i>Required for fold-down push handles.</i>	E2611	\$ 369.00
ZBUP5	Tension Adjustable By Straps, AIR Breathable	<i>Required for fold-down push handles.</i>	E2611	\$ 369.00
ZUPC3	Ultrasuede Backrest Stripe, Black	<i>Not available with AIR breathable upholstery.</i>		\$ 115.00
ZBUP0	Omit Back Upholstery			No charge

Backrest

Back Type

ZSB1	Folding Aluminum Adjustable Height Back	<i>Required for heavy duty frame.</i>		Standard
ZSB2	Folding Titanium Adjustable Height Back, Rounded Rigidizer Bar	<i>Titanium back will match frame finish.</i>		\$ 277.00
ZSB15	Non-Folding Fixed Height Titanium Back, Rounded Rigidizer Bar	<i>Not available with ZRA, heavy duty frame, transport brackets or push handles. Titanium back will match frame finish.</i>		\$ 315.00

Back Options

ZSB11	Backrest Locks When Folded	<i>Available on all backs. Not available with ergonomic seat.</i>		Standard
ZSB12	4" Deep Backrest Rigidizer Bar	<i>Not a rounded style rigidizer bar. Required for heavy duty frame. Not available with non-folding fixed height titanium back type.</i>		No charge
ZSB16	Depth Adjustable Back Frame	<i>Seat depths of 10"-18" allow 2" of rearward adjustment. Height adjustable push handles have a limited height range. Not with 19" or 20" seat depths, heavy duty frame, tension adjustable by straps seat upholstery or carbon fiber side guards.</i>		No charge
ZSB3	TiShaft Back Release Bar, Aluminum	<i>Not available with non-folding fixed height titanium backs.</i>		Standard
ZSB5	TiShaft Back Release Bar, 2x2 Twill Weave Carbon Fiber	<i>Not available with non-folding fixed height titanium backs.</i>		\$ 82.00

Push Handles

ZSB6	Integrated Push Handles			No charge
ZSB7	Integrated Fold-Down Push Handles	<i>Requires tension adjustable by straps back upholstery.</i>		\$ 125.00
ZSB8	Bolt-On Push Handles	<i>Includes 4" deep backrest rigidizer bar. Not available with titanium back.</i>		\$ 217.00
ZSB9	Bolt-On Height-Adjustable Push Handles	<i>Includes 4" deep backrest rigidizer bar.</i>		\$ 363.00

Folding Back Height 'F' | ZSB1/ZSB2. Set at 0.5" increments.

Push Handles	None / Bolt-On	9.5"-11.5"	12"-16.5"	16"-20.5"	Set at:
	Integrated	11"-12.5"	13.5"-16.5"	17.5"-20.5"	Set at:
	Fold-Down	13.5"-15"	14.5"-18.5"	18.5"-21"	Set at:
		Short	Medium	Tall	Specify
Folding Back Height					

Fixed Back Height 'F'

Select height:	8"	8.5"	9"	9.5"	10"	10.5"	11"	11.5"	12"	12.5"	13"	13.5"
	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"
	20"											

Required for fixed side guards. Not with ZRA, swing-away tubular or swing-away flip-back tubular armrests, integrated push handles, integrated fold-down push handles or bolt-on height adjustable push handles.

Back Angle 'G' | ZSBA1 Measure angle from front of back post to floor.

Folding Back:	80°	83°	86°	89°	92°	95°	98°	101°				
Fixed Back:	80°	81°	82°	83°	84°	85°	86°	87°	88°	89°	90°	91°
	93°	94°	95°	96°	97°	98°	99°	100°	101°			92°

Rear Wheel Positioning

Center of Gravity 'H' | Measure from front of back post at the seat tube to the center of the rear axle. If no center of gravity is selected the chair will ship with the center of gravity in the most rearward position. Ergo seat selection may limit available center of gravity. 0° not available with swing-away push handles.

ZCOG1:	0"	0.25"	0.5"	0.75"	1"	1.25"	1.5"	1.75"	2"	2.25"	2.5"	2.75"
	3"	3.25"	3.5"	3.75"	4"	4.25"	4.5"	4.75"	5"			

Rear Wheel Spacing 'J' | Measure from outside of back post to inside of rear tire. 0.75" rear spacing not available with 0° camber, and not available with 2° camber for ZRA. Ergonomic seat selection may limit available center of gravity.

ZRWS1	Fixed Rear Wheel Spacing											Standard
		0.75"	1"	1.25"								
ZCBR12	Adjustable Rear Wheel Spacing Camber Plugs											No charge
ZCBR9	Additional Spacers to Increase Rear Wheel Spacing											No charge
ZCBR11	Zero Play Aluminum Camber Plugs	Not available with power adaptable frame or more than 6° camber.										No charge

Camber 'K'

ZCBR1:	0°	2°	4°	6°								E2601
	More than 0° camber not with power adaptable frame. 8° camber may be available through our Creative Mobility Division.											
ZCBR15	Aluminum Camber Tube											Standard
ZCBR3	Titanium Camber Tube											\$ 217.00
ZCBR4	Carbon Fiber Camber Tube	Only available on ZR and only with 0°, 2°, 4° zero play plugs. Not available with anti-tips, heavy duty frame.										#N/A

Seat Upholstery

ZSUP1	Tension Adjustable Bolt-On Upholstery, Black	Not available with ergonomic seat or tapered seat.										Standard
ZSUP2	Tension Adjustable By Straps, Black	Standard with ergonomic seat and seat taper of 1"-3".										\$ 201.00
ZSUP8	Aluminum Solid Seat	Ergonomic seat requires CMD #16. Not with 10" or 11" rear seat width.										E2231 \$ 309.00
ZSUP0	Omit Seat Upholstery											No charge
CMD018	Custom iSleeve	Please provide sleeve width and depth (inches) in Notes section at the end of this order form. Add 1" to device depth and 0.25" to width. Best when used with tension adjustable bolt-on upholstery. Some limitations may apply.										CMD

Seat Cushions

ZCSH2	Comfort Curve Cushion	Stretch-Air cover.										E2601 \$ 122.00
ZCSH23	Comfort Acta-Embrace Comfort-Tek Zero Elevation Cushion											E2607 \$ 428.00
ZCSH11	ROHO Single Compartment	Fits seat widths 14"-19".										E2622 \$ 534.00
	Low Profile	Mid Profile	High Profile									
ZCSH14	ROHO QUADTRO SELECT Cushion	Fits seat widths 14"-19".										E2624 \$ 534.00
	Low Profile	Mid Profile	High Profile									
ZCSH16	ROHO Hybrid Select Cushion	Fits seat widths 14"-18", and 20"										E2624 \$ 571.00
ZCSH17	Vicar Vector O2 with Comfair Cover	Fits seat widths 14"-19". Non-stock upcharges may apply.										E2624 \$ 542.00

Front Wheels

ZFW24	3"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire	\$ 206.00
	Silver Hub Black Hub	
ZFW14	4"x0.75" Plastic Wheel w/Poly Tire <i>Not available with Slipstream fork.</i>	No charge
ZFW3	4"x0.75" Light-Up Micro Caster <i>Not available with Slipstream fork.</i>	\$ 88.00
ZFW26	4"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire <i>Not available with power assist or heavy duty frame.</i>	No charge
ZFW27	4"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire	\$ 217.00
	Silver Hub Black Hub	
ZFW6	5"x1" Plastic Wheel w/Poly Tire <i>Not available with 90° front angle or Slipstream fork.</i>	No charge
ZFW28	5"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire <i>Not available with 90° front angle.</i>	No charge
ZFW30	5"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire <i>Not available with 90° front angle.</i>	\$ 233.00
	Silver Hub Black Hub	
ZFW9	6"x1" Plastic Wheel w/Poly Tire <i>Not available with 85° or 90° front angle or Slipstream fork.</i>	No charge
ZFW32	6"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire <i>Not available with 85° or 90° front angle or Slipstream fork.</i>	No charge
ZFW33	6"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire <i>Not available with 85° or 90° front angle.</i>	\$ 244.00
	Silver Hub Black Hub	

Precision Lock Front Forks

ZFK1	TiLite Standard Fork	Standard
ZFK2	TiLite Slipstream Single-Sided Fork <i>Not available with heavy duty frame.</i>	\$ 185.00
ZGL1	Glide Suspension Fork <i>Black anodize only. Not available with heavy duty frame. Select rider weight.</i>	E1015 \$ 526.00
	125lbs 125-180lbs 180-265lbs	
ZTFS1	Titanium Fork Stem <i>Comes standard on ZR.</i>	\$ 66.00

Rear Wheels

See rear seat height ranges chart for rear seat height 'D' limitations.

ZRW14	TiLite SHADOW Rear Wheel <i>Select size. 25" and 26" wheels require 15" minimum seat depth.</i>	Standard
	20" (451) 22" (501) 24" (540) 25" (559) 26" (590)	
ZRW7	Spinergy SPOX Rear Wheel <i>Select size and spoke color. 25" and 26" wheels require 15" minimum seat depth. Orange, pink, purple and green spokes only available in 20" size.</i>	\$ 1,195
	20" (451) 22" (501) 24" (540) 25" (559) 26" (590)	
	Black Red Blue White Orange Pink Purple Green	
ZRW8	Spinergy LX Rear Wheel <i>Select size and spoke color. 25" and 26" wheels require 15" minimum seat depth.</i>	\$ 1,295
	22" (501) 24" (540) 25" (559) 26" (590)	
	Black Red Blue White Orange Pink Purple Green	
ZRW9	Spinergy Design Lab LX Rear Wheel <i>Black spokes only. Select size and anodizing color. If same color is selected, anodized hubs and rims match colour anodized package. 25" and 26" wheels require 15" minimum seat depth.</i>	\$ 1,995
	22" (501) 24" (540) 25" (559) 26" (590)	
	Gold Green Blue Red Purple Orange	
ZRW16	Spinergy Carbon Blade Rear Wheel <i>2x2 twill weave carbon fiber. Not available with heavy duty frame. 25" wheels require 15" minimum seat depth.</i>	\$ 1,795
	24" (540) 25" (559)	
ZRW17	X-Core Mag Rear Wheel <i>Select size. Not available with heavy duty frame, short tabs, quad release axles, spoke guards or omit handrims options. Requires aluminum silver anodized handrims.</i>	\$ 225
	20" (451) 22" (501) 24" (540)	

Axles

ZAXL1	Stainless Quick Release Axle	Standard
ZAXL2	Titanium Quick Release Axle	\$ 277.00
ZAXL3	Quad Release Axle <i>Not available with Mag or 26" (590) Spinergy Spox or LX rear wheels.</i>	\$ 169.00

Pneumatic Rear Tires

ZRTR1	Treaded Tire 22"(501), 24"(540), 26"(590)		Standard
ZRTR2	Airless Inserts <i>Not available with Spinergy Carbon Blade rear wheels.</i>	E2213	\$ 98.00
ZRTR31	MBL SpeedLite Tire 20"(451), 22"(501), 24"(540), 25"(559), 26"(590)		No charge
ZRTR10	MBL Gazelle Tire 24"(540), 25"(559) <i>Kevlar bead, lightest weight option.</i>		\$ 99
ZRTR30	Primo Sentinel Tire 24"(540), 25"(559), 26"(590)		\$ 136.00
ZRTR8	Schwalbe Marathon Plus Evolution Tire 24"(540), 25"(559), 26"(590)		\$ 169
ZRTR32	Primo C-764 Knobby Tire 24"(540)		\$ 104.00

Solid Rear Tires

ZRTR22	Primo Orion Tire 20"(451), 22"(501), 24"(540), 25"(559) <i>Not available with Spinergy Carbon Blade rear wheels.</i>		\$ 120.00
ZRTR23	Primo Active Tire 24"(540), 25"(559) <i>Not available with Spinergy Carbon Blade rear wheels.</i>	E2213	\$ 147.00

Handrims

ZRIM1	Aluminum Handrim, Silver Anodized		Standard
ZRIM21	Aluminum Handrim, Black Anodized		No charge
ZRIM12	Titanium Handrim, Satin Finish <i>Not available with 20"(451) or 22"(501) rear wheels or 24"(540) knobby tires.</i>		\$ 509.00
ZRIM20	Q-Grip Handrim <i>Exclusive high-friction, anti-slip coating. Not available with 20"(451) rear wheels.</i>		\$ 358.00
ZRIM5	Plastic Coated Handrim, Black		\$ 174.00
ZRIM10	The Surge Handrim <i>Large profile. All-In-One oval with Gription strip.</i>		\$ 433.00
ZRIM11	The Surge LT Handrim <i>Small profile. All-In-One oval with Gription strip.</i>		\$ 466.00
ZRIM14	The Natural-Fit Handrim, Standard Grip <i>Original, large profile. Not available with short tabs or Spinergy Carbon Blade rear wheels.</i>		\$ 428.00
ZRIM15	The Natural-Fit Handrim, Super Grip <i>Original, large profile. Not available with short tabs or Spinergy Carbon Blade rear wheels.</i>		\$ 466.00
ZRIM16	The Natural-Fit Handrim, Without Thumb Grip <i>Original, large profile.</i>		\$ 390.00
ZRIM17	The Natural-Fit LT Handrim, Standard Grip <i>Reduced weight, small profile. Not available with short tabs or Spinergy Carbon Blade rear wheels.</i>		\$ 440.00
ZRIM18	The Natural-Fit LT Handrim, Super Grip <i>Reduced weight, small profile. Not available with short tabs or Spinergy Carbon Blade rear wheels.</i>		\$ 514.00
ZRIM19	The Natural-Fit LT Handrim, Without Thumb Grip <i>Reduced weight, small profile.</i>		\$ 423.00
ZRIM22	Omit Handrims		No charge

Handrim Tab Length

ZRIM7	Long Tabs		Standard
ZRIM8	Short Tabs <i>Not available with Mag rear wheels, The Natural-Fit handrims with thumb grip, or Natural-Fit handrims on Shadow wheels.</i>		\$ 28.00

Wheel Locks

ZWLK7	Compact Composite Push to Lock <i>Available with 0.75"-1.5" rear wheel spacing only.</i>		No charge
ZWLK8	Compact Composite Pull to Lock <i>Available with 0.75"-1.5" rear wheel spacing only.</i>		No charge
ZWLK10	Extension Handles for Compact Composite Push/Pull to Lock	E0961	\$ 93.00
ZWLK11	MBL Composite Scissor Lock		\$ 158.00
ZWLK12	MBL Performance Aluminum Scissor Lock		\$ 175
ZWLK3	Uni-Lock		\$ 179.00
	Push Pull		
ZWLK4	Extension Handles for Uni-Lock	E0961	\$ 93.00
ZWLK15	Grade Aids Push to Lock <i>Not available with solid rear tires.</i>	E0974	\$ 190.00

Side Guards

TSDG13	Side Guards Custom Cut to Profile of Rear Wheels Available for aluminum and carbon non-fendered side guards.		\$ 342.00
CMD007	Custom Offset Profiled Side Guards (Removable or Fold-Down) Please specify desired offset to nearest 0.5" and side guards type (rigid or fold-down) in the Notes section of this order form. Side guard material must be selected. Some limitations may apply.	CMD	
CMD	Fixed Side Guards CMD33 (fendered carbon fiber) or CMD73 (flat panel aluminum or carbon fiber) available. ZRA also requires CMD19 (fixed back). Please visit our CMD catalogue, pricing will vary.	CMD	

Aluminum

ZSDG2	Aluminum Side Guards, Rigid Removable	E0956, E1034	\$ 217.00
	Small Standard		
ZSDG14	Aluminum LT Side Guards, Rigid Removable	E0956, E1034	\$ 217
	Small Standard		
ZSDG3	Aluminum Side Guards, Quad Removable w/Hole	E0956, E1034	\$ 233.00
	Small Standard		
ZSDG15	Aluminum LT Side Guards, Quad Removable w/Hole	E0956, E1034	\$ 233
	Small Standard		
ZSDG4	Aluminum Fold Down Side Guards, Rigid Not available with ergonomic seats.	E0956, E1034	\$ 363.00
	Small Standard		
ZSDG6	Aluminum Fendered Side Guards, Adjustable 24"-25" rear wheels require rear seat height 14.5"-19.5". 26" rear wheels require rear seat height 15.5"-18.5".	E0956, E1034	\$ 363.00

2x2 Twill Weave Carbon Fiber

ZSDG7	Carbon Side Guards, Rigid Removable	E0956, E1034	\$ 336.00
	Small Standard		
ZSDG16	Carbon LT Side Guards, Rigid Removable	E0956, E1034	\$ 336
	Small Standard		
ZSDG8	Carbon Side Guards, Quad Removable w/Hole	E0956, E1034	\$ 347.00
	Small Standard		
ZSDG17	Carbon LT Side Guards, Quad Removable w/Hole	E0956, E1034	\$ 347
	Small Standard		
ZSDG9	Carbon Fold Down Side Guards, Rigid Not available with ergonomic seats.	E0956, E1034	\$ 499.00
	Small Standard		
ZSDG12	Carbon Fendered Side Guards, Removable	E0956, E1034	\$ 493.00

Armrests

ZARM1	Swing Away Tubular Armrests Select range. Not available with non-folding fixed height titanium back or height adjustable push handles.	E0973	\$ 244.00
	6.5"-9.5" 9.5"-12.5"		
ZARM6	Swing Away Fixed Height Flip Back Tubular Armrests Not available with heavy duty frame.		No charge
ZARM4	Swing Away Flip Back Tubular Armrests, 8.5"-11.5" Not available with non-folding fixed height titanium back.	E0973	\$ 342.00
ZARM2	Removable Height Adjustable Desk Arm with Rigid Side Guard Select range and pad size. Requires 14" minimum seat depth. Not available with ergonomic seat.	E0973	\$ 363.00
	9"-12" 12"-15"		
	Desk Length (10") Full Length (14")		
ZARM3	Removable Height Adjustable Desk Arm with Rigid Side Guard and Transfer Loop Select range and pad size. Requires 15" minimum seat depth. Not available with ergonomic seat.	E0973	\$ 406.00
	9"-12" 12"-15"		
	Desk Length (10") Full Length (14")		
ZARM5	Armrests Covered in Black Ultrasuede Not available with desk arms.		\$ 120.00

Positioning Belts

ZBLT1	Velcro Adjustable Belt, Black Nylon Webbing	E0978	\$ 61.00
ZBLT2	Bodypoint 2" Belt, Push-Button Buckle (Auto Style)	E0978	\$ 103.00
ZBLT4	Bodypoint 1" Padded Belt, Push-Button Buckle (Auto Style) <i>Not available with 19" or 20" rear seat widths.</i>	E0978	\$ 196.00
ZBLT5	Bodypoint 1.5" Padded Belt, Push-Button Buckle (Auto Style) <i>Not available with 19" or 20" rear seat widths.</i>	E0978	\$ 185.00
ZBLT8	Q'Straint WC-19 Wheelchair-Anchored Pelvic Belt <i>Recommended for use with Transit Tie Down. Not available with ZR.</i>	E0978	\$ 217.00
ZBLT10	Bodypoint Optional Band Clamp Mounting Hardware	K0108	\$ 22.00
ZBLT16	Bodypoint EvoFlex Pelvic Stabilizer, Push Button	E0978	\$ 201.00

Calf Straps

ZCLF2	Velcro Adjustable	K0038	\$ 44.00
	Small (9"-12" footrest widths) Medium (11"-15" footrest widths) Large (14"-17" footrest widths) Xtra Large (16"-19" footrest widths)		
ZCLF3	Bodypoint Padded Velcro Adjustable	K0038	\$ 109.00
	Small (8"-11" footrest widths) Medium (11"-15" footrest widths) Large (16"-19" footrest widths)		
ZCLF1	Buckle Adjustable, Black Nylon Webbing	K0038	\$ 34.00

Residual Limb Supports

ZRSL1	Residual Limb Support, 8"W x 10"L, Single	E1020	\$ 428.00
	Left Right		
ZRSL2	Residual Limb Support, 8"W x 10"L, Pair	E1020	\$ 855.00

Anti-Tips

ZTIP1	Rear Aluminum, Black	E0971	\$ 174.00
ZTIP2	User Friendly Flip-Up	E0971	\$ 363.00

Accessories

ZGRD1	Clear Lexan Spoke Guards <i>Not available with Mag rear wheels.</i>	K0065	\$ 147.00
ZBCK3	Permobil Medical Necessities Backpack		\$ 103.00
ZPCH3	Permobil Medical Necessities Seat Pouch		\$ 66.00
ZPCH4	Bodypoint® Mobility Bag		\$ 109.00
ZUEST1	Upper Extremity Support Tray <i>Requires armrests.</i>	E0950	\$ 315.00
	12"-14" seat widths 13"-15" seat widths		
ZLUG1	Luggage Carriers		\$ 158.00
ZMPCT2	Black Neoprene Impact Guards with Ultrasuede Stripe		\$ 115.00
ZCH1	Crutch Holder	E2207	\$ 244.00
ZAFP1	4" Military Service Patch		\$ 115.00
	Army Navy Air Force Marines Coast Guard		

General

Permobil manufactures TiLite manual wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. HCPCS Codes are for reference purposes only. Actual coding and coverage varies by payer. Proper HCPCS coding is not a guarantee of payment. Please verify the appropriate payer's coverage policy. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

Notes

TiFit Reference Diagrams



TiLite ZRA



TiLite ZR

