



Order form

# F3 Corpus - AADL W688



Dealer Name\*:

Account #\*:

City, Province, Postal Code\*:

Phone:

Email:

ATP:

Rep Name:

P.O. Number:

Quote Number:

Client Reference\*:

Client Weight\*:

*\*Required field*

Special considerations that need to be addressed here (e.g., diagnosis):

## Chair Model and Wheels

Front-wheel drive F3 power base includes: Corpus seat frame, ComfortRide suspension, 4-point tie-down hardware, VoltPro charger, and R-net 120 amp electronics.

### Choose Base

11331	<b>F3 Base Corpus - F3 NPO</b> <i>FWD F3 power base set up for no power seat functions, max speed: 6 mph, weight capacity: 300 lb.</i>	\$8,274.00
11332	<b>F3 Base Corpus - F3 SPO</b> <i>FWD F3 power base set up for one power seat function, max speed: 6 mph, weight capacity: 300 lb.</i>	\$8,274.00
108398-99-0	<b>F3 Base Corpus - F3 MPO</b> <i>FWD F3 power base set up for multiple power seat functions, max speed: 6 mph, weight capacity: 300 lb.</i>	\$8,274.00
108852-99-0	<b>Permobil Connect</b>	No charge

### Support Wheel Options | Choose one

107330-99-0	<b>Support Wheels for F3</b>	No charge
11313	<b>Omit Front Support Wheels</b> <i>User Weight must not exceed 220lbs. Max speed will be limited to 5mph on the F3 Corpus and 6mph on the F5 Corpus. Additional restrictions may apply based on chair configuration.</i>	No charge

## Chair Model and Wheels

### Choose Drive Tires

108414-99-0	<b>Flat Free Tire 3"W x 14"T Black - Standard Tread</b> <i>Includes flat free standard tread caster wheels.</i>	No charge
108411-99-0	<b>Pneumatic Tire 3"W x 14"T Black - Standard Tread</b> <i>Includes flat free standard tread caster wheels.</i>	No charge
11333	<b>Flat Free Tire 3"W x 14"T Black - Aggressive Tread</b> <i>Includes flat free aggressive tread caster wheels and mud guards.</i>	\$596.00
108412-99-0	<b>Pneumatic/Snow Tire 3"W x 14"T Black - Aggressive Tread (Client Cost)</b> <i>Not recommended for year-long use and must be ordered along with standard drive wheels. Does not include caster wheels.</i>	\$596.00
11334	<b>Mud Guard Kit F3</b>	\$67.00

## Base Options

### Choose Base Color | Choose both shroud and rim color

<b>Standard Shroud Colors</b>	Sky	Midnight	No charge
<b>Standard Rim Colors</b>	Sky	Midnight	No charge

### Batteries | Choose one

108999-99-0	<b>Permobil Batteries, Grp 24 AGM GEL - Pair (AADL Contract Price)</b>	\$1,064.00
109183-99-0	<b>Permobil Batteries, Grp 34 AGM GEL - Pair (AADL Contract Price)</b>	\$1,052.00

### Optional Items

ESP UNIT	<b>Enhanced Steering Perform Unit (ESP), R-net</b>	\$1,349.00
109170-99-0	<b>Essentials Box with Reflectors (Client Cost)</b>	\$609.00
109169-99-0	<b>Essentials Box with Lights (Client Cost)</b> <i>Must also select LED Light Kit.</i>	\$1,014.00
109379-99-0	<b>LED Light Kit for F3 (Client Cost)</b> <i>No charge with 4 power seat function configurations.</i>	\$1,318.00
10715	<b>Compact Ventilator Tray - Corpus Seat</b> <i>Custom option. 13"T, adjustable width between 8" - 13", and adjustable depth between 2" - 6.5". Additional drive and seating restrictions may apply.</i>	\$1,910.00
10787	<b>Battery Tray for Laptop Ventilator Tray</b> <i>Custom option. 11 3/8" W x 6 1/2" H x 1 5/8" D.</i>	\$306.00
10788	<b>O2 Holder for Laptop Ventilator Tray</b> <i>Custom option. Mounts to left or right side of compact laptop ventilator tray. 7"T, adjustable width between 3.5" - 7", and adjustable depth between 3.75" - 5".</i>	\$432.00
109620-99-0	<b>Programming Key for R-net - Permobil Version</b>	\$596.00

## Joystick Options

Choose joystick type, handle, and mount. NOTE: For additional input devices, please refer to the Alternative Drive Systems section.

### Joystick Type

107779-99-0	<b>Color Joystick Module (CJSM)</b> <i>High-resolution 3.5" VGA color display.</i>	\$1,517.00
11186	<b>Permobil Joystick Module (PJSM), NC</b> <i>2.8" color display.</i>	No charge
105700-99-0	<b>R-net Remote Color Joystick w/ Lights, Heavy Duty Gimbal</b> <i>Color LCD joystick, includes 3" knob with 3/4" diameter. Not compatible with alternative joystick handles.</i>	\$2,319.00
107829-99-0	<b>3.5mm Y-Splitter Cable for PJSM/CJSM</b> <i>Allows 2 mechanical switches to control up to 4 separate functions through assignable buttons programming.</i>	\$48.00

## Joystick Options

### Joystick Handle

10045	<b>Cone Shaped Handle - Standard</b>	No charge
100960-99-0	<b>Chin Cup Handle</b>	\$136.00
100961-99-0	<b>T Handle</b>	\$200.00
100962-99-0	<b>Large Ball</b>	\$136.00
100963-99-0	<b>Softball</b>	\$182.00
100964-99-0	<b>Mushroom Handle</b>	\$182.00
100965-99-0	<b>Stick Handle</b>	\$182.00
105442-99-0	<b>1.75" Black Foam Ball</b>	\$136.00
105846-99-0	<b>Bodypoint 3" U-Shaped Handle w/ Flex-Shaft</b>	\$175.00
105847-99-0	<b>Bodypoint 4" U-Shaped Handle w/ Flex-Shaft</b>	\$175.00
105848-99-0	<b>Bodypoint Rubber Dome Handle</b>	\$175.00

### Joystick Mount

<b>Fixed</b>	Mount on Left	Mount on Right	No charge
<b>Slimline Retractable</b>	Mount on Left	Mount on Right	\$514.00
11232	<b>Adjustable Height Panel Bracket</b> <i>Provides 1.77" of vertical adjustment and 25° of angle adjustment.</i>		\$398.00
105210-99-0	<b>Mount Standard Joystick on Corpus Back Hinge</b> <i>Can be chosen instead of attendant control.</i>		\$356.00
	Mount on Left	Mount on Right	
108268-99-0	<b>Bodypoint Tri-Lock Midline Mount Assembly</b> <i>Must select specific joystick attachment separately.</i>		\$785.00
	Mount on Left	Mount on Right	
108280-99-0	<b>Bodypoint Tri-Lock Mount for R-net Joystick</b>		\$322.00
108279-99-0	<b>Bodypoint Tri-Lock Mount for Compact Joystick</b>		\$322.00
10842	<b>PERMOFix Mount for mo-vis/Compact Joystick</b> <i>PermoFix mount for attaching Compact and mo-vis joysticks to Tri-Lock midline mount.</i>		\$322.00

### Other Options

10411	<b>Expandable Controller - R-net</b> <i>Required when ordering power tilt, power recline, and power elevating legrests together OR when ordering an alternative drive control.</i>	No charge
10513	<b>Harness for Expandable Controller</b> <i>Must be ordered with any expandable controller.</i>	No charge
10582	<b>Output Module (IOM) for R-net Control System (Client Cost)</b> <i>Enables user to control an external device (e.g., communication device, door opener, etc.) through the drive input (e.g., JSM, head array, etc.). May require specialty cable or wiring to be ordered separately.</i>	\$1,986.00

## Power Seating

*Please Note: Power seat adjustability may be limited to prevent contact with base and/or backrest-mounted items.*

**Adjustable Seat Height** | *If ordering power adjustable seat height, elevator travel is used to set seat heights greater than 17.5".*

11110	<b>ActiveHeight Power Adjustable Seat Height - 12" Travel (Client Cost)</b> <i>Must also select power tilt when selecting this option.</i>	\$3,588.00
11117	<b>Seat Tube Fixed</b>	No charge

## Power Seating

Seat-to-Floor Height	17.5"	18.5"	19.5"	No charge
----------------------	-------	-------	-------	-----------

### Posterior Tilt/Recline Packages | Choose one

11211	Corpus Seat w/ Manual Recline	\$2,751.00
11131	Corpus Seat w/ No Tilt & 180° Power Recline (AADL W993 for Power Tilt)	\$10,404.00
11132	Corpus Seat w/ 50° Tilt (AADL W993 for Power Tilt)	\$9,329.00
11133	Corpus Seat w/ 50° Tilt & 180° Power Recline (AADL W993 for Power Tilt)	\$15,370.00

**ActiveReach Packages** | Chest support belt is strongly recommended. Power legrest recommended for full functionality. Seating and drive capabilities may be limited in certain configurations.

11141	5° ActiveReach	No charge
107090-99-0	10° ActiveReach <i>ActiveHeight and power tilt required.</i>	No charge
107091-99-0	20° ActiveReach (Client Cost) <i>ActiveHeight, power tilt/recline and Corpus VS elevating legrests required.</i>	\$3,359.00
11100	30° ActiveReach w/Transfer (Client Cost) <i>ActiveHeight, power tilt/recline and Corpus VS elevating legrests required.</i>	\$3,777.00

### Seat Function Control Options | ICS switchbox availability may be limited to certain configurations.

10418	Seat Function Control Switch for R-net <i>Does not include operation of power seat functions through drive control.</i>	No charge
10419	Single Seat Function Control Kit for R-net <i>Allows operation of a single power seat function through drive control.</i>	\$1,887.00
10585	Multiple Seat Function Control Kit for R-net <i>Allows operation of a two or more power seat functions through drive control.</i>	\$3,144.00
102798-99-0	Push Buttons on ICS Switchbox	No charge
102799-99-0	Toggle Switches on ICS Switchbox	\$336.00

### Switchbox Mounting

10626	Mount ICS Switchbox onto R-net Joystick Bracket	No charge
104651-99-0	Mount ICS Switchbox onto Armrest Mount on Left      Mount on Right	No charge
106256-99-0	Mount ICS Switchbox on Back (Client Cost) <i>No charge when ordering Omni2 Control Module.</i> Mount on Left      Mount on Right	\$189.00
105409-99-0	Do Not Install ICS Switchbox, Send as Part	No charge

### Power Seating Options

104196-99-0	ICS Alternative Switchbox for R-net Electronics <i>Accommodates up to 8 momentary switches.</i>	\$1,484.00
10750	Memory Seat Program, Corpus w/ R-net	\$138.00
106209-99-0	Legrest Movement Linked to Recline Electronically (Client Cost)	\$56.00
10898	Enable Independent Repositioning Mode (IRM) <i>IRM allows for sequential movement of power seat functions into a saved posterior position and returns to a saved upright position.</i>	No charge

## Backrest Options

Installed backrest heights are measured from the seat pan to the top of the backrest shell with the backrest set at 90°. Accessible seat depth will be dependent upon backrest selection.

### Corpus Backrest Kits

105067-99-0	<b>Sliding Backrest Kit - 5" Mechanical Shear</b>	No charge
105066-99-0	<b>Fixed Backrest Kit - No Mechanical Shear</b>	\$61.00

**Corpus Backrests** | *Corpus Ergo Backs include lateral wedges and foam lumbar pads. Solid Backs do not include foam lumbar pads.*

Solid Corpus Ergo	20" Tall <i>Power Recline will be limited to 150°</i>				23" - 28" Tall <i>Height adjustable in 1" increments</i>				
	16" W	18" W	20" W		16 x 23-28	18 x 23-28	20 x 23-28		
16" W	16 x 20								No charge
18" W	18 x 20								No charge
20" W	20 x 20								No charge
<b>Set Back Height:</b>	20"	23"	24"	25"	26"	27"	28"		

### Corpus Backrest Options

11335	<b>Corpus Ergo Back - Leatherette</b>	No charge
11336	<b>Corpus Ergo Back - Stretch-Air</b>	No charge
10725	<b>Split Backrest Cushion - Corpus</b> <i>Custom option.</i>	\$404.00
10786	<b>Tailored Backrest Cushion - Corpus</b> <i>Custom option. Sides of cushion are trimmed to provide additional clearance for lateral trunk supports.</i>	\$126.00
109029-99-0	<b>BodiLink Premium 3" x 5" Lateral Supports - Pair (Seating Benefits)</b> <i>Must also select hardware.</i>	\$470.00
109030-99-0	<b>BodiLink Premium 4" x 6" Lateral Supports - Pair (Seating Benefits)</b> <i>Must also select hardware.</i>	\$470.00
109031-99-0	<b>Swing-Away Adj BodiLink Lateral Hardware Corpus - Pair (Seating Benefits)</b>	\$725.00

**ROHO Backrests** | *Includes backrest shell, cushion, air inserts, UT rail for headrest and Direct Mount adapter. (Seating Benefits)*

ROHO AGILITY <i>Max Contour</i>	16" Tall <i>Installed height between 22" - 25"</i>				18" Tall <i>Installed height between 22" - 26"</i>				20" Tall <i>Installed height between 24" - 27"</i>									
	14" W	16" W	18" W	20" W	14 x 16	16 x 16	18 x 16	20 x 16	14 x 18	16 x 18	18 x 18	20 x 18		14 x 20	16 x 20	18 x 20	20 x 20	
14" W	14 x 16																	\$1,226.00
16" W	16 x 16																	\$1,226.00
18" W	18 x 16																	\$1,226.00
20" W	20 x 16																	\$1,226.00
<b>Set Back Height:</b>	22"	23"	24"	25"	26"	27"												

## Backrest Options

**Comfort Backrests** | Includes backrest shell, cushion and Compass Power Mount.

User weight may be limited based upon the size of the backrest selected. ([Seating Benefits](#))

	<b>14" Tall</b> <i>Installed height between 17" – 23"</i>	<b>16" Tall</b> <i>Installed height between 19" – 25"</i>	<b>18" Tall</b> <i>Installed height between 21" – 27"</i>	<b>20" Tall</b> <i>Installed height between 23" – 29"</i>	
<b>Acta-Back</b>					
<i>Stretch-Air</i>					
14" W		14 x 16	14 x 18		\$896.00
16" W		16 x 16	16 x 18	16 x 20	\$896.00
18" W		18 x 16	18 x 18	18 x 20	\$896.00
20" W			20 x 18	20 x 20	\$896.00
<b>Acta-Embrace</b>					
<i>Comfort-Tek</i>					
12" W	12 x 14				\$1,242.00
14" W		14 x 16	14 x 18		\$1,242.00

### Comfort Backrest Options | ([Seating Benefits](#))

11714	<b>BodiLink Premium 3" x 5" Lateral Supports - Pair</b> <i>Must also select hardware.</i>	\$470.00
11715	<b>BodiLink Premium 4" x 6" Lateral Supports - Pair</b> <i>Must also select hardware.</i>	\$470.00
109903-99-0	<b>BodiLink GT2 Swing-Away Lateral Hardware - Pair</b> <i>Must select with BodiLink lateral trunk supports on Comfort backrests.</i>	\$725.00

### Aftermarket Backrest Options

105782-99-0	<b>Tubular Backrest Frame - Corpus (1" Canes)</b> <i>Provides 1" diameter, vertical canes for mounting aftermarket backrests. Frame can be adjusted to support backrest widths from 14" to 20". Up to several inches of seat depth may be lost when mounting this system.</i>	\$542.00
106072-99-0	<b>Direct Backrest Frame - Corpus</b> <i>Used to mount backrests ranging in width from 14" to 20". Some seat depth may be lost when mounting this system.</i>	\$542.00

## Seat Sizing

Please select both width and depth options below.

### Seat Pan Width | Choose one

<b>Set at:</b>	17" W	19" W	21" W	23" W	No charge
----------------	-------	-------	-------	-------	-----------

### Seat Pan Depth | Choose one

<b>Kit 1, Set at:</b>	14"	15"	16"	No charge
<b>Kit 2, Set at:</b>	17"	18"	19"	No charge
<b>Kit 3, Set at:</b>	20"	21"	22"	No charge

## Seat Cushion Options

**Corpus Ergo Seat Cushions** | *Corpus Ergo Air cushion features ROHO Dry Floatation technology.* [\(Seating Benefits\)](#)

	15" D	17" D	19" D	21" D	23" D	
<b>Corpus Ergo</b>						
<i>Leatherette</i>						
17" W	17 x 15	17 x 17	17 x 19	17 x 21		\$687.00
19" W	19 x 15	19 x 17	19 x 19	19 x 21	19 x 23	\$687.00
21" W		21 x 17	21 x 19	21 x 21	21 x 23	\$687.00
23" W			23 x 19	23 x 21	23 x 23	\$687.00
<b>Corpus Ergo</b>						
<i>Stretch-Air</i>						
17" W		17 x 17	17 x 19	17 x 21		\$687.00
19" W		19 x 17	19 x 19	19 x 21	19 x 23	\$687.00
21" W		21 x 17	21 x 19	21 x 21	21 x 23	\$687.00
23" W			23 x 19	23 x 21	23 x 23	\$687.00
<b>Corpus Ergo Air</b>						
<i>C3GRHOSCA</i>						
17" W		17 x 17	17 x 19	17 x 21	17 x 23	\$729.00
19" W		19 x 17	19 x 19	19 x 21		\$729.00
21" W			21 x 19	21 x 21		\$729.00

## ROHO Neoprene Seat Cushions | [\(Seating Benefits\)](#)

	8 cells 14.75" D	9 cells 16.5" D	10 cells 18.25" D	11 cells 20" D	12 cells 21.75" D	13 cells 23.5" D	
<b>Quadro Select</b>							
<i>High Profile</i>							
9 cells / 16.5" W	9 x 8	9 x 9	9 x 10	9 x 11	9 x 12		\$816.00
10 cells / 18.25" W		10 x 9	10 x 10	10 x 11	10 x 12	10 x 13	\$816.00
11 cells / 20" W		11 x 9	11 x 10	11 x 11	11 x 12	11 x 13	\$816.00
12 cells / 21.75" W			12 x 10	12 x 11	12 x 12	12 x 13	\$816.00
13 cells / 23.5" W		13 x 9	13 x 10	13 x 11	13 x 12	13 x 13	\$816.00
<b>Single Compartment</b>							
<i>High Profile, Sensor-Ready (* = non-SR)</i>							
9 cells / 16.5" W	9 x 8	9 x 9	9 x 10	9 x 11	9 x 12*		\$668.00
10 cells / 18.25" W		10 x 9	10 x 10	10 x 11	10 x 12	10 x 13*	\$668.00
11 cells / 20" W		11 x 9	11 x 10	11 x 11	11 x 12	11 x 13*	\$668.00
12 cells / 21.75" W			12 x 10	12 x 11	12 x 12	12 x 13*	\$668.00
13 cells / 23.5" W			13 x 10	13 x 11			\$668.00

## ROHO Hybrid Seat Cushions | [\(Seating Benefits\)](#)

	14" D	15" D	16" D	17" D	18" D	20" D	
<b>Hybrid Select</b>							
14" W	14 x 14		14 x 16				\$1,106.00
15" W		15 x 15	15 x 16	15 x 17			\$1,106.00
16" W		16 x 15	16 x 16	16 x 17	16 x 18	16 x 20	\$1,106.00
17" W			17 x 16	17 x 17	17 x 18	17 x 20	\$1,106.00
18" W			18 x 16	18 x 17	18 x 18	18 x 20	\$1,106.00
20" W			20 x 16	20 x 17	20 x 18	20 x 20	\$1,106.00

## Seat Cushion Options

**Comfort Seat Cushions** | *User weight may be limited based upon the size of the cushion selected. (Seating Benefits)*

	14" D	16" D	18" D	20" D	22" D	
<b>Embrace Anti-Thrust</b> <i>Comfort-Tek</i>						
12" W	12 x 14					\$508.00
14" W	14 x 14	14 x 14				\$508.00
<b>M2 Anti-Thrust</b> <i>Comfort-Tek</i>						
14" W	14 x 14	14 x 16				\$590.00
16" W		16 x 16	16 x 18	16 x 20		\$590.00
18" W		18 x 16	18 x 18	18 x 20	18 x 22	\$590.00
20" W		20 x 16	20 x 18	20 x 20	20 x 22	\$590.00
22" W			22 x 18	22 x 20	22 x 22	\$762.00

## Cushion Options | (Seating Benefits)

107505-99-0	<b>ROHO Smart Check</b>	\$313.00
11113	<b>Incontinence Cover for Ergo Cushion</b> <i>Sent as a part.</i>	\$219.00
11173	<b>Incontinence Cover for ROHO Cushion</b> <i>Sent as a part.</i>	\$153.00
105118-99-0	<b>Omit Seat Cushion</b>	No charge

## Armrest Options

*Corpus armrests are infinitely adjustable from 8" - 15", measured from the seat pan to the top of the armrest pad.*

**Armrest Set-Up** | *Width must be wider than backrest to ensure proper fit.*

106543-99-0	<b>Narrow Armbar Kit (13"W - 19"W)</b> <i>For 13" width, armrests will be mounted to inside of bracket. Set Width Between Armrests to:</i> 13"    15"    17"    19"	No charge
106544-99-0	<b>Wide Armbar Kit (19"W - 23"W)</b> <i>Set Width Between Armrests to:</i> 19"    21"    23"	No charge
105069-99-0	<b>Height Adjustable Armrests - NC</b> <i>No charge when ordering power tilt and/or recline. Set Armrest Height to:</i> 8"    9"    10"    11"    12"    13"    14"    15"	No charge
105071-99-0	<b>Medial Swivel</b> <i>Medial swivel cannot be enabled when selecting fixed option. Must select fixed option when ordering arm trough adapters.</i> Left Fixed    Left Enabled    Right Fixed    Right Enabled	No charge

## Armrest Options

Choose Armrest Pads | *Must select left and right. Priced per side.*

CL Pads	<b>Corpus Leatherette</b>	13" Left	16" Left	18" Left	13" Right	16" Right	18" Right	No charge
CG Pads	<b>Corpus Gel (Seating Benefits)</b>	13" Left	16" Left	18" Left	13" Right	16" Right	18" Right	\$229.00
MFS	<b>Moldable Forearm Support (Seating Benefits)</b>	Left	Right					\$378.00
MCHS	<b>Contoured Hand Support - Medium (Seating Benefits)</b> <i>Used with Moldable Forearm Support.</i>	Left Medium	Right Medium					\$188.00
MLHS	<b>Contoured Hand Support - Large (Seating Benefits)</b> <i>Used with Moldable Forearm Support.</i>	Left Large	Right Large					\$199.00
ATA	<b>Arm Trough Adapters (Seating Benefits)</b>	1" Pair	7/8" Pair	Install		Send as Parts		\$394.00

Optional Items | *Priced per side/support. (Seating Benefits)*

<b>Armrest Pouch</b>	Left	Right	\$287.00	
<b>Armrest Supports</b>	Left Outer Right Inner	Left Inner	Right Outer	\$411.00
<b>Elbow Supports (Must also select hardware)</b>	Left Elbow Support	Right Elbow Support	\$221.00	
<b>Elbow Support Hardware</b>	Left	Right	\$322.00	

## Legrest Options

Legrest Type | *Complete column for legrest type chosen*

	<b>Corpus</b> <i>Includes mechanical articulation Footplate range is 6" - 17"</i>	<b>Corpus VS</b> <i>Only available with &gt; or = 20° ActiveReach; Includes power articulation Footplate range is 9" - 17"</i>	
<b>Type</b>	Power <b>(Client Cost)</b> Manual - NC	Power <b>(Client Cost)</b>	\$4,030.00 No charge
<b>Footplate Style</b>	2pc Wide - 9"D x 7.5"W 2pc Narrow - 9"D x 6"W 1pc Narrow - 9.75"D x 12.5"W	2pc Wide - 9"D x 7.5"W 2pc Narrow - 9"D x 6"W 1pc Wide - 9"D x 16.5"W	No charge No charge \$431.00
<b>Footplate Height</b> <i>Height measured from seat pan with legrest at 90°</i>	Set footplate height at:	Set footplate height at:	
<b>Calf Supports</b>			
B Pads	3.5"W x 4.5"T	3.5"W x 4.5"T	No charge
J Pads	5"W x 6.5"T	5"W x 6.5"T	No charge
A Pads	5.5"W x 6"T	5.5"W x 6"T	No charge
H Pads (curved)	7"W x 7.5"T	7"W x 7.5"T	\$446.00
D Pads	8"W x 8"T	8"W x 8"T	\$446.00
<b>Knee Support</b> <i>Both items required for Corpus VS Legrest</i>		Single Post Knee Support/Abductor Adj/Removable Knee Support Hardware	\$238.00 \$347.00

## Legrest Options

### Footplate Options

106530-99-0	<b>Stealth Gel Pads - Narrow Footplates - Pair (Seating Benefits)</b>	\$454.00
106225-99-0	<b>Stealth Gel Pads - Wide Footplates - Pair (Seating Benefits)</b>	\$467.00

## Accessories

### Belt Options | *Must select one of the following*

104518-99-0	<b>Bodypoint Non-Padded Hip Belt - Medium</b>	No charge
104520-99-0	<b>Bodypoint Padded Hip Belt - Medium (Seating Benefits)</b>	\$163.00
104514-99-0	<b>Bodypoint Non-Padded Hip Belt - Large (Seating Benefits)</b>	\$132.00
104516-99-0	<b>Bodypoint Padded Hip Belt - Large (Seating Benefits)</b>	\$224.00
100519-99-0	<b>Positioning Belt Retractable (Client Cost)</b>	\$339.00
100600-99-0	<b>Positioning Belt Retractable Long (Client Cost)</b>	\$368.00
108267-99-0	<b>Bodypoint Padded Hip Belt "Airplane" Latch - Medium (Seating Benefits)</b>	\$287.00
108266-99-0	<b>Bodypoint Padded Hip Belt "Airplane" Latch - Large (Seating Benefits)</b>	\$304.00

### Head Support Options | (Seating Benefits)

107470-99-0	<b>Permobil 7"W Headrest</b> <i>Must also select headrest mounting hardware.</i>	\$310.00
107469-99-0	<b>Permobil 10"W Headrest</b> <i>Must also select headrest mounting hardware.</i>	\$310.00
107468-99-0	<b>Permobil 14"W Headrest</b> <i>Must also select headrest mounting hardware.</i>	\$310.00
107480-99-0	<b>Adj/Removable Headrest Link Hardware - Standard</b>	\$331.00
109178-99-0	<b>Extra Link for Link Headrest Hardware</b> <i>Sent as a part.</i>	\$67.00
108721-99-0	<b>Adj/Removable BodiLink Head Support Hardware - MAX 8</b> <i>Required with more than 2 swing-away headrest accessories or chin boom.</i>	\$342.00
109624-99-0	<b>Extra 4" Link for BodiLink MAX 8 Hardware</b> <i>Sent as a part.</i>	\$102.00
11305	<b>Swing-Away Egg Switch Assembly</b> Mount on Upper Left      Mount on Upper Right	\$678.00
11306	<b>Swing-Away Microlight Assembly</b> Mount on Upper Left      Mount on Upper Right	\$777.00
11307	<b>Swing-Away Facial Pad Assembly 4"L x 2"T</b> Mount on Upper Left      Mount on Upper Right	\$552.00
11308	<b>Swing-Away Facial Pad Assembly 6"L x 2"T</b> Mount on Upper Left      Mount on Upper Right	\$552.00
11276	<b>Provide Headrest Receiver for AGILITY Backrest</b> <i>Select for AGILITY backs and third party backs ordered off chair.</i>	No charge
106331-99-0	<b>Permobil 10"W Headrest, Leatherette</b> <i>Only available with Corpus backrests. Must also select hardware.</i>	\$706.00
10552	<b>Adj/Removable Headrest Hardware</b>	\$152.00
107499-99-0	<b>Universal Headrest Adapter</b>	\$285.00

### Leg Support Options | (Seating Benefits)

103070-99-0	<b>UT Amp. Support Pad (H), Mesh Cover - 1 Each</b> <i>7"W x 7"T curved pad. Must also select hardware.</i>	\$424.00
103071-99-0	<b>UT Amp. Support Pad (D), Mesh Cover - 1 Each</b> <i>8"W x 8"T flat pad. Must also select hardware.</i>	\$424.00

## Accessories

10553	<b>Adj Removable - Amp Support Hardware - 1 Each</b> Mount on Left      Mount on Right	\$152.00
109040-99-0	<b>BodiLink Premium LPTS Pad (B), Stretch-Air - Pair</b> <i>4.5"W x 3.5"T pad. Must also select hardware.</i>	\$299.00
109043-99-0	<b>BodiLink Premium LPTS Pad (A), Stretch-Air - Pair</b> <i>6"W x 5.5"T pad. Must also select hardware.</i>	\$299.00
109041-99-0	<b>BodiLink Premium LPTS Pad (C), Stretch-Air - Pair</b> <i>8"W x 3.5"T pad. Must also select hardware.</i>	\$299.00
109042-99-0	<b>BodiLink Premium LPTS Pad (G), Stretch-Air - Pair</b> <i>12"W x 3.5"T pad. Must also select hardware.</i>	\$299.00
109044-99-0	<b>Adj/Removable BodiLink LPTS Hardware Med Corpus - Pair</b> <i>Height adjustable.</i>	\$353.00
109194-99-0	<b>Adj/Removable Permobil Lateral Thigh/Hip Hardware - Pair</b>	\$152.00
106762-99-0	<b>Abductor Kit Corpus - Leatherette</b> <i>Not compatible with Corpus VS Legrest. Must also select hardware.</i>	\$419.00
10555	<b>Adj Removable - Abductor Hardware - 1 Each</b>	\$152.00

### Chest Support Options | (Seating Benefits)

105840-99-0	<b>Bodypoint Monoflex Chest Support Belt Kit - Medium</b> <i>18.25"L x 4"W Pad Size.</i>	\$214.00
105839-99-0	<b>Bodypoint Monoflex Chest Support Belt Kit - Large</b> <i>21.5"L x 4.5"W Pad Size.</i>	\$214.00
105850-99-0	<b>Bodypoint Monoflex Chest Support Belt Kit - X-Large</b> <i>24.5"L x 5"W Pad Size.</i>	\$214.00

### Other Options

105260-99-0	<b>Transfer Handles (Seating Benefits)</b>	\$582.00
105292-99-0	<b>Push Handles, Corpus Seat</b>	No charge
106210-99-0	<b>RAM Self-Leveling Cup Holder (Client Cost)</b>	\$101.00
106806-99-0	<b>RAM X-Grip Phone Holder (Client Cost)</b>	\$101.00
107100-99-0	<b>RAM X-Grip 10" Tablet Holder w/ Swing-Away Mount (Client Cost)</b>	\$556.00
107189-99-0	<b>USB Charger (Client Cost)</b> <i>Provides up to 1.5A (at 5V) to charge most smartphones and tablets. No charge with 4 power seat function configurations.</i> Mount on Left Armrest      Mount on Left Seat Rail      Mount on Right Armrest Mount on Right Seat Rail	\$238.00
105805-99-0	<b>Medical Necessities Bag Hooks (Client Cost)</b>	\$110.00
102059-99-0	<b>Medical Necessities Bag (Client Cost)</b>	\$309.00
105262-99-0	<b>Joystick Cutout Tray - 16" W (Seating Benefits)</b> <i>Must also select hardware.</i>	\$517.00
105263-99-0	<b>Joystick Cutout Tray - 18" W (Seating Benefits)</b> <i>Must also select hardware.</i>	\$517.00
100578-99-0	<b>Small Angle Adjustable Tray (Seating Benefits)</b> <i>Must also select hardware.</i> Mount on Left      Mount on Right	\$517.00
10556	<b>Adj Removable - Tray Hardware - 1 Each (Seating Benefits)</b>	\$152.00
108230-99-0	<b>Permobil Adjustable Tray - 13"W (Seating Benefits)</b> <i>Must also select hardware.</i>	\$381.00
108228-99-0	<b>Permobil Adjustable Tray - 15"W (Seating Benefits)</b> <i>Must also select hardware.</i>	\$381.00
108227-99-0	<b>Permobil Adjustable Tray - 17"W (Seating Benefits)</b> <i>Must also select hardware.</i>	\$381.00
108229-99-0	<b>Permobil Adjustable Tray - 19"W (Seating Benefits)</b> <i>Must also select hardware.</i>	\$381.00

## Accessories

<b>Permobil Adjustable Tray Hardware</b>	Left Mount ( <a href="#">Seating Benefits</a> )	\$277.00
	Right Mount ( <a href="#">Seating Benefits</a> )	
105224-99-0	<b>Oxygen Holder, Corpus</b> ( <a href="#">Seating Benefits</a> )	\$376.00
105223-99-0	<b>Crutch Holder, Corpus</b> ( <a href="#">Seating Benefits</a> )	\$391.00

## Alternative Drive Systems

*Omni2 Control Module required for all alternative drive systems.*

### Alt Proportional Joysticks

109525-99-0	<b>mo-vis Micro Joystick for Omni</b> <i>Operating Force = 8 g (0.28 oz). Neutral to Full Deflection = 0.13" (3.3 mm).</i>	\$3,550.00
109527-99-0	<b>mo-vis Multi Joystick for Omni</b> <i>Operating Force = 50 g (1.8 oz). Neutral to Full Deflection = 0.31" (8 mm).</i>	\$3,550.00
103894-99-0	<b>Compact Joystick Lite for Omni - (9-pin connector)</b> <i>Operating Force = 50 g (1.8 oz). Neutral to Full Deflection = 0.7" (17mm).</i>	\$1,922.00
103893-99-0	<b>Compact Joystick for Omni - (9-pin connector)</b> <i>Operating Force = 220 g (7.8 oz). Neutral to Full Deflection = 1.1" (28mm).</i>	\$1,679.00

**Alt Proportional Joystick Mounts** | *Please see Alt Drive Accessories section for additional mounting options.*

107857-99-0	<b>Swing-Away Chin Boom for Micro/Multi Joystick - Left</b>	\$623.00
107947-99-0	<b>Swing-Away Chin Boom for Micro/Multi Joystick - Right</b>	\$623.00
107946-99-0	<b>Swing-Away Chin Boom for Compact Lite Joystick - Left</b>	\$623.00
107948-99-0	<b>Swing-Away Chin Boom for Compact Lite Joystick - Right</b>	\$623.00
104690-99-0	<b>Alt Joystick Swing-Away Mount Kit - Armrest</b> Mount on Left      Mount on Right	\$327.00
107135-99-0	<b>Alt Joystick / Omni Combo Swing-Away Mount Kit - Armrest</b> Mount on Left      Mount on Right	\$661.00
104143-99-0	<b>Alt Joystick Adjustable Mount Kit - Chin Boom</b> <i>Swing-Away / Removable Chin Boom mounted on right of backrest shell. Custom charge may be applied if ordering lower backrest heights.</i>	\$1,058.00
10752	<b>Chin Control Harness - Large</b> <i>For neck width up to 7½".</i>	\$755.00
10753	<b>Chin Control Harness - Small</b> <i>For neck width up to 5½".</i>	\$755.00

### Alt Proportional Joystick Accessories

104904-99-0	<b>Dual Switches for Compact Joystick - Fixed Mount</b> Mount Left of JS      Mount Right of JS      Mount front of JS      Mount back of JS	\$313.00
104641-99-0	<b>Dual Switches for Compact Joystick - Gooseneck Mount</b>	\$571.00

## Alternative Drive Systems

### Total Control Head Array

109558-99-0	<b>Permobil Total Control Head Array</b> <i>Features an occipital pad sensor and two lateral arms configured with proximity sensors. Must select swing away hardware and mounting hardware below to complete assembly.</i>	\$6,537.00
11072	<b>Program Head Array as a 5-Switch System</b> <i>Select additional arms/sensors below and indicate the intended use of each sensor in the notes section.</i>	No charge
109563-99-0	<b>Adj/Removable BodiLink HA Hardware - MAX 8</b>	\$347.00
109564-99-0	<b>Adj/Removable BodiLink HA Hardware - MAX 8 Flip Down</b> <i>Flip-down joint that allows head array to swing back 120°.</i>	\$604.00
10976	<b>(Optional) Add Left Lateral Arm w Proximity Sensor</b>	\$472.00
10977	<b>(Optional) Add Left Lateral Arm w Egg Switch</b>	\$314.00
10978	<b>(Optional) Add Left Lateral Arm w Pad (No Sensor or Switch)</b>	\$193.00
10979	<b>(Optional) Add Right Lateral Arm w Proximity Sensor</b>	\$472.00
10980	<b>(Optional) Add Right Lateral Arm w Egg Switch</b>	\$314.00
10981	<b>(Optional) Add Right Lateral Arm w Pad (No Sensor or Switch)</b>	\$193.00

### ASL FUSION Head Array

109582-99-0	<b>ASL FUSION Proportional and Digital Head Array</b> <i>Includes occipital (ASL 935) and lateral (ASL 933) pads, and pulse mechanical switch (ASL 320BK) for reset/mode changes. Must select hardware listed below to complete assembly.</i>	\$8,941.00
11580	<b>Swing-Away Hardware for Lateral Sensor - Left</b>	\$447.00
11581	<b>Swing-Away Hardware for Lateral Sensor - Right</b>	\$447.00
11582	<b>Swing-Away Hardware for Reset/Mode Switch - Right</b>	\$447.00
11583	<b>Adj/Removable Headrest Mounting Hardware</b> <i>Includes ASL 606 and BodiLink receiver.</i>	\$447.00
109583-99-0	<b>ASL FUSION Display, Programmer and Attendant Control w/ Mount</b> <i>Includes ASL 165 and quick release mount.</i>	\$951.00

### Sip & Puff Systems

10573	<b>Sip &amp; Puff System for Omni</b>	\$2,020.00
10574	<b>Sip &amp; Puff System with Fiber Optic Opti-Stop</b>	\$2,687.00
10762	<b>Adjustable Mounting Hardware - Sip &amp; Puff</b> <i>Required with any sip &amp; puff system.</i>	\$327.00

### Proportional Foot Control

10767	<b>Proportional Foot Control</b> <i>Not available with ActiveReach packages.</i> Mount on Left      Mount on Right	\$3,545.00
10766	<b>Adjustable Mounting Hardware - Foot Control</b> <i>Required with proportional foot control.</i>	\$327.00

## Alt Drive System Setup

### Omni2 Options | *Required for all alternative drive systems.*

108500-99-0	<b>Omni2 Control Module for R-net - Corpus</b> <i>3.5" VGA high resolution display, USB charger port, built in IR control and integrated Bluetooth technology. Must select mounting hardware below. Specify alternative input device being used:</i>	\$2,200.00
103476-99-0	<b>Adjustable Swing-Away Hardware - Omni/Tablet</b> Mount on Left      Mount on Right	\$321.00

## Alt Drive Accessories

### PERMOfix Mounting Solutions

106655-99-0	<b>PERMOfix UniTrack Fixed Mount Kit</b>	\$351.00
106656-99-0	<b>PERMOfix UniTrack Removable Mount Kit</b>	\$617.00
106657-99-0	<b>PERMOfix UniTrack Swing-Away Mount Kit</b>	\$889.00
105907-99-0	<b>Proximity Sensor + Swivel Assembly</b> <i>Swivel assembly allows sensor to be mounted with PERMOfix. Since Proximity Sensors require power, CANNOT be used as On/Off switches.</i>	\$459.00
105919-99-0	<b>Egg Switch + Swivel Mount Assembly</b> <i>Swivel assembly allows switch to be mounted with PERMOfix.</i>	\$306.00
105912-99-0	<b>Solid Oval Pad + Swivel Assembly</b> <i>Swivel assembly allows switch to be mounted with PERMOfix.</i>	\$186.00
106675-99-0	<b>Oval Plate + Swivel Assembly - No Switch</b> <i>Blank plate provides mounting option for Microlight, Buddy Button or other switch options. Swivel assembly allows pad to be mounted with PERMOfix.</i>	\$161.00

### Switches and Cables

101769-99-0	<b>Microlight Switch   Qty:</b>	\$202.00
101767-99-0	<b>Buddy Button   Qty:</b>	\$153.00
106293-99-0	<b>Egg Switch Black   Qty:</b>	\$168.00
101771-99-0	<b>Mini Cup Switch   Qty:</b>	\$150.00
107500-99-0	<b>Proximity Switch   Qty:</b>	\$706.00
10772	<b>Dome Switch, Backrest Mount w/ Gooseneck   Qty:</b> <i>Mounts at top of backrest shell. Custom charge may be applied if ordering with short backrest heights.</i> Mount on Left      Mount on Right	\$629.00
10846	<b>Dome Switch, Backrest Swing-Away Mount w/ Gooseneck   Qty:</b> <i>Swing-away, mounts at top of backrest shell. Custom charge may be applied if ordered with short backrest heights.</i> Mount on Left      Mount on Right	\$943.00
105935-99-0	<b>3.5mm Ext. Cable for Sensor/Switch, 0.3m Long</b>	\$44.00
105931-99-0	<b>3.5mm Ext. Cable for Sensor/Switch, 1m Long</b>	\$48.00
105932-99-0	<b>3.5mm Ext. Cable for Sensor/Switch, 2m Long</b>	\$52.00
105933-99-0	<b>3.5mm Ext. Cable for Sensor/Switch, 3m Long</b>	\$58.00
105934-99-0	<b>3.5mm Y-splitter Cable for Sensor/Switch</b> <i>Allows two Proximity Sensors/Mechanical Switches to connect to the same direction jack.</i>	\$48.00

## Additional Order Instructions (for Permobil):

## Notes & Comments (not for Permobil):

## Chair Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- The client's name or code and the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

## Order Acknowledgement:

I, \_\_\_\_\_, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: \_\_\_\_\_

Printed Name: \_\_\_\_\_

Title: \_\_\_\_\_

Date: \_\_\_\_\_

*NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty.*