


Tricycle dimensions

User dimensions (inches)	X320 small	X330 medium	X340 large
Inside leg length	12 – 22	19 – 28	25 – 38
 Key user dimension: inside leg length The extended leg should reach from seat to pedal comfortably when both feet and torso are strapped into place. Important: User's weight must not exceed the maximum working load. Rifton Tricycles are not intended for clients with strong reflexes or poorly developed trunk balance and head control. Consult your client's therapist before ordering.			
Item dimensions (inches)	X320 small	X330 medium	X340 large
Overall width	25½	27½	31½
Overall length	44	53	65
Overall height (at lowest position with wheels / front end removed for transport)	21¼	21¼	30
Overall length with no rear handle (wheels and front end removed)	36	41	49
Overall length with rear handle (wheels and front end removed)	48	51	60
Wheel diameter	12	15	18
Seat center to extended pedal	12 – 22	19 – 28	25 – 36
Seat to top of trunk support (max. height)	16	24	24
Gear ratio	1:1	1:1	1:1
Item weight (lb)	47	56	65
Adjustment range of laterals	small trunk support: 6–12	small trunk support: 6–12 large trunk support: 9½–15½	large trunk support: 9½–15½
Max. working load (lb)	125	250	310

Small and medium Trikes come in your choice of five colors.




Blue




Red



Teal



Lime



Pink

Large Trike comes in your choice of five colors.



Blue



Red



Teal



White



Black



Retail order form - January 2025

All Rifton products are made to order and are non-returnable and non-exchangeable.



		* Required	X320 small	X330 medium	X340 large
*Tricycle frame includes pedals, parking brake with variable resistance, and a rear basket.			X321 \$3,562 ☐	X331 \$3,969 ☐	X341 \$5,050 ☐
*Color (note that the choices are different between sizes)			CHOOSE Blue ☐ Red ☐ Teal ☐ Lime ☐ Pink ☐ Blue ☐ Red ☐ Teal ☐ Lime ☐ Pink ☐ Blue ☐ Red ☐ Teal ☐ Black ☐ White ☐		
*Handlebar	CHOOSE	Loop	X322 \$422 ☐	X332 \$444 ☐	X342 \$449 ☐
		Conventional	X323 \$422 ☐	X333 \$444 ☐	X343 \$449 ☐
*Back support (trunk support comes with an appropriately-sized butterfly harness for anterior support)	CHOOSE	Small backrest pad	X326 \$304 ☐	X326 \$304 ☐	N/A
		Large backrest pad	N/A	X336 \$360 ☐	X336 \$360 ☐
		Small trunk support	X327 \$657 ☐	X327 \$657 ☐	N/A
		Large trunk support	N/A	X337 \$764 ☐	X337 \$764 ☐
*Seat (seat belt is included)	CHOOSE	Small	X324 \$612 ☐	N/A	N/A
		Medium	X325 \$702 ☐	X334 \$702 ☐	N/A
		Large	N/A	X335 \$764 ☐	(included)
Headrest (requires trunk support)	CHOOSE	Contoured	X376 \$506 ☐	X376 \$506 ☐	X376 \$506 ☐
		Winged	X378 \$562 ☐	X378 \$562 ☐	X378 \$562 ☐
Laterals (require trunk support)	CHOOSE	For small trunk support	X328 \$432 ☐	X328 \$432 ☐	N/A
		For large trunk support	N/A	X338 \$461 ☐	X338 \$461 ☐
Rear handle (rear steer must be factory installed)	CHOOSE	Push handle	X312 \$202 ☐	X312 \$202 ☐	X313 \$202 ☐
		Push handle w/ caregiver brake	X314 \$354 ☐	X314 \$354 ☐	X315 \$360 ☐
		Rear steer	X352 \$708 ☐	X353 \$708 ☐	X354 \$708 ☐
		Rear steer w/ caregiver brake	X362 \$854 ☐	X363 \$854 ☐	X364 \$866 ☐
Hand anchors (pair)	CHOOSE	Small	X317 \$253 ☐	X317 \$253 ☐	N/A
		Large	N/A	X318 \$259 ☐	X318 \$259 ☐
Front handbrake (must be factory installed)			X329 \$304 ☐	X339 \$304 ☐	X349 \$304 ☐
Front guide bar			X319 \$360 ☐	X319 \$360 ☐	X319 \$360 ☐
Abductor			X355 \$422 ☐	X356 \$427 ☐	X357 \$461 ☐
Abductor with adduction straps			X365 \$562 ☐	X366 \$584 ☐	X367 \$607 ☐
Pelvic harness			X372 \$365 ☐	X373 \$382 ☐	X374 \$404 ☐
Pedal raise kit			X344 \$107 ☐	X344 \$107 ☐	X344 \$107 ☐
Stationary stand			X345 \$208 ☐	X345 \$208 ☐	X345 \$208 ☐
Mounting bar (for device holder)			X346 \$107 ☐	X346 \$107 ☐	X346 \$107 ☐
Front basket			X347 \$107 ☐	X347 \$107 ☐	X347 \$107 ☐
Outdoor cover			X348 \$714 ☐	X358 \$753 ☐	X368 \$764 ☐
Total			X320 \$ _____	X330 \$ _____	X340 \$ _____