

Titanium adjustable wheelchairs

TiLite TRA | TiLite ZRA



permobil

TiLite TRA

Our most advanced titanium adjustable chair, the TRA's dual-tube frame minimises weight while maximising performance.

The geometry of its titanium frame stiffens the wheelchair, transferring more of your push stroke directly to the ground resulting in greater stamina and a more effortless ride.



TiLite ZRA

Distinguished by a clean, minimal look, the mono-tube frame's 1.25" tubing gives the ZRA strength and durability.

In a class of its own, the ZRA has a lightweight and fully adjustable titanium frame. The minimal design results in excellent transportability and easy placement on a car seat.



Tattoo options


- 1 Mosaic
- 2 Flower Garden
- 3 Need for Speed



Frame colours



Component colours



Specifications

TiLite TRA

Frame features	Titanium Dual-Tube Rigid, TiFit, Full adjustability	Front seat width (seat taper)	0", 1", 2"
Total chair width*	23"	Folding seat back	Adj. height 9.5" – 22" Adj. angle 80° – 101°
Seat width	12" – 20"	Ergonomic seat	6", 7", 8"
Seat depth	12" – 20"	Transit options	Q-Straint lock down approved
Front seat height	15" – 21"	Transport weight*	5 kg
Rear seat height	13" – 21"	Weight limit	120 kg Standard 159 kg Heavy duty
Front frame angle	70°, 75°, 80°, 85°, or 90°	Warranty	Five years on frame One year on TiLite components
Camber	0°, 2°, 4°, 6°, or 8°		
Centre of gravity adjustment	5"		

TiLite ZRA

Frame features	Titanium Mono-Tube Rigid, TiFit, Full adjustability	Front seat width (seat taper)	0", 1.5"
Total chair width*	23"	Folding seat back	Adj. height 9.5" – 22" Adj. angle 80° – 101°
Seat width	10" – 20"	Ergonomic seat	6", 7", 8"
Seat depth	10" – 20"	Transit options	Transit tie-down
Front seat height	16" – 21"	Transport weight*	4.8 kg
Rear seat height	12.5" – 21"	Weight limit	120 kg Standard 136 kg Heavy duty
Front frame angle	70°, 75°, 80°, 85°, or 90°	Warranty	Five years on frame One year on TiLite components
Camber	0°, 2°, 4°, 6°, or 8°		
Centre of gravity adjustment	6"		

* Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 3.3 kg.