

Order form

Vicair® Vector O2

HCPCS coding: E2624/E2625 | Adjustable skin protection & positioning

Vector O2



Company Name/ACCT #: _____

P.O. Number: _____

Requested By: _____

Phone: _____ Fax: _____

Email: _____

Ship To: _____

Client Reference: _____

For best results, do not fill out in your browser. Interactive form should be completed using Adobe Reader after saving to your local drive. Then email or print and fax to Customer Support at orders.comfort@permobil.com or 1.406.522.8563.

DO NOT SEND PROTECTED HEALTH INFORMATION. IT IS NOT NEEDED TO MAKE THE PRODUCT YOU ARE REQUESTING.

QUESTIONS? PLEASE CONTACT CUSTOMER SUPPORT FOR ASSISTANCE. | 800.736.0925

1 Size & style

Choose the fabric type and size for the Vector O2 Standard Profile. If a Kwik Strap is needed, choose a size on page 2.

VECTOR O2

HCPCS CODE: E2624/E2625

 = Nonstock items available at extended lead times and a \$40.00 upcharge

INCOTEC									
\$515.00					\$567.00				
22"					OVRs-CI-1822	OVRs-CI-1922	OVRs-CI-2022	OVRs-CI-2222	OVRs-CI-2422
21"						OVRs-CI-1921			
20"					OVRs-CI-1620	OVRs-CI-1720	OVRs-CI-1820	OVRs-CI-1920	OVRs-CI-2020
19"					OVRs-CI-1619	OVRs-CI-1719	OVRs-CI-1819	OVRs-CI-1919	OVRs-CI-2019
18"	OVRs-CI-1418	OVRs-CI-1518	OVRs-CI-1618	OVRs-CI-1718	OVRs-CI-1818	OVRs-CI-1918	OVRs-CI-2018	OVRs-CI-2218	OVRs-CI-2418
17"					OVRs-CI-1617	OVRs-CI-1717	OVRs-CI-1817	OVRs-CI-2017	
16"	OVRs-CI-1416	OVRs-CI-1516	OVRs-CI-1616	OVRs-CI-1716	OVRs-CI-1816		OVRs-CI-2016		
15"					OVRs-CI-1515	OVRs-CI-1615			
14"	OVRs-CI-1414								
	14"	15"	16"	17"	18"	19"	20"	22"	24"

INCOTEC

For fluid protection & an easily cleaned surface.

Or

COMFAIR

Fore user comfort & heat dissipation.

COMFAIR									
\$515.00					\$567.00				
22"					OVRs-CC-1822	OVRs-CC-1922	OVRs-CC-2022	OVRs-CC-2222	OVRs-CC-2422
21"						OVRs-CC-1921			
20"					OVRs-CC-1620	OVRs-CC-1720	OVRs-CC-1820	OVRs-CC-1920	OVRs-CC-2020
19"					OVRs-CC-1619	OVRs-CC-1719	OVRs-CC-1819	OVRs-CC-1919	OVRs-CC-2019
18"	OVRs-CC-1418	OVRs-CC-1518	OVRs-CC-1618	OVRs-CC-1718	OVRs-CC-1818	OVRs-CC-1918	OVRs-CC-2018	OVRs-CC-2218	OVRs-CC-2418
17"					OVRs-CC-1617	OVRs-CC-1717	OVRs-CC-1817	OVRs-CC-2017	
16"	OVRs-CC-1416	OVRs-CC-1516	OVRs-CC-1616	OVRs-CC-1716	OVRs-CC-1816		OVRs-CC-2016		
15"					OVRs-CC-1515	OVRs-CC-1615			
14"	OVRs-CC-1414								
	14"	15"	16"	17"	18"	19"	20"	22"	24"

Choose the size needed for the Vector O2 with Kwik Strap. If no Kwik Strap is needed, choose a size and fabric for Vector O2 Standard Profile on page 1.


Vector O2
with Kwik Strap

INCOTEC					
Depth	\$578.00			\$525.00	
	20"	OVRS-CI-1620-KS	OVRS-CI-1820-KS	OVRS-CI-2020-KS	OVRS-CI-2220-KS
	18"	OVRS-CI-1618-KS	OVRS-CI-1818-KS	OVRS-CI-2018-KS	OVRS-CI-2218-KS
	16"	OVRS-CI-1616-KS	OVRS-CI-1816-KS	OVRS-CI-2016-KS	
		16"	18"	20"	22"
				24"	
		Width			

2 Active O2 accessories

2.1 Removable solid seat pan
& hardware kit
HCPCS code: E2231 MSRP \$440.00

The solid seat pan kit includes a slotted aluminum pan and attaching hardware to accommodate 7/8" or 1" tubing. The kit also comes with two different cross bars to accommodate various wheelchair frame types. This will fit both folding and non-folding seat rails. Lateral thigh support hardware can be easily attached.



Solid seat pan & hardware kit
(RSSP-1-___)

Extra Vicair® Air Cells
MSRP \$50.00 / 50 air cell pack



Pack of 50 Vicair® Air Cells (AIRCELL50)


3 BodiLink® Lateral Pelvic/Thigh support hardware & pads

Please select sizes in the appropriate fabric & style charts below. Contact Customer Support for custom sizes that are not listed below 800.736.0925.

3.1 Premium pad shape & size
HCPCS coding: E0953 MSRP \$72.00

- You may make one selection per column (Left, Extra Left, Right, and Extra Right) in section 3.1.
- Make a fabric selection in 3.2.

*Length (L) refers to the actual size dimension of the support from top to bottom edge. Depth (D) refers to the actual size dimension from anterior to posterior edge.

Style	Size	Left	Extra left	Right	Extra right
	3.5"L x 4"D	BL-LPTSP2Z-3L4D-LH	BL-LPTSP2Z-3L4D-LH	BL-LPTSP2Z-3L4D-RH	BL-LPTSP2Z-3L4D-RH
	3.5"L x 8"D	BL-LPTSP2Z-3L8D-LH	BL-LPTSP2Z-3L8D-LH	BL-LPTSP2Z-3L8D-RH	BL-LPTSP2Z-3L8D-RH
	3.5"L x 12"D	BL-LPTSP2Z-3L12D-LH	BL-LPTSP2Z-3L12D-LH	BL-LPTSP2Z-3L12D-RH	BL-LPTSP2Z-3L12D-RH
	5.5"L x 6"D	BL-LPTSP2Z-5L6D-LH	BL-LPTSP2Z-5L6D-LH	BL-LPTSP2Z-5L6D-RH	BL-LPTSP2Z-5L6D-RH

3.2 Premium pad cover with foam insert

Cover options	MSRP	Left		Extra left		Right		Extra right	
COMFORT-TEK	\$0.00/ea.		LPTS-P2-COMFORT-FM		LPTS-P2-COMFORT-FM		LPTS-P2-COMFORT-FM		LPTS-P2-COMFORT-FM
STRETCH-AIR	\$0.00/ea.		LPTS-P2-STRETCH-FM		LPTS-P2-STRETCH-FM		LPTS-P2-STRETCH-FM		LPTS-P2-STRETCH-FM
GLIDEWEAR™	\$15.00/ea.		LPTS-P2-GLIDE-FM		LPTS-P2-GLIDE-FM		LPTS-P2-GLIDE-FM		LPTS-P2-GLIDE-FM



COMFORT-TEK

For fluid protection & an easily cleaned surface.



STRETCH-AIR

For comfort & heat dissipation.



For skin protection & shear reduction.
Not fluid resistant.

3.3 BodiLink® Lateral Pelvic/Thigh Support hardware

- You may make one selection per column (Left, Extra Left, Right, and Extra Right) in section 3.3 from group A, B, -OR- C.
- TT long extension arms (section 3.4) will only work with TT style hardware.

A. TT hardware, slot mount


TT hardware slot mount - Fixed



TT hardware slot mount - Removable



Hardware size VS. cushion thickness




Hardware size	Max cushion thickness clearance	Recommended cushion thickness
Small	3.0"	2.0"- 3.0"
Medium	4.5"	3.0"- 4.5"
Large	6.5"	4.0"- 6.5"


		Size	MSRP	Left		Extra left		Right		Extra right	
Fixed		Small	\$161.00		BL-LPTS-TT1FXSL1-LH		BL-LPTS-TT1FXSL1-LH		BL-LPTS-TT1FXSL1-RH		BL-LPTS-TT1FXSL1-RH
		Medium	\$161.00		BL-LPTS-TT1FXSL2-LH		BL-LPTS-TT1FXSL2-LH		BL-LPTS-TT1FXSL2-RH		BL-LPTS-TT1FXSL2-RH
		Large	\$161.00		BL-LPTS-TT1FXSL3-LH		BL-LPTS-TT1FXSL3-LH		BL-LPTS-TT1FXSL3-RH		BL-LPTS-TT1FXSL3-RH
Or Removable E1028		Small	\$282.00		BL-LPTS-TT1RMSL1-LH		BL-LPTS-TT1RMSL1-LH		BL-LPTS-TT1RMSL1-RH		BL-LPTS-TT1RMSL1-RH
		Medium	\$282.00		BL-LPTS-TT1RMSL2-LH		BL-LPTS-TT1RMSL2-LH		BL-LPTS-TT1RMSL2-RH		BL-LPTS-TT1RMSL2-RH
		Large	\$282.00		BL-LPTS-TT1RMSL3-LH		BL-LPTS-TT1RMSL3-LH		BL-LPTS-TT1RMSL3-RH		BL-LPTS-TT1RMSL3-RH

B. TT hardware, power mount


TT hardware power mount - Fixed



TT hardware power mount - Removable



Hardware size VS. cushion thickness



Hardware size	Max cushion thickness clearance	Recommended cushion thickness
Small	3.5"	2.0"- 3.5"
Medium	5.5"	3.0"- 5.5"
Large	7.5"	4.0"- 7.5"

		Size	MSRP	Left		Extra left		Right		Extra right	
Fixed		Small	\$211.00		BL-LPTS-TT1FXPWL1-LH		BL-LPTS-TT1FXPWL1-LH		BL-LPTS-TT1FXPWL1-RH		BL-LPTS-TT1FXPWL1-RH
		Medium	\$211.00		BL-LPTS-TT1FXPWL2-LH		BL-LPTS-TT1FXPWL2-LH		BL-LPTS-TT1FXPWL2-RH		BL-LPTS-TT1FXPWL2-RH
		Large	\$211.00		BL-LPTS-TT1FXPWL3-LH		BL-LPTS-TT1FXPWL3-LH		BL-LPTS-TT1FXPWL3-RH		BL-LPTS-TT1FXPWL3-RH
Or Removable E1028		Small	\$324.00		BL-LPTS-TT1RMPWL1-LH		BL-LPTS-TT1RMPWL1-LH		BL-LPTS-TT1RMPWL1-RH		BL-LPTS-TT1RMPWL1-RH
		Medium	\$324.00		BL-LPTS-TT1RMPWL2-LH		BL-LPTS-TT1RMPWL2-LH		BL-LPTS-TT1RMPWL2-RH		BL-LPTS-TT1RMPWL2-RH
		Large	\$324.00		BL-LPTS-TT1RMPWL3-LH		BL-LPTS-TT1RMPWL3-LH		BL-LPTS-TT1RMPWL3-RH		BL-LPTS-TT1RMPWL3-RH


TT hardware, power mount option

Note: If choosing the Power Mount Hardware, choose the appropriate power mount type below.

Power mount	MSRP	Part number
QUANTUM	add \$0.00	LPTS-TT1-PW1
PERMOBIL	add \$0.00	LPTS-TT1-PW2
QUICKIE	add \$0.00	LPTS-TT1-PW3
ROVI	add \$0.00	LPTS-TT1-PW4
AVID REHAB	add \$0.00	LPTS-TT1-PW5

3.4 TT Long Extension Arm 4"

MSRP \$43.00/arm





Left	Extra left	Right	Extra right
BL-LPTS-TT1LARM4	BL-LPTS-TT1LARM4	BL-LPTS-TT1LARM4	BL-LPTS-TT1LARM4

4 BodiLink® Medial Knee/Thigh support hardware & pads

- You may make one selection per column.
- Choose a cover material in 4.2.
- Choose between wedge or oval in 4.1.
- Choose hardware size in 4.3.

4.1 Pad size & shape

HCPCS coding: E0957 MSRP \$148.00

Style	Size	Part number
 Wedge	2.5"W x 3.5"D	BL-MKTSP1-3W4D
	3.5"W x 5"D	BL-MKTSP1-4W5D
	4"W x 6"D	BL-MKTSP1-4W6D
 Oval	2.5"W x 3.5"D	BL-MKTSP2-3W4D
	3.5"W x 5"D	BL-MKTSP2-4W5D
	5"W x 6"D	BL-MKTSP2-5W6D


4.2 Pad cover material

 COMFORT-TEK 	For skin protection & shear reduction. Not fluid resistant.
STRETCH-AIR 	For fluid protection & an easily cleaned surface.
	For patient comfort & heat dissipation.

	GlideWear® (MKTSP-GLIDE) (add \$15.00)
	Comfort-Tek® (MKTSP-COMFORT)
	Stretch-Air® (MKTSP-STRETCH)



4.3 Medial Knee/Thigh support Swing Away hardware

HCPCS coding: E1028 MSRP \$205.00



Size	Part number
Small	BL-MKTS-ST1SASL1
Medium	BL-MKTS-ST1SASL2
Large	BL-MKTS-ST1SASL3

*Swing Away Hardware only, not available in Fixed Hardware
 **MKTS Hardware must be mounted to Removable Solid Seat Pan in section 2.4
 ***If Leg Length Discrepancy is needed on Removable Solid Seat Pan, Offset Orientation must be selected



Size	Recommended inferior thigh thickness
Small	1.00" - 3.00"
Medium	2.00" - 5.00"
Large	3.00" - 7.00"