



# PDG Fuze T50 - AADL W653



Dealer Name\*: \_\_\_\_\_ Account #\*: \_\_\_\_\_

City, State, Zip: \_\_\_\_\_

Phone: \_\_\_\_\_ Email: \_\_\_\_\_ ATP: \_\_\_\_\_

Rep Name: \_\_\_\_\_ P.O. Number: \_\_\_\_\_ Quote Number: \_\_\_\_\_

Client Reference\*: \_\_\_\_\_ Client Weight: \_\_\_\_\_

Special considerations that need to be addressed here (e.g., diagnosis): \_\_\_\_\_ \*Required field

## October 2024, Revised January 2026. Prices are in Canadian Dollars (CAD)

### Wheelchair & Frame

**Fuze T50 50° Manual Tilt Wheelchair** Also available: the narrow Fuze T50n with an overall chair width of only seat width + 4". See separate order form. W653 \$ 2,685.00

#### Frame Type

- Fixed Width Frame, 300lbs Weight Capacity** Select a fixed seat width below. No charge
- Fixed Width Frame, 400lbs Weight Capacity** Not compatible with transit tie down system. Select a fixed seat width below. Prior Approval \$ 300.00
- Adjustable Width Frame. 250lbs Weight Capacity** Select an adjustable seat width below. \$ 200.00

#### Frame Options

**Transit Tie Down System (TTDS)** Includes 4 factory installed transport brackets. TTDS has been tested according to ANSI/RESNA WC19. See Owners Manual for complete details. Tested to maximum weight capacity of 300lbs. No charge

### Seat Measurements

**Seat Width** | Overall chair width is +7.5" with 12/16" rear wheels and +8.5" with 20/22/24" rear wheels.

<b>Fixed:</b>	15"	16"	17"	18"	19"		No charge
	20"	21"	22"			W111	\$ 250.00
	23"	24"				W111	\$ 500.00
<b>Adjustable 15"-20":</b>	15"	16"	17"	18"	19"		No charge
	20"					W111	\$ 250.00

#### Seat Depth

<b>Select Depth:</b>	15"	16"	17"	18"	19"		No charge
	20"	21"	22"			W111	\$ 250.00

*21"-22" seat depths are not adjustable and require straight stroller handle back. Not compatible with reclining or fold-down back.*

### Frame Finishes

Select Color: High Gloss Black      Hyper Blue

## Armrests

### Armrest Type

No Armrests

No charge

Single Post, Height Adjustable Armrests

No charge

10-14" Range    12"-16" Range

Armrest Position:    Rear    Front

*If left blank default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions PDG may not be able to accommodate stated preference.*

### Armpads

Regular Armpads

No charge

Desk Length 2" x 10"    Full Length 2" x 14"

## Seat Height

Seat to floor height is adjustable within the range of the chosen axle plate (Low/Standard/High). Chairs with swing-away front rigging are set with 3° of pre-tilt to avoid interference. This increases the actual seat to floor height at the front of the seat by up to 1" compared to the selected value. Easy-Out and standard plug-in front rigging are not affected.

RW

	5" Casters					6" Casters & 6" Forks					6" Casters & 8" Forks					8" Casters & 8" Forks				
24"			Low	Low	Low		Low	Low	Low	Std	Low	Low	Low	Std	Std	Low	Low	Std	Std	Std
22"		Low	Low	Low	Std	Low	Low	Low	Std	Std	Low	Low	Std	Std	Std	Low	Std	Std	Std	
20"	Low	Low	Low	Std	Std	Low	Low	Std	Std	Std	Low	Std	Std	Std		Std	Std	Std		High
16"		Std	Std	Std		Std	Std	Std		High	Std	Std		High	High	Std		High	High	High
12"	Std	Std		High	High	Std		High	High	High		High	High	High		High	High	High	High	High
	13"	14"	15"	16"	17"	14"	15"	16"	17"	18"	15"	16"	17"	18"	19"	16"	17"	18"	19"	20"
	Seat to Floor Height																			

## Seat Pan

### Seat Pan

No Seat Pan

No charge

Solid Aluminium Seat Pan

No charge

## Front Casters

### Urethane Casters

5" x 1.25" Caster

No charge

6" x 1.25" Caster

No charge

8" x 1.25" Caster

No charge

## Rear Wheels

Composite Rear Wheel

No charge

12"    16"    20"    22"    24"

Heavy Duty Spoke Wheel

\$ 350.00

24"

Rear Wheel Position:    Rear (More Stable)    Mid (Default)    Front (Less Stable)

*If left blank, default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.*

### Axles

Quick Release Rear Axles

No charge

Fixed Rear Axles

No charge

### Rear Tires

Urethane Tires

No charge

**Handrims** | Handrims are only available with 20", 22" and 24" rear wheels.

No Handrims	No charge
Standard Aluminium Handrims	No charge
Plastic Coated Handrims <i>Not compatible with drum brakes.</i>	\$ 175.00

## Wheel Locks

User Wheel Locks	No charge
Push to Lock    Pull to Lock	

## Wheel Lock Options

Wheel Lock Extensions <i>Not compatible with attendant wheel lock.</i>	No charge
--	-----------

## Anti-Tippers

Rear Anti-Tippers <i>Required safety option.</i>	No charge
--	-----------

## Backrest

Back Post Style:	Stroller Handle Back
	<i>Cane handle back not compatible with power tilt with manual override.</i>
Back Post Type:	Straight    8° Bend
Back Post Height:	Short (17")    Tall (23")
	<i>Measured from seat pan to handle.</i>
Folding Style:	Fold-Down
	<i>Fold-down back not compatible with reclining backrest or dynamic rocker back and will be substituted with a non-folding back.</i>
Seat to Back Angle:	Set at 90°    Set at 98°    Set at 106°    Set at 114°    Set at 122°

## Backrest Options

Angle/Height Adjustable Stroller Handle: <i>Not compatible with cane handle back.</i>	No charge
Button on Outside    Button on Inside	

## Tilt Activation

### Attendant Tilt Only

Attendant Tilt Only, Right Back Post	No charge
Attendant Tilt Only, Left Back Post	No charge

**Tilt Limit** | *Optional. If no tilt limit is selected the full tilt range will be available.*

Tilt Limit:	Set at 20°    Set at 30°    Set at 40°
-------------	--

## Front Rigging

No Front Rigging	No charge
------------------	-----------

**Easy-Out Front Rigging** | *Effortless insertion and removal. No hanger brackets remain on the frame.*

Fixed 70°	\$ 95.00
-----------	----------

**Swing Away Front Rigging** | *Legacy design. Not compatible with Easy-Out.*

Swing Away 70°	No charge
Swing Away 90° <i>Includes 11"-16" extension tubes. Not compatible with 8" casters.</i>	\$ 0.00
Swing Away Elevating Legrests (ELR) <i>Includes calf pads.</i>	W106 \$ 375.00
1" Offset Hangers for Swing Away Front Rigging	\$ 95.00

**Plug In Front Rigging** | *Legacy design. Not compatible with Easy-Out.*

Plug In 70°	No charge
-------------	-----------

## Footplates

Composite Footplates		\$ 0.00
Aluminium Angle Adjustable Footplates	W107	\$ 165.00

### Extension Tubes | Consider the height of your seat cushion.

Standard Extension Tubes, Short	No charge
Standard Extension Tubes, Long	\$ 75.00

#### Footplate Height:

Specify:

Ext	Short	11" - 16"	12" - 16"	4" - 9"	13" - 16"
	Long	14" - 22"	14" - 19"	12" - 16"	14" - 17"
		Standard Easy-Out	60°, 70°, ELR Swing-Away	90° Swing-Away	Standard Plug-In
Legrest					

## Front Rigging Accessories

Length Adjustable Heel Loops	No charge
Standard Calf Strap	\$ 50.00

## General

Permobil manufactures PDG wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

## Notes