

Customer and Order Information

Purchase Order *	Tag	v.0910
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in
		*required
Order Notes/Instructions		

Base Chair & Frame Options

☒ 1901000	Stellar GL 45° Manual Tilt Wheelchair	W321	\$ 2,685
Chair Weight Capacity			
☐ 1902000	250 lbs Capacity Available for 14" to 22" Seat Width.		
☐ 1902050	300 lbs Capacity Available for 14" to 22" Seat Width.	Prior Approval	\$ 200
☐ 1905006	450 lbs Capacity Available for 20" to 32" Seat Width. Not compatible with Transit Tie Down System.	Prior Approval	\$ 500
User Weight Range Required for setup of tilt uplift force.			
☐ 1902430	Up to 75 lbs	☐ 1902431	75 lbs to 100 lbs
☐ 1902432	100 lbs to 115 lbs	☐ 1902433	115 lbs to 150 lbs
☐ 1902434	150 lbs to 200 lbs	☐ 1902435	200 lbs to 250 lbs
☐ 1902436	250 lbs to 300 lbs	☐ 1902437	300 lbs to 375 lbs
☐ 1902438	375 lbs to 450 lbs		
☐ 1902300	Transit Tie Down System (TTDS) TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. Not compatible with Seat Widths greater than 22". Tested to max. Weight Capacity of 300 lbs.		
PDG Premium Colours			
☐ 1905001	Black Shimmer	☐ 1905015	Sky Blue
☐	Other Standard PDG Premium Colour	Specify Colour _____	pdgmobility.com/colours
			\$ 125

Seat Width & Depth

Overall chair width is +9.5" for 12/16" RW and +10.5" for 20/22/24" RW

Standard Seat Widths			
☐ 1903014	14" Seat Width	☐ 1903015	15" Seat Width
☐ 1903016	16" Seat Width	☐ 1903017	17" Seat Width
☐ 1903018	18" Seat Width	☐ 1903019	19" Seat Width
☐ 1903020	20" Seat Width	☐ 1903021	21" Seat Width
☐ 1903022	22" Seat Width		
			W111 \$ 250
Custom Seat Widths Require 450 lbs Weight Capacity Frame.			
☐ 1903023	23" Seat Width	☐ 1903024	24" Seat Width
☐ 1903025	25" Seat Width	☐ 1903026	26" Seat Width
☐ 1903027	27" Seat Width	☐ 1903028	28" Seat Width
☐ 1903029	29" Seat Width	☐ 1903030	30" Seat Width
☐ 1903031	31" Seat Width		
☐ 1903032	32" Seat Width		
			W111 \$ 500
			W111 \$ 500
			W111 \$ 750
			W111 \$ 750
Seat Depths			
☐ 1904016	16" Seat Depth	☐ 1904017	17" Seat Depth
☐ 1904018	18" Seat Depth	☐ 1904019	19" Seat Depth
☐ 1904020	20" Seat Depth	☐ 1904021	21" Seat Depth
☐ 1904022	22" Seat Depth		
			W111 \$ 250
☒ 1908001	Solid Aluminum Seat Pan		

Armrests

Single Post Height Adjustable Armrests			
☐ 1906111	8"-11" Range	☐ 1906114	10"-14" Range
☐ 1906116	12"-16" Range		
Dual Post Height Adjustable Armrests			
☐ 1906211	8"-11" Range	☐ 1906214	10"-14" Range
☐ 1906216	12"-16" Range		
Regular Armpads			
☐ 1907100	Desk Length - 2" x 10"	☐ 1907101	Full Length - 2" x 14"
☐ 1907310	Clothing Guards		
			\$ 125

Seat Height

➤

Don't forget

Select desired Seat-to-Floor Height

STF Height is adjustable within the range of the chosen Subframe Plate (Low/High).

6" Casters

	Low							High						
RW	13"	14"	15"	16"	17"	18"	RW	15"	16"	17"	18"	19"	20"	
12"	☐	☐	☐	☐			12"	☐	☐	☐	☐			
16"				☐	☐	☐	16"			☐	☐	☐	☐	
20"	☐	☐	☐	☐			20"	☐	☐	☐	☐			
22"		☐	☐	☐	☐		22"		☐	☐	☐	☐		
24"			☐	☐	☐	☐	24"			☐	☐	☐	☐	

8" Casters

	Low							High						
RW	13"	14"	15"	16"	17"	18"	RW	15"	16"	17"	18"	19"	20"	
16"			☐	☐	☐	☐	16"			☐	☐	☐	☐	
24"			☐	☐	☐	☐	24"			☐	☐	☐	☐	

Casters & Rear Wheels

Urethane Casters

☐ 1911061 6" x 1.25" Caster size

☐ 1911081 8" x 1.25" Caster size

1914000 Composite Wheels with Urethane Tires

☐ 1912012 12" Wheel size

☐ 1912016 16" Wheel size

☐ 1912020 20" Wheel size

☐ 1912022 22" Wheel size

☐ 1912024 24" Wheel size

☒ 1913000 Quick Release Rear Axles

Handdrims

Handdrims are only available with 20", 22", and 24" Rear Wheels.

☐ 1915010 Standard Aluminum

☐ 1915000 No Handdrims

☐ 1915020 Plastic Coated

\$ 175

Wheel Locks

☒ 1916010 Push to Lock Wheel Locks

☐ 1916100 Wheel Lock Extensions

Anti-Tippers

☒ 1917010 Rear Anti-Tippers

Required Safety Option.

Back Post Style

☒ 1918010 1919010 Angle/Height Adjustable Stroller Handle

Button on Outside.

Back Post Height

Measured from the Seat Pan to the Handle.

☐ 1918425 Short (22")

☐ 1918431 Tall (25")

☒ 1919090 90° Seat to Back Angle

☒ 1919220 Fold Down Back

Tilt/Recline Options

☒ 1923011 Attendant Tilt only

Trigger on Right Back Post.

Front Rigging

Plug In Front Rigging

☐ 1925410 Easy-Out Fixed 70°

Effortless insertion and removal. No Hanger Bracket remains on frame.

\$ 95

☐ 1925110 Standard Plug In 70°

Legacy Design. Not compatible with Easy-Out.

Swing Away Front Rigging

☐ 1925140 Swing Away 70°

☐ 1925200 Swing Away Elevating Legrests

☐ 1925300 1" Offset Hangers for Swing Away

Includes Calf Pads.

W106

\$ 375

\$ 95

Permobil | 318 East Kent Ave South | Vancouver BC V5X 4N6 Canada | Phone 888 858 4422 | Email sales.pdg@permobil.com

2 / 3

Foot Plates & Options				
Footplates				
☐	1926010	Composite		
☐	1926020	Aluminum Angle Adjustable	<i>W107</i>	\$ 165
☒	1927040	Length Adjustable Heel Loops		
Front Rigging Extension Tubes				
☐	1928010	Short 14"-18" for Easy-Out Front Rigging 13"-17" for Swing Away and Standard Plug In Front Rigging		
☐	1928020	Long 18"-24" for Easy-Out Front Rigging 15"-19" for Swing Away and Standard Plug In Front Rigging		\$ 75

Please send your completed order form by email to sales.pdg@permobil.com.