

Customer and Order Information			
Purchase Order *		Tag	v.0910
Dealer/Provider Name *			
Purchaser Name and Contact Information *			
Shipping Address *			
User Information (if available)		Weight _____ lbs	Height _____ ft/in
			*required
Order Notes/Instructions			

Base Chair & Frame Options						
☒	2701000	Bentley LT-R 20° Manual Tilt Wheelchair	Reverse Wheel Layout enhances user access for self-propulsion. pdgmobility.com/blt-r	W340	\$ 2,685	
☒	2401009	Lightweight Aluminum Frame and Reverse Wheel Layout Upgrade Required Option for the Bentley LT-R.			\$ 300	
Frame Type All frames offer 300 lbs Weight Capacity.						
☐	2702010	Fixed Width Frame Select a Fixed Seat Width below.				
☐	2702025	Adjustable Width Frame Select an Adjustable Seat Width below.				
					\$ 200	
User Weight Range Required for setup of tilt uplift force.						
☐	2702431	75 lbs to 100 lbs	☐	2702432	100 lbs to 115 lbs	
☐	2702434	150 lbs to 200 lbs	☐	2702435	200 lbs to 250 lbs	
☐	2702433	115 lbs to 150 lbs		☐	2702436	250 lbs to 300 lbs
☐	2702300	Transit Tie Down System (TTDS)				
TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details.						
PDG Premium Colours						
☐	2705001	Black Shimmer		☐	2705015	Sky Blue
☐		Other Standard PDG Premium Colour		Specify Colour _____	pdgmobility.com/colours	\$ 125

Seat Width & Depth				Overall chair width is +10.5"	
Fixed Seat Widths Order with Fixed Frame.					
☐	2703015	15" Seat Width	☐	2703016	16" Seat Width
☐	2703018	18" Seat Width	☐	2703019	19" Seat Width
☐	2703020	20" Seat Width	☐	2703021	21" Seat Width
☐	2703023	23" Seat Width	☐	2703024	24" Seat Width
					W111 \$ 250
					W111 \$ 500
Adjustable Seat Widths Order with Adjustable Frame.					
☐	2703215	15"-20" Range - Set at 15"	☐	2703216	15"-20" Range - Set at 16"
☐	2703218	15"-20" Range - Set at 18"	☐	2703219	15"-20" Range - Set at 19"
☐	2703220	15"-20" Range - Set at 20"			W111 \$ 250
Seat Depths					
☐	2704016	16" Seat Depth	☐	2704017	17" Seat Depth
☐	2704019	19" Seat Depth	☐	2704018	18" Seat Depth
☐	2704020	20" Seat Depth	☐	2704021	21" Seath Depth
☐	2704022	22" Seat Depth			W111 \$ 250
☒	2708001	Solid Aluminum Seat Pan			

Armrests					
Single Post Height Adjustable Armrests					
☐	2706111	8"-11" Range	☐	2706114	10"-14" Range
☐	2706116	12"-16" Range			
Regular Armpads					
☐	2707100	Desk Length - 2" x 10"	☐	2707101	Full Length - 2" x 14"
☐	2707310	Clothing Guards			\$ 125

Seat Height



Select desired Seat-to-Floor Height

STF Height is adjustable within the range of the chosen Subframe Plate (Low/High).

6" Casters	Low							High						
	RW	12"	13"	14"	15"	16"	17"	RW	14"	15"	16"	17"	18"	19"
	20"	☐	☐	☐	☐	☐	☐	20"	☐	☐	☐	☐	☐	☐
	22"	☐	☐	☐	☐	☐	☐	22"	☐	☐	☐	☐	☐	☐
	24"	☐	☐	☐	☐	☐	☐	24"	☐	☐	☐	☐	☐	☐

8" Casters	Low							High						
	RW	12"	13"	14"	15"	16"	17"	RW	14"	15"	16"	17"	18"	19"
	24"	☐	☐	☐	☐	☐	☐	24"	☐	☐	☐	☐	☐	☐

Casters & Large Wheels

Urethane Casters

8" Casters are compatible only with 24" Large Wheels.

☐ 2711061 6" x 1.25" Caster size ☐ 2711081 8" x 1.25" Caster size

Composite Wheels with Urethane Tires

☐ 2712020 20" Wheel size ☐ 2712022 22" Wheel size ☐ 2712024 24" Wheel size

☒ 2713000 Quick Release Axles

Handrim

☐ 2715010 Standard Aluminum ☐ 2715000 No Handrim

☐ 2715020 Plastic Coated

\$ 175

Wheel Locks

☒ 2716010 Push to Lock Wheel Locks

☐ 2716100 Wheel Lock Extensions

Anti-Tippers

Front Anti-Tippers

Required Safety Option. Must choose one type.

☐ 2717010 Standard

☐ 2717015 Spring-loaded Dampens the forward-tip movement of the chair during transfers and allows for improved ground clearance. Not removable.

\$ 200

Back Post

☒ 2718010 2719010 Angle/Height Adjustable Stroller Handle Button on Outside.

Back Post Height

Measured from the Seat Pan to the Handle.

☐ 2718425 Short (22") ☐ 2418431 Tall (25")

☒ 2719090 90° Seat to Back Angle

☒ 2719220 Fold Down Back

Tilt/Recline Options

☒ 2723011 Attendant Tilt only Trigger on Right Back Post.

Front Rigging

Plug In Front Rigging

☐ 2725410 Easy-Out Fixed 70° Effortless insertion and removal. No Hanger Bracket remains on frame.

\$ 95

☐ 2725110 Standard Plug In 70° Legacy Design. Not compatible with Easy-Out.

Footplates & Options				
Footplates				
☐	2726010	Composite		
☐	2726020	Aluminum Angle Adjustable	W107	\$ 165
☒	2727040	Length Adjustable Heel Loops		
Front Rigging Extension Tubes				
☐	2728010	Short 14"-18" for Easy-Out Front Rigging 13"-17" for Standard Plug In Front Rigging		
☐	2728020	Long 18"-24" for Easy-Out Front Rigging 15"-19" for Standard Plug In Front Rigging		\$ 75

Please send your completed order form by email to sales.pdg@permobil.com.