

TiLite TWIST

Play. Giggle. Grow.



Designed from the ground up for kiddos, the TWIST encourages exploration, independence, and, yes, playfulness. There are no parts or growth kits to buy. Instead, the TWIST grows with them.

permobil



Other chairs say they can grow too, but parts are required. With the TWIST, growth is built in. Sized from 8" x 8" up to 15" x 16", the TWIST offers 2" of growth in seat width and 3" of growth in seat depth.

Custom configuration allows you to add options and components to your child's chair to meet their needs and to optimize their riding experience. This includes choices of wheels, handrims, tires, wheel locks, color, and more.

Specifications

Frame features	Aluminum mono-tube rigid Growable Full adjustability	Camber	0°, 2°, 4°, 6°, 8°, or 12°
Total chair width*	19"	Center of gravity adjustment	6"
Seat width	8" – 15" (2" of built-in Seat Width Adj.)	Transport options	Transit tie-down ANSI RESNA WC-19 tested
Seat depth	8" – 16" (3" of built-in Seat Depth Adj.)	Transport weight*	12 lb.
Front seat height	13" – 20" (up to 2.25" of built-in front seat height adj.)	Weight limit	165 lb.
Rear seat height	12" – 19.5" (up to 3.5" of built-in rear seat height adj.)	HCPCS code	E1235, E1237
Front frame angle	70°, 80°, or 90°	Warranty	Five years on frame One year on TiLite components
Seat back type	Folding seat back Adj. height 8.5" – 21" Adj. angle 80° – 101°	* Weights and widths based on 12" x 12" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weights. For reference, a pair of standard 22" rear wheels weighs 7 lb.	


Frame colors

High gloss paint



TiLite red Neon pink Super wet black

Metallic paint



Tangerine metallic Acid green Metallic Ocean blue metallic

Pearlescent paint



Midnight blue pearl White pearl

Anodized



Red Pepper Burnt orange Sapphire

Candy paint




Ruby red Azure blue Electric plum

Matte textured paint



Flame orange OD green Hyper blue

Matte paint



Titanium grey Black lava

Matte paint



Ultraviolet Fuchsia Black