

TiLite T-Series Canadian Order Form - AADL #W395

Effective October 1, 2021 (New Items in Red)

DUAL-TUBE RIGID FRAME

Page 1 of 7

DEALER INFORMATION		RTS INFORMATION	
Dealer:	Account #:	RTS:	Phone:
PO:	Ordered By:	Email:	
Email:		SHIP TO INFORMATION	
Mark For:		Name:	Ship Via:
TiLite Sales Rep:		Address:	
TiLite Information:	Promo Code:	City:	Province: Postal Code:
Ship Date:	Order #:	Phone:	Fax:

Please Email or Fax to the Permobil Manual Customer Service Team at tilitesales@permobil.com or 866.586.2413 / 509.586.2413



FRAME STYLE

- ATFS1** ☐ **Aero T Aluminum Frame** (Full Adjustability. 265 Pound Weight Limit. For Rear Seat Width **A** and Seat Depth **B** of 15" or greater, use K0005. For Rear Seat Width **A** and Seat Depth **B** of less than 15", use E1235.) **\$3,660**
- THDR1** ☐ **Heavy Duty Frame** (Not Available with TR Titanium Frame. 350 Pound Weight Limit. No charge when 20" Rear Seat Width **A** and 20" Seat Depth **B** are selected. Includes Tension Adjustable by Straps Back Upholstery. Requires Folding Aluminum Adjustable Height Backrest, Billet Aluminum Front Wheels, Shadow or Spinergy SPOX Rear Wheels. LX Rear Wheels only available with 0 or 2 degree camber. Not available with SuperLite Edition, Slipstream Fork, Velcro Adjustable Upholstery, Open Loop Footrest, V-Front Taper or Glide™ Forks.) **\$445**
- TQST1** ☐ **Q'Straight Adaptable Frame** (350 Pound Weight Limit. No charge when 20" Rear Seat Width **A** and 20" Seat Depth **B** are selected. Includes Tension Adjustable by Straps Back Upholstery. Requires Folding Aluminum Adjustable Height Backrest, Billet Aluminum Front Wheels, Shadow or Spinergy SPOX Rear Wheels. LX Rear Wheels only available with 0 or 2 degree camber. Not available with SuperLite Edition, Slipstream Fork, Velcro Adjustable Upholstery, Open Loop Footrest, V-Front Taper or Glide™ Forks.) **\$445**
- TRFR1** ☐ **Reinforced Frame** (No charge when 20" Rear Seat Width **A** and 20" **B** Seat Depth are selected.) **\$400**

FRAME FINISHES

- High Gloss Paint** ☐ TFF8 Super Wet Black
- Metallic Paint** ☐ TFF14 Ocean Blue Metallic

COLOR ANODIZE PACKAGE

(Color Anodize Package includes: Aluminum Camber Tube (Color Anodize not available with Power Adaptable Camber Tube), Camber Mount, Footrest Clamp, Caster Mount, Footrest Clamp, Caster Fork, Caster Cap and Agility Carbon Back clamps)

- ☐ TANO1 Black **STD**



TiFit: TAILORED INDIVIDUALIZED FRAME

REAR SEAT WIDTH A Measure from outside of seat tube at back post to the outside of the opposite seat tube at back post.

- TSW1 ☐ 12" ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19" **STD**
- TSW2 ☐ 20" **\$495**

FRONT SEAT WIDTH (SEAT TAPER) ① Measure from inside of front frame tube to inside of opposite front frame tube. (Required with V-Front End ①. Not Available with Standard Front End ①.)


- TFSW1 ☐ No Taper **STD**
- TFSW2 ☐ 1" (2.5cm) Taper (3" Narrower than Rear Seat Width **A**. Includes Tension Adjustable by Straps Seat Upholstery.) **\$280**
- TFSW4 ☐ 2" (5.1cm) Taper (4" Narrower than Rear Seat Width **A**. Includes Tension Adjustable by Straps Seat Upholstery.) **\$280**
- TFSW6 ☐ 3" (7.6cm) Taper (5" Narrower than Rear Seat Width **A**. Includes Tension Adjustable by Straps Seat Upholstery. Not available with 12" Seat Width.) **\$280**

SEAT DEPTH B Measure from front of back post to the front edge of seat sling. Frame bend will start .5" in front of seat sling unless Custom Frame Depth is specified, or the combination of specifications affect the rolling dynamics of the chair.

- TSD1 ☐ 12" ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19" **STD**
- TSD2 ☐ 20" **\$540**
- TSD3 ☐ Custom Frame Depth: _____" Frame Depth Longer or Shorter Than Seat Depth (Measure from front of seat sling to beginning of frame bend. Ranges from 0.25" to 3" and can be given as specific as 1/4" increments. A negative Custom Frame Depth or specified above 3" may be available through our Creative Mobility Division. A Custom Frame Depth of 20" will have an additional charge) **NCO**

ERGONOMIC SEAT ① Measure from front of back post to beginning of upward bend of seat tube.

- TERG2** ☐ **AeroT Ergonomic Seat** (Center of Gravity adjustment may be affected by this option. Not available with V-Front or Seat Taper, Seat Back that Locks When Folded, Desk Arms, Depth Adjustable Backrests, Fold Down Rigid Sideguards, Fendered Adjustable Sideguards, Solid Seat, Tension Adjustable Bolt On Seat Upholstery, or the Curve cushion.)
- ☐ 6" (15.2cm) (Minimum Seat Depth **B** 12") ☐ 7" (17.8cm) (Minimum Seat Depth **B** 12")
- ☐ 8" (20.3cm) (Minimum Seat Depth **B** 13") **\$520**



Permobil.com
2701 West Court Street | Pasco, Washington 99301
509.586.6117 Tel | 800.545.2266 Toll-Free
509.586.2413 Fax | 866.586.2413 Toll-Free Fax

OF0029c T Series 0719RevAADL

TiLite T-Series Canadian Order Form - AADL #W395

Effective October 1, 2021 (New Items in Red)

DUAL-TUBE RIGID FRAME

Page 2 of 7

TiFit: TAILORED INDIVIDUALIZED FRAME

FRONT SEAT HEIGHT ● Measure from floor to top of seat tube at front edge of seat sling.

TFSH1 ☐ 15" (38.1cm) ☐ 16" (40.6cm) ☐ 16.5" (41.9cm) ☐ 17" (43.2cm) ☐ 17.5" (44.5cm)
☐ 18" (45.7cm) ☐ 18.5" (47cm) ☐ 19" (48.3cm) ☐ 19.5" (49.5cm) ☐ 20" (50.8cm)
☐ 20.5" (52.1cm) ☐ 21" (53.3cm) **STD**

REAR SEAT HEIGHT ● Measure from floor to top of seat tube at back post.

Rear Wheel Size	Rear Seat Height
20"	12.5"-21" (31.8-53.3 cm)
22"	13.5"-21" (34.3-53.3 cm)
24"	14.5"-21" (15" on AeroT/TRA) (36.8-53.3 cm)
25"	15.0"-21" (38.1-53.3 cm)
26"	15.5"-21" (38.1-53.3 cm)

TRSH1 ☐ 12.5" (31.8cm) ☐ 13" (33cm) ☐ 13.5" (34.3cm) ☐ 14" (35.6cm) ☐ 14.5" (36.8cm)
☐ 15" (38.1cm) ☐ 15.5" (39.4cm) ☐ 16" (40.6cm) ☐ 16.5" (41.9cm) ☐ 17" (43.2cm)
☐ 17.5" (44.5cm) ☐ 18" (45.7cm) ☐ 18.5" (47cm) ☐ 19" (48.3cm) ☐ 19.5" (49.5cm)
☐ 20" (50.8cm) **STD**

FRONT ANGLE ● Measure from frame front to floor.

TFA1 ☐ 70° ☐ 75° ☐ 80° ☐ 85° (Not available with 6" Front Wheels.) ☐ 90° (Not available with 5" or 6" Front Wheels.) **STD**

SEAT TO FOOTREST ● /FOOTRESTS Measure from front edge of seat sling to top rear of footrest.

Standard Footrests Measurement must be at least 2" shorter than Front Seat Height ●.

☐ 13.5" (34.3cm) ☐ 14" (35.6cm) ☐ 14.5" (36.8cm) ☐ 15" (38.1cm) ☐ 15.5" (39.4cm) ☐ 16" (40.6cm)
☐ 16.5" (41.9cm) ☐ 17" (43.2cm) ☐ 17.5" (44.5cm) ☐ 18" (45.7cm) ☐ 18.5" (47cm) ☐ 19" (48.3cm)

(Due to current supplier issues, titanium footrest options will use aluminum tubing instead of titanium and will be powder coated Black Lava. During this time, these options will be available in aluminum at no charge to our customers. If a titanium footrest is desired at a later point, the option will be offered at on-chair pricing when available.)

TFTR3 ☐ Angle Adjustable Footrest **STD**
 TFTR1 ☐ Titanium Footrest w/Flat ABS Cover **\$200**
 TFTR2 ☐ Titanium Open Loop Footrest **\$255**
 TFTR12 ☐ Flip-Back Footrest (Not available with ABS Cover with 2" Sides or Front Roller.) **\$235**
High Mount Footrest
☐ 7"-9" ☐ 9"-11" ☐ 11"-13" Set at _____" (Use 1/2" increments)
 TFTR4 ☐ Angle Adjustable High Mount Footrest **\$145**

FOOTREST WIDTH ● or ▼ Select either a Standard Front End ● if no Seat Taper ● is selected or a V-Front End ▼ if a Seat Taper ● is selected.





STANDARD FRONT END ● Find your Rear Seat Width ● in the Chart below and select one of the options to the right.
 Measure from inside of front frame tube to inside of opposite front frame tube.

TFE1 ☐ NO V-Front Taper **STD**

REAR SEAT WIDTH ● (Outside Measurement)	FOOTREST WIDTH ●		
	No Taper (2" Narrower than Rear Seat Width)	2" Taper (4" Narrower than Rear Seat Width)	4" Taper (6" Narrower than Rear Seat Width)
12"	☐ 10"	☐ 8"	
13"	☐ 11"	☐ 9"	
14"	☐ 12"	☐ 10"	☐ 8"
15"	☐ 13"	☐ 11"	☐ 9"
16"	☐ 14"	☐ 12"	☐ 10"
17"	☐ 15"	☐ 13"	☐ 11"
18"	☐ 16"	☐ 14"	☐ 12"
19"	☐ 17"	☐ 15"	☐ 13"
20"	☐ 18"	☐ 16"	☐ 14"

V-FRONT END ▼ Measure from inside of front frame tube to inside of opposite front frame tube 2-1/2" above footrest. (Titanium Footrests have 1" of adjustment up and down. Not available with Ergonomic Seat frames.)

TFE2 ☐ V-Front Taper **NCO**
☐ 9" (22.9cm) ☐ 10" (25.4cm) ☐ 11" (27.9cm)
☐ 12" (30.5cm) ☐ 13" (33cm) ☐ 14" (35.6cm)
☐ 15" (38.1cm)



TiLite T-Series Canadian Order Form - AADL #W395

Effective October 1, 2021 (New Items in Red)

DUAL-TUBE RIGID FRAME

Page 3 of 7

TiFit: TAILORED INDIVIDUALIZED FRAME

SEAT BACK TYPE/PUSH HANDLES/SEAT BACK HEIGHT ① Measure from top of back post to top of seat tube at rear of frame.

- TSB1 ☐ Folding Aluminum Adjustable Height w/Precision Lock Back **STD**
- TSB2 ☐ Folding Titanium Adjustable Height w/Precision Lock Back and Rounded Rigidizer Bar (Titanium back will match Frame Finish. Not available with Bolt-On Push Handles.) **\$295**
- TSB15 ☐ Non-Folding, Fixed Height Titanium with New Rounded Rigidizer Bar (Only available with TR. Titanium back will match Frame Finish. Not available with Push Handles.) **\$325**

TISHAFT™ BACK RELEASE BARS (Not available with Non-Folding, Fixed Height Titanium Back.)

- TSB3 ☐ Aluminum **STD**
- TSB5 ☐ Carbon Fiber **\$85**

PUSH HANDLES

- TSB6 ☐ Integrated Push Handles **NCO**
- TSB7 ☐ Integrated Fold-Down Push Handles (Requires Tension Adjustable by Straps Back Upholstery.) **\$115**
- TSB8 ☐ Bolt-On Push Handles (Includes 4" Deep Backrest Rigidizer Bar. Not available with Ti Back.) **\$200**
- TSB9 ☐ Bolt-On Height-Adjustable Push Handles (Includes 4" Deep Backrest Rigidizer Bar. Not available with Ti Back.) **\$410**

FOLDING SEAT BACK HEIGHT ② (TSB1 and TSB2)

Push Handles (Selected Above)	Short	Medium	Tall	Set at: (Use 1/2" increments)
None/ Bolt-On	☐ 9.5"-12.5"	☐ 13"-17.5"	☐ 17"-21.5"	Set at: _____
Integrated	☐ 11.5"-13.5"	☐ 14.5"-17.5"	☐ 18.5"-21.5"	Set at: _____
Integrated Fold-Down	☐ 14"-16.5"	☐ 15"-19.5"	☐ 19.5"-22"	Set at: _____

NON-FOLDING, FIXED-HEIGHT ③ TITANIUM SEAT BACK with New Rounded Rigidizer Bar (TSB15) (Titanium back will match Frame Finish. Required for fixed Side Guards. Not available with Swing Away Tubular or Swing Away - Flip Back Tubular Armrests and Integrated Push Handles, Integrated Fold-Down Push Handles or Bolt-On Height Adjustable Push Handles.)

- ☐ 5" (12.7cm) ☐ 5.5" (14cm) ☐ 6" (15.2cm) ☐ 6.5" (16.5cm) ☐ 7" (17.8cm) ☐ 7.5" (19.1cm) ☐ 8" (20.3cm)
- ☐ 8.5" (21.6cm) ☐ 9" (22.9cm) ☐ 9.5" (24.1cm) ☐ 10" (25.4cm) ☐ 10.5" (26.7cm) ☐ 11" (27.9cm) ☐ 11.5" (29.2cm)
- ☐ 12" (30.5cm) ☐ 12.5" (31.8cm) ☐ 13" (33cm) ☐ 13.5" (34.3cm) ☐ 14" (35.6cm) ☐ 14.5" (36.8cm) ☐ 15" (38.1cm)
- ☐ 15.5" (39.4cm) ☐ 16" (40.6cm) ☐ 16.5" (41.9cm) ☐ 17" (43.2cm) ☐ 17.5" (44.5cm) ☐ 18" (45.7cm) ☐ 18.5" (47cm)
- ☐ 19" (48.3cm) ☐ 19.5" (49.5cm) ☐ 20" (50.8cm)

AVAILABLE ON ALL BACKS

- TSB11 ☐ Backrest that Locks when Folded (Not available with Ergonomic Seat.) **STD**
- TSB12 ☐ 4" Deep Backrest Rigidizer Bar (Currently available with Aluminum backrest option only. Not available with a Non-Folding, Fixed Height Titanium Seat Back Type. 4" Deep Rigidizer Bar is NOT a rounded style rigidizer bar.) **NCO**

SEAT BACK ANGLE ④ Measure angle from front of back post to floor.




- TSBA1 **Folding Back**
- ☐ 80° ☐ 83° ☐ 86° ☐ 89° ☐ 92° ☐ 95° ☐ 98° ☐ 101° **STD**
- TSBA1 **Non-Folding Back**
- ☐ 80° ☐ 81° ☐ 82° ☐ 83° ☐ 84° ☐ 85° ☐ 86° ☐ 87° ☐ 88° ☐ 89° ☐ 90°
- ☐ 91° ☐ 92° ☐ 93° ☐ 94° ☐ 95° ☐ 96° ☐ 97° ☐ 98° ☐ 99° ☐ 100° ☐ 101° **NCO**

CENTER OF GRAVITY ⑤ Measure from front of back post at the seat tube to the center of the rear axle. If no center of gravity is selected, chair will ship with the center of gravity in the rearward-most position.

- TCOG1 ☐ 1" (2.5cm) ☐ 1.25" (3.2cm) ☐ 1.5" (3.8cm) ☐ 1.75" (4.4cm) ☐ 2" (5.1cm)
- ☐ 2.25" (5.7cm) ☐ 2.5" (6.4cm) ☐ 2.75" (7cm) ☐ 3" (7.6cm) ☐ 3.25" (8.3cm)
- ☐ 3.5" (8.9cm) ☐ 3.75" (9.5cm) ☐ 4" (10.2cm) ☐ 4.25" (10.8cm) ☐ 4.5" (11.4cm)
- ☐ 4.75" (12.1cm) ☐ 5" (12.7cm) **STD**

REAR WHEEL SPACING ⑥ /PLUGS Measure from outside of seat back post to inside of rear tire. 0.5" Rear Wheel Spacing on TR not available with 0° or 2° Camber. TRA 0.75" Rear Wheel Spacing not available with 0° or 2° Camber

- TRWS1 ☐ .75" ☐ 1" ☐ 1.25" **STD**
- TCBR12 ☐ Adjustable Rear Wheel Spacing Camber Plugs (Not available with Power Adaptable Frame.) **STD**
- TCBR9 ☐ Additional Spacers to Increase Rear Wheel Spacing **NCO**
- TCBR11 ☐ Zero Play Aluminum Camber Plugs (Not available with Power Adaptable Frame or more than 6° of Camber ⑦.) **NCO**



TiLite T-Series Canadian Order Form - AADL #W395

Effective October 1, 2021 (New Items in Red)

DUAL-TUBE RIGID FRAME

Page 4 of 7

TiFit: TAILORED INDIVIDUALIZED FRAME

CAMBER ☉ /CAMBER TUBES

TCBR1	☐ 0° ☐ 2° ☐ 4° ☐ 6° (8° of Camber ☉ may be available through our Creative Mobility Division.)	STD
TCBR15	☐ Aluminum Camber Tube	STD
TCBR3	☐ Titanium Camber Tube (Not available with Power Adaptable Frame.)	\$225

CUSTOM CONFIGURATION OPTIONS

FRONT WHEELS

TFW24	☐ 3"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire ☐ Silver Hub OR ☐ Black Hub	\$180
TFW14	☐ 4"x.75" Plastic Wheel w/Poly Tire (Not available with Slipstream Fork, Glide™ Single-Sided Fork or Power Adaptable Frame.)	NCO
TFW3	☐ 4"x.75" Light-Up Micro-Caster (Not available with Slipstream Fork, Glide™ Single-Sided Fork or Power Adaptable Frame.)	\$80
TFW25	☐ 4"x.75" LiteSpeed Billet Aluminum Wheel w/Poly Tire ☐ Silver Hub OR ☐ Black Hub	\$150
TFW26	☐ 4"x1.5" LiteSpeed Plastic Wheel w/Soft Roll Tire (Not available with Power Assist, Power Add-on System or Heavy Duty Frame.)	\$95
TFW27	☐ 4"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire ☐ Silver Hub OR ☐ Black Hub	\$190
TFW6	☐ 5"x1" Plastic Wheel w/Poly Tire (Not available with 90° Front Angle ☉, Slipstream Fork, Glide™ Single-Sided Fork or Power Adaptable Frame.)	NCO
TFW28	☐ 5"x1.5" LiteSpeed Plastic Wheel w/Soft Roll Tire (Not available with 90° Front Angle ☉ or Power Adaptable Frame.)	\$95
TFW30	☐ 5"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire ☐ Silver Hub OR ☐ Black Hub (Not available with 90° Front Angle ☉.)	\$250
TFW9	☐ 6"x1" Plastic Wheel w/Poly Tire (Not available with 85° or 90° Front Angle ☉, Slipstream Fork, Glide™ Single-Sided Fork or Power Adaptable Frame.)	NCO
TFW32	☐ 6"x1.5" LiteSpeed Plastic Wheel w/Soft Roll Tire (Not available with 85° or 90° Front Angle ☉, Slipstream Fork, Glide™ Single-Sided Fork or Power Adaptable Frame.)	\$100
TFW33	☐ 6"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire ☐ Silver Hub OR ☐ Black Hub (Not available with 85° or 90° Front Angle ☉, Slipstream Fork, Glide™ Single-Sided Fork.)	\$250

PRECISION LOCK FRONT FORKS

TFK1	☐ TiLite Standard Fork	STD
TFK2	☐ TiLite Slipstream Single-Sided Fork (Not available with Power Adaptable Frame.)	\$185
TGL1	☐ Out-Front Glide™ Suspension Fork (Black anodize only.)	\$545
	Rider's Weight: ☐ 125 lbs. ☐ 125-180 lbs. ☐ 180-265 lbs.	
TGL2	☐ Out-Front Glide™ Single-Sided Suspension Fork (Black anodize only. Not available with 6" Front Wheels.)	\$680
	Rider's Weight: ☐ 125 lbs. ☐ 125-180 lbs. ☐ 180-265 lbs.	
TTFS1	☐ Titanium Fork Stem	\$60

STANDARD REAR WHEELS (See Rear Seat Height Ranges Chart for Rear Seat Height ☉ limitations.)

TRW14	TiLite SHADOW™	STD
	☐ 20"(451)	
	☐ 22"(501)	
	☐ 24"(540)	
	☐ 25"(559) (Requires 15" minimum Seat Depth ☉.)	
	☐ 26"(590) (Requires 15" minimum Seat Depth ☉.)	

PERFORMANCE REAR WHEELS (See Rear Seat Height Ranges Chart for Rear Seat Height ☉ limitations.)

TRW7	Spinergy SPOX ☐ Black ☐ Red ☐ Yellow ☐ Blue ☐ White	\$900
	☐ 22"(501) - Std. Hub	
	☐ 24"(540) - Std. Hub	
	☐ 25"(559) - Std. Hub (Requires 15" minimum Seat Depth ☉.)	
	☐ 26"(590) - Sport Hub (Requires 15" minimum Seat Depth ☉. Not available with Quad Release Axles.)	
TRW8	Spinergy LX ☐ Black ☐ Red ☐ Yellow ☐ Blue ☐ White ☐ Go Green ☐ Bright Orange ☐ Flamingo Pink	\$1,055
	☐ 22"(501) - Std. Hub	
	☐ 24"(540) - Std. Hub	
	☐ 25"(559) - Std. Hub (Requires 15" minimum Seat Depth ☉.)	
	☐ 26"(590) - Sport Hub (Requires 15" minimum Seat Depth ☉. Not available with Quad Release Axles.)	
TRW16	Spinergy Carbon Blade (Not available with Heavy Duty Frame, Natural-Fit or Natural-Fit LT with Standard Grip or Super Grip.)	\$1,540
	☐ 24"(540)	
	☐ 25"(559) (Requires 15" minimum Seat Depth ☉.)	
TRW17	X-Core Mag (Not available with Heavy Duty Frame, Short Tabs, Quad Release Axles, Spoke Guards or Omit Handrims Option. Requires Aluminum Silver Anodized Handrims.)	\$110
	☐ 20"(451)	
	☐ 22"(501)	
	☐ 24"(540)	

AXLES

TAXL1	☐ Stainless Quick Release	STD
TAXL2	☐ Titanium Quick Release	\$290
TAXL3	☐ Quad Release (Not available with Mag or 26"(590) Spinergy Spox or LX Rear Wheels.)	\$170

TiLite T-Series Canadian Order Form - AADL #W395

Effective October 1, 2021 (New Items in Red)

DUAL-TUBE RIGID FRAME

Page 5 of 7

CUSTOM CONFIGURATION OPTIONS

PNEUMATIC REAR TIRES

TRTR1	☐ Treaded (Not available with 25"(559) or 20" (451) Rear Wheels.)	STD
TRTR2	☐ Airless Inserts (Not available with Spinergy Carbon Blade Rear Wheels.)	\$90
TRTR3	☐ 25"(559) Primo Treaded (Standard with 25"(559) Rear Wheels.)	STD
TRTR4	☐ Primo Slick	NCO
TRTR30	☐ 24"(540), 25"(559) and 26"(590) Primo Sentinel	\$160
TRTR8	☐ 24"(540), 25"(559) and 26"(590) Schwalbe® Marathon Plus Evolution	\$155
TRTR9	☐ 24"(540) and 25"(559) Schwalbe® Slick	\$150
TRTR24	☐ 24"(540) Kenda Kobra Knobby	\$95

SOLID REAR TIRES (Not available with Spinergy Carbon Blade Rear Wheels.)

TRTR22	☐ 20"(451), 22"(501), 24"(540) and 25"(559) Primo Express	\$110
TRTR23	☐ 24"(540) and 25"(559) Primo Active	\$135

HANDRIMS

Handrims by OUT-FRONT®

TRIM1	☐ Aluminum – Silver Anodized	STD
TRIM21	☐ Aluminum – Black Anodized (Not available with 20"(451) Rear Wheels.)	\$60
TRIM12	☐ Titanium – Satin Finish (Not available with 20"(451) or 22"(501) Rear Wheels, or 24"(540) Knobby tires.)	\$490
TRIM10	☐ The Surge (All-In-One Oval with Gription™ Strip, Larger Profile) (Not available with 20"(451) Rear Wheels.)	\$435
TRIM11	☐ The Surge LT (All-In-One Oval with Gription™ Strip, Smaller Profile) (Not available with 20"(451) Rear Wheels.)	\$470
TRIM20	☐ Q-Grip (Exclusive High-Friction, Anti-Slip Coating) (Not available with 20"(451) Rear Wheels.)	\$375

The Natural-Fit® (Original, Larger Profile) (Not available with 20"(451) Rear Wheels.)

TRIM14	☐ Standard Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.)	\$415
TRIM15	☐ Super Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.)	\$470
TRIM16	☐ Without Thumb Grip	\$380

The Natural-Fit® LT (Reduced Weight, Smaller Profile) (Not available with 20"(451) Rear Wheels.)

TRIM17	☐ Standard Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.)	\$500
TRIM18	☐ Super Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.)	\$535
TRIM19	☐ Without Thumb Grip	\$460

Additional Handrim Options

TRIM5	☐ Plastic Coated – Black	\$160
TRIM22	☐ Omit Handrims	NCO

Tab Length

TRIM7	☐ Long Tabs	STD
TRIM8	☐ Short Tabs (Not available with Mag Rear Wheels or The Natural-Fit® Handrims with Thumb Grip.)	\$35

WHEEL LOCKS

Wheel Locks by OUT-FRONT®

TWLK7	☐ Compact Composite Push to Lock (Available with .75" to 1.5" Rear Wheel Spacing ① only.)	NCO
TWLK8	☐ Compact Composite Pull to Lock (Available with .75" to 1.5" Rear Wheel Spacing ① only.)	NCO
TWLK10	☐ Extension Handles for Compact Composite Push/Pull to Lock	\$85
TWLK9	☐ Composite Scissor Lock	\$175

Additional Wheel Locks

TWLK3	Uni-Lock Set At: ☐ Push OR ☐ Pull	\$165
TWLK4	☐ Extension Handles for Uni-Lock	\$85
TWLK15	☐ Grade Aids Push to Lock (Not available with solid rear tires)	\$230

BACK UPHOLSTERY (Want more seating options? View the TiLite Seating Order Form on our website.)

TBUP2	☐ Velcro Adjustable (Not available with Fold-Down Push Handles.)	STD
-------	---	-----

Tension Adjustable by Straps (Required with Fold-Down Push Handles.)

TBUP3	☐ Padded Nylon	\$415
TBUP5	☐ AIR Breathable Upholstery (Not available with Ultrasuede Backrest Stripe.)	\$400
TUPC3	Ultrasuede Backrest Stripe ☐ Black	\$105
TBUP0	☐ Omit Back Upholstery	NCO

SEAT UPHOLSTERY

TSUP1	☐ Tension Adjustable Bolt On (Not available with Ergonomic Seat or Tapered Seat ①. Black only.)	STD
TSUP2	☐ Tension Adjustable by Straps (Standard with Ergonomic Seat and FSW (Seat Taper) ① of 1", 2" or 3". Black upholstery only.)	\$185
TSUP8	☐ Aluminum Solid Seat (Not available with 10" or 11" Rear Seat Width ② or Ergonomic Seat ②.)	\$285
TSUP0	☐ Omit Seat Upholstery	NCO

SEAT CUSHION (Want more seating options? View the TiLite Seating Order Form on our website.)

TCSH2	☐ Comfort Curve Cushion	\$155
-------	--	-------

SIDE GUARDS

TSDG1	☐ Cloth (Not available with Ride Designs Corbac Back Support.)	\$120
-------	---	-------

TiLite T-Series Canadian Order Form - AADL #W395

Effective October 1, 2021 (New Items in Red)

DUAL-TUBE RIGID FRAME

Page 6 of 7

CUSTOM CONFIGURATION OPTIONS

SIDE GUARDS (CONTINUED)

Aluminum

TSDG2	☐ Rigid Removable	☐ Small	OR	☐ Std		\$290
TSDG3	☐ Quad Removable w/Hole	☐ Small	OR	☐ Std		\$315
TSDG4	☐ Fold Down Rigid	☐ Small	OR	☐ Std	(Not available with Ergonomic Seat.)	\$375
TSDG6	☐ Fendered Adjustable	(Available with 24" and 25" Rear Wheels. Requires Rear Seat Height ① of 14.5" to 18.5".)				\$400

Carbon Fiber

TSDG7	☐ Rigid Removable	☐ Small	OR	☐ Std		\$340
TSDG8	☐ Quad Removable w/Hole	☐ Small	OR	☐ Std		\$375
TSDG9	☐ Fold Down Rigid	☐ Small	OR	☐ Std	(Not available with Ergonomic Seat or Carbon Edition.)	\$510
TSDG11	☐ Fendered Adjustable	(Available with 24" and 25" Rear Wheels. Requires Rear Seat Height ① of 14.5" to 18.5".)				\$565
TSDG12	☐ Fendered Removable					\$510

Available on Aluminum and Carbon Fiber Non-Fendered Side Guards

TSDG13	☐ Side Guards Custom Cut to Profile of Rear Wheels	(Not available with Cloth.)		\$315
--------	---	-----------------------------	-------	--------------

ARMRESTS (Desk Arms come with built-in Side Guards.)

TARM1	Swing Away Tubular ☐ 6.5" to 9.5" OR ☐ 9.5" to 12.5" (Not Available with Non-Folding, Fixed Height Titanium Seat Back or Height Adjustable Push Handles.)	\$260
TARM4	☐ Swing Away - Flip Back Tubular 8.5" to 11.5" (Not Available with Non-Folding, Fixed Height Titanium Seat Back.)	\$315
TARM2	Removable Height Adjustable Desk Arm with Rigid Side Guard (Requires 14" minimum Seat Depth ①. Not available with Ergonomic Seat.)	\$360
	☐ 9" to 12" OR ☐ 12" to 15"	
	☐ Standard Length Pad (10") OR ☐ Full Length Pad (14")	
TARM3	Removable Height Adjustable Desk Arm with Rigid Side Guard and Transfer Loop (Requires 15" minimum Seat Depth ①. Not Available with Ergonomic Seat.)	\$375
	☐ 9" to 12" OR ☐ 12" to 15"	
	☐ Standard Length Pad (10") OR ☐ Full Length Pad (14")	

Available Armrest Option

TARM5	☐ Covered in Black Ultrasuede	(Not available with Desk Arms)		\$115
-------	--	--------------------------------	-------	--------------

SEAT BELTS

TBLT1	☐ Velcro Adjustable, Black Nylon Webbing		\$55
TBLT2	☐ Bodypoint® 2" Push-Button Buckle (Auto-Style.)		\$95
TBLT4	☐ Bodypoint® 1" Wide Push-Button Buckle Padded (Auto-Style.)	(Not available with 19" or 20" Rear Seat Widths ②.)	\$170
TBLT5	☐ Bodypoint® 1.5" Wide Push-Button Buckle Padded (Auto-Style.)	(Not available with 19" or 20" Rear Seat Widths ②.)	\$170
TBLT10	☐ Bodypoint® Optional Band Clamp Mounting Hardware		\$20
TBLT15	☐ Bodypoint® EvoFlex® Pelvic Stabilizer Swivel		\$185
TBLT16	☐ Bodypoint® EvoFlex® Pelvic Stabilizer Push Button		\$185
TBLT17	☐ Bodypoint® EvoFlex® Pelvic Stabilizer Rehab Latch		\$185

CALF STRAP

TCLF2	Velcro Adjustable		\$45
	☐ Xtra Small (Recommended for 6" - 9" Footrest Widths.)		
	☐ Small (Recommended for 9" - 12" Footrest Widths.)		
	☐ Medium (Recommended for 11" - 15" Footrest Widths.)		
	☐ Large (Recommended for 14" - 17" Footrest Widths.)		
	☐ Xtra Large (Recommended for 16" - 19" Footrest Widths.)		
TCLF3	Bodypoint® Padded Velcro Adjustable		\$80
	☐ Medium (Recommended for 11" - 15" Footrest Widths.)		
	☐ Large (Recommended for 16" - 19" Footrest Widths.)		
TCLF1	☐ Buckle Adjustable, Black Nylon Webbing		\$30

RESIDUAL LIMB SUPPORT

TRSL1	☐ Residual Limb Support		\$515
TRSL2	☐ Residual Limb Support Pair (Both Left and Right Pad information required)		\$1,030
	Left Pad Width ☐ 8" ☐ 9" ☐ 10"	Right Pad Width ☐ 8" ☐ 9" ☐ 10"	
	Left Pad Length ☐ 10" ☐ 14"	Right Pad Length ☐ 10" ☐ 14"	

BACKPACK

TBCK3	☐ Black Permobil Backpack		\$95
-------	--	-------	-------------

SEAT POUCH

TPCH3	☐ Black Permobil Seat Pouch		\$60
TPCH4	☐ Bodypoint® Mobility Bag		\$130

TiLite T-Series Canadian Order Form - AADL #W395

Effective October 1, 2021 (New Items in Red)

DUAL-TUBE RIGID FRAME

Page 7 of 7

CUSTOM CONFIGURATION OPTIONS

SPOKE GUARDS

TGRD1 ☐ Clear Lexan with TiMotion Logo (Not available with Mag Rear Wheels.) \$135

ANTI-TIPS

TTIP1 ☐ Rear Aluminum - Black \$195

TTIP7 ☐ Black Aluminum Single \$100

Mount on ☐ Left Side or ☐ Right Side NCO

TTIP2 ☐ User Friendly Flip-up (Not available with Frog Legs® Rear Suspension.) \$360

ZTIP8 ☐ User Friendly Flip-up Single \$170

Mount on ☐ Left Side or ☐ Right Side NCO

ACCESSORIES

TUEST1 Upper Extremity Support Tray ☐ 12" to 14" OR ☐ 13" to 15" Seat Widths (Requires Armrests.) \$290


TLUG1 ☐ Luggage Carriers \$145

TMPCT2 Neoprene Impact Guards with Ultrasuede Stripes ☐ Black \$105

TCH1 ☐ Crutch Holder \$145

FRAME NOTES

Aero T will always be a Rounded Rear Frame.



A. Round Rear Frame

GENERAL

Permobil manufactures TiLite manual wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. HCPCS Codes are for reference purposes only. Actual coding and coverage varies by payer. Proper HCPCS coding is not a guarantee of payment. Please verify the appropriate payer's coverage policy.

Get the latest order forms at permobilus.com. Prices, standard features, options and accessories are subject to change without notice. All rights reserved. Trademarks are identified by symbols™ and ®. TiLite, TiFit, the TiLite logo, Ti, the Ti logo, The Titanium Experts, and TiShaft are registered trademarks of TiSport, LLC and are licensed to Permobil. Permobil® is a registered trademark of Permobil. Adobe®, Acrobat®, Acrobat Professional® and Reader® are registered trademarks of Adobe Systems Inc. Bodypoint® is a registered trademark of Bodypoint Inc. Kenda® is a registered trademark of Kenda Rubber Industrial Co. Ltd. Ride Designs® and Corbac® are registered trademarks of Aspen Seating, LLC. ROHO®, Low Profile® and Quadro Select® are registered trademarks of ROHO Inc. Schwalbe® is a registered trademark of Ralf Bohle GmbH. Supracor®, Stimulite® and Corbee® are registered trademarks of Supracor, Inc. The Natural-Fit® and Out-Front® are registered trademarks of Three Rivers Holdings, LLC. Varilite® is a registered trademark of Cascade Designs, Inc.

NOTES