

# Tricycle features



**Variable pedal resistance**  
The adjustable brake can provide drag, or lock the wheels for safe transfer.



**Backrest angle**  
Position the backrest vertically or at a 7½° or 15° recline.



**Carbon fiber drive belt**  
The drive belt will not rust or stretch, and will last at least five times longer than a bicycle chain.



**Seat depth**  
The seat depth adjusts easily with a knob underneath.

# What's new?



New colors are available for the large Trike, including teal, black and white.



## Unique and new features

- 310 lb weight limit—increased from 200 lb
- The seat height adjusts on a low-friction track, and can be adjusted while in use, depending on the rider's weight and ability.
- Pedals adjust forward and backward on pedal crank to position rider's feet optimally.
- Variable pedal resistance to help riders moderate their speed, and for stationary cycling.
- The sandals can be removed to reveal standard bike pedals underneath.
- Durable, maintenance-free carbon fiber drive belt.
- Disassembles easily for transport.

