

Permobil

M4 Nexa[®] Jr



permobil

Made to grow with you

M4 Nexa Jr brings the next generation of Permobil power mobility to children, young adults, and people who need a smaller fit. It gives users with smaller body dimensions the access, support, and confidence they need in everyday environments, without compromising power functions, seating capability, or performance.

Small in size. Big in capability.

With a compact, maneuverable base, a seat width starting at just 30 cm, and a seat-to-floor height down to 39 cm, M4 Nexa Jr helps users move more easily through classrooms, corridors, doorways, bathrooms, and busy social spaces. Its lower seat height and tuck-under footplate design support closer access to tables, vehicles, people, and the activities that matter. Quieter driving and braking help create a smoother, more discreet mobility experience throughout the day.

Built to fit you today — and adapt with you tomorrow

The Nexa biomechanical seating system is designed to maintain a supported position with minimal displacement as power seat functions are used. With adjustable seat width and depth, M4 Nexa Jr can adapt as the user grows or as posture, function, or lifestyle changes over time. Armrest swivel and side/out functions also support individual positioning needs, including a midline joystick position for children and easier transfers.

Easier fitting

With fewer components, protected cabling, improved battery access, and tool-efficient adjustments, M4 Nexa Jr helps simplify fitting, reconfiguration, and service. Seat depth can be adjusted while the user remains seated, helping reduce interruptions for users and families as needs change.

Connectivity for more control

With Permobil +Connect, M4 Nexa supports remote insights, proactive service, and wireless programming. MyPermobil app provides useful information such as battery range, distance, seating behavior, and the ability to personalise LED lights with animated themes.





Specifications

Max. user weight	136 kg
Max. speed	10 km/h
Battery	52 Ah/68 Ah
Range (52 Ah battery)	33 km*
Weight incl. batteries 52 Ah/68 Ah	144 kg (52 Ah) 160.5 kg (68 Ah)
Obstacle climbing	70 mm
Ground clearance	74 mm
Base width	540 mm
Base length	850 mm
Seat depth	250–560 mm
Seat width	300–560 mm
Seat-to-floor height	390/430 mm
Back frame height	Low 300–410 mm** Standard 460–560 mm
R-net electronics Power Platform	120 A
Crash test ISO 7176-19:2022 tie-down	Yes–max 136 kg user weight
Tested according to EN12184:2022 CE	Yes

*Driving range (Actual driving range will vary based on driving, battery conditions and terrain)

**To fit back support sizes

Measurements are approximates and subject to change without notice.

Accent colours

BLACK	DARK BLUE	RED	GRAY	LIGHT BLUE	ORANGE
WHITE	DARK PURPLE	LIGHT PURPLE	YELLOW	PINK	GREEN

Power functions



Anterior tilt
Active Reach 10°



Seat elevator
Active Height
300 mm



Power leg rest
85–175°



Posterior tilt
50°



Recline
150°

