

Fitting / Order form


F3 Corpus (Model 2019)



F3 Corpus

Front wheel drive chair with independent 4 wheel suspension. Standard equipped with airfilled tires. Corpus seating system with ergonomic seat and backrest cushions. De F3 Corpus is prepared to be equipped with powered backrest, legrest, seat elevation and tilt.

Dealer information	
Customer (name / code)	/
Contact	
Address	
Post code / City	/
Phone / Fax	/
Client Information	
Name	/
Address	
Post code / City	/
Phone	
E-mail	
Diagnosis	



Anatomic Dimension



Weight:	kg	3:	cm	7:	cm	
Height:	cm	4:	cm	8:	cm	
	1:	cm	5:	cm	9:	cm
	2:	cm	6:	cm	10:	cm

Fitting date							<input checked="" type="checkbox"/> Required choice <input type="checkbox"/> / <input type="radio"/> Optional choice <input type="checkbox"/> / <input checked="" type="radio"/> Default required option		
Article code	Position	Price	Funding	Additional					
	Left Right	excl VAT	incl VAT	code	Info				



Chassis						
F3 CORPUS 3G R-NET	108398-99-0	<input checked="" type="checkbox"/>				
MAX SPEED 10KMH	107305-99-0	<input checked="" type="radio"/>				
MAX SPEED 8KMH	107628-99-0	<input type="radio"/>				
Language user manuals						
LANGUAGE						
ESP						
ESP	107699-99-0	<input checked="" type="checkbox"/>				
Connect-Me						
CONNECTME 4G EU, ACTIVATE D	109465-99-0	<input checked="" type="checkbox"/>				
Mud guards rear						
MUDGUARD KIT	108445-99-0	<input checked="" type="checkbox"/>				
Accent colour						
Chassis / Wheels (no charge)						
COVER ACCENT COLOUR AQUA	108419-99-0	<input type="radio"/>	<input type="radio"/>			
COVER ACCENT COLOUR SKY	108420-99-0	<input type="radio"/>	<input type="radio"/>			
COVER ACCENT COLOUR VIOLET	108421-99-0	<input type="radio"/>	<input type="radio"/>			
COVER ACCENT COLOUR FLAMINGO	108422-99-0	<input type="radio"/>	<input type="radio"/>			
COVER ACCENT COLOUR CHILI	109328-99-0	<input type="radio"/>	<input type="radio"/>			

COVER ACCENT COLOUR CANARY	108424-99-0	<input type="radio"/>	<input type="radio"/>			
COVER ACCENT COLOUR SPRING	108425-99-0	<input type="radio"/>	<input type="radio"/>			
COVER ACCENT COLOUR ARCTIC	108426-99-0	<input type="radio"/>	<input type="radio"/>			
COVER ACCENT COLOUR SUNSET	108427-99-0	<input type="radio"/>	<input type="radio"/>			
COVER ACCENT COLOUR MIDNIGHT	108428-99-0	<input type="radio"/>	<input type="radio"/>			
COVER ACCENT COLOUR GOLD	108429-99-0	<input type="radio"/>	<input type="radio"/>			
COVER ACCENT COLOUR BRONZE	108430-99-0	<input type="radio"/>	<input type="radio"/>			
COVER ACCENT COLOUR PLATINUM	108431-99-0	<input type="radio"/>	<input type="radio"/>			
Drive wheels						
WHEELS 3.00x8" AIR FILLED TIRES BLACK	108411-99-0	<input checked="" type="radio"/>				 108412-99-0  108411-99-0
WHEELS 3.00x8" SOLID TIRES BLACK	108414-99-0	<input type="radio"/>				108414-99-0
WHEELS 3.00x8" AIR FILLED TIRES/SNOW AND MUD BLACK	108412-99-0	<input type="radio"/>				
Caster wheels						
WHEELS 210X65 AIR FILLED TIRES BLACK	108408-99-0	<input checked="" type="radio"/>				
WHEELS 210X65 SOLID TIRES BLACK	108409-99-0	<input type="radio"/>				
Batteries						
BATTERIES, 65Ah, SET OF 2 PERMOBIL	109183-99-0	<input checked="" type="radio"/>				
BATTERIES, 85Ah, SET OF 2 PERMOBIL	108999-99-0	<input type="radio"/>				
BATTERIES, 85Ah, SET OF 2 PERMOBIL UPCHARGE	109307-99-0	<input type="radio"/>				
WITHOUT BATTERIES	100661-99-0	<input type="radio"/>				
Charger						
CHARGER VOLTPRO 24/10 PER MOBIL GEL-AGM IP44 EU	109166-99-0	<input checked="" type="radio"/>				
WITHOUT CHARGER	100662-99-0	<input type="radio"/>				
Light kit						
LIGHT KIT	108418-99-0	<input checked="" type="radio"/>				
REFLECTOR KIT	108415-99-0	<input type="radio"/>				
Anti-tippers						
ANTI TIPPER WHEELS	106997-99-0	<input type="checkbox"/>				
Accessories chassis						
LUGGAGE CARRIER WITH REFLECTORS	109170-99-0	<input type="radio"/>				
LUGGAGE CARRIER WITH TAIL LIGHTS	109169-99-0	<input type="radio"/>				
WITHOUT LUGGAGE CARRIER	109284-99-0	<input checked="" type="radio"/>				
PREPARED FOR DAHL CARLOCK	108254-99-0	<input type="radio"/>				101888-99-0 101897-99-0
MOUNTDAHL LOCKING PLATE	109019-99-0	<input type="radio"/>				
REMOTE STOP	108526-99-0	<input type="checkbox"/>				
TARPAULIN/RAINCOVER FOR ALL MODLES (NYLON)	101888-99-0	<input type="radio"/>				
TARPAULIN/RAINCOVER YELLOW (WATERPROOF)	101897-99-0	<input type="radio"/>				109019-99-0 108254-99-0
Seat assembly						
SEAT HEIGHT FIXED	107097-99-0	<input checked="" type="radio"/>				
POWERED SEAT LIFT	107098-99-0	<input type="radio"/>				
FIXED SEAT ANGEL	107096-99-0	<input checked="" type="radio"/>				
POWERED ADJ.SEAT TILT 50°	107095-99-0	<input type="radio"/>				$\alpha > -10^\circ$: 107090-99-0 $\alpha > -20^\circ$: 107091-99-0 $\alpha > -30^\circ$: 107627-99-0
WITHOUT ANTERIOR/NEGATIVE SEAT TILT	107093-99-0	<input checked="" type="radio"/>				107098-99-0 107095-99-0 107751-99-0 101959-99-0
5° ACTIVE SEATING	107099-99-0	<input type="radio"/>				
10° ACTIVE SEATING	107090-99-0	<input type="radio"/>				107099-99-0
20° ACTIVE SEATING (REQ. LEGREST(X)+KNEESUPPORT)	107091-99-0	<input type="radio"/>				
30° ACTIVE SEATING (REQ. LEGREST(X)+KNEESUPPORT)	107627-99-0	<input type="radio"/>				
Seat						
SEAT CORPUS	107022-99-0	<input checked="" type="checkbox"/>				
Seat sizes						

		Depth													
Size (mm)		370	395	420	445	470	495	520	545	570					
Width	420	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>					
	470	<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>					
	520	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>					
	570	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>					
Seat cushion															
SEAT CUSHION VINYL		10XXXX-99-0									<input checked="" type="radio"/>				
SEAT CUSHION STRETCH AIR		10XXXX-99-0									<input type="radio"/>				
WITHOUT SEAT CUSHION		105118-99-0									<input type="radio"/>				
ERGO SUPPORTS SEAT CUSHION		105219-99-0									<input checked="" type="radio"/>				
Adjustment backrest															
POWERED ADJ. BACKREST 150° (LA28)		107825-99-0									<input type="radio"/>				
POWERED ADJ. BACKREST 180° (LA28)		107755-99-0									<input type="radio"/>				
MANUAL ADJUST. BACKREST		108195-99-0									<input checked="" type="radio"/>				
BACKREST FIXED		105066-99-0									<input checked="" type="radio"/>				
BACKREST BIO-MECH		105067-99-0									<input type="radio"/>				
Backrest type															
SHORT		107262-99-0									<input type="radio"/>				
STANDARD		107264-99-0									<input checked="" type="radio"/>				
DIRECT BACKREST ADAPTER		107261-99-0									<input type="radio"/>				
TUBELAR BACKREST ADAPTER		107265-99-0									<input type="radio"/>				
WITHOUT BACKREST SHELL		107272-99-0									<input type="radio"/>				
Backrest sizes															
		Short	Standard												
Size (mm)		470	545	570	595	620	645	670							
Width	360	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>							
	410	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>							
	460	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>							
	510	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>							
Backrest cushion															
BACKREST CUSHION VINYL		10XXXX-99-0									<input checked="" type="radio"/>				
BACKREST CUSHION STRETCH AIR		10XXXX-99-0									<input type="radio"/>				
WITHOUT BACKREST CUSHION		105117-99-0									<input type="radio"/>				
ERGO SUPPORTS BACK SMALL (5 CM DEEP)		105569-99-0									<input checked="" type="radio"/>				
ERGO SUPPORTS BACK LARGE (7 CM DEEP)		105570-99-0									<input type="radio"/>				
Agility Max Contour															
		(H)	Height			(MH)	Mounting height								
		(W)	Size ("/mm)	16"	18"	20"		16"	18"	20"					
Width	14"/355mm	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>		470mm	<input type="radio"/>							
	15"/380mm	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>		495mm	<input type="radio"/>							
	16"/405mm	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>		520mm	<input type="radio"/>	<input type="radio"/>						
	17"/430mm	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>		545mm	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>					
	18"/455mm	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>		570mm	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>					
	19"/485mm	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>		595mm	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>					
	20"/510mm	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>		620mm		<input type="radio"/>	<input type="radio"/>					
							645mm			<input type="radio"/>					
Height Agility lateral supports															
AGILITY LATERAL SUPPORTS LOW POSITION (STANDARD)		108538-99-0									<input checked="" type="radio"/>				
AGILITY LATERAL SUPPORTS HIGH POSITION		108537-99-0									<input type="radio"/>				



Mounting armrests										
WIDTH BETWEEN ARMRESTS: NARROW (38-48CM)	107251-99-0	<input checked="" type="checkbox"/>								
WIDTH BETWEEN ARMRESTS: WIDE (48-58CM)	107252-99-0	<input type="checkbox"/>								
WIDTH BETWEEN ARMRESTS SET TO 380MM	105083-99-0	<input type="checkbox"/>								
WIDTH BETWEEN ARMRESTS SET TO 430MM	105084-99-0	<input checked="" type="checkbox"/>								
WIDTH BETWEEN ARMRESTS SET TO 480MM	105085-99-0	<input type="checkbox"/>								
WIDTH BETWEEN ARMRESTS SET TO 530MM	105086-99-0	<input type="checkbox"/>								
WIDTH BETWEEN ARMRESTS SET TO 580MM	105082-99-0	<input type="checkbox"/>								
ARMREST HEIGHT POS 1: 230MM (FROM SEAT PAN)	106631-99-0	<input type="checkbox"/>								
ARMREST HEIGHT POS 2: 255MM (FROM SEAT PAN)	106632-99-0	<input type="checkbox"/>								
ARMREST HEIGHT POS 3: 280MM (FROM SEAT PAN)	106633-99-0	<input checked="" type="checkbox"/>								
ARMREST HEIGHT POS 4: 305MM (FROM SEAT PAN)	106634-99-0	<input type="checkbox"/>								
ARMREST HEIGHT POS 5: 330MM (FROM SEAT PAN)	106635-99-0	<input type="checkbox"/>								
ARMREST HOLDERS FIXED	105070-99-0	<input checked="" type="checkbox"/>								
ARMREST HOLDERS SWIVEL	105072-99-0	<input type="checkbox"/>								
Armrests										
		Left		Right						
		260	335	410	460	260	335	410	460	
WITHOUT ARMREST		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
TEXTILE ARMREST		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
VINYL ARMREST		<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
GEL ARMREST		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
MOLDED ARMREST LEFT	108860-99-0					<input type="checkbox"/>				
MOLDED ARMREST RIGHT	108861-99-0						<input type="checkbox"/>			
MOLDED HAND PAD MEDIUM LEFT	108862-99-0					<input type="checkbox"/>				
MOLDED HAND PAD MEDIUM RIGHT	108863-99-0						<input type="checkbox"/>			
MOLDED HAND PAD LARGE LEFT	108864-99-0					<input type="checkbox"/>				
MOLDED HAND PAD LARGE RIGHT	108865-99-0						<input type="checkbox"/>			
Adjustment legrest										
POWERED ADJ. LEGREST CORPUS (LA28)	107757-99-0	<input type="checkbox"/>								
POWERED ADJ. LEGREST CORPUS POW EXT LEG (LA28)	107758-99-0	<input type="checkbox"/>								
MANUAL ADJUSTMENT LEGREST (BRACKET LA28)	107763-99-0	<input checked="" type="checkbox"/>								
Legrest										
LEGREST W. 205MM LENGTH COMPENSATION	107012-99-0	<input checked="" type="checkbox"/>								
LEGREST W. 185MM LENGTH COMPENSATION	107184-99-0	<input type="checkbox"/>								
LEGREST WITH POW. LENGTH COMPENSATION 20CM (X)	107185-99-0	<input type="checkbox"/>								
LEGREST W. 32-41(+10)CM HEIGHT ADJ. FOOTPL.	107219-99-0	<input type="checkbox"/>								
LEGREST W. 35-44(+10)CM HEIGHT ADJ. FOOTPL.	107220-99-0	<input type="checkbox"/>								
WITHOUT LEGREST	105760-99-0	<input type="checkbox"/>								
Footplates										
DUAL FOOTPLATE 16CM LEG LENGTH 29-50CM AL.	105436-99-0	<input checked="" type="checkbox"/>								
DUAL FOOTPLATE 20CM LEG LENGTH 29-50CM	105437-99-0	<input type="checkbox"/>								
FOOTPLATE 1 PIECE LEG LENGTH 29-50CM	106760-99-0	<input type="checkbox"/>								
DUAL FOOTPLATE 16CM LEG LENGTH 14-30CM AL.	105438-99-0	<input type="checkbox"/>								
DUAL FOOTPLATE 20CM LEG LENGTH 14-30CM	105439-99-0	<input type="checkbox"/>								
DUAL FOOTPLATE 16CM LEG LENGTH 25-50CM (PWR)	106967-99-0	<input type="checkbox"/>								
FOOTPLATE 1P. 42CM LEG LENGTH 25-50CM (PWR FOOT)	106969-99-0	<input type="checkbox"/>								
DUAL FOOTPLATE 20CM LEG LENGTH 25-50CM (PWR FOOT)	106971-99-0	<input type="checkbox"/>								
DUAL FOOTPLATE 16CM LEG LENGTH <25CM (PWR)	106966-99-0	<input type="checkbox"/>								
DUAL FOOTPLATE 20CM LEG LENGTH <25CM (PWR FOOT)	106970-99-0	<input type="checkbox"/>								



FOOTPLATE 1P. 42CM LEG LENGTH <25CM (PWR FOOT)	106968-99-0	<input type="radio"/>				106969-99-0 1x420mm	106968-99-0 1x420mm
DUAL FOOTPLATE 16CM	107453-99-0	<input type="radio"/>				 104646-99-0 2x160mm 106008-99-0 2x200mm	
DUAL FOOTPLATE 20CM	106008-99-0	<input type="radio"/>					
FOOTPLATE CHOICE NOT REQUIRED	106154-99-0	<input type="radio"/>					

Accessories seat

Headrest

HEAD SUPPORT, ADJUSTABLE TEXTILE(BLACK)	106336-99-0	<input type="radio"/>				 106331-99-0 106336-99-0 106332-99-0 106333-99-0 106637-99-0 106638-99-0 106634-99-0 106635-99-0
HEAD SUPPORT, ADJUSTABLE TEXTILE(BLACK)/SHORT BACK	106337-99-0	<input type="radio"/>				
HEAD SUPPORT, ADJUSTABLE VINYL	106331-99-0	<input type="radio"/>				
HEAD SUPPORT, ADJUSTABLE VINYL/SHORT BACK	106338-99-0	<input type="radio"/>				
HEAD SUPPORT "ERGO" TEXTILE(BLACK)	106333-99-0	<input type="radio"/>				
HEAD SUPPORT "ERGO" TEXT. (BLACK) SHORT BACK	106335-99-0	<input type="radio"/>				
HEAD SUPPORT "ERGO"	106332-99-0	<input type="radio"/>				
HEAD SUPPORT "ERGO" SHORT BACK	106334-99-0	<input type="radio"/>				
HEADREST ORIGINAL 340x130 LONG HOLDER, STD BACK	107716-99-0	<input type="radio"/>				 340x130mm 250x125mm 180x110mm 107712-99-0 107711-99-0 107714-99-0 107713-99-0 107708-99-0 107710-99-0 107719-99-0 107719-99-0 107709-99-0 107716-99-0 107715-99-0 107718-99-0 107715-99-0 108580-99-0 107717-99-0 108582-99-0 108582-99-0 108579-99-0 108584-99-0 108581-99-0 108583-99-0 STD LOW BACK / AGILITY 107714-99-0 107713-99-0 108579-99-0 107712-99-0 107711-99-0 108580-99-0 107710-99-0 107709-99-0 108581-99-0 107708-99-0 107719-99-0 108582-99-0 107718-99-0 107717-99-0 108583-99-0 107716-99-0 107715-99-0 108584-99-0
HEADREST ORIGINAL 340x130 LONG HOLDER, SHORT BACK	107715-99-0	<input type="radio"/>				
HEADREST ORIGINAL 340x130 STD HOLDER, STD BACK	107718-99-0	<input type="radio"/>				
HEADREST ORIGINAL 340x130 STD HOLDER, SHORT BACK	107717-99-0	<input type="radio"/>				
HEADREST ORIGINAL 250x125 LONG HOLDER, STD BACK	107708-99-0	<input type="radio"/>				
HEADREST ORIGINAL 250x125 LONG HOLDER, SHORT BACK	107719-99-0	<input type="radio"/>				
HEADREST ORIGINAL 250x125 STD HOLDER, STD BACK	107710-99-0	<input type="radio"/>				
HEADREST ORIGINAL 250x125 STD HOLDER, SHORT BACK	107709-99-0	<input type="radio"/>				
HEADREST ORIGINAL 180x110 LONG HOLDER, STD BACK	107712-99-0	<input type="radio"/>				
HEADREST ORIGINAL 180x110 LONG HOLDER, SHORT BACK	107711-99-0	<input type="radio"/>				
HEADREST ORIGINAL 180x110 STD HOLDER, STD BACK	107714-99-0	<input type="radio"/>				
HEADREST ORIGINAL 180x110 STD HOLDER, SHORT BACK	107713-99-0	<input type="radio"/>				
HEADREST ORIGINAL 180x110 STD HOLDER, AGILITY	108579-99-0	<input type="radio"/>				
HEADREST ORIGINAL 180x110 LONG HOLDER, AGILITY	108580-99-0	<input type="radio"/>				
HEADREST ORIGINAL 250x125 STD HOLDER, AGILITY	108581-99-0	<input type="radio"/>				
HEADREST ORIGINAL 250x125 LONG HOLDER, AGILITY	108582-99-0	<input type="radio"/>				
HEADREST ORIGINAL 340x130 STD HOLDER, AGILITY	108583-99-0	<input type="radio"/>				
HEADREST ORIGINAL 340x130 LONG HOLDER, AGILITY	108584-99-0	<input type="radio"/>				
BODILINK HEADREST L MAX 8 STD BACK	108719-99-0	<input type="radio"/>				 340x130mm Large 250x125mm Medium 180x110mm Small
BODILINK HEADREST L MAX 8 LOW BACK	108720-99-0	<input type="radio"/>				
BODILINK HEADREST M MAX 8 STD BACK	108717-99-0	<input type="radio"/>				
BODILINK HEADREST M MAX 8 LOW BACK	108718-99-0	<input type="radio"/>				
BODILINK HEADREST S MAX 8 STD BACK	109003-99-0	<input type="radio"/>				
BODILINK HEADREST S MAX 8 LOW BACK	109004-99-0	<input type="radio"/>				

Facial support pad / switch Left

SWING AWAY CHIN BOOM ASSY SUPPORT STD LEFT/HIGH	107860-99-0	<input type="radio"/>				107860-99-0	107941-99-0
SWING AWAY CHIN BOOM ASSY SUPPORT LONG LEFT/HIGH	107861-99-0	<input type="radio"/>					
SWING AWAY CHIN BOOM ASSY F. HEADREST LEFT/HIGH	107858-99-0	<input type="radio"/>				107861-99-0	107944-99-0
SWING AWAY ASSY W. TASH MICRO LIGHT LEFT/HIGH	107859-99-0	<input type="radio"/>					

Facial support pad / switch Right

SWING AWAY ASSY W. FACIAL SUPPORT STD RIGHT/HIGH	107941-99-0	<input type="radio"/>				107858-99-0	107931-99-0
SWING AWAY ASSY W. FACIAL SUPPORT LONG RIGHT/HIGH	107944-99-0	<input type="radio"/>					
SWING AWAY EGG SW. ASSY F. HEADREST RIGHT/HIGH	107931-99-0	<input type="radio"/>				107859-99-0	107935-99-0
SWING AWAY ASSY W. TASH MICRO LIGHT RIGHT/HIGH	107935-99-0	<input type="radio"/>					

Extension cable

EXTENSION CABLE 3MM STEREO, L=1M	105931-99-0	<input type="radio"/>					
----------------------------------	-------------	-----------------------	--	--	--	--	--

EXTENSION CABLE 3MM STEREO, L=2M	105932-99-0	<input type="radio"/>					
Lateral support							
BODILINK LATERAL SUPP. L. STRETCH-AIR 75X140	109025-99-0	<input type="radio"/>					
BODILINK LATERAL SUPP. L. STRETCH-AIR 100X165	109027-99-0	<input type="radio"/>					
BODILINK LATERAL SUPP. R. STRETCH-AIR 75X140	109026-99-0	<input type="radio"/>					
BODILINK LATERAL SUPP. R. STRETCH-AIR 100X165	109028-99-0	<input type="radio"/>					
LATERAL SUPPORT ORIGINAL LEFT 180X130	107720-99-0	<input type="radio"/>					
LATERAL SUPPORT ORIGINAL LEFT 130X80	107722-99-0	<input type="radio"/>					
LATERAL SUPPORT ORIGINAL RIGHT 180X130	107721-99-0	<input type="radio"/>					
LATERAL SUPPORT ORIGINAL RIGHT 130X80	107723-99-0	<input type="radio"/>					
Thigh supports							
BODILINK THIGH SUPPORT L. STRETCH-AIR 90X110 S(B)	109179-99-0	<input type="radio"/>					
BODILINK THIGH SUPPORT L. STRETCH-AIR 200X90 L(C)	109181-99-0	<input type="radio"/>					
BODILINK THIGH SUPPORT L. STRETCH-AIR 300X90 XL(G)	109036-99-0	<input type="radio"/>					
BODILINK THIGH SUPPORT R. STRETCH-AIR 90X110 S(B)	109180-99-0	<input type="radio"/>					
BODILINK THIGH SUPPORT R. STRETCH-AIR 200X90 L(C)	109182-99-0	<input type="radio"/>					
BODILINK THIGH SUPPORT R. STRETCH-AIR 300X90 XL(G)	109037-99-0	<input type="radio"/>					
THIGH SUPPORT LEFT TEXT. (MESH) 110x90 SMALL (B)	107392-99-0	<input type="radio"/>					
THIGH SUPPORT LEFT TEXT. (MESH) 165x120 MEDIUM (J)	107396-99-0	<input type="radio"/>					
THIGH SUPPORT LEFT TEXT. (MESH) 200x90 LARGE (C)	107394-99-0	<input type="radio"/>					
THIGH SUPPORT LEFT TEXT. (MESH) 300x90 XLARGE (G)	107390-99-0	<input type="radio"/>					
THIGH SUPPORT RIGHT TEXT. (MESH) 110x90 SMALL (B)	107393-99-0	<input type="radio"/>					
THIGH SUPPORT RIGHT TEXT. (MESH) 165x120 MEDIUM (J)	107397-99-0	<input type="radio"/>					
THIGH SUPPORT RIGHT TEXT. (MESH) 200x90 LARGE (C)	107395-99-0	<input type="radio"/>					
THIGH SUPPORT RIGHT TEXT. (MESH) 300x90 XLARGE (G)	107391-99-0	<input type="radio"/>					
Hip supports							
BODILINK HIP SUPPORT L. STRETCH-AIR 90X110 S(B)	109032-99-0	<input type="radio"/>					
BODILINK HIP SUPPORT L. STRETCH-AIR 200X90 L(C)	109034-99-0	<input type="radio"/>					
BODILINK HIP SUPPORT L. STRETCH-AIR 150X140 M(A)	109038-99-0	<input type="radio"/>					
BODILINK HIP SUPPORT R. STRETCH-AIR 90X110 S(B)	109033-99-0	<input type="radio"/>					
BODILINK HIP SUPPORT R. STRETCH-AIR 200X90 L(C)	109035-99-0	<input type="radio"/>					
BODILINK HIP SUPPORT R. STRETCH-AIR 150X140 M(A)	109039-99-0	<input type="radio"/>					
HIP SUPPORT LEFT TEXT. (MESH) 110x90 SMALL (B)	107400-99-0	<input type="radio"/>					
HIP SUPPORT LEFT TEXT. (MESH) 150x140 MEDIUM (A)	107398-99-0	<input type="radio"/>					
HIP SUPPORT LEFT TEXT. (MESH) 200x90 LARGE (C)	107402-99-0	<input type="radio"/>					
HIP SUPPORT LEFT TEXT. (MESH) 200x200 XLARGE (D)	107404-99-0	<input type="radio"/>					
HIP SUPPORT RIGHT TEXT. (MESH) 110x90 SMALL (B)	107401-99-0	<input type="radio"/>					
HIP SUPPORT RIGHT TEXT. (MESH) 150x140 MEDIUM (A)	107399-99-0	<input type="radio"/>					
HIP SUPPORT RIGHT TEXT. (MESH) 200x90 LARGE (C)	107403-99-0	<input type="radio"/>					
HIP SUPPORT RIGHT TEXT. (MESH) 200x200 XLARGE (D)	107405-99-0	<input checked="" type="radio"/>					
Calf supports							
CALF SUPPORTS S (B) 110x90 (MESH)	105252-99-0	<input type="radio"/>					
CALF SUPPORTS M (J) 165x120 (MESH)	105258-99-0	<input type="radio"/>					
CALF SUPPORTS L (A) 150x140 (MESH)	105250-99-0	<input type="radio"/>					
CALF SUPPORTS XL (H) 170x175 (MESH)	105256-99-0	<input type="radio"/>					
CALF SUPPORTS XXL (D) 200x200 (MESH)	105254-99-0	<input type="radio"/>					
CALF SUPPORTS S (B) 110x90 (MESH) (VS/HD)	107210-99-0	<input type="radio"/>					
CALF SUPPORTS M (J) 165x120 (MESH) (VS/HD)	107211-99-0	<input type="radio"/>					
CALF SUPPORTS L (A) 150x140 (MESH) (VS/HD)	106190-99-0	<input type="radio"/>					
CALF SUPPORTS XL (H) 170x175 (MESH) (VS/HD)	106188-99-0	<input type="radio"/>					

CALF SUPPORTS XXL (D) 200x200 (MESH) (VS/HD)	106186-99-0	<input type="radio"/>				
Elbow and side supports armrests						
ELBOW SUPPORT LEFT	107838-99-0	<input type="checkbox"/>				 107840-99-0 107863-99-0 107864-99-0 107865-99-0  107838-99-0 107839-99-0
ELBOW SUPPORT RIGHT	107839-99-0	<input type="checkbox"/>				
ARMREST SIDE SUPPORT LEFT OUTSIDE	107863-99-0	<input type="checkbox"/>				
ARMREST SIDE SUPPORT RIGHT OUTSIDE	107840-99-0	<input type="checkbox"/>				
ARMREST SIDE SUPPORT LEFT INSIDE	107865-99-0	<input type="checkbox"/>				
ARMREST SIDE SUPPORT RIGHT INSIDE	107864-99-0	<input type="checkbox"/>				
Gel pads footplates						
STEALTH GEL PADS 16CM WIDE FOOTPATES	106530-99-0	<input type="checkbox"/>				
STEALTH GEL PADS 20CM WIDE FOOTPATES	106225-99-0	<input type="checkbox"/>				
Abduction wedge						
ABDUCTION WEDGE	106762-99-0	<input type="checkbox"/>				
Shoulder support						
SHOULDER SUPPORT FOR 36CM WIDE BACKREST	105410-99-0	<input type="radio"/>				 36cm: 105410-99-0 41cm: 105411-99-0 46cm: 105412-99-0 51cm: 105854-99-0
SHOULDER SUPPORT FOR 41CM WIDE BACKREST	105411-99-0	<input type="radio"/>				
SHOULDER SUPPORT FOR 46CM WIDE BACKREST	105412-99-0	<input type="radio"/>				
SHOULDER SUPPORT FOR 51CM WIDE BACKREST	105854-99-0	<input type="radio"/>				
Stump supports						
STUMP SUPPORT SMALL (A) 150x140 TEXTILE (MESH)	106497-99-0	<input type="radio"/>	<input type="radio"/>			 A: 150x140mm 103067-99-0 J: 165x120mm 103068-99-0 D: 200x200mm 103071-99-0 H: 170x175mm 103070-99-0
STUMP SUPPORT MEDIUM (J) 165x120 TEXTILE (MESH)	106500-99-0	<input type="radio"/>	<input type="radio"/>			
STUMP SUPPORT LARGE (H) 170x175 TEXTILE (MESH)	106499-99-0	<input type="radio"/>	<input type="radio"/>			
STUMP SUPPORT XLARGE (D) 200x200 TEXTILE (MESH)	106498-99-0	<input type="radio"/>	<input type="radio"/>			
Knee support						
KNEE SUPPORT WITH TEXTIEL UPHOLSTERY	107285-99-0	<input type="checkbox"/>				 107285-99-0  107028-99-0
HOLDER KNEE SUPPORT UNITRACK	107028-99-0	<input type="checkbox"/>				
Lumbar support						
LUMBAR SUPPORT,AIR	101972-99-0	<input type="checkbox"/>				
Seatbelts						
CHEST STRAP BODYPOINT ELASTIC LARGE	105843-99-0	<input type="radio"/>				 105845-99-0 105843-99-0
CHEST STRAP BODYPOINT ELASTIC SMALL	105845-99-0	<input type="radio"/>				
SEATBELT, BODYPOINT 2- POINT, MEDIUM, NON-PADDED	104518-99-0	<input type="radio"/>				 Hip belts <u>non padded</u>  Hip belts padded <u>Rehab Latch</u>  Hip belts <u>padded</u>  <u>4-Point R pull</u>  100518-99-0  SAFETY BELT,BUCKLE 100519-99-0  SAFETY BELT,AUTOTENSION 100600-99-0  <u>Monoflex™</u> <u>Chest Support</u>  <u>PivotFit™</u> <u>Shoulder Harness</u>  
SEATBELT, BODYPOINT 2- POINT, LARGE, NON-PADDED	104514-99-0	<input type="radio"/>				
SEATBELT, BODYPOINT 2- POINT, MEDIUM, PADDED	104520-99-0	<input type="radio"/>				
SEATBELT, BODYPOINT 2- POINT, LARGE, PADDED	104516-99-0	<input type="radio"/>				
BP HIP BELT, 2P C PULL LATCH 46CM PAD	108267-99-0	<input type="radio"/>				
BP HIP BELT, 2P C PULL LATCH 62CM PAD	108266-99-0	<input type="radio"/>				
BP HIP BELT, 4P R PULL PUSHBUTT. 46CM PAD	108265-99-0	<input type="radio"/>				
SAFETY BELT,BUCKLE	100518-99-0	<input type="radio"/>				
SAFETY BELT,AUTOTENSION	100519-99-0	<input type="radio"/>	<input type="radio"/>			
SAFETY BELT,AUTOTENSION LONG VERSION	100600-99-0	<input type="radio"/>	<input type="radio"/>			
SAFETY BELT,4-POINT	100529-99-0	<input type="radio"/>				
SEATBELT, BODYPOINT MONOFLEX MEDIUM	105840-99-0	<input type="radio"/>				
SEATBELT, BODYPOINT MONOFLEX LARGE	105839-99-0	<input type="radio"/>				
SEATBELT, BODYPOINT MONOFLEX X-LARGE	105850-99-0	<input type="radio"/>				
BP HARNESS PIVOTFIT S SHOULDER SUPPORT FP	108275-99-0	<input type="radio"/>				
BP HARNESS PIVOTFIT M SHOULDER SUPPORT FP	108276-99-0	<input type="radio"/>				
BP HARNESS PIVOTFIT L SHOULDER SUPPORT FP	108277-99-0	<input type="radio"/>				
BP HARNESS PIVOTFIT XL SHOULDER SUPPORT FP	108278-99-0	<input type="radio"/>				
BP HARNESS STAYFLEX XS CHEST SUPPORT FP	108269-99-0	<input type="radio"/>				

BP HARNESS STAYFLEX S CHEST SUPPORT FP	108270-99-0	<input type="radio"/>				 <p>Stayflex™ Chest Support 100529-99-0</p>
BP HARNESS STAYFLEX M CHEST SUPPORT FP	108271-99-0	<input type="radio"/>				
BP HARNESS STAYFLEX L CHEST SUPPORT FP	108272-99-0	<input type="radio"/>				
BP HARNESS STAYFLEX XL CHEST SUPPORT FP	108273-99-0	<input type="radio"/>				
FOOTSTRAPS WITH PADS (SET OF 2)	103104-99-0	<input type="radio"/>				 <p>103104-99-0 103427-99-0 108242-99-0 108243-99-0</p>
FOOTSTRAPS WITHOUT PADS (SET OF 2)	103427-99-0	<input type="radio"/>				
BP ANKEL HUGGERS MEDIUM	108242-99-0	<input type="radio"/>				
BP ANKEL HUGGERS LARGE	108243-99-0	<input type="radio"/>				
Alternative seat cushions / options						
INCONTINENCE PROTECTIVE COVER	0	<input type="radio"/>				
SEAT CUSHION ERGO AIR		<input type="radio"/>				
Tray tables						
TABLE, SWING AWAY 550MM FOR JOYSTICK L/R	109325-99-0	<input type="radio"/>	<input type="radio"/>			 <p>109325-99-0 105262-99-0 105282-99-0 105281-99-0 105283-99-0 106611-99-0 106809-99-0 106612-99-0 106810-99-0</p>
TABLE, SWING AWAY 600MM FOR JOYSTICK L/R	105262-99-0	<input type="radio"/>	<input type="radio"/>			
TABLE, SWING AWAY 650MM FOR JOYSTICK L/R	105263-99-0	<input type="radio"/>	<input type="radio"/>			
TABLE, SWING AWAY 600MM FOR JOYSTICK CENTER	105264-99-0	<input type="radio"/>	<input type="radio"/>			
TABLE, SWING AWAY 650MM FOR JOYSTICK CENTER	105261-99-0	<input type="radio"/>	<input type="radio"/>			
TABLE, REMOVABLE ANGLE ADJ LEFT, WIDTH 330MM	108510-99-0	<input type="radio"/>				 <p>108509-99-0 108504-99-0 108510-99-0 108505-99-0 108511-99-0 108506-99-0 108512-99-0 108507-99-0 108513-99-0 108508-99-0</p> <p>106863-99-0 107295-99-0</p> <p>105686-99-0 106618-99-0 105831-99-0 106619-99-0</p>
TABLE, REMOVABLE ANGLE ADJ RIGHT, WIDTH 330MM	108505-99-0	<input type="radio"/>				
TABLE, REMOVABLE ANGLE ADJ LEFT, WIDTH 381MM	108511-99-0	<input type="radio"/>				
TABLE, REMOVABLE ANGLE ADJ RIGHT, WIDTH 381MM	108506-99-0	<input type="radio"/>				
TABLE, REMOVABLE ANGLE ADJ LEFT, WIDTH 432MM	108512-99-0	<input type="radio"/>				
TABLE, REMOVABLE ANGLE ADJ RIGHT, WIDTH 432MM	108507-99-0	<input type="radio"/>				
TABLE, REMOVABLE ANGLE ADJ LEFT, WIDTH 482MM	108513-99-0	<input type="radio"/>				
TABLE, REMOVABLE ANGLE ADJ RIGHT, WIDTH 482MM	108508-99-0	<input type="radio"/>				
TABLE, ANGLE ADJUST. 1-SIDE FIXATION	106863-99-0	<input type="radio"/>				
BRACKET ADDITIONAL POSITION TABLE	107295-99-0					
COLLEGE TABLE, LEFT	105686-99-0	<input type="radio"/>				
COLLEGE TABLE RIGHT	105831-99-0	<input type="radio"/>				
Pouches/Bags						
ARMREST POUCH, RIGHT BLACK	107339-99-0	<input type="checkbox"/>	<input type="checkbox"/>			 <p>107338-99-0 107339-99-0 107340-99-0 102059-99-0</p>
ARMREST POUCH, LEFT BLACK	107338-99-0	<input type="checkbox"/>				
SEATFRAME POUCH BLACK	107340-99-0	<input type="checkbox"/>				
BACK PACK	102059-99-0	<input type="checkbox"/>				
BACKPACK HOOKS, SET UNITRACK	105563-99-0	<input type="radio"/>				 <p>105563-99-0 Pos. 1 109186-99-0 Pos. 2 109187-99-0 Pos. 3</p>
BACKPACK HOOK BACKREST PROFILE HIGH	109186-99-0	<input type="radio"/>				
BACKPACK HOOK BACKREST PROFILE LOW	109187-99-0	<input type="radio"/>				
POSITION LEFT BACKPACK HOOK ONLY	109189-99-0	<input type="radio"/>				
POSITION RIGHT BACKPACK HOOK ONLY	109188-99-0	<input type="radio"/>				
Other						
SUPPORT HANDLE ARMREST	106058-99-0	<input type="radio"/>	<input type="radio"/>			 <p>106058-99-0 105260-99-0 105292-99-0 105795-99-0</p>
TRANSFER HANDLES, SET	105260-99-0	<input type="checkbox"/>				
PUSH HANDLES	105292-99-0	<input type="radio"/>				
PUSH HANDLES (BACKREST 360)	105795-99-0	<input type="radio"/>				
CRUTCH HOLDER	105223-99-0	<input type="radio"/>	<input type="radio"/>			 <p>105223-99-0 106972-99-0 106806-99-0</p> <p>105224-99-0 106210-99-0</p>
HOLDER OXYGEN BOTTLE C3G	105224-99-0	<input type="checkbox"/>				
CUP HOLDER	106210-99-0	<input type="checkbox"/>				
PHONE HOLDER X-GRIP	106806-99-0	<input type="checkbox"/>				
TABLET HOLDER X-GRIP	106972-99-0	<input type="checkbox"/>				
Controls						

User control						
JOYSTICK R-NET WITH 3.5" COLOR DISPLAY	107159-99-0	<input checked="" type="radio"/>				<p>107159-99-0 107779-99-0 108138-99-0 108139-99-0</p> <p>105246-99-0</p> <p>108323-99-0 108324-99-0 105293-99-0 105294-99-0</p> <p>103814-99-0 109141-99-0</p> <p>105700-99-0 108140-99-0</p>
JOYSTICK R-NET WITH 3.5" COLOR DISPLAY WITH BT	107779-99-0	<input type="radio"/>				
JOYSTICK R-NET "EL" LED	105246-99-0	<input type="radio"/>				
JOYSTICK R-NET HEAVY DUTY WITH COLOUR DISPLAY	105700-99-0	<input type="radio"/>				
MDT/RNET: JS COMPACT ADVANCED	103814-99-0	<input type="radio"/>				
TABLE STEERING, FLIP-OVER 450MM, BACK W. 350	109323-99-0	<input type="radio"/>	<input type="radio"/>			
TABLE STEERING, FLIP-OVER 500MM, BACK W. 350/400	109324-99-0	<input type="radio"/>	<input type="radio"/>			
TABLE STEERING, FLIP-OVER 600MM, BACK W. 400/450	105293-99-0	<input type="radio"/>	<input type="radio"/>			
TABLE STEERING, FLIP-OVER 650MM, BACK W. 450	105294-99-0	<input type="radio"/>	<input type="radio"/>			
NON STANDARD STEERING (SEE SPECIAL CONTROLS)	101352-99-0	<input type="radio"/>				
Joystick holder <input type="checkbox"/>						
PANELHOLDER PARALLEL MOTION UNITRACK	107754-99-0	<input type="radio"/>	<input checked="" type="radio"/>			<p>107754-99-0 109154-99-0</p> <p>103473-99-0</p> <p>103474-99-0</p> <p>108019-99-0 108020-99-0</p> <p>107788-99-0</p> <p>108019-99-0 109243-99-0</p> <p>108020-99-0 109144-99-0</p>
PANELHOLDER SWING AWAY L=335MM	103473-99-0	<input type="radio"/>	<input type="radio"/>			
PANELHOLDER SWING AWAY L=435MM	103474-99-0	<input type="radio"/>	<input type="radio"/>			
BODYPOINT MIDLINE MOUNT ASSEMBLY KIT R-NET	108019-99-0	<input type="radio"/>	<input type="radio"/>			
BODYPOINT MIDLINE MOUNT ASSEMBLY KIT CJS	108020-99-0	<input type="radio"/>	<input type="radio"/>			
PANEL HOLDER ATTENDANT LONG	107788-99-0	<input type="radio"/>	<input type="radio"/>			
WITHOUT PANELHOLDER	102835-99-0	<input type="radio"/>				
Height adjustable holder						
HEIGHT ADJUSTABLE RIGHT JOYSTICK HOLDER	107960-99-0	<input type="radio"/>				<p>107960-99-0 107963-99-0 108005-99-0</p>
HEIGHT ADJUSTABLE LEFT JOYSTICK HOLDER	107963-99-0	<input type="radio"/>				
HEIGHT ADJUSTABLE JOYSTICK HOLDER	108005-99-0	<input type="radio"/>				
Switch box						
SWITCH BOX WITH TOGGLE SWITCHES	102799-99-0	<input type="radio"/>				<p>102799-99-0 109148-99-0</p> <p>102798-99-0 109147-99-0</p>
SWITCH BOX WITH PUSH BUTTONS	102798-99-0	<input checked="" type="radio"/>				
WITHOUT SWITCH BOX	104621-99-0	<input type="radio"/>				
Bracket switch box						
BRACKET SWITCH BOX JOYSTICK	103469-99-0	<input checked="" type="radio"/>				<p>103469-99-0 104651-99-0 106256-99-0</p>
BRACKET SWITCH BOX UNITRACK ARMREST	104651-99-0	<input type="radio"/>	<input type="radio"/>			
BRACKET SWITCH BOX BACKREST	106256-99-0	<input type="radio"/>	<input type="radio"/>			
SWITCH BOX NOT MOUNTED	105409-99-0	<input type="radio"/>				
Layout switch box						
SWITCHBOX(1) LAYOUT "A": SINGLE LAYER (STD)	103744-99-0	<input checked="" type="radio"/>				<p>(A) 103744-99-0 (B) 103786-99-0 (C) 103787-99-0 (D) 103745-99-0 (E) 103788-99-0 (F) 103789-99-0 (L) 105832-99-0</p>
SWITCHBOX(1) LAYOUT "B": SINGLE LAYER, MEM	103786-99-0	<input type="radio"/>				
SWITCHBOX(1) LAYOUT "C": SWAP LAYERS, MEM	103787-99-0	<input type="radio"/>				
SWITCHBOX(1) LAYOUT "D": SINGLE LAYER (ALT SEQ)	103745-99-0	<input type="radio"/>				
SWITCHBOX(1) LAYOUT "E": SINGLE LAYER, NO TILT	103788-99-0	<input type="radio"/>				
SWITCHBOX(1) LAYOUT "F": SINGLE STAND, SWAP, MEM	103789-99-0	<input type="radio"/>				
SWITCHBOX(1) LAYOUT "L": SWAP, MEM, 5 OPTIONS	105832-99-0	<input type="radio"/>				
Indication level switch box						
SWITCHBOX(1) IND. LEVEL 1: NONE	103746-99-0	<input type="radio"/>				
SWITCHBOX(1) IND. LEVEL 2: SEAT FUNCTION	103747-99-0	<input type="radio"/>				
SWITCHBOX(1) IND. LEVEL 3: SPEED LIMITATION	103748-99-0	<input type="radio"/>				
SWITCHBOX(1) IND. LEVEL 4: COMPLETE	103749-99-0	<input checked="" type="radio"/>				
External charging socket (required with table- or special control)						
EXTRA CHARGER CONNECTOR MOUNTED ON THE ARMREST	104163-99-0	<input type="radio"/>	<input type="radio"/>			<p>104163-99-0 104164-99-0 104165-99-0</p>
EXTRA CHARGER CONNECTOR MOUNTED ON ACC. RAIL	104164-99-0	<input type="radio"/>	<input type="radio"/>			
EXTRA CHARGER CONNECTOR MOUNTED UNDER SEAT FRA	104165-99-0	<input type="radio"/>	<input type="radio"/>			
Lock Key						

LOCKING KEY R-NET	104620-99-0	<input type="checkbox"/>				
Alternative Joystick knobs						
JOYSTICK KNOB CHIN CUP	100960-99-0	<input type="radio"/>				
JOYSTICK KNOB T-GRIP	100961-99-0	<input type="radio"/>				
JOYSTICK KNOB BALL BIG	100962-99-0	<input type="radio"/>				
JOYSTICK KNOB SOFT BALL	100963-99-0	<input type="radio"/>				
JOYSTICK KNOB MUSHROOM	100964-99-0	<input type="radio"/>				
JOYSTICK KNOB STICK	100965-99-0	<input type="radio"/>				
JOYSTICK KNOB BALL SMALL	101165-99-0	<input type="radio"/>				
JOYSTICK KNOB STANDARD	101166-99-0	<input type="radio"/>				
JOYSTICK KNOB SOFT BALL BLACK (CONE)	105442-99-0	<input type="radio"/>				
JOYSTICK KNOB BODYPOINT 7.5CM U-SHAPED	105846-99-0	<input type="radio"/>				
JOYSTICK KNOB BODYPOINT 10CM U-SHAPED	105847-99-0	<input type="radio"/>				
JOYSTICK KNOB BODYPOINT RUBBER DOME	105848-99-0	<input type="radio"/>				
Attendant control						
CO-PILOT	105225-99-0	<input type="radio"/>				
JOYSTICK ATTENDANT RNET WITH 3.5" COLOR DISPLAY	107318-99-0	<input type="radio"/>	<input type="radio"/>			
ATTENDANT CONTROL JOYSTICK LIGHT LED	107320-99-0	<input type="radio"/>	<input type="radio"/>			
JOYSTICK ATTENDANT RNET SPEED/MODE	107321-99-0	<input type="radio"/>	<input type="radio"/>			
Emergency switch						
EMERGENCY STOP SWITCH	101792-99-0	<input type="checkbox"/>				
Supplementary accessoires						
(Mounting costs can apply)						
MDT/RNET: JS COMPACT (LOW FORCE)	103892-99-0	<input type="radio"/>				
MDT/RNET: JS COMPACT	103817-99-0	<input type="radio"/>				
MDT/RNET: MINI JOYSTICK	105289-99-0	<input type="radio"/>				
SWITCH KIT FOR COMPACT JS (CHIN)	104641-99-0	<input type="checkbox"/>				
PANEL HOLDER ATTENDANT	105210-99-0	<input type="checkbox"/>				
CHEST SUPPORT FOR TABLE CII/VII VINYL	101303-99-0	<input type="checkbox"/>				
BLUETOOTH IDEVICE	107350-99-0	<input type="checkbox"/>				
ADAPT. CABLE BLUETOOTH MOUSE	107359-99-0	0	▼			
CHARGER USB	107189-99-0	<input type="radio"/>	<input type="radio"/>			
BUDDY BUTTON BLACK 63 MM	101766-99-0	0	▼			
BUDDY BUTTON RED 63 MM	101767-99-0	0	▼			
BUDDY BUTTON YELLOW 63 MM	101768-99-0	0	▼			
MICRO LIGHT SWITCH	101769-99-0	0	▼			
CUP SWITCH 45 MM	101770-99-0	0	▼			
MINI CUP SWITCH 24 MM	101771-99-0	0	▼			
GRASP SWITCH	101772-99-0	0	▼			
CONNECTION BOX EXTERNAL SWITCHES SEATING	104196-99-0	<input type="checkbox"/>				
CONNECTION BOX EXTERNAL SWITCHES R-NET (RBI)	104696-99-0	<input type="checkbox"/>				
Notes						