



ROHO® Canadian Product List



Table of Contents

Back Systems

Back Systems Comparison, including Mounting Hardware.....	3
AGILITY Min Contour Back with Foam and Lumbar Pad.....	4
AGILITY Min Contour Back with Air and Lumbar Pad.....	4
AGILITY Mid Contour Back with Air and Lumbar Pad.....	5
AGILITY Max Contour Back with Air and Lumbar Pad.....	5
Custom AGILITY Backs.....	6
AGILITY Back Replacement Air Cell Pad.....	6
AGILITY Back Replacement Covers - Standard.....	6-7
AGILITY Min Contour Replacement Foam - with Cut Out.....	8
AGILITY Mid Contour Replacement Foam - with Cut Out.....	8
AGILITY Max Contour Replacement Foam - with Cut Out.....	9
AGILITY Min Contour Replacement Foam - no Cut Out.....	9
AGILITY Headrest Accessories.....	9
AGILITY Lateral Support - Swing-Away and Fixed.....	10
AGILITY Accessories, including Mounting Hardware.....	10
AGILITY Active Aluminum Back.....	11
AGILITY Active Aluminum Replacement Covers.....	12
AGILITY Active Aluminum Replacement Foam.....	12
AGILITY Active Aluminum Replacement Components.....	12
AGILITY Active Carbon Back with Air only.....	13
AGILITY Active Carbon Replacement Covers.....	14
AGILITY Active Carbon Replacement Foam.....	14
AGILITY Active Carbon Replacement Components.....	15

Seat Cushions

NEW ROHO® Hybrid Select Cushions with Cover.....	16
Therapeutic Seat Cushions with Cover: "HIGH PROFILE", "LOW PROFILE" and "MID PROFILE".....	17-18
Therapeutic Cushions with Incontinent Cover: "HIGH PROFILE" and "LOW PROFILE".....	18-19
QUADTRO SELECT® Cushions: "HIGH PROFILE", "LOW PROFILE" and "MID PROFILE".....	19-20
QUADTRO SELECT® Cushions with Incontinent Cover: "HIGH PROFILE" and "LOW PROFILE".....	20-21
ENHANCER® Cushions (dual valve) with Cover.....	21
CONTOUR SELECT® Cushions with Cover.....	21
ROHO® "Smart Check" Cushions with Electronic Cushion Monitoring System.....	22
Sensor Ready Single Valve ROHO® HYBRID ELITE® Cushions with Cover.....	23
Dual Valve ROHO® HYBRID ELITE® Cushions with Cover.....	24
AirLITE® Cushions with Cover.....	25
MOSAIC® Cushions with Cover.....	25
LTV Office Cushions with Cover.....	25
Solid Seat Inserts.....	26
MINI-MAX® Cushions with Cover.....	27
Special Application Cushions.....	27
Accessories/Replacement Pumps/Repair Kits.....	28
ROHO® Mattress Sections & Therapeutic Surfaces.....	28
Replacement Cushion Covers.....	29
Incontinent Cushion Covers.....	29
ROHO® HYBRID ELITE® Cushion Covers.....	30
NEW ROHO® Hybrid Select Cushion Covers.....	30
Size Table for ROHO® Cushions.....	31
Dimensions for ROHO® products.....	32
Notes regarding Cushion Sizing.....	32
Terms and Policies.....	32
Customer Design-a-ROHO® Grid.....	33

AGILITY Back Systems

AGILITY is offered in a Standard Series and a new Active Series. All backs include: mounting hardware, foam with air cell cavity, cover and user manual.



Standard Series
AGILITY Standard Aluminum

OR



Active Series
AGILITY Active Aluminum



Active Series
AGILITY Carbon
Note: Carbon weave pattern was updated in 2021

The three backs have different contours: Min = 3" depth, Mid = 4" depth, Max = 7" depth



Min
Code: AG2



Mid
Code: AG3



Max
Code: AG6



Min
Code: AG2



Mid
Code: AG3

For exact dimensions of the AGILITY Contour Backs, please see the sizing matrix on our website.

Warning: The measurements are given for information only. ROHO strongly recommends that clinicians perform a thorough evaluation of the client and equipment to ensure optimal fit of the AGILITY back.



Standard Series Width

	14"			15"			16"			17"			18"			19"			20"		
Length	Min	Mid	Max	Min	Mid	Max	Min	Mid	Max	Min	Mid	Max	Min	Mid	Max	Min	Mid	Max	Min	Mid	Max
10"																					
13"																					
16"																					
18"																					
20"																					

AGILITY backs have the following mounting hardware available:



Quick Release
Default hardware used for AGILITY Backs, indicated by a suffix in the part number as (-1).



Fixed
Please use (-2) in the part number for that option.



Direct Mount*
Please use (-3) in the part number for that option.
**For Permobil power wheelchairs*

AGILITY Min Contour Back with Foam and Lumbar Pad

Recommended for people using wheelchairs who need basic postural support and minimal trunk support.

Item Number	Width	Height
AG2-1410FL-1	14" / 35.5 cm	10" / 25.5 cm
AG2-1413FL-1	14" / 35.5 cm	13" / 33.0 cm
AG2-1416FL-1	14" / 35.5 cm	16" / 40.5 cm
AG2-1418FL-1	14" / 35.5 cm	18" / 45.5 cm
AG2-1510FL-1	15" / 38.1 cm	10" / 25.5 cm
AG2-1513FL-1	15" / 38.1 cm	13" / 33.0 cm
AG2-1516FL-1	15" / 38.1 cm	16" / 40.5 cm
AG2-1518FL-1	15" / 38.1 cm	18" / 45.5 cm
AG2-1610FL-1	16" / 40.5 cm	10" / 25.5 cm
AG2-1613FL-1	16" / 40.5 cm	13" / 33.0 cm
AG2-1616FL-1	16" / 40.5 cm	16" / 40.5 cm
AG2-1618FL-1	16" / 40.5 cm	18" / 45.5 cm
AG2-1710FL-1	17" / 43.2 cm	10" / 25.5 cm
AG2-1713FL-1	17" / 43.2 cm	13" / 33.0 cm
AG2-1716FL-1	17" / 43.2 cm	16" / 40.5 cm
AG2-1718FL-1	17" / 43.2 cm	18" / 45.5 cm



Item Number	Width	Height
AG2-1810FL-1	18" / 45.5 cm	10" / 25.5 cm
AG2-1813FL-1	18" / 45.5 cm	13" / 33.0 cm
AG2-1816FL-1	18" / 45.5 cm	16" / 40.5 cm
AG2-1818FL-1	18" / 45.5 cm	18" / 45.5 cm
AG2-1910FL-1	19" / 48.3 cm	10" / 25.5 cm
AG2-1913FL-1	19" / 48.3 cm	13" / 33.0 cm
AG2-1916FL-1	19" / 48.3 cm	16" / 40.5 cm
AG2-1918FL-1	19" / 48.3 cm	18" / 45.5 cm
AG2-2010FL-1	20" / 51.0 cm	10" / 25.5 cm
AG2-2013FL-1	20" / 51.0 cm	13" / 33.0 cm
AG2-2016FL-1	20" / 51.0 cm	16" / 40.5 cm
AG2-2018FL-1	20" / 51.0 cm	18" / 45.5 cm

AGILITY Min Contour Back with Air and Lumbar Pad

Recommended for people using wheelchairs who need basic postural support and minimal trunk support.

Item Number	Width	Height
AG2-1410AL-1	14" / 35.5 cm	10" / 25.5 cm
AG2-1413AL-1	14" / 35.5 cm	13" / 33.0 cm
AG2-1416AL-1	14" / 35.5 cm	16" / 40.5 cm
AG2-1418AL-1	14" / 35.5 cm	18" / 45.5 cm
AG2-1510AL-1	15" / 38.1 cm	10" / 25.5 cm
AG2-1513AL-1	15" / 38.1 cm	13" / 33.0 cm
AG2-1516AL-1	15" / 38.1 cm	16" / 40.5 cm
AG2-1518AL-1	15" / 38.1 cm	18" / 45.5 cm
AG2-1610AL-1	16" / 40.5 cm	10" / 25.5 cm
AG2-1613AL-1	16" / 40.5 cm	13" / 33.0 cm
AG2-1616AL-1	16" / 40.5 cm	16" / 40.5 cm
AG2-1618AL-1	16" / 40.5 cm	18" / 45.5 cm
AG2-1710AL-1	17" / 43.2 cm	10" / 25.5 cm
AG2-1713AL-1	17" / 43.2 cm	13" / 33.0 cm
AG2-1716AL-1	17" / 43.2 cm	16" / 40.5 cm
AG2-1718AL-1	17" / 43.2 cm	18" / 45.5 cm



Item Number	Width	Height
AG2-1810AL-1	18" / 45.5 cm	10" / 25.5 cm
AG2-1813AL-1	18" / 45.5 cm	13" / 33.0 cm
AG2-1816AL-1	18" / 45.5 cm	16" / 40.5 cm
AG2-1818AL-1	18" / 45.5 cm	18" / 45.5 cm
AG2-1910AL-1	19" / 48.3 cm	10" / 25.5 cm
AG2-1913AL-1	19" / 48.3 cm	13" / 33.0 cm
AG2-1916AL-1	19" / 48.3 cm	16" / 40.5 cm
AG2-1918AL-1	19" / 48.3 cm	18" / 45.5 cm
AG2-2010AL-1	20" / 51.0 cm	10" / 25.5 cm
AG2-2013AL-1	20" / 51.0 cm	13" / 33.0 cm
AG2-2016AL-1	20" / 51.0 cm	16" / 40.5 cm
AG2-2018AL-1	20" / 51.0 cm	18" / 45.5 cm

AGILITY Mid Contour Back with Air and Lumbar Pad

Recommended for people using wheelchairs who need postural support and mid thoracic trunk support.



Item Number	Width	Height
AG3-1410AL-1	14" / 35.5 cm	10" / 25.5 cm
AG3-1413AL-1	14" / 35.5 cm	13" / 33.0 cm
AG3-1416AL-1	14" / 35.5 cm	16" / 40.5 cm
AG3-1418AL-1	14" / 35.5 cm	18" / 45.5 cm
AG3-1420AL-1	14" / 35.5 cm	20" / 51.0 cm
AG3-1510AL-1	15" / 38.1 cm	10" / 25.5 cm
AG3-1513AL-1	15" / 38.1 cm	13" / 33.0 cm
AG3-1516AL-1	15" / 38.1 cm	16" / 40.5 cm
AG3-1518AL-1	15" / 38.1 cm	18" / 45.5 cm
AG3-1520AL-1	15" / 38.1 cm	20" / 51.0 cm
AG3-1610AL-1	16" / 40.5 cm	10" / 25.5 cm
AG3-1613AL-1	16" / 40.5 cm	13" / 33.0 cm
AG3-1616AL-1	16" / 40.5 cm	16" / 40.5 cm
AG3-1618AL-1	16" / 40.5 cm	18" / 45.5 cm
AG3-1620AL-1	16" / 40.5 cm	20" / 51.0 cm
AG3-1710AL-1	17" / 43.2 cm	10" / 25.5 cm
AG3-1713AL-1	17" / 43.2 cm	13" / 33.0 cm
AG3-1716AL-1	17" / 43.2 cm	16" / 40.5 cm
AG3-1718AL-1	17" / 43.2 cm	18" / 45.5 cm
AG3-1720AL-1	17" / 43.2 cm	20" / 51.0 cm

Item Number	Width	Height
AG3-1810AL-1	18" / 45.5 cm	10" / 25.5 cm
AG3-1813AL-1	18" / 45.5 cm	13" / 33.0 cm
AG3-1816AL-1	18" / 45.5 cm	16" / 40.5 cm
AG3-1818AL-1	18" / 45.5 cm	18" / 45.5 cm
AG3-1820AL-1	18" / 45.5 cm	20" / 51.0 cm
AG3-1910AL-1	19" / 48.3 cm	10" / 25.5 cm
AG3-1913AL-1	19" / 48.3 cm	13" / 33.0 cm
AG3-1916AL-1	19" / 48.3 cm	16" / 40.5 cm
AG3-1918AL-1	19" / 48.3 cm	18" / 45.5 cm
AG3-1920AL-1	19" / 48.3 cm	20" / 51.0 cm
AG3-2010AL-1	20" / 51.0 cm	10" / 25.5 cm
AG3-2013AL-1	20" / 51.0 cm	13" / 33.0 cm
AG3-2016AL-1	20" / 51.0 cm	16" / 40.5 cm
AG3-2018AL-1	20" / 51.0 cm	18" / 45.5 cm
AG3-2020AL-1	20" / 51.0 cm	20" / 51.0 cm

AGILITY Max Contour Back with Air and Lumbar Pad

Recommended for people using wheelchairs who need a deep back that conforms to their individual needs.



Item Number	Width	Height
AG6-1416AL-1	14" / 35.5 cm	16" / 40.5 cm
AG6-1418AL-1	14" / 35.5 cm	18" / 45.5 cm
AG6-1420AL-1	14" / 35.5 cm	20" / 51.0 cm
AG6-1516AL-1	15" / 38.1 cm	16" / 40.5 cm
AG6-1518AL-1	15" / 38.1 cm	18" / 45.5 cm
AG6-1520AL-1	15" / 38.1 cm	20" / 51.0 cm
AG6-1616AL-1	16" / 40.5 cm	16" / 40.5 cm
AG6-1618AL-1	16" / 40.5 cm	18" / 45.5 cm
AG6-1620AL-1	16" / 40.5 cm	20" / 51.0 cm
AG6-1716AL-1	17" / 43.2 cm	16" / 40.5 cm
AG6-1718AL-1	17" / 43.2 cm	18" / 45.5 cm
AG6-1720AL-1	17" / 43.2 cm	20" / 51.0 cm

Item Number	Width	Height
AG6-1816AL-1	18" / 45.5 cm	16" / 40.5 cm
AG6-1818AL-1	18" / 45.5 cm	18" / 45.5 cm
AG6-1820AL-1	18" / 45.5 cm	20" / 51.0 cm
AG6-1916AL-1	19" / 48.3 cm	16" / 40.5 cm
AG6-1918AL-1	19" / 48.3 cm	18" / 45.5 cm
AG6-1920AL-1	19" / 48.3 cm	20" / 51.0 cm
AG6-2016AL-1	20" / 51.0 cm	16" / 40.5 cm
AG6-2018AL-1	20" / 51.0 cm	18" / 45.5 cm
AG6-2020AL-1	20" / 51.0 cm	20" / 51.0 cm

Custom AGILITY Backs

Item Number

AG2CUSTC

AG3CUSTC

AG6CUSTC

Custom AGILITY Back order forms available upon request. Please contact Permobil for more information. Measurements of the width (L) are taken from the outer part of the backrest.

**There will be a loss of 0.5" (1.5 cm) seat depth if the AGILITY is recessed to the smallest part of the width of the wheelchair.

Width	Adjust to fit wheelchair widths	
	Smallest**	Largest
14" / 35.5 cm	13" / 33 cm	14" / 35.5 cm
15" / 38 cm	14" / 35.5 cm	15" / 38 cm
16" / 40.5 cm	15" / 38 cm	16" / 40.5 cm
17" / 43 cm	16" / 40.5 cm	17" / 43 cm
18" / 45.5 cm	17" / 43 cm	18" / 45.5 cm
19" / 48.5 cm	18" / 45.5 cm	19" / 48.5 cm
20" / 51.0 cm	19" / 48.5 cm	20" / 51.0 cm

AGILITY Back Replacement Air Cell Pad

Item Number	Description	Height
AGAIR-10	Air cell pad, AGILITY back	10" / 25.5 cm
AGAIR-13	Air cell pad, AGILITY back	13" / 33.0 cm
AGAIR-16	Air cell pad, AGILITY back	16" / 40.5 cm
AGAIR-18	Air cell pad, AGILITY back	18" / 45.5 cm
AGAIR-20	Air cell pad, AGILITY back	20" / 51.0 cm
AG6-AIR	Air cell pad, AGILITY Max (AG6) side bolster	



AGILITY Back Replacement Covers - Standard

Durable, two-way stretch breathable, vapor-permeable polyester.

Item Number	Width	Height
COV-AG2-1410	14" / 35.5 cm	10" / 25.5 cm
COV-AG2-1413	14" / 35.5 cm	13" / 33.0 cm
COV-AG2-1416	14" / 35.5 cm	16" / 40.5 cm
COV-AG2-1418	14" / 35.5 cm	18" / 45.5 cm
COV-AG2-1510	15" / 38.1 cm	10" / 25.5 cm
COV-AG2-1513	15" / 38.1 cm	13" / 33.0 cm
COV-AG2-1516	15" / 38.1 cm	16" / 40.5 cm
COV-AG2-1518	15" / 38.1 cm	18" / 45.5 cm
COV-AG2-1610	16" / 40.5 cm	10" / 25.5 cm
COV-AG2-1613	16" / 40.5 cm	13" / 33.0 cm
COV-AG2-1616	16" / 40.5 cm	16" / 40.5 cm
COV-AG2-1618	16" / 40.5 cm	18" / 45.5 cm
COV-AG2-1710	17" / 43.2 cm	10" / 25.5 cm
COV-AG2-1713	17" / 43.2 cm	13" / 33.0 cm

Item Number	Width	Height
COV-AG2-1716	17" / 43.2 cm	16" / 40.5 cm
COV-AG2-1718	17" / 43.2 cm	18" / 45.5 cm
COV-AG2-1810	18" / 45.5 cm	10" / 25.5 cm
COV-AG2-1813	18" / 45.5 cm	13" / 33.0 cm
COV-AG2-1816	18" / 45.5 cm	16" / 40.5 cm
COV-AG2-1818	18" / 45.5 cm	18" / 45.5 cm
COV-AG2-1910	19" / 48.3 cm	10" / 25.5 cm
COV-AG2-1913	19" / 48.3 cm	13" / 33.0 cm
COV-AG2-1916	19" / 48.3 cm	16" / 40.5 cm
COV-AG2-1918	19" / 48.3 cm	18" / 45.5 cm
COV-AG2-2010	20" / 51.0 cm	10" / 25.5 cm
COV-AG2-2013	20" / 51.0 cm	13" / 33.0 cm
COV-AG2-2016	20" / 51.0 cm	16" / 40.5 cm
COV-AG2-2018	20" / 51.0 cm	18" / 45.5 cm

AGILITY Back Replacement Covers - Standard (cont.)

Durable, two-way stretch breathable, vapor-permeable polyester.

Item Number	Width	Height
COV-AG3-1410	14" / 35.5 cm	10" / 25.5 cm
COV-AG3-1413	14" / 35.5 cm	13" / 33.0 cm
COV-AG3-1416	14" / 35.5 cm	16" / 40.5 cm
COV-AG3-1418	14" / 35.5 cm	18" / 45.5 cm
COV-AG3-1420	14" / 35.5 cm	20" / 51.0 cm
COV-AG3-1510	15" / 38.1 cm	10" / 25.5 cm
COV-AG3-1513	15" / 38.1 cm	13" / 33.0 cm
COV-AG3-1516	15" / 38.1 cm	16" / 40.5 cm
COV-AG3-1518	15" / 38.1 cm	18" / 45.5 cm
COV-AG3-1520	15" / 38.1 cm	20" / 51.0 cm
COV-AG3-1610	16" / 40.5 cm	10" / 25.5 cm
COV-AG3-1613	16" / 40.5 cm	13" / 33.0 cm
COV-AG3-1616	16" / 40.5 cm	16" / 40.5 cm
COV-AG3-1618	16" / 40.5 cm	18" / 45.5 cm
COV-AG3-1620	16" / 40.5 cm	20" / 51.0 cm
COV-AG3-1710	17" / 43.2 cm	10" / 25.5 cm
COV-AG3-1713	17" / 43.2 cm	13" / 33.0 cm
COV-AG3-1716	17" / 43.2 cm	16" / 40.5 cm
COV-AG3-1718	17" / 43.2 cm	18" / 45.5 cm
COV-AG3-1720	17" / 43.2 cm	20" / 51.0 cm
COV-AG3-1810	18" / 45.5 cm	10" / 25.5 cm
COV-AG3-1813	18" / 45.5 cm	13" / 33.0 cm
COV-AG3-1816	18" / 45.5 cm	16" / 40.5 cm
COV-AG3-1818	18" / 45.5 cm	18" / 45.5 cm
COV-AG3-1820	18" / 45.5 cm	20" / 51.0 cm
COV-AG3-1910	19" / 48.3 cm	10" / 25.5 cm
COV-AG3-1913	19" / 48.3 cm	13" / 33.0 cm
COV-AG3-1916	19" / 48.3 cm	16" / 40.5 cm

Item Number	Width	Height
COV-AG3-1918	19" / 48.3 cm	18" / 45.5 cm
COV-AG3-1920	19" / 48.3 cm	20" / 51.0 cm
COV-AG3-2010	20" / 51.0 cm	10" / 25.5 cm
COV-AG3-2013	20" / 51.0 cm	13" / 33.0 cm
COV-AG3-2016	20" / 51.0 cm	16" / 40.5 cm
COV-AG3-2018	20" / 51.0 cm	18" / 45.5 cm
COV-AG3-2020	20" / 51.0 cm	20" / 51.0 cm
COV-AG6-1416	14" / 35.5 cm	16" / 40.5 cm
COV-AG6-1418	14" / 35.5 cm	18" / 45.5 cm
COV-AG6-1420	14" / 35.5 cm	20" / 51.0 cm
COV-AG6-1516	15" / 38.1 cm	16" / 40.5 cm
COV-AG6-1518	15" / 38.1 cm	18" / 45.5 cm
COV-AG6-1520	15" / 38.1 cm	20" / 51.0 cm
COV-AG6-1616	16" / 40.5 cm	16" / 40.5 cm
COV-AG6-1618	16" / 40.5 cm	18" / 45.5 cm
COV-AG6-1620	16" / 40.5 cm	20" / 51.0 cm
COV-AG6-1716	17" / 43.2 cm	16" / 40.5 cm
COV-AG6-1718	17" / 43.2 cm	18" / 45.5 cm
COV-AG6-1720	17" / 43.2 cm	20" / 51.0 cm
COV-AG6-1816	18" / 45.5 cm	16" / 40.5 cm
COV-AG6-1818	18" / 45.5 cm	18" / 45.5 cm
COV-AG6-1820	18" / 45.5 cm	20" / 51.0 cm
COV-AG6-1916	19" / 48.3 cm	16" / 40.5 cm
COV-AG6-1918	19" / 48.3 cm	18" / 45.5 cm
COV-AG6-1920	19" / 48.3 cm	20" / 51.0 cm
COV-AG6-2016	20" / 51.0 cm	16" / 40.5 cm
COV-AG6-2018	20" / 51.0 cm	18" / 45.5 cm
COV-AG6-2020	20" / 51.0 cm	20" / 51.0 cm



AGILITY Min Contour Replacement Foam - With Cut Out

Item Number	Width	Height
AG2FC-1410	14" / 35.5 cm	10" / 25.5 cm
AG2FC-1413	14" / 35.5 cm	13" / 33.0 cm
AG2FC-1416	14" / 35.5 cm	16" / 40.5 cm
AG2FC-1418	14" / 35.5 cm	18" / 45.5 cm
AG2FC-1510	15" / 38.1 cm	10" / 25.5 cm
AG2FC-1513	15" / 38.1 cm	13" / 33.0 cm
AG2FC-1516	15" / 38.1 cm	16" / 40.5 cm
AG2FC-1518	15" / 38.1 cm	18" / 45.5 cm
AG2FC-1610	16" / 40.5 cm	10" / 25.5 cm
AG2FC-1613	16" / 40.5 cm	13" / 33.0 cm
AG2FC-1616	16" / 40.5 cm	16" / 40.5 cm
AG2FC-1618	16" / 40.5 cm	18" / 45.5 cm
AG2FC-1710	17" / 43.2 cm	10" / 25.5 cm
AG2FC-1713	17" / 43.2 cm	13" / 33.0 cm
AG2FC-1716	17" / 43.2 cm	16" / 40.5 cm
AG2FC-1718	17" / 43.2 cm	18" / 45.5 cm
AG2FC-1810	18" / 45.5 cm	10" / 25.5 cm
AG2FC-1813	18" / 45.5 cm	13" / 33.0 cm



Item Number	Width	Height
AG2FC-1816	18" / 45.5 cm	16" / 40.5 cm
AG2FC-1818	18" / 45.5 cm	18" / 45.5 cm
AG2FC-1910	19" / 48.3 cm	10" / 25.5 cm
AG2FC-1913	19" / 48.3 cm	13" / 33.0 cm
AG2FC-1916	19" / 48.3 cm	16" / 40.5 cm
AG2FC-1918	19" / 48.3 cm	18" / 45.5 cm
AG2FC-2010	20" / 51.0 cm	10" / 25.5 cm
AG2FC-2013	20" / 51.0 cm	13" / 33.0 cm
AG2FC-2016	20" / 51.0 cm	16" / 40.5 cm
AG2FC-2018	20" / 51.0 cm	18" / 45.5 cm

AGILITY Mid Contour Replacement Foam - With Cut Out

Item Number	Width	Height
AG3FC-1410	14" / 35.5 cm	10" / 25.5 cm
AG3FC-1413	14" / 35.5 cm	13" / 33.0 cm
AG3FC-1416	14" / 35.5 cm	16" / 40.5 cm
AG3FC-1418	14" / 35.5 cm	18" / 45.5 cm
AG3FC-1420	14" / 35.5 cm	20" / 51.0 cm
AG3FC-1510	15" / 38.1 cm	10" / 25.5 cm
AG3FC-1513	15" / 38.1 cm	13" / 33.0 cm
AG3FC-1516	15" / 38.1 cm	16" / 40.5 cm
AG3FC-1518	15" / 38.1 cm	18" / 45.5 cm
AG3FC-1520	15" / 38.1 cm	20" / 51.0 cm
AG3FC-1610	16" / 40.5 cm	10" / 25.5 cm
AG3FC-1613	16" / 40.5 cm	13" / 33.0 cm
AG3FC-1616	16" / 40.5 cm	16" / 40.5 cm
AG3FC-1618	16" / 40.5 cm	18" / 45.5 cm
AG3FC-1620	16" / 40.5 cm	20" / 51.0 cm
AG3FC-1710	17" / 43.2 cm	10" / 25.5 cm
AG3FC-1713	17" / 43.2 cm	13" / 33.0 cm
AG3FC-1716	17" / 43.2 cm	16" / 40.5 cm

Item Number	Width	Height
AG3FC-1718	17" / 43.2 cm	18" / 45.5 cm
AG3FC-1720	17" / 43.2 cm	20" / 51.0 cm
AG3FC-1810	18" / 45.5 cm	10" / 25.5 cm
AG3FC-1813	18" / 45.5 cm	13" / 33.0 cm
AG3FC-1816	18" / 45.5 cm	16" / 40.5 cm
AG3FC-1818	18" / 45.5 cm	18" / 45.5 cm
AG3FC-1820	18" / 45.5 cm	20" / 51.0 cm
AG3FC-1910	19" / 48.3 cm	10" / 25.5 cm
AG3FC-1913	19" / 48.3 cm	13" / 33.0 cm
AG3FC-1916	19" / 48.3 cm	16" / 40.5 cm
AG3FC-1918	19" / 48.3 cm	18" / 45.5 cm
AG3FC-1920	19" / 48.3 cm	20" / 51.0 cm
AG3FC-2010	20" / 51.0 cm	10" / 25.5 cm
AG3FC-2013	20" / 51.0 cm	13" / 33.0 cm
AG3FC-2016	20" / 51.0 cm	16" / 40.5 cm
AG3FC-2018	20" / 51.0 cm	18" / 45.5 cm
AG3FC-2020	20" / 51.0 cm	20" / 51.0 cm

AGILITY Max Contour Replacement Foam - With Cut Out

Item Number	Width	Height
AG6FC-1416	14" / 35.5 cm	16" / 40.5 cm
AG6FC-1418	14" / 35.5 cm	18" / 45.5 cm
AG6FC-1420	14" / 35.5 cm	20" / 51.0 cm
AG6FC-1516	15" / 38.1 cm	16" / 40.5 cm
AG6FC-1518	15" / 38.1 cm	18" / 45.5 cm
AG6FC-1520	15" / 38.1 cm	20" / 51.0 cm
AG6FC-1616	16" / 40.5 cm	16" / 40.5 cm
AG6FC-1618	16" / 40.5 cm	18" / 45.5 cm
AG6FC-1620	16" / 40.5 cm	20" / 51.0 cm
AG6FC-1716	17" / 43.2 cm	16" / 40.5 cm
AG6FC-1718	17" / 43.2 cm	18" / 45.5 cm

Item Number	Width	Height
AG6FC-1720	17" / 43.2 cm	20" / 51.0 cm
AG6FC-1816	18" / 45.5 cm	16" / 40.5 cm
AG6FC-1818	18" / 45.5 cm	18" / 45.5 cm
AG6FC-1820	18" / 45.5 cm	20" / 51.0 cm
AG6FC-1916	19" / 48.3 cm	16" / 40.5 cm
AG6FC-1918	19" / 48.3 cm	18" / 45.5 cm
AG6FC-1920	19" / 48.3 cm	20" / 51.0 cm
AG6FC-2016	20" / 51.0 cm	16" / 40.5 cm
AG6FC-2018	20" / 51.0 cm	18" / 45.5 cm
AG6FC-2020	20" / 51.0 cm	20" / 51.0 cm

AGILITY Min Contour Replacement Foam - No Cut Out

Item Number	Width	Height
AG2F-1410	14" / 35.5 cm	10" / 25.5 cm
AG2F-1413	14" / 35.5 cm	13" / 33.0 cm
AG2F-1416	14" / 35.5 cm	16" / 40.5 cm
AG2F-1418	14" / 35.5 cm	18" / 45.5 cm
AG2F-1510	15" / 38.1 cm	10" / 25.5 cm
AG2F-1513	15" / 38.1 cm	13" / 33.0 cm
AG2F-1516	15" / 38.1 cm	16" / 40.5 cm
AG2F-1518	15" / 38.1 cm	18" / 45.5 cm
AG2F-1610	16" / 40.5 cm	10" / 25.5 cm
AG2F-1613	16" / 40.5 cm	13" / 33.0 cm
AG2F-1616	16" / 40.5 cm	16" / 40.5 cm
AG2F-1618	16" / 40.5 cm	18" / 45.5 cm
AG2F-1710	17" / 43.2 cm	10" / 25.5 cm
AG2F-1713	17" / 43.2 cm	13" / 33.0 cm
AG2F-1716	17" / 43.2 cm	16" / 40.5 cm

Item Number	Width	Height
AG2F-1718	17" / 43.2 cm	18" / 45.5 cm
AG2F-1810	18" / 45.5 cm	10" / 25.5 cm
AG2F-1813	18" / 45.5 cm	13" / 33.0 cm
AG2F-1816	18" / 45.5 cm	16" / 40.5 cm
AG2F-1818	18" / 45.5 cm	18" / 45.5 cm
AG2F-1910	19" / 48.3 cm	10" / 25.5 cm
AG2F-1913	19" / 48.3 cm	13" / 33.0 cm
AG2F-1916	19" / 48.3 cm	16" / 40.5 cm
AG2F-1918	19" / 48.3 cm	18" / 45.5 cm
AG2F-2010	20" / 51.0 cm	10" / 25.5 cm
AG2F-2013	20" / 51.0 cm	13" / 33.0 cm
AG2F-2016	20" / 51.0 cm	16" / 40.5 cm
AG2F-2018	20" / 51.0 cm	18" / 45.5 cm

AGILITY Headrest Accessories

Item Number	Description
AG-HR-BRKT	Headrest, AGILITY universal adapter plate
AG-HR-LG	AGILITY Headrest Bracket with large pad and cover
AG-HR-MED	AGILITY Headrest Bracket with medium pad and cover
AG-HR-SM	AGILITY Headrest Bracket with small pad and cover



AGILITY Lateral Support - Swing-Away and Fixed

Please note that Lateral Supports can be installed on all heights and are sold individually.

Swing-Away with curved pad	Pad size	Offset type
AG-LATSA-M-FL	Medium; 4" T x 4" L	Flush
AG-LATSA-M-OF	Medium; 4" T x 4" L	1/2" offset
AG-LATSA-M-OF1	Medium; 4" T x 4" L	1" offset
AG-LATSA-M-OF2	Medium; 4" T x 4" L	2" offset
AG-LATSA-L-FL	Large; 4" T x 5" L	Flush
AG-LATSA-L-OF	Large; 4" T x 5" L	1/2" offset
AG-LATSA-L-OF1	Large; 4" T x 5" L	1" offset
AG-LATSA-L-OF2	Large; 4" T x 5" L	2" offset



Swing-Away with curved pad	Pad size	Offset type
AG-LATSA-X-FL	X-Large; 5" T x 5" L	Flush
AG-LATSA-X-OF	X-Large; 5" T x 5" L	1/2" offset
AG-LATSA-X-OF1	X-Large; 5" T x 5" L	1" offset
AG-LATSA-X-OF2	X-Large; 5" T x 5" L	2" offset
AG-LATSACUST	Price upon request	

Fixed with flat pad	Pad size	Offset type
AG-LATFX-M-FL	Medium; 4" T x 4" L	Flush
AG-LATFX-M-OF	Medium; 4" T x 4" L	1/2" offset
AG-LATFX-L-FL	Large; 4" T x 5" L	Flush
AG-LATFX-L-OF	Large; 4" T x 5" L	1/2" offset
AG-LATFX-X-FL	X-Large; 5" T x 5" L	Flush
AG-LATFX-X-OF	X-Large; 5" T x 5" L	1/2" offset
AG-LATFXCUST	Price upon request	



AGILITY Accessories

Item Number	Description
AGLUM	Lumbar pad, AGILITY Back
AG-WRENCH	Allen wrench, 4 mm AGILITY Back
AG-SPACKIT	Spacing plates
AG-PS14	AGILITY Back privacy shield fitting chair size 14/15
AG-PS16	AGILITY Back privacy shield fitting chair size 16/17
AG-PS18	AGILITY Back privacy shield fitting chair size 18/19
AG-PS20	AGILITY Back privacy shield fitting chair size 20/21
AG-HDWR-FX-SET	Fixed hardware for AGILITY Back (pair)
AG-HDWR-QR-SET	"Quick Release" hardware for AGILITY Back (pair)
AG-HDWR-DM14	Direct mount hardware for AGILITY; 14" wide
AG-HDWR-DM15	Direct mount hardware for AGILITY; 15" wide
AG-HDWR-DM16	Direct mount hardware for AGILITY; 16" wide
AG-HDWR-DM17	Direct mount hardware for AGILITY; 17" wide
AG-HDWR-DM18	Direct mount hardware for AGILITY; 18" wide
AG-HDWR-DM19	Direct mount hardware for AGILITY; 19" wide
AG-HDWR-DM20	Direct mount hardware for AGILITY; 20" wide



AGILITY Active Series is available with an aluminum or carbon shell, and in both Min and Mid contours.

AGILITY Active Aluminum comes standard with **Fixed Adjustable Hardware**.

Please use (-2A) in the part number.



Length	14"				15"				16"				17"				18"			
	Aluminum Min	Aluminum Mid	Carbon Min	Carbon Mid	Aluminum Min	Aluminum Mid	Carbon Min	Carbon Mid	Aluminum Min	Aluminum Mid	Carbon Min	Carbon Mid	Aluminum Min	Aluminum Mid	Carbon Min	Carbon Mid	Aluminum Min	Aluminum Mid	Carbon Min	Carbon Mid
8"																				
10"																				
13"																				
16"																				

AGILITY Active Aluminum Back with Foam

Recommended for people using wheelchairs who need basic postural support and minimal or medium trunk support.

AGILITY Active Aluminum Min Contour		
Item Number	Width	Height
AG2-1408F-2A	14" / 35.5 cm	8" / 20 cm
AG2-1410F-2A	14" / 35.5 cm	10" / 25.5 cm
AG2-1508F-2A	15" / 38.1 cm	8" / 20 cm
AG2-1510F-2A	15" / 38.1 cm	10" / 25.5 cm
AG2-1608F-2A	16" / 40.5 cm	8" / 20 cm
AG2-1610F-2A	16" / 40.5 cm	10" / 25.5 cm
AG2-1708F-2A	17" / 43.2 cm	8" / 20 cm
AG2-1710F-2A	17" / 43.2 cm	10" / 25.5 cm
AG2-1808F-2A	18" / 45.5 cm	8" / 20 cm
AG2-1810F-2A	18" / 45.5 cm	10" / 25.5 cm

AGILITY Active Aluminum Mid Contour		
Item Number	Width	Height
AG3-1408F-2A	14" / 35.5 cm	8" / 20 cm
AG3-1410F-2A	14" / 35.5 cm	10" / 25.5 cm
AG3-1508F-2A	15" / 38.1 cm	8" / 20 cm
AG3-1510F-2A	15" / 38.1 cm	10" / 25.5 cm
AG3-1608F-2A	16" / 40.5 cm	8" / 20 cm
AG3-1610F-2A	16" / 40.5 cm	10" / 25.5 cm
AG3-1708F-2A	17" / 43.2 cm	8" / 20 cm
AG3-1710F-2A	17" / 43.2 cm	10" / 25.5 cm
AG3-1808F-2A	18" / 45.5 cm	8" / 20 cm
AG3-1810F-2A	18" / 45.5 cm	10" / 25.5 cm



AGILITY Active Aluminum Replacement Covers

AGILITY Active Aluminum Min Contour		
Item Number	Width	Height
COV-AG2-1408CB	14" / 35.5 cm	8" / 20 cm
COV-AG2-1410CB	14" / 35.5 cm	10" / 25.5 cm
COV-AG2-1508CB	15" / 38.1 cm	8" / 20 cm
COV-AG2-1510CB	15" / 38.1 cm	10" / 25.5 cm
COV-AG2-1608CB	16" / 40.5 cm	8" / 20 cm
COV-AG2-1610CB	16" / 40.5 cm	10" / 25.5 cm
COV-AG2-1708CB	17" / 43.2 cm	8" / 20 cm
COV-AG2-1710CB	17" / 43.2 cm	10" / 25.5 cm
COV-AG2-1808CB	18" / 45.5 cm	8" / 20 cm
COV-AG2-1810CB	18" / 45.5 cm	10" / 25.5 cm

AGILITY Active Aluminum Mid Contour		
Item Number	Width	Height
COV-AG3-1408CB	14" / 35.5 cm	8" / 20 cm
COV-AG3-1410CB	14" / 35.5 cm	10" / 25.5 cm
COV-AG3-1508CB	15" / 38.1 cm	8" / 20 cm
COV-AG3-1510CB	15" / 38.1 cm	10" / 25.5 cm
COV-AG3-1608CB	16" / 40.5 cm	8" / 20 cm
COV-AG3-1610CB	16" / 40.5 cm	10" / 25.5 cm
COV-AG3-1708CB	17" / 43.2 cm	8" / 20 cm
COV-AG3-1710CB	17" / 43.2 cm	10" / 25.5 cm
COV-AG3-1808CB	18" / 45.5 cm	8" / 20 cm
COV-AG3-1810CB	18" / 45.5 cm	10" / 25.5 cm

AGILITY Active Aluminum Replacement Foam

AGILITY Active Aluminum Min Contour		
Item Number	Width	Height
AG2FO-1408AC	14" / 35.5 cm	8" / 20 cm
AG2FO-1410AC	14" / 35.5 cm	10" / 25.5 cm
AG2FO-1508AC	15" / 38.1 cm	8" / 20 cm
AG2FO-1510AC	15" / 38.1 cm	10" / 25.5 cm
AG2FO-1608AC	16" / 40.5 cm	8" / 20 cm
AG2FO-1610AC	16" / 40.5 cm	10" / 25.5 cm
AG2FO-1708AC	17" / 43.2 cm	8" / 20 cm
AG2FO-1710AC	17" / 43.2 cm	10" / 25.5 cm
AG2FO-1808AC	18" / 45.5 cm	8" / 20 cm
AG2FO-1810AC	18" / 45.5 cm	10" / 25.5 cm

AGILITY Active Aluminum Mid Contour		
Item Number	Width	Height
AG3FO-1408AC	14" / 35.5 cm	8" / 20 cm
AG3FO-1410AC	14" / 35.5 cm	10" / 25.5 cm
AG3FO-1508AC	15" / 38.1 cm	8" / 20 cm
AG3FO-1510AC	15" / 38.1 cm	10" / 25.5 cm
AG3FO-1608AC	16" / 40.5 cm	8" / 20 cm
AG3FO-1610AC	16" / 40.5 cm	10" / 25.5 cm
AG3FO-1708AC	17" / 43.2 cm	8" / 20 cm
AG3FO-1710AC	17" / 43.2 cm	10" / 25.5 cm
AG3FO-1808AC	18" / 45.5 cm	8" / 20 cm
AG3FO-1810AC	18" / 45.5 cm	10" / 25.5 cm

AGILITY Active Aluminum Replacement Components

Item Number	Description
AG-HDWR-2A	Hardware Active Aluminum Back - Fixed
AG-CLAMP-FXA	Clamp Active Aluminum - Fixed
AG-BOLTKIT-A	Bolt Kit Active Aluminum - Fixed



AGILITY Active Aluminum
Fixed Adjustable Hardware

AGILITY Active Carbon Back with Air

Recommended for people using wheelchairs who need basic postural support and minimal or medium trunk support.

AGILITY Active Carbon Min Contour		
Item Number	Width	Height
AG2-1410A-2CB	14" / 35.5 cm	10" / 25.5 cm
AG2-1413A-2CB	14" / 35.5 cm	13" / 33.0 cm
AG2-1416A-2CB	14" / 35.5 cm	16" / 40.5 cm
AG2-1510A-2CB	15" / 38.1 cm	10" / 25.5 cm
AG2-1513A-2CB	15" / 38.1 cm	13" / 33.0 cm
AG2-1516A-2CB	15" / 38.1 cm	16" / 40.5 cm
AG2-1610A-2CB	16" / 40.5 cm	10" / 25.5 cm
AG2-1613A-2CB	16" / 40.5 cm	13" / 33.0 cm
AG2-1616A-2CB	16" / 40.5 cm	16" / 40.5 cm
AG2-1710A-2CB	17" / 43.2 cm	10" / 25.5 cm
AG2-1713A-2CB	17" / 43.2 cm	13" / 33.0 cm
AG2-1716A-2CB	17" / 43.2 cm	16" / 40.5 cm
AG2-1810A-2CB	18" / 45.5 cm	10" / 25.5 cm
AG2-1813A-2CB	18" / 45.5 cm	13" / 33.0 cm
AG2-1816A-2CB	18" / 45.5 cm	16" / 40.5 cm

AGILITY Active Carbon Mid Contour		
Item Number	Width	Height
AG3-1410A-2CB	14" / 35.5 cm	10" / 25.5 cm
AG3-1413A-2CB	14" / 35.5 cm	13" / 33.0 cm
AG3-1416A-2CB	14" / 35.5 cm	16" / 40.5 cm
AG3-1510A-2CB	15" / 38.1 cm	10" / 25.5 cm
AG3-1513A-2CB	15" / 38.1 cm	13" / 33.0 cm
AG3-1516A-2CB	15" / 38.1 cm	16" / 40.5 cm
AG3-1610A-2CB	16" / 40.5 cm	10" / 25.5 cm
AG3-1613A-2CB	16" / 40.5 cm	13" / 33.0 cm
AG3-1616A-2CB	16" / 40.5 cm	16" / 40.5 cm
AG3-1710A-2CB	17" / 43.2 cm	10" / 25.5 cm
AG3-1713A-2CB	17" / 43.2 cm	13" / 33.0 cm
AG3-1716A-2CB	17" / 43.2 cm	16" / 40.5 cm
AG3-1810A-2CB	18" / 45.5 cm	10" / 25.5 cm
AG3-1813A-2CB	18" / 45.5 cm	13" / 33.0 cm
AG3-1816A-2CB	18" / 45.5 cm	16" / 40.5 cm



AGILITY Active Carbon weave pattern was updated to the above style in 2021

AGILITY Active Carbon Replacement Covers

AGILITY Active Carbon Min Contour		
Item Number	Width	Height
COV-AG2-1410CB	14" / 35.5 cm	10" / 25.5 cm
COV-AG2-1413CB	14" / 35.5 cm	13" / 33.0 cm
COV-AG2-1416CB	14" / 35.5 cm	16" / 40.5 cm
COV-AG2-1510CB	15" / 38.1 cm	10" / 25.5 cm
COV-AG2-1513CB	15" / 38.1 cm	13" / 33.0 cm
COV-AG2-1516CB	15" / 38.1 cm	16" / 40.5 cm
COV-AG2-1610CB	16" / 40.5 cm	10" / 25.5 cm
COV-AG2-1613CB	16" / 40.5 cm	13" / 33.0 cm
COV-AG2-1616CB	16" / 40.5 cm	16" / 40.5 cm
COV-AG2-1710CB	17" / 43.2 cm	10" / 25.5 cm
COV-AG2-1713CB	17" / 43.2 cm	13" / 33.0 cm
COV-AG2-1716CB	17" / 43.2 cm	16" / 40.5 cm
COV-AG2-1810CB	18" / 45.5 cm	10" / 25.5 cm
COV-AG2-1813CB	18" / 45.5 cm	13" / 33.0 cm
COV-AG2-1816CB	18" / 45.5 cm	16" / 40.5 cm

AGILITY Active Carbon Mid Contour		
Item Number	Width	Height
COV-AG3-1410CB	14" / 35.5 cm	10" / 25.5 cm
COV-AG3-1413CB	14" / 35.5 cm	13" / 33.0 cm
COV-AG3-1416CB	14" / 35.5 cm	16" / 40.5 cm
COV-AG3-1510CB	15" / 38.1 cm	10" / 25.5 cm
COV-AG3-1513CB	15" / 38.1 cm	13" / 33.0 cm
COV-AG3-1516CB	15" / 38.1 cm	16" / 40.5 cm
COV-AG3-1610CB	16" / 40.5 cm	10" / 25.5 cm
COV-AG3-1613CB	16" / 40.5 cm	13" / 33.0 cm
COV-AG3-1616CB	16" / 40.5 cm	16" / 40.5 cm
COV-AG3-1710CB	17" / 43.2 cm	10" / 25.5 cm
COV-AG3-1713CB	17" / 43.2 cm	13" / 33.0 cm
COV-AG3-1716CB	17" / 43.2 cm	16" / 40.5 cm
COV-AG3-1810CB	18" / 45.5 cm	10" / 25.5 cm
COV-AG3-1813CB	18" / 45.5 cm	13" / 33.0 cm
COV-AG3-1816CB	18" / 45.5 cm	16" / 40.5 cm

AGILITY Active Carbon Replacement Foam

AGILITY Active Carbon Min Contour		
Item Number	Width	Height
AG2FC-1410CB	14" / 35.5 cm	10" / 25.5 cm
AG2FC-1413CB	14" / 35.5 cm	13" / 33.0 cm
AG2FC-1416CB	14" / 35.5 cm	16" / 40.5 cm
AG2FC-1510CB	15" / 38.1 cm	10" / 25.5 cm
AG2FC-1513CB	15" / 38.1 cm	13" / 33.0 cm
AG2FC-1516CB	15" / 38.1 cm	16" / 40.5 cm
AG2FC-1610CB	16" / 40.5 cm	10" / 25.5 cm
AG2FC-1613CB	16" / 40.5 cm	13" / 33.0 cm
AG2FC-1616CB	16" / 40.5 cm	16" / 40.5 cm
AG2FC-1710CB	17" / 43.2 cm	10" / 25.5 cm
AG2FC-1713CB	17" / 43.2 cm	13" / 33.0 cm
AG2FC-1716CB	17" / 43.2 cm	16" / 40.5 cm
AG2FC-1810CB	18" / 45.5 cm	10" / 25.5 cm
AG2FC-1813CB	18" / 45.5 cm	13" / 33.0 cm
AG2FC-1816CB	18" / 45.5 cm	16" / 40.5 cm

AGILITY Active Carbon Mid Contour		
Item Number	Width	Height
AG3FC-1410CB	14" / 35.5 cm	10" / 25.5 cm
AG3FC-1413CB	14" / 35.5 cm	13" / 33.0 cm
AG3FC-1416CB	14" / 35.5 cm	16" / 40.5 cm
AG3FC-1510CB	15" / 38.1 cm	10" / 25.5 cm
AG3FC-1513CB	15" / 38.1 cm	13" / 33.0 cm
AG3FC-1516CB	15" / 38.1 cm	16" / 40.5 cm
AG3FC-1610CB	16" / 40.5 cm	10" / 25.5 cm
AG3FC-1613CB	16" / 40.5 cm	13" / 33.0 cm
AG3FC-1616CB	16" / 40.5 cm	16" / 40.5 cm
AG3FC-1710CB	17" / 43.2 cm	10" / 25.5 cm
AG3FC-1713CB	17" / 43.2 cm	13" / 33.0 cm
AG3FC-1716CB	17" / 43.2 cm	16" / 40.5 cm
AG3FC-1810CB	18" / 45.5 cm	10" / 25.5 cm
AG3FC-1813CB	18" / 45.5 cm	13" / 33.0 cm
AG3FC-1816CB	18" / 45.5 cm	16" / 40.5 cm

AGILITY Active Carbon Replacement Components

Item Number	Description
AGAIR-10	AGILITY back cell pad, 10" Height
AGAIR-13	AGILITY back cell pad, 13" Height
AGAIR-16	AGILITY back cell pad, 16" Height
AG-HDWR-2CB	Hardware AGILITY CARBON Back - 2CB (two required)
AG-CLAMP-FXCB	Clamp AG Carbon Inner/Outer set (one side of carbon back)
AG-DISC-CB	Angle adjustable disc for AGILITY Carbon
BOLTKIT-CB	Bolt Kit AG Carbon Fiber Back



AGILITY Active Carbon Back Hardware



Seat Cushions

Please refer to the "Size Table" on page 31 to help you choose the right cushion size.

The ISOFLO® Memory Control opening has been eliminated from ROHO® Universal Covers that come with HIGH PROFILE®, LOW PROFILE®, MID PROFILE™, MINI-MAX® and ENHANCER® Cushions.

Universal Covers ordered separately come standard with the ISOFLO® Memory Control opening unless otherwise indicated. Please consult page 29 of this price list for all ROHO® cover options.

NEW! ROHO® Hybrid Select Cushions with Cover

Advanced hybrid seating with superior positioning. The ROHO overlay and removable IT insert for offloading utilizes ISOFLO® Memory Control. The contoured foam base has deep leg troughs and a coccyx relief cutout.




Item Number	Width	Depth
HS1414C	14" / 35.5 cm	14" / 35.5 cm
HS1416C	14" / 35.5 cm	16" / 40.5 cm
HS1515C	15" / 38 cm	15" / 38 cm
HS1516C	15" / 38 cm	16" / 40.5 cm
HS1517C	15" / 38 cm	17" / 43 cm
HS1615C	16" / 40.5 cm	15" / 38 cm
HS1616C	16" / 40.5 cm	16" / 40.5 cm
HS1617C	16" / 40.5 cm	17" / 43 cm
HS1618C	16" / 40.5 cm	18" / 45.5 cm
HS1620C	16" / 40.5 cm	20" / 51 cm
HS1716C	17" / 43 cm	16" / 40.5 cm

Item Number	Width	Depth
HS1717C	17" / 43 cm	17" / 43 cm
HS1718C	17" / 43 cm	18" / 45.5 cm
HS1720C	17" / 43 cm	20" / 51 cm
HS1816C	18" / 45.5 cm	16" / 40.5 cm
HS1817C	18" / 45.5 cm	17" / 43 cm
HS1818C	18" / 45.5 cm	18" / 45.5 cm
HS1820C	18" / 45.5 cm	20" / 51 cm
HS2016C	20" / 51 cm	16" / 40.5 cm
HS2017C	20" / 51 cm	17" / 43 cm
HS2018C	20" / 51 cm	18" / 45.5 cm
HS2020C	20" / 51 cm	20" / 51 cm

ROHO has created an all-in-one seating solution that brings foam and air technology together with optional offloading.

Above photo shows cushion with a see-through cover.

Below photos show the Hybrid Select without the cover and with lines indicating the air flow and contours; and another with the IT air insert removed.



HIGH PROFILE® Therapeutic Cushions with Cover

Recommended for individuals with a history of skin breakdown; 4" cell height.

The Dual Valve cushion is divided into two separate air compartments. Can be used for positioning, amputees, scoliosis, and more. (Contact sales representative for more info).



SINGLE VALVE - STANDARD SIZES		DUAL VALVE - STANDARD SIZES	
Adult Sizes		Adult Sizes	
Item Number	Description	Item Number	Description
1R88C	15"w x 15"d x 4"h (38 cm x 38 cm x 10 cm)	2R88C	15"w x 15"d x 4"h (38 cm x 38 cm x 10 cm)
1R89C	15"w x 17"d x 4"h (38 cm x 43 cm x 10 cm)	2R89C	15"w x 17"d x 4"h (38 cm x 43 cm x 10 cm)
1R99C	17"w x 17"d x 4"h (43 cm x 43 cm x 10 cm)	2R99C	17"w x 17"d x 4"h (43 cm x 43 cm x 10 cm)
1R910C	17"w x 19"d x 4"h (43 cm x 47 cm x 10 cm)	2R910C	17"w x 19"d x 4"h (43 cm x 47 cm x 10 cm)
1R109C	19"w x 17"d x 4"h (47 cm x 43 cm x 10 cm)	2R109C	19"w x 17"d x 4"h (47 cm x 43 cm x 10 cm)
1R1010C	19"w x 19"d x 4"h (47 cm x 47 cm x 10 cm)	2R1010C	19"w x 19"d x 4"h (47 cm x 47 cm x 10 cm)
Pediatric Sizes		Pediatric Sizes	
1R66C	12"w x 12"d x 4"h (30 cm x 30 cm x 10 cm)	2R66C	12"w x 12"d x 4"h (30 cm x 30 cm x 10 cm)
1R67C	12"w x 13"d x 4"h (30 cm x 33 cm x 10 cm)	2R67C	12"w x 13"d x 4"h (30 cm x 33 cm x 10 cm)
Item Number	Description		
1RXXXXCUST	Single valve custom cushion "High profile" (specify width and depth) up to 160 cells		
2RXXXXCUST	Dual valve custom cushion "High profile" (specify width and depth) up to 160 cells		
C-ROHO	Custom design cushion (provide diagram) up to 160 cells		

*Largest size available is 20 x 15 cells. Smallest size available is 3 x 4 or 4 x 3 cells.

LOW PROFILE® Therapeutic Cushions with Cover

Recommended for active individuals without skin problems; 2" cell height.

The Dual Valve cushion is divided into two separate air compartments. Can be used for positioning, amputees, scoliosis, and more. (Contact sales representative for more info).



SINGLE VALVE - STANDARD SIZES		DUAL VALVE - STANDARD SIZES	
Adult Sizes		Adult Sizes	
Item Number	Description	Item Number	Description
1R88LPC	15"w x 15"d x 2"h (38 cm x 38 cm x 5 cm)	2R88LPC	15"w x 15"d x 2"h (38 cm x 38 cm x 5 cm)
1R89LPC	15"w x 17"d x 2"h (38 cm x 43 cm x 5 cm)	2R89LPC	15"w x 17"d x 2"h (38 cm x 43 cm x 5 cm)
1R99LPC	17"w x 17"d x 2"h (43 cm x 43 cm x 5 cm)	2R99LPC	17"w x 17"d x 2"h (43 cm x 43 cm x 5 cm)
1R910LPC	17"w x 19"d x 2"h (43 cm x 47 cm x 5 cm)	2R910LPC	17"w x 19"d x 2"h (43 cm x 47 cm x 5 cm)
1R109LPC	19"w x 17"d x 2"h (47 cm x 43 cm x 5 cm)	2R109LPC	19"w x 17"d x 2"h (47 cm x 43 cm x 5 cm)
1R1010LPC	19"w x 19"d x 2"h (47 cm x 47 cm x 5 cm)	2R1010LPC	19"w x 19"d x 2"h (47 cm x 47 cm x 5 cm)
Pediatric Sizes		Pediatric Sizes	
1R66LPC	12"w x 12"d x 2"h (30 cm x 30 cm x 5 cm)	2R66LPC	12"w x 12"d x 2"h (30 cm x 30 cm x 5 cm)
1R67LPC	12"w x 13"d x 2"h (30 cm x 33 cm x 5 cm)	2R67LPC	12"w x 13"d x 2"h (30 cm x 33 cm x 5 cm)
Item Number	Description		
1RXXXXLPCCUST	Single valve custom cushion "Low profile" (specify width and depth) up to 160 cells		
2RXXXXLPCCUST	Dual valve custom cushion "Low profile" (specify width and depth) up to 160 cells		
C-ROHO	Custom design cushion (provide diagram) up to 160 cells		

*Largest size available is 20 x 15 cells. Smallest size available is 3 x 4 or 4 x 3 cells.

MID PROFILE™ Therapeutic Cushions with Cover

Recommended for active individuals at increased risk for skin breakdown. 3" cell height.

The Dual Valve cushion is divided into two separate air compartments. Can be used for positioning, amputees, scoliosis, and more. (Contact sales representative for more info).



SINGLE VALVE - STANDARD SIZES	
Adult Sizes	
Item Number	Description
1R77MPC	14"w x 14"d x 3"h (35 cm x 35 cm x 7 cm)
1R87MPC	15"w x 14"d x 3"h (38 cm x 35 cm x 7 cm)
1R88MPC	15"w x 15"d x 3"h (38 cm x 38 cm x 7 cm)
1R89MPC	15"w x 17"d x 3"h (38 cm x 43 cm x 7 cm)
1R99MPC	17"w x 17"d x 3"h (43 cm x 43 cm x 7 cm)
1R910MPC	17"w x 19"d x 3"h (43 cm x 47 cm x 7 cm)
1R108MPC	19"w x 15"d x 3"h (47 cm x 38 cm x 7 cm)
1R109MPC	19"w x 17"d x 3"h (47 cm x 43 cm x 7 cm)
1R1010MPC	19"w x 19"d x 3"h (47 cm x 47 cm x 7 cm)

DUAL VALVE - STANDARD SIZES	
Adult Sizes	
Item Number	Description
2R88MPC	15"w x 15"d x 3"h (38 cm x 38 cm x 7 cm)
2R89MPC	15"w x 17"d x 3"h (38 cm x 43 cm x 7 cm)
2R99MPC	17"w x 17"d x 3"h (43 cm x 43 cm x 7 cm)
2R910MPC	17"w x 19"d x 3"h (43 cm x 47 cm x 7 cm)

Item Number	Description
1RXXXXMPCCUST	Single valve custom cushion "Mid profile" (specify width and depth) up to 160 cells
2RXXXXMPCCUST	Dual valve custom cushion "Mid profile" (specify width and depth) up to 160 cells
C-ROHO	Custom design cushion (provide diagram) up to 160 cells

*Largest size available is 20 x 15 cells. Smallest size available is 3 x 4 or 4 x 3 cells.

HIGH PROFILE® Therapeutic Cushions with (HD) Incontinent Cover

Recommended for individuals requiring ease of cleaning and heavy usage. Ideal in a Long Term Care setting; 4" cell height. See page 29 for the description of the waterproof cover.



Adult Standard Sizes	
Item Number	Description
1R88HD	15"w x 15"d x 4"h (38 cm x 38 cm x 10 cm)
1R89HD	15"w x 17"d x 4"h (38 cm x 43 cm x 10 cm)
1R99HD	17"w x 17"d x 4"h (43 cm x 43 cm x 10 cm)

Adult Standard Sizes	
Item Number	Description
1R910HD	17"w x 19"d x 4"h (43 cm x 47 cm x 10 cm)
1R109HD	19"w x 17"d x 4"h (47 cm x 43 cm x 10 cm)
1R1010HD	19"w x 19"d x 4"h (47 cm x 47 cm x 10 cm)

Item Number	Description
1RXXXXHDCUST	Other dimension "High profile" HD (specify width and depth)

LOW PROFILE® Therapeutic Cushions with (HD) Incontinent Cover

Recommended for individuals without existing skin integrity issues requiring ease of cleaning and heavy usage. Ideal in a Long Term Care setting; 2" cell height. See page 29 for the description of the waterproof cover.



Adult Standard Sizes	
Item Number	Description
1R88LPHD	15"w x 15"d x 2"h (38 cm x 38 cm x 5 cm)
1R89LPHD	15"w x 17"d x 2"h (38 cm x 43 cm x 5 cm)
1R99LPHD	17"w x 17"d x 2"h (43 cm x 43 cm x 5 cm)

Adult Standard Sizes	
Item Number	Description
1R910LPHD	17"w x 19"d x 2"h (43 cm x 47 cm x 5 cm)
1R109LPHD	19"w x 17"d x 2"h (47 cm x 43 cm x 5 cm)
1R1010LPHD	19"w x 19"d x 2"h (47 cm x 47 cm x 5 cm)

Item Number	Description
1RXXXXLPHDCUST	Other dimension "Low profile" HD (specify width and depth)

QUADTRO SELECT® HIGH PROFILE Cushions with Cover

Featuring ISOFLOR® Memory Control. The unique four separate sections of the cushion allows for progressive positioning for short and long term changes; 4" cell height.



Adult Standard Sizes	
Item Number	Description
QS88C	15"w x 15"d x 4"h (38 cm x 38 cm x 10 cm)
QS89C	15"w x 17"d x 4"h (38 cm x 43 cm x 10 cm)
QS99C	17"w x 17"d x 4"h (43 cm x 43 cm x 10 cm)
QS910C	17"w x 19"d x 4"h (43 cm x 47 cm x 10 cm)
QS109C	19"w x 17"d x 4"h (47 cm x 43 cm x 10 cm)
QS1010C	19"w x 19"d x 4"h (47 cm x 47 cm x 10 cm)

Pediatric Sizes	
QS66C	12"w x 12"d x 4"h (30 cm x 30 cm x 10 cm)
QS67C	12"w x 13"d x 4"h (30 cm x 33 cm x 10 cm)

Item Number	Description
QSXXXXCUST	Other dimension "High profile" (specify width and depth) up to 160 cells
C-ROHO	Custom design cushion (provide diagram) up to 160 cells

*Largest size available is 13 x 18 and 14 x 13 cells. Smallest size available is 6 x 4 cells.



QUADTRO SELECT® LOW PROFILE Cushions with Cover

Featuring ISOFLO® Memory Control. Ideal for active users; 2" cell height and four valves promote lateral, interior and posterior stability.



Adult Standard Sizes	
Item Number	Description
QS88LPC	15"w x 15"d x 2"h (38 cm x 38 cm x 5 cm)
QS89LPC	15"w x 17"d x 2"h (38 cm x 43 cm x 5 cm)
QS99LPC	17"w x 17"d x 2"h (43 cm x 43 cm x 5 cm)
QS910LPC	17"w x 19"d x 2"h (43 cm x 47 cm x 5 cm)
QS109LPC	19"w x 17"d x 2"h (47 cm x 43 cm x 5 cm)
QS1010LPC	19"w x 19"d x 2"h (47 cm x 47 cm x 5 cm)

Pediatric Sizes	
QS66LPC	12"w x 12"d x 2"h (30 cm x 30 cm x 5 cm)
QS67LPC	12"w x 13"d x 2"h (30 cm x 33 cm x 5 cm)

Item Number	Description
QSXXXXLPCUST	Other dimension "Low profile" (specify width and depth) up to 160 cells
C-ROHO	Custom design cushion (provide diagram) up to 160 cells

*Largest size available is 13 x 18 cells. Smallest size available is 6 x 4 cells.

QUADTRO SELECT® MID PROFILE Cushions with Cover

Featuring ISOFLO® Memory Control. Offers shape fitting capabilities while the user is seated, allowing quick and easy on-demand adjustment to maximize function. 3" cell height.



Adult Standard Sizes	
Item Number	Description
QS88MPC	15"w x 15"d x 3"h (38 cm x 38 cm x 7 cm)
QS89MPC	15"w x 17"d x 3"h (38 cm x 43 cm x 7 cm)
QS99MPC	17"w x 17"d x 3"h (43 cm x 43 cm x 7 cm)

Adult Standard Sizes	
Item Number	Description
QS910MPC	17"w x 19"d x 3"h (43 cm x 47 cm x 7 cm)
QS109MPC	19"w x 17"d x 3"h (47 cm x 43 cm x 7 cm)
QS1010MPC	19"w x 19"d x 3"h (47 cm x 47 cm x 7 cm)

Item Number	Description
QSXXXXMPCUST	Other dimension "Mid profile" (specify width and depth)
C-ROHO	Custom design cushion (provide diagram)

*Largest size available is 13 x 18 cells. Smallest size available is 6 x 4 cells.

QUADTRO SELECT® HIGH PROFILE Cushions with (HD) Incontinent Cover

Featuring ISOFLO® Memory Control. The unique four separate sections of the cushion allows for progressive positioning for short and long term changes. 4" cell height. See page 29 for the description of the waterproof cover.



Adult Standard Sizes	
Item Number	Description
QS88HD	15"w x 15"d x 4"h (38 cm x 38 cm x 10 cm)
QS89HD	15"w x 17"d x 4"h (38 cm x 43 cm x 10 cm)
QS99HD	17"w x 17"d x 4"h (43 cm x 43 cm x 10 cm)

Adult Standard Sizes	
Item Number	Description
QS910HD	17"w x 19"d x 4"h (43 cm x 47 cm x 10 cm)
QS109HD	19"w x 17"d x 4"h (47 cm x 43 cm x 10 cm)
QS1010HD	19"w x 19"d x 4"h (47 cm x 47 cm x 10 cm)

Item Number	Description
QSXXXXHD	Other dimension "High profile" (specify width and depth)

*Largest size available is 13 x 18 cells. Smallest size available is 6 x 4 cells.

QUADTRO SELECT® LOW PROFILE Cushions with (HD) Incontinent Cover

Featuring ISOFLO® Memory Control. Ideal for active users; 2" cell height and four valves promote lateral, interior and posterior stability. See page 29 for the description of the waterproof cover.



Adult Standard Sizes	
Item Number	Description
QS88LPHD	15"w x 15"d x 2"h (38 cm x 38 cm x 5 cm)
QS89LPHD	15"w x 17"d x 2"h (38 cm x 43 cm x 5 cm)
QS99LPHD	17"w x 17"d x 2"h (43 cm x 43 cm x 5 cm)

Adult Standard Sizes	
Item Number	Description
QS910LPHD	17"w x 19"d x 2"h (43 cm x 47 cm x 5 cm)
QS109LPHD	19"w x 17"d x 2"h (47 cm x 43 cm x 5 cm)
QS1010LPHD	19"w x 19"d x 2"h (47 cm x 47 cm x 5 cm)

Item Number	Description
QSXXXXLPHD	Other dimension "Low profile" (specify width and depth)

*Largest size available is 13 x 18 cells. Smallest size available is 6 x 4 cells.

ENHANCER® Cushions (dual valve) with Cover

Recommended for enhanced pressure distribution, positioning and posture.



Adult Standard Sizes	
Item Number	Description
ENH88C	15"w x 15"d (38 cm x 38 cm)
ENH89C	15"w x 17"d (38 cm x 43 cm)
ENH99C	17"w x 17"d (43 cm x 43 cm)
ENH910C	17"w x 19"d (43 cm x 47 cm)
ENH109C	19"w x 17"d (47 cm x 43 cm)
ENH1010C	19"w x 19"d (47 cm x 47 cm)

Pediatric Sizes	
ENH66C	12"w x 12"d (30 cm x 30 cm)
ENH67C	12"w x 13"d (30 cm x 33 cm)

Item Number	Description
ENHXXXXC	ENHANCER custom design cushion (specify width and depth)

CONTOUR SELECT® Cushions with Cover

Recommended for enhanced pressure distribution and posture with the advantage of ISOFLO® Memory Control.



Adult Standard Sizes	
Item Number	Description
CS88C	15"w x 15"d (38 cm x 38 cm)
CS89C	15"w x 17"d (38 cm x 43 cm)
CS99C	17"w x 17"d (43 cm x 43 cm)

Adult Standard Sizes	
Item Number	Description
CS910C	17"w x 19"d (43 cm x 47 cm)
CS109C	19"w x 17"d (47 cm x 43 cm)
CS1010C	19"w x 19"d (47 cm x 47 cm)

Item Number	Description
CSXXXXCUST	CONTOUR SELECT custom design cushion (specify width and depth)
C-ROHO	Custom design cushion (provide diagram)

ROHO® “Smart Check” Cushions with Electronic Cushion Monitoring System

This cushion is equipped with a pressure monitor that facilitates proper inflation and allows users to achieve optimal skin and soft tissue protection. Size available for “HIGH PROFILE” (cell height 4”): Suitable for chairs with a width between 11.5” to 26.75”/29 cm to 68 cm; Size available for “MID PROFILE” (cell height 3”): Suitable for chairs with a width between 13.25” to 20”/33.5 cm to 50.5 cm.

SMARTCHECK-01 “Smart Check” pressure monitor for replacement only.



HIGH PROFILE (4”) with standard cover	
Item Number	Item Number
1R66H-CA-SC	1R1010H-CA-SC
1R67H-CA-SC	1R1011H-CA-SC
1R77H-CA-SC	1R1012H-CA-SC
1R78H-CA-SC	1R119H-CA-SC
1R79H-CA-SC	1R1110H-CA-SC
1R87H-CA-SC	1R1111H-CA-SC
1R88H-CA-SC	1R1112H-CA-SC
1R89H-CA-SC	1R129H-CA-SC
1R810H-CA-SC	1R1210H-CA-SC
1R98H-CA-SC	1R1211H-CA-SC
1R99H-CA-SC	1R1212H-CA-SC
1R910H-CA-SC	1R1310H-CA-SC
1R911H-CA-SC	1R1311H-CA-SC
1R108H-CA-SC	1R1411H-CA-SC
1R109H-CA-SC	1R1511H-CA-SC

MID PROFILE (3”) with standard cover	
Item Number	Item Number
1R77M-CA-SC	1R1010M-CA-SC
1R78M-CA-SC	1R1011M-CA-SC
1R87M-CA-SC	1R1012M-CA-SC
1R88M-CA-SC	1R119M-CA-SC
1R89M-CA-SC	1R1110M-CA-SC
1R810M-CA-SC	1R1111M-CA-SC
1R98M-CA-SC	1R1112M-CA-SC
1R99M-CA-SC	1R1210M-CA-SC
1R910M-CA-SC	1R1211M-CA-SC
1R911M-CA-SC	1R1212M-CA-SC
1R108M-CA-SC	1R1310M-CA-SC
1R109M-CA-SC	1R1311M-CA-SC

HIGH PROFILE (4”) with heavy duty (HD) cover	
Item Number	Item Number
1R66H-HD-SC	1R1010H-HD-SC
1R67H-HD-SC	1R1011H-HD-SC
1R77H-HD-SC	1R1012H-HD-SC
1R78H-HD-SC	1R119H-HD-SC
1R79H-HD-SC	1R1110H-HD-SC
1R87H-HD-SC	1R1111H-HD-SC
1R88H-HD-SC	1R1112H-HD-SC
1R89H-HD-SC	1R129H-HD-SC
1R810H-HD-SC	1R1210H-HD-SC
1R98H-HD-SC	1R1211H-HD-SC
1R99H-HD-SC	1R1212H-HD-SC
1R910H-HD-SC	1R1310H-HD-SC
1R911H-HD-SC	1R1311H-HD-SC
1R108H-HD-SC	1R1411H-HD-SC
1R109H-HD-SC	1R1511H-HD-SC

MID PROFILE (3”) with heavy duty (HD) cover	
Item Number	Item Number
1R77M-HD-SC	1R1010M-HD-SC
1R78M-HD-SC	1R1011M-HD-SC
1R87M-HD-SC	1R1012M-HD-SC
1R88M-HD-SC	1R119M-HD-SC
1R89M-HD-SC	1R1110M-HD-SC
1R810M-HD-SC	1R1111M-HD-SC
1R98M-HD-SC	1R1112M-HD-SC
1R99M-HD-SC	1R1210M-HD-SC
1R910M-HD-SC	1R1211M-HD-SC
1R911M-HD-SC	1R1212M-HD-SC
1R108M-HD-SC	1R1310M-HD-SC
1R109M-HD-SC	1R1311M-HD-SC

Sensor Ready Single Valve ROHO® HYBRID ELITE® Cushions with Cover

Combines the skin protection and shape matching ability of ROHO® DRY FLOATATION® with the stability of Jay® base technology. This new cushion design allows for improved pelvic alignment and stability for improved position and function, while maintaining optimal immersion and envelopment of both ischial tuberosities and greater trochanters. Note: Table on the right is for heavy duty (HD) incontinent covers. The single valve cushion comes Sensor Ready (SR) and can be added at any time. Pressure monitor is sold separately.



Single Valve - SENSOR READY		
Item Number	Width	Height
1RHE1414C-SR	14" / 36 cm	14" / 36 cm
1RHE1415C-SR	14" / 36 cm	15" / 38 cm
1RHE1416C-SR	14" / 36 cm	16" / 41 cm
1RHE1515C-SR	15" / 38 cm	15" / 38 cm
1RHE1516C-SR	15" / 38 cm	16" / 41 cm
1RHE1517C-SR	15" / 38 cm	17" / 43 cm
1RHE1615C-SR	16" / 41 cm	15" / 38 cm
1RHE1616C-SR	16" / 41 cm	16" / 41 cm
1RHE1617C-SR	16" / 41 cm	17" / 43 cm
1RHE1618C-SR	16" / 41 cm	18" / 46 cm
1RHE1620C-SR	16" / 41 cm	20" / 51 cm
1RHE1716C-SR	17" / 43 cm	16" / 41 cm
1RHE1717C-SR	17" / 43 cm	17" / 43 cm
1RHE1718C-SR	17" / 43 cm	18" / 46 cm
1RHE1816C-SR	18" / 46 cm	16" / 41 cm
1RHE1817C-SR	18" / 46 cm	17" / 43 cm
1RHE1818C-SR	18" / 46 cm	18" / 46 cm
1RHE1820C-SR	18" / 46 cm	20" / 51 cm
1RHE2017C-SR	20" / 51 cm	17" / 43 cm
1RHE2018C-SR	20" / 51 cm	18" / 46 cm
1RHE2020C-SR	20" / 51 cm	20" / 51 cm
1RHE2218C-SR	22" / 56 cm	18" / 46 cm
1RHE2220C-SR	22" / 56 cm	20" / 51 cm
1RHE2222C-SR	22" / 56 cm	22" / 56 cm
1RHE2418C-SR	24" / 61 cm	18" / 46 cm
1RHE2420C-SR	24" / 61 cm	20" / 51 cm

Single Valve - SENSOR READY with HD cover		
Item Number	Width	Height
1RHEHD1414C-SR	14" / 36 cm	14" / 36 cm
1RHEHD1415C-SR	14" / 36 cm	15" / 38 cm
1RHEHD1416C-SR	14" / 36 cm	16" / 41 cm
1RHEHD1515C-SR	15" / 38 cm	15" / 38 cm
1RHEHD1516C-SR	15" / 38 cm	16" / 41 cm
1RHEHD1517C-SR	15" / 38 cm	17" / 43 cm
1RHEHD1615C-SR	16" / 41 cm	15" / 38 cm
1RHEHD1616C-SR	16" / 41 cm	16" / 41 cm
1RHEHD1617C-SR	16" / 41 cm	17" / 43 cm
1RHEHD1618C-SR	16" / 41 cm	18" / 46 cm
1RHEHD1620C-SR	16" / 41 cm	20" / 51 cm
1RHEHD1716C-SR	17" / 43 cm	16" / 41 cm
1RHEHD1717C-SR	17" / 43 cm	17" / 43 cm
1RHEHD1718C-SR	17" / 43 cm	18" / 46 cm
1RHEHD1816C-SR	18" / 46 cm	16" / 41 cm
1RHEHD1817C-SR	18" / 46 cm	17" / 43 cm
1RHEHD1818C-SR	18" / 46 cm	18" / 46 cm
1RHEHD1820C-SR	18" / 46 cm	20" / 51 cm
1RHEHD2017C-SR	20" / 51 cm	17" / 43 cm
1RHEHD2018C-SR	20" / 51 cm	18" / 46 cm
1RHEHD2020C-SR	20" / 51 cm	20" / 51 cm
1RHEHD2218C-SR	22" / 56 cm	18" / 46 cm
1RHEHD2220C-SR	22" / 56 cm	20" / 51 cm
1RHEHD2222C-SR	22" / 56 cm	22" / 56 cm
1RHEHD2418C-SR	24" / 61 cm	18" / 46 cm
1RHEHD2420C-SR	24" / 61 cm	20" / 51 cm

Dual Valve ROHO® HYBRID ELITE® Cushions with Cover

Combines the skin protection and shape matching ability of ROHO® DRY FLOATATION® with the stability of Jay® base technology. This new cushion design allows for improved pelvic alignment and stability for improved position and function, while maintaining optimal immersion and envelopment of both ischial tuberosities and greater trochanters. Note: Table on the right is for heavy duty (HD) incontinent covers.



Dual Valve		
Item Number	Width	Height
2RHE1414C	14" / 36 cm	14" / 36 cm
2RHE1415C	14" / 36 cm	15" / 38 cm
2RHE1416C	14" / 36 cm	16" / 41 cm
2RHE1515C	15" / 38 cm	15" / 38 cm
2RHE1516C	15" / 38 cm	16" / 41 cm
2RHE1517C	15" / 38 cm	17" / 43 cm
2RHE1615C	16" / 41 cm	15" / 38 cm
2RHE1616C	16" / 41 cm	16" / 41 cm
2RHE1617C	16" / 41 cm	17" / 43 cm
2RHE1618C	16" / 41 cm	18" / 46 cm
2RHE1620C	16" / 41 cm	20" / 51 cm
2RHE1716C	17" / 43 cm	16" / 41 cm
2RHE1717C	17" / 43 cm	17" / 43 cm
2RHE1718C	17" / 43 cm	18" / 46 cm
2RHE1816C	18" / 46 cm	16" / 41 cm
2RHE1817C	18" / 46 cm	17" / 43 cm
2RHE1818C	18" / 46 cm	18" / 46 cm
2RHE1820C	18" / 46 cm	20" / 51 cm
2RHE2017C	20" / 51 cm	17" / 43 cm
2RHE2018C	20" / 51 cm	18" / 46 cm
2RHE2020C	20" / 51 cm	20" / 51 cm
2RHE2218C	22" / 56 cm	18" / 46 cm
2RHE2220C	22" / 56 cm	20" / 51 cm
2RHE2222C	22" / 56 cm	22" / 56 cm
2RHE2418C	24" / 61 cm	18" / 46 cm
2RHE2420C	24" / 61 cm	20" / 51 cm

Dual Valve with HD cover		
Item Number	Width	Height
2RHEHD1414C	14" / 36 cm	14" / 36 cm
2RHEHD1415C	14" / 36 cm	15" / 38 cm
2RHEHD1416C	14" / 36 cm	16" / 41 cm
2RHEHD1515C	15" / 38 cm	15" / 38 cm
2RHEHD1516C	15" / 38 cm	16" / 41 cm
2RHEHD1517C	15" / 38 cm	17" / 43 cm
2RHEHD1615C	16" / 41 cm	15" / 38 cm
2RHEHD1616C	16" / 41 cm	16" / 41 cm
2RHEHD1617C	16" / 41 cm	17" / 43 cm
2RHEHD1618C	16" / 41 cm	18" / 46 cm
2RHEHD1620C	16" / 41 cm	20" / 51 cm
2RHEHD1716C	17" / 43 cm	16" / 41 cm
2RHEHD1717C	17" / 43 cm	17" / 43 cm
2RHEHD1718C	17" / 43 cm	18" / 46 cm
2RHEHD1816C	18" / 46 cm	16" / 41 cm
2RHEHD1817C	18" / 46 cm	17" / 43 cm
2RHEHD1818C	18" / 46 cm	18" / 46 cm
2RHEHD1820C	18" / 46 cm	20" / 51 cm
2RHEHD2017C	20" / 51 cm	17" / 43 cm
2RHEHD2018C	20" / 51 cm	18" / 46 cm
2RHEHD2020C	20" / 51 cm	20" / 51 cm
2RHEHD2218C	22" / 56 cm	18" / 46 cm
2RHEHD2220C	22" / 56 cm	20" / 51 cm
2RHEHD2222C	22" / 56 cm	22" / 56 cm
2RHEHD2418C	24" / 61 cm	18" / 46 cm
2RHEHD2420C	24" / 61 cm	20" / 51 cm

AirLITE® Cushions with Cover

Combines the benefits of sealed air support pad and durable contoured foam base for comfort, enhanced positioning and pressure relief. Designed for client with tactile sensation and who are low risk for pressure ulcers.



Item Number	Width	Depth
AL1313HD	13" / 33 cm	13" / 33 cm
AL1315HD	13" / 33 cm	15" / 38 cm
AL1414HD	14" / 35.5 cm	14" / 35.5 cm
AL1416HD	14" / 35.5 cm	16" / 40.5 cm
AL1515HD	15" / 38 cm	15" / 38 cm
AL1517HD	15" / 38 cm	17" / 43 cm
AL1616HD	16" / 40.5 cm	16" / 40.5 cm
AL1618HD	16" / 40.5 cm	18" / 45.5 cm
AL1717HD	17" / 43 cm	17" / 43 cm

Item Number	Width	Depth
AL1718HD	17" / 43 cm	18" / 45.5 cm
AL1720HD	17" / 43 cm	20" / 51 cm
AL1816HD	18" / 45.5 cm	16" / 40.5 cm
AL1818HD	18" / 45.5 cm	18" / 45.5 cm
AL1820HD	18" / 45.5 cm	20" / 51 cm
AL1917HD	19" / 48.5 cm	17" / 43 cm
AL1919HD	19" / 48.5 cm	19" / 48.5 cm
AL2017HD	20" / 51 cm	17" / 43 cm
AL2018HD	20" / 51 cm	18" / 45.5 cm

ROHO has improved the AirLITE with you in mind!

- Lower anterior height for improved transfers
- Same great positioning, with even deeper leg troughs
- Built-in protection in a non-adjustable positioning cushion
- Softer, more compliant foam

MOSAIC® Cushions with Cover

Provides basic preventative skin and soft tissue protection. New design offers a more comfortable fit.

Item Number	Description
MOSAIC1616	16"w x 16"d (41 cm x 41 cm)
MOSAIC1618	16"w x 18"d (41 cm x 46 cm)
MOSAIC1717	17"w x 17"d (43 cm x 43 cm)
MOSAIC1816	18"w x 16"d (46 cm x 41 cm)
MOSAIC1818	18"w x 18"d (46 cm x 46 cm)
MOSAIC2018	20"w x 18"d (51 cm x 46 cm)



LTV Cushions with Cover

The valve allows you to customize the seat to your specific comfort level.

LTV 18" x 16" Office cushion with standard cover

LTV2B 18" x 16" Office cushion with ultra leather black cover



Solid Seat Inserts

Custom sizes available. Not recommended for individuals over 340 lb.

Contour Base (used with sling seat)	
Item Number	Description
BASE1515	15"w x 15"d (38 cm x 38 cm)
BASE1715	17"w x 15"d (43 cm x 38 cm)
BASE1717	17"w x 17"d (43 cm x 43 cm)

Planar Solid Seat Insert (used with sling seat)	
Item Number	Description
SI1416	14"w x 16"d (36 cm x 41 cm)
SI1616	16"w x 16"d (41 cm x 31 cm)
SI1618	16"w x 18"d (41 cm x 46 cm)
SI1818	18"w x 18"d (46 cm x 46 cm)
SI1820	18"w x 20"d (46 cm x 51 cm)
SI2020	20"w x 20"d (51 cm x 51 cm)
SI CUSTOM	Customized, please specify the size



MINI-MAX® Cushions with Cover (New expanded standard sizes)

Recommended for active individuals with low risk of skin breakdown & individuals seeking protection during recreational activities.

Item Number	Description
MINI1515C	15" w x 15" d x 1" (38 cm x 38 cm x 3 cm)
MINI1616C	16" w x 16" d x 1" (41 cm x 41 cm x 3 cm)
MINI1816C	18" w x 16" d x 1" (46 cm x 41 cm x 3 cm)
MINI1818C	18" w x 18" d x 1" (46 cm x 46 cm x 3 cm)
MINI2018C	20" w x 18" d x 1" (51 cm x 46 cm x 3 cm)
MINI2020C	20" w x 20" d x 1" (51 cm x 51 cm x 3 cm)



Special Application Cushions * Does not include cover

Unique design solutions for any seat and positioning problem.

ADAPTOR 9" x 13" x 1" Sealed 1 square inch air cell padding (can be cut to fit virtually any surface).

SACRAL82 Sacral support cushion (4"h x 15"w) (10 cm x 38 cm)

SACRAL92 Sacral support cushion (4"h x 17"w) (10 cm x 42 cm)

SACRAL102 Sacral support cushion (4"h x 19"w) (10 cm x 47 cm)

Shower Commode Cushion (low profile - 2" cell height only)

SEATCMD88LP (15" x 15") (38 cm x 38 cm)

SEATCMD89LP (15" x 17") (38 cm x 43 cm)

SEATCMD99LP (17" x 17") (43 cm x 43 cm)

SEATCMD910LP (17" x 18") (43 cm x 47 cm)

SEATCMD CUSTOM

SEATTOILET Toilet seat cushion; fits standard toilet seat and no cover available (no custom size available).

SEATTOILET Adjustable straps.

RECLINER-C Standard cover included free of charge. Recliner wheelchair cushion (made to fit per order, up to 320 cells overall).

CUST RECLINER-C "CUSTOM" Standard cover included free of charge. Recliner wheelchair cushion.

PACK-IT®-B with cover, designed to suit short term needs such as bathing and traveling (16" w x 9" d x 1" h) (41 cm x 23 cm x 3 cm).

HEAD/NECK PILLOW Cushion the head while minimizing the risk of pressure ulcer development with this DRY FLOATATION specialized pillow.



Accessories/Replacement Pumps/Repair Kits

Item Number	Description
HP-20	Inflation pump for ROHO mattresses, double action
HP-45	Inflation pump for ROHO cushions, NEW double action
A5	Rubber patch repair kit
MRK-1	Master repair kit
CV-3	Complete valve assembly (metal valve and 1" neoprene nose)
Factory labour (per repair)	

ROHO® Mattress Sections & Therapeutic Surfaces

MATTSECT3400 ROHO® neoprene section (19" x 34" / 48 cm x 85 cm)

MATTSECT1610 ROHO® neoprene section (19" x 30" / 48 cm x 76 cm)

MATTSECTCUST ROHO® neoprene section custom (up to 19" x 34" / 48 cm x 85 cm)

MATTSECTBARI neoprene section with snaps for BARI add-on (19" x 34" / 48 cm x 85 cm)

BARI310EA Bariatric mattress kit single section (3 x 10 cells)

BARI410EA Bariatric mattress kit single section (4 x 10 cells)

BARI310SET Bariatric kit system (3 x 10 cells) - 4 sections

BARI410SET Bariatric kit system (4 x 10 cells) - 4 sections

FOOTBOARDCUST Add-on length for neoprene sections (2, 3 or 4 cells by 18 cells)

Custom Neoprene Section size not listed above

3800-1 Prodigy® Mattress Overlay single section (27" x 35" / 68 cm x 89 cm)

SFLX2EA SOFFLEX 2® standard size section (1 section)

SFLX2NEA SOFFLEX 2® narrow system (1 section)



From left to right:

- ROHO® Prodigy® Mattress Overlay
- ROHO® SOFFLEX 2® Mattress Overlay System
- ROHO® Dry Floatation® Mattress Overlay System

Replacement Cushion Covers

Assist in slide transfers and minimizes wear on the cushion. Six-way stretch cover: top material made up of 80% nylon and 20% spandex; side material 100% polyester and bottom is 100% PVC coated polyester non-skid material. Universal Cover standard sizes (for single, dual valve, QUADTRO SELECT, ENHANCER and CONTOUR SELECT). Applicable to covers for cushions up to 160 cells.

**To order a standard replacement cover without ISO FLO opening please change the COV-SSA code to COV-A followed by the size. There is no additional charge.

HIGH or MID PROFILE	LOW PROFILE
COV-SSA88	COV-SSA88LP
COV-SSA89	COV-SSA89LP
COV-SSA99	COV-SSA99LP
COV-SSA910	COV-SSA910LP
COV-SSA109	COV-SSA109LP
COV-SSA1010	COV-SSA1010LP

LTV
COV-LTV1816 - Standard
COV-LTVLTHR1816 - Black ultra leather

MINI MAX
COV-MINI15
COV-MINI16
COV-MINI18
COV-MINI20

AirLITE	
COV-HDAL1313	COV-HDAL1718
COV-HDAL1315	COV-HDAL1720
COV-HDAL1414	COV-HDAL1816
COV-HDAL1416	COV-HDAL1818
COV-HDAL1515	COV-HDAL1820
COV-HDAL1517	COV-HDAL1917
COV-HDAL1616	COV-HDAL1919
COV-HDAL1618	COV-HDAL2017
COV-HDAL1717	COV-HDAL2018

Incontinent Cushion Covers

Protect your ROHO cushion from incontinent incidents. Cover material protects against micro-fungi and bacteria. Allows for easy cleaning to ensure a clean, sterile sitting environment for the user (for single, dual valve, QUADTRO SELECT, ENHANCER and CONTOUR SELECT).

STANDARD SIZES (match with cushion size)

CUSTOM SIZES (not listed) up to 160 cells

LTC waterproof		
COV-HD66	COV-HD89	COV-HD910
COV-HD67	COV-HD98	COV-HD911
COV-HD77	COV-HD99	COV-HD1010
COV-HD78	COV-HD108	COV-HD1011
COV-HD87	COV-HD109	COV-HD1110
COV-HD88	COV-HD810	

LTC waterproof "LOW PROFILE"		
COV-HD66LP	COV-HD89LP	COV-HD910LP
COV-HD67LP	COV-HD98LP	COV-HD1010LP
COV-HD77LP	COV-HD99LP	
COV-HD78LP	COV-HD108LP	
COV-HD87LP	COV-HD109LP	
COV-HD88LP	COV-HD810LP	

ROHO® HYBRID ELITE® Cushion Covers

STANDARD SIZES Two-way stretch cover with encased comfort foam cells, conforms to the individual's body shape. Solid, non-skid base material offers a durable easy to clean surface.

HEAVY DUTY Protect your HYBRID ELITE ROHO cushion from incontinent incidents. Cover material protects against micro-fungi and bacteria. Allow for easy cleaning to ensure a clean, sterile sitting environment.



STANDARD SIZES		
Item Number	Width	Height
COV-HE1414	14" / 36 cm	14" / 36 cm
COV-HE1415	14" / 36 cm	15" / 38 cm
COV-HE1416	14" / 36 cm	16" / 41 cm
COV-HE1515	15" / 38 cm	15" / 38 cm
COV-HE1516	15" / 38 cm	16" / 41 cm
COV-HE1517	15" / 38 cm	17" / 43 cm
COV-HE1615	16" / 41 cm	15" / 38 cm
COV-HE1616	16" / 41 cm	16" / 41 cm
COV-HE1617	16" / 41 cm	17" / 43 cm
COV-HE1618	16" / 41 cm	18" / 46 cm
COV-HE1620	16" / 41 cm	20" / 51 cm
COV-HE1716	17" / 43 cm	16" / 41 cm
COV-HE1717	17" / 43 cm	17" / 43 cm
COV-HE1718	17" / 43 cm	18" / 46 cm
COV-HE1816	18" / 46 cm	16" / 41 cm
COV-HE1817	18" / 46 cm	17" / 43 cm
COV-HE1818	18" / 46 cm	18" / 46 cm
COV-HE1820	18" / 46 cm	20" / 51 cm
COV-HE2017	20" / 51 cm	17" / 43 cm
COV-HE2018	20" / 51 cm	18" / 46 cm
COV-HE2020	20" / 51 cm	20" / 51 cm
COV-HE2218	22" / 56 cm	18" / 46 cm
COV-HE2220	22" / 56 cm	20" / 51 cm
COV-HE2222	22" / 56 cm	22" / 56 cm
COV-HE2418	24" / 61 cm	18" / 46 cm
COV-HE2420	24" / 61 cm	20" / 51 cm

WATERPROOF COVERS	
Item Number	Description
COV-HD77	14"w x 14"d (36 cm x 36 cm)
COV-HD78	14"w x 15"d (36 cm x 38 cm)
COV-HD79	14"w x 16"d (36 cm x 41 cm)
COV-HD88	15"w x 15"d (38 cm x 38 cm)
COV-HD89	15"w x 16-17"d (38 cm x 41-43 cm)
COV-HD98	16"w x 15"d (41 cm x 38 cm)
COV-HD99	16"w x 16-17"d (41 cm x 41-43 cm)
COV-HD910	16"w x 18"d (41 cm x 46 cm)
COV-HD911	16"w x 20"d (41 cm x 51 cm)
COV-HD109	17-18"w x 16-17"d (43-46 cm x 41-43 cm)
COV-HD1010	17-18"w x 18"d (43-46 cm x 46 cm)
COV-HD1011	17-18"w x 20"d (43-46 cm x 51 cm)
COV-HD119	20"w x 17"d (51 cm x 43 cm)
COV-HD1110	20"w x 18"d (51 cm x 46 cm)
COV-HD1111	20"w x 20"d (51 cm x 51 cm)
COV-HD1210	22"w x 18"d (56 cm x 46 cm)
COV-HD1211	22"w x 20"d (56 cm x 51 cm)
COV-HD1212	22"w x 22"d (56 cm x 56 cm)
COV-HD1310	24"w x 18"d (61 cm x 46 cm)
COV-HD1311	24"w x 20"d (61 cm x 51 cm)

ROHO® Hybrid Select Cushion Covers

The 4-way stretch, breathable, Stretch-Air fabric ensures greater comfort for end users while offering an easy-to-clean moisture barrier that protects the underlying foam (match with cushion size listed below).



Hybrid Select Stretch-Air cover			
COV-HS1414	COV-HS1616	COV-HS1718	COV-HS2016
COV-HS1416	COV-HS1617	COV-HS1720	COV-HS2017
COV-HS1515	COV-HS1618	COV-HS1816	COV-HS2018
COV-HS1516	COV-HS1620	COV-HS1817	COV-HS2020
COV-HS1517	COV-HS1716	COV-HS1818	
COV-HS1615	COV-HS1717	COV-HS1820	

Size Table for ROHO® Cushions

To select the correct cushion for your wheelchair:


STEP 1 - Measure the wheelchair seat width and depth. Note: Wheelchair manufacturers' published seat dimensions may vary.

STEP 2 - Use your wheelchair seat dimensions and the table below to help you determine the number of cells in cushion width and in depth.


STEP 3 - Select cushion height: LOW PROFILE (2.25"/5.5 cm), MID PROFILE (3.25"/8.5 cm), HIGH PROFILE (4.25"/10.5 cm), ENHANCER or CONTOUR SELECT (various cell heights).



Standard DRY FLOATATION® Single (1R) or Dual (2R) Valve and ENHANCER® (ENH)				
# of Cells	A & C cushion cell to cell	B & D cushion base	Measured chair range (narrowest to widest)	
5	8.25" / 21 cm	10" / 25 cm	8.25" / 21 cm	10" / 25.5 cm
6	10" / 25.5 cm	11.75" / 29 cm	10" / 25.5 cm	11.75" / 29.5 cm
7	11.75" / 29.5 cm	13.25" / 33.5 cm	11.75" / 29.5 cm	13.25" / 34 cm
8	13.25" / 34 cm	15" / 37.5 cm	13.25" / 34 cm	15" / 38 cm
9	15" / 38 cm	16.75" / 42 cm	15" / 38 cm	16.75" / 42.5 cm
10	16.75" / 42.5 cm	18.5" / 46.5 cm	16.75" / 42.5 cm	18.5" / 47 cm
11	18.5" / 47 cm	20.25" / 50.5 cm	18.5" / 47 cm	20.25" / 51 cm
12	20.25" / 51 cm	21.75" / 55 cm	20.25" / 51 cm	21.75" / 55.5 cm
13	21.75" / 55.5 cm	23.5" / 59.5 cm	21.75" / 55.5 cm	23.5" / 60 cm
14	23.5" / 60 cm	25.25" / 63.5 cm	23.5" / 60 cm	25.25" / 64 cm
15	25.25" / 64 cm	27" / 68 cm	25.25" / 64 cm	27" / 68 cm



SELECT® Series QUADTRO SELECT (QS) and CONTOUR SELECT (CS)			
# of Cells	A & C cushion cell to cell	B & D cushion base	Measured chair range (narrowest to widest)
5	8.25" / 21 cm	11" / 28 cm	8.25" / 21 cm 10" / 25.5 cm
6	10" / 25.5 cm	12.75" / 32.5 cm	10" / 25.5 cm 11.75" / 29.5 cm
7	11.75" / 29.5 cm	14.5" / 36.5 cm	11.75" / 29.5 cm 13.25" / 34 cm
8	13.25" / 34 cm	16" / 41 cm	13.25" / 34 cm 15" / 38 cm
9	15" / 38 cm	17.75" / 45 cm	15" / 38 cm 16.75" / 42.5 cm
10	16.75" / 42.5 cm	19.5" / 49.5 cm	16.75" / 42.5 cm 18.5" / 47 cm
11	18.5" / 47 cm	21.25" / 54 cm	18.5" / 47 cm 20.25" / 51 cm
12	20.25" / 51 cm	23" / 58 cm	20.25" / 51 cm 21.75" / 55.5 cm
13	21.75" / 55.5 cm	24.5" / 62.5 cm	21.75" / 55.5 cm 23.5" / 60 cm
14	23.5" / 60 cm	26.25" / 67 cm	23.5" / 60 cm 25.25" / 64 cm



Note: Cell to cell and edge to edge dimensions are rounded to the nearest .5 cm. Please check price list to confirm that a size is available. Sizes vary per model.

Sizing Exceptions - The above sizes are compatible with all ROHO seat cushions **except AirLITE®** cushions. The number of cells is compatible with the model numbers listed next to each product. The wheelchair seat format should be used to determine the proper number of cells required to suit your seat. For the AirLITE® cushions, the model number indicates the dimensions of the wheelchair seat in inches, width X length.

Dimensions for ROHO® Products

Read the measurements in this guide and in the product documentation at www.permobil.com.



Cushion Height – The height of the cushion without load and without possible edging or seat shell. The point of departure is a normally inflated cushion. These measurements may also be found in the product descriptions and on product labels.



Cushion Width – Cushion width measured from edge to edge without load. These measurements may also be found in the product descriptions and on product labels.



Cushion Depth – Cushion depth measured from edge to edge without load. These measurements may also be found in the product descriptions and on product labels.



Seat Width – The appropriate seat width for using the cushion. See also “appropriate width” in the list of prices.



Seat Depth – The appropriate seat width for using the cushion. See also “appropriate depth” in the list of prices.



Effective height of DRY FLOATATION® / AIR FLOATATION™ – Spacing between the lower portion of the body and the top of the seat cushion or mattress. The effective height of DRY FLOATATION® (height for AIR FLOATATION products with AIR FLOATATION™ technology) makes assumptions regarding the depth to which the body sinks into the cushion. The deeper, the greater the contact surface and the lower the pressure per cm² of skin.

Notes

NOTE 1: Ventilation/drainage openings through base of cushion available. See Custom Price.

NOTE 2: All cushion prices include a cover, inflation pump, repair kit, and operating instructions, except where noted by a *.

NOTE 3: All measurements for cushions are approximate from edge to edge.

NOTE 4: Custom pricing is for cushions up to 160 cells. Cushions over 160 cells require a quote.

NOTE 5: Custom items are non-refundable and non-exchangeable.

Terms and Policies

1. Return Policy: All returns require prior authorization and are subject to a restocking fee.
2. All custom items are non-refundable and non-exchangeable.
3. All cushion and mattress prices include cushion or mattress, inflation pump, rubber patch repair kit, and operating instructions. Cushion prices also include cover unless noted by *.

Permobil reserves the right to make changes to prices and specifications without notice.



For Internal Use	
Customer Order #	
Line #	
CSR	

		Front of Cushion													
		1	2	3	4	5	6	7	8	9	10	11	12	13	14
Valve		A	B	C	D	E	F	G	H	I	J	K	L	M	N
Deleted 	1														
	2														
Cut-Out 	3														
	4														
Valve & 	5														
	6														
	7														
	8														
Compartment 	9														
	10														
	11														
VELCRO 	12														
	13														
	14														

LENGTH

1	3	3	3	3	3	3
2	3	3	3	3	3	3
3	3	3	3	3	3	3
4	3	3	3	3	3	3
5	3	3	X	X	3	3
6	3	3	X	X	3	3

33

