

Permobil

Progeo Noir 2.0

Performance and design in perfect harmony



Carbon



Rigid



Minimized
weight



Made-to-
measure

permobil

The perfect fit

Progeo Noir 2.0 combines the benefits of carbon with the superior manoeuvrability of a made to measure wheelchair. This together gives the user an optimal balanced ride allowing for energy saving throughout the day.



Progeo Noir 2.0

Seat width	33 - 36 - 39 - 42 - 45 cm
Seat depth	35 - 37,5 - 40 - 42,5 - 45 - 47,5 cm
Seat height	Made-to-measure
Footplate distance	Made-to-measure (standard) Adjustable (on request)
Backrest height	From 24,5 to 47 cm (adjustable)
Front frame angle	85° - 90° - 95°
Total length	85 cm approx. (standard configuration)
Total width	Seat width + 17 cm (with 0° camber)
Backrest angle	From 74° to 94° to the seat (adjustable)
Camber	0° - 2° - 4°
Weight without rear wheels	Starting from 4,0 kg approx.
Wheelchair weight	Starting from 6,6 kg approx.
Maximum load	125 kg

Personalise your ride

CARBON
ESSENCE

LIME
ESSENCE

RED
ESSENCE

BLUE
ESSENCE

OPAQUE
ESSENCE

Actual product colors may vary.

Permobil reserves the right to change specifications and design without prior notice.