




Retail order form - February 2023

Use dimension chart to select appropriate size.

All Rifton products are made to order and are non-returnable and non-exchangeable.

Compass Chair dimensions

User dimensions (inches)	R310 size 1	R320 size 2	R330 size 3	R340 size 4	R350 size 5
Height	36 – 43	43 – 48	48 – 55	55 – 63	63+
Knee to heel	9 – 11	10½ – 12½	12½ – 14½	14½ – 16½	16½ – 18½
 <p>Key user dimension: height / knee to heel The user's overall height and measurement from knee to heel is a general guide to help you select the appropriate chair. Choose the model that allows for growth.</p> <p>Important: Make sure that seat width, depth and height are adequate for each user, and that the user's weight does not exceed the maximum working load.</p>					
Item dimensions (inches)	R310	R320	R330	R340	R350
Seat height	9, 10, 11	10½, 11½, 12½	12½, 13½, 14½	14½, 15½, 16½	16½, 17½, 18½
Seat depth	9½	11	12¾	14¾	16½
Seat width	10½	12¾	14	16¾	18¾
Backrest height	9½	11	12½	14¾	16
Armrest height above the seat	5½	6	7	8¾	8¾
Armrest length	6½	7½	8½	9¾	11
Max. working load (lb)	100	120	175	220	250

	R310 size 1		R320 size 2		R330 size 3		R340 size 4		R350 size 5	
* Compass chair	R311	\$706 <input type="checkbox"/>	R321	\$764 <input type="checkbox"/>	R331	\$827 <input type="checkbox"/>	R341	\$885 <input type="checkbox"/>	R351	\$941 <input type="checkbox"/>
Seat belt	R312	\$95 <input type="checkbox"/>	R312	\$95 <input type="checkbox"/>	R332	\$95 <input type="checkbox"/>	R342	\$95 <input type="checkbox"/>	R342	\$95 <input type="checkbox"/>
Stability feet	R316	\$142 <input type="checkbox"/>	R316	\$142 <input type="checkbox"/>	R316	\$142 <input type="checkbox"/>	R316	\$142 <input type="checkbox"/>	R316	\$142 <input type="checkbox"/>
New! Pelvic harness	R313	\$324 <input type="checkbox"/>	R323	\$343 <input type="checkbox"/>	R333	\$361 <input type="checkbox"/>	R343	\$385 <input type="checkbox"/>	R353	\$404 <input type="checkbox"/>
Total	R310		R320		R330		R340		R350	

* Required