

TiLite X Order Form (Replaces Aero X)
Folding Frame

Effective October 2024 (New items in red)

PROVIDER INFORMATION		RTS INFORMATION			
Dealer:	Account #:	RTS:	Phone:		
PO:	Ordered By:	Email:			
Email:		SHIP TO INFORMATION			
Mark For:	User Weight:	Name:	Ship Via:		
TiLite Sales Rep:		Address:			
TiLite Information:	Promo Code:	City:	State:	Zip:	
Ship Date:	Order #:	Phone:	Fax:		

Please contact the Permobil Manual Customer Service Team at chairs.tilite@permobil.com - Phone: 800.736.0925 or Fax: 866.586.2413

FRAME STYLE

- Standard Aluminum Frame (Full Adjustability)..... \$2,725 K0005 E1236
 - 250 lb. weight limit..... STD
 - 300 lb. weight limit..... NCO
- Standard Aluminum Frame with Transit Tie-Down..... \$2,995 K0005 E1236
 - 250 lb. weight limit..... STD
 - 300 lb. weight limit..... NCO
- Heavy Duty Aluminum Frame (350 Pound (158 Kg) Weight Limit. Available with 18" to 22" Seat Widths and 16" to 20" Seat Depths) Requires Aluminum Casters – Provided as NCO)..... \$575

SMARTDRIVE POWER ASSIST SYSTEM



- Permobil Smart Drive MX2+* (N/A with rear seat heights less than 13.5")..... \$7,529 E0986
- SpeedControl Dial* Extended Mount **OR** SwitchControl Buttons..... NCO
*SpeedControl Dial mounting may be limited based on chair and wheel lock configuration. Please refer to table on page 8 for more information.
- SwitchControl Button with Mono Jack Option \$159
 - Buddy Button \$262
- Carrying Bag \$84
- Universal Padded Handle..... \$22
*SmartDrive apps available on Apple App Store and Google Play Store for Apple Watch Series 5 or newer and Samsung Galaxy Watch 4 or newer.

TRANSIT OPTIONS

- Transit Tie-Down (Consists of four Tie-Down Brackets. N/A with HD Frame, Fixed Height Non-Removable Armrest or Swing Away Tubular Armrest. (Please refer to the user manual for important information regarding the Transit Tie-Down Option)..... \$363
- Transport Bracket Option (Consists of four Tie-Down Brackets. Unoccupied Use Only. N/A with HD Frame, Fixed Height Non-Removable Armrest or Swing Away Tubular Armrest..... \$379

FRAME FINISHES

- | | | |
|---|--|---|
| <input type="checkbox"/> High Gloss Black | <input type="checkbox"/> TiLite Red | <input type="checkbox"/> Neon Pink |
| <input type="checkbox"/> Ocean Blue Metallic | <input type="checkbox"/> Tangerine Metallic | <input type="checkbox"/> Acid Green Metallic |
| <input type="checkbox"/> Midnight Blue Pearl | <input type="checkbox"/> White Pearl | <input type="checkbox"/> Black Lava |
| <input type="checkbox"/> Titanium Grey | <input checked="" type="checkbox"/> Pomegranate Red | <input checked="" type="checkbox"/> Blue Mist |
| <input checked="" type="checkbox"/> Indigo Purple | <input type="checkbox"/> Hyper Blue | <input type="checkbox"/> Flame Orange |
| <input type="checkbox"/> OD Green | <input checked="" type="checkbox"/> Silver Streak | <input checked="" type="checkbox"/> Lavender Dusk |
| <input checked="" type="checkbox"/> Black Galaxy | <input checked="" type="checkbox"/> Seafoam Green..... | STD |

FRAME MEASUREMENTS (Please see information on page 8 regarding frame measurements)

For odd size seat depths, frame depth will be 1" longer than seat sling.

		Seat Width A									
		14	15	16	17	18	19	20	21	22	
Seat Depth B	14	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>				
	15	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>				
	16	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
	17	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
	18	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	19		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	20			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

- Note: Configurations including Seat Widths 20"-22" \$559 E2201
- Configurations including 20" Seat Depth..... \$675 E2203



Permobil.com
2701 West Court Street | Pasco, Washington 99301
509.352.2245 Tel | 866.586.2413 Toll-Free
chairs.tilite@permobil.com

FRONT SEAT HEIGHT RANGES

Front Wheel Size	Front Seat Heights C
3" Front Wheel	13.5" - 18"
4" Front Wheel	14.5" - 18.5"
5" Front Wheel	15.5" - 19"
6" Front Wheel	16.5" - 19.5"
8" Front Wheel	18.5" - 20.5"

FRONT SEAT HEIGHT c Measure from floor to top of seat tube at front edge of seat sling.

- 13.5" 14" 14.5" 15" 15.5" 16" 16.5" 17" 17.5" 18" 18.5" 19" 19.5" 20" 20.5"STD

REAR SEAT HEIGHT RANGES

Rear Wheel Size	Rear Seat Heights D
22" (501)	13" - 18.5"
24" (540)	14" - 19.5"
25" (559)	14.5" - 20"
26" (590)	15" - 20"

REAR SEAT HEIGHT D Measure from floor to top of seat tube at back post.

- 13" 13.5" 14" 14.5" 15" 15.5" 16" 16.5" 17" 17.5" 18" 18.5" 19" 19.5" 20"STD

CENTER OF GRAVITY E Measure from front of back post at the seat tube to the center of the rear axle.

- 1" 1.5" 2" 2.5" 3" 3.5" 4"STD

Amputee Axle Plate (Center of Gravity positioned 1" behind back canes. N/A with User Friendly Anti-tips. Limitations will apply)\$225 **E0959**

FRONT HANGERS

Composite Footrest with Std Fork						Depth Adj Footrest with Std Fork					
	3"	4"	5"	6"	8"		3"	4"	5"	6"	8"
85°	85°
80°	80°
70°	70°

Composite Footrest with Glide Fork						Depth Adj Footrest with Glide Fork					
	3"	4"	5"	6"	8"		3"	4"	5"	6"	8"
85°	85°
80°	80°
70°	70°

- 70° 80° 85°STD
 Omit HangersNCO

SEAT TO FOOTREST F Measure from front edge of seat sling to top rear of footrest. Measurement must be at least 2" shorter than Front Seat Height c.

Footrest Type

- Composite Flip-up FootrestSTD
 10.5"- 14" (High Mount Hanger) Set at\$202 **K0037**
 13.5"- 18.5" (Standard Hanger) Set atNCO
 Depth and Angle Adjustable Footrest (7"-8.5" will receive 3" extension tube - 9"-10" will receive 5" extension tube)\$293 **K0040**
 7"- 10" (Clamp-on High Mount) Set at\$202 **K0037**
 10.5"- 14" (High Mount Hanger) Set at\$202 **K0037**
 13.5"- 18.5" (Standard Hanger) Set atNCO
 Extra-long footplate (9" long)\$189 **K0041**

Elevating Leg Rest

Elevating Leg Rest, Single (When choosing a left or right ELR, please select the opposite hanger above. N/A with 14" Seat Width A)

- Install on Left **OR** Right\$250 **E0990**
 Elevating Leg Rest, Pair (N/A with 14" Seat Width A)\$425 **E0990**

Seat to Footrest/ Elevating Leg Rests F Measure from front edge of seat sling to top rear of ELR footrest.

- 14" 15" 16" 17" 18" 19" 20" 21"

INTEGRATED RESIDUAL LIMB SUPPORT Mounts in place of front hangers.

- Residual Limb Support.....\$449 E1020
 - 8" x 10" pad **OR** 8" x 14" pad Install on L **OR** R
- Residual Limb Support Pair.....\$850 E1020
 - 8" x 10" pad **OR** 8" x 14" pad.

FOOTREST ACCESSORIES

- Safe Tread Non-Skid Adhesive Patch.....\$20
- Heel Loops.....\$55 E0951

CALF STRAPS

- Velcro Adjustable.....\$44 K0038
 - Large (Recommended for 14"- 17" seat widths.)
 - Extra Large (Recommended for 16"- 22" seat widths.)
- Bodypoint® Padded Velcro Adjustable.....\$109 K0038
 - Medium (Recommended for 16"- 20" seat widths.)
 - Large (Recommended for 20"- 22" seat widths.)
- Buckle Adjustable
 - Black Nylon Webbing\$50 K0038

SEAT BACK CONFIGURATIONS

SEAT BACK TYPE

- Fixed Angle, Adjustable Height Seat Back.....STD
- Fixed Angle, Adjustable Height Seat Back w/ 8° Bend.....NCO
- Folding Adjustable-Angle and Height Seat Back (N/A with HD).....\$250
- Folding Stabilizer Bar.....\$423

PUSH HANDLES

- None.....NCO
- Integrated Push Handles.....NCO
- Integrated Fold-Down Push Handles (N/A with HD).....\$125
- Bolt-on height adjustable Push Handles (N/A with HD).....\$363

SEAT BACK HEIGHT G Measure from top of seat tube at rear of frame to top of backrest upholstery.

Fixed Angle, Adjustable Height Seat Back

Push Handles	Short	Tall	Set at: (Use 1" increments)
None	<input type="checkbox"/> 12.5"- 15.5"	<input type="checkbox"/> 16.5"- 20.5"	Set at: _____
Integrated	<input type="checkbox"/> 12"- 15"	<input type="checkbox"/> 16"- 20"	Set at: _____
Fold-Down	<input type="checkbox"/> 12"- 15"	<input type="checkbox"/> 16"- 20"	Set at: _____
Bolt-on Height Adj.	<input type="checkbox"/> 12.5"- 15.5"	<input type="checkbox"/> 16.5"- 20.5"	Set at: _____

Folding Angle, Adjustable Height Seat Back

Push Handles	Short	Medium	Tall	Set at: (Use .5" increments.)
None	12.5"- 14.5"	15"- 17"	17.5"- 20"	Set at: _____
Integrated	14"- 15.5"	16"- 17"	17.5"- 20"	Set at: _____
Fold-Down	14.5"-15.5"	16"- 17.5"	18"- 20.5"	Set at: _____

BACK UPHOLSTERY

- Standard Back Upholstery (Tension Adjustable. Includes Privacy Flap).....STD
- AIR Breathable Back Upholstery (Tension Adjustable. Includes Privacy Flap).....\$369 Pending E2611/E2612
 - Removable Lumbar Contour Pack Removable Lateral Contour Pack
- Omit Back Upholstery.....NCO
- 4" Military Service Patch.....\$115
 - Army Navy Air Force Marines Coast Guard Purple Heart

PERMOBIL BACK SUPPORTS (Please see information on page 8 regarding product shipment)

COMFORT ACTA-BACK® Back Support (Backs over 20" wide or tall subject to increase pricing)

- Shell Height 10" 12" 14" 16" 18" 20" 22".....\$623 E2613/E2614
 - Comfort-Tek Cover Stretch-Air Cover.....NCO
 - Comfort Privacy Flap.....\$44

TiLite X Order Form
FOLDING FRAME

COMFORT ACTA-BACK® Deep Back Support (Backs over 20" wide or tall subject to increase pricing)

Shell Height	<input type="checkbox"/> 10" <input type="checkbox"/> 12" <input type="checkbox"/> 14" <input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20" <input type="checkbox"/> 22"	\$738	E2615/E2616
	<input type="checkbox"/> Comfort-Tek Cover <input type="checkbox"/> Stretch-Air Cover	NCO	
	<input type="checkbox"/> 3" Deep Wing <input type="checkbox"/> 5" Deep Wing	NCO	
	<input type="checkbox"/> Comfort Privacy Flap	\$44	

COMFORT ACTA-EMBRACE® (Backs over 20" wide or tall subject to increase pricing)

Shell Height	<input type="checkbox"/> 10" <input type="checkbox"/> 12" <input type="checkbox"/> 14" <input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20" <input type="checkbox"/> 22"	\$649	E2613/E2614
	<input type="checkbox"/> Comfort-Tek Cover <input type="checkbox"/> Stretch-Air Cover	NCO	
	<input type="checkbox"/> Comfort Privacy Flap	\$44	

COMFORT ACTA-RELIEF

Shell Height	<input type="checkbox"/> 10" <input type="checkbox"/> 12" <input type="checkbox"/> 14" <input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20" <input type="checkbox"/> 22"	\$886	E2615
	<input type="checkbox"/> Comfort-Tek Cover <input type="checkbox"/> Stretch-Air Cover	NCO	
	<input type="checkbox"/> Comfort Privacy Flap	\$44	

SEAT UPHOLSTERY AND SEATING

<input type="checkbox"/> Tension Adjustable	STD
<input type="checkbox"/> Omit Seat Upholstery	NCO

PERMOBIL SEAT CUSHIONS (Please see information on page 8 regarding product shipment)

Comfort Cushions (Cushions over 20" wide subject to increase pricing)

<input type="checkbox"/> Comfort Curve Comfort-Tek Cushion	\$122	E2601/E2602
<input type="checkbox"/> Comfort Elements® Comfort-Tek Cushion (Odd size cushions subject to increased pricing)	\$112	E2601/E2602
<input type="checkbox"/> Comfort Acta-Embrace™ Comfort-Tek Zero Elevation Cushion	\$428	E2607/E2608
<input type="checkbox"/> Comfort M2 Zero Elevation Cushion with Comfort-Tek cover	\$474	E2607/E260

ROHO® Cushions (Cushions over 20" wide subject to increase pricing)

<input type="checkbox"/> ROHO Single Compartment <input type="checkbox"/> Low Profile <input type="checkbox"/> Mid Profile <input type="checkbox"/> High Profile	\$534	E2622/E2623
<input type="checkbox"/> ROHO Quadro Select <input type="checkbox"/> Low Profile <input type="checkbox"/> Mid Profile <input type="checkbox"/> High Profile	\$534	E2624/E2625
<input type="checkbox"/> ROHO Hybrid Elite <input type="checkbox"/> Single Compartment <input type="checkbox"/> Dual Compartment	\$552	E2622/E2623
<input type="checkbox"/> ROHO Hybrid Select	\$571	E2624

CUSTOM CONFIGURATION OPTIONS

FRONT WHEELS (Please reference table in front hanger section for compatibility)

3" Casters

<input type="checkbox"/> 3"x1.5" LiteSpeed Billet Aluminum Wheel w/ Soft Roll Tire <input type="checkbox"/> Silver Hub OR <input type="checkbox"/> Black Hub	\$205
---	-------

4" Casters

<input type="checkbox"/> 4"x1.5" LiteSpeed Plastic Wheel w/Soft Roll Tire	STD
<input type="checkbox"/> 4"x.75" Plastic Wheel w/Poly Tire	\$65
<input type="checkbox"/> 4"x.75" Light-Up Micro-Caster	\$125
<input type="checkbox"/> 4"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire <input type="checkbox"/> Silver Hub OR <input type="checkbox"/> Black Hub	\$217

5" Casters

<input type="checkbox"/> 5"x1.5" LiteSpeed Plastic Wheel w/Soft Roll Tire	STD
<input type="checkbox"/> 5"x1" Plastic Wheel w/Poly Tire	\$65
<input type="checkbox"/> 5"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire <input type="checkbox"/> Silver Hub OR <input type="checkbox"/> Black Hub	\$233

6" Casters

<input type="checkbox"/> 6"x1.5" LiteSpeed Plastic Wheel w/Soft Roll Tire	STD
<input type="checkbox"/> 6"x1" Plastic Wheel w/Poly Tire	\$65
<input type="checkbox"/> 6"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire <input type="checkbox"/> Silver Hub OR <input type="checkbox"/> Black Hub	\$244

8" Front Wheels

<input type="checkbox"/> 8"x1" Plastic Wheel w/Poly Tire	NCO
--	-----

FRONT FORKS

<input type="checkbox"/> TiLite Standard Fork	STD	
<input type="checkbox"/> Glide Suspension Fork (Max user weight 265 lbs)	\$526	E1015
Rider's Weight: <input type="checkbox"/> 125 lbs. <input type="checkbox"/> 125-180 lbs. <input type="checkbox"/> 180-265 lbs.		

TiLite X Order Form
FOLDING FRAME

REAR WHEELS (See Rear Seat Height Ranges Chart for Rear Seat Height D limitations)

ONE ARM DRIVE

- One Arm Drive (Requires 24" Rear Wheel and 0° Camber, N/A with HD.)\$1,298 E0958
Install on Left **OR** Right

STANDARD REAR WHEELS (See Rear Seat Height Ranges Chart for Rear Seat Height D limitations)

	22" (501)	24" (540)	25" (559)	26" (590)
TiLite Wheels	<input type="checkbox"/> STD	<input type="checkbox"/> STD	<input type="checkbox"/> STD	<input type="checkbox"/> STD

When ordering 24" wheel with treaded pneumatic or treaded tire with airless inserts and aluminum black or silver anodized handrims and long tabs, the chair will ship with TiLite standard rear wheel. Otherwise, the chair will ship TiLite Shadow rear wheel.

ALTERNATIVE REAR WHEELS (See Rear Seat Height Ranges Chart for Rear Seat Height D limitations.)

	22" (501)	24" (540)	25" (559)	26" (590)
X-Core Mag®	<input type="checkbox"/> \$225	<input type="checkbox"/> \$225		
Spinerge Spox <input type="checkbox"/> Black <input type="checkbox"/> Red <input type="checkbox"/> Blue <input type="checkbox"/> White	<input type="checkbox"/> \$1,195	<input type="checkbox"/> \$1,195	<input type="checkbox"/> \$1,195	<input type="checkbox"/> \$1,195
Spinerge LX <input type="checkbox"/> Black <input type="checkbox"/> Red <input type="checkbox"/> Blue <input type="checkbox"/> White <input type="checkbox"/> Bright Orange <input type="checkbox"/> Flamingo Pink <input type="checkbox"/> Purple	<input type="checkbox"/> \$1,295	<input type="checkbox"/> \$1,295	<input type="checkbox"/> \$1,295	<input type="checkbox"/> \$1,295
Spinerge Design Lab LX (Available with black spokes only.) <input type="checkbox"/> Gold Anodize <input type="checkbox"/> Green Anodize <input type="checkbox"/> Blue Anodize <input type="checkbox"/> Red Anodize <input type="checkbox"/> Purple Anodize <input type="checkbox"/> Orange Anodize	<input type="checkbox"/> \$1,995	<input type="checkbox"/> \$1,995	<input type="checkbox"/> \$1,995	<input type="checkbox"/> \$1,995

X-Core Mags requires Aluminum Silver Anodized Handrims with long tabs. N/A with Heavy Duty Frame, Quad Release Axles, Spoke Guards

PNEUMATIC REAR TIRES

- Treaded pneumatic (N/A with 25"(559) rear wheels).....STD
- Treaded with Airless Inserts (N/A with 25"(559) rear wheels).....\$145 E2213
- MBL SpeedLite.....STD
- Primo Sentinel.....\$159
- Schwalbe® Marathon Plus Evolution (N/A with Mag or 22" (501) wheels).....\$169
- Primo C-763 Knobby (Available on 24"(540) rear wheels).....\$129

SOLID REAR TIRES

- Primo Active (Available with 24"(540) and 25"(559) Rear Wheels).....\$160
- Primo Orion (Replaces Primo Express) (N/A with 26"(590) Rear Wheels).....\$160

CAMBER (Rear wheel spacing will be set to minimum value allowed by your configuration)

- 0° 2° 4°.....STD

AXLES

- Stainless Quick Release.....STD
- Quad Release (N/A with Mag, 24"(540)/25" (559) Spinerge Spox or Spinerge LX Rear Wheels).....\$169
- Fixed Non-removable.....\$179

HANDRIMS

- Aluminum – Silver Ano **OR** Black Ano.....STD
- Plastic Coated – Black.....\$174
- Omit Handrims.....NCO

The Natural-Fit® – Original, Larger Profile

- Standard Grip (Requires Long Tabs).....\$428
- Super Grip (Requires Long Tabs).....\$466
- Without Thumb Grip.....\$390

The Natural-Fit® LT – Reduced Weight, Smaller Profile

- Standard Grip (Requires Long Tabs).....\$450
- Super Grip (Requires Long Tabs).....\$514
- Without Thumb Grip.....\$423

The Surge (All-In-One Oval with Gription™ Strip)

- The Surge (Larger Profile).....\$495
- The Surge LT (Smaller Profile).....\$466

TiLite X Order Form
FOLDING FRAME

Effective October 2024 (New items in red)

Q Grip (High-Friction, Anti-Slip Coating)

Q-Grip\$495

TAB LENGTH

Long Tabs.....STD

Short Tabs (N/A with Mag Rear Wheels).....\$39

SPOKE GUARDS

Clear Lexan (N/A with Mag Rear Wheels.).....\$179 K0065

WHEEL LOCKS

Aluminum Wheel Locks

Push to LockSTD

Pull to LockSTD

Extension Handles for Aluminum Push/Pull to Lock.....\$75 E0961

Composite Wheel Locks

Compact Composite Push to Lock.....\$129

Compact Composite Pull to Lock\$129

Extension Handles for Compact Composite Push/Pull to Lock.....\$93 E0961

Compact Composite Scissor Lock.....\$158

Additional Wheel Locks

Grade Aids, Push to Lock (Requires pneumatic tires).....\$395 E0974

Hemi Wheel Lock – Left **OR** Right.....\$330

SIDE GUARDS

ALUMINUM SIDE GUARDS

Rigid Removable Small Standard.....\$217

Quad Removable w/ Hole Small Standard.....\$233

ARMRESTS

Omit.....NCO

Standard Armrest

Fixed Height Non-Removable Armrest (N/A with Heavy Duty).....STD

Swing Away Tubular Armrests (N/A HD, Bolt on push handles or Transit/Transport option)

Height Adjustable\$244 E0973

Range: 9" to 12" **OR** 12" to 15"

Height Adjustable with Sideguards (Side guard will be NCO) (N/A with Heavy Duty).....\$425 E0973

9" to 12" **OR** 12" to 15"

T-Style Armrests (Requires minimum seat depth of 15")

2-Point Flip Back Height Adjustable T-Style Armrest (N/A with Folding Angle, Adjustable Height Seat Back).....\$395 E0973

8" to 12" **OR** 10" to 14"

Desk Length Pad (10") **OR** Full Length Pad (14")

Removeable Height Adjustable T-Style Armrest\$375 E0973

8" to 12" **OR** 10" to 14"

Desk Length Pad (10") **OR** Full Length Pad (14")

POSITIONING BELTS

Velcro Adjustable, Black Nylon Webbing.....\$61 E0978

Bodypoint® 2" Wide Push-Button Buckle (Auto-Style).....\$103 E0978

Bodypoint® 1.5" Wide Push-Button Buckle Padded (Auto-Style, N/A with 19" -20" Rear Widths).....\$185 E0978

Bodypoint® 1" Wide Push-Button Buckle Padded (Auto-Style, N/A with 19" -20" Rear Widths).....\$196 E0978

Bodypoint® EvoFlex® Pelvic Stabilizer Push Button.....\$201 E0978

Q'Straint® WC-19 Wheelchair-Anchored Pelvic Belt (Recommended for use with Transit Tie Down).....\$217 E0978

Bodypoint® Optional Band Clamp Mounting Hardware.....\$70

ACCESSORIES

ANTI-TIPS

- Standard Rear Aluminum - Black.....\$174 E0971
- User Friendly Flip-up.....\$363 E0971

FUNCTIONAL ACCESSORIES

- Upper Extremity Support Tray
 - 12" to 14" OR 13" to 15" Seat Widths A (Requires T Style Armrests.).....\$315 E0950
- Luggage Carriers (N/A with High Mount Footrests.).....\$158
- Crutch Holder.....\$244 E2207
- Oxygen Tank Holder – Left OR Right.....\$330 E2208

SOFT GOODS ACCESSORIES

- Permobil Medical Necessities Bag.....\$103
- Permobil Medical Necessities Seat Pouch.....\$66
- Bodypoint® Mobility Bag.....\$109
- Neoprene Impact Guards, Black.....\$195

GENERAL

Permobil manufactures TiLite manual wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. HCPCS Codes are for reference purposes only. Actual coding and coverage varies by payer. Proper HCPCS coding is not a guarantee of payment. Please verify the appropriate payer's coverage policy.

Get the latest order forms at permobil.com. Prices, standard features, options and accessories are subject to change without notice. All rights reserved. Trademarks are identified by symbols™ and ®. TiLite, TiFit, the TiLite logo, Ti, the Ti logo, and TiShaft are registered trademarks of TiSport, LLC and are licensed to Permobil. Permobil® is a registered trademark of Permobil. Adobe®, Acrobat®, Acrobat Professional® and Reader® are registered trademarks of Adobe Systems Inc. Bodypoint® is a registered trademark of Bodypoint Inc. Kenda® is a registered trademark of Kenda Rubber Industrial Co. Ltd. ROHO®, AGILITY®, Low Profile® and Quadro Select® are registered trademarks of ROHO Inc. Schwalbe® is a registered trademark of Ralf Bohle GmbH. The Natural-Fit® and Out-Front® are registered trademarks of Three Rivers Holdings, LLC. The ACTA-BACK®, ACTA-EMBRACE™, COMFORT-TEK™, COMPASS®, CURVE®, ELEMENTS®, and STRETCH AIR™ are registered trademarks of THE COMFORT COMPANY® and are licensed to Permobil. Q'Straint® is a registered trademark of Q'Straint Inc.

NOTES


QUOTE THIS CHAIR

If you are using Adobe Acrobat, this order form may be converted into a quote request or into an order by clicking one of the buttons (this feature is available exclusively to Permobil dealers and distributors). A Quote or an Order Confirmation will be sent to you within 24 hours.

ORDER THIS CHAIR

Measurement Guide

- A** Seat Width - Measure from outside of seat tube at back post to the outside of the opposite seat tube at back post.
- B** Seat Depth - Measure from front of back post to the front edge of seat sling.
- C** Front Seat Height - Measure from floor to top of seat tube at front edge of seat sling.
- D** Rear Seat Height - Measure from floor to top of seat tube at back post.
- E** Center of Gravity - Measure from front of back post at the seat tube to the center of the rear axle.
- F** Seat to Footrest - Measure from front edge of seat sling to top rear of footrest.
- G** Seat Back Height - Measure from top of back post to top of seat tube at rear of frame



Seating and Positioning Shipping Information

The following Seating and Positioning products will ship with the chair. Cushions are stocked with Comfort Tek Cover and backs are stocked with Stretch Air Cover. All other Seating and Positioning products will be drop-shipped from the manufacturer. (Width x Depth)

Curve Cushion	Elements Cushion	Acta Embrace Cushion	M2 Cushion Zero Elevation	ROHO Single High Profile Cushion	ROHO Quadtro Cushion
18" x 18"	18" x 18"	18" x 18"	18" x 18"	18" x 18"	18" x 18"
16" x 16"	16" x 18"	16" x 18"	18" x 16"	16" x 18"	16" x 16"
16" x 18"	18" x 16"	16" x 16"	20" x 16"	16" x 16"	20" x 16"
			20" x 20"		20" x 18"
			16" x 16"		

Acta Back	Acta Embrace Back	Acta Relief Back	ROHO Agility Mid Active	ROHO Agility Min Active
18" x 12"	18" x 16"	18" x 16"	18" x 10"	18" x 10"
18" x 16"	16" x 16"	16" x 16"	16" x 10"	16" x 10"
16" x 16"		18" x 12"		
16" x 12"		16" x 12"		

SpeedControl Dial Compatibility Information

SpeedControl Dial Compatibility:

SpeedControl Dial mounting will require Aluminum or Composite push/pull wheel locks and a seat depth greater than 14". The following configurations are compatible with the SpeedControl Dial. Both wheel locks work with all 19" and 20" seat depth configurations. Please consider this when selecting your free wired controller for SmartDrive.

- : Aluminum push/pull wheel lock compatible
- : Composite push/pull wheel lock compatible

Seat depth	Center of gravity	Rear wheel sizes				Seat depth	Center of gravity	Rear wheel sizes			
		22"	24"	25"	26"			22"	24"	25"	26"
15" - 16"	1"	●	●	●	●	17" - 18"	1"	●○	●○	●○	●○
	1.5"	●	●	●	●		1.5"	●○	●○	●○	●○
	2"	●	●				2"	●○	●○	●○	●
	2.5"	●					2.5"	●○	●○	●	●
	3"	●					3"	●○	●	●	●
	3.5"						3.5"	●○	●		
	4"						4"	●	●		