



PDG Eclipse - AADL W320



Dealer Name*:		Account #*:
City, Province, Post Code:		
Phone:	Email:	ATP:
Rep Name:	P.O. Number:	Quote Number:
Client Reference*:	Client Weight:	
Special considerations that need to be addressed here (e.g., diagnosis):		

*Required field

October 1, 2026. Prices are in Canadian Dollars (CAD)

Wheelchair & Frame

☒ Eclipse Bariatric Manual Wheelchair	W320	\$ 3,500.00
---	------	-------------

Frame Type

☐ 600lbs Weight Capacity	No charge
☐ 1000lbs Weight Capacity <i>Not compatible with Transit Tie Down System</i>	\$ 300.00

Frame Options

☐ Transit Tie Down System (TTDS) <i>TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. Tested to max. weight capacity of 600lbs. Not compatible with O2 holder or IV pole.</i>	\$ 395.00
--	-----------

Seat Measurements

Seat Width | Overall chair width is +7.5"

Select Width:	☐ 20" Seat Width	☐ 22" Seat Width	☐ 24" Seat Width	☐ 26" Seat Width	\$ 675.00
	☐ 28" Seat Width	☐ 30" Seat Width			\$ 875.00
	☐ 32" Seat Width	☐ 36" Seat Width			\$ 1,100.00
	☐ 40" Seat Width	☐ 44" Seat Width	☐ 48" Seat Width		\$ 1,300.00

Seat Depth (Adjustable)

Select Depth:	☐ 16" Seat Depth	☐ 17" Seat Depth	☐ 18" Seat Depth	☐ 19" Seat Depth	☐ 20" Seat Depth	No charge
	☐ 21" Seat Depth	☐ 22" Seat Depth				\$ 400.00
	☐ 23" Seat Depth	☐ 24" Seat Depth				\$ 400.00

24" seat depth not available with TTDS.

Frame Finishes

Select Color:	☐ High Gloss Black	☐ Hyper Blue
---------------	---	-------------------------------------

Armrests

Armrest Type

☐ No Armrests	No charge
☐ Dual Post, Height Adjustable Armrests	\$ 295.00
☐ 10-14" Range ☐ 12"-16" Range	
☐ Dual Post, Fixed Height Armrests (12")	No charge

Armpads

☐ Regular Armpads	No charge
☐ Desk Length 2" x 10" ☐ Full Length 2" x 14"	

- ☐ **Comfort Armpads with Interlocking Hardware** *May interfere with wheels in some tilt configurations.* A1 \$ 496.00
☐ Desk Length 3.75" x 10" ☐ Full Length 3.75" x 14"
- ☐ **Gel Armpads** *May interfere with wheels in some tilt configurations.* A1 \$ 305.00
☐ Narrow 2" x 12" ☐ Wide 3.5" x 12"

Armrest Options

- ☐ **Forward Guards** \$ 200.00
- ☐ **Armrest Extensions - 3' Forward** *Not compatible with desk length armpads.* \$ 75.00

Seat Height

Seat Height | *Front and rear height can be maximum 1" apart.*

- Front Seat Height:** ☐ 14" STF Front ☐ 15" STF Front ☐ 16" STF Front ☐ 17" STF Front ☐ 18" STF Front
☐ 19" STF Front ☐ 20" STF Front \$ 200.00
- Rear Seat Height:** ☐ 14" STF Rear ☐ 15" STF Rear ☐ 16" STF Rear ☐ 17" STF Rear ☐ 18" STF Rear
☐ 19" STF Rear ☐ 20" STF Rear \$ 200.00

Seat Pan

- ☒ **Solid Aluminium Seat Pan** \$ 165.00

Seat Cushions

Cushions over 20" wide subject to increased pricing.

- ☐ **Comfort M2 Cushion Zero Elevation** \$ 580.00
☐ Comfort-Tek Cover ☐ Stretch Air Cover
- ☐ **Comfort M2 Cushion Anti-Thrust** \$ 580.00
☐ Comfort-Tek Cover ☐ Stretch Air Cover
- ☐ **Comfort Embrace Cushion Zero Elevation** \$ 508.00
☐ Comfort-Tek Cover ☐ Stretch Air Cover
- ☐ **Comfort Embrace Cushion Anti-Thrust** \$ 508.00
☐ Comfort-Tek Cover ☐ Stretch Air Cover
- ☐ **Comfort Saddle 7 Series Foam Cushion Zero Elevation** *Solace cover.* \$ 818.00
- ☐ **Comfort Saddle 7 Series Gel Cushion Zero Elevation** *Solace cover.* \$ 1,133.00
- ☐ **Comfort Support Pro 7 Series Foam Cushion Zero Elevation** *Solace cover.* \$ 445.00
- ☐ **Comfort Support Pro 7 Series Gel Cushion Zero Elevation** *Solace cover.* \$ 693.00
- ☐ **ROHO Quadtro Select High Profile Cushion** \$ 777.00
- ☐ **ROHO Sensor Ready High Profile Cushion** \$ 637.00

Front Casters

- ☒ **5" Casters with Heavy Duty Steel Fork** No charge

Rear Wheels

- ☒ **24" Heavy Duty Spoke Wheels** No charge

Rear Wheel Position: ☐ Position 1 (Rear) ☐ Position 2 ☐ Position 3 ☐ Position 4 (Front)

Positions 3 and 4 not compatible with TTDS. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

Axles

- ☒ **Quick Release Rear Axles** No charge

Rear Tires

- ☒ **Urethane Tires** No charge

Handrims

- ☐ **No Handrims** No charge
- ☐ **Standard Aluminium Handrims** No charge
- ☐ **Plastic Coated Handrims** No charge

Wheel Locks

- ☐ **User Wheel Locks** No charge
☐ Push to Lock ☐ Pull to Lock
- ☐ **Wheel Lock Extensions** No charge

Anti-Tippers

- ☐ **Rear Anti-Tippers** Required safety option. \$ 75.00

Backrest

- Back Post Style:** ☐ Stroller Handle Back ☐ Cane Handle Back
Cane handle not compatible with headrest mounting plate.
- Back Post Height:** ☐ Short (22") ☐ Tall (25")
Measured from seat pan to handle.
- Seat to Back Angle:** ☐ Set at 85° ☐ Set at 90° ☐ Set at 95° ☐ Set at 100°

Backrest Options

- ☐ **Endomorph Back** Backrest will be narrower than seat by total of 2x offset selected. Minimum backrest width is 12". Not available with reclining back. \$ 395.00
☐ 1" Offset Per Side ☐ 2" Offset Per Side \$ 395.00
- ☐ **Removable Backrest** No charge

Back Supports

Back Upholstery

- ☐ **No Back Upholstery** No charge
- ☐ **Tension Adjustable Back Upholstery** \$ 375.00

Position of Back Upholstery

Uph Type	22" Back Post						25" Back Post									
	Regular	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	
	Tall				☐	☐	☐			☐	☐	☐	☐	☐	☐	
		17"	18"	19"	20"	21"	22"	17"	18"	19"	20"	21"	22"	23"	24"	25"
		Upholstery Height														

Comfort Back Supports | Backs over 20" wide or tall subject to increased pricing.

- ☐ **Comfort Unite 1 Back** Includes 2-point quick release Compass hardware. Select shell height and cover. \$ 961.00
☐ 12" ☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22"
☐ Comfort-Tek Cover ☐ Diamond Spacer Cover
- ☐ **Comfort Unite 1 Back - 5" Wing** Includes 2-point quick release Compass hardware. Select shell height and cover. \$ 1,245.00
☐ 12" ☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22"
☐ Comfort-Tek Cover ☐ Diamond Spacer Cover
- ☐ **Comfort Unite 1 Back w/ROHO Insert - 5" Wing** Includes 2-point quick release Compass hardware. Select shell height and cover. \$ 1,499.00
☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22"
☐ Comfort-Tek Cover ☐ Diamond Spacer Cover
- ☐ **Comfort Unite 1 Back - 7" Wing** Includes 2-point quick release Compass hardware. Select shell height and cover. \$ 1,245.00
☐ 12" ☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22"
☐ Comfort-Tek Cover ☐ Diamond Spacer Cover
- ☐ **Comfort Unite 1 Back w/ROHO Insert - 7" Wing** Includes 2-point quick release Compass hardware. Select shell height and cover. \$ 1,499.00
☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22"
☐ Comfort-Tek Cover ☐ Diamond Spacer Cover
- ☐ **Comfort Unite 2 Back** Includes 2-point quick release Compass hardware. Select shell height and cover. \$ 1,106.00
☐ 12" ☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22"
☐ Comfort-Tek Cover ☐ Diamond Spacer Cover
- ☐ **Comfort Unite 2 Back w/ROHO Insert** Includes 2-point quick release Compass hardware. Select shell height and cover. \$ 1,254.00
☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22"
☐ Comfort-Tek Cover ☐ Diamond Spacer Cover

☐ Comfort Unite 3 Back <i>Includes 2-point quick release Compass hardware. Select shell height and cover.</i>	\$ 1,287.00
☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22"	
☐ Comfort-Tek Cover ☐ Diamond Spacer Cover	
☐ Comfort Unite 3 Back w/ROHO Insert <i>Includes 2-point quick release Compass hardware. Select shell height and cover.</i>	\$ 1,499.00
☐ 16" ☐ 18" ☐ 20" ☐ 22"	
☐ Comfort-Tek Cover ☐ Diamond Spacer Cover	
☐ Comfort Embrace Back <i>Includes 2-point quick release Compass hardware. Select shell height and cover.</i>	\$ 1,254.00
☐ 10" ☐ 12" ☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22"	
☐ Comfort-Tek Cover ☐ Diamond Spacer Cover	
☐ Comfort Relief Back <i>Includes Comfort-Tek cover and 2-point quick release Compass hardware. Select shell height.</i>	\$ 1,329.00
☐ 10" ☐ 12" ☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22"	

Positioning Accessories

Lateral Supports

☐ BodiLink Lateral Supports <i>Select size, cover and side.</i>	
☐ 3" x 5.5" ☐ 4" x 6.5" ☐ 6" x 6.5"	
☐ Comfort-Tek Cover ☐ Stretch Air Cover	
☐ Left ☐ Right	\$ 220.00
☐ Pair	\$ 440.00
☐ BodiLink Lateral Support Mounting Hardware	
☐ Individual	\$ 368.00
☐ Pair	\$ 736.00

Headrests

☐ BodiLink Headrest Pad <i>Select size and cover.</i>	\$ 260.00
☐ 10" Wide ☐ 12" Wide ☐ 14" Wide ☐ 18" Wide	
☐ Comfort-Tek Cover ☐ Stretch Air Cover	
☐ BodiLink Headrest Hardware, Removable	
☐ Mid 6 ☐ Mid 8	\$ 373.00
☐ Lat 6 ☐ Lat 8	\$ 489.00
☐ Max 8	\$ 545.00
☐ BodiLink Headrest Hardware, Flip-Down	
☐ Mid 6 ☐ Mid 8	\$ 419.00
☐ Lat 6 ☐ Lat 8	\$ 530.00
☐ Max 8	\$ 591.00

Recline

☐ Reclining Backrest (Manual) <i>30° maximum recline. Not compatible with endomorph back or cane handle back. Reclining backrest is 1" narrower than seat width. Back frame is not removable. Finished back post height is 3" taller than selected back post height. Not recommended for user weight above 450lbs. Maximum seat depth is 23".</i>	\$ 895.00
---	-----------

Front Rigging

☐ No Front Rigging	No charge
☐ Center Mount Front Rigging Bracket	\$ 300.00
Easy-Out Front Rigging Effortless insertion and removal. No hanger brackets remain on the frame.	
☐ Fixed 60°	\$ 95.00
☐ Universal Adjustable Angle (0°-60°) <i>Select Angle. Angle adjustment by loosening one bolt. Tool-free leg length adjustment.</i>	\$ 295.00
☐ 0° ☐ 10° ☐ 20° ☐ 30° ☐ 40° ☐ 50° ☐ 60°	
☐ AA High Mount Option <i>For shorter lower leg lengths of 8"-12.5".</i>	\$ 295.00
Swing Away Front Rigging Legacy design. Not compatible with Easy-Out.	
☐ Swing Away 60°	No charge
☐ Swing Away Elevating Legrests (ELR) <i>Includes calf pads.</i>	W106 \$ 375.00
☐ 1" Offset Hangers for Swing Away Front Rigging	\$ 75.00

Footplates

Aluminium Angle Adjustable Footplates

W107 \$ 165.00

Extension Tubes | *No selection required for Easy-Out AA High-Mount option.*

☐ **Standard Extension Tubes, Short**

No charge

☐ **Standard Extension Tubes, Long**

\$ 75.00

Footplate Height: ☐ Specify:

Ext	Short	8" - 13"	11" - 16"	12" - 16"
	Long	N/A	14" - 22"	14" - 19"
		High Mount Easy-Out	Standard Easy-Out	60°, 70°, ELR Swing-Away
		Legrest		

Front Rigging Accessories

☐ Oversized 9"x10" Calf Pad <i>Not compatible with seat widths less than 25"</i>	Basic Seating	\$ 695.00
☐ Calf Pads <i>Not compatible with Easy-Out front rigging.</i>	Basic Seating	\$ 305.00
☐ Gel Knee Pads, Easy-Out Angle Adjustable Front Rigging	Basic Seating	\$ 400.00
☐ Gel Knee Pads, Swing Away ELR Front Rigging	Basic Seating	\$ 250.00
☐ Length Adjustable Heel Loops	No charge	
☐ Standard Calf Strap	Basic Seating	\$ 50.00
☐ Bodypoint Padded Calf Strap <i>Available up to a maximum size of 32"</i>	Basic Seating	\$ 117.00
☐ Comfort Foot Box, Double <i>Without separator. Select seat width</i>		
☐ 16 ☐ 18"		\$ 454.00
☐ 20" ☐ 22"		\$ 644.00
☐ 24"		\$ 776.00
☐ Comfort Foot Box, Single		
☐ Left ☐ Right		\$ 285.00
☐ Pair		\$ 570.00

Positioning Belts

☐ Bodypoint Non-Padded Center Pull Push Button Belt, Large <i>Not compatible with gel padding. Includes reinforced premium calf pad mount.</i>		\$ 125.00
☐ Bodypoint 2 Point Padded Center Pull Push Button Belt <i>Not compatible with gel padding.</i>		\$ 215.00
☐ Bodypoint 4 Point Padded Center Pull Push Button Belt <i>Not compatible with calf pads or Easy-Out front rigging.</i>		\$ 175.00
☐ Bodypoint Non-Padded Chest Belt 46"		\$ 165.00
☐ Q'Straint Lap Belt <i>Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" belt length.</i>	Basic Seating	\$ 470.00

Miscellaneous Options

☐ IV Pole <i>Not compatible with TTDS.</i>	A1	\$ 545.00
☐ O2 Holder <i>Not compatible with TTDS.</i>	W110	\$ 695.00
☐ Medical Necessities Backpack <i>Black with Permobil logo.</i>	A1	\$ 245.00
☐ Adhesive Back Hook Strip <i>For cushion mounting, 2"x30".</i>	A1	\$ 15.00

General

Permobil manufactures PDG wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

Notes