

Customer and Order Information

Purchase Order *	Tag	v.0101
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in

Base Chair

☒ 1701000	Fuze JR 50° Pediatric Manual Tilt Wheelchair	\$ 4,300
---	---	----------

Frame Type & Weight Capacity

☒ 1702025	Adjustable/Growing Width Frame	250 lbs Weight Capacity.
---	---------------------------------------	--------------------------

Frame Options

☐ 1702222	Take-Apart Frame	The Take-Apart frame allows attendants to take the Fuze apart in seconds without tools. pdgmobility.com/take-apart	\$ 500
This feature reduces the weight of the heaviest piece of the wheelchair to 29 lbs or less. Perfect for transporting it in the trunk of a car or SUV. Take-Apart rear crossbar is not width adjustable. PDG will provide one no-charge crossbar in a different width over the lifetime of each chair ordered with adjustable width frame (can be ordered with the chair or later).			
			
☐ 1702300	Transit Tie Down System (TTDS)	TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details.	

Seat Width

Overall chair width is +7.5" for 12/16" RW and +8.5" for 20/22/24" RW

Adjustable Seat Widths

☐ 1703112	12"-15" Range - Set at 12"	☐ 1703113	12"-15" Range - Set at 13"	☐ 1703114	12"-15" Range - Set at 14"
☐ 1703115	12"-15" Range - Set at 15"				
☐ 1703215	15"-18" Range - Set at 15"	☐ 1703216	15"-18" Range - Set at 16"	☐ 1703217	15"-18" Range - Set at 17"
☐ 1703218	15"-18" Range - Set at 18"				

Seat Depth

Standard Seat Depths (Adjustable)

☐ 1704013	13" Seat Depth	☐ 1704014	14" Seat Depth	☐ 1704015	15" Seat Depth
☐ 1704016	16" Seat Depth	☐ 1704017	17" Seat Depth	☐ 1704018	18" Seat Depth

Frame Color

PDG Premium Colors

View color swatches on PDG's website.

pdgmobility.com/colors

☐ 1705001	Black Shimmer	☐ 1705015	Sky Blue	☐ 1705016	Poppy Red
☐ 1705003	Fine Texture Black	☐ 1705010	Lollipop Blue	☐ 1705011	Lollipop Red
☐ 1705017	Arctic Grey	☐ 1705018	Tropical Teal	☐ 1705019	Tangerine Orange

See next page for more options >



☐ 1705020 Olive Green	☐ 1705021 Lush Lilac	☐ 1705012 Bubble Gum Pink
☐ 1705022 Leaf Green	☐ 1705023 Boysenberry Purple	☐ 1705014 Daffodil Yellow
☐ 1705200 Custom Color Contact PDG for specifications.		
\$ 300		

Armrests

☐ 1706000 No Armrests		
☐ 1706111 Height Adjustable Armrests (Single Post) 8"-11" Range.		\$ 295
Armrest Height If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used. Seat Pan to top of Armrest: _____ Does not include height of Armpad		
☐ Locking Flip-up Armrests Not recommended for users with excessive strength or tone. Not compatible with Fold-down Back.		\$ 295
☐ 1706308 Set at 8"	☐ 1706309 Set at 9"	
☐ 1706310 Set at 10"	☐ 1706311 Set at 11"	

Armpads

Regular Armpads		
☐ 1707100 Desk Length - 2" x 10"	☐ 1707101 Full Length - 2" x 14"	
Waterfall Armpads		\$ 115
☐ 1707105 Desk Length - 2" x 10"	☐ 1707106 Full Length - 2" x 14"	
Wide Armpads May interfere with wheels in some tilt configurations.		\$ 200
☐ 1707110 Desk Length - 3" x 10"	☐ 1707111 Full Length - 3" x 12"	
Gel Armpads May interfere with wheels in some tilt configurations.		\$ 135
☐ 1707120 Full Length - 2" x 12"	☐ 1707121 Full Length - 3.5" x 12"	

Armrest Options

Armrest Wideners Move Armrests outward by up to 1" on each side. Seat size does not change.		
☐ 1707210 1" Total Additional Width		\$ 95
☐ 1707211 2" Total Additional Width		\$ 195

Seat Pan

☐ 1708000 No Seat Pan	
☐ 1708001 Solid Aluminum Seat Pan	\$ 135

Seat Height

Don't forget

Select desired Seat-to-Floor Height STF Height is adjustable.

5" Casters									
RW	13"	14"	15"	16"	17"	18"	19"	20"	
12"	☐	☐							
16"		☐	☐	☐					
20"				☐	☐				
22"					☐				
24"									

6" Casters & 8" Forks									
RW	13"	14"	15"	16"	17"	18"	19"	20"	
12"						☐	☐		
16"			☐	☐					
20"				☐	☐	☐			
22"					☐	☐	☐		
24"						☐	☐		

6" Casters & 6" Forks									
RW	13"	14"	15"	16"	17"	18"	19"	20"	
12"		☐				☐			
16"		☐	☐	☐					
20"				☐	☐	☐			
22"					☐	☐			
24"						☐			

8" Casters & 8" Forks									
RW	13"	14"	15"	16"	17"	18"	19"	20"	
12"						☐	☐	☐	
16"				☐				☐	
20"				☐	☐	☐			
22"					☐	☐	☐		
24"						☐	☐	☐	

Casters

Urethane Casters

☐ 1711051	5" x 1.25" Caster size	☐ 1711061	6" x 1.25" Caster size	☐ 1711081	8" x 1.25" Caster size	
☐ 1711063	6" x 1.50" Caster size					\$ 85
☐ 1711084	8" x 2.00" Caster size					\$ 135

Pneumatic Casters

☐ 1711161	6" x 1.25" Caster size					\$ 110
☐ 1711181	8" x 1.25" Caster size					\$ 120
☐ 1711184	8" x 2.00" Caster size					\$ 150

Rear Wheels

Composite Wheels

☐ 1712012	12" Wheel size	☐ 1712016	16" Wheel size	☐ 1712020	20" Wheel size	
☐ 1712022	22" Wheel size	☐ 1712024	24" Wheel size			

Spoke Wheels

☐ 1712120	20" Wheel size	☐ 1712122	22" Wheel size	☐ 1712124	24" Wheel size	\$ 75
----------------------------------	----------------	----------------------------------	----------------	----------------------------------	----------------	-------



Rear Wheel Position

Indicate Rear Wheel positioning preference. If left blank, the default setting will be used.

In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

☐ 1712401	Rear More stable	☐ 1712402	Mid	☐ 1712403	Front Less stable	
----------------------------------	------------------	----------------------------------	-----	----------------------------------	-------------------	--

Axle Type

☐ 1713000	Quick Release Rear Axles	
☐ 1713010	Fixed Rear Axles	

Rear Tires

☐ 1714000	Urethane Tires	
☐ 1714010	Pneumatic Tires	May result in a finished Seat Height 1/2" greater than stated. \$ 60
☐ 1714020	Pneumatic Tires with Flat-Free Inserts	May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels. \$ 95

Handrims

Handrims are available only with 20", 22", and 24" Rear Wheels.

☐ 1715000	No Handrims	
☐ 1715010	Standard Aluminum Handrims	
☐ 1715020	Plastic Coated Handrims	Not compatible with Drum Brakes. \$ 150
☐ 1715030	Projection Handrims	Composite Wheels only. 10 Projections, 30°, 1-1.8" Length. \$ 495

Wheel Locks

User Wheel Locks

☐ 1716010	Push to Lock	☐ 1716020	Pull to Lock	
☐ 1716030	Attendant Wheel Locks (Push to Lock)	Not compatible with Wheel Lock Extensions.		\$ 85
☐ 1716005	Foot Operated Wheel Lock	Rear-mounted unilateral single-activation lock. Compatible with 19" and 20" STF Height only.		\$ 595
☐ 1716300	Drum Brakes	Available with 12" Composite or 20" Spoke Rear Wheels. Urethane Tires only. 12" Wheels compatible with 13", 14", 17.5", 18.5", 19.5" STF Height only. 20" Wheels compatible with 16", 17", 18" STF Height only.	12" WHEELS COMING SOON	\$ 795
	User-activated Unilateral Wheel Lock	Pull to Lock.		\$ 795
☐ 1716310	Left Activation	☐ 1716311	Right Activation	

Wheel Lock Options

☐ 1716100	Wheel Lock Extensions	Not compatible with Attendant Wheel Lock. \$ 60
----------------------------------	-----------------------	---

Anti-Tippers

☒ 1717010	Rear Anti-Tippers	Required Safety Option. \$ 75
---	-------------------	-------------------------------

Back Posts

Back Post Style			
☐ 1718100	Stroller Handle Back	☐ 1718105	Cane Handle Back Not compatible with Headrest Mounting Plate.
Back Post Type			
☐ 1718200	Straight	☐ 1718205	8° Bend
Back Post Height Measured from the Seat Pan to the Handle.			
☐ 1718420	Short (17")	☐ 1718425	Tall (23")
Back Post Folding Fold-down Back is not compatible with Reclining Backrest.			
☐ 1719220	Fold-down	☐ 1719221	Non-folding
Seat to Back Angle			
☐ 1719090	Set at 90°	☐ 1719096	Set at 98°
☐ 1719102	Set at 106°		
☐ 1719108	Set at 114°	☐ 1719114	Set at 122°

Back Options

Angle/Height Adjustable Stroller Handle		Not compatible with Headrest Mounting Plate or Cane Handle Back.	\$ 150
☐ 1719010	Button on Outside	☐ 1719011	Button on Inside
☐ 1719050	Removable Cane Handle Push Bar	Not compatible with Stroller Handle Back.	\$ 140
☐ 1719212	Headrest Mounting Plate		\$ 95
Not compatible with Angle/Height Adjustable Stroller Handle, Cane Handle Back. Not required if mounting to a solid Back.			

Tilt Activation

Attendant Tilt only Select trigger mounting location below.			
☐ 1723011	Right Back Post	☐ 1723012	Left Back Post
☐ 1723400	Power Tilt 40° tilt range. Includes 2 batteries for off-chair charging. pdgmobility.com/power-tilt		\$ 1,795
	Toggle Switch	☐ 1723421	Installed on Right side
	Button Switch	☐ 1723431	Installed on Right side
	Stealth Egg Switch	☐ 1723440	Not installed
	Optional for Egg Switch:	☐ 1723445	Swing Away Mounting Hardware Not compatible with all Headrest makes and models.
			\$ 625
Specify Headrest to be used: _____			
Tilt Limit Optional. If no tilt limit is selected, the full tilt range will be available.			
☐ 1723220	Set at 20°	☐ 1723230	Set at 30°
☐ 1723240	Set at 40°		

Recline

☐ 1723310	Reclining Backrest (Manual)	\$ 895
30° maximum recline. Not compatible with Cane Handle Back, Fold-down Back, Requires 23" Angle/Height Adjustable Stroller Handle with 8° Bend. Not available with 15" or 16" Seat Depth.		

Front Rigging

☐ 1725000	No Front Rigging	☐ 1725001	Sector Blocks only
Easy-Out™ Front Rigging Effortless insertion and removal. No Hanger Bracket remains on frame. pdgmobility.com/easy-out			
☐ 1725410	Fixed 70°		\$ 95
☐ 1725420	Universal Angle Adjustable (0°-110°)	Set at angle: _____ (use increments of 10°)	\$ 295
Adjust angle by loosening one bolt. Tool-free leg-length adjustment. Front Rigging can interfere with casters in the 90-110° settings.			
Swing Away Front Rigging Legacy Design. Not compatible with Easy-Out.			
☐ 1725130	Swing Away 60°	☐ 1725140	Swing Away 70°
☐ 1725160	Swing Away 90°	Adjustable down to 3.5" Leg Length. Includes 11"-16" Extension Tubes. Not compatible with 8" Casters.	
☐ 1725200	Swing Away Elevating Legrests	Includes Calf Pads. Not compatible with 12" or 13" Seat Width.	
☐ 1725300	1" Offset Hangers for Swing Away		\$ 75
☐ 1725110	Standard Plug In 70° Front Rigging	Legacy Design. Not compatible with Easy-Out.	

Footplates

☐	1726010	Composite Footplates	Not available on 12" Seat Width.	
☐	1726020	Aluminum Angle Adjustable Footplates		\$ 165
☐	1726030	One Piece Fixed Footplate	Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.	\$ 350
☐	1726040	One Piece Flip-up Footplate	Requires selection of Composite or Aluminum Angle Adjustable.	\$ 395
☐	1726050	Interlocking Angle Adjustable Footplates		\$ 195

Footrest/Legrest Accessories

		Calf Pad	Not compatible with Gel Padding.	
☐	1727012	For Easy-Out™ Front Rigging	Includes reinforced premium Calf Pad mount.	\$ 200
☐	1727010	For Swing Away Front Rigging		\$ 135
☐	1727020	Gel Padding	Not compatible with Calf Pad or Easy-Out Front Rigging.	\$ 125
		Gel Knee Pad		
☐	1727032	For Easy-Out™ Angle Adjustable Front Rigging		\$ 200
☐	1727030	For Swing Away ELR Front Rigging		\$ 100
☐	1727040	Length Adjustable Heel Loops		\$ 50
☐	1727050	Bodypoint Padded Calf Strap		\$ 195

Front Rigging Extension Tubes

		Regular Extension Tubes	No selection required for 90° Front Rigging.	
☐	1728010	Short	14"-18" for Easy-Out Front Rigging 13"-17" for Swing Away and Standard Plug In Front Rigging	
☐	1728020	Long	18"-24" for Easy-Out Front Rigging 15"-19" for Swing Away and Standard Plug In Front Rigging	\$ 75
		Lower Leg Length	If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.	
		Seat Pan to Foot Plates:	_____	Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.

Miscellaneous Options

		Spoke Guards	Available for 20"-24" Composite Rear Wheels. Made from durable and lightweight aluminum.	\$ 195
☐	1730020	Black	☐ 1730021 Frame Color	
		Stealth Flip Down Abductor	Requires Solid Seat Pan.	\$ 650
☐	1730060	Small Pad	3"x4"x3"x2" size ☐ 1730062 Large Pad 5"x5"x4"x3" size	
☐	1730090	Hip Belt Mount	Provides improved mounting options for Bodypoint Hip Belt or other third-party belts.	\$ 55
☐	1730160	Clear Lap Tray (Slide On)	Available for Seat Widths of 17"-18" only. Not compatible with Contour Arm Trough.	\$ 490
☐	1730200	Q'Straint Lap Belt	Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" Belt Length.	\$ 245
☐	1730222	Wheelchair Backpack	Black with PDG Logo. 17"x14" size.	\$ 145
☐	1730170	Adhesive Back Hook Strip, 2" x 30"	For Cushion mounting.	\$ 10

On-Chair Seating & Positioning

 Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. pdgmobility.com/seating

Order Notes

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.