

TiLite T-Series Canadian Order Form

Effective August 1, 2020 (New Items in Red)

DUAL-TUBE RIGID FRAME

Page 1 of 8

DEALER INFORMATION		RTS INFORMATION	
Dealer:	Account #:	RTS:	Phone:
PO:	Ordered By:	Email:	
Email:		SHIP TO INFORMATION	
Mark For:	User Weight:	Name:	Ship Via:
TiLite Sales Rep:		Address:	
TiLite Information:	Promo Code:	City:	Province: Postal Code:
Ship Date:	Order #:	Phone:	Fax:

More options and frame modifications may be available through our Creative Mobility Division! Please contact customer service for details or [click here](#) to view options.

Please Email or Fax to the Permobil Manual Customer Service Team at tilitesales@permobil.com or 800.736.0925, option 2.



FRAME STYLE

- | | | |
|-------|---|---------|
| ATFS1 | ☐ Aero T Aluminum Frame (Full Adjustability. 265 Pound Weight Limit. For Rear Seat Width A and Seat Depth B of 15" or greater, use K0005. For Rear Seat Width A and Seat Depth B of less than 15", use E1235.)..... | \$3,660 |
| TAFS1 | ☐ TRA Titanium Frame (Full Adjustability. 265 Pound Weight Limit. For Rear Seat Width A and Seat Depth B of 15" or greater, use K0005. For Rear Seat Width A and Seat Depth B of less than 15", use E1235.)..... | \$4,715 |
| TRFS1 | ☐ TR Titanium Frame Ultimate TiFIT (Minimal Adjustability. 265 Pound Weight Limit. For Rear Seat Width A and Seat Depth B of 15" or greater, use K0005. For Rear Seat Width A and Seat Depth B of less than 15", use E1235.)..... | \$4,920 |
| THDR1 | ☐ Heavy Duty Frame (Not Available with TR Titanium Frame. 350 Pound Weight Limit. No charge when 20" Rear Seat Width A and 20" Seat Depth B are selected. Includes Tension Adjustable by Straps Back Upholstery. Requires Folding Aluminum Adjustable Height Backrest, Billet Aluminum Front Wheels, Shadow or Spinergy SPOX Rear Wheels. LX Rear Wheels only available with 0 or 2 degree camber. Not available with Slipstream Fork, Velcro Adjustable Upholstery, Open Loop Footrest, V-Front Taper or Glide™ Forks.)..... | \$445 |
| TQST1 | ☐ Q'Straint Adaptable Frame (350 Pound Weight Limit. No charge when 20" Rear Seat Width A and 20" Seat Depth B are selected. Includes Tension Adjustable by Straps Back Upholstery. Requires Folding Aluminum Adjustable Height Backrest, Billet Aluminum Front Wheels, Shadow or Spinergy SPOX Rear Wheels. LX Rear Wheels only available with 0 or 2 degree camber. Not available with Slipstream Fork, Velcro Adjustable Upholstery, Open Loop Footrest, V-Front Taper or Glide™ Forks.)..... | \$445 |
| TRFR1 | ☐ Reinforced Frame (No charge when 20" Rear Seat Width A and 20" B Seat Depth are selected.)..... | \$400 |

SMARTDRIVE POWER ASSIST SYSTEM



- | | | |
|---------|--|---------|
| TMAXSD1 | ☐ Permobil Smart Drive MX2+ (SwitchControl comes standard. N/A with rear seat heights less than 13.5", or 20" rear wheels).... | \$7,495 |
| TMAXSD6 | ☐ NEW PushTracker E2 OR ☐ ZMAXSD7 ☐ PushTracker | NCO |
| TMAXSD8 | ☐ Upgrade to SwitchControl with Mono Jack Option | \$194 |
| TMAXSD4 | ☐ Carrying Bag | \$102 |

FRAME FINISHES

- | | | |
|-------|---|-------|
| TFF1 | ☐ Satin (Available with Titanium Frame only)..... | STD |
| SCBT1 | ☐ Satin Finish Buffing Pads | \$20 |
| TFF3 | ☐ Tattooed Titanium™ (Tattoos on front frame tubes. Titanium frame finished in Satin. Tattoos will be trimmed and adjusted to fit your frame.)..... | \$360 |
| TFF45 | ☐ Need for Speed | |
| | | |
| TFF43 | ☐ Mosaic | |
| | | |
| TFF44 | ☐ Flower Garden | |
| | | |

FRAME FINISHES

- | | | | |
|-----------------------------|--|---|--|
| High Gloss Paint | ☐ TFF8 Super Wet Black | ☐ TFF10 TiLite Red | ☐ TFF29 Neon Pink |
| Metallic Paint | ☐ TFF14 Ocean Blue Metallic | ☐ TFF15 Tangerine Metallic | ☐ TFF30 Acid Green Metallic |
| Pearlescent Paint | ☐ TFF17 Midnight Blue Pearl | ☐ TFF18 White Pearl | |
| Matte Paint | ☐ TFF19 Black Lava | ☐ TFF27 Titanium Grey | |
| Candy Paint | ☐ TFF22 Ruby Red | ☐ TFF23 Azure Blue | ☐ TFF25 Electric Plum |
| Matte Textured Paint | ☐ TFF40 Hyper Blue | ☐ TFF41 Flame Orange | ☐ TFF42 OD Green..... |
- STD

COLOR ANODIZE PACKAGE

(Color Anodize Package includes: Aluminum Camber Tube (Color Anodize not available with Power Adaptable Camber Tube), Camber Mount, Footrest Clamp, Caster Mount, Footrest Clamp, Caster Fork, Caster Cap and Agility Carbon Back clamps. If no anodize color is selected, the default will be black.)

- | | |
|--|-------|
| ☐ TANO1 Black..... | STD |
| ☐ TANO5 Sapphire ☐ TANO3 Red Pepper ☐ TANO6 Burnt Orange ☐ TANO8 Ultraviolet ☐ TANO9 Fuchsia..... | \$200 |



Permobil.com
2701 West Court Street | Pasco, Washington 99301
509.586.6117 Tel | 800.545.2266 Toll-Free
509.586.2413 Fax | 866.586.2413 Toll-Free Fax

OF0029c T Series 0419RevP

TiLite T-Series Canadian Order Form

Effective August 1, 2020 (New Items in Red)

DUAL-TUBE RIGID FRAME

Page 2 of 8

TiFit: TAILORED INDIVIDUALIZED FRAME

REAR SEAT WIDTH ④ Measure from outside of seat tube at back post to the outside of the opposite seat tube at back post.

TSW1 ☐ 12" ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19".....**STD**
TSW2 ☐ 20".....**\$495**

FRONT SEAT WIDTH (SEAT TAPER) ① Measure from inside of front frame tube to inside of opposite front frame tube. (Required with V-Front End ⑤. Not Available with Standard Front End ⑥.)

TFSW1 ☐ No Taper.....**STD**
TFSW2 ☐ 1" (2.5cm) Taper (3" Narrower than Rear Seat Width ④. Includes Tension Adjustable by Straps Seat Upholstery.).....**\$280**
TFSW4 ☐ 2" (5.1cm) Taper (4" Narrower than Rear Seat Width ④. Includes Tension Adjustable by Straps Seat Upholstery.).....**\$280**

SEAT DEPTH ⑤ Measure from front of back post to the front edge of seat sling. Frame bend will start .5" in front of seat sling unless Custom Frame Depth is specified, or the combination of specifications affect the rolling dynamics of the chair.

TSD1 ☐ 12" ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19".....**STD**
TSD2 ☐ 20".....**\$540**
TSD3 ☐ Custom Frame Depth: _____" Frame Depth Longer or Shorter Than Seat Depth (Measure from front of seat sling to beginning of frame bend. Ranges from 0.25" to 3" and can be given as specific as 1/4" increments. A negative Custom Frame Depth or specified above 3" may be available through our Creative Mobility Division. A Custom Frame Depth of 20" will have an additional charge).....**NCO**

ERGONOMIC SEAT ① Measure from front of back post to beginning of upward bend of seat tube.

TERG1 ☐ TR Ergonomic Seat (Not available with Tension Adjustable Bolt On Seat Upholstery)
☐ 5" (12.7cm) ☐ 6" (15.2cm) ☐ 7" (17.8cm) ☐ 8" (20.3cm) ☐ 9" (22.9cm).....**\$425**
TERG2 ☐ AeroT and TRA Ergonomic Seat (Center of Gravity adjustment may be affected by this option. Not available with V-Front or Seat Taper, Seat Back that Locks When Folded, Desk Arms, Depth Adjustable Backrests, Fold Down Rigid Sideguards, Fendered Adjustable Sideguards, Solid Seat, Tension Adjustable Bolt On Seat Upholstery, or the Curve cushion.)
☐ 6" (15.2cm) (Minimum Seat Depth ⑤ 12") ☐ 7" (17.8cm) (Minimum Seat Depth ⑤ 12")
☐ 8" (20.3cm) (Minimum Seat Depth ⑤ 13")**\$520**

FRONT SEAT HEIGHT ⑥ Measure from floor to top of seat tube at front edge of seat sling.

TFSH1 ☐ 15" (38.1cm) ☐ 16" (40.6cm) ☐ 16.5" (41.9cm) ☐ 17" (43.2cm) ☐ 17.5" (44.5cm)
☐ 18" (45.7cm) ☐ 18.5" (47cm) ☐ 19" (48.3cm) ☐ 19.5" (49.5cm) ☐ 20" (50.8cm)
☐ 20.5" (52.1cm) ☐ 21" (53.3cm).....**STD**

REAR SEAT HEIGHT ⑥ Measure from floor to top of seat tube at back post.

Rear Wheel Size	Rear Seat Height
20"	12.5"-21" (31.8-53.3 cm)
22"	13.5"-21" (34.3-53.3 cm)
24"	14.5"-21" (15" on AeroT/TRA) (36.8-53.3 cm)
25"	15.0"-21" (38.1-53.3 cm)
26"	15.5"-21" (38.1-53.3 cm)

TRSH1 ☐ 12.5" (31.8cm) ☐ 13" (33cm) ☐ 13.5" (34.3cm) ☐ 14" (35.6cm) ☐ 14.5" (36.8cm)
☐ 15" (38.1cm) ☐ 15.5" (39.4cm) ☐ 16" (40.6cm) ☐ 16.5" (41.9cm) ☐ 17" (43.2cm)
☐ 17.5" (44.5cm) ☐ 18" (45.7cm) ☐ 18.5" (47cm) ☐ 19" (48.3cm) ☐ 19.5" (49.5cm)
☐ 20" (50.8cm)**STD**

FRONT ANGLE ① Measure from frame front to floor.





TFA1 ☐ 70° ☐ 75° ☐ 80° ☐ 85° (Not available with 6" Front Wheels.) ☐ 90° (Not available with 5" or 6" Front Wheels.).....**STD**

SEAT TO FOOTEST ⑥/FOOTESTS Measure from front edge of seat sling to top rear of footrest.

Standard Footrests Measurement must be at least 2" shorter than Front Seat Height ⑥.

☐ 13.5" (34.3cm) ☐ 14" (35.6cm) ☐ 14.5" (36.8cm) ☐ 15" (38.1cm) ☐ 15.5" (39.4cm) ☐ 16" (40.6cm)
☐ 16.5" (41.9cm) ☐ 17" (43.2cm) ☐ 17.5" (44.5cm) ☐ 18" (45.7cm) ☐ 18.5" (47cm) ☐ 19" (48.3cm)

TFTR3 ☐ Angle Adjustable Footrest.....**STD**
TFTR1 ☐ Titanium Footrest w/Flat ABS Cover.....**\$200**
TFTR2 ☐ Titanium Open Loop Footrest.....**\$255**
TFTR12 ☐ Flip-Back Footrest (Not available with ABS Cover with 2" Sides or Front Roller.).....**\$235**
High Mount Footrest
☐ 7"-9" ☐ 9"-11" ☐ 11"-13" Set at _____(Use 1/2" increments)
TFTR4 ☐ Angle Adjustable High Mount Footrest.....**\$145**



TiLite T-Series Canadian Order Form

Effective August 1, 2020 (New Items in Red)

DUAL-TUBE RIGID FRAME

Page 3 of 8

TiFit: TAILORED INDIVIDUALIZED FRAME

FOOTREST WIDTH ❶ or **❷** Select either a Standard Front End ❶ if no Seat Taper ❶ is selected or a V-Front End ❷ if a Seat Taper ❶ is selected.

STANDARD FRONT END ❶ Find your Rear Seat Width **A** in the Chart below and select one of the options to the right. Measure from inside of front frame tube to inside of opposite front frame tube.



TFE1 ☐ NO V-Front Taper.....STD

REAR SEAT WIDTH A (Outside Measurement)	FOOTREST WIDTH ❶		
	No Taper (2" Narrower than Rear Seat Width)	2" Taper (4" Narrower than Rear Seat Width)	4" Taper (6" Narrower than Rear Seat Width)
12"	☐ 10"	☐ 8"	
13"	☐ 11"	☐ 9"	
14"	☐ 12"	☐ 10"	☐ 8"
15"	☐ 13"	☐ 11"	☐ 9"
16"	☐ 14"	☐ 12"	☐ 10"
17"	☐ 15"	☐ 13"	☐ 11"
18"	☐ 16"	☐ 14"	☐ 12"
19"	☐ 17"	☐ 15"	☐ 13"
20"	☐ 18"	☐ 16"	☐ 14"

V-FRONT END ❷ Measure from inside of front frame tube to inside of opposite front frame tube 2-1/2" above footrest. (Titanium Footrests have 1" of adjustment up and down. Not available with Ergonomic Seat frames.)

TFE2 ☐ V-Front Taper.....NCO

- ☐ 9" (22.9cm) ☐ 10" (25.54cm) ☐ 11" (27.9cm)
☐ 12" (30.5cm) ☐ 13" (33cm) ☐ 14" (35.6cm)
☐ 15" (38.1cm)



BACK SUPPORTS

SSH3 Back Support Installation ☐ Install Back OR ☐ Do Not Install Back.....NCO

ROHO® AGILITY® ACTIVE ALUMINUM BACK SUPPORT

Shell Height ☐ 8" ☐ 10"
 BRS14 ☐ Min Contour AGILITY Active with Foam.....\$895
 BRS15 ☐ Mid Contour AGILITY Active with Foam.....\$895

ROHO® AGILITY® ACTIVE CARBON BACK SUPPORT

Shell Height ☐ 10" ☐ 13" ☐ 16"
☐ Min Contour ROHO AGILITY Carbon Back with Foam w/ Air Insert
 BRS7 Foam w/ Air Insert.....\$1095
 ACB1 AGILITY™ Active Carbon Upgrade (Not available with Bolt-on Push Handles).....\$230
 Fixed Hardware.....NCO
☐ Mid Contour ROHO AGILITY Carbon Back with Foam w/ Air Insert
 BRS9 Foam w/ Air Insert.....\$1095
 ACB2 AGILITY Active Carbon Upgrade (Not available with Bolt-on Push Handles).....\$230
 Fixed Hardware.....NCO

COMFORT ACTA-EMBRACE™

Shell Height ☐ 10" ☐ 12" ☐ 14" ☐ 16"
 SBRS16 ☐ Acta-Embrace Back with 2-PT Quick Release Compass Hardware.....\$775
 CPSB1 ☐ Comfort Privacy Shield.....\$55

COMFORT ACTA-RELIEF

Shell Height ☐ 10" ☐ 12" ☐ 14" ☐ 16"
 SBRS13 ☐ Acta-Relief Back with 2-PT Quick Release Compass Hardware.....\$945
 CPSB1 ☐ Comfort Privacy Shield.....\$55

BACK UPHOLSTERY (Want more seating options? [Click here](#) to view the new TiLite Seating Order Form.)

TBUP2 ☐ Velcro Adjustable (Not available with Fold-Down Push Handles.).....STD

Tension Adjustable by Straps (Required with Fold-Down Push Handles.)

TBUP3 ☐ Padded Nylon.....\$415
 TBUP5 ☐ AIR Breathable Upholstery (Not available with Ultrasuede Backrest Stripe.).....\$400
 TUPC3 ☐ Ultrasuede Backrest Stripe ☐ Black.....\$105
 TBUP0 ☐ Omit Back Upholstery.....NCO

TiLite T-Series Canadian Order Form

Effective August 1, 2020 (New Items in Red)

DUAL-TUBE RIGID FRAME

Page 4 of 8

TiFit: TAILORED INDIVIDUALIZED FRAME

SEAT BACK TYPE/PUSH HANDLES/SEAT BACK HEIGHT  Measure from top of back post to top of seat tube at rear of frame.

- TSB1 ☐ Folding Aluminum Adjustable Height w/Precision Lock Back.....**STD**
- TSB2 ☐ Folding Titanium Adjustable Height w/Precision Lock Back and Rounded Rigidizer Bar (Titanium back will match Frame Finish. Not available with Bolt-On Push Handles.).....**\$295**
- TSB15 ☐ Non-Folding, Fixed Height Titanium with New Rounded Rigidizer Bar (Only available with TR. Titanium back will match Frame Finish. Not available with Push Handles.).....**\$325**

TISHAFT™ BACK RELEASE BARS (Not available with Non-Folding, Fixed Height Titanium Back.)

- TSB3 ☐ Aluminum.....**STD**
- TSB5 ☐ Carbon Fiber.....**\$85**

PUSH HANDLES

- TSB6 ☐ Integrated Push Handles.....**NCO**
- TSB7 ☐ Integrated Fold-Down Push Handles (Requires Tension Adjustable by Straps Back Upholstery.).....**\$115**
- TSB8 ☐ Bolt-On Push Handles (Includes 4" Deep Backrest Rigidizer Bar. Not available with Ti Back.).....**\$200**
- TSB9 ☐ Bolt-On Height-Adjustable Push Handles (Includes 4" Deep Backrest Rigidizer Bar. Not available with Ti Back.).....**\$410**

FOLDING SEAT BACK HEIGHT (TSB1 and TSB2)

Push Handles (Selected Above)	Short	Medium	Tall	Set at: (Use 1/2" increments)
None/ Bolt-On	☐ 9.5"-12.5"	☐ 13"-17.5"	☐ 17"-21.5"	Set at: _____
Integrated	☐ 11.5"-13.5"	☐ 14.5"-17.5"	☐ 18.5"-21.5"	Set at: _____
Integrated Fold-Down	☐ 14"-16.5"	☐ 15"-19.5"	☐ 19.5"-22"	Set at: _____

NON-FOLDING, FIXED-HEIGHT  TITANIUM SEAT BACK with New Rounded Rigidizer Bar (TSB15) (Titanium back will match Frame Finish. Required for fixed Side Guards. Not available with Swing Away Tubular or Swing Away - Flip Back Tubular Armrests and Integrated Push Handles, Integrated Fold-Down Push Handles or Bolt-On Height Adjustable Push Handles.)

- ☐ 5" (12.7cm) ☐ 5.5" (14cm) ☐ 6" (15.2cm) ☐ 6.5" (16.5cm) ☐ 7" (17.8cm) ☐ 7.5" (19.1cm) ☐ 8" (20.3cm)
- ☐ 8.5" (21.6cm) ☐ 9" (22.9cm) ☐ 9.5" (24.1cm) ☐ 10" (25.4cm) ☐ 10.5" (26.7cm) ☐ 11" (27.9cm) ☐ 11.5" (29.2cm)
- ☐ 12" (30.5cm) ☐ 12.5" (31.8cm) ☐ 13" (33cm) ☐ 13.5" (34.3cm) ☐ 14" (35.6cm) ☐ 14.5" (36.8cm) ☐ 15" (38.1cm)
- ☐ 15.5" (39.4cm) ☐ 16" (40.6cm) ☐ 16.5" (41.9cm) ☐ 17" (43.2cm) ☐ 17.5" (44.5cm) ☐ 18" (45.7cm) ☐ 18.5" (47cm)
- ☐ 19" (48.3cm) ☐ 19.5" (49.5cm) ☐ 20" (50.8cm)

AVAILABLE ON ALL BACKS

- TSB11 ☐ Backrest that Locks when Folded (Not available with Ergonomic Seat.).....**STD**
- TSB12 ☐ 4" Deep Backrest Rigidizer Bar (Not available with a Non-Folding, Fixed Height Titanium Seat Back Type. 4" Deep Rigidizer Bar is NOT a rounded style rigidizer bar.).....**NCO**

SEAT BACK ANGLE  Measure angle from front of back post to floor.




- TSBA1 **Folding Back**
- ☐ 80° ☐ 83° ☐ 86° ☐ 89° ☐ 92° ☐ 95° ☐ 98° ☐ 101°.....**STD**
- TSBA1 **Non-Folding Back**
- ☐ 80° ☐ 81° ☐ 82° ☐ 83° ☐ 84° ☐ 85° ☐ 86° ☐ 87° ☐ 88° ☐ 89° ☐ 90°
- ☐ 91° ☐ 92° ☐ 93° ☐ 94° ☐ 95° ☐ 96° ☐ 97° ☐ 98° ☐ 99° ☐ 100° ☐ 101°.....**NCO**

CENTER OF GRAVITY  Measure from front of back post at the seat tube to the center of the rear axle. If no center of gravity is selected, chair will ship with the center of gravity in the rearward-most position.

- TCOG1 ☐ 1" (2.5cm) ☐ 1.25" (3.2cm) ☐ 1.5" (3.8cm) ☐ 1.75" (4.4cm) ☐ 2" (5.1cm)
- ☐ 2.25" (5.7cm) ☐ 2.5" (6.4cm) ☐ 2.75" (7cm) ☐ 3" (7.6cm) ☐ 3.25" (8.3cm)
- ☐ 3.5" (8.9cm) ☐ 3.75" (9.5cm) ☐ 4" (10.2cm) ☐ 4.25" (10.8cm) ☐ 4.5" (11.4cm)
- ☐ 4.75" (12.1cm) ☐ 5" (12.7cm).....**STD**

REAR WHEEL SPACING /PLUGS Measure from outside of seat back post to inside of rear tire. 0.5" Rear Wheel Spacing on TR not available with 0° or 2° Camber. TRA 0.75" Rear Wheel Spacing not available with 0° or 2° Camber

- K** Ergonomic Seat selection may limit available center of gravity.
- TRWS1 ☐ .75" ☐ 1" ☐ 1.25".....**STD**
- TCBR12 ☐ Adjustable Rear Wheel Spacing Camber Plugs (Not available with Power Adaptable Frame.).....**STD**
- TCBR9 ☐ Additional Spacers to Increase Rear Wheel Spacing.....**NCO**
- TCBR11 ☐ Zero Play Aluminum Camber Plugs (Not available with Power Adaptable Frame or more than 6° of Camber **K**.).....**NCO**



TiLite T-Series Canadian Order Form

Effective August 1, 2020 (New Items in Red)

DUAL-TUBE RIGID FRAME

Page 5 of 8

TiFit: TAILORED INDIVIDUALIZED FRAME

SEAT UPHOLSTERY

TSUP1	☐ Tension Adjustable Bolt On (Not available with Ergonomic Seat or Tapered Seat ①. Black only.)	STD
TSUP2	☐ Tension Adjustable by Straps (Standard with Ergonomic Seat and FSW (Seat Taper) ① of 1", 2" or 3". Black upholstery only.)	\$185
TSUP8	☐ Aluminum Solid Seat (Not available with 10" or 11" Rear Seat Width ④ or Ergonomic Seat ①.)	\$285
TSUP0	☐ Omit Seat Upholstery	NCO

SEAT CUSHION (Want more seating options? [Click here](#) to view the new TiLite Seating Order Form.)

TCSH2	☐ Comfort Curve Cushion	\$155
TCSH23	☐ Comfort Acta-Embrace™ Comfort-Tek Zero Elevation Cushion	\$470
TCSH11	☐ ROHO Single Compartment (Fits seat widths 14" to 19") Profile: ☐ Low ☐ Mid ☐ High	\$720
TCSH14	☐ ROHO Quadtro Select (Fits seat widths 14" to 19") Profile: ☐ Low ☐ Mid ☐ High	\$720

CAMBER ⑤/CAMBER TUBES

TCBR1	☐ 0° ☐ 2° ☐ 4° ☐ 6° (8° of Camber ⑤ may be available through our Creative Mobility Division.)	STD
TCBR15	☐ Aluminum Camber Tube	STD
TCBR6	☐ NEW Carbon Fiber Camber Tube (Not available with 6° of Camber ⑤ and seat widths greater than 18")	\$310
TCBR3	☐ Titanium Camber Tube (Not available with Power Adaptable Frame.)	\$225

CUSTOM CONFIGURATION OPTIONS

FRONT WHEELS

TFW24	☐ 3"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire ☐ Silver Hub OR ☐ Black Hub	\$180
TFW14	☐ 4"x.75" Plastic Wheel w/Poly Tire (Not available with Slipstream Fork, Glide™ Single-Sided Fork or Power Adaptable Frame.)	NCO
TFW3	☐ 4"x.75" Light-Up Micro-Caster (Not available with Slipstream Fork, Glide™ Single-Sided Fork or Power Adaptable Frame.)	\$80
TFW26	☐ 4"x1.5" LiteSpeed Plastic Wheel w/Soft Roll Tire (Not available with Power Assist, Power Add-on System or Heavy Duty Frame.)	\$95
TFW27	☐ 4"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire ☐ Silver Hub OR ☐ Black Hub	\$190
TFW6	☐ 5"x1" Plastic Wheel w/Poly Tire (Not available with 90° Front Angle ①, Slipstream Fork, Glide™ Single-Sided Fork or Power Adaptable Frame.)	NCO
TFW28	☐ 5"x1.5" LiteSpeed Plastic Wheel w/Soft Roll Tire (Not available with 90° Front Angle ① or Power Adaptable Frame.)	\$95
TFW30	☐ 5"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire ☐ Silver Hub OR ☐ Black Hub (Not available with 90° Front Angle ①)	\$250
TFW9	☐ 6"x1" Plastic Wheel w/Poly Tire (Not available with 85° or 90° Front Angle ①, Slipstream Fork, Glide™ Single-Sided Fork or Power Adaptable Frame.)	NCO
TFW32	☐ 6"x1.5" LiteSpeed Plastic Wheel w/Soft Roll Tire (Not available with 85° or 90° Front Angle ①, Slipstream Fork, Glide™ Single-Sided Fork or Power Adaptable Frame.)	\$100
TFW33	☐ 6"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire ☐ Silver Hub OR ☐ Black Hub (Not available with 85° or 90° Front Angle ①, Slipstream Fork, Glide™ Single-Sided Fork.)	\$250

PRECISION LOCK FRONT FORKS

TFK1	☐ TiLite Standard Fork	STD
TFK2	☐ TiLite Slipstream Single-Sided Fork (Not available with Power Adaptable Frame.)	\$185
TGL1	☐ Out-Front Glide™ Suspension Fork (Black anodize only.)	\$545
	Rider's Weight: ☐ 125 lbs. ☐ 125-180 lbs. ☐ 180-265 lbs.	
TTFS1	☐ Titanium Fork Stem	\$60

STANDARD REAR WHEELS (See Rear Seat Height Ranges Chart for Rear Seat Height ⑥ limitations.)

TRW14	TiLite SHADOW™	STD
	☐ 20"(451)	
	☐ 22"(501)	
	☐ 24"(540)	
	☐ 25"(559) (Requires 15" minimum Seat Depth ⑥.)	
	☐ 26"(590) (Requires 15" minimum Seat Depth ⑥.)	

PERFORMANCE REAR WHEELS (See Rear Seat Height Ranges Chart for Rear Seat Height ⑥ limitations.)

TRW7	Spinergy SPOX ☐ Black ☐ Red ☐ Blue ☐ White	\$900
	☐ 22"(501) - Std. Hub	
	☐ 24"(540) - Std. Hub	
	☐ 25"(559) - Std. Hub (Requires 15" minimum Seat Depth ⑥.)	
	☐ 26"(590) - Sport Hub (Requires 15" minimum Seat Depth ⑥. Not available with Quad Release Axles.)	
TRW8	Spinergy LX ☐ Black ☐ Red ☐ Blue ☐ White ☐ Bright Orange ☐ Flamingo Pink	\$1,055
	☐ 22"(501) - Std. Hub	
	☐ 24"(540) - Std. Hub	
	☐ 25"(559) - Std. Hub (Requires 15" minimum Seat Depth ⑥.)	
	☐ 26"(590) - Sport Hub (Requires 15" minimum Seat Depth ⑥. Not available with Quad Release Axles.)	

TiLite T-Series Canadian Order Form

Effective August 1, 2020 (New Items in Red)

DUAL-TUBE RIGID FRAME

Page 6 of 8

CUSTOM CONFIGURATION OPTIONS

TRW16	Spinergy Carbon Blade (Not available with Heavy Duty Frame, Natural-Fit or Natural-Fit LT with Standard Grip or Super Grip.).....	\$1,540
	☐ 24"(540)	
	☐ 25"(559) (Requires 15" minimum Seat Depth ③.)	
TRW17	X-Core Mag (Not available with Heavy Duty Frame, Short Tabs, Quad Release Axles, Spoke Guards or Omit Handrims Option. Requires Aluminum Silver Anodized Handrims.).....	\$110
	☐ 20"(451)	
	☐ 22"(501)	
	☐ 24"(540)	

AXLES

TAXL1	☐ Stainless Quick Release.....	STD
TAXL2	☐ Titanium Quick Release.....	\$290
TAXL3	☐ Quad Release (Not available with Mag or 26"(590) Spinergy Spox or LX Rear Wheels.).....	\$170

PNEUMATIC REAR TIRES

TRTR1	☐ Treaded (Not available with 25"(559) or 20" (451) Rear Wheels.).....	STD
TRTR2	☐ Airless Inserts (Not available with Spinergy Carbon Blade Rear Wheels.).....	\$90
TRTR3	☐ 25"(559) Primo Treaded (Standard with 25"(559) Rear Wheels.).....	STD
TRTR4	☐ Primo Slick.....	NCO
TRTR30	☐ 24"(540), 25"(559) and 26"(590) Primo Sentinel.....	\$160
TRTR8	☐ 24"(540), 25"(559) and 26"(590) Schwalbe® Marathon Plus Evolution.....	\$155
TRTR9	☐ 24"(540) and 25"(559) Schwalbe® Slick.....	\$150
TRTR24	☐ 24"(540) Kenda Kobra Knobby.....	\$95

SOLID REAR TIRES (Not available with Spinergy Carbon Blade Rear Wheels.)

TRTR22	☐ 20"(451), 22"(501), 24"(540) and 25"(559) Primo Express.....	\$110
TRTR23	☐ 24"(540) and 25"(559) Primo Active.....	\$135

HANDRIMS

Handrims by OUT-FRONT®

TRIM1	☐ Aluminum – Silver Anodized.....	STD
TRIM21	☐ Aluminum – Black Anodized (Not available with 20"(451) Rear Wheels.).....	\$60
TRIM12	☐ Titanium – Satin Finish (Not available with 20"(451) or 22"(501) Rear Wheels, or 24"(540) Knobby tires.).....	\$490
TRIM10	☐ The Surge (All-In-One Oval with Gription™ Strip, Larger Profile) (Not available with 20"(451) Rear Wheels.).....	\$435
TRIM11	☐ The Surge LT (All-In-One Oval with Gription™ Strip, Smaller Profile) (Not available with 20"(451) Rear Wheels.).....	\$470
TRIM20	☐ Q-Grip (Exclusive High-Friction, Anti-Slip Coating) (Not available with 20"(451) Rear Wheels.).....	\$375

The Natural-Fit® (Original, Larger Profile) (Not available with 20"(451) Rear Wheels.)

TRIM14	☐ Standard Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.).....	\$415
TRIM15	☐ Super Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.).....	\$470
TRIM16	☐ Without Thumb Grip.....	\$380

The Natural-Fit® LT (Reduced Weight, Smaller Profile) (Not available with 20"(451) Rear Wheels.)

TRIM17	☐ Standard Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.).....	\$500
TRIM18	☐ Super Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.).....	\$535
TRIM19	☐ Without Thumb Grip.....	\$460

Additional Handrim Options

TRIM5	☐ Plastic Coated – Black.....	\$160
TRIM22	☐ Omit Handrims.....	NCO

Tab Length

TRIM7	☐ Long Tabs.....	STD
TRIM8	☐ Short Tabs (Not available with Mag Rear Wheels or The Natural-Fit® Handrims with Thumb Grip.).....	\$35

WHEEL LOCKS

Wheel Locks by OUT-FRONT®

TWLK7	☐ Compact Composite Push to Lock (Available with .75" to 1.5" Rear Wheel Spacing ① only.).....	NCO
TWLK8	☐ Compact Composite Pull to Lock (Available with .75" to 1.5" Rear Wheel Spacing ① only.).....	NCO
TWLK10	☐ Extension Handles for Compact Composite Push/Pull to Lock.....	\$85
TWLK9	☐ Composite Scissor Lock.....	\$175

Additional Wheel Locks

TWLK3	Uni-Lock Set At: ☐ Push OR ☐ Pull.....	\$165
TWLK4	☐ Extension Handles for Uni-Lock.....	\$85
TWLK15	☐ Grade Aids Push to Lock (Not available with solid rear tires).....	\$230

TiLite T-Series Canadian Order Form

Effective August 1, 2020 (New Items in Red)

DUAL-TUBE RIGID FRAME

Page 7 of 8

CUSTOM CONFIGURATION OPTIONS

SIDE GUARDS

Aluminum

TSDG2	☐ Rigid Removable	☐ Small OR ☐ Std	\$290
TSDG3	☐ Quad Removable w/Hole	☐ Small OR ☐ Std	\$315
TSDG4	☐ Fold Down Rigid	☐ Small OR ☐ Std (Not available with Ergonomic Seat.)	\$375
TSDG6	☐ Fendered Adjustable (Available with 24" and 25" Rear Wheels. Requires Rear Seat Height ⑩ of 14.5" to 18.5").		\$400

Carbon Fiber

TSDG7	☐ Rigid Removable	☐ Small OR ☐ Std	\$340
TSDG8	☐ Quad Removable w/Hole	☐ Small OR ☐ Std	\$375
TSDG9	☐ Fold Down Rigid	☐ Small OR ☐ Std (Not available with Ergonomic Seat or Carbon Edition.)	\$510
TSDG11	☐ Fendered Adjustable (Available with 24" and 25" Rear Wheels. Requires Rear Seat Height ⑩ of 14.5" to 18.5").		\$565
TSDG12	☐ NEW Fendered Removable		\$595

Available on Aluminum and Carbon Fiber Non-Fendered Side Guards

TSDG13	☐ Side Guards Custom Cut to Profile of Rear Wheels	\$315
--------	---	-------

ARMRESTS (Desk Arms come with built-in Side Guards.)

TARM1	Swing Away Tubular ☐ 6.5" to 9.5" OR ☐ 9.5" to 12.5" (Not Available with Non-Folding, Fixed Height Titanium Seat Back or Height Adjustable Push Handles.)	\$260
TARM4	☐ Swing Away - Flip Back Tubular 8.5" to 11.5" (Not Available with Non-Folding, Fixed Height Titanium Seat Back.)	\$315
TARM2	Removable Height Adjustable Desk Arm with Rigid Side Guard (Requires 14" minimum Seat Depth ⑩. Not available with Ergonomic Seat.)	\$360
	☐ 9" to 12" OR ☐ 12" to 15"	
	☐ Standard Length Pad (10") OR ☐ Full Length Pad (14")	
TARM3	Removable Height Adjustable Desk Arm with Rigid Side Guard and Transfer Loop (Requires 15" minimum Seat Depth ⑩. Not Available with Ergonomic Seat.)	\$375
	☐ 9" to 12" OR ☐ 12" to 15"	
	☐ Standard Length Pad (10") OR ☐ Full Length Pad (14")	

Available Armrest Option

TARM5	☐ Covered in Black Ultrasuede (Not available with Desk Arms)	\$115
-------	---	-------

SEAT BELTS

TBLT1	☐ Velcro Adjustable, Black Nylon Webbing	\$55
TBLT2	☐ Bodypoint® 2" Push-Button Buckle (Auto-Style.)	\$95
TBLT4	☐ Bodypoint® 1" Wide Push-Button Buckle Padded (Auto-Style.) (Not available with 19" or 20" Rear Seat Widths ①.)	\$170
TBLT5	☐ Bodypoint® 1.5" Wide Push-Button Buckle Padded (Auto-Style.) (Not available with 19" or 20" Rear Seat Widths ①.)	\$170
TBLT10	☐ Bodypoint® Optional Band Clamp Mounting Hardware	\$20
TBLT16	☐ Bodypoint® EvoFlex® Pelvic Stabilizer Push Button	\$185

CALF STRAP

TCLF2	Velcro Adjustable	\$45
	☐ Xtra Small (Recommended for 6" - 9" Footrest Widths.)	
	☐ Small (Recommended for 9" - 12" Footrest Widths.)	
	☐ Medium (Recommended for 11" - 15" Footrest Widths.)	
	☐ Large (Recommended for 14" - 17" Footrest Widths.)	
	☐ Xtra Large (Recommended for 16" - 19" Footrest Widths.)	
TCLF3	Bodypoint® Padded Velcro Adjustable	\$80
	☐ Medium (Recommended for 11" - 15" Footrest Widths.)	
	☐ Large (Recommended for 16" - 19" Footrest Widths.)	
TCLF1	☐ Buckle Adjustable, Black Nylon Webbing	\$30

RESIDUAL LIMB SUPPORT

TRSL1	☐ Residual Limb Support ☐ Left OR ☐ Right (8" width x 10" length)	\$515
TRSL2	☐ Residual Limb Support Pair (8" width x 10" length)	\$1,030

BACKPACK

TBCK3	☐ Black Permobil Backpack	\$95
-------	--	------

SEAT POUCH

TPCH3	☐ Black Permobil Seat Pouch	\$60
TPCH4	☐ Bodypoint® Mobility Bag	\$130

SPOKE GUARDS

TGRD1	☐ Clear Lexan with TiMotion Logo (Not available with Mag Rear Wheels.)	\$135
-------	---	-------

TiLite T-Series Canadian Order Form

Effective August 1, 2020 (New Items in Red)

DUAL-TUBE RIGID FRAME

Page 8 of 8

CUSTOM CONFIGURATION OPTIONS

ANTI-TIPS


TTIP1	☐ Rear Aluminum - Black.....	\$195
TTIP7	☐ Black Aluminum Single.....	\$100
	Mount on ☐ Left Side or ☐ Right Side.....	NCO
TTIP2	☐ User Friendly Flip-up (Not available with Frog Legs® Rear Suspension.).....	\$360
ZTIP8	☐ User Friendly Flip-up Single.....	\$170
	Mount on ☐ Left Side or ☐ Right Side.....	NCO

ACCESSORIES


TUEST1	Upper Extremity Support Tray ☐ 12" to 14" OR ☐ 13" to 15" Seat Widths [Ⓐ] (Requires Armrests.).....	\$290
TLUG1	☐ Luggage Carriers.....	\$145
TMPCT2	Neoprene Impact Guards with Ultrasuede Stripe ☐ Black.....	\$105
TCH1	☐ Crutch Holder.....	\$145

FRAME NOTES

The rear frame for a TR will be square when the Center of Gravity [Ⓐ] measurement is less than 3.5". Aero T and TRA will always be a Rounded Rear Frame.



Ⓐ Round Rear Frame



Ⓑ Square Rear Frame

GENERAL

More options and frame modifications may be available through our Creative Mobility Division! Please contact customer service for details or [click here](#) to view options.

Permobil manufactures TiLite manual wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. HCPCS Codes are for reference purposes only. Actual coding and coverage varies by payer. Proper HCPCS coding is not a guarantee of payment. Please verify the appropriate payer's coverage policy.

Get the latest order forms at permobilus.com. Prices, standard features, options and accessories are subject to change without notice. All rights reserved. Trademarks are identified by symbols™ and ®. TiLite, TiFit, the TiLite logo, Ti, the Ti logo, and TiShaft are registered trademarks of TiSport, LLC and are licensed to Permobil. Permobil® is a registered trademark of Permobil. Adobe®, Acrobat®, Acrobat Professional® and Reader® are registered trademarks of Adobe Systems Inc. Bodypoint® is a registered trademark of Bodypoint Inc. Kenda® is a registered trademark of Kenda Rubber Industrial Co. Ltd. ROHO®, AGILITY®, Low Profile® and Quadtro Select® are registered trademarks of ROHO Inc. Schwalbe® is a registered trademark of Ralf Bohle GmbH. The Natural-Fit® and Out-Front® are registered trademarks of Three Rivers Holdings, LLC. The ACTA-BACK®, ACTA-EMBRACE™, COMFORT-TEK™, COMPASS®, CURVE®, ELEMENTS®, and STRETCH AIR™ are registered trademarks of THE COMFORT COMPANY® and are licensed to Permobil. Q'Straint® is a registered trademark of Q'Straint Inc.

NOTES

QUOTE THIS CHAIR

If you are using Adobe Acrobat, this order form may be converted into a quote request or into an order by clicking one of the buttons (this feature is available exclusively to Permobil dealers and distributors). A Quote or an Order Confirmation will be sent to you within 24 hours.

ORDER THIS CHAIR