

TiLite 2GX Order Form

Effective Aug 2025 (New Items in Red)

FOLDING FRAME

Page 1 of 7

DEALER INFORMATION		RTS INFORMATION		
Dealer:	Account #:	RTS:	Phone:	
PO:	Ordered By:	Email:		
Email:		SHIP TO INFORMATION		
Mark For:	User Weight:	Name:	Ship Via:	
TiLite Sales Rep:		Address:		
TiLite Information:	Promo Code:	City:	State:	Zip:
Ship Date:	Order #:	Phone:	Fax:	
More options and frame modifications may be available through our Creative Mobility Division! Please contact customer service for details or click here to view options.				
Please Email or Fax to the Permobil Manual Customer Experience Team at tilitesales@permobil.com or 509.586.2413				



FRAME STYLE

- G2FS1 ☐ **2GX Titanium Frame** (Full Adjustability. 265 Pound Weight Limit. Light Weight and High Strength. Not Available with Frame Transit Tie-Down.).....\$4,692 K0005
- XFS1 ☐ **Swing Away**.....STD
- XFS3 ☐ **Hemi Height Frame** (Available with a Front Seat Height of 13.5" - 20" (34-51cm). Not Available with Depth Adjustable Back/Frame, 12" or 13" Seat Depth or Heavy Duty Frame.).....NCO
- XFS2 ☐ **Fixed Front**.....NCO K0005
- XHDR1 ☐ **Heavy Duty Frame** (350 Pound Weight Limit. Includes Tension Adjustable by Straps Back Upholstery. Available with 18" to 22" Seat Widths and 16" to 20" Seat Depths. Not Available with 90° front end).....NCO
- XPWR1 ☐ **Power Adaptable Frame** (Requires LiteSpeed Billet Aluminum Front Wheels and 0 degree Camber. Not available with Single-Sided Forks).....NCO
- XRFR1 ☐ **Reinforced Frame**.....NCO

SMARTDRIVE POWER ASSIST SYSTEM



- XMAXSD1 ☐ **Permobil SmartDrive MX2+*** (Not available with rear seat heights less than 13.5").....\$7,529 E0986
- XMAXSD18 ☐ **SwitchControl Buttons**.....NCO
- XMAXSD8 ☐ **SwitchControl Button with Mono Jack Option**.....NCO
- XMAXSD19 ☐ **NEW Buddy Button**.....NCO
- XMAXSD4 ☐ **Carrying Bag**.....\$84
- XMAXSD10 ☐ **Universal Padded Handle**.....\$22

*SmartDrive apps available on Apple App Store and Google Play Store for Apple Watch Series 5 or newer and Samsung Galaxy Watch 4 or newer.

TRANSIT OPTIONS

- XWTB ☐ **Transport Bracket Option** (Consists of four Tie-Down Brackets. Unoccupied Use Only.) (Not Available when Swing Away frame and Heavy Duty Frame are both selected. Not available with Swing Away Tubular Armrest or Swing Away -Flip Back Tubular Armrests. Requires a minimum seat depth of 14" or 15" seat Depth with Hemi Height Frame.).....\$379 E1022

FRAME FINISHES

Natural Finishes

- XFF1 ☐ **Satin** (Available with Titanium Frame only.).....STD
- XSCBT1 ☐ **Satin Finish Buffing Pads**.....\$22
- XFF3 ☐ **Tattoos** (Tattoos on front frame tubes. Titanium frame finished in Satin. Aluminum frame tattooed in blasted aluminum with clear coat. Tattoos will be trimmed and adjusted to fit your frame.).....\$363
- ☐ **XFF45 Need for Speed**
- ☐ **XFF44 Flower Garden**



XFF43 Mosaic



- | | | |
|--|---|--|
| ☐ XFF8 High Gloss Black | ☐ XFF10 Red Pepper | ☐ XFF29 Neon Pink |
| ☐ XFF14 Ocean Blue Metallic | ☐ XFF15 Tangerine Metallic | ☐ XFF30 Acid Green Metallic |
| ☐ XFF17 Midnight Blue Pearl | ☐ XFF18 White Pearl | ☐ XFF19 Black Lava |
| ☐ XFF27 Titanium Grey | ☐ XFF59 Pomegranate Red | ☐ XFF62 Blue Mist |
| ☐ XFF57 Indigo Purple | ☐ XFF40 Hyper Blue | ☐ XFF41 Flame Orange |
| ☐ XFF42 OD Green | ☐ XFF60 Silver Streak | ☐ XFF63 Lavender Dusk |
| ☐ XFF58 Black Galaxy | ☐ XFF61 Seafoam Green | ☐ XFF64 Daffodil Yellow.....STD |



Permobil.com
2701 West Court Street / Pasco, Washington 99301
Tel 800 736 0925 or Email tilitesales@permobil.com

OF0003 2GX 07132021RevF

TiLite 2GX Order Form

Effective Aug 2025 (New Items in Red)

FOLDING FRAME

Page 2 of 7

FRAME MEASUREMENTS

Measurements shown below are our most popular selections. **Not seeing what you want?** Contact customer service or visit our Creative Mobility Division (CMD) catalog [here](#) to view additional options.



SEAT WIDTH A Measure from outside of seat tube at back post to the outside of the opposite seat tube at back post.

XSW1 ☐ 12" ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19".....**STD**
XSW2 ☐ 20" ☐ 21" ☐ 22" (21" and 22" Seat Widths A are only available with Heavy Duty Frame.).....**\$536 E2201**

SEAT DEPTH B Measure from 1/2" in front of back post to the front edge of seat sling. (12"-14" Seat Depths B not available with 25" (559) and 26"

(590) Rear Wheels)

XSD1 ☐ 12" ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19".....**STD**
XSD2 ☐ 20".....**\$568 E2203**
XSD3 ☐ CustomFrame Depth: _____" Frame Depth Longer Than Seat Depth (Use 1" increments with a maximum of 2". A Custom Frame Depth with a maximum of 20" is \$546.).....**NCO**



FRONT SEAT HEIGHT C Measure from floor to top of seat tube at front edge of seat sling. Adjustable in 1/2" increments. (13.5" - 14" Front Seat Height C requires Hemi Height Frame Style. 20.5" - 21.5" Front Seat Height C not available with Hemi Height Frame Style.).

FRONT SEAT HEIGHT RANGES

	Front Wheel Size	Front Seat Height C	Hemi Front Seat Height C
	3" Front Wheel	14.5" – 18.5"	13.5" – 17.5"
	3" Front Wheel w/ Glide™	16.5" – 17.5"	15.0" – 17.0"
	4" Front Wheel	15.5" – 19.0"	14.5" – 18.0"
	4" Front Wheel w/ Glide™	17.0" – 18.0"	15.5" – 17.0"
	5" Front Wheel	16.5" – 19.5"	16.0" – 18.5"
	5" Front Wheel w/ Glide™	17.5" – 18.5"	16.5" – 18.0"
	6" Front Wheel	17.5" – 20.0"	16.5" – 19.0"
	6" Front Wheel w/ Glide™	19.0"	17.5" – 18.0"
	8" Front Wheel	19.5" – 21.0"	18.5" – 19.5"
XFSH1	☐ 13.5" ☐ 14" ☐ 14.5" ☐ 15" ☐ 15.5" ☐ 16" ☐ 16.5" ☐ 17.0" ☐ 16.5" ☐ 16.5" ☐ 17.0"	☐ 17.5" ☐ 18.0" ☐ 18.5" ☐ 19.0" ☐ 19.5" ☐ 20.0" ☐ 20.5" ☐ 21.0"	 STD

REAR SEAT HEIGHT D Measure from floor to top of seat tube at back post. Adjustable in 3/8" increments.

Rear Wheel Size	Rear Seat Height
22"	13.0" – 18.0"
24"	14.0" – 19.0"
25"	14.0" – 19.0"
26"	15.0" – 20.0"

XRSH1 ☐ 13" ☐ 13.5" ☐ 14" ☐ 14.5" ☐ 15" ☐ 15.5" ☐ 16" ☐ 16.5" ☐ 17" ☐ 17.5" ☐ 18" ☐ 18.5"
☐ 19" ☐ 19.5" ☐ 20".....**STD**

BACK SUPPORTS

COMFORT ACTA-BACK® Back Support

☐ Acta-Back with 2-PT Quick Release Compass Hardware.....**\$616 E2613**
Shell Height ☐ 10" ☐ 12" ☐ 14" ☐ 16" ☐ 18"**STD**
☐ 20".....**\$83**
☐ Comfort-Tek Cover ☐ Stretch-Air Cover.....**NCO**
☐ Color Accent (Colored piping along the edge of the cover on the backside of the backrest) ☐ Silver.....**STD**
☐ Red ☐ Purple ☐ Pink ☐ Blue ☐ Green ☐ Orange ☐ Remove (no color).....**\$30**
☐ Comfort Privacy Shield.....**\$44**

COMFORT ACTA-BACK® Deep Back Support

☐ Acta-Back Deep with 2-PT Quick Release Compass Hardware.....**\$730 E2615**
Shell Height ☐ 10" ☐ 12" ☐ 14" ☐ 16" ☐ 18"**STD**
☐ 20".....**\$93**

TiLite 2GX Order Form

Effective Aug 2025 (New Items in Red)

FOLDING FRAME

Page 3 of 7

☐ Comfort-Tek Cover	☐ Stretch-Air Cover	NCO					
☐ 3" Deep Wing	☐ 5" Deep Wing	NCO					
☐ Color Accent (Colored piping along the edge of the cover on the backside of the backrest)	☐ Silver	STD					
☐ Red	☐ Purple	☐ Pink	☐ Blue	☐ Green	☐ Orange	☐ Remove (no color)	\$30
☐ Comfort Privacy Shield							\$44

COMFORT ACTA-EMBRACE® Back Support

Shell Height ☐ 10" ☐ 12" ☐ 14" ☐ 16"

SBRS16	☐ Acta-Embrace Back with 2-PT Quick Release Compass Hardware	\$642	E2613
CPSB1	☐ Comfort Privacy Shield	\$44	

COMFORT ACTA-RELIEF Back Support

Shell Height ☐ 10" ☐ 12" ☐ 14" ☐ 16"

SBRS13	☐ Acta-Relief Back with 2-PT Quick Release Compass Hardware	\$846	E2615
CPSB1	☐ Comfort Privacy Shield	\$44	

SEAT BACK HEIGHT E/SEAT BACK/PUSH HANDLES Measure from top of back post to top of seat tube at rear of frame.

SEAT BACK TYPE

XSB13	☐ Fixed Angle, Adjustable-Height Seat Back	STD
XLB1	☐ Folding Stabilizer Bar (Required with Transit Tie-Down Option unless Depth Adjustable Backrest is selected. Includes Tension Adjustable by Strap Back Upholstery.)	\$423
XSB14	☐ Fixed Angle, Adjustable Height Seat Back w/8-Degree Bend. (Not available with Heavy Duty Frame.)	NCO
XSB18	☐ Folding Adjustable-Angle and Height Seat Back (Not available with Heavy Duty Frame.)	\$250
XSB16	☐ Depth Adjustable Back/Frame (Not Available with Hemi Height Frame, Heavy Duty Frame, 19" or 20" Seat Depths B, Swing Away Tubular or Swing Away Flip-Back Armrest. Seat Depths B of 12"-17" allow 3" of rearward adjustment. Seat Depth B of 18" allow 2" of rearward adjustment. Requires Tension Adjustable Seat Upholstery)	\$244

PUSH HANDLES

XSB6	☐ Integrated Push Handles	NCO
XSB7	☐ Integrated Fold-Down Push Handles (Requires Tension Adjustable by Straps Back Upholstery. Not available with Heavy Duty Frame.)	\$125
XSB9	☐ Bolt-On Height Adjustable Push Handles (Requires Tension Adjustable by Straps Back Upholstery. Not available with Heavy Duty Frame, Folding Adjustable-Angle and Height Seat Back, Swing Away Tubular Armrests or Swing Away - Flip Back Tubular Armrests.)	\$363

XSB20/XSB14 **FIXED ANGLE ADJUSTABLE SEAT BACK HEIGHT E** (11"- 15" Fixed Seat Back Height not available with Heavy Duty Frame, Swing Away Tubular Arms 12"-15", Bolt-On Height Adjustable Push Handles or Folding Stabilizer Bar.)

Push Handles (Selected on Pg. 2)	Short	Tall	Set at: (Use 1" Increments)
None	☐ 11.0" - 15.0"	☐ 16" - 20"	Set at: _____
Integrated	☐ 11.0" - 15.0"	☐ 16" - 20"	Set at: _____
Integrated Fold-Down	☐ 11.0" - 15.0"	☐ 16" - 20"	Set at: _____
Bolt-On Height Adj.	-----	☐ 16" - 20"	Set at: _____

XSB18 **FOLDING ADJUSTABLE- ANGLE AND HEIGHT SEAT BACK E** (Not available with a Folding Stabilizer Bar or Heavy Duty Frame)

Push Handles (Selected on Pg. 2)	Short	Medium	Tall	Set at: (Use 1" Increments)
None	☐ 11.0" - 13.0"	☐ 13.5" - 16.5"	☐ 16.5" - 19.5"	Set at: _____
Integrated	☐ 12.5" - 14.5"	☐ 13.5" - 16.5"	☐ 17.0" - 19.5"	Set at: _____
Integrated Fold-Down	☐ 15.0" - 17.0"	☐ 16.0" - 19.0"	☐ 18.5" - 20.5"	Set at: _____

FRAME MEASUREMENTS

Options shown below are our most popular selections. **Not seeing what you want?** Contact customer service or visit our Creative Mobility Division (CMD) catalog [here](#) to view additional options.

CAMBER F

XCBR8	☐ 0° ☐ 2° ☐ 4° ☐ 6°	STD
-------	---	-----

FRONT ANGLE H Front Angle H is measured from inside of front frame tube to top of bottom frame tube.

Swing Away Front End

XFA3	☐ 70° ☐ 80° ☐ 85° (Not available with 6" Front Wheels.) ☐ 90° (Requires aluminum footplate. Not available with Heavy Duty Frame, 5" or 6" Front Wheels.)	STD
XFA0	☐ Omit Hangers	NCO

Fixed Front End

XFA6	☐ 70° ☐ 80° ☐ 85° (85° Not available with 6" Front Wheels.)	STD
------	--	-----

SEAT TO FOOTREST I Measure from front edge of seat sling to top rear of footrest.

If your Front Angle H is 70°, 80°, OR 85° select one Seat to Footrest I length and one Footrest type:

Standard Footrests Measurement must be at least 2" shorter than Front Seat Height C. 70° Fixed Front has a minimum seat to footrest of 14.5"

☐ 14" ☐ 14.5" ☐ 15" ☐ 15.5" ☐ 16" ☐ 16.5" ☐ 17" ☐ 17.5" ☐ 18" ☐ 18.5" ☐ 19" ☐ 19.5" ☐ 20"

TiLite 2GX Order Form

Effective Aug 2025 (New Items in Red)

FOLDING FRAME

Page 4 of 7

XFTR5	☐ Flip-up Footrest.....	STD	
XFTR9	☐ Depth Adjustable Flip-up Footrest.....	\$120	K0045
XFTR10	☐ Depth & Angle Adjustable Flip-up Footrest (Available for 15"- 20" Seat to Footrest I only.).....	\$293	K0040, K0045
XFTR11	☐ Angle Adjustable Flip-up Footrest (Available for 15"- 20" Seat to Footrest I only.).....	\$185	K0040, K0045

High Mount Footrests (Not Available with Luggage Carriers, or Impact Guards.)

	☐ 7" ☐ 7.5" ☐ 8" ☐ 8.5" ☐ 9" ☐ 9.5" ☐ 10" ☐ 10.5" ☐ 11" ☐ 11.5" ☐ 12" ☐ 12.5" ☐ 13" ☐ 13.5" ☐ 14" ☐ 14.5"		
XFTR20	☐ High Mount Flip-up.....	\$202	K0037
XFTR15	☐ Angle Adjustable High Mount Flip-up Footrest.....	\$352	K0040, K0037
XFTR16	☐ Depth & Angle Adjustable High Mount Flip-up Footrest.....	\$406	K0040, K0037

If your Front Angle H is 90°, select one Seat to Footrest I length and one Footrest type:

	☐ 7" ☐ 8" ☐ 9" ☐ 10" ☐ 11" ☐ 12" ☐ 13"		
XFTR5	☐ Flip-up Footrest.....	STD	
XFTR9	☐ Depth Adjustable Flip-up Footrest.....	\$120	K0037
XFTR10	☐ Depth & Angle Adjustable Flip-up Footrest.....	\$293	K0040, K0037
XHLP1	☐ Heel Loops (Not available with 13" or narrower Seat Width A or 1-Piece Footrests.).....	\$55	E0951

Elevating Leg Rest

	☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19" ☐ 20" ☐ 21"		
XFA5	☐ Single Depth Adjustable Flip-up Footrest - ☐ Left OR ☐ Right.....	\$185	E0990
XFTR17	☐ Depth Adjustable Flip-up Footrest - Pair.....	\$363	

CUSTOM CONFIGURATION OPTIONS

FRONT WHEELS (See Front Seat Height Ranges Chart on Page 2 for Front Seat Height C limitations.)

XFW24	☐ 3"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire ☐ Silver Hub OR ☐ Black Hub.....	\$206	
XFW14	☐ 4"x.75" Plastic Wheel w/Poly Tire (Not available with Power Adaptable Frame or Heavy Duty Frame).....	NCO	
XFW3	☐ 4"x.75" Light-Up Micro-Caster (Not available with Power Adaptable Frame or Heavy Duty Frame.).....	\$88	
XFW26	☐ 4"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire (Not available with Power Assist, Power Add-on System or Heavy Duty Frame.).....	\$112	
XFW27	☐ 4"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire ☐ Silver Hub OR ☐ Black Hub.....	\$217	
XFW6	☐ 5"x1" Plastic Wheel w/Poly Tire (Not available with Power Adaptable Frame, Heavy Duty Frame or with 90° Front Angle H.).....	NCO	
XFW28	☐ 5"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire (Not available with Power Adaptable Frame, Heavy Duty Frame or 90° Front Angle H.).....	\$103	
XFW30	☐ 5"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire ☐ Silver Hub OR ☐ Black Hub (Not available with 90° Front Angle H.).....	\$233	
XFW9	☐ 6"x1" Plastic Wheel w/Poly Tire (Not available with Power Adaptable Frame, Heavy Duty Frame or 85° or 90° Front Angle H.).....	NCO	
XFW32	☐ 6"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire (Not available with Power Adaptable Frame, Heavy Duty Frame, 85° or 90° Front Angle H.).....	\$119	
XFW33	☐ 6"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire ☐ Silver Hub OR ☐ Black Hub (Not available with 85° or 90° Front Angle H.).....	\$244	
XFW16	☐ 8"x1" Plastic Wheel w/Poly Tire (Not available with Power Assist, Power Add-on, Heavy Duty Frame or 12"-14" Seat Depths. Available with 70° Front Angle H only.).....	NCO	

FRONT FORKS

XFK1	☐ TiLite Standard Fork.....	STD	
XGL1	☐ Glide™ Suspension Fork (Available in Black anodize only. Not available with 8" Front Wheels or Heavy Duty Frame).....	\$526	E1015
	Rider's Weight: ☐ 125 lbs. ☐ 125-180 lbs. ☐ 180-265 lbs.		

ONE ARM DRIVE

XOAD1	☐ One Arm Drive (Available with 24" Mag Wheels only. Rear Seat Height, Camber or axle restrictions may apply.).....	\$1298	E0958
	Drive Handrim Side (Available with standard aluminum handrim only.) ☐ Install on Left ☐ Install on Right		

STANDARD REAR WHEELS (See Rear Seat Height Ranges Chart on Page 2 for Rear Seat Height D limitations.)

XRW14	TiLite SHADOW™.....	STD	
	☐ 22"(501)		
	☐ 24"(540)		
	☐ 25"(559) (Requires 15" minimum Seat Depth B.)		
	☐ 26"(590) (Requires 15" minimum Seat Depth B.)		

PERFORMANCE REAR WHEELS (See Rear Seat Height Ranges Chart on Page 2 for Rear Seat Height D limitations.)

XRW7	Spinergy® SPOX ☐ Black ☐ Red ☐ Blue ☐ White ☐ Green.....	\$915	
	☐ 22"(501) - Std. Hub		
	☐ 24"(540) - Std. Hub		
	☐ 25"(559) - Std. Hub (Requires 15" minimum Seat Depth B.)		
	☐ 26"(590) - Std. Hub (Requires 15" minimum Seat Depth B.)		
XRW8	Spinergy® LX (Available with Heavy Duty Frame only when 0 or 2 degree camber is selected.)		
	☐ Black ☐ Red ☐ Blue ☐ White ☐ Bright Orange ☐ Flamingo Pink ☐ Purple ☐ Green.....	\$1148	
	☐ 22"(501) - Std. Hub		
	☐ 24"(540) - Std. Hub		
	☐ 25"(559) - Std. Hub (Requires 15" minimum Seat Depth B.)		
	☐ 26"(590) - Std. Hub (Requires 15" minimum Seat Depth B.)		

TiLite 2GX Order Form

Effective Aug 2025 (New Items in Red)

FOLDING FRAME

Page 5 of 7

ZRW9 Spinergy® Design Lab LX Wheels (Available with black spokes only.)
☐ Gold Anodize ☐ Green Anodize ☐ Blue Anodize ☐ Red Anodize ☐ Purple Anodize ☐ Orange Anodize.....\$1,364
☐ 22"(501) - Std. Hub
☐ 24"(540) - Std. Hub
☐ 25"(559) - Std. Hub (Requires 15" minimum Seat Depth B.)
☐ 26"(590) - Std. Hub (Requires 15" minimum Seat Depth B.)

XRW16 Spinergy® Carbon Blade (Not available with Heavy Duty Frame.).....\$1,602
☐ 24"(540)
☐ 25"(559) (Requires 15" minimum Seat Depth B.)

XRW17 X-Core Mag (Not available with Heavy Duty Frame, Short Tabs, Quad Release Axles, Spoke Guards or Omit Handrims Option. Requires Aluminum Sliver Anodized Handrims.).....NCO
☐ 22"(501)
☐ 24"(540)

AXLES

XAXL1 ☐ Stainless Quick Release.....STD
XAXL2 ☐ Titanium Quick Release.....\$277
XAXL3 ☐ Quad Release (Not available with Mag or 26"(590) Spinergy Spox or Spinergy LX Rear Wheels.).....\$169

PNEUMATIC REAR TIRES

XRTR1 ☐ Treaded 22"(501), 24"(540), and 26"(590) (Not available with 25"(559) Rear Wheels.).....STD
XRTR2 ☐ Airless Inserts.....\$98
XRTR31 ☐ New MBL SpeedLite 20"(451), 22"(501), 24"(540), 25"(559) and 26"(590).....STD
XRTR30 ☐ 24"(540), 25"(559) and 26"(590) Primo Sentinel.....\$136
XRTR8 ☐ 24"(540), 25"(559) and 26"(590) Schwalbe® Marathon Plus Evolution.....\$152
XRTR32 ☐ 24"(540) Primo C-763 Knobby.....\$104

E2213

SOLID REAR TIRES

XRTR22 ☐ 22"(501), 24"(540) and 25"(559) Primo Orion (Not available with Spinergy Carbon Blade Rear Wheels.).....\$120
XRTR23 ☐ 24"(540) and 25"(559) Primo Active.....\$147

HANDRIMS

XRIM1 ☐ Aluminum – Silver Anodized.....STD
XRIM21 ☐ Aluminum – Black Anodized.....NCO
XRIM12 ☐ Titanium – Satin Finish (Not available with 20"(451) or 22"(501) Rear Wheels, or 24"(540) Knobby tires.).....\$509
XRIM10 ☐ The Surge (All-In-One Oval with Gription™ Strip, Larger Profile).....\$433
XRIM11 ☐ The Surge LT (All-In-One Oval with Gription™ Strip, Smaller Profile).....\$466
XRIM20 ☐ Q-Grip (Exclusive High-Friction, Anti-Slip Coating).....\$358

The Natural-Fit® (Original, Larger Profile) (Not available with 20"(451) Rear Wheels.)

TRIM14 ☐ Standard Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.).....\$428
TRIM15 ☐ Super Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.).....\$466
TRIM16 ☐ Without Thumb Grip.....\$390

The Natural-Fit® LT (Reduced Weight, Smaller Profile) (Not available with 20"(451) Rear Wheels.)

TRIM17 ☐ Standard Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.).....\$440
TRIM18 ☐ Super Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.).....\$514
TRIM19 ☐ Without Thumb Grip.....\$423

Additional Handrim Options

XRIM5 ☐ Plastic Coated – Black.....\$174
XRIM22 ☐ Omit Handrims.....NCO

Tab Length

XRIM7 ☐ Long Tabs.....STD
XRIM8 ☐ Short Tabs (Not available with Mag Rear Wheels or The Natural-Fit® Handrims with Thumb Grip.).....\$28

WHEEL LOCKS (Composite wheel locks not always compatible with 12"-13" Seat Depth. Uni-locks will be substituted at no charge.)

XWLK7 ☐ Compact Composite Push to Lock.....NCO
XWLK8 ☐ Compact Composite Pull to Lock.....NCO
XWLK10 ☐ Extension Handles for Compact Composite Push/Pull to Lock.....\$93

E0961

Additional Wheel Locks

XWLK11 ☐ MBL Composite Scissor Lock (Not available on 2GX Swing Away.).....\$158
XWLK3 Uni-Lock Set At: ☐ Push **OR** ☐ Pull.....\$179
XWLK4 ☐ Extension Handles for Uni-Lock.....\$93
XWLK15 ☐ Grade Aids Push to Lock (Not available with solid rear tires).....\$190

E0961

E0974

TiLite 2GX Order Form

Effective Aug 2025 (New Items in Red)

FOLDING FRAME

Page 6 of 7

BACK UPHOLSTERY (Want more seating options? [Click here](#) to view the new TiLite Seating Order Form.)

XBUP1 ☐ Padded Nylon (Black Only. Not available with Heavy Duty Frame, Folding Back, Folding Stabilizer Bar, Depth Adjustable Frame, Fold-Down Push Handles or Ultrasuede Backrest Stripe.)..... **STD**

Tension Adjustable by Straps (Included with Heavy Duty Frame. Required with Aluminum Fendered Side Guards and Fold-Down Push Handles.)

XBUP3 ☐ Padded Nylon..... **\$369 E2611**
XBUP5 ☐ AIR Breathable Upholstery (Not available with Ultrasuede Backrest Stripe.)..... **\$369 E2611**
XUPC3 Ultrasuede Backrest Stripe ☐ Black (Not available with Padded Nylon or AIR Breathable Back Upholstery.)..... **\$115**
XBUP0 ☐ Omit Back Upholstery..... **NCO**

SEAT UPHOLSTERY

XSUP4 ☐ Tension Adjustable (Black only)..... **STD**
XSUP0 ☐ Omit Seat Upholstery..... **NCO**

SEAT CUSHION (Want more seating options? [Click here](#) to view the new TiLite Seating Order Form.)

XCSH2 ☐ Comfort Curve Cushion (Stretch-Air cover)..... **\$122 E2601**
ZCSH23 ☐ Comfort Acta-Embrace™ Comfort-Tek Zero Elevation Cushion..... **\$456 E2607**
XCSH15 ☐ ROHO® Hybrid Elite SR™ Single Compartment Cushion (Fits seat widths 14" to 19")..... **\$552 E2622**
XCSH14 ☐ ROHO® QUADTRO SELECT® Cushion (Fits seat widths 14" to 19") Profile: ☐ Low ☐ Mid ☐ High..... **\$594 E2624**
XCSH16 ☐ ROHO® Hybrid Select Cushion (Fits seat widths 14" to 18", and 20")..... **\$571 E2624**
XCSH17 ☐ Vicair Active O2 Low (Fits even seat widths A 14" to 20". Non-stock upcharges may apply)..... **\$515 E2624**

CUSTOM CONFIGURATION OPTIONS

SIDE GUARDS

Aluminum

XSDG2 ☐ Rigid Removable ☐ Small **OR** ☐ Std..... **\$217**
XSDG3 ☐ Quad Removable w/Hole ☐ Small **OR** ☐ Std..... **\$233**
XSDG6 ☐ Fendered Adjustable (Requires Tension Adjustable by Straps Back Upholstery. Not available with Angle-Adjustable Fold-Down Seat Backs.)..... **\$363**

2x2 Twill Weave Carbon Fiber

XSDG7 ☐ Rigid Removable ☐ Small **OR** ☐ Std..... **\$336**
XSDG8 ☐ Quad Removable w/Hole ☐ Small **OR** ☐ Std..... **\$347**

☐ **Custom Offset Profiled Sideguards** CMD 007. Please specify desired offset to nearest $\pm\frac{1}{2}$ " (e.g., ± 0.50 ") and sideguard type (rigid or fold-down) in the Notes section of the order form. Sideguard material must be selected. Some limitations may apply.

ARMRESTS (Desk Arms come with built-in Side Guards.)

XARM1 Swing Away Tubular (Not available with Heavy Duty Frame, Folding Stabilizer Bar with Short or Medium Seat Back Height ranges or Depth Adjustable Back.)..... **\$244 E0973**
☐ 9" - 12" **OR** ☐ 12" - 15"
XARM4 ☐ Swing Away - Flip Back Tubular (Not available with Heavy Duty Frame, Folding Stabilizer Bar or Depth Adjustable Back.)..... **\$342 E0973**
XARM2 Removable Height Adjustable DeskArm with Rigid Side Guard (Requires 14" minimum Seat Depth B.)..... **\$363 E0973**
☐ 9" to 12" **OR** ☐ 12" to 15"
☐ Standard Length Pad (10") **OR** ☐ Full Length Pad (14")
XARM3 Removable Height Adjustable Desk Arm with Rigid Side Guard and Transfer Loop (Requires 16" minimum Seat Depth B.)..... **\$406 E0973**
☐ 9" to 12" **OR** ☐ 12" to 15"
☐ Standard Length Pad (10") **OR** ☐ Full Length Pad (14")

Available Armrest Option

XARM5 ☐ Covered in Black Ultrasuede (Not available with Desk Arms)..... **\$120**

SEAT BELTS

XBLT1 ☐ Velcro Adjustable, Black Nylon Webbing..... **\$61 E0978**
XBLT2 ☐ Bodypoint® 2" Push-Button Buckle (Auto-Style.)..... **\$103 E0978**
XBLT4 ☐ Bodypoint® 1" Wide Push-Button Buckle Padded (Auto-Style.) (Not available with 19" - 22" Rear Seat Widths A.)..... **\$196 E0978**
XBLT5 ☐ Bodypoint® 1.5" Wide Push-Button Buckle Padded (Auto-Style.) (Not available with 19" - 22" Rear Seat Widths A.)..... **\$185 E0978**
XBL10 ☐ Bodypoint® Optional Band Clamp Mounting Hardware..... **\$22 K0108**
XBLT16 ☐ Bodypoint® EvoFlex® Pelvic Stabilizer Push Button..... **\$201 E0978**

CALF STRAP

XCLF2 Velcro Adjustable..... **\$44 K0038**
☐ Small (Recommended for 9"- 12" Footrest Widths)
☐ Medium (Recommended for 11"- 15" Footrest Widths)
☐ Large (Recommended for 14"- 17" Footrest Widths)
☐ Xtra Large (Recommended for 16"- 19" Footrest Widths)
XCLF3 Bodypoint® Padded Velcro Adjustable..... **\$109 K0038**
☐ Medium (Recommended for 12"- 15" Footrest Widths)
☐ Large (Recommended for 16"- 20" Footrest Widths)
XCLF1 ☐ Buckle Adjustable, Black Nylon Webbing..... **\$34 K0038**

TiLite 2GX Order Form

Effective Aug 2025 (New Items in Red)

FOLDING FRAME

Page 7 of 7

RESIDUAL LIMB SUPPORT

XRSL1	☐ Residual Limb Support (8" width x 10" length).....	\$428	E1020
XRSL2	☐ Residual Limb Support Pair (8" width x 10" length).....	\$855	E1020

BACKPACK

XBCK3	☐ Permobil Medical Necessities Backpack.....	\$103	
-------	---	-------	--

SEAT POUCH

XPCH3	☐ Permobil Medical Necessities Seat Pouch.....	\$66	
XPCH4	☐ Bodypoint® Mobility Bag.....	\$109	

CUSTOM CONFIGURATION OPTIONS

SPOKE GUARDS

XGRD1	☐ Clear (Not available with Mag Rear Wheels).....	\$147	K0065
-------	--	-------	-------

ANTI-TIPS

XTIP1	☐ Rear Aluminum - Black.....	\$174	E0971
XTIP2	☐ User Friendly Flip-up.....	\$363	E0971
XTIP9	☐ Tip Assist.....	\$115	E0971

ACCESSORIES

XUEST1	Upper Extremity Support Tray ☐ 12" to 14" OR ☐ 13" to 15" Seat Widths A (Requires Armrests.).....	\$315	E0950
XLUG1	☐ Luggage Carriers (Not available with High Mount Swing Away Footrests or 1-Piece Flip-up Footrest.).....	\$158	
XMPCT2	Neoprene Impact Guards with Ultrasuede Stripe (Not available with High Mount Footrests, 90° Hangers or Elevating Leg Rest.) ☐ Black.....	\$115	
XCH1	☐ Crutch Holder.....	\$244	
XAFP1	4" Logo Patch ☐ Army ☐ Navy ☐ Air Force ☐ Marines.....	\$115	

GENERAL

More options and frame modifications may be available through our Creative Mobility Division! Please contact customer service for details or [click here](#) to view options.

Permobil manufactures TiLite manual wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. HCPCS Codes are for reference purposes only. Actual coding and coverage varies by payer. Proper HCPCS coding is not a guarantee of payment. Please verify the appropriate payer's coverage policy.

Get the latest order forms at permobilus.com. Prices, standard features, options and accessories are subject to change without notice. All rights reserved. Trademarks are identified by symbols™ and ®. TiLite, TiFit, the TiLite logo, Ti, the Ti logo, and TiShaft are registered trademarks of TiSport, LLC and are licensed to Permobil. Permobil® is a registered trademark of Permobil. Adobe®, Acrobat®, Acrobat Professional® and Reader® are registered trademarks of Adobe Systems Inc. Bodypoint® is a registered trademark of Bodypoint Inc. Kenda® is a registered trademark of Kenda Rubber Industrial Co. Ltd. ROHO®, AGILITY®, Low Profile® and Quadro Select® are registered trademarks of ROHO Inc. Schwalbe® is a registered trademark of Ralf Bohle GmbH. The Natural-Fit® and Out-Front® are registered trademarks of Three Rivers Holdings, LLC. The ACTA-BACK®, ACTA-EMBRACE™, COMFORT-TEK™, COMPASS®, CURVE®, ELEMENTS®, and STRETCH AIR™ are registered trademarks of THE COMFORT COMPANY® and are licensed to Permobil. Q'Strain® is a registered trademark of Q'Strain Inc.

NOTES

QUOTE THIS CHAIR

If you are using Adobe Acrobat, this order form may be converted into a quote request or into an order by clicking one of the buttons (this feature is available exclusively to Permobil dealers and distributors). A Quote or an Order Confirmation will be sent to you within 24 hours.

ORDER THIS CHAIR