



TiLite Z - Order Form



Dealer Name*:		Account #*:
City, State, Zip:		
Phone:	Email:	ATP:
Rep Name:	P.O. Number:	Quote Number:
Client Reference*:	Client Weight:	
Special considerations that need to be addressed here (e.g., diagnosis):		

*Required field

Frame Style

- | | | |
|---|--------------|-------------|
| ☐ Aluminum Rigid Frame Full adjustability. | K0005, E1235 | \$ 2,750.00 |
| ☐ 250lb Weight Limit | | Standard |
| ☐ 275lb Weight Limit | | No charge |
| ☐ Aluminum Rigid Frame with Transit Tie-Downs Full adjustability. | K0005, E1235 | \$ 2,995.00 |
| ☐ 250lb Weight Limit | | Standard |
| ☐ Heavy Duty Frame 350lb weight limit. Not compatible with Transit Tie-Down. Available with 18"-20" seat widths and 16"-20" seat depths. Requires aluminum casters, provided at no charge. | | \$ 575.00 |
| ☐ 350lb Weight Limit | | |

Transit Options

- | | | |
|---|-------|-----------|
| ☐ Transit Tie-Down Tested to WC-19. Consists of four tie-down brackets and two pelvic belt mounting brackets. Not compatible with Heavy Duty Frame. Please refer to user manual for important information regarding Transit Tie-Down option. | E1023 | \$ 363.00 |
| ☐ Transport Bracket Option Consists of four tie-down brackets. Unoccupied use only. | E1022 | \$ 379.00 |

Frame Finishes

- | | | | | | |
|-----------------|---|--|--|--|-----------|
| Frame Finishes: | ☐ High Gloss Black | ☐ Red Pepper | ☐ Neon Pink | ☐ Ocean Blue Metallic | No charge |
| | ☐ Tangerine Metallic | ☐ Acid Green Metallic | ☐ Midnight Blue Pearl | ☐ White Pearl | |
| | ☐ Black Lava | ☐ Titanium Grey | ☐ Pomegranate Red | ☐ Blue Mist | |
| | ☐ Indigo Purple | ☐ Hyper Blue | ☐ Flame Orange | ☐ OD Green | |
| | ☐ Silver Streak | ☐ Lavendar Dusk | ☐ Black Galaxy | ☐ Seafoam Green | |
| | ☐ Daffodil Yellow | | | | |

Color Anodize Package

- | | | |
|------------------|---|-----------|
| Anodize Package: | ☐ Black | No charge |
| | ☐ Gold ☐ Green ☐ Purple ☐ Orange ☐ Blue ☐ Red | \$ 195.00 |
- Includes caster mount, caster fork and caster caps. If no anodize colour is selected, default will be black.

TiFit: Tailored Individualized Frame

Please refer to Reference Diagrams at end of form for more information.

Seat Width 'A' | Measure from outside of seat tube at back post to the outside of the opposite seat tube at back post.

- | | | |
|-------------|---|-----------------|
| Seat Width: | ☐ 12" ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19" | No charge |
| | ☐ 20" | E2201 \$ 559.00 |

Seat Depth 'B' | Measure from front of back post to the front edge of seat sling.

- | | | |
|-------------|---|-----------------|
| Seat Depth: | ☐ 12" ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19" | No charge |
| | ☐ 20" | E2203 \$ 675.00 |

Seat Width 'A'	12"	☐ 12"x12"	☐ 12"x13"	☐ 12"x14"	☐ 12"x15"	☐ 12"x16"	☐ 12"x17"	☐ 12"x18"		
	13"	☐ 13"x12"	☐ 13"x13"	☐ 13"x14"	☐ 13"x15"	☐ 13"x16"	☐ 13"x17"	☐ 13"x18"		
	14"		☐ 14"x13"	☐ 14"x14"	☐ 14"x15"	☐ 14"x16"	☐ 14"x17"	☐ 14"x18"	☐ 14"x19"	☐ 14"x20**
	15"		☐ 15"x13"	☐ 15"x14"	☐ 15"x15"	☐ 15"x16"	☐ 15"x17"	☐ 15"x18"	☐ 15"x19"	☐ 15"x20**
	16"		☐ 16"x13"	☐ 16"x14"	☐ 16"x15"	☐ 16"x16"	☐ 16"x17"	☐ 16"x18"	☐ 16"x19"	☐ 16"x20**
	17"			☐ 17"x14"	☐ 17"x15"	☐ 17"x16"	☐ 17"x17"	☐ 17"x18"	☐ 17"x19"	☐ 17"x20**
	18"			☐ 18"x14"	☐ 18"x15"	☐ 18"x16"	☐ 18"x17"	☐ 18"x18"	☐ 18"x19"	☐ 18"x20**
	19"			☐ 19"x14"	☐ 19"x15"	☐ 19"x16"	☐ 19"x17"	☐ 19"x18"	☐ 19"x19"	☐ 19"x20**
	20"					☐ 20"x16**	☐ 20"x17**	☐ 20"x18**	☐ 20"x19**	☐ 20"x20**
		12"	13"	14"	15"	16"	17"	18"	19"	20"
Seat Depth 'B'										

Custom Frame Depth | Measure from front of seat sling to beginning of frame bend. Frame bend will start 1" in front of seat sling unless custom frame depth is specified, or the combination of specifications affect the rolling dynamics of the chair.

Additional Frame Depth: ☐ 1" Standard
☐ 0.5" ☐ 1.5" ☐ 2" No charge

Front Seat Height 'C' | Measure from floor to top of seat tube at front edge of seat sling.

Front Seat Height: ☐ 16" ☐ 16.5" ☐ 17" ☐ 17.5" ☐ 18" ☐ 18.5" ☐ 19" ☐ 19.5" ☐ 20" ☐ 20.5" ☐ 21"

Rear Seat Height 'D' | Measure from floor to top of seat tube at back post. Will be optimized to achieve desired rear seat height allowing +/- 1" adjustment.

Rear Seat Height: ☐ 13" ☐ 13.5" ☐ 14" ☐ 14.5" ☐ 15" ☐ 15.5" ☐ 16" ☐ 16.5" ☐ 17* ☐ 17.5" ☐ 18" ☐ 18.5"
☐ 19" ☐ 19.5" ☐ 20"

Rear Wheel	20" (451)	13" - 18.5"
	22" (501)	14" - 19.5"
	24" (540)	15" - 20.5"
	25" (559)	15.5" - 20.5"
	26" (590)	16" - 20.5"
Rear Seat Height		

Front Angle 'E' | Measure from frame front to floor.

Front Angle: ☐ 75° ☐ 80° ☐ 85° ☐ 90°
85° not available with 6" front wheels. 90° not available with 5" or 6" front wheels.

Center of Gravity 'F' | Measure from front of back post at the seat tube to the centre of the rear axle.

Center of Gravity: ☐ 1" ☐ 1.5" ☐ 2" ☐ 2.5" ☐ 3" ☐ 3.5" ☐ 4" ☐ 4.5" ☐ 5"
☐ Amputee Axle Plate E0959 \$ 225.00
Amputee axle plate is set to 1" behind back canes.

Footrest

Footrest Width 'G'

Seat Width 'A'	12"	☐ 9.5"	☐ 8"	
	13"	☐ 10.5"	☐ 9"	
	14"	☐ 11.5"	☐ 10"	
	15"	☐ 12.5"	☐ 11"	☐ 9"
	16"	☐ 13.5"	☐ 12"	☐ 10"
	17"	☐ 14.5"	☐ 13"	☐ 11"
	18"	☐ 15.5"	☐ 14"	☐ 12"
	19"	☐ 16.5"	☐ 15"	☐ 13"
	20"	☐ 17.5"	☐ 16"	☐ 14"
		No Taper	1.5" Taper	3.5" Taper
Footrest Width 'G'				

Standard Footrests

Seat to Footrest 'H': ☐ 13.5" ☐ 14" ☐ 14.5" ☐ 15" ☐ 15.5" ☐ 16" ☐ 16.5" ☐ 17" ☐ 17.5" ☐ 18" ☐ 18.5" ☐ 19"
Measure from front edge of seat sling to top rear of footrest. Must allow for at least 2" of ground clearance.

☐ Aluminum Open Loop Footrest No charge

☐ Angle Adjustable Footrest K0040 \$ 185.00

☐ Angle Adjustable Flip-Back Footrest K0040 \$ 244.00

High Mount Footrests | *Taper will affect displayed ranges.*

Seat to Footrest 'H': ☐ 7"-9" (4" risers) ☐ 9"-11" (2" risers) ☐ 11"-13"
☐ Set at

Measure from front edge of seat sling to top rear of footrest. Measurement must be at least 2" shorter than front seat height.

☐ Angle Adjustable High Mount Footrest K0037, K0040 \$ 155.00

Footrest Accessories

☐ Safe Tread Non-Skid Adhesive Patches \$ 20.00
☐ ABS Footrest Cover *Only for open loop footrest.* \$ 30.00

Clamp On Residual Limb Support

☐ Residual Limb Support, 8"W x 10"L, Single \$ 495.00
☐ Install on Left ☐ Install on Right E1020
☐ Residual Limb Support, 8"W x 10"L, Pair E1020 \$ 895.00

Calf Straps

☐ Velcro Adjustable Calf Strap *Select footrest width range.* \$ 44.00
☐ Small (9"-12") ☐ Medium (11"-15") ☐ Large (14"-17") ☐ Extra Large (16"-19") K0038
☐ Bodypoint Padded Calf Strap, Velcro Adjustable *Select footrest width range.* \$ 109.00
☐ Small (8"-11") ☐ Medium (11"-15") ☐ Large (16"-19") K0038
☐ Buckle Adjustable Calf Strap, Black Nylon Webbing K0038 \$ 50.00

Backrest

Back Type

☐ Folding Aluminum Adjustable Height Back No charge

Back Options

☐ Locking Backrest *Locks in up and down positions.* No charge
☐ Non-Locking Backrest *Locks in up position only.* No charge
☐ Depth Adjustable Back Frame *Not available with 19"-20" seat depths. Allows 2" of rearward adjustment in 1" increments. Bolt on height adjustable push handles have a limited height range with this option.* \$ 185.00
☐ Rigidizer Bar No charge
☐ 3" Deep (Standard) ☐ 4" Deep

Push Handles

☐ Integrated Push Handles No charge
☐ Integrated Fold-Down Push Handles *Requires tension adjustable by straps back upholstery.* \$ 125.00
☐ Bolt-On Height-Adjustable Push Handles *Requires 4" deep backrest rigidizer bar.* \$ 363.00
☐ Omit Push Handles No charge

Back Height 'I' | Set at 0.5" increments.

Push Handles	None / Bolt-On	☐ 9" - 11.5"	☐ 12" - 15.5"	☐ 16" - 20.5"	Set at:
	Integrated	☐ 10" - 12.5"	☐ 13.5" - 16.5"	☐ 17.5" - 20.5"	Set at:
	Fold-Down	☐ 12" - 13.5"	☐ 14" - 17.5"	☐ 18" - 20.5"	Set at:
		Short	Medium	Tall	Specify
Back Height 'I'					

Back Angle 'J' | Measure angle from front of back post to floor.

Back Angle: ☐ 80° ☐ 83° ☐ 86° ☐ 89° ☐ 92° ☐ 95° ☐ 98° ☐ 101°

Back Supports

Back Upholstery

- ☐ **Standard Back Upholstery, Tension Adjustable** *Includes privacy flap.*
- ☐ **AIR Breathable Back Upholstery, Tension Adjustable** *Includes privacy flap.* E2611, E2612 \$ 369.00
- ☐ Removable Lumbar Contour Pack ☐ Removable Lateral Contour Pack
- ☐ **Omit Back Upholstery** No charge
- ☐ **4" Military Service Patch** \$ 115.00
- ☐ Army ☐ Navy ☐ Air Force ☐ Marines ☐ Coast Guard ☐ Purple Heart

ROHO Back Supports | See information at end of form regarding product shipment.

- ☐ **ROHO Agility Active Aluminium Back, Min Contour with Foam** *Not compatible with Heavy Duty Frame. Inc fixed hardware. Select shell height.* E2620 \$ 738.00
- ☐ **ROHO Agility Active Aluminium Back, Mid Contour with Foam** *Not compatible with Heavy Duty Frame. Inc fixed hardware. Select shell height.* E2620 \$ 738.00
- ☐ 8" ☐ 10"
- ☐ ROHO Privacy Shield \$ 65.00
- ☐ **ROHO Agility Carbon Back, Min Contour with Foam/Air Insert** *Not compatible with Heavy Duty Frame. Inc fixed hardware. Select shell height.* E2620 \$ 1,202.00
- ☐ **ROHO Agility Carbon Back, Mid Contour with Foam/Air Insert** *Not compatible with Heavy Duty Frame. Inc fixed hardware. Select shell height.* E2620 \$ 1,202.00
- ☐ 10" ☐ 13" ☐ 16"
- ☐ ROHO Privacy Shield \$ 65.00

Comfort Back Supports | Backs over 20" wide or tall subject to increased pricing. See information at end of form regarding product shipment.

- ☐ **Comfort Unite 1 Back** *Replaces Comfort Acta-Back. Inc 2-Point Quick Release Compass Hardware. Select shell height.* E2620, E2621 \$ 712.00
- ☐ 12" ☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22"
- ☐ Comfort-Tek Cover ☐ Diamond Spacer Cover
- ☐ Comfort Privacy Flap \$ 52.00
- ☐ **Comfort Embrace Back** *Replaces Comfort Acta-Embrace. Inc 2-Point Quick Release Compass Hardware. Select shell height.* E2613, E2614 \$ 700.00
- ☐ 10" ☐ 12" ☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22"
- ☐ Comfort-Tek Cover ☐ Diamond Spacer Cover
- ☐ Comfort Privacy Flap \$ 52.00
- ☐ **Comfort Relief Back** *Replaces Comfort Acta-Relief. Inc 2-Point Quick Release Compass Hardware. Select shell height.* E2620 \$ 886.00
- ☐ 10" ☐ 12" ☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22"
- ☐ Comfort-Tek Cover ☐ Diamond Spacer Cover
- ☐ Comfort Privacy Flap \$ 52.00

Seat Options

Seat Upholstery

- ☐ **Tension Adjustable Seat Upholstery, Black**
- ☐ **Aluminum Solid Seat Pan** *Limited sizes available, visit the Permobil Store for more information.* E2231 \$ 309.00
- ☐ **Omit Seat Upholstery**

Seat Cushions | Cushions over 20" wide subject to increased pricing. See information at end of form regarding product shipment.

- ☐ **Comfort Curve Cushion** *Comfort-Tek cover.* E2601, E2602 \$ 122.00
- ☐ **Comfort Acta-Embrace™ Zero Elevation Cushion** *Comfort-Tek cover.* E2607, E2608 \$ 428.00
- ☐ **Comfort M2 Zero Elevation Cushion** *Comfort-Tek cover.* E2607, E2608 \$ 474.00

☐ ROHO Single Compartment Cushion	E2622, E2623	\$ 534.00
☐ Low Profile ☐ Mid Profile ☐ High Profile		
☐ ROHO Quadtro Select Cushion	E2624, E2625	\$ 534.00
☐ Low Profile ☐ Mid Profile ☐ High Profile		
☐ ROHO Hybrid Elite Cushion	E2622, E2623	\$ 552.00
☐ Single Compartment ☐ Dual Compartment		
☐ ROHO Hybrid Select Cushion	E2624	\$ 571.00

Front Forks

☐ TiLite Standard Fork		Standard
☐ TiLite Slipstream Single-Sided Fork		\$ 185.00
☐ Glide Suspension Fork <i>Black anodize only. Select user weight (maximum 265lb). Not compatible with Heavy Duty Frame.</i>	E1015	\$ 526.00
☐ 125lbs ☐ 125-180lbs ☐ 180-265lbs		

Front Wheels

☐ 3"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire		\$ 205.00
☐ Silver Hub ☐ Black Hub		
☐ 4"x0.75" Plastic Wheel w/Poly Tire		\$ 65.00
☐ 4"x0.75" Light-Up Micro Caster		\$ 125.00
☐ 4"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire		No charge
☐ 4"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire		\$ 217.00
☐ Silver Hub ☐ Black Hub		
☐ 5"x1" Plastic Wheel w/Poly Tire <i>Not available with 90° front frame or slipstream fork.</i>		\$ 65.00
☐ 5"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire <i>Not available with 90° front frame.</i>		No charge
☐ 5"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire <i>Not available with 90° front frame.</i>		\$ 233.00
☐ Silver Hub ☐ Black Hub		
☐ 6"x1" Plastic Wheel w/Poly Tire <i>Not available with 85°-90° front frame or slipstream fork.</i>		\$ 65.00
☐ 6"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire <i>Not available with 85°-90° front frame or slipstream fork.</i>		No charge
☐ 6"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire <i>Not available with 85°-90° front frame or slipstream fork.</i>		\$ 244.00
☐ Silver Hub ☐ Black Hub		

Rear Wheels

See rear seat height ranges chart for rear seat height limitations.

☐ TiLite Standard Rear Wheel <i>Select size. When ordering 24" wheels with treaded pneumatic or treaded tire with airless inserts and aluminum silver or aluminum black handrims with long tabs, the chair will ship with TiLite Standard rear wheels. Otherwise, the chair will ship TiLite Shadow™ rear wheels.</i>		No charge
☐ 20" (451) ☐ 22" (501) ☐ 24" (540) ☐ 25" (559) ☐ 26" (590)		
☐ Spinergy SPOX Rear Wheel <i>Select size and spoke color.</i>		\$ 1,195.00
☐ 20" (451) ☐ 22" (501) ☐ 24" (540) ☐ 25" (559) ☐ 26" (590)		
☐ Black ☐ Red ☐ Blue ☐ White		
☐ Spinergy LX Rear Wheel <i>Select size and spoke color.</i>		\$ 1,295.00
☐ 22" (501) ☐ 24" (540) ☐ 25" (559) ☐ 26" (590)		
☐ Black ☐ Red ☐ Blue ☐ White ☐ Orange ☐ Pink ☐ Purple		
☐ Spinergy Design Lab LX Rear Wheel <i>Black spokes only. Select size and anodizing color.</i>		\$ 1,995.00
☐ 22" (501) ☐ 24" (540) ☐ 25" (559) ☐ 26" (590)		
☐ Gold ☐ Green ☐ Blue ☐ Red ☐ Purple ☐ Orange		
☐ X-Core Mag Rear Wheel <i>Select size. Not available with short tabs, quad release axles, spoke guard or omit handrims options. Requires aluminum silver anodized handrims.</i>		\$ 225.00
☐ 20" (451) ☐ 22" (501) ☐ 24" (540)		

Camber 'K'

Camber Angle: ☐ 0° ☐ 2° ☐ 4° ☐ 6°

Rear wheel spacing will be set to minimum value allowed by your configuration.

Axles

- | | |
|--|-----------|
| ☐ Stainless Quick Release Axle | No charge |
| ☐ Quad Release Axle <i>Not available with Mag or 24"/25" Spinergy Spox or LX rear wheels.</i> | \$ 169.00 |

Pneumatic Rear Tires

- | | |
|---|-----------------|
| ☐ Treaded Tire 20" (451), 22"(501), 24"(540), 26"(590) <i>Not available with 25" rear wheels.</i> | No charge |
| ☐ Treaded Tire w/Airless Inserts <i>Not available with 25" rear wheels.</i> | E2213 \$ 145.00 |
| ☐ MBL SpeedLite Tire 20" (451), 22"(501), 24"(540), 25"(559), 26"(590) | No charge |
| ☐ Primo Sentinel Tire 20" (451), 22"(501), 24"(540), 25"(559), 26"(590) | \$ 159.00 |
| ☐ Schwalbe Marathon Plus Evolution Tire 24"(540), 25"(559), 26"(590) <i>Not available with Mag or 20"-22" rear wheels.</i> | \$ 169.00 |
| ☐ Primo C-764 Knobby Tire 24"(540) <i>24" rear wheels only.</i> | \$ 129.00 |

Solid Rear Tires

- | | |
|--|-----------|
| ☐ Primo Orion Tire 20" (451), 22"(501), 24"(540), 25"(559) <i>Not available with 26" rear wheels.</i> | \$ 160.00 |
| ☐ Primo Active Tire 24"(540), 25"(559) <i>Available with 24" and 25" rear wheels only.</i> | \$ 160.00 |

Handrims

- | | |
|--|-----------|
| ☐ Aluminum Handrim | No charge |
| ☐ Silver Anodized ☐ Black Anodized | |
| ☐ Q-Grip Handrim <i>Exclusive high-friction, anti-slip coating. Not available with 20" rear wheels.</i> | \$ 495.00 |
| ☐ Plastic Coated Handrim, Black | \$ 174.00 |
| ☐ The Surge Handrim <i>Large profile. All-In-One oval with Gription strip. Not available with 20" rear wheels.</i> | \$ 495.00 |
| ☐ The Surge LT Handrim <i>Small profile. All-In-One oval with Gription strip. Not available with 20" rear wheels.</i> | \$ 466.00 |
| ☐ The Natural-Fit Handrim, Standard Grip <i>Original, large profile. Not with short tabs.</i> | \$ 428.00 |
| ☐ The Natural-Fit Handrim, Super Grip <i>Original, large profile. Not with short tabs.</i> | \$ 466.00 |
| ☐ The Natural-Fit Handrim, Without Thumb Grip <i>Original, large profile.</i> | \$ 390.00 |
| ☐ The Natural-Fit LT Handrim, Standard Grip <i>Reduced weight, small profile. Not with short tabs.</i> | \$ 450.00 |
| ☐ The Natural-Fit LT Handrim, Super Grip <i>Reduced weight, small profile. Not with short tabs.</i> | \$ 514.00 |
| ☐ The Natural-Fit LT Handrim, Without Thumb Grip <i>Reduced weight, small profile.</i> | \$ 423.00 |
| ☐ Omit Handrims | No charge |

Handrim Tab Length

- | | |
|--|----------|
| ☐ Long Tabs | Standard |
| ☐ Short Tabs <i>Not available with Mag rear wheels or The Natural-Fit handrims with thumb grip.</i> | \$ 39.00 |

Spoke Guards

- | | |
|---|-----------------|
| ☐ Clear Lexan Spoke Guards with TiMotion Logo <i>Not available with Mag rear wheels.</i> | K0065 \$ 179.00 |
|---|-----------------|

Wheel Locks

- | | |
|--|-----------------|
| ☐ Aluminum Push to Lock | No charge |
| ☐ Aluminum Pull to Lock | No charge |
| ☐ Extension Handles for Aluminum Push/Pull to Lock | E0961 \$ 75.00 |
| ☐ Compact Composite Push to Lock | \$ 129.00 |
| ☐ Compact Composite Pull to Lock | \$ 129.00 |
| ☐ Extension Handles for Compact Composite Push/Pull to Lock | E0961 \$ 93.00 |
| ☐ MBL Composite Scissor Lock | \$ 158.00 |
| ☐ MBL Performance Scissor Lock | \$ 175.00 |
| ☐ Grade Aids, Push to Lock <i>Requires pneumatic tires.</i> | E0974 \$ 395.00 |

Side Guards

☐ Aluminum Rigid Side Guards, Removable	Select size.	E0956, E1034	\$ 217.00
☐ Small ☐ Standard			
☐ Aluminum Quad Side Guards, Removable w/Hole	Select size.	E0956, E1034	\$ 233.00
☐ Small ☐ Standard			
☐ Aluminum Fold Down Side Guards, Rigid	Select size.	E0956, E1034	\$ 363.00
☐ Small ☐ Standard			

Armrests

☐ Omit Armrests			No charge
☐ Fixed Height Non-Removable Armrests			No charge
☐ Short ☐ Tall			
☐ Height Adjustable Tubular Armrests		E0973	\$ 244.00
☐ 6.5"-9.5" ☐ 9.5"-12.5"			
☐ Height Adjustable Flip Back Tubular Armrests	Height adjustable between 8.5"-11.5".	E0973	\$ 342.00
☐ Height Adjustable Tubular Armrests with Side Guards		E0973	\$ 425.00
☐ 6.5"-9.5" ☐ 9.5"-12.5"			
☐ Removable Height Adjustable T-Style Armrests	Requires 15" minimum seat depth.	E0973	\$ 375.00
☐ 8"-12" ☐ 10"-14"			
☐ Desk (10") ☐ Full (14")			

Positioning Belts

☐ Velcro Adjustable Belt, Black Nylon Webbing		E0978	\$ 61.00
☐ Bodypoint 1" Belt, Push-Button Padded (Auto Style)	Not available with 19" or 20" rear seat widths.	E0978	\$ 185.00
☐ Bodypoint 1.5" Belt, Push-Button Padded (Auto Style)	Not available with 19" or 20" rear seat widths.	E0978	\$ 196.00
☐ Bodypoint 2" Belt, Push-Button (Auto Style)		E0978	\$ 103.00
☐ Bodypoint EvoFlex Pelvic Stabilizer, Push Button		E0978	\$ 201.00
☐ Q'Straint WC-19 Wheelchair Anchored Pelvic Belt	Recommended for use with transit tie down.	E0978	\$ 217.00
☐ Bodypoint Optional Band Clamp Mounting Hardware			\$ 70.00

Accessories

Anti-Tips

☐ Rear Anti-Tips, Aluminum Black	E0971	\$ 174.00
☐ Rear Anti-Tips, User Friendly Flip-Up	E0971	\$ 363.00

Soft Good Accessories

☐ Permobil Medical Necessities Backpack		\$ 103.00
☐ Permobil Medical Necessities Seat Pouch		\$ 66.00
☐ Bodypoint Mobility Bag		\$ 109.00
☐ Neoprene Impact Guards, Black	Not available with high mount footrests.	\$ 195.00

Functional Accessories

☐ Upper Extremity Support Tray	Requires T-style armrests.		\$ 315.00
☐ 12"-14" seat widths ☐ 13"-15" seat widths		E0950	
☐ Luggage Carriers			\$ 158.00
☐ Crutch Holder		E2207	\$ 244.00
☐ Left ☐ Right			

SmartDrive Power Assist System

☐ Permobil SmartDrive MX2+ <i>Not available with rear seat heights less than 13.5". Maximum user weight capacity 331 lbs.</i>	E0986	\$ 7,529.00
☐ Switch Control Buttons		No charge
☐ SwitchControl Button with Mono Jack Option		No charge
☐ Buddy Button <i>For use with SwitchControl button with mono jack option.</i>		No charge
☐ Carrying Bag		\$ 84.00
☐ Universal Padded Handle		\$ 22.00


General

Permobil manufactures TiLite manual wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. HCPCS Codes are for reference purposes only. Actual coding and coverage varies by payer. Proper HCPCS coding is not a guarantee of payment. Please verify the appropriate payer's coverage policy. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

Notes

Measurement Guide

- **A Seat Width** *Measure from outside of seat tube at back post to outside of opposite seat tube at back post.*
- **B Seat Depth** *Measure from front of back post to front edge of seat sling. Frame bend will start 1" in front of seat sling unless Custom Frame Depth is specified. Custom Frame Depth is measured from front of seat sling to beginning of frame bend.*
- **C Front Seat Height** *Measure from floor to top of seat tube at front edge of seat sling.*
- **D Rear Seat Height** *Measure from floor to top of seat tube at back post.*
- **E Front Angle** *Measure from front frame to floor.*
- **F Center of Gravity** *Measure from front of back post at seat tube to center of rear axle.*
- **G Footrest Width** *Measure from inside of front frame tube to inside of opposite front frame tube.*
- **H Seat to Footrest** *Measure from front edge of seat sling to top rear of footrest.*
- **I Back Height** *Measure from top of back post to top of seat tube at rear of frame.*
- **J Back Angle** *Measure angle from front of back post to floor.*
- **K Camber** *This will affect overall width of chair.*



Seating and Positioning Shipping Information

The following Seating and Positioning products will ship with the chair. Cushions are stocked with Comfort-Tek covers and backs are stocked with Stretch Air covers. All other Seating and Positioning products will be drop-shipped from the manufacturer. (Width x Depth)

Shipped With
Chair

16"x16"	16"x18"	16"x16"	16"x16"	16"x16"	16"x16"
16"x18"	18"x16"	16"x18"	18"x16"	16"x18"	18"x18"
18"x18"	18"x18"	18"x18"	18"x18"	18"x18"	20"x16"
			20"x16"		20"x18"
			20"x20"		

Curve

Elements

Acta Embrace

M2 Zero Elevation

ROHO Single HP

ROHO Quadtro

Seat Cushions

Shipped
With Chair

16"x12"	16"x16"	16"x12"	16"x10"	16"x10"
16"x16"	18"x16"	16"x16"	18"x10"	18"x10"
18"x12"		18"x12"		
18"x16"		18"x16"		

Comfort Unite 1

Comfort Embrace

Comfort Relief

ROHO Agility Min

ROHO Agility Mid

Back Supports