

Pacer dimensions

Key user dimensions (inches)		K610 mini	K620 small	K630 medium	K640 large	K650 XL
Elbow height		15½–20½	18½–27½	24–35	32–47	34–49
 Key user dimension: elbow height Measure the vertical distance from the bent elbow to the floor while the user is standing upright. Choose the size that allows for growth.						
Item dimensions (inches)		K610 mini	K620 small	K630 medium	K640 large	K650 XL
Floor to top of arm prompt		15½–20½	18½–27½	24–35	32–47	34–49
Floor to top of MPS (not tilted)*		N/A	10–20	9½–23½	17–35	19½–37½
Overall width	Standard base:	20½	23	26	28	31½
	Utility base:	N/A	N/A	30	32	N/A
	Treadmill/stability base:	N/A	35	35	35	N/A
Overall length	Standard base:	22½	26	30	36	39¾
	Utility base:	N/A	N/A	36	40½	N/A
	Treadmill/stability base:	N/A	41	41	41	N/A
Frame height	Standard & utility bases:	12½	16–21	21½–27	28½–38½	31–41
	Treadmill/stability base:	N/A	22–34	23½–36	29½–46½	N/A
Base height (without upper frame)	Standard & utility bases:	N/A	11	14½	16	18½
	Treadmill/stability base:	N/A	21½–28	21½–28	21½–28	N/A
Frame weight (lb)	Standard base:	10¼	11½	15	15½	22½
	Utility base:	N/A	N/A	18½	19¼	N/A
	Treadmill/stability base:	N/A	27½	27½	27½	N/A
	Standard upper:	N/A	6	7½	9½	11
	Dynamic upper:	N/A	7½	11½	13½	15
Dynamic upper movement	Vertical:	N/A	1¼	2	2	2
	Horizontal:	N/A	1½	2	2	2
Maximum treadmill width	Treadmill/stability base:	N/A	29	29	29	N/A
Maximum treadmill height	Treadmill/stability base:	N/A	11	11	11	N/A
Maximum working load (lb)		50	75	150	200	250

*This measurement is taken at the lowest frame and MPS height, and the highest frame and MPS height.

Pelvic support

User dimensions (inches)


Indicated by the red arrow below

Key user dimension:

half of thigh circumference


measured at the groin.

K556 small	K557 med.	K558 large
4 – 8	6 – 10	8 – 14




Product dimensions (inches)			
	K556 small	K557 med.	K558 large
a	6½	8½	12
b	12	16	20
c	7	9	12

Hip positioner

Dimensions (inches) – without pad		
		
	K513 small	K533 large
a Crotch length	2½	3½
b Width	1¼	1¾
c Length	9	11

MPS saddle

Dimensions (inches)		
		
	K672/K673 small/med	K674 large
a Width at rear	3½	5¼
b Width at front	3	4½
c Length	9½	10¾



All five sizes are available in five colors: red, blue, gray, lime and pink.



Retail order form - January 2023

Use dimension chart to select appropriate size.

All Rifton products are made to order and are non-returnable and non-exchangeable.

		* Required	K610 mini	K620 small	K630 medium	K640 large	K650 XL		
* Color	CHOOSE	Red	<input type="checkbox"/>	Red	<input type="checkbox"/>	Red	<input type="checkbox"/>	Red	<input type="checkbox"/>
		Blue	<input type="checkbox"/>	Blue	<input type="checkbox"/>	Blue	<input type="checkbox"/>	Blue	<input type="checkbox"/>
		Gray	<input type="checkbox"/>	Gray	<input type="checkbox"/>	Gray	<input type="checkbox"/>	Gray	<input type="checkbox"/>
		Lime	<input type="checkbox"/>	Lime	<input type="checkbox"/>	Lime	<input type="checkbox"/>	Lime	<input type="checkbox"/>
		Pink	<input type="checkbox"/>	Pink	<input type="checkbox"/>	Pink	<input type="checkbox"/>	Pink	<input type="checkbox"/>
* Upper frame	CHOOSE	Dynamic	N/A	K621 \$2,096 <input type="checkbox"/>	K631 \$2,776 <input type="checkbox"/>	K641 \$3,124 <input type="checkbox"/>	K651 \$3,515 <input type="checkbox"/>		
		Standard	N/A	K622 \$1,116 <input type="checkbox"/>	K632 \$1,460 <input type="checkbox"/>	K642 \$1,721 <input type="checkbox"/>	K652 \$2,511 <input type="checkbox"/>		
Base *	Standard base	CHOOSE	With odometer	K613 \$1,622 <input type="checkbox"/>	K623 \$1,314 <input type="checkbox"/>	K633 \$1,440 <input type="checkbox"/>	K643 \$1,720 <input type="checkbox"/>	K653 \$2,175 <input type="checkbox"/>	
			Without odometer	K614 \$1,342 <input type="checkbox"/>	K624 \$701 <input type="checkbox"/>	K634 \$844 <input type="checkbox"/>	K644 \$1,197 <input type="checkbox"/>	K654 \$1,877 <input type="checkbox"/>	
	Utility base (larger wheels)	CHOOSE	With odometer	N/A	N/A	K635 \$1,897 <input type="checkbox"/>	K645 \$2,247 <input type="checkbox"/>	N/A	
			Without odometer	N/A	N/A	K636 \$1,616 <input type="checkbox"/>	K646 \$1,897 <input type="checkbox"/>	N/A	
	Treadmill/stability base (comes in gray only)			N/A	K676 \$2,340 <input type="checkbox"/>	K676 \$2,340 <input type="checkbox"/>	K676 \$2,340 <input type="checkbox"/>	N/A	
Arm supports (pair)	CHOOSE	Small arm prompts	K512 \$845 <input type="checkbox"/>	K512 \$845 <input type="checkbox"/>	K512 \$845 <input type="checkbox"/>	K512 \$845 <input type="checkbox"/>	K512 \$845 <input type="checkbox"/>		
		Large arm prompts	N/A	N/A	K532 \$968 <input type="checkbox"/>	K532 \$968 <input type="checkbox"/>	K532 \$968 <input type="checkbox"/>		
		Arm platforms	N/A	N/A	K697 \$738 <input type="checkbox"/>	K697 \$738 <input type="checkbox"/>	K697 \$738 <input type="checkbox"/>		
		Arm platforms w/ handgrips	N/A	N/A	K696 \$849 <input type="checkbox"/>	K696 \$849 <input type="checkbox"/>	K696 \$849 <input type="checkbox"/>		
Handbrake for utility base (requires arm supports)			N/A	N/A	K687 \$570 <input type="checkbox"/>	K688 \$570 <input type="checkbox"/>	N/A		
Handlebars (pair)	CHOOSE	Small	K625 \$426 <input type="checkbox"/>	K625 \$426 <input type="checkbox"/>	K625 \$426 <input type="checkbox"/>	K625 \$426 <input type="checkbox"/>	K625 \$426 <input type="checkbox"/>		
		Large	N/A	N/A	K675 \$448 <input type="checkbox"/>	K675 \$448 <input type="checkbox"/>	K675 \$448 <input type="checkbox"/>		
Hand loops (pair)			K579 \$331 <input type="checkbox"/>	K579 \$331 <input type="checkbox"/>	K579 \$331 <input type="checkbox"/>	K579 \$331 <input type="checkbox"/>	K579 \$331 <input type="checkbox"/>		
Chest prompt	CHOOSE	Small	K514 \$977 <input type="checkbox"/>	K514 \$977 <input type="checkbox"/>	K514 \$977 <input type="checkbox"/>	N/A	N/A		
		Medium	N/A	K524 \$1,037 <input type="checkbox"/>	K524 \$1,037 <input type="checkbox"/>	K524 \$1,037 <input type="checkbox"/>	K524 \$1,037 <input type="checkbox"/>		
		Large	N/A	N/A	K534 \$1,100 <input type="checkbox"/>	K534 \$1,100 <input type="checkbox"/>	K534 \$1,100 <input type="checkbox"/>		
Chest pad			N/A	K629 \$638 <input type="checkbox"/>	K639 \$686 <input type="checkbox"/>	K649 \$734 <input type="checkbox"/>	N/A		
Multi-Position Saddle	CHOOSE	Small	N/A	K672 \$1,966 <input type="checkbox"/>	N/A	N/A	N/A		
		Medium	N/A	N/A	K673 \$1,889 <input type="checkbox"/>	N/A	N/A		
		Large	N/A	N/A	K674 \$1,966 <input type="checkbox"/>	K674 \$1,966 <input type="checkbox"/>	K674 \$1,966 <input type="checkbox"/>		
MPS seat cover	CHOOSE	Small/Medium	N/A	K627 \$152 <input type="checkbox"/>	K627 \$152 <input type="checkbox"/>	N/A	N/A		
		Large	N/A	N/A	K628 \$152 <input type="checkbox"/>	K628 \$152 <input type="checkbox"/>	K628 \$152 <input type="checkbox"/>		
Pelvic support (includes handholds)	CHOOSE	Small	K556 \$699 <input type="checkbox"/>	K556 \$699 <input type="checkbox"/>	K556 \$699 <input type="checkbox"/>	N/A	N/A		
		Medium	N/A	K557 \$742 <input type="checkbox"/>	K557 \$742 <input type="checkbox"/>	K557 \$742 <input type="checkbox"/>	K557 \$742 <input type="checkbox"/>		
		Large	N/A	N/A	K558 \$781 <input type="checkbox"/>	K558 \$781 <input type="checkbox"/>	K558 \$781 <input type="checkbox"/>		
Hip positioner (includes handholds)	CHOOSE	Small	K513 \$699 <input type="checkbox"/>	K513 \$699 <input type="checkbox"/>	K513 \$699 <input type="checkbox"/>	K513 \$699 <input type="checkbox"/>	K513 \$699 <input type="checkbox"/>		
		Large	N/A	N/A	K533 \$781 <input type="checkbox"/>	K533 \$781 <input type="checkbox"/>	K533 \$781 <input type="checkbox"/>		
Hip positioner pad	CHOOSE	Small	K541 \$180 <input type="checkbox"/>	K541 \$180 <input type="checkbox"/>	K541 \$180 <input type="checkbox"/>	K541 \$180 <input type="checkbox"/>	K541 \$180 <input type="checkbox"/>		
		Large	N/A	N/A	K542 \$237 <input type="checkbox"/>	K542 \$237 <input type="checkbox"/>	K542 \$237 <input type="checkbox"/>		
Thigh prompts (pair)	CHOOSE	Small	N/A	K515 \$699 <input type="checkbox"/>	K515 \$699 <input type="checkbox"/>	K515 \$699 <input type="checkbox"/>	N/A		
		Large	N/A	N/A	K535 \$845 <input type="checkbox"/>	K535 \$845 <input type="checkbox"/>	K535 \$845 <input type="checkbox"/>		
Ankle prompts for standard base (pair)			K681 \$318 <input type="checkbox"/>	K682 \$336 <input type="checkbox"/>	K683 \$366 <input type="checkbox"/>	K684 \$423 <input type="checkbox"/>	K684 \$423 <input type="checkbox"/>		
Ankle prompts for utility base (pair)			N/A	N/A	K693 \$366 <input type="checkbox"/>	K694 \$423 <input type="checkbox"/>	N/A		
Ankle prompts for treadmill/stability base (pair)			N/A	K664 \$423 <input type="checkbox"/>	K664 \$423 <input type="checkbox"/>	K664 \$423 <input type="checkbox"/>	N/A		
Additional handholds (pair)			K519 \$191 <input type="checkbox"/>	K519 \$191 <input type="checkbox"/>	K519 \$191 <input type="checkbox"/>	K519 \$191 <input type="checkbox"/>	K519 \$191 <input type="checkbox"/>		
Communication tray			K518 \$647 <input type="checkbox"/>	K698 \$675 <input type="checkbox"/>	K698 \$675 <input type="checkbox"/>	K698 \$675 <input type="checkbox"/>	K698 \$675 <input type="checkbox"/>		
Attendant guide bar			K589 \$237 <input type="checkbox"/>	K679 \$276 <input type="checkbox"/>	K689 \$276 <input type="checkbox"/>	K689 \$276 <input type="checkbox"/>	N/A		
Rifton accessories tote			K522 \$139 <input type="checkbox"/>	K522 \$139 <input type="checkbox"/>	K522 \$139 <input type="checkbox"/>	K522 \$139 <input type="checkbox"/>	K522 \$139 <input type="checkbox"/>		
Total			K610 \$ _____	K620 \$ _____	K630 \$ _____	K640 \$ _____	K650 \$ _____		