

Permobil

Panthera S3 Swing Series



permobil

Adapts to the active life you lead

The Panthera S3 Swing offers the same comfort and low transport weight as the Panthera S3 series, the S3 Swing provides the added flexibility of a swing-away foot rest for easier transfers. Also available in S3 Swing Short and Long option to suit smaller or taller builds.

- Key benefits:
- Made to maneuver
 - Low weight and robust
 - Adjusts for comfort
 - Integrated anti-tip device
 - Easy to maintain



Personalise your ride*

ANODIC BLACK

PEARL WHITE

GREY METALLIC

BLUE CLASSIC

PINE GREEN

Actual product colors may vary.

Specifications

	S3 Swing	S3 Swing Short
Seat width	36 - 39 - 42 - 45 - 50 cm	33 - 36 - 39 cm
Seat depth	40 cm	35 cm
Seat height	Front 45 cm (24" rear wheels) Back 43 cm (24" rear wheels)	Front 45 cm (24" rear wheels) Back 43 cm (24" rear wheels)
Seat angle	4°	4°
Total length	78 - 90 cm	88 cm
Total width	Seat width +21 cm	Seat width +21 cm
Front frame angle	105°	105°
Backrest height	20 - 45 cm	20 - 45 cm
Backrest angle	13° - (-8°)	13° - (-8°)
Footplate distance	38 - 46 cm	38 - 46 cm
Camber	2.2°, 4° (seat width 36 - 45 cm) 8° (seat width 36 - 39 cm)	2.2°, 4°, 6 (seat width 33 cm) 8°
Weight without rear wheels	4.1 - 4.9 kg	4.0 - 4.1 kg
Wheelchair weight	9.0 - 10.0 kg	8.8 -9.0 kg
Maximum load	100 kg, 150 kg (seat width 45 - 50 cm)	100 kg
Material (chassis)	Chrome molybdenum tubing	Chrome molybdenum tubing
Material (upholstery)	Polyurethane covered polyester	Polyurethane covered polyester

S3 Swing Low 0°	S3 Swing Long
36 - 39 - 42 - 45 - 50 cm	36 - 39 - 42 - 45 - 50 cm
35 - 37.5 - 40 cm	37.5 - 40 - 42.5 - 45 cm
Front 39 cm (24" rear wheels) Back 39 cm (24" rear wheels)	Front 45 cm (24" rear wheels) Back 43 cm (24" rear wheels)
4°	4°
78 - 90cm	83 - 95 cm
Seat width +21 cm	Seat width +21 cm
105°	105°
20 - 45 cm	20 - 45 cm
13° - (-8°)	13° - (-8°)
38 - 46 cm	38 - 46cm
2.2° 4° (seat width 33-45cm) 8° (seat width 33-39cm)	2.2° 4° (seat width 33-45cm) 8° (seat width 33-39cm)
4.6 - 5.4 kg	4.9 - 5.7 kg
9.4 - 10.2 kg	9.7 - 10.5 kg
100 kg, 150 kg (seat width 45- 50 cm)	100 kg, 150 kg (seat width 45- 50 cm)
Chrome molybdenum tubing	Chrome molybdenum tubing
Polyurethane covered polyester	Polyurethane covered polyester

Innovating for individuals