



# PDG Eclipse - AADL W320



Dealer Name\*: \_\_\_\_\_ Account #\*: \_\_\_\_\_  
 City, State, Zip: \_\_\_\_\_  
 Phone: \_\_\_\_\_ Email: \_\_\_\_\_ ATP: \_\_\_\_\_  
 Rep Name: \_\_\_\_\_ P.O. Number: \_\_\_\_\_ Quote Number: \_\_\_\_\_  
 Client Reference\*: \_\_\_\_\_ Client Weight: \_\_\_\_\_  
 Special considerations that need to be addressed here (e.g., diagnosis): \_\_\_\_\_ \*Required field

## October 2024, Revised January 2026. Prices are in Canadian Dollars (CAD)

### Wheelchair & Frame

Eclipse Bariatric Manual Wheelchair W320 \$ 3,100.00

#### Frame Type

600lbs Weight Capacity No charge  
 1000lbs Weight Capacity *Not compatible with Transit Tie Down System* \$ 300.00

#### Frame Options

**Transit Tie Down System (TTDS)** *TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. Tested to max. weight capacity of 600lbs. Not compatible with O2 holder or IV pole.* \$ 395.00

### Seat Measurements

#### Seat Width

<b>Select Width:</b>	20" Seat Width	22" Seat Width	24" Seat Width	26" Seat Width	\$ 675.00
	28" Seat Width	30" Seat Width			\$ 875.00
	32" Seat Width	36" Seat Width			\$ 1,100.00
	40" Seat Width	44" Seat Width	48" Seat Width		\$ 1,300.00

#### Seat Depth (Adjustable)

<b>Select Depth:</b>	16" Seat Depth	17" Seat Depth	18" Seat Depth	19" Seat Depth	20" Seat Depth	No charge
	21" Seat Depth	22" Seat depth				
	23" Seat Depth	24" Seat Depth				\$ 400.00
	<i>24" seat depth not available with TTDS.</i>					

### Frame Finishes

Select Color: High Gloss Black Hyper Blue

## Armrests

### Armrest Type

No Armrests	No charge
Dual Post, Height Adjustable Armrests	\$ 295.00
10-14" Range    12"-16" Range	
Dual Post, Fixed Height Armrests (12")	No charge

### Armpads

Regular Armpads	No charge
Desk Length 2" x 10"    Full Length 2" x 14"	
Comfort Armpads with Interlocking Hardware <i>May interfere with wheels in some tilt configurations.</i>	A1 \$ 200.00
Desk Length 3.75" x 10"    Full Length 3.75" x 14"	
Gel Armpads <i>May interfere with wheels in some tilt configurations.</i>	A1 \$ 135.00
Narrow 2" x 12"    Wide 3.5" x 12"	

### Armrest Options

Forward Guards	\$ 200.00
Armrest Extensions - 3' Forward <i>Not compatible with desk length armpads.</i>	\$ 75.00

## Seat Height

**Seat Height** | *Front and rear height can be maximum 1" apart.*

Front Seat Height:	14" STF Front	15" STF Front	16" STF Front	17" STF Front	18" STF Front	
	19" STF Front	20" STF Front				\$ 200.00
Rear Seat Height:	14" STF Rear	15" STF Rear	16" STF Rear	17" STF Rear	18" STF Rear	
	19" STF Rear	20" STF Rear				\$ 200.00

## Seat Pan

Solid Aluminium Seat Pan	\$ 165.00
--------------------------	-----------

## Front Casters

5" Casters with Heavy Duty Steel Fork	No charge
---------------------------------------	-----------

## Rear Wheels

24" Heavy Duty Spoke Wheels	No charge
-----------------------------	-----------

Rear Wheel Position:	Position 1 (Rear)	Position 2	Position 3	Position 4 (Front)	
	<i>Positions 3 and 4 not compatible with TTDS. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.</i>				

### Axles

Quick Release Rear Axles	No charge
--------------------------	-----------

### Rear Tires

Urethane Tires	No charge
----------------	-----------

### Handrims

No Handrims	No charge
Standard Aluminium Handrims	No charge
Plastic Coated Handrims	\$ 175.00

## Wheel Locks

User Wheel Locks	No charge
Push to Lock    Pull to Lock	
Wheel Lock Extensions	No charge

## Anti-Tippers

Rear Anti-Tippers *Required safety option.*

\$ 75.00

## Backrest

**Back Post Style:** Stroller Handle Back Cane Handle Back  
*Cane handle not compatible with headrest mounting plate.*

**Back Post Height:** Short (22") Tall (25")  
*Measured from seat pan to handle.*

**Seat to Back Angle:** Set at 85° Set at 90° Set at 95° Set at 100°

## Backrest Options

**Endomorph Back** *Backrest will be narrower than seat by total of 2x offset selected. Minimum backrest width is 12".*

\$ 395.00

1" Offset Per Side 2" Offset Per Side

**Removable Backrest**

\$ 0.00

## Back Supports

### Back Upholstery

**No Back Upholstery**

No charge

**Tension Adjustable Back Upholstery**

\$ 375.00

Regular (18"H) Tall (22"H)

### Position of Back Upholstery

Uph Type	22" Back Post						25" Back Post								
	17"	18"	19"	20"	21"	22"	17"	18"	19"	20"	21"	22"	23"	24"	25"
Regular															
Tall															
	Upholstery Height														

## Recline

**Reclining Backrest (Manual)** 30° maximum recline. Not compatible with cane handle back or backrest extension. Reclining backrest is 1" narrower than seat width. Back frame is not removable. Finished back post height is 3" taller than selected back post height. Not recommended for user weight above 450lbs. Maximum seat depth is 23".

\$ 895.00

## Front Rigging

**No Front Rigging**

No charge

**Center Mount Front Rigging Bracket**

\$ 300.00

**Easy-Out Front Rigging** | Effortless insertion and removal. No hanger brackets remain on the frame.

**Fixed 60°**

\$ 95.00

**Universal Adjustable Angle (0°-60°)** Select Angle. Angle adjustment by loosening one bolt. Tool-free leg length adjustment.

\$ 295.00

0° 10° 20° 30° 40° 50° 60°

**AA High Mount Option** For shorter lower leg lengths of 8"-12.5".

\$ 295.00

**Swing Away Front Rigging** | Legacy design. Not compatible with Easy-Out.

**Swing Away 60°**

No charge

**Swing Away Elevating Legrests (ELR)** Includes calf pads.

W106

\$ 375.00

**1" Offset Hangers for Swing Away Front Rigging**

\$ 75.00

## Footplates

**Aluminium Angle Adjustable Footplates**

W107

\$ 165.00

**Extension Tubes** | No selection required for Easy-Out AA High-Mount option.

**Standard Extension Tubes, Short**

No charge

**Standard Extension Tubes, Long**

\$ 75.00

**Footplate Height:** Specify:

Ext	Short	8" - 13"	11" - 16"	12" - 16"
	Long	N/A	14" - 22"	14" - 19"
		High Mount Easy-Out	Standard Easy-Out	60°, 70°, ELR Swing-Away
		Legrest		

## Front Rigging Accessories

<b>Oversized 9"x10" Calf Pad</b> <i>Not compatible with seat widths less than 25"</i>	Basic Seating	\$ 695.00
<b>Calf Pads</b> <i>Not compatible with Easy-Out front rigging.</i>	Basic Seating	\$ 135.00
<b>Gel Knee Pads, Easy-Out Angle Adjustable Front Rigging</b>	Basic Seating	\$ 200.00
<b>Gel Knee Pads, Swing Away ELR Front Rigging</b>	Basic Seating	\$ 100.00
<b>Length Adjustable Heel Loops</b>		No charge
<b>Standard Calf Strap</b>	Basic Seating	\$ 50.00
<b>Bodypoint Padded Calf Strap</b> <i>Available up to a maximum size of 32"</i>	Basic Seating	\$ 0.00

## Miscellaneous Options

<b>Q'Straint Lap Belt</b> <i>Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" belt length.</i>	Basic Seating	\$ 245.00
<b>IV Pole</b> <i>Not compatible with TTDS.</i>	A1	\$ 295.00
<b>O2 Holder</b> <i>Not compatible with TTDS.</i>	W110	\$ 395.00
<b>Medical Necessities Backpack</b> <i>Black with Permobil logo.</i>	A1	\$ 95.00
<b>Adhesive Back Hook Strip</b> <i>For cushion mounting, 2"x30".</i>	A1	\$ 10.00

## General

Permobil manufactures PDG wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

## Notes