

Order form

BodiLink[®] Lateral Pelvic/Thigh Support

Pricing in Canadian dollars



Lateral Pelvic Thigh Support with TT Power Mount Removable Hardware



Removable Solid Seat Pan

Company Name/ACCT #: _____
 P.O. Number: _____
 Requested By: _____
 Phone: _____ Fax: _____
 Email: _____
 Ship To: _____
 Client Reference: _____

For best results, do not fill out in your browser. Interactive form should be completed using Adobe Reader after saving to your local drive. Then email or print and fax to Customer Support at sales.ca@permobil.com or 1.877.636.8944.

Please select sizes in the appropriate fabric & style charts below. Contact Customer Support for custom sizes that are not listed below at 800.265.9830.

1 Support pads

1.1 Premium pad shape & size

MSRP \$150.00

- You may make one selection per column (Left, Extra Left, Right, and Extra Right) in section 1.1 for Premium Pads, or 1.3 for Basic Pads.
- Make a pad cover fabric selection in 1.2 or 1.4.

*Length (L) refers to the actual size dimension of the support from top to bottom edge. Depth (D) refers to the actual size dimension from anterior to posterior edge.

Style	Size	Left		Extra left		Right		Extra right	
		Code	Code	Code	Code	Code	Code	Code	Code
	3.5" L x 4" D	BL-LPTSP2Z-3L4D-LH	BL-LPTSP2Z-3L4D-LH	BL-LPTSP2Z-3L4D-RH	BL-LPTSP2Z-3L4D-RH				
	3.5" L x 8" D	BL-LPTSP2Z-3L8D-LH	BL-LPTSP2Z-3L8D-LH	BL-LPTSP2Z-3L8D-RH	BL-LPTSP2Z-3L8D-RH				
	3.5" L x 12" D	BL-LPTSP2Z-3L12D-LH	BL-LPTSP2Z-3L12D-LH	BL-LPTSP2Z-3L12D-RH	BL-LPTSP2Z-3L12D-RH				
	5.5" L x 6" D	BL-LPTSP2Z-5L6D-LH	BL-LPTSP2Z-5L6D-LH	BL-LPTSP2Z-5L6D-RH	BL-LPTSP2Z-5L6D-RH				

1.2 Premium pad cover with foam insert

Cover options	MSRP	Left		Extra left		Right		Extra right	
		Code	Code	Code	Code	Code	Code	Code	Code
COMFORT-TEK	\$0.00/ea.	LPTS-P2-COMFORT-FM	LPTS-P2-COMFORT-FM	LPTS-P2-COMFORT-FM	LPTS-P2-COMFORT-FM	LPTS-P2-COMFORT-FM	LPTS-P2-COMFORT-FM	LPTS-P2-COMFORT-FM	LPTS-P2-COMFORT-FM
STRETCH-AIR	\$0.00/ea.	LPTS-P2-STRETCH-FM	LPTS-P2-STRETCH-FM	LPTS-P2-STRETCH-FM	LPTS-P2-STRETCH-FM	LPTS-P2-STRETCH-FM	LPTS-P2-STRETCH-FM	LPTS-P2-STRETCH-FM	LPTS-P2-STRETCH-FM
GLIDEWEAR™	\$36.00/ea.	LPTS-P2-GLIDE-FM	LPTS-P2-GLIDE-FM	LPTS-P2-GLIDE-FM	LPTS-P2-GLIDE-FM	LPTS-P2-GLIDE-FM	LPTS-P2-GLIDE-FM	LPTS-P2-GLIDE-FM	LPTS-P2-GLIDE-FM

COMFORT-TEK
 For fluid protection & an easily cleaned surface.


STRETCH-AIR
 For comfort & heat dissipation.

GlideWear
 For skin protection & shear reduction. Not fluid resistant.

1.3 Basic pad

MSRP \$139.00/ea.

*Length (L) refers to the actual size dimension of the support from top to bottom edge. Depth (D) refers to the actual size dimension from anterior to posterior edge.

Style	Size	Left	Extra left	Right	Extra right
	3"L x 4"D	BL-LPTSP1Z-3L4D-LH	BL-LPTSP1Z-3L4D-LH	BL-LPTSP1Z-3L4D-RH	BL-LPTSP1Z-3L4D-RH
	3"L x 5"D	BL-LPTSP1Z-3L5D-LH	BL-LPTSP1Z-3L5D-LH	BL-LPTSP1Z-3L5D-RH	BL-LPTSP1Z-3L5D-RH
	4"L x 4"D	BL-LPTSP1Z-4L4D-LH	BL-LPTSP1Z-4L4D-LH	BL-LPTSP1Z-4L4D-RH	BL-LPTSP1Z-4L4D-RH
	4"L x 6"D	BL-LPTSP1Z-4L6D-LH	BL-LPTSP1Z-4L6D-LH	BL-LPTSP1Z-4L6D-RH	BL-LPTSP1Z-4L6D-RH
	4"L x 8"D	BL-LPTSP1Z-4L8D-LH	BL-LPTSP1Z-4L8D-LH	BL-LPTSP1Z-4L8D-RH	BL-LPTSP1Z-4L8D-RH
	4"L x 10"D	BL-LPTSP1Z-4L10D-LH	BL-LPTSP1Z-4L10D-LH	BL-LPTSP1Z-4L10D-RH	BL-LPTSP1Z-4L10D-RH
	4"L x 12" D	BL-LPTSP1Z-4L12D-LH	BL-LPTSP1Z-4L12D-LH	BL-LPTSP1Z-4L12D-RH	BL-LPTSP1Z-4L12D-RH
	4"L x 14" D	BL-LPTSP1Z-4L14D-LH	BL-LPTSP1Z-4L14D-LH	BL-LPTSP1Z-4L14D-RH	BL-LPTSP1Z-4L14D-RH
	4"L x 16" D	BL-LPTSP1Z-4L16D-LH	BL-LPTSP1Z-4L16D-LH	BL-LPTSP1Z-4L16D-RH	BL-LPTSP1Z-4L16D-RH
5"L x 7" D	BL-LPTSP1Z-5L7D-LH	BL-LPTSP1Z-5L7D-LH	BL-LPTSP1Z-5L7D-RH	BL-LPTSP1Z-5L7D-RH	

1.4 Basic pad cover

COMFORT-TEK
For fluid protection & an easily cleaned surface.

STRETCH-AIR
For comfort & heat dissipation.

GlideWear For skin protection & shear reduction. **Not fluid resistant.**



Cover options	MSRP	Left	Extra left	Right	Extra right
COMFORT-TEK	\$0.00/ea.	LPTS-P1-COMFORT	LPTS-P1-COMFORT	LPTS-P1-COMFORT	LPTS-P1-COMFORT
STRETCH-AIR	\$0.00/ea.	LPTS-P1-STRETCH	LPTS-P1-STRETCH	LPTS-P1-STRETCH	LPTS-P1-STRETCH
GLIDEWEAR	\$36.00/ea.	LPTS-P1-GLIDE	LPTS-P1-GLIDE	LPTS-P1-GLIDE	LPTS-P1-GLIDE

2 Mounting hardware

2.1 BodiLink® Lateral Pelvic/Thigh Support hardware

- You may make one selection per column (Left, Extra Left, Right, and Extra Right) in section 2.1 from group A, B, -OR- C.
- TT long extension arms will only work with TT style hardware.
- GT hardware extra links will only work with GT style hardware.

A. TT hardware, slot mount




Hardware size	Max cushion thickness clearance	Recommended cushion thickness
Small	3.0"	2.0"- 3.0"
Medium	4.5"	3.0"- 4.5"
Large	6.5"	4.0"- 6.5"

	Size	MSRP	Left	Extra left	Right	Extra right
Fixed	Small	\$218.00/ea.	BL-LPTS-TT1FXSL1-LH	BL-LPTS-TT1FXSL1-LH	BL-LPTS-TT1FXSL1-RH	BL-LPTS-TT1FXSL1-RH
	Medium	\$218.00/ea.	BL-LPTS-TT1FXSL2-LH	BL-LPTS-TT1FXSL2-LH	BL-LPTS-TT1FXSL2-RH	BL-LPTS-TT1FXSL2-RH
	Large	\$218.00/ea.	BL-LPTS-TT1FXSL3-LH	BL-LPTS-TT1FXSL3-LH	BL-LPTS-TT1FXSL3-RH	BL-LPTS-TT1FXSL3-RH
Or						
	Small	\$311.00/ea.	BL-LPTS-TT1RMSL1-LH	BL-LPTS-TT1RMSL1-LH	BL-LPTS-TT1RMSL1-RH	BL-LPTS-TT1RMSL1-RH
	Medium	\$311.00/ea.	BL-LPTS-TT1RMSL2-LH	BL-LPTS-TT1RMSL2-LH	BL-LPTS-TT1RMSL2-RH	BL-LPTS-TT1RMSL2-RH
	Large	\$311.00/ea.	BL-LPTS-TT1RMSL3-LH	BL-LPTS-TT1RMSL3-LH	BL-LPTS-TT1RMSL3-RH	BL-LPTS-TT1RMSL3-RH

B. TT hardware, power mount


TT hardware power mount - Fixed



TT hardware power mount - Removable



Hardware size VS. cushion thickness



Hardware size	Max cushion thickness clearance	Recommended cushion thickness
Small	3.5"	2.0" - 3.5"
Medium	5.5"	3.0" - 5.5"
Large	7.5"	4.0" - 7.5"

	Size	MSRP	Left		Extra left		Right		Extra right	
	Fixed	Small	\$282.00/ea.	BL-LPTS-TT1FXPWL1-LH	BL-LPTS-TT1FXPWL1-LH	BL-LPTS-TT1FXPWL1-RH	BL-LPTS-TT1FXPWL1-RH	BL-LPTS-TT1FXPWL2-LH	BL-LPTS-TT1FXPWL2-LH	BL-LPTS-TT1FXPWL2-RH
	Medium	\$282.00/ea.	BL-LPTS-TT1FXPWL2-LH	BL-LPTS-TT1FXPWL2-LH	BL-LPTS-TT1FXPWL2-RH	BL-LPTS-TT1FXPWL2-RH	BL-LPTS-TT1FXPWL3-LH	BL-LPTS-TT1FXPWL3-LH	BL-LPTS-TT1FXPWL3-RH	BL-LPTS-TT1FXPWL3-RH
	Large	\$282.00/ea.	BL-LPTS-TT1FXPWL3-LH	BL-LPTS-TT1FXPWL3-LH	BL-LPTS-TT1FXPWL3-RH	BL-LPTS-TT1FXPWL3-RH				
Or										
Removable	Small	\$368.00/ea.	BL-LPTS-TT1RMPWL1-LH	BL-LPTS-TT1RMPWL1-LH	BL-LPTS-TT1RMPWL1-RH	BL-LPTS-TT1RMPWL1-RH	BL-LPTS-TT1RMPWL2-LH	BL-LPTS-TT1RMPWL2-LH	BL-LPTS-TT1RMPWL2-RH	BL-LPTS-TT1RMPWL2-RH
	Medium	\$368.00/ea.	BL-LPTS-TT1RMPWL2-LH	BL-LPTS-TT1RMPWL2-LH	BL-LPTS-TT1RMPWL2-RH	BL-LPTS-TT1RMPWL2-RH	BL-LPTS-TT1RMPWL3-LH	BL-LPTS-TT1RMPWL3-LH	BL-LPTS-TT1RMPWL3-RH	BL-LPTS-TT1RMPWL3-RH
	Large	\$368.00/ea.	BL-LPTS-TT1RMPWL3-LH	BL-LPTS-TT1RMPWL3-LH	BL-LPTS-TT1RMPWL3-RH	BL-LPTS-TT1RMPWL3-RH				


TT hardware, power mount option

Note: If choosing the Power Mount Hardware, choose the appropriate power mount type below.

Power mount	MSRP	Part number	
QUANTUM	add \$0.00	LPTS-TT1-PW1	
PERMOBIL	add \$0.00	LPTS-TT1-PW2	
QUICKIE	add \$0.00	LPTS-TT1-PW3	
ROVI	add \$0.00	LPTS-TT1-PW4	
AVID REHAB	add \$0.00	LPTS-TT1-PW5	

TT long extension arm 4"

MSRP \$55.00/arm




Left		Extra left		Right		Extra right	
BL-LPTS-TT1LARM4		BL-LPTS-TT1LARM4		BL-LPTS-TT1LARM4		BL-LPTS-TT1LARM4	


QUESTIONS? PLEASE CONTACT CUSTOMER SUPPORT FOR ASSISTANCE. | 800.265.9830

C. GT hardware, only slot mount


TT hardware slot mount - Fixed



TT hardware slot mount - Removable



No medial



1" Medial




Hardware size	Max cushion thickness clearance	Recommended cushion thickness
No medial	*5.0"	2.5" - 5.0"
1" Medial	*4.0"	2.25" - 4.0"

*If you need to increase the height or medial reach, add additional links in section 2.2.

Size	MSRP	Left		Extra left		Right		Extra right	
Fixed	\$218.00	BL-LPTS-GT2FXS-LH	BL-LPTS-GT2FXS-LH	BL-LPTS-GT2FXS-LH	BL-LPTS-GT2FXS-LH	BL-LPTS-GT2FXS-RH	BL-LPTS-GT2FXS-RH	BL-LPTS-GT2FXS-RH	BL-LPTS-GT2FXS-RH
Removable	\$260.00	BL-LPTS-GT2RMS-LH	BL-LPTS-GT2RMS-LH	BL-LPTS-GT2RMS-LH	BL-LPTS-GT2RMS-LH	BL-LPTS-GT2RMS-RH	BL-LPTS-GT2RMS-RH	BL-LPTS-GT2RMS-RH	BL-LPTS-GT2RMS-RH

2.2 Additional options



GT additional hardware links

MSRP \$58.00/link
(2 Additional 1.25" Links Allowed Per Hardware Ordered)

Left		Extra left		Right		Extra right	
BL-LPTS-GT2LINK		BL-LPTS-GT2LINK		BL-LPTS-GT2LINK		BL-LPTS-GT2LINK	
BL-LPTS-GT2LINK		BL-LPTS-GT2LINK		BL-LPTS-GT2LINK		BL-LPTS-GT2LINK	

3 Solid seat pans


3.1 Removable solid seat pan & hardware kit

MSRP \$690.00

The solid seat pan kit includes a slotted aluminum pan and attaching hardware to accommodate 7/8" or 1" tubing. The kit also comes with two different cross bars to accommodate various wheelchair frame types. This will fit both folding and non-folding seat rails. Lateral thigh support hardware can be easily attached.

Please fill in the below with dimensions. Example: RRSP-1-1416 for a 14"W x 16"D solid seat pan.

RSSP-1- _____



Size range:
Width: 10" - 28"
Depth: 10" - 22"