

Permobil

# Progeo Tekna Tilt 2.0



Aluminium



Foldable



Adjustable



Tilt

permobil

# A lightweight and compact folding wheelchair

A versatile tilting wheelchair with a range of options for seating and positioning, plus compact folding for transportation and storage

For wheelchair users with different seating and positioning needs. The Tekna Tilt 2.0 can be fitted with small or large rear wheels, and also with a drum brake.



## Progeo Tekna Tilt 2.0

Seat width	33 - 36 - 39 - 42 - 45 - 48 cm
Seat depth	35 - 37.5 - 40 - 42.5 - 45 - 47.5 - 50 cm
Seat height	42 - 48 cm (adjustable)
Footsupport distance	From 27 to 46 (adjustable)
Backrest height	From 35 to 55 cm (adjustable)
Front frame angle	Standard 95° - 110° - 120° Elevating legsupport 100° - 180° (adjustable)
Total length	95 cm approx.
Total width	Seat width + 21 cm
Total width folded	33 cm approx.
Backrest angle	From 78° to 90° (adjustable)
Tilt-in-space	0° - 7.5° - 15° - 22.5° - 30° - 37.5° - 45° (adjustable)
Weight without rear wheels	Starting from 15 kg approx.
Wheelchair weight	Starting from 22 kg approx.
Maximum load	125 kg

Personalise your ride\*

BLACK

BRIGHT  
LIME

BRIGHT  
FUCHSIA

BRIGHT  
BLUE

BRIGHT  
SILVER

Permobil reserves the right to change specifications and design without prior notice. \*More color options available.

Actual product colors may vary.

permobil

permobil.com



GR3354-1UK230303