

Adjusted by:

Tel:

Save the form to the desktop and open it with Adobe Reader.

E-mail:

Date:

The price changes depending on the choices made in the form £:

Standard configuration includes; chassis, seat, backrest lock, colour , backrest, carbon rear axle, castor and fork, foot rest, push to lock Bambino 3 brake, Panthera crosslaced steel spoke rear wheel, titanium pushrims and Schwalbe right Run tires, anti-tip device

Choice required:

Max. user weight : Seat width **27-33 cm 60 kg**, Click on "**blue underlined text**" for more info or image

Optional **Upcharge /Add-price** **Marked yellow**

Original design

Alternative options

Other options

1) Chassis

Bambino 3

Seat width 24-33cm, Seat depth 25-30 cm

Bambino 3 Short

Seat width 24-27cm, Seat depth 18-25 cm

Clear

2) Chassis colour

Blue - Racing Blue

Optional colour - A delivery time of 20 days for optional colour. Fill in art.no. in the square

Black - Anodic Black

Clear

3) Seat width

24 cm

27 cm

30 cm

33 cm

Clear

4) Backrest adjustable height

20-28 cm blue

523BB17

27-35 cm blue

523BB18

20-28 cm black

523BB17BL

27-35 cm black

523BB18BL

Clear

5) Rear axle

Carbon 6dgr

33030BB

Carbon 6dgr +10 mm

33031BB

Other rear axle or seat height

Increase the seat height

Clear

With 20" rearwheels +120 mm castors

6) Castor and fork

art.no.

90 mm

Other castor fork

120 mm

Standard with 20" rearwheels

Excel guide art.no:

Clear

7) Footrest

Footrest B3

42341BB

Clear

8) Brake

Brake B3

4630200

Clear

Brake blade B3 PU tyre SET

4630234

Assistant brake

art.no.

All prices for rearwheels are for pair

9) Rear wheel - without pushrim

[Rear wheel sizes](#)

[22" standard](#)

20" [standard](#)

[One arm drive](#)

Clear

24" [standard](#)

Art.nr.

[24" Spinergy X, black spokes](#)

[24" Spinergy X, white spokes](#)

10) Tyre - choose tyre

All prices for tyres are for pair

[Schwalbe Right Run](#) 20", 22", 24"

[Rugged 24"](#)

[PU 20", 22", 24" \(massive\)](#)

[Schwalbe Marathon Plus, 22", 24"](#)

reTyre does not fit PU tyres

[Schwalbe One 24"](#)

[reTyre skins 24"](#)

Clear

11) Pushrim

All prices for pushrims are for pair

[Aluminum](#)

[Paragrip](#)

[Maxgrip 20", 22", 24"](#)

1902430 (x2 / wheelchair)

19051DDBT x 2

1903070 x 2 (70=20", 80=22", 60=24"

(contains latex)

Clear

[Tetragrip](#)

[Titanium](#)

19052DDBT x 2

190DD11 (x2 / wheelchair)

DD = Wheel diameter

Accessories

Side guard

[Side guard B3](#)

4235000

Clear

Arm support

[Arm support B3](#) with side guard

Adjustable in 3 heights

Seat width adaptor

[Seat width adaptor](#) - pair

[Seat width adaptor](#) - left

[Seat width adaptor](#) - right

Limits the back rest angle forward to 16 mm.
Narrows the seat width 7,5mm on each side

Clear

Lateral support

[Lateral support B3 adjustable](#)

4236000 / pcs

Left Right

Std Std

Narrow-9mm Narrow-9mm

Push bar

[Push bar B3](#)

4233000

Height, angle adjustable
with Quick release

Head support

[Head support B3 Blue](#)

4710100

[Head support B3 Black](#)

4710100BL

Clear

Low position.

High position

Cushion

[2,5 cm](#)

572BB34

[5 cm](#)

575BB34

Clear

Wedge cushion

[Wedge cushion](#)

57710BB

Not for width 24 cm

Attachment pelvic belt

[Attachment plate pelvic belt cpl](#)

9900247 backrest 2 point attachment

[Attachment clamp pelvic belt pair](#)

9900248 2 point attachment

[Attachment pelvic belt 4-p plate + clamp](#)

9900250 4 point attachment

[Attachment pelvic belt 4-p with clamps](#)

9900248 x 2 4 point attachment

Spoke protectors

[20"](#) [22"](#) [24"](#) [LOGO](#)

4500021 4500100 4500000

[22"](#) [24"](#) [NO LOGO](#)

4500011 4500111

Clear

Spoke protectors with theme

[20"](#) [22"](#) [24"](#) [Splash](#)

4500320 4500322 4500324

[20"](#) [22"](#) [24"](#) [Car rim](#)

4500120 4500122 4500124

[20"](#) [22"](#) [24"](#) [Wild animals](#)

4500420 4500422 4500424

[20"](#) [22"](#) [24"](#) [Smiley](#)

4500220 4500222 4500224

Clear

SPLASH CAR RIM WILD ANIMALS SMILEY

Other

[Tube protection, front *](#)

9900251 (2 pcs)

[Tube protection, rear](#)

9900251 (1 pcs)

[Tray-table](#)

6810000

[Magnetic holder](#)

6810010

[Backrest wedge](#)

4710012

[Torso strap](#)

4380238

[Freewheel *](#)

Freewheel-KIT assembly parts

42342BB

* Tube protection front;
Can be too long if footrest is in its
high position. Fits when rear axle is
at std balance and footrest at low
pos. 16 cm or lower, high pos. 28 cm
or lower, from floor to top of
footplate.

* Choose art.no. for Freewheel according to the footrest
height of the wheelchair. Read more about Freewheel on
the [homepage](#)

Set-up adjustments

Back rest angle

Note measure (A) for desired back rest angle.
14-28 mm. Standard = 18 mm

A= _____ mm

Rear axle adjustment

Note measure (A) for desired balance
Standard = 77 mm.

A= _____ mm

Footrest position

Note the height measure from the floor, (A) or (B)
for the footrest.

[Bambino footplate with Freewheel - info](#)

A= _____ mm B= _____ mm

At 22" rear wheels

Measure high pos. 210-310 mm

Measure low pos. 150-220 mm

Note depth position

Note height position

Std position

High position

Moved forward 30 mm position

Low position

Extra