

Aluminum dual-tube rigid wheelchair

TiLite[®] Aero T

Adjustable seating angle
to maximize comfort

Find your perfect
ride balance with
our adjustable
canter of gravity

Custom designed
and built to your
specifications, with
available ergo seating

Dual-tube geometry
provides better load
distribution and less
frame flex

Adjustable caster
fork angle to ensure
a precise ride




permobil

Tattoo options




Frame colors


High gloss paint




Candy paint



Metallic paint




Matte paint



Pearlescent paint




Matte textured paint



Component colors

Anodized



Specifications

Frame features	Aluminum Dual-Tube Rigid, TiFit, full adjustability
Total chair width*	23"
Seat width	12" – 20"
Seat depth	12" – 20"
Front seat height	15" – 21"
Rear seat height	13" – 21"
Front frame angle	70°, 75°, 80°, 85°, or 90°
Camber	0°, 2°, 4°, or 6°
Center of gravity adjustment	5"

Front seat width (seat taper)	0", 1", 2", 3"
Folding seat back	Adj. height 9.5" – 22" Adj. angle 80° – 101°
Ergonomic seat	6", 7", 8"
Transit options	Q-Straint lock down approved
Transport weight*	12 lb
Weight limit	265 lb Standard 350 lb Heavy duty
HCPSC Code	K0005, E1235
Warranty	Five years on frame One year on TiLite components

* Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lb.