



PDG Bentley - AADL W322



Dealer Name*:

Account #*:

City, State, Zip:

Phone:

Email:

ATP:

Rep Name:

P.O. Number:

Quote Number:

Client Reference*:

Client Weight:

*Required field

Special considerations that need to be addressed here (e.g., diagnosis):

October 1, 2026. Prices are in Canadian Dollars (CAD)

Wheelchair & Frame

Bentley 20° Manual Tilt Wheelchair Also available: Bentley LT with lightweight aluminium frame. See separate order form.

W322

C\$ 2,955.00

Frame Type

- ☐ 250lbs Weight Capacity Available for 14" - 22" seat widths. No charge
- ☐ 300lbs Weight Capacity Available for 14" - 22" seat widths. Prior Approval C\$ 400.00
- ☐ 450lbs Weight Capacity Available for 20" - 32" seat widths. Prior Approval C\$ 850.00

User Weight | Required for setup of tilt uplift force.

- User Weight Range: ☐ Up to 75lbs ☐ 75lbs - 100lbs ☐ 100lbs - 115lbs ☐ 115lbs - 150lbs ☐ 150lbs - 200lbs ☐ 200lbs - 250lbs
- ☐ 250lbs - 300lbs ☐ 300lbs - 375lbs ☐ 375lbs - 450lbs

Frame Options

- ☐ Extended Base with Rear Floor Jacks Lifts rear wheels off the ground for increased stability. Limited ground clearance. Not compatible with attendant wheel lock, HD anti-tipper, IV pole, O2 holder. C\$ 1,450.00
- ☐ Transit Tie Down System (TTDS) Includes 4 factory installed transport brackets. TTDS has been tested according to ANSI/RESNA WC19. See owners manual for complete details. Not compatible with seat widths > 22". Tested to maximum weight capacity of 300lbs. No charge

Seat Measurements

Seat Width | Overall chair width is +9.5" with 12"/16" rear wheels, +10.5" with 20"/22"/24" rear wheels.

- Select Width: ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19" No charge
- ☐ 20" ☐ 21" ☐ 22" W111 C\$ 200.00
- ☐ 23" ☐ 24" ☐ 25" ☐ 26" ☐ 27" W111 C\$ 500.00
- ☐ 28" ☐ 29" ☐ 30" ☐ 31" ☐ 32" W111 C\$ 1,150.00

Seat widths 23"-32" require 450lb Frame. Extended delivery lead times may apply.

Seat Depth

- Select Depth: ☐ 16" ☐ 17" ☐ 18" ☐ 19" No charge
- ☐ 20" W111 C\$ 100.00
- ☐ 21" ☐ 22" C\$ 700.00

Frame Finishes

- Select Color: ☐ High Gloss Black ☐ Hyper Blue

Armrests

Armrest Type

- ☐ No Armrests No charge
- ☐ Single Post, Height Adjustable Armrests No charge
- ☐ 8"-11" Range ☐ 10-14" Range ☐ 12"-16" Range

- ☐ **Dual Post, Height Adjustable Armrests** C\$ 545.00
☐ 8"-11" Range ☐ 10-14" Range ☐ 12"-16" Range
- ☐ **Fixed Height Armrests** 12" height, full length armpads only. No charge
☐ Single Post ☐ Dual Post
- Armrest Position:** ☐ Rear ☐ Front
If left blank default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions PDG may not be able to accommodate stated preference.
- ☐ **Locking Flip-Up Armrests** Not recommended for users with excessive strength or tone. Not compatible with fold-down back, or clothing guards. Only with regular armpads. C\$ 545.00
☐ Set at 10" ☐ Set at 11" ☐ Set at 12" ☐ Set at 13" ☐ Set at 14"

Armpads

- ☐ **Regular Armpads** No charge
☐ Desk Length 2" x 10" ☐ Full Length 2" x 14"
- ☐ **Waterfall Armpads** C\$ 275.00
☐ Desk Length 2" x 10" ☐ Full Length 2" x 14"
- ☐ **Comfort Armpads with Interlocking Hardware** May interfere with wheels in some tilt configurations. C\$ 496.00
☐ Desk Length 3.75" x 10" ☐ Full Length 3.75" x 14"
- ☐ **Gel Armpads** May interfere with wheels in some tilt configurations. C\$ 305.00
☐ Narrow 2" x 12" ☐ Wide 3.5" x 12"
- ☐ **Comfort Standard Arm Trough with Interlocking Hardware, Right** Not compatible with lap tray or locking flip-up armrests. May interfere with wheels in some tilt configurations. C\$ 562.00
☐ 6" x 17" ☐ 6" x 22"
- ☐ **Comfort Standard Arm Trough with Interlocking Hardware, Left** Not compatible with lap tray or locking flip-up armrests. May interfere with wheels in some tilt configurations. C\$ 562.00
☐ 6" x 17" ☐ 6" x 22"

Armrest Wideners | Moves armrests outwards by up to 1" each side. Seat size does not change. May cause tilt interference.

- ☐ **Armrest Wideners, 1" Total Additional Width** C\$ 245.00
- ☐ **Armrest Wideners, 2" Total Additional Width** C\$ 395.00

Clothing Guards | Not compatible with flip-up armrests.

- ☐ **Clothing Guards, Single Post Armrests** Not compatible with 8"-11" height range. C\$ 125.00
- ☐ **Clothing Guards, Dual Post Armrests** Not compatible with 8"-11" height range. C\$ 350.00

Armrest Extensions | Not compatible with flip-up armrests or desk length armpads.

- ☐ **Armrest Extensions - 3" Forward** C\$ 215.00
- ☐ **Armrest Extensions - 3" Rearward** C\$ 215.00

Seat Height

Seat Height - 6" Casters

		Low					High					
RW	24"		☐	☐	☐			☐	☐	☐	☐	
	22"	☐	☐	☐			☐	☐	☐	☐		
	20"	☐	☐	☐			☐	☐	☐			
	16"		☐	☐	☐			☐	☐	☐	☐	
	12"	☐	☐	☐			☐	☐	☐			
		13"	14"	15"	16"	17"	15"	16"	17"	18"	19"	20"
Seat to Floor Height												

Seat Height - 8" Casters

RW	Low			High			
	24"	☐	☐	☐	☐	☐	☐
	16"	☐	☐	☐	☐	☐	☐
	15"	16"	17"	17"	18"	19"	20"
Seat to Floor Height							

Seat Pan

- ☐ **No Seat Pan** No charge
- ☐ **Solid Aluminium Seat Pan** No charge

Seat Cushions

Cushions over 20" wide subject to increased pricing.

☐ Comfort M2 Cushion Zero Elevation	C\$ 580.00
☐ Comfort-Tek Cover ☐ Stretch Air Cover	
☐ Comfort M2 Cushion Anti-Thrust	C\$ 580.00
☐ Comfort-Tek Cover ☐ Stretch Air Cover	
☐ Comfort Embrace Cushion Zero Elevation	C\$ 508.00
☐ Comfort-Tek Cover ☐ Stretch Air Cover	
☐ Comfort Embrace Cushion Anti-Thrust	C\$ 508.00
☐ Comfort-Tek Cover ☐ Stretch Air Cover	
☐ Comfort Saddle 7 Series Foam Cushion Zero Elevation Solace cover.	C\$ 818.00
☐ Comfort Saddle 7 Series Gel Cushion Zero Elevation Solace cover.	C\$ 1,133.00
☐ Comfort Support Pro 7 Series Foam Cushion Zero Elevation Solace cover.	C\$ 445.00
☐ Comfort Support Pro 7 Series Gel Cushion Zero Elevation Solace cover.	C\$ 693.00
☐ ROHO Quadtro Select High Profile Cushion	C\$ 777.00
☐ ROHO Sensor Ready High Profile Cushion	C\$ 637.00

Front Casters

Urethane Casters

☐ 6" x 1.25" Caster	No charge
☐ 8" x 1.25" Caster Only compatible with 16" and 24" rear wheels.	No charge
☐ 6" x 1.50" Caster	C\$ 230.00
☐ 8" x 2.00" Caster Only compatible with 16" and 24" rear wheels.	C\$ 305.00

Pneumatic Casters

☐ 6" x 1.25" Pneumatic Caster	C\$ 265.00
☐ 8" x 1.25" Pneumatic Caster Only compatible with 16" and 24" rear wheels.	C\$ 280.00
☐ 8" x 2.00" Pneumatic Caster Only compatible with 16" and 24" rear wheels.	C\$ 325.00

Soft-Roll Casters

☐ 6" x 1.5" Soft-Roll Caster Urethane tire and silver anodized aluminum hub.	C\$ 515.00
--	------------

Frog Legs Suspension Forks

☐ FrogLegs Suspension Forks Only compatible with 16" and 24" rear wheels. Requires 6" soft-roll caster.	C\$ 995.00
---	------------

Rear Wheels

☐ Composite Rear Wheel	No charge
☐ 12" ☐ 16" ☐ 20" ☐ 22" ☐ 24"	
☐ Spoke Rear Wheel	C\$ 215.00
☐ 20" ☐ 22" ☐ 24"	
☐ Heavy Duty Spoke Wheel	C\$ 625.00
☐ 24"	
Rear Wheel Position: ☐ Rear ☐ Front (Default)	

If left blank, default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

Axles

☐ Quick Release Rear Axles	No charge
☐ Fixed Rear Axles	No charge

Rear Tires | Pneumatic tires may result in a finished seat height 1/2" greater than stated.

☐ Urethane Tires	No charge
☐ Pneumatic Tires	C\$ 190.00
☐ Pneumatic Tires with Flat-Free Inserts Not available with spoke wheels.	C\$ 245.00

Handrims

- | | |
|--|------------|
| ☐ No Handrims | No charge |
| ☐ Standard Aluminium Handrims | No charge |
| ☐ Plastic Coated Handrims <i>Not compatible with drum brakes.</i> | No charge |
| ☐ Projection Handrims <i>Composite wheels only. 10 projections, 30°, 1-1.8" length.</i> | C\$ 845.00 |

Wheel Locks

- | | |
|--|--------------|
| ☐ User Wheel Locks | No charge |
| ☐ Push to Lock ☐ Pull to Lock | |
| ☐ Attendant Wheel Locks <i>Not compatible with wheel lock extensions.</i> | C\$ 230.00 |
| ☐ Push to Lock ☐ Pull to Lock | |
| ☐ User and Attendant Wheel Locks <i>Not compatible with drum brakes on 20" rear wheels.</i> | C\$ 290.00 |
| ☐ Push to Lock ☐ Pull to Lock | |
| ☐ Drum Brakes <i>Available with 12" composite or 20" spoke rear wheels only. Urethane tires only. Not compatible with plastic coated handrims or two trigger tilt system.</i> | C\$ 795.00 |
| ☐ User Activated Unilateral Wheel Lock <i>Pull to lock. Available up to 20" seat widths only. Not compatible with hill holders.</i> | C\$ 1,295.00 |
| ☐ Left Activation ☐ Right Activation | |

Wheel Lock Options

- | | |
|---|------------|
| ☐ Wheel Lock Extensions <i>Not compatible with attendant wheel lock.</i> | No charge |
| ☐ Hill Holders <i>Requires pneumatic rear tires with flat-free inserts. Not compatible with unilateral wheel lock.</i> | C\$ 775.00 |

Anti-Tippers

- | | |
|---|--------------|
| ☒ Rear Anti-Tippers <i>Required safety option.</i> | No charge |
| ☐ Front Anti-Tippers | C\$ 625.00 |
| ☐ Heavy Duty Rear Anti-Tippers <i>Provides additional rearward stability for users with excessive movement. Not compatible with attendant wheel locks, rear floor jacks, IV pole or O2 holder.</i> | C\$ 1,445.00 |

Backrest

- | | |
|--|--|
| Back Post Style: | ☐ Stroller Handle Back ☐ Cane Handle Back |
| <i>Seat widths greater than 22" require stroller handle back. Cane handle back may reduce tilt range. Cane handle back not compatible with power tilt with manual override, dynamic rocker back or endomorph back.</i> | |
| Back Post Height: | ☐ Short (22") ☐ Tall (25") |
| <i>Measured from seat pan to handle.</i> | |
| Folding Style: | ☐ Fold-Down ☐ Non-Folding |
| <i>Fold-down back not compatible with reclining backrest or dynamic rocker back and will be substituted with a non-folding back.</i> | |
| Seat to Back Angle: | ☐ Set at 90° ☐ Set at 98° ☐ Set at 106° ☐ Set at 114° ☐ Set at 122° |
| <i>Back angles >98° in combination with seat depths >19" may limit range of tilt.</i> | |

Backrest Options

- | | |
|--|--------------|
| ☐ Angle/Height Adjustable Stroller Handle: <i>Not compatible with cane handle back.</i> | No charge |
| ☐ Button on Outside ☐ Button on Inside | |
| ☐ Endomorph Back <i>Backrest will be narrower than seat by total of 2x offset selected. Minimum backrest width is 12". Not compatible with flip-up armrests, dynamic rocker back, reclining backrest or cane handle back. Requires straight back post.</i> | C\$ 695.00 |
| ☐ 1" Offset Per Side ☐ 1.5" Offset Per Side ☐ 2" Offset Per Side | |
| ☐ Spring Back Suspension <i>Not compatible with reclining back or dynamic rocker back.</i> | No charge |
| ☐ Seating Dynamics Dynamic Rocker Back Interface <i>Not compatible with fold-down back, spring back suspension, endomorph back, 30° recline or cane handle. Not available with 20", 21" and 22" seat depths. Finished back post height is 4" taller than selected back post height.</i> | C\$ 2,995.00 |

Back Supports

Comfort Back Supports | Backs over 20" wide or tall subject to increased pricing.

- | | |
|---|--------------|
| ☐ Comfort Unite 1 Back <i>Includes 2-point quick release Compass hardware. Select shell height and cover.</i> | C\$ 961.00 |
| ☐ 12" ☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22" | |
| ☐ Comfort-Tek Cover ☐ Diamond Spacer Cover | |
| ☐ Comfort Unite 1 Back - 5" Wing <i>Includes 2-point quick release Compass hardware. Select shell height and cover.</i> | C\$ 1,245.00 |
| ☐ 12" ☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22" | |
| ☐ Comfort-Tek Cover ☐ Diamond Spacer Cover | |

- ☐ **Comfort Unite 1 Back w/ROHO Insert - 5" Wing** Includes 2-point quick release Compass hardware. Select shell height and cover. C\$ 1,499.00
☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22"
☐ Comfort-Tek Cover ☐ Diamond Spacer Cover
- ☐ **Comfort Unite 1 Back - 7" Wing** Includes 2-point quick release Compass hardware. Select shell height and cover. C\$ 1,245.00
☐ 12" ☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22"
☐ Comfort-Tek Cover ☐ Diamond Spacer Cover
- ☐ **Comfort Unite 1 Back w/ROHO Insert - 7" Wing** Includes 2-point quick release Compass hardware. Select shell height and cover. C\$ 1,499.00
☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22"
☐ Comfort-Tek Cover ☐ Diamond Spacer Cover
- ☐ **Comfort Unite 2 Back** Includes 2-point quick release Compass hardware. Select shell height and cover. C\$ 1,106.00
☐ 12" ☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22"
☐ Comfort-Tek Cover ☐ Diamond Spacer Cover
- ☐ **Comfort Unite 2 Back w/ROHO Insert** Includes 2-point quick release Compass hardware. Select shell height and cover. C\$ 1,254.00
☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22"
☐ Comfort-Tek Cover ☐ Diamond Spacer Cover
- ☐ **Comfort Unite 3 Back** Includes 2-point quick release Compass hardware. Select shell height and cover. C\$ 1,287.00
☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22"
☐ Comfort-Tek Cover ☐ Diamond Spacer Cover
- ☐ **Comfort Unite 3 Back w/ROHO Insert** Includes 2-point quick release Compass hardware. Select shell height and cover. C\$ 1,499.00
☐ 16" ☐ 18" ☐ 20" ☐ 22"
☐ Comfort-Tek Cover ☐ Diamond Spacer Cover
- ☐ **Comfort Embrace Back** Includes 2-point quick release Compass hardware. Select shell height and cover. C\$ 1,254.00
☐ 10" ☐ 12" ☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22"
☐ Comfort-Tek Cover ☐ Diamond Spacer Cover
- ☐ **Comfort Relief Back** Includes Comfort-Tek cover and 2-point quick release Compass hardware. Select shell height. C\$ 1,329.00
☐ 10" ☐ 12" ☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22"

Positioning Accessories

Lateral Supports

- ☐ **BodiLink Lateral Supports** Select size, cover and side. C\$ 220.00
☐ 3" x 5.5" ☐ 4" x 6.5" ☐ 6" x 6.5"
☐ Comfort-Tek Cover ☐ Stretch Air Cover
☐ Left ☐ Right
☐ Pair C\$ 440.00
- ☐ **BodiLink Lateral Support Mounting Hardware**
☐ Individual C\$ 368.00
☐ Pair C\$ 736.00

Headrests

- ☐ **BodiLink Headrest Pad** Select size and cover. C\$ 260.00
☐ 10" Wide ☐ 12" Wide ☐ 14" Wide ☐ 18" Wide
☐ Comfort-Tek Cover ☐ Stretch Air Cover
- ☐ **BodiLink Headrest Hardware, Removable**
☐ Mid 6 ☐ Mid 8 C\$ 373.00
☐ Lat 6 ☐ Lat 8 C\$ 489.00
☐ Max 8 C\$ 545.00
- ☐ **BodiLink Headrest Hardware, Flip-Down**
☐ Mid 6 ☐ Mid 8 C\$ 419.00
☐ Lat 6 ☐ Lat 8 C\$ 530.00
☐ Max 8 C\$ 591.00

Recline

- ☐ **Reclining Backrest (Manual)** 30° maximum recline. Not compatible with cane handle back, spring back suspension, fold-down back, dynamic rocker back or endomorph back. Requires longer frame. Maximum seat depth is 19". Finished back post height is 3" taller than selected back post height. C\$ 1,445.00

Tilt Activation

Attendant Tilt Only | Chairs with 450lbs weight capacity require two tilt triggers for tilt activation.

- ☐ **Attendant Tilt Only, Right Back Post** No charge
- ☐ **Attendant Tilt Only, Left Back Post** No charge
- ☐ **Two Trigger Setup** A trigger mounted on each side of the stroller handle. Reduces grip force required to activate tilt. Not compatible with reclining backrest or drum brakes. No charge

User Tilt Only | Not available for 450lbs weight capacity. Full length armrests only.

- ☐ **User Tilt Only, Right Armrest** C\$ 250.00
- ☐ **User Tilt Only, Left Armrest** C\$ 250.00

Attendant & User Tilt | Not available for 450lbs weight capacity. Not compatible with tall (25") back posts or reclining backrest. Full length armrests only.

- ☐ **Attendant & User Tilt Operation** C\$ 395.00

Attendant Trigger Position: ☐ Right Back Post ☐ Left Back Post

User Trigger Position: ☐ Right Armrest ☐ Left Armrest

Power Tilt | 250lbs weight capacity. 20° tilt range. Includes 2 batteries for off-chair charging. Power tilt is compatible with either O2 holder or IV pole, not both. Not available with 21" and 22" seat depths. Requires stroller handle back.

- ☐ **Standard Power Tilt** C\$ 1,795.00
- ☐ **Power Tilt with Manual Tilt Override** Allows for instant manual release of tilt in case of power tilt failure. Not compatible with cane handle back. C\$ 2,095.00
- ☐ **Toggle Switch Control** No charge
- ☐ Installed on Right ☐ Installed on Left No charge
- ☐ **Button Switch Control** No charge
- ☐ Installed on Right ☐ Installed on Left No charge
- ☐ **Stealth Egg Switch Control** Not installed. No charge
- ☐ **Stealth Swing-Away Mounting Hardware** Only for Stealth egg switch. Not compatible with all headrest makes and models. C\$ 1,040.00

Front Rigging

- ☐ **No Front Rigging** No charge

Easy-Out Front Rigging | Effortless insertion and removal. No hanger brackets remain on the frame.

- ☐ **Fixed 70°** C\$ 95.00
- ☐ **Universal Adjustable Angle (0°-100°)** Select angle. Angle adjustment by loosening one bolt in 10° increments. Tool-free leg length adjustment. Front rigging may interfere with castors in the 90°-100° settings. C\$ 545.00
- ☐ 0° ☐ 10° ☐ 20° ☐ 30° ☐ 40° ☐ 50° ☐ 60° ☐ 70° ☐ 80° ☐ 90° ☐ 100°
- ☐ **High Mount Option for Easy-Out Adjustable Angle** For shorter lower leg lengths of 8"-12.5". C\$ 545.00

Swing Away Front Rigging | Legacy design. Not compatible with Easy-Out.

- ☐ **Swing Away 60°** No charge
- ☐ **Swing Away 70°** No charge
- ☐ **Swing Away Elevating Legrests (ELR)** Includes calf pads. W106 C\$ 375.00
- ☐ **1" Offset Hangers for Swing Away Front Rigging** C\$ 95.00

Plug In Front Rigging | Legacy design. Not compatible with Easy-Out.

- ☐ **Plug In 70°** No charge

Footplates

- ☐ **Composite Footplates** No charge
- ☐ **Aluminium Angle Adjustable Footplates** W107 C\$ 165.00
- ☐ **One Piece Fixed Footplate** Available on plug in front rigging only. Requires selection of composite or aluminum angle adjustable footplates. C\$ 625.00
- ☐ **One Piece Flip-Up Footplate** Requires selection of composite or aluminum angle adjustable footplates. C\$ 695.00
- ☐ **Interlocking Angle Adjustable Footplates** Available for seat widths up to 21" C\$ 395.00

Extension Tubes | Consider the height of your seat cushion.

- | | |
|--|------------|
| ☐ Standard Extension Tubes, Short | No charge |
| ☐ Standard Extension Tubes, Long | C\$ 75.00 |
| ☐ Heavy Duty Extension Tubes, Short | C\$ 250.00 |
| ☐ Heavy Duty Extension Tubes, Long | C\$ 250.00 |

Footplate Height: ☐ Specify:

Ext Tube	Std Short	8" - 13"	11" - 16"	12" - 16"	13" - 16"
	Std Long	N/A	14" - 22"	14" - 19"	14" - 17"
	HD Short	N/A	N/A	N/A	13" - 17"
	HD Long	N/A	N/A	N/A	15" - 19"
		High Mount Easy-Out	Standard Easy-Out	60°, 70°, ELR Swing-Away	Standard Plug-In
Legrest					

Front Rigging Accessories

- | | | |
|---|--|------------|
| ☐ Calf Pads, Easy-Out Front Rigging | Not compatible with gel padding. Includes reinforced premium calf pad mount. | C\$ 400.00 |
| ☐ Calf Pads, Swing Away Front Rigging | Not compatible with gel padding. | C\$ 305.00 |
| ☐ Gel Padding | Not compatible with calf pads or Easy-Out front rigging. | C\$ 250.00 |
| ☐ Gel Knee Pads, Easy-Out Angle Adjustable Front Rigging | | C\$ 400.00 |
| ☐ Gel Knee Pads, Swing Away ELR Front Rigging | | C\$ 250.00 |
| ☐ Length Adjustable Heel Loops | | No charge |
| ☐ Standard Calf Strap | | C\$ 50.00 |
| ☐ Bodypoint Padded Calf Strap | Sizes greater than 22" are not padded. | C\$ 117.00 |
| ☐ Miller's Dynamic Footrest Mounts, Plug In Front Rigging | Not compatible with Easy-Out front rigging. | C\$ 770.00 |
| ☐ Miller's Dynamic Footrest Mounts, Swing Away Front Rigging | Not compatible with Easy-Out front rigging. Not available with 13" STF height. | C\$ 965.00 |
| ☐ Comfort Foot Box, Double | Without separator. Select seat width | |
| ☐ 16 | ☐ 18" | C\$ 454.00 |
| ☐ 20" | ☐ 22" | C\$ 644.00 |
| ☐ 24" | | C\$ 776.00 |
| ☐ Comfort Foot Box, Single | | |
| ☐ Left | ☐ Right | C\$ 285.00 |
| ☐ Pair | | C\$ 570.00 |

Positioning Belts

- | | | |
|---|--|------------|
| ☐ Bodypoint Non-Padded Center Pull Push Button Belt, Large | Not compatible with gel padding. Includes reinforced premium calf pad mount. | C\$ 125.00 |
| ☐ Bodypoint 2 Point Padded Center Pull Push Button Belt | Not compatible with gel padding. | C\$ 215.00 |
| ☐ Bodypoint 4 Point Padded Center Pull Push Button Belt | Not compatible with calf pads or Easy-Out front rigging. | C\$ 175.00 |
| ☐ Bodypoint Non-Padded Chest Belt 46" | | C\$ 165.00 |
| ☐ Q'Straint Lap Belt | Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" belt length. | C\$ 470.00 |

Miscellaneous Options

- | | | |
|---|---|--------------|
| ☐ One-Arm Drive | Available with 22" and 24" composite rear wheels. Minimum 17" STF and 18" STF respectively. Large hand setup includes 24" rear wheels with 22" driving handrim. The maximum achievable tilt could be as low as 5°-10° and be further reduced to 0° when transit tie down system is installed. Check with PDG prior to ordering. Not compatible with user and attendant wheel locks, drum brakes, quick release rear axles or pneumatic tires. | C\$ 2,800.00 |
| Installed On: | ☐ Right ☐ Left | |
| Handrim Setup: | ☐ Standard ☐ Large Hand | |
| | Large hand setup only available with 24" rear wheels. | |
| ☐ Spoke Guards | Available for 20"-24" composite rear wheels. Made from durable and lightweight aluminium. | C\$ 395.00 |
| ☐ Black | ☐ Frame color | |
| ☐ Gel Hip Guides | Reduces available seat width by 2" | C\$ 845.00 |
| ☐ IV Pole | Not compatible with heavy duty rear anti tippers, rear floor jacks or attendant wheel locks. | C\$ 545.00 |
| ☐ Installed on Right | ☐ Installed on Left | |
| ☐ O2 Holder | Not compatible with heavy duty rear anti tippers, rear floor jacks or attendant wheel locks. | C\$ 695.00 |
| ☐ Installed on Right | ☐ Installed on Left | |
| ☐ Medical Necessities Backpack | Black with Permobil logo. | C\$ 245.00 |
| ☐ Adhesive Back Hook Strip | For cushion mounting, 2"x30". | C\$ 15.00 |

General

Permobil manufactures PDG wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

Notes