


TRAM dimensions


User dimensions (inches)		K310 and K320
Girth		22 min – 55 max
 Key user dimension: girth Important: User's weight or height must not exceed 350 lb or 76 inches		
Item dimensions (inches)	K310 TRAM	K320 TRAM
Overall length	45½	Same
Overall width	27½–40	Same
Overall height	43½–58	Same
TRAM weight (lb)	70	Same
Turning diameter	50	Same
Floor to top of base leg for furniture clearance	6¾	4¾
Min. user armpit height (when standing or ambulating)	30½	Same
Max. user height (when standing or ambulating)	76	Same
Max. working load (lb)	350	Same

Pelvic support

User dimensions (inches)
Indicated by the red arrow below

Key user dimension: half of thigh circumference measured at the groin.


K306 small	K316 med.	K326 large
6 – 10	8 – 14	10 – 18



The diagram shows a top-down view of a person's pelvis and upper thighs. A red curved arrow points from the center of the pelvis towards the left thigh, indicating the measurement point for dimension 'a'. Dimension 'a' is the horizontal distance between the center of the pelvis and the center of the left thigh. Dimension 'b' is the horizontal distance between the center of the pelvis and the center of the right thigh. Dimension 'c' is the vertical distance from the center of the pelvis to the center of the right thigh.

Product dimensions (inches)			
	K306 small	K316 med.	K326 large
a	8½	12	13½
b	16	20	25
c	9	12	14

Hip positioner

Product dimensions (inches) – without pad		
		
	K667 small	K668 large
a Crotch length	2½	3½
b Width	1¼	1¾
c Length	9	11



Retail order form - January 2023

Use dimension chart to select appropriate size.

All Rifton products are made to order and are non-returnable and non-exchangeable.



* Required

		K310 Standard (better for gait training)	K320 Low-base (fits under low beds)
* Frame includes one battery, one battery charger, accessories tote and charger power cord. (International power cords available on request.)	CHOOSE	Frame with scale	K301 \$9,970 <input type="checkbox"/> K321 \$9,970 <input type="checkbox"/>
		Frame	K302 \$8,175 <input type="checkbox"/> K322 \$8,175 <input type="checkbox"/>
* Thigh straps (pair)	CHOOSE	Narrow	K308 \$130 <input type="checkbox"/> K308 \$130 <input type="checkbox"/>
		Wide	K338 \$130 <input type="checkbox"/> K338 \$130 <input type="checkbox"/>
Switch pole		K666 \$384 <input type="checkbox"/>	K666 \$384 <input type="checkbox"/>
Arm supports (pair)	CHOOSE	Small arm prompts	K305 \$846 <input type="checkbox"/> K305 \$846 <input type="checkbox"/>
		Large arm prompts	K315 \$975 <input type="checkbox"/> K315 \$975 <input type="checkbox"/>
		Arm platforms	K397 \$738 <input type="checkbox"/> K397 \$738 <input type="checkbox"/>
		Arm platforms w/ handgrips	K396 \$845 <input type="checkbox"/> K396 \$845 <input type="checkbox"/>
Pelvic support	CHOOSE	Small	K306 \$424 <input type="checkbox"/> K306 \$424 <input type="checkbox"/>
		Medium	K316 \$477 <input type="checkbox"/> K316 \$477 <input type="checkbox"/>
		Large	K326 \$508 <input type="checkbox"/> K326 \$508 <input type="checkbox"/>
Hip positioner	CHOOSE	Small	K667 \$520 <input type="checkbox"/> K667 \$520 <input type="checkbox"/>
		Large	K668 \$575 <input type="checkbox"/> K668 \$575 <input type="checkbox"/>
Hip positioner pad	CHOOSE	Small	K541 \$180 <input type="checkbox"/> K541 \$180 <input type="checkbox"/>
		Large	K542 \$237 <input type="checkbox"/> K542 \$237 <input type="checkbox"/>
Additional battery (recommended for institutional use)		K307 \$235 <input type="checkbox"/>	K307 \$235 <input type="checkbox"/>
Additional thigh straps	CHOOSE	Narrow single	K318 \$66 <input type="checkbox"/> K318 \$66 <input type="checkbox"/>
		Narrow pair	K308 \$130 <input type="checkbox"/> K308 \$130 <input type="checkbox"/>
		Wide single	K328 \$66 <input type="checkbox"/> K328 \$66 <input type="checkbox"/>
		Wide pair	K338 \$130 <input type="checkbox"/> K338 \$130 <input type="checkbox"/>
Swivel locks (pair)		K319 \$318 <input type="checkbox"/>	K324 \$318 <input type="checkbox"/>
Total		K310 \$ _____	K320 \$ _____