







M3 Corpus



M5 Corpus





Accent options







Mossy Oak[®] Breakup[®] Country™



Nater Lily



Fireball





Please note: Actual product colours may vary from colours shown

Midnight Arctic Violet Sky Aqua Spring Canary **Sunset** Chili

Flamingo

Platinum

Bronze

Gold



The redesigned essentials carrier provides the ultimate storage solution. A compact shell and smooth interior maximises storage space while minimising overall length.

~ • • •		
Specif	CLCST	inne.
JUECII	ILAL	IUI IO

	M3 Corpus	M5 Corpus
Maximum user weight	136kg	136kg
Maximum speed	10km/h	10km/h
Range	20.9km (Grp 34s) 32.2km (Grp 24s)	32.2km (Grp 24s)*
Base width	24"	25"
Base length	35.5"	37.5"
Base turning radius	20.5"	21.5"
Ground clearance	3"	3"
Recommended battery type	Group 34 or 24	Group 24
Weight including batteries	175kg (Grp 34s) 184kg (Grp 24s)	192kg (Grp 24s)
Drive electronics	R-Net 120A	R-Net 120A
Seat-to-floor height	17.5", 18.5", 19.5"	17.5", 18.5", 19.5"
Seat elevation	Optional 12"	Optional 14"
Drive wheel	14"	14"
Casters	7"	8"
Tilt options	Posterior — 0° – 50° ActiveReach — 5°, 10° & 20°	Posterior — 0° – 50° ActiveReach — 5°, 10° & 20°
Legrest angle - power	85° – 180°	85° – 180°
Backrest angle - power	85° – 180°	85° – 180°
Backrest height	20", 23" - 28" (1" incr.)	20", 23" – 28" (1" incr.)
Seat depth	14" – 22" (1" incr.)	14" – 22" (1" incr.)
Seat width	17" – 23" (2" incr.)	17" – 23" (2" incr.)
Arm pad lengths	10", 13", 16", 18"	10", 13", 16", 18"

*Actual driving range will vary based on driving and battery conditions. | Specifications may vary based upon configuration.
Please refer to your owners manual for warnings and instructions for the safe operation of your wheelchair in a variety of driving conditions.

