



# Retail order form - January 2023

Use dimension chart to select appropriate size.

All Rifton products are made to order and are non-returnable and non-exchangeable.



The **conversion kit** converts the K110 mini stander to the K130 small. It includes the deck and complete body support of the K130. Request from customer service.

**Conversion kit** K116




**\* Required**

**\* Mobile Stander** frame includes body support with padding, seat pad, straps and fleece.

**Specify color:**  
pink, blue or lime

CHOOSE

User dimensions (inches)	K110 mini	K130 small	K150 med.	K170 large
Height	25–32	30–40	37–52	48–65
 <b>Key user dimension: height</b> Select the appropriate Mobile Stander by the user's overall height. Choose the model that allows for growth.  <b>Important:</b> User's weight must not exceed the maximum working load.				
Item dimensions (inches)	K110 mini	K130 small	K150 med.	K170 large
Overall length and width	25 x 25	25 x 25½	32 x 29	39 x 33½
Overall height	22½–28	23½–30	31–39	37½–49
Width without large wheels	19½	19½	24½	28¾
Wheel diameter	20	20	27	35
Width of body support	7–10	7–10	9–12	10–14
Height of body support	16½–22	22–29	29–38	36–48
Item weight (lb)	28	30	40	65
Max. working load (lb)	50	50	120	175

		K110 mini		K130 small		K150 medium		K170 large	
<b>Wheels</b> (pair) (include brakes)		K101	\$3,212 <input type="checkbox"/>	K103	\$3,585 <input type="checkbox"/>	K105	\$4,356 <input type="checkbox"/>	K107	\$5,308 <input type="checkbox"/>
		Pink <input type="checkbox"/>		Pink <input type="checkbox"/>		Pink <input type="checkbox"/>		Pink <input type="checkbox"/>	
		Blue <input type="checkbox"/>		Blue <input type="checkbox"/>		Blue <input type="checkbox"/>		Blue <input type="checkbox"/>	
		Lime <input type="checkbox"/>		Lime <input type="checkbox"/>		Lime <input type="checkbox"/>		Lime <input type="checkbox"/>	
		K113	\$1,901 <input type="checkbox"/>	K113	\$1,901 <input type="checkbox"/>	K115	\$2,135 <input type="checkbox"/>	K117	\$2,576 <input type="checkbox"/>
<b>Mounting bracket</b>		K133	\$101 <input type="checkbox"/>	K133	\$101 <input type="checkbox"/>	K173	\$113 <input type="checkbox"/>	K173	\$113 <input type="checkbox"/>
<b>Arm prompts</b> (pair) (require mounting bracket)		K512	\$845 <input type="checkbox"/>	K512	\$845 <input type="checkbox"/>	K512	\$845 <input type="checkbox"/>	K512	\$845 <input type="checkbox"/>
<b>Hand loops</b> (require mounting bracket)	Single	K578	\$172 <input type="checkbox"/>	K578	\$172 <input type="checkbox"/>	K578	\$172 <input type="checkbox"/>	K578	\$172 <input type="checkbox"/>
	Pair	K579	\$331 <input type="checkbox"/>	K579	\$331 <input type="checkbox"/>	K579	\$331 <input type="checkbox"/>	K579	\$331 <input type="checkbox"/>
<b>Communication tray</b> (requires mounting bracket)		K518	\$647 <input type="checkbox"/>	K518	\$647 <input type="checkbox"/>	K518	\$647 <input type="checkbox"/>	K518	\$647 <input type="checkbox"/>
<b>Additional seat pad</b>	Small	N/A		N/A		K135	\$228 <input type="checkbox"/>	K135	\$228 <input type="checkbox"/>
	Large	N/A		N/A		K155	\$251 <input type="checkbox"/>	K155	\$251 <input type="checkbox"/>
<b>Sandals with wedges</b> (pair)	Small	K114	\$176 <input type="checkbox"/>	K114	\$176 <input type="checkbox"/>	N/A		N/A	
	Medium	N/A		N/A		K154	\$187 <input type="checkbox"/>	N/A	
	Large	N/A		N/A		N/A		K174	\$199 <input type="checkbox"/>
<b>Total</b>		<b>K110</b>	<b>\$ _____</b>	<b>K130</b>	<b>\$ _____</b>	<b>K150</b>	<b>\$ _____</b>	<b>K170</b>	<b>\$ _____</b>