

New product

Dreamline STX2 Cushion

Learn more about the
Dreamline STX2 cushion



The Dreamline STX2 cushion is a highly supportive, deep contoured, Gel Fluid cushion with abundant positioning options, for the less active user with mild to severe asymmetrical positioning requirements.

Key Features

- **Adjustability.** With an array of positioning tools for on-the-spot adjustability and adaptability for various postures.
- **Stability.** A firm positioning base with a deep well gives high lateral stability accentuated with its deep contouring.
- **Pressure Protection.** Large dual Cluster Gel Fluid packs nestled in 3cm of Enduro foam provides effective protection from skin breakdown, all secured within the FLO2™ cover for extra protection during transfers.



Unparalleled Positioning

The Dreamline STX2 cushion has a deep contour profile and incorporates advanced positioning options for symmetrical to severely asymmetrical postures. Pressure redistribution is effective through its patented Cluster Gel Fluid packs, which are neatly hidden in a unique secure zip-up compartment. A soft Comfort Wrap together with the FLO2™ fabric facilitates immersion and long-term comfort.



Technical Specifications

Weight Limit	136kg (300lbs) - Up to 20"			
Weight of Product	3.8kg* (8.38lbs)			
Width	12" - 20" (30.5cm - 50.8cm)			
Depth	12" - 22" (30.5cm - 55.9cm)			
Pelvic Well	1.8" (4.5cm)			
Positioning Base Height	12"/14"	16"/18"	20"	22"
	4.5cm	5.5cm	6.0cm	6.5cm
Wrap Thickness	2.5cm	3.0cm	3.0cm	3.0cm
Gel Packs	2 x Gel Fluid Packs			
Base Foam	EVA Closed Cell Foam			
Comfort Wrap	Enduro Foam with FLO2™ Cover			
Covers	FLO2™ Incontinence and Outer Cover			
	3D Spacer Side Panel			
	Neoprene Front & Side Panel Non-Slip Rubberised Urethane Base Panel			
Adjustable Positioning	Pelvic Obliquity / Pelvic Rotation			
	Abduction/Adduction			
	Windswept			
	Hip Flexion			
Warranty	18 Months			

*Weight Based on 18" x 18" (45.8cm x 45.8cm)

Product Code DX2C

