

			
M3 Corpus	M5 Corpus	M Corpus VS	M300 Corpus HD
150kg	150kg	120kg	204kg
10km/h	12 km/h	10km/h	8km/h
32.2km*	32.2km*	32.2km*	32.2km*
24"	25"	25"	26"
35.5"	37.5"	36.5"	38"
20.5"	21.5"	22"	22"
3"	3"	3"	3"
184kg (Grp 24)	192kg (Grp 24)	193kg (Grp 24)	203kg (Grp 24)
17.5", 18.5", 19.5"	17.5", 18.5", 19.5"	17.5", 18.5", 19.5"	18", 19", 20" (elevator) 18", 19", 20" (fixed)
Optional 12"	Optional 14"	12"	Optional 8"
Posterior: 0° – 50° ActiveReach: 5°, 10° & 20°	Posterior: 0° – 50° ActiveReach: 5°, 10° & 20°	Posterior: 0° – 50° Standing	Posterior: 0° – 45°
85° – 120° (manual) 85° – 180° (power)	85° – 120° (manual) 85° – 180° (power)	85° – 180° (power)	85° – 120° (manual) 85° – 150° (power)
85° – 180° (manual & power)	85° – 180° (manual & power)	85° – 180° (power)	85° – 175° (manual & power)
14" – 22" (by 1" incr.)	14" – 22" (by 1" incr.)	14" – 22" (by 1" incr.)	19" – 25" (by 1" incr.)
17", 19", 21", 23"	17", 19", 21", 23"	17", 19", 21", 23"	19", 21", 23", 25"
20", 23" – 28" (by 1" incr.)	20", 23" – 28" (by 1" incr.)	23" – 28" (by 1" incr.)	21", 24" – 29" (by 1" incr.)

Unmatched power positioning

Corpus is the most ergonomic and capable wheelchair seating system in the world. It is thoughtfully designed and configured to empower functional movements while supporting your optimal body position.



MyPermobility connectivity is an available option on all new M & F-Series Corpus & Corpus VS models.



App Store is a trademark of Apple Inc. Google Play is a trademark of Google LLC.

Power wheelchairs

Wheelchair matrix



permobil

Front-wheel

Mid-wheel

				
Rehab Series	F3 Corpus	F5 Corpus	F5 Corpus VS	M1
Max user weight	150kg	150kg	136kg	136kg
Max speed	10km/h 8km/h without support wheels	12km/h 10km/h without support wheels	12km/h	8km/h
Range	32.2km*	32.2km*	32.2km*	25.8km*
Base width	24"	25.5"	25.5"	24"
Base length	40" with support wheels 36" without support wheels	42" with support wheels 38" without support wheels	43"	35.5"
Base turning radius	26.5"	30"	30"	20.5"
Ground clearance	3"	3"	3"	3"
Weight w/ batteries	185kg (Grp 24)	187kg (Grp 24)	196kg (Grp 24)	145kg (Grp 34)
Seat-to-floor height	17.5", 18.5", 19.5"	17.5", 18.5", 19.5"	17.5", 18.5", 19.5"	16", 17", 18"
Seat elevator	Optional 12"	Optional 14"	14"	—
Power tilt options	Posterior: 0° – 50° ActiveReach: 5°, 10°, 20° & 30°	Posterior: 0° – 50° ActiveReach: 5°, 10°, 20° & 45°	Posterior: 0° – 50° Standing	Posterior: 0° – 50°
Recline options	85° – 120° (manual) 85° – 180° (power)	85° – 120° (manual) 85° – 180° (power)	85° – 180° (power)	—
Legrest options	85° – 180° (manual & power)	85° – 180° (manual & power)	85° – 180° (power)	85° – 140° (manual) 65° – 80° (adjust. swing-away)
Seat depth	14" – 22" (by 1" incr.)	14" – 22" (by 1" incr.)	14" – 22" (by 1" incr.)	14" – 22" (by 1" incr.)
Seat width	17", 19", 21", 23"	17", 19", 21", 23"	17", 19", 21", 23"	14" – 22" (by 2" incr.)
Backrest height	20", 23" – 28" (by 1" incr.)	20", 23" – 28" (by 1" incr.)	23" – 28" (by 1" incr.)	18" – 29"

All specifications may vary based upon configuration. | Chair weights are approximate and include base, batteries, power seat functions and seating system without accessories. | * Actual driving range will vary based on driving and battery conditions.



Front-wheel

Mid-wheel

Rear-wheel

					
Paediatric Series	Explorer Mini	Koala	K300 PS Jr.	M300 PS Jr.	K450 MX
Intended age range	12 months – 3 years	—	—	—	—
Max user weight	16kg	35kg	75kg	75kg	60kg
Max speed	2.4km/h	6.5km/h	8km/h	10km/h	8km/h
Range	5.6km*	24.1km*	20.9km*	20.9km*	27.4km*
Base width	19"	23"	24"	24"	25"
Base length	25"	32"	40"	36.5"	35"
Min turning radius	—	21"	25"	20"	28"
Ground clearance	1"	2.5"	3"	3"	2.5"
Weight w/ batteries	23.6kg	102kg (M 40)	123kg (Grp 34)	142kg (Grp 34)	153kg (Grp 34)
Seat-to-floor height	12.1" – 21.7"	15" (elevator)	16.5" (elevator) 16.5", 17.25", 18.5" (fixed)	17" (elevator) 17", 17.75", 18.5" (fixed)	3" – 26" (elevator)
Seat elevator	—	Up to 8"	Optional 8"	Optional 8"	8"
Power tilt options	—	0° – 25°	0° – 45°	0° – 45°	0° – 45°
Recline options	—	85° – 120° (manual)	85° – 120° (manual)	85° – 120° (manual)	85° – 105° (manual)
Legrest options	—	90° – 150° (manual)	70° (swing-away) 90 – 150° (center mount)	70° or 90° (swing-away) 90° – 150° (center mount)	90° – 125° (manual)
Seat depth	—	10", 12", 14"	10" – 18"	10" – 18"	10" – 18" (by 1" incr.)
Seat width	—	10", 12", 14"	11" – 16" (by 1" incr.)	11" – 16" (by 1" incr.)	10", 12", 14", 16"
Backrest height	Total chair height 29" – 37"	14", 16", 20"	16.5" – 22.5" (by 1" incr.)	16.5" – 22.5" (by 1" incr.)	14", 16", 20" (MX) 16.5" – 22.5" (PS Jr)

All specifications may vary based upon configuration. | Chair weights are approximate and include base, batteries, power seat functions and seating system without accessories. *Actual driving range will vary based on driving and battery conditions.