

TiLite X-Series Canadian Order Form

Effective Sept 2023 (New Items in Red)

FOLDING FRAME

Page 1 of 7

DEALER INFORMATION		RTS INFORMATION	
Dealer:	Account #:	RTS:	Phone:
PO:	Ordered By:	Email:	
Email:		SHIP TO INFORMATION	
Mark For:	User Weight:	Name:	Ship Via:
TiLite Sales Rep:		Address:	
TiLite Information:	Promo Code:	City:	Province: Postal Code:
Ship Date:	Order #:	Phone:	Fax:

More options and frame modifications may be available through our Creative Mobility Division! Please contact customer service for details or [click here](#) to view options.

Please Email or Fax to the Permobil Manual Customer Experience Team at tilitesales@permobil.com or 866.586.2413

FRAME STYLE

- AXFS1 ☐ **Aero X Aluminum Frame** (Full Adjustability. 265 (120.2kg) Pound Weight Limit.).....\$3270
- GXFS4 ☐ **2GX Titanium Frame** (Full Adjustability. 265 (120.2kg) Pound Weight Limit. Light Weight and High Strength. Not Available with Frame Transit Tie-Down.).....\$5478
- XFS1 ☐ **Swing Away**.....STD
- XFS3 ☐ **Hemi Height Frame** (Available with a Front Seat Height of 13.5" - 20" (34.3-50.8cm). Not Available with Depth Adjustable Back/Frame, 12" or 13" (30.5cm or 33cm) Seat Depth or Heavy Duty Frame.).....NCO
- XFS2 ☐ **Fixed Front**.....NCO
- XHDR1 ☐ **Heavy Duty Frame (350 (158.75kg) Pound Weight Limit.** Includes Tension Adjustable by Straps Back Upholstery. Available with 18" to 22" (45.7 to 55.9cm) Seat Widths & 16" to 20" (40.6 to 50.8cm) Seat Depths. Not Available with 90° front end).....\$602
- XRFR1 ☐ **Reinforced Frame** (No charge when 20" (50.8cm) Rear Seat Width and 20" (50.8cm) Seat Depth are selected.).....\$437

SMARTDRIVE POWER ASSIST SYSTEM



- XMAXSD1 ☐ **Permobil Smart Drive MX2+*** (Not available with rear seat heights less than 13.5", or 20" rear wheels).....\$7989
- XMAXSD16 ☐ **NEW** SpeedControl Dial XMAXSD17 ☐ Extended Mount **OR** XMAXSD18 ☐ SwitchControl Buttons.....NCO
- XMAXSD8 ☐ SwitchControl with Mono Jack Option.....\$208
- XMAXSD19 ☐ **NEW** Buddy Button.....\$368
- XMAXSD4 ☐ Carrying Bag.....\$110
- XMAXSD10 ☐ Universal Padded Handle.....\$31

*SmartDrive apps available on Apple App Store and Google Play Store for Apple Watch Series 5 or newer and Samsung Galaxy Watch 4 or newer.

TRANSIT OPTIONS

- XTTDO1 ☐ **Aluminum Frame Transit Tie-Down** (Consists of four Tie-Down Brackets. Not Available with 2GX2, Heavy Duty Frame, Hemi Height Frame, Swing Away Tubular Armrests or Swing Away- Flip Back Tubular Armrests. Requires a Folding Stabilizer Bar unless a Depth Adjustable Backrest is selected. Requires minimum seat depth of 14" (35.6cm)) (Please refer to the Aero X2 Transit Tie-Down Option Brochure for important information regarding the Transit Tie-Down Options.).....\$487
- XWTB ☐ **Transport Bracket Option** (Consists of four Tie-Down Brackets. Unoccupied Use Only.) (Not Available when Swing Away frame and Heavy Duty Frame are both selected. Not available with Swing Away Tubular Armrest or Swing Away -Flip Back Tubular Armrests. Requires a minimum seat depth of 14" (35.6cm) or 15" (38.1cm) seat Depth with Hemi Height Frame.).....\$422

FRAME FINISHES

Natural Finishes

- XFF1 ☐ **Satin** (Available with Titanium Frame only).....STD
- XSCBT1 ☐ **Satin Finish Buffing Pads**.....\$22
- XFF3 ☐ **Tattoos** (Tattoos on front frame tubes. Titanium frame finished in Satin Aluminum frame finished in blasted aluminum with clear coat. Tattoos will be trimmed and adjusted to fit your frame.).....\$394

XFF45 ☐ **Need For Speed**



XFF44 ☐ **Flower Garden**



XFF43 ☐ **Mosaic**



- | | | | |
|-----------------------------|--|---|--|
| High Gloss Paint | <input type="checkbox"/> XFF8 High Gloss Black | <input type="checkbox"/> XFF10 TiLite Red | <input type="checkbox"/> XFF29 Neon Pink |
| Metallic Paint | <input type="checkbox"/> XFF14 Ocean Blue Metallic | <input type="checkbox"/> XFF15 Tangerine Metallic | <input type="checkbox"/> XFF30 Acid Green Metallic |
| Pearlescent Paint | <input type="checkbox"/> XFF17 Midnight Blue Pearl | <input type="checkbox"/> XFF18 White Pearl | |
| Matte Paint | <input type="checkbox"/> XFF19 Black Lava | <input type="checkbox"/> XFF27 Titanium Grey | |
| Candy Paint | <input type="checkbox"/> XFF22 Ruby Red | <input type="checkbox"/> XFF23 Azure Blue | <input type="checkbox"/> XFF25 Electric Plum |
| Matte Textured Paint | <input type="checkbox"/> XFF40 Hyper Blue | <input type="checkbox"/> XFF41 Flame Orange | <input type="checkbox"/> XFF42 OD Green.....STD |

TiLite X-Series Canadian Order Form

Effective Sept 2023 (New Items in Red)

FOLDING FRAME



Page 2 of 7

FRAME MEASUREMENTS

Measurements shown below are our most popular selections. Not seeing what you want? Contact customer service or visit our Creative Mobility Division (CMD) catalog [here](#) to view additional options.

SEAT WIDTH A Measure from outside of seat tube at back post to the outside of the opposite seat tube at back post.

XSW1 ☐ 12" ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19"..... **STD**
XSW2 ☐ 20" ☐ 21" ☐ 22" (21" and 22" Seat Widths A are only available with Heavy Duty Frame.)..... **\$541**



SEAT DEPTH B Measure from 1/2" (1cm) in front of back post to the front edge of seat sling. (12"-14" (30-35.6cm) Seat Depths B not available

With 25" (559) and 26" (560) Rear Wheels.)

XSD1 ☐ 12" (30.5cm) ☐ 13" (33cm) ☐ 14" (35.6cm) ☐ 15" (38.1cm) ☐ 16" (40.6cm) ☐ 17" (43.2cm)..... **STD**
XSD2 ☐ 18" (45.7cm) ☐ 19" (48.3cm)..... **\$590**
XSD3 ☐ 20" (50.8cm)..... **NCO**
☐ Custom Frame Depth: _____" Frame Depth Longer Than Seat Depth (Use 1" (2.5cm) increments with a maximum of 2" (5.1cm). A Custom Frame Depth with a maximum of 20" (50.8cm) is additional charge.)

FRONT SEAT HEIGHT C Measure from floor to top of seat tube at front edge of seat sling Adjustable in 1/2" (1.3cm) increments. (13.5" - 14" (34.3-35.6cm) Front Seat Height C requires Hemi Height Frame Style. 20.5" - 21" (52.1-53.3cm) Front Seat Height C not available with Hemi Height Frame Style.

FRONT SEAT HEIGHT RANGES

Front Wheel Size	Front Seat Heights C	Hemi Front Seat Heights C
3" Front Wheel	14.5" - 18.5" (36.8-47cm)	13.5" - 17.5" (34.3-44.5cm)
3" Front Wheel w/Glide™	16.0" - 17.5" (40.6-44.5cm)	15.0" - 17.0" (38.1-43.2cm)
4" Front Wheel	15.5" - 19.0" (39.4-48.3cm)	14.5" - 18.0" (36.8-45.7cm)
4" Front Wheel w/Glide™	16.5" - 18.0" (41.9-45.7cm)	15.5" - 17" (39.4-43.2cm)
5" Front Wheel	16.5" - 19.5" (41.9-49.5cm)	16.0" - 18.5" (40.6-47cm)
5" Front Wheel w/Glide™	17.5" - 18.5" (44.5-47cm)	16.5" - 18.0" (41.9-45.7cm)
6" Front Wheel	17.5" - 20.0" (44.5-50.8cm)	16.5" - 19.0" (41.9-48.3cm)
6" Front Wheel w/Glide™	19.0" (48.3cm)	17.5" - 18" (44.5-45.7cm)
8" Front Wheel	19.5" - 21.0" (49.5-53.3cm)	18.5" - 19.5" (47-50.8cm)

XFSH1 ☐ 13.5" (34.3cm) ☐ 14" (35.6cm) ☐ 14.5" (36.8cm)
☐ 15" (38.1cm) ☐ 15.5" (39.4cm) ☐ 16" (40.6cm)
☐ 16.5" (41.9cm) ☐ 17" (43.2cm) ☐ 17.5" (44.5cm)
☐ 18" (45.7cm) ☐ 18.5" (47cm) ☐ 19" (48.3cm)
☐ 19.5" (49.5cm) ☐ 20" (50.8cm) ☐ 20.5" (52.1cm)
☐ 21" (53.3cm)..... **STD**

REAR SEAT HEIGHT D Measure from floor to top of seat tube at back post. Adjustable in 3/8" increments.

REAR SEAT HEIGHT RANGES

Rear Wheel Size	Rear Seat Height D
22"	13" - 18" (33-45.7cm)
24"	14" - 19" (35.6-48.3cm)
25"	14" - 19" (35.6-48.3cm)
26"	15" - 20" (38.1-50.8cm)

XRSH1 ☐ 13" (33cm) ☐ 13.5" (34.3cm) ☐ 14" (35.6cm)
☐ 14.5" (36.8cm) ☐ 15" (38.1cm) ☐ 15.5" (39.4cm)
☐ 16" (40.6cm) ☐ 16.5" (41.9cm) ☐ 17" (43.2cm)
☐ 17.5" (44.5cm) ☐ 18" (45.7cm) ☐ 18.5" (47cm)
☐ 19" (48.3cm) ☐ 19.5" (49.5cm) ☐ 20" (50.8cm) **STD**

BACK SUPPORTS

COMFORT ACTA-BACK® Back Support

Shell Height ☐ 10" ☐ 12" ☐ 14" ☐ 16" ☐ 18"..... **\$896**
☐ 20"..... **\$896**
☐ Acta-Back with 2-PT Quick Release Compass Hardware..... **NCO**
☐ Comfort-Tek Cover ☐ Stretch-Air Cover..... **NCO**
☐ Color Accent (Colored piping along the edge of the cover on the backside of the backrest) ☐ Silver..... **STD**
☐ Red ☐ Purple ☐ Pink ☐ Blue ☐ Green ☐ Orange ☐ Remove (no color)..... **\$58**
☐ Comfort Privacy Shield..... **\$80**

COMFORT ACTA-BACK® Deep Back Support

Shell Height ☐ 10" ☐ 12" ☐ 14" ☐ 16" ☐ 18"..... **\$896**
☐ 20"..... **\$896**
☐ Acta-Back with 2-PT Quick Release Compass Hardware..... **NCO**

TiLite X-Series Canadian Order Form

Effective Sept 2023 (New Items in Red)

FOLDING FRAME

Page 3 of 7

- ☐ Comfort-Tek Cover ☐ Stretch-Air Cover..... **NCO**
☐ 3" Deep Wing ☐ 5" Deep Wing..... **NCO**
☐ Color Accent (Colored piping along the edge of the cover on the backside of the backrest) ☐ Silver..... **STD**
☐ Red ☐ Purple ☐ Pink ☐ Blue ☐ Green ☐ Orange ☐ Remove (no color)..... **\$58**
☐ Comfort Privacy Shield **\$80**

BACK SUPPORTS

COMFORT ACTA-EMBRACE™

Shell Height ☐ 10" ☐ 12" ☐ 14" ☐ 16"

SBRS16 ☐ Acta-Embrace Back with 2-PT Quick Release Compass Hardware..... **\$840**

CPSB1 ☐ Comfort Privacy Shield..... **\$58**

COMFORT ACTA-RELIEF

Shell Height ☐ 10" ☐ 12" ☐ 14" ☐ 16"

SBRS13 ☐ Acta-Relief Back with 2-PT Quick Release Compass Hardware..... **\$1061**

CPSB1 ☐ Comfort Privacy Shield..... **\$58**

SEAT BACK HEIGHT E/SEAT BACK/PUSH HANDLES Measure from top of back post to top of seat tube at rear of frame.

SEAT BACK TYPE

XSB20 ☐ Fixed Angle, Adjustable-Height Seat Back..... **STD**

XLB1 ☐ Folding Stabilizer Bar (Required with Transit Tie-Down Option unless Depth Adjustable Backrest is selected. Includes Tension Adjustable by Strap Back Upholstery.)..... **\$427**

XSB14 ☐ Fixed Angle, Adjustable Height Seat Back w/8-Degree Bend. (Not available with Heavy Duty Frame.)..... **NCO**

XSB18 ☐ Folding Adjustable-Angle and Height Seat Back Not Available with Hemi Height Frame, Heavy Duty Frame, 19" or 20" (48.3 or 50.8cm) Seat Depths B or Swing Away Tubular Armrest. Seat Depths B of 12"-17" (30.5-43.2cm) allow 3" (7.6cm) of rearward adjustment. Seat Depth B of 18" (45.7cm) allow 2" (5.1cm) of rearward adjustment. Requires Tension Adjustable Seat Upholstery)..... **\$252**

XSB16 ☐ Depth Adjustable Back/Frame (Not Available with Hemi Height Frame, Heavy Duty Frame, 19" or 20" (48.3 or 50.8cm) Seat Depths B, Swing Away Tubular or Swing Away Flip-Back Armrest. Seat Depths B of 12"-17" (30.5-43.2cm) allow 3" (7.6cm) of rearward adjustment. Seat Depth B of 18" (46cm) allow 2" (5cm) of rearward adjustment. Requires Tension Adjustable Seat Upholstery)..... **NCO**

PUSH HANDLES

XSB6 ☐ Integrated Push Handles..... **NCO**

XSB7 ☐ Integrated Fold-Down Push Handles (Requires Tension Adjustable by Straps Back Upholstery. Not available with Heavy Duty Frame.)..... **\$126**

XSB9 ☐ Bolt-On Height Adjustable Push Handles (Requires Tension Adjustable by Straps Back Upholstery. Not available with Heavy Duty Frame, Folding Adjustable-Angle and Height Seat Back, Swing Away Tubular Armrests, Swing Away - Flip Back Tubular Armrests.)..... **\$449**

FRAME MEASUREMENTS

Options shown below are our most popular selections. **Not seeing what you want?** Contact customer service or visit our Creative Mobility Division (CMD) catalog [here](#) to view additional options.

XSB20/XSB14 **FIXED ANGLE ADJUSTABLE SEAT BACK HEIGHT E** (11"-15" Fixed Seat Back Height not available with Heavy Duty Frame, Swing Away Tubular Arms 12"-15", Bolt-On Adjustable Height Push Handles or Folding Stabilizer Bar.)

Push Handles

(Selected on Pg. 2)	Short	Tall	Set at: (Use 1" Increments)
None	<input type="checkbox"/> 11.0" - 15.0"	<input type="checkbox"/> 16" - 20"	Set at: _____
Integrated	<input type="checkbox"/> 11.0" - 15.0"	<input type="checkbox"/> 16" - 20"	Set at: _____
Integrated Fold-Down	<input type="checkbox"/> 11.0" - 15.0"	<input type="checkbox"/> 16" - 20"	Set at: _____
Bolt-On Height Adj.	-----	<input type="checkbox"/> 16" - 20"	Set at: _____

XSB18 **FOLDING ADJUSTABLE- ANGLE AND HEIGHT SEAT BACK E** (Not available with a Folding Stabilizer Bar or Heavy Duty Frame)

Push Handles

(Selected on Pg. 2)	Short	Medium	Tall	Set at: (Use 1" Increments)
None	<input type="checkbox"/> 11.0" - 13.0"	<input type="checkbox"/> 13.5" - 16.5"	<input type="checkbox"/> 16.5" - 19.5"	Set at: _____
Integrated	<input type="checkbox"/> 12.5" - 14.5"	<input type="checkbox"/> 13.5" - 16.5"	<input type="checkbox"/> 17.0" - 19.5"	Set at: _____
Integrated Fold-Down	<input type="checkbox"/> 15.0" - 17.0"	<input type="checkbox"/> 16.0" - 19.0"	<input type="checkbox"/> 18.5" - 20.5"	Set at: _____

CAMBER F

XCBR8 ☐ 0° ☐ 2° ☐ 4° ☐ 6°..... **STD**

FRONT ANGLE H Front Angle H is measured from inside of front frame tube to top of bottom frame tube.

Swing Away Front End

XFA3 ☐ 70° ☐ 80° ☐ 85° (Not available with 6" Front Wheels.) ☐ 90° (Requires aluminum footplate. Not available with Heavy Duty Frame, 5" or 6" Front Wheels.)..... **STD**

XFA0 ☐ Omit Hangers..... **NCO**

Fixed Front End

XFA6 ☐ 70° ☐ 80° ☐ 85° (85° Not available with 6" Front Wheels)..... **STD**

TiLite X-Series Canadian Order Form

Effective Sept 2023 (New Items in Red)

FOLDING FRAME

Page 4 of 7

SEAT TO FOOTREST //FOOTRESTS Measure from front edge of seat sling to top rear of footrest.

If your Front Angle H is 70°, 80° OR 85° select one Seat to Footrest I length and one Footrest type:

Standard Footrests Measurement must be at least 2" shorter than Front Seat Height C. 70° Fixed Front has a minimum seat to footrest of 14.5"

☐ 14" (35.6cm) ☐ 14.5" (36.8cm) ☐ 15" (38.1cm) ☐ 15.5" (39.4cm) ☐ 16" (40.6cm) ☐ 16.5" (41.9cm)
☐ 17" (43.2cm) ☐ 17.5" (44.5cm) ☐ 18" (45.7cm) ☐ 18.5" (47cm) ☐ 19" (48.3cm) ☐ 19.5" (49.5cm) ☐ 20" (50.8cm)

XFTR5	<input type="checkbox"/> Flip-up Footrest.....	STD
XFTR9	<input type="checkbox"/> Depth Adjustable Flip-up Footrest.....	\$121
XFTR10	<input type="checkbox"/> Depth & Angle Adjustable Flip-up Footrest (Available for 15"- 20" (38.1-50.8cm) Seat to Footrest I only.).....	\$296
XFTR11	<input type="checkbox"/> Angle Adjustable Flip-up Footrest (Available for 15"- 20" (38.1-50.8cm) Seat to Footrest I only.).....	\$257

High Mount Footrests (Not Available with Luggage Carriers, or Impact Guards)

☐ 7" (17.8cm) ☐ 7.5" (19.1cm) ☐ 8" (20.3cm) ☐ 8.5" (21.6cm) ☐ 9" (22.9cm) ☐ 9.5" (24.1cm)
☐ 10" (25.4cm) ☐ 10.5" (26.7cm) ☐ 11" (27.9cm) ☐ 11.5" (29.2cm) ☐ 12" (30.5cm) ☐ 12.5" (31.8cm)
☐ 13" (33cm) ☐ 13.5" (34.3cm) ☐ 14" (35.6cm) ☐ 14.5" (36.8cm)

XFTR20	<input type="checkbox"/> High Mount Flip-up.....	\$208
XFTR15	<input type="checkbox"/> Angle Adjustable High Mount Flip-up Footrest.....	\$356
XFTR16	<input type="checkbox"/> Depth & Angle Adjustable High Mount Flip-up Footrest.....	\$410

Elevating Leg Rest

☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19" ☐ 20" ☐ 21"

XFA5	<input type="checkbox"/> Single Depth Adjustable Flip-up Footrest - <input type="checkbox"/> Left OR <input type="checkbox"/> Right.....	\$187
XFTR17	<input type="checkbox"/> Depth Adjustable Flip-up Footrest - Pair.....	\$394

If your Front Angle H is 90°, select one Seat to Footrest I length and one Footrest type:

☐ 7" (17.8cm) ☐ 8" (20.3cm) ☐ 9" (22.9cm) ☐ 10" (25.4cm) ☐ 11" (27.9cm) ☐ 12" (30.5cm) ☐ 13" (33cm)

XFTR5	<input type="checkbox"/> Flip-up Footrest.....	STD
XFTR9	<input type="checkbox"/> Depth Adjustable Flip-up Footrest.....	\$121
XFTR10	<input type="checkbox"/> Depth & Angle Adjustable Flip-up Footrest.....	\$296
XHLP1	<input type="checkbox"/> Heel Loops (Not available with 13" (33cm) or narrower Seat Width A or 1-Piece Footrests.).....	\$66

CUSTOM CONFIGURATION OPTIONS

FRONT WHEELS (See Front Seat Height Ranges Chart on Page 2 for Front Seat Height C limitations.)

XFW24	<input type="checkbox"/> 3"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire <input type="checkbox"/> Silver Hub OR <input type="checkbox"/> Black Hub.....	\$197
XFW14	<input type="checkbox"/> 4"x.75" Plastic Wheel w/Poly Tire (Not available with Power Adaptable Frame, Heavy Duty Frame or Glide™ Single-Sided Suspension Fork.).....	NCO
XFW3	<input type="checkbox"/> 4"x.75" Light-Up Micro-Caster (Not available with Power Adaptable Frame, Heavy Duty Frame or Glide™ Single-Sided Suspension Fork.).....	\$88
XFW26	<input type="checkbox"/> 4"x1.5" LiteSpeed Plastic Wheel w/Soft Roll Tire (Not available with Power Assist, Power Add-on System or Heavy Duty Frame.).....	\$104
XFW27	<input type="checkbox"/> 4"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire <input type="checkbox"/> Silver Hub OR <input type="checkbox"/> Black Hub.....	\$208
XFW6	<input type="checkbox"/> 5"x1" Plastic Wheel w/Poly Tire (Not available with Heavy Duty Frame or with 90° Front Angle H or Glide™ Single-Sided Suspension Fork.).....	NCO
XFW28	<input type="checkbox"/> 5"x1.5" LiteSpeed Plastic Wheel w/Soft Roll Tire (Not available with Heavy Duty Frame or 90° Front Angle H.).....	\$100
XFW30	<input type="checkbox"/> 5"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire <input type="checkbox"/> Silver Hub OR <input type="checkbox"/> Black Hub (Not available with +90° Front Angle H.).....	\$274
XFW9	<input type="checkbox"/> 6"x1" Plastic Wheel w/Poly Tire (Not available with Power Adaptable Frame, Heavy Duty Frame or 85° or 90° Front Angle H or Glide™ Single-Sided Suspension Fork.).....	NCO
XFW32	<input type="checkbox"/> 6"x1.5" LiteSpeed Plastic Wheel w/Soft Roll Tire (Not available with Power Adaptable Frame, Heavy Duty Frame, 85° or 90° Front Angle H, or Glide™ Single-Sided Fork.).....	\$110
XFW33	<input type="checkbox"/> 6"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire <input type="checkbox"/> Silver Hub OR <input type="checkbox"/> Black Hub (Not available with 85° or 90° Front Angle H or Glide™ Single-Sided Suspension Fork.).....	\$274
X2FW16	<input type="checkbox"/> 8"x1" Plastic Wheel w/Poly Tire (Not available with Power Assist, Power Add-on, Heavy Duty Frame or Glide™ Single-Sided Available with 70° Front Angle H only.) (Not available with 12"-14" Seat Depths).....	NCO

FRONT FORKS

XFk1	<input type="checkbox"/> TiLite Standard Fork.....	STD
XGL1	<input type="checkbox"/> Glide™ Suspension Fork (Available in Black anodize only. Not available with 8" Front Wheels or Heavy Duty Frame).....	\$596

Rider's Weight: ☐ 125 lbs. ☐ 125-180 lbs. ☐ 180-265 lbs.

ONE ARM DRIVE

XOAD1	<input type="checkbox"/> One Arm Drive (Available with 24" Mag Wheels only. Rear Seat Height, Camber or axle restrictions may apply.).....	\$1704
	Drive Handrim Side (Available with standard aluminum handrim only.)	
	<input type="checkbox"/> Install on Left <input type="checkbox"/> Install on Right	

STANDARD REAR WHEELS (See Rear Seat Height Ranges Chart on Page 2 for Rear Seat Height D limitations.)

XRW14	TiLite SHADOW™.....	STD
	<input type="checkbox"/> 22"(501)	
	<input type="checkbox"/> 24"(540)	
	<input type="checkbox"/> 25"(559) (Requires 15" (38.1cm) minimum Seat Depth B.)	
	<input type="checkbox"/> 26"(590) (Requires 15" (38.1cm) minimum Seat Depth B.)	

TiLite X-Series Canadian Order Form

Effective Sept 2023 (New Items in Red)

FOLDING FRAME

Page 5 of 7

XRW7	Spinergy SPOX <input type="checkbox"/> Black <input type="checkbox"/> Red <input type="checkbox"/> Blue <input type="checkbox"/> White.....	\$983
	<input type="checkbox"/> 22"(501) - Std. Hub	
	<input type="checkbox"/> 24"(540) - Std. Hub	
	<input type="checkbox"/> 25"(559) - Std. Hub (Requires 15" (38.1cm) minimum Seat Depth B.)	
	<input type="checkbox"/> 26"(590) - Std. Hub (Requires 15" minimum Seat Depth B.)	
XRW8	Spinergy LX (Only available with Heavy Duty Frame when 0 or 2 degree camber is selected.) <input type="checkbox"/> Black <input type="checkbox"/> Red <input type="checkbox"/> Blue <input type="checkbox"/> White <input type="checkbox"/> Bright Orange <input type="checkbox"/> Flamingo Pink <input type="checkbox"/> Purple <input type="checkbox"/> Green.....	\$1153
	<input type="checkbox"/> 22"(501) - Std. Hub	
	<input type="checkbox"/> 24"(540) - Std. Hub	
	<input type="checkbox"/> 25"(559) - Std. Hub (Requires 15" (38.1cm) minimum Seat Depth B.)	
	<input type="checkbox"/> 26"(590) - Std. Hub (Requires 15" (38.1cm) minimum Seat Depth B. Not available with Quad Release Axles.)	
XRW9	Spinergy Design Lab LX Wheels (If same color is selected, anodized hubs and rims match color anodized package on page 1. Available with black spokes only.) <input type="checkbox"/> Gold Anodize <input type="checkbox"/> Green Anodize <input type="checkbox"/> Blue Anodize <input type="checkbox"/> Red Anodize <input type="checkbox"/> Purple Anodize <input type="checkbox"/> Orange Anodize.....	\$1928
	<input type="checkbox"/> 22"(501) - Std. Hub	
	<input type="checkbox"/> 24"(540) - Std. Hub	
	<input type="checkbox"/> 25"(559) - Std. Hub (Requires 15" minimum Seat Depth B.)	
	<input type="checkbox"/> 26"(590) - Std. Hub (Requires 15" minimum Seat Depth B.)	
XRW16	Spinergy® Carbon Blade with 3k Weave Carbon Fiber (Not available with Heavy Duty Frame.).....	\$1682
	<input type="checkbox"/> 24"(540)	
	<input type="checkbox"/> 25"(559) (Requires 15" minimum Seat Depth B.)	
XRW17	X-Core Mag (Not available with Heavy Duty Frame, Short Tabs, Quad Release Axles, Spoke Guards or Omit Handrims Option Requires Aluminum Silver Anodized Handrims)	\$121
	<input type="checkbox"/> 22"(501)	
	<input type="checkbox"/> 24"(540)	

AXLES

XAXL1	<input type="checkbox"/> Stainless Quick Release.....	STD
XAXL2	<input type="checkbox"/> Titanium Quick Release.....	\$318
XAXL3	<input type="checkbox"/> Quad Release (Not available with Mag or 26"(590) Spinergy Spox or Spinergy LX Rear Wheels.).....	\$187

PNEUMATIC REAR TIRES

XRTR1	<input type="checkbox"/> Treaded 22"(501), 24"(540), and 26"(590) (Not available with 25"(559) Rear Wheels.).....	STD
XRTR2	<input type="checkbox"/> Airless Inserts.....	\$99
XRTR3	<input type="checkbox"/> 25"(559) Primo Treaded (Standard with 25"(559) Rear Wheels.).....	STD
XRTR4	<input type="checkbox"/> Primo.....	NCO
XRTR31	<input type="checkbox"/> MBL SpeedLite 20"(451), 22"(501), 24"(540), 25"(559) and 26"(590).....	STD
XRTR30	<input type="checkbox"/> 24"(540), 25"(559) and 26"(590) Primo Sentinel.....	\$175
XRTR8	<input type="checkbox"/> 24"(540), 25"(559) and 26"(590) Schwalbe® Marathon Plus Evolution.....	\$170
XRTR32	<input type="checkbox"/> 24"(540) Primo C-763 Knobby Tires.....	\$106

SOLID REAR TIRES

XRTR22	<input type="checkbox"/> 22"(501), 24"(540) and 25"(559) Primo Express (Not available with Spinergy Carbon Blade Rear Wheels.).....	\$121
XRTR23	<input type="checkbox"/> 24"(540) and 25"(559) Primo Active.....	\$148

HANDRIMS

XRIM1	<input type="checkbox"/> Aluminum – Silver Anodized.....	STD
XRIM21	<input type="checkbox"/> Aluminum – Black Anodized.....	\$66
XRIM12	<input type="checkbox"/> Titanium – Satin Finish (Not available with 20"(451) or 22"(501) Rear Wheels, or 24"(540) Knobby tires.).....	\$536
The Natural-Fit® (Original, Larger Profile)		
XRIM14	<input type="checkbox"/> Standard Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.).....	\$454
XRIM15	<input type="checkbox"/> Super Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.).....	\$514
XRIM16	<input type="checkbox"/> Without Thumb Grip.....	\$415
The Natural-Fit® LT (Reduced Weight, Smaller Profile)		
XRIM17	<input type="checkbox"/> Standard Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.).....	\$546
XRIM18	<input type="checkbox"/> Super Grip (Not available with Short Tabs or Spinergy Carbon Blade Rear Wheels.).....	\$585
XRIM19	<input type="checkbox"/> Without Thumb Grip.....	\$503
XRIM10	<input type="checkbox"/> The Surge (All-In-One Oval with Gription™ Strip, Larger Profile).....	\$476
XRIM11	<input type="checkbox"/> The Surge LT (All-In-One Oval with Gription™ Strip, Smaller Profile).....	\$514
XRIM20	<input type="checkbox"/> Q-Grip (Exclusive High-Friction, Anti-Slip Coating).....	\$410

Additional Handrim Options

XRIM5	<input type="checkbox"/> Plastic Coated – Black.....	\$175
XRIM22	<input type="checkbox"/> Omit Handrims.....	NCO

Tab Length

XRIM7	<input type="checkbox"/> Long Tabs.....	STD
XRIM8	<input type="checkbox"/> Short Tabs (Not available with Mag Rear Wheels or The Natural-Fit® Handrims with Thumb Grip.).....	\$39

WHEEL LOCKS (Composite wheel locks not always compatible with 12"-13" (30-33cm) Seat Depth. Uni-locks will be substituted at no charge.)

XWLK7	<input type="checkbox"/> Compact Composite <u>Push</u> to Lock.....	NCO
XWLK8	<input type="checkbox"/> Compact Composite <u>Pull</u> to Lock.....	NCO

TiLite X-Series Canadian Order Form

Effective Sept 2023 (New Items in Red)

FOLDING FRAME

Page 6 of 7

XWLK10	<input type="checkbox"/> Extension Handles for Compact Composite Push/Pull to Lock.....	\$94
Additional Wheel Locks		
XWLK11	<input type="checkbox"/> MBL Composite Scissor Lock (Not available on Swing Away.).....	\$225
XWLK3	Uni-Lock Set At: <input type="checkbox"/> Push OR <input type="checkbox"/> Pull.....	\$181
XWLK4	<input type="checkbox"/> Extension Handles for Uni-Lock.....	\$94
XWLK15	<input type="checkbox"/> Grade Aids Push to Lock (Not available with solid rear tires).....	\$252

BACK UPHOLSTERY (Want more seating options? [Click here](#) to view the new TiLite Seating Order Form.)

XBUP1	<input type="checkbox"/> Padded Nylon (Black Only. Not available with Heavy Duty Frame, Folding Back, Folding Stabilizer Bar, Depth Adjustable Frame, Fold-Down Push Handles or Ultrasuede Backrest Stripe.).....	STD
Tension Adjustable by Straps (Included with Heavy Duty Frame. Required with Aluminum Fendered Side Guards and Fold-Down Push Handles.)		
XBUP3	<input type="checkbox"/> Padded Nylon.....	\$454
XBUP5	<input type="checkbox"/> AIR Breathable Upholstery (Not available with Ultrasuede Backrest Stripe.).....	\$437
XUPC3	Ultrasuede Backrest Stripe <input type="checkbox"/> Black (Not available with Padded Nylon or AIR Breathable Back Upholstery.).....	\$116
XBUP0	<input type="checkbox"/> Omit Back Upholstery.....	NCO

SEAT UPHOLSTERY

XSUP4	<input type="checkbox"/> Tension Adjustable (Black only).....	STD
XSUP0	<input type="checkbox"/> Omit Seat Upholstery.....	NCO

SEAT CUSHION (Want more seating options? [Click here](#) to view the new TiLite Seating Order Form.)

XCSH2	<input type="checkbox"/> Comfort Curve Cushion (Stretch-Air cover).....	\$150
XCSH23	<input type="checkbox"/> Comfort Embrace Zero Elevation Cushion (Comfort-Tek cover).....	\$522
XCSH15	<input type="checkbox"/> ROHO® Hybrid Elite SR™ Single Compartment Cushion (Fits seat widths 14" to 19").....	\$739
XCSH14	<input type="checkbox"/> ROHO® QUADTRO SELECT® Cushion (Fits seat widths 14" to 19") Profile: <input type="checkbox"/> Low <input type="checkbox"/> Mid <input type="checkbox"/> High.....	\$711
XCSH16	<input type="checkbox"/> ROHO® Hybrid Select Cushion (Fits seat widths 14" to 19").....	\$756

CUSTOM CONFIGURATION OPTIONS

SIDE GUARDS

Aluminum

XSDG2	<input type="checkbox"/> Rigid Removable <input type="checkbox"/> Small OR <input type="checkbox"/> Std.....	\$318
XSDG3	<input type="checkbox"/> Quad Removable w/Hole <input type="checkbox"/> Small OR <input type="checkbox"/> Std.....	\$345
XSDG6	<input type="checkbox"/> Fendered Adjustable (Requires Tension Adjustable by Straps Back Upholstery. Not available with Angle-Adjustable Fold-Down Seat Backs.).....	\$437

Carbon Fiber (2x2 Twill Weave)

XSDG7	<input type="checkbox"/> Rigid Removable <input type="checkbox"/> Small OR <input type="checkbox"/> Std.....	\$372
XSDG8	<input type="checkbox"/> Quad Removable w/Hole <input type="checkbox"/> Small OR <input type="checkbox"/> Std.....	\$410

☐ **Custom Offset Profiled Sideguards** CMD 007. Please specify desired offset to nearest ±½" (e.g., ±0.50") and sideguard type (rigid or fold-down) in the Notes section of the order form. Sideguard material must be selected. Some limitations may apply.

ARMRESTS (Desk Arms come with built-in Side Guards.)

XARM1	Swing Away Tubular (Not available with Heavy Duty Frame, Folding Stabilizer Bar with Short or Medium Seat Back Height ranges or Depth Adjustable Back.).....	\$284
	<input type="checkbox"/> 9" - 12" (22.9 to 30.5cm) OR <input type="checkbox"/> 12" - 15" (30.5 to 38.1cm)	
XARM4	<input type="checkbox"/> Swing Away - Flip Back Tubular (Not available with Heavy Duty Frame, Folding Stabilizer Bar or Depth Adjustable Back.).....	\$345
XARM2	Removable Height Adjustable Desk Arm with Rigid Side Guard (Requires 14" minimum Seat Depth B.).....	\$394
	<input type="checkbox"/> 9" to 12" (22.9 to 30.5cm) OR <input type="checkbox"/> 12" to 15" (30.5 to 38.1cm)	
	<input type="checkbox"/> Standard Length Pad (10") OR <input type="checkbox"/> Full Length Pad (14")	
XARM3	Removable Height Adjustable Desk Arm with Rigid Side Guard and Transfer Loop (Requires 16" minimum Seat Depth B.).....	\$410
	<input type="checkbox"/> 9" to 12" (22.9 to 30.5cm) OR <input type="checkbox"/> 12" to 15" (30.5 to 38.1cm)	
	<input type="checkbox"/> Standard Length Pad (10") OR <input type="checkbox"/> Full Length Pad (14")	

Available Armrest Option

XARM5	<input type="checkbox"/> Covered in Black Ultrasuede (Not available with Desk Arms).....	\$126
-------	--	-------

SEAT BELTS

XBLT1	<input type="checkbox"/> Velcro Adjustable, Black Nylon Webbing.....	\$61
XBLT2	<input type="checkbox"/> Bodypoint® 2" Push-Button Buckle (Auto-Style.).....	\$104
XBLT4	<input type="checkbox"/> Bodypoint® 1" Wide Push-Button Buckle Padded (Auto-Style.) (Not available with 19" - 22" Rear Seat Widths A.).....	\$187
XBLT5	<input type="checkbox"/> Bodypoint® 1.5" Wide Push-Button Buckle Padded (Auto-Style.) (Not available with 19" - 22" Rear Seat Widths A.).....	\$187
XBL10	<input type="checkbox"/> Bodypoint® Optional Band Clamp Mounting Hardware.....	\$22
XBL16	<input type="checkbox"/> Bodypoint® EvoFlex® Pelvic Stabilizer Push Button.....	\$263

TiLite X-Series Canadian Order Form

Effective Sept 2023 (New Items in Red)

FOLDING FRAME

Page 7 of 7

CALF STRAP

XCLF2	Velcro Adjustable.....	\$50
	<input type="checkbox"/> Small (Recommended for 9"- 12" (22.9-30.5cm) Footrest Widths)	
	<input type="checkbox"/> Medium (Recommended for 11"- 15" (27.9-38.1cm) Footrest Widths)	
	<input type="checkbox"/> Large (Recommended for 14"- 17" (35.6-43.2cm) Footrest Widths)	
	<input type="checkbox"/> Xtra Large (Recommended for 16"- 19" (40.6-48.3cm) Footrest Widths)	
XCLF3	Bodypoint® Padded Velcro Adjustable.....	\$88
	<input type="checkbox"/> Medium (Recommended for 12"- 15" (30.5-38.1cm) Footrest Widths)	
	<input type="checkbox"/> Large (Recommended for 16"- 20" (40.6-50.8cm) Footrest Widths)	
XCLF1	<input type="checkbox"/> Buckle Adjustable, Black Nylon Webbing.....	\$34

RESIDUAL LIMB SUPPORT

XRSL1	<input type="checkbox"/> Residual Limb Support (8" width x 10" length).....	\$563
XRSL2	<input type="checkbox"/> Residual Limb Support Pair (8" width x 10" length).....	\$1126

BACKPACK

XBCK3	<input type="checkbox"/> Black Permobil Backpack.....	\$104
-------	---	-------

SEAT POUCH

XPCH3	<input type="checkbox"/> Black Permobil Seat Pouch.....	\$66
XPCH4	<input type="checkbox"/> Bodypoint® Mobility Bag.....	\$143

SPOKE GUARDS

XGRD1	<input type="checkbox"/> Clear with TiMotion Logo (Not available with Mag Rear Wheels.).....	\$148
-------	--	-------

CUSTOM CONFIGURATION OPTIONS

ANTI-TIPS

XTIP1	<input type="checkbox"/> Rear Aluminum - Black.....	\$214
XTIP7	<input type="checkbox"/> Black Aluminum Single.....	\$110
	Mount on <input type="checkbox"/> Left Side or <input type="checkbox"/> Right Side.....	NCO
XTIP2	<input type="checkbox"/> User Friendly Flip-up.....	\$394
XTIP8	<input type="checkbox"/> User Friendly Flip-up Single.....	\$187
	Mount on <input type="checkbox"/> Left Side or <input type="checkbox"/> Right Side.....	NCO
XTIP9	<input type="checkbox"/> Tip Assist.....	\$111

ACCESSORIES

XUEST1	Upper Extremity Support Tray	
	<input type="checkbox"/> 12" to 14" (30.5-35.6cm) OR <input type="checkbox"/> 13" to 15" (33-38.1cm) Seat Widths A (Requires Armrests.).....	\$305
XLUG1	Luggage Carriers	
	(Not available with High Mount Swing Away Footrests or 1-Piece Flip-up Footrest.).....	\$153
XMPCT2	Neoprene Impact Guards with Ultrasuede Stripe (Not available with High Mount Footrests, 90° Hangers or Elevating Leg Rest.)	
	<input type="checkbox"/> Black.....	\$116
XCH1	<input type="checkbox"/> Crutch Holder.....	\$160

GENERAL

More options and frame modifications may be available through our Creative Mobility Division! Please contact customer service for details or [click here](#) to view options.

Permobil manufactures TiLite manual wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional.

Get the latest order forms at permobilus.com. Prices, standard features, options and accessories are subject to change without notice. All rights reserved. Trademarks are identified by symbols™ and ®. TiLite, TiFit, the TiLite logo, Ti, the Ti logo, and TiShaft are registered trademarks of TiSport, LLC and are licensed to Permobil. Permobil® is a registered trademark of Permobil. Adobe®, Acrobat®, Acrobat Professional® and Reader® are registered trademarks of Adobe Systems Inc. Bodypoint® is a registered trademark of Bodypoint Inc. Kenda® is a registered trademark of Kenda Rubber Industrial Co. Ltd. ROHO®, AGILITY®, Low Profile® and Quadro Select® are registered trademarks of ROHO Inc. Schwalbe® is a registered trademark of Ralf Bohle GmbH. The Natural-Fit® and Out-Front® are registered trademarks of Three Rivers Holdings, LLC. The ACTA-BACK®, ACTA-EMBRACE™, COMFORT-TEK™, COMPASS®, CURVE®, ELEMENTS®, and STRETCH AIR™ are registered trademarks of THE COMFORT COMPANY® and are licensed to Permobil. Q Straint® is a registered trademark of Q Straint Inc.

NOTES

QUOTE THIS CHAIR

If you are using Adobe Acrobat, this order form may be converted into a quote request or into an order by clicking one of the buttons (this feature is available exclusively to TiLite dealers and distributors). A TiLite Quote or a TiLite Order Confirmation will be sent to you within 24 hours.

ORDER THIS CHAIR