

# Power Wheelchairs

Independence in motion.



permobil

Front wheel



Mid wheel



Rehab Series	F3 Corpus	F5 Corpus	F5 Corpus VS	M1	M3 Corpus	M5 Corpus	M300 HD
<b>Max user weight</b>	300 lb	300 lb	300 lb	300 lb	300 lb	300 lb	450 lb
<b>Max speed</b>	6 mph   5 mph <i>without anti-tips</i>	7.5 mph   6 mph <i>without anti-tips</i>	7.5 mph	5 mph	6 mph	7.5 mph	5 mph
<b>Range</b>	13 mi (Grp 34s)   20 mi (Grp 24s)*	20 mi*	20 mi*	13 mi*	13 mi (Grp 34s)   20 mi (Grp 24s)*	20 mi*	13 mi (Grp 34s)   20 mi (Grp 24s)*
<b>Base width</b>	24"	25.5"	25.5"	24"	24"	25"	26"
<b>Base length</b>	40" <i>with anti-tips</i> 36" <i>without anti-tips</i>	42" <i>with anti-tips</i> 38" <i>without anti-tips</i>	43"	35.5"	35.5"	37.5"	38"
<b>Base turning radius</b>	26.5"	30"	30"	20.5"	20.5"	21.5"	22"
<b>Ground clearance</b>	3"	3"	3"	3"	3"	3"	3"
<b>Weight w/ batteries</b>	386 lb (Gr 34), 407 lb (Gr 24)	412 lb (Grp 24s)	432 lb (Grp 24s)	319 lb (Gr 34)	385 lb (Grp 34s)   405 lb (Grp 24s)	423 lb (Grp 24s)	425 lb (Gr 34)   447 lb (Gr 24)
<b>Seat-to-floor height</b>	17.5", 18.5", 19.5"	17.5", 18.5", 19.5"	17.5", 18.5", 19.5"	16", 17", 18"	17.5", 18.5", 19.5"	17.5", 18.5", 19.5"	18", 19", 20" (elevator) 18", 19", 20" (fixed)
<b>Seat elevator</b>	<i>Optional 12"</i>	<i>Optional 14"</i>	14"	—	<i>Optional 12"</i>	<i>Optional 14"</i>	<i>Optional 8"</i>
<b>Power tilt options</b>	0° – 50° <i>Posterior</i> 5°, 10°, 20° & 30° <i>Anterior</i>	0° – 50° <i>Posterior</i> 5°, 10°, 20° & 45° <i>Anterior</i>	0° – 50° <i>Posterior</i> 45° <i>Anterior</i>	0° – 50° <i>Posterior</i>	0° – 50° <i>Posterior</i> 5°, 10° & 20° <i>Anterior</i>	0° – 50° <i>Posterior</i> 5°, 10° & 20° <i>Anterior</i>	0° – 45° <i>Posterior</i>
<b>Recline options</b>	85° – 120° (manual) 85° – 180° (power)	85° – 120° (manual) 85° – 180° (power)	85° – 180° (power)	—	85° – 120° (manual) 85° – 180° (power)	85° – 120° (manual) 85° – 180° (power)	85° – 120° (manual) 85° – 150° (power)
<b>Legrest options</b>	85° – 180°	85° – 180°	85° – 180°	85° – 140° 65° – 80° (adjustable swing-away)	85° – 180°	85° – 180°	85° – 175°
<b>Seat depth</b>	14" – 22" (by 1" incr.)	14" – 22" (by 1" incr.)	14" – 22" (by 1" incr.)	14" – 22" (by 1" incr.)	14" – 22" (by 1" incr.)	14" – 22" (by 1" incr.)	19" – 25" (by 1" incr.)
<b>Seat width</b>	17", 19", 21", 23"	17", 19", 21", 23"	17", 19", 21", 23"	14" – 22" (by 2" incr.)	17", 19", 21", 23"	17", 19", 21", 23"	19", 21", 23", 25"
<b>Backrest height</b>	20", 23" – 28" (by 1" incr.)	20", 23" – 28" (by 1" incr.)	23" – 28" (by 1" incr.)	18" – 29"	20", 23" – 28" (by 1" incr.)	20", 23" – 28" (by 1" incr.)	21", 24" – 29" (by 1" incr.)

*All specifications may vary based upon configuration. | Chair weights are approximate and include base, batteries, power seat functions and seating system without accessories. | \*Actual driving range will vary based on driving and battery conditions.*

## Front wheel

## Mid wheel

## Rear wheel



Pediatric Series	Koala	K300 PS Jr.	M300 PS Jr.	K450 MX
<b>Max user weight</b>	80 lb	165 lb	165 lb	125 lb
<b>Max speed</b>	4 mph	5 mph   6.5 mph	6 mph   7.5 mph	5 mph
<b>Range</b>	15 mi*	16 mi*	13 mi (Grp 34s) 20 mi (Grp 24s)*	17 mi*
<b>Base width</b>	23"	24"	24"	25"
<b>Base length</b>	32"	40"	36.5"	35"
<b>Min turning radius</b>	21"	25"	20"	28"
<b>Ground clearance</b>	2.5"	3"	3"	2.5"
<b>Weight w/ batteries</b>	226 lb (M 40)	272 lb (Gr 34)	313 lb (Gr 34) 335 lb (Gr 24)	337 lb (Gr 34)
<b>Seat-to-floor height</b>	15" (elevator)	16.5" (elevator) 16.5", 17.25", 18.5" (fixed)	17" (elevator) 17", 17.75", 18.5" (fixed)	3" - 26" (elevator)
<b>Seat elevator</b>	4"	Optional 8"	Optional 8"	8"
<b>Power tilt options</b>	0° - 25°	0° - 45° (power or manual)	0° - 45° (power or manual)	0° - 45°
<b>Recline options</b>	85° - 120° (manual)	85° - 120° (manual)	85° - 120° (manual)	85° - 105° (manual)
<b>Legrest options</b>	90° - 150° (manual)	70° (swing away) 90 - 150° (center mount)	70° or 90° (swing away) 90° - 150° (center mount)	90° - 125° (manual)
<b>Seat depth</b>	10", 12", 14"	10" - 18"	10" - 18"	10" - 18" (by 1" incr.)
<b>Seat width</b>	10", 12", 14"	11" - 16" (by 1" incr.)	11" - 16" (by 1" incr.)	10", 12", 14", 16"
<b>Backrest height</b>	14", 16", 20"	16.5" - 22.5" (by 1" incr.)	16.5" - 22.5" (by 1" incr.)	14", 16", 20" (MX) 16.5" - 22.5" (PS Jr)

All specifications may vary based upon configuration. | Chair weights are approximate and include base, batteries, power seat functions and seating system without accessories.

\*Actual driving range will vary based on driving and battery conditions.