

PDG Fuze JR - Atlantic Canada Full Feature Order Form



Dealer Name*:	Account #*:	
City, Province, Post Code:		
Phone:	Email:	ATP:
Rep Name:	P.O. Number:	Quote Number:
Client Reference*:	Client Weight:	
Special considerations that need to be addressed here (e.g., diagnosis):		<i>*Required field</i>

Wheelchair & Frame

- Fuze JR Pediatric 50° Manual Tilt Wheelchair, Standard One-Piece Frame** \$ 4,300.00
- Fuze JR Pediatric 50° Manual Tilt Wheelchair, Take-Apart Frame** *The Take-Apart frame allows attendants to take the Fuze apart in seconds without tools. This feature reduces the weight of the heaviest piece of the wheelchair to 29lbs or less, perfect for transporting the chair in the trunk of a car or SUV. Not compatible with power tilt or drum brakes.* \$ 4,800.00

Frame Type

- Adjustable/Growable Width Frame. 250lbs Weight Capacity** No charge

Frame Options

- Transit Tie Down System (TTDS)** *Includes 4 factory installed transport brackets. TTDS has been tested according to ANSI/RESNA WC19. See Owners Manual for complete details.* No charge

Seat Measurements

Seat Width | Overall chair width is +7.5" with 12/16" rear wheels and +8.5" with 20/22/24" rear wheels.

- | | | | | | | | |
|---------------------|-----|-----|-----|-----|-----|--|-----------|
| Adjustable 12"-16": | 12" | 13" | 14" | 15" | 16" | | No charge |
| Adjustable 15"-18": | 15" | 16" | 17" | 18" | | | No charge |

Seat Depth (Adjustable)

- | | | | | | | | | |
|---------------|-----|-----|-----|-----|-----|-----|--|-----------|
| Select Depth: | 13" | 14" | 15" | 16" | 17" | 18" | | No charge |
|---------------|-----|-----|-----|-----|-----|-----|--|-----------|

Frame Finishes

- | | | | | | |
|----------------------|---------------------|---------------------|---------------|---------------------|--------------------|
| Select Color: | High Gloss Black | Red Pepper | Neon Pink | Ocean Blue Metallic | Tangerine Metallic |
| | Acid Green Metallic | Midnight Blue Pearl | White Pearl | Black Lava | Titanium Grey |
| | Pomegranate Red | Blue Mist | Indigo Purple | Hyper Blue | Flame Orange |
| | OD Green | Silver Streak | Lavendar Dusk | Black Galaxy | Seafoam Green |
| | Daffodil Yellow | | | | |

Armrests

Armrest Type

- No Armrests** No charge
 - Single Post, Height Adjustable Armrests** *8"-11" Range* \$ 295.00
- Armrest Position:** Rear (Default) Front

If left blank default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions PDG may not be able to accommodate stated preference.

- Locking Flip-Up Armrests** *Not recommended for users with excessive strength or tone. Not compatible with fold-down back. Only with regular armpads.* \$ 295.00
- Set at 8" Set at 9" Set at 10" Set at 11"

Armpads

Regular Armpads

Desk Length 2" x 10" Full Length 2" x 14"

No charge

Waterfall Armpads

Desk Length 2" x 10" Full Length 2" x 14"

\$ 115.00

Comfort Armpads with Interlocking Hardware *May interfere with wheels in some tilt configurations.*

Desk Length 3.75" x 10" Full Length 3.75" x 14"

\$ 200.00

Gel Armpads *May interfere with wheels in some tilt configurations.*

Narrow 2" x 12" Wide 3.5" x 12"

\$ 135.00

Armrest Wideners | *Moves armrests outwards by up to 1" each side. Seat size does not change.*

Armrest Wideners, 1" Total Additional Width

\$ 95.00

Armrest Wideners, 2" Total Additional Width

\$ 195.00

Seat Height

Chairs with swing-away front rigging are set up with 3° of pre-tilt to avoid interference. This increases the actual seat to floor height at the front of the seat by up to 1" compared to the selected value. Easy-Out and standard plug-in front rigging are not affected.

RW	5" Casters					6" Casters & 6" Forks					6" Casters & 8" Forks					8" Casters & 8" Forks				
	13"	14"	15"	16"	17"	14"	15"	16"	17"	18"	15"	16"	17"	18"	19"	16"	17"	18"	19"	20"
24"																				
22"																				
20"																				
16"																				
12"																				

Seat to Floor Height

Seat Pan

Seat Pan

No Seat Pan

No charge

Solid Aluminium Seat Pan

\$ 135.00

Seat Cushions

Cushions over 20" wide subject to increased pricing.

Comfort M2 Cushion Zero Elevation

\$ 580.00

Comfort-Tek Cover Stretch Air Cover

Comfort M2 Cushion Anti-Thrust

\$ 580.00

Comfort-Tek Cover Stretch Air Cover

Comfort Embrace Cushion Zero Elevation

\$ 508.00

Comfort-Tek Cover Stretch Air Cover

Comfort Embrace Cushion Anti-Thrust

\$ 508.00

Comfort-Tek Cover Stretch Air Cover

ROHO Quadro Select High Profile Cushion

\$ 777.00

ROHO Sensor Ready High Profile Cushion

\$ 637.00

Front Casters

Urethane Casters

5" x 1.25" Caster	No charge
6" x 1.25" Caster	No charge
8" x 1.25" Caster	No charge
6" x 1.50" Caster	\$ 85.00
8" x 2.00" Caster	\$ 135.00

Pneumatic Casters

6" x 1.25" Pneumatic Caster	\$ 110.00
8" x 1.25" Pneumatic Caster	\$ 120.00
8" x 2.00" Pneumatic Caster	\$ 150.00

Rear Wheels

Composite Rear Wheel	No charge
12" 16" 20" 22" 24"	

Spoke Rear Wheel	\$ 75.00
20" 22" 24"	

Rear Wheel Position: Rear (More Stable) Mid (Default) Front (Less Stable)

If left blank, default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

Axles

Quick Release Rear Axles	No charge
Fixed Rear Axles	No charge

Rear Tires

Urethane Tires	No charge
Pneumatic Tires <i>May result in a finished seat height 1/2" greater than stated.</i>	\$ 60.00
Pneumatic Tires with Flat-Free Inserts <i>May result in a finished seat height 1/2" greater than stated. Not available with spoke wheels.</i>	\$ 95.00

Handrims | Handrims are only available with 20", 22" and 24" rear wheels.

No Handrims	No charge
Standard Aluminium Handrims	No charge
Plastic Coated Handrims <i>Not compatible with drum brakes.</i>	\$ 150.00
Projection Handrims <i>Composite wheels only. 10 projections, 30°, 1-1.8" length.</i>	\$ 495.00

Wheel Locks

User Wheel Locks	No charge
Push to Lock Pull to Lock	
Attendant Wheel Locks <i>Not compatible with wheel lock extensions.</i>	\$ 85.00
Push to Lock Pull to Lock	
Foot Operated Wheel Lock <i>Rear-mounted unilateral single-activation lock. Not available with power tilt in combination with 12" or 16" rear wheels. Compatible with 19" and 20" seat to floor heights only.</i>	\$ 595.00
Drum Brakes <i>Available with 12" composite or 20" spoke rear wheels only. Urethane tires only. Not with plastic coated handrims. 12" wheels add 1" to overall chair width. 12" compatible with 13", 14", 17.5", 18.5", 19.5" STF heights only. 20" wheels compatible with 16", 17", 18" STF heights only.</i>	\$ 795.00
User Activated Unilateral Wheel Lock <i>Pull to lock.</i>	\$ 795.00
Left Activation Right Activation	

Wheel Lock Options

Wheel Lock Extensions <i>Not compatible with attendant wheel lock.</i>	\$ 60.00
---	----------

Anti-Tippers

Rear Anti-Tippers <i>Required safety option.</i>	\$ 75.00
---	----------

Backrest

Back Post Style:	Stroller Handle Back	Cane Handle Back			
	<i>Cane handle back not compatible with headrest mounting plate.</i>				
Back Post Type:	Straight	8° Bend			
Back Post Height:	Short (17")	Tall (23")			
	<i>Measured from seat pan to handle.</i>				
Folding Style:	Fold-Down	Non-Folding			
	<i>Fold-down back not compatible with reclining backrest.</i>				
Seat to Back Angle:	Set at 90°	Set at 98°	Set at 106°	Set at 114°	Set at 122°

Backrest Options

Angle/Height Adjustable Stroller Handle:	<i>Not compatible with cane handle back.</i>	\$ 150.00
Button on Outside	Button on Inside	

Back Supports

Comfort Back Supports | Backs over 20" wide or tall subject to increased pricing.

Comfort Unite 1 Back	<i>Includes 2-point quick release Compass hardware. Select shell height and cover.</i>	\$ 961.00
12" 14" 16" 18" 20" 22"		
Comfort-Tek Cover	Diamond Spacer Cover	
Comfort Unite 1 Back - 5" Wing	<i>Includes 2-point quick release Compass hardware. Select shell height and cover.</i>	\$ 1,245.00
12" 14" 16" 18" 20" 22"		
Comfort-Tek Cover	Diamond Spacer Cover	
Comfort Unite 1 Back w/ROHO Insert - 5" Wing	<i>Includes 2-point quick release Compass hardware. Select shell height and cover.</i>	\$ 1,499.00
14" 16" 18" 20" 22"		
Comfort-Tek Cover	Diamond Spacer Cover	
Comfort Unite 1 Back - 7" Wing	<i>Includes 2-point quick release Compass hardware. Select shell height and cover.</i>	\$ 1,245.00
12" 14" 16" 18" 20" 22"		
Comfort-Tek Cover	Diamond Spacer Cover	
Comfort Unite 1 Back w/ROHO Insert - 7" Wing	<i>Includes 2-point quick release Compass hardware. Select shell height and cover.</i>	\$ 1,499.00
14" 16" 18" 20" 22"		
Comfort-Tek Cover	Diamond Spacer Cover	
Comfort Unite 2 Back	<i>Includes 2-point quick release Compass hardware. Select shell height and cover.</i>	\$ 1,106.00
12" 14" 16" 18" 20" 22"		
Comfort-Tek Cover	Diamond Spacer Cover	
Comfort Unite 2 Back w/ROHO Insert	<i>Includes 2-point quick release Compass hardware. Select shell height and cover.</i>	\$ 1,254.00
14" 16" 18" 20" 22"		
Comfort-Tek Cover	Diamond Spacer Cover	
Comfort Unite 3 Back	<i>Includes 2-point quick release Compass hardware. Select shell height and cover.</i>	\$ 1,287.00
14" 16" 18" 20" 22"		
Comfort-Tek Cover	Diamond Spacer Cover	
Comfort Unite 3 Back w/ROHO Insert	<i>Includes 2-point quick release Compass hardware. Select shell height and cover.</i>	\$ 1,499.00
16" 18" 20" 22"		
Comfort-Tek Cover	Diamond Spacer Cover	
Comfort Embrace Back	<i>Includes 2-point quick release Compass hardware. Select shell height and cover.</i>	\$ 1,099.00
10" 12" 14" 16" 18" 20" 22"		
Comfort-Tek Cover	Diamond Spacer Cover	
Comfort Relief Back	<i>Includes Comfort-Tek cover and 2-point quick release Compass hardware. Select shell height.</i>	\$ 1,329.00
10" 12" 14" 16" 18" 20" 22"		

Positioning Accessories

Lateral Supports

BodiLink Lateral Supports *Select size, cover and side.*

3" x 5.5" 4" x 6.5" 6" x 6.5"

Comfort-Tek Cover Stretch Air Cover

Left Right

\$ 220.00

Pair

\$ 440.00

BodiLink Lateral Support Mounting Hardware

Individual

\$ 368.00

Pair

\$ 736.00

Headrests

BodiLink Headrest Pad *Select size and cover.*

\$ 260.00

10" Wide 12" Wide 14" Wide 18" Wide

Comfort-Tek Cover Stretch Air Cover

BodiLink Headrest Hardware, Removable

Mid 6 Mid 8

\$ 373.00

Lat 6 Lat 8

\$ 489.00

Max 8

\$ 545.00

BodiLink Headrest Hardware, Flip-Down

Mid 6 Mid 8

\$ 419.00

Lat 6 Lat 8

\$ 530.00

Max 8

\$ 591.00

Recline

Reclining Backrest (Manual) *30° maximum recline. Not compatible with cane handle back or fold-down back. Requires 23" angle/height adjustable stroller handle with 8° bend. Not available with 15" or 16" seat depth.* \$ 895.00

Tilt Activation

Attendant Tilt Only

Attendant Tilt Only, Right Back Post

No charge

Attendant Tilt Only, Left Back Post

No charge

Power Tilt | 40° tilt range. Includes 2 batteries for off-chair charging.

Power Tilt

\$ 1,795.00

Toggle Switch Control

No charge

Installed on Right Installed on Left

Button Switch Control

No charge

Installed in Right Installed in Left

Stealth Egg Switch Control *Not installed.*

No charge

Swing-Away Mounting Hardware *Only for Stealth egg switch. Not compatible with all headrest makes and models.*

\$ 625.00

Tilt Limit | *Optional. If no tilt limit is selected the full tilt range will be available.*

Tilt Limit:

Set at 20°

Set at 30°

Set at 40°

Front Rigging

No Front Rigging No charge

Easy-Out Front Rigging | Effortless insertion and removal. No hanger brackets remain on the frame.

Fixed 70° \$ 95.00

Universal Adjustable Angle (0°-100°) Select angle. Angle adjustment by loosening one bolt, in 10° increments. Tool-free leg length adjustment. Front rigging may interfere with castors in the 90°-100° settings. \$ 295.00

0° 10° 20° 30° 40° 50° 60° 70° 80° 90° 100°

High Mount Option for Easy-Out Adjustable Angle For shorter lower leg lengths of 8"-12.5". \$ 295.00

Swing Away Front Rigging | Legacy design. Not compatible with Easy-Out.

Swing Away 60° No charge

Swing Away 70° No charge

Swing Away 90° Adjustable up to 3.5" leg length. Includes 11"-16" extension tubes. Not compatible with 8" casters. \$ 75.00

Swing Away Elevating Legrests (ELR) Includes calf pads. Not compatible with 12" and 13" seat widths. \$ 375.00

1" Offset Hangers for Swing Away Front Rigging \$ 75.00

Plug In Front Rigging | Legacy design. Not compatible with Easy-Out.

Plug In 70° No charge

Footplates

Composite Footplates Not available with 12" seat width. No charge

Aluminium Angle Adjustable Footplates \$ 165.00

One Piece Fixed Footplate Available on plug in front rigging only. Requires selection of composite or aluminum angle adjustable footplates. \$ 350.00

One Piece Flip-Up Footplate Requires selection of composite or aluminum angle adjustable footplates. \$ 395.00

Interlocking Angle Adjustable Footplates \$ 195.00

Extension Tubes | Consider the height of your seat cushion.

Standard Extension Tubes, Short No charge

Standard Extension Tubes, Medium Only with 90° swing-away front rigging. No charge

Standard Extension Tubes, Long \$ 75.00

Footplate Height:

Specify:

Ext Tube	Short	8" - 13"	11" - 16"	12" - 16"	4" - 9"	13" - 16"
	Medium	N/A	N/A	N/A	10" - 12"	N/A
	Long	N/A	14" - 22"	14" - 19"	12" - 16"	14" - 17"
		High Mount Easy-Out	Standard Easy-Out	60°, 70°, ELR Swing-Away	90° Swing-Away	Standard Plug-In
		Legrest				

Front Rigging Accessories

Calf Pads, Easy-Out Front Rigging Not compatible with gel padding. Includes reinforced premium calf pad mount. \$ 200.00

Calf Pads, Swing Away Front Rigging Not compatible with gel padding. \$ 135.00

Gel Padding Not compatible with calf pads or Easy-Out front rigging. \$ 125.00

Gel Knee Pads, Easy-Out Angle Adjustable Front Rigging \$ 200.00

Gel Knee Pads, Swing Away ELR Front Rigging \$ 100.00

Length Adjustable Heel Loops \$ 50.00

Standard Calf Strap \$ 50.00

Bodypoint Padded Calf Strap \$ 195.00

Comfort Foot Box, Double Without separator. Select seat width

16 18"

Comfort Foot Box, Single

Left Right

\$ 285.00

Pair

\$ 570.00

Positioning Belts

Bodypoint Non-Padded Center Pull Push Button Belt, Large <i>Not compatible with gel padding. Includes reinforced premium calf pad mount.</i>	\$ 115.00
Bodypoint 2 Point Padded Center Pull Push Button Belt <i>Not compatible with gel padding.</i>	\$ 215.00
Bodypoint 4 Point Padded Center Pull Push Button Belt <i>Not compatible with calf pads or Easy-Out front rigging.</i>	\$ 255.00
Bodypoint Non-Padded Chest Belt 46"	\$ 165.00
Q'Straight Lap Belt <i>Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" belt length.</i>	\$ 245.00

Miscellaneous Options

Spoke Guards <i>Available for 20"-24" composite rear wheels. Made from durable and lightweight aluminium.</i>	\$ 195.00
Black Frame color	
Hip Belt Mount <i>Provides improved mounting options for Bodypoint hip belt or other third-party belts.</i>	\$ 55.00
Clear Lap Tray (Slide On) <i>Available for seat widths 17"-18" only.</i>	\$ 490.00
Medical Necessities Backpack <i>Black with Permobil logo.</i>	\$ 95.00
Adhesive Back Hook Strip <i>For cushion mounting, 2"x30".</i>	\$ 10.00

General

Permobil manufactures PDG wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

Notes