


Wave dimensions

User dimensions (inches)	Z210 small	Z220 medium	Z230 large
Height	up to 46	40–56	50–74
 <p>Key user dimension: height Select the appropriate bath chair by the user's overall height. Choose the size that allows for growth.</p> <p>Important: Make sure that seat depth, width and back height are adequate. The user's weight must not exceed the maximum working load.</p>			
Item dimensions (inches)	Z210 small	Z220 medium	Z230 large
Seat depth	10½	13	17
Seat width	15	15	17½
Back height	23	26½	32
Seat height	2–7	2–7	2–7
Seat height (incl. tub stand)	12½ – 18½	12½ – 18½	12½ – 18½
Seat height (incl. shower stand)	27½	27½	27½
Seat height (incl. tub transfer base) with long leg inserts	21–26 24½–29½	21–26 24½–29½	21–26 24½–29½
Length of mesh area when flat (without calf rest)	35	41	50½
Length of mesh area when flat (with calf rest)	42½	49	61
Base width	18	18	21
Base length (legs down) with legs flat	24 29½	24 29½	24 29½
Tub stand height (legs flat)	8	8	8
Tub stand length (legs flat)	39	39	39
Tub stand width (legs flat)	17½	17½	17½
Shower stand height	22	22	22
Shower stand width	27	27	27
Shower stand length	34	34	34
Tub transfer base footprint	28 x 30	28 x 30	28 x 30
Max. tub height for transfer base with long leg inserts	20 24	20 24	20 24
Maximum working load (lb)	75	150	250
Back and seat angle adjustments (degrees)	0°, 15°, 30°, 45°, 60°, 75°		



All three sizes are available in blue and pink.

permobil

Distributed in Canada by Permobil www.permobil.com
Vancouver Toronto Montreal Tel: 800-265-9830 Fax: 877-636-8944



Retail order form - January 2023

Use dimension chart to select appropriate size.

All Rifton products are made to order and are non-returnable and non-exchangeable.

Conversion kits

The tool-free conversion kit converts your small chair to medium as your child grows.

Kit without calf rest

Blue Z268

Pink Z269

(Kit includes seat, back, cover and chest strap)



Kit with calf rest

Blue Z266

Pink Z267

(Kit includes seat, back, cover and chest strap)



* Required

		Z21 small	Z220 medium	Z230 large
* Chair		Z211 \$854 <input type="checkbox"/>	Z221 \$922 <input type="checkbox"/>	Z231 \$1,055 <input type="checkbox"/>
* Fabric cover	CHOOSE	Blue Z244 N/C <input type="checkbox"/>	Z246 N/C <input type="checkbox"/>	Z248 N/C <input type="checkbox"/>
	Pink Z245 N/C <input type="checkbox"/>	Z247 N/C <input type="checkbox"/>	Z249 N/C <input type="checkbox"/>	
* Chest strap	CHOOSE	With lateral positioning Z212 \$117 <input type="checkbox"/>	Z222 \$135 <input type="checkbox"/>	Z232 \$154 <input type="checkbox"/>
	Without lateral positioning Z213 \$78 <input type="checkbox"/>	Z223 \$97 <input type="checkbox"/>	Z233 \$117 <input type="checkbox"/>	
Head blocks		Z214 \$166 <input type="checkbox"/>	Z214 \$166 <input type="checkbox"/>	Z234 \$166 <input type="checkbox"/>
Leg straps		Z215 \$78 <input type="checkbox"/>	Z225 \$97 <input type="checkbox"/>	Z235 \$117 <input type="checkbox"/>
Calf rest (Calf rest must be chosen at the time of the original purchase. It cannot be added later.)		Z216 \$214 <input type="checkbox"/>	Z226 \$242 <input type="checkbox"/>	Z236 \$269 <input type="checkbox"/>
Tub stand		Z217 \$481 <input type="checkbox"/>	Z217 \$481 <input type="checkbox"/>	Z217 \$481 <input type="checkbox"/>
Shower stand		Z218 \$1,315 <input type="checkbox"/>	Z218 \$1,315 <input type="checkbox"/>	Z218 \$1,315 <input type="checkbox"/>
Tub transfer base (legs required, calf rest recommended)	CHOOSE	Counterclockwise into tub Z227 \$1,630 <input type="checkbox"/>	Z227 \$1,630 <input type="checkbox"/>	Z227 \$1,630 <input type="checkbox"/>
	Clockwise into tub Z228 \$1,630 <input type="checkbox"/>	Z228 \$1,630 <input type="checkbox"/>	Z228 \$1,630 <input type="checkbox"/>	
Legs for tub transfer base	CHOOSE	Standard Z224 \$99 <input type="checkbox"/>	Z224 \$99 <input type="checkbox"/>	Z224 \$99 <input type="checkbox"/>
	Long (tub edge over 20" high) Z238 \$99 <input type="checkbox"/>	Z238 \$99 <input type="checkbox"/>	Z238 \$99 <input type="checkbox"/>	
Additional head blocks		Z214 \$166 <input type="checkbox"/>	Z214 \$166 <input type="checkbox"/>	Z234 \$166 <input type="checkbox"/>
Additional fabric cover		Z219 \$355 <input type="checkbox"/>	Z229 \$372 <input type="checkbox"/>	Z239 \$393 <input type="checkbox"/>
Additional chest strap	CHOOSE	With lateral positioning Z212 \$117 <input type="checkbox"/>	Z222 \$135 <input type="checkbox"/>	Z232 \$154 <input type="checkbox"/>
	Without lateral positioning Z213 \$78 <input type="checkbox"/>	Z223 \$97 <input type="checkbox"/>	Z233 \$117 <input type="checkbox"/>	
Total		Z210 \$ _____	Z220 \$ _____	Z230 \$ _____