

Permobil

# Progeo Exelle Vario

Everyday chair designed with the user in mind



Aluminium



Foldable



Adjustable

permobil

# Compact and elegant with numerous adjustments

Progeo Exelle Vario is a lightweight, folding wheelchair designed for the active user. It's contoured folding frame and the swing-in/out footrests enable easy transfers and storage. A variety of options to support seating and positioning such as the three legrest frame angles can address knee flexion/extension and spasticity management. The adaptability built into the frame makes this a versatile solution for many different users.



## Progeo Exelle Vario

|                            |  |
|----------------------------|--|
| Seat width                 | 33 - 36 - 39 - 42 - 45 - 48 cm   |
| Seat depth                 | 35 - 37,5 - 40 - 42,5 - 45 - 47,5 cm   |
| Seat height                | Front from 41 to 56 cm (adjustable)<br>Rear from 36 to 53 cm (adjustable)                                      |
| Footplate distance         | From 32 to 51 cm (adjustable)  |
| Backrest height            | From 30 to 54 cm (adjustable)  |
| Front frame angle          | 95° - 110° - 120°  |
| Total length               | 90 cm approx. (std. config.)<br>70 cm approx. (without footplates) (std. config.)                              |
| Total width                | Seat width + 17 cm (with 0° camber)  |
| Encumbrance                | 26 cm approx. (with 0° camber)   |
| Backrest angle             | 90° to the seat (with standard rear frame)<br>Adjustable from 81° to 99° to the seat (with Dynamic rear frame) |
| Camber                     | 0° - 1,5° - 3° (with standard rear frame)<br>0° - 1,5° - 2° - 4° (with Dynamic rear frame)                     |
| Weight without rear wheels | Starting from 8,0 kg approx.   |
| Wheelchair weight          | Starting from 10,6 kg approx.  |
| Maximum load               | 125 kg   |

Permobil reserves the right to change specifications and design without prior notice. \*More colour options available.

Personalise your ride\*

- BLACK
- BRIGHT LIME
- BRIGHT FUCHSIA
- BRIGHT BLUE
- BRIGHT SILVER

Actual product colours may vary.