



TiLite Pilot - Order Form



Dealer Name*:

Account #*:

City, State, Zip:

Phone:

Email:

ATP:

Rep Name:

P.O. Number:

Quote Number:

Client Reference*:

Client Weight:

Special considerations that need to be addressed here (e.g., diagnosis):

*Required field

Frame Style

- ☐ TPFS1 **TiLite Pilot INCLUDING Growth Kit** Full adjustability. 120lb weight limit. Seat width 'A' grows 2" and seat depth 'B' grows 3". Price includes one-time growth kit and Transit Tie-Down option. Additional kits sold separately. E1235, E1237 \$ 3,640.00
- ☐ TPFS2 **TiLite Pilot WITHOUT Growth Kit** E1235, E1237 \$ 3,380.00

Transit Options

- ☐ TPWTB **Transport Brackets** Consists of four tie-down brackets. E1022 No charge
- ☐ TPTDO **Transit Tie-Downs** Tested to WC-19. Refer to the Transit Tie-Down Option Brochure for important information regarding transit tie-down options. WC-19 tie-down brackets require the TiLite back and head support. 95lb weight limit. Not available for user weight less than 51lbs. E1023 \$ 271.00

SmartDrive Power Assist System

- ☐ TMAXSD1 **Permobil SmartDrive MX2+** Not available with rear seat heights less than 13.5" or 20" rear wheels. E0986 \$ 7,529.00
- ☐ TMAXSD18 **Switch Control Buttons** No charge
- ☐ TMAXSD8 **SwitchControl Button with Mono Jack Option** No charge
- ☐ TMAXSD19 **Buddy Button** For use with SwitchControl button with mono jack option. No charge
- ☐ TMAXSD4 **Carrying Bag** \$ 84.00
- ☐ TMAXSD10 **Universal Padded Handle** \$ 22.00

Frame Finishes

Frame Color:

- | | | |
|---|--|---|
| ☐ TPFF8 High Gloss Black | ☐ TPFF10 Red Pepper | ☐ TPFF29 Neon Pink |
| ☐ TPFF14 Ocean Blue Metallic | ☐ TPFF15 Tangerine Metallic | ☐ TPFF30 Acid Green Metallic |
| ☐ TPFF17 Midnight Blue Pearl | ☐ TPFF18 White Pearl | ☐ TPFF19 Black Lava |
| ☐ TPFF27 Titanium Grey | ☐ TPFF59 Pomegranate Red | ☐ TPFF62 Blue Mist |
| ☐ TPFF57 Indigo Purple | ☐ TPFF40 Hyper Blue | ☐ TPFF41 Flame Orange |
| ☐ TPFF42 OD Green | ☐ TPFF60 Silver Streak | ☐ TPFF63 Lavendar Dusk |
| ☐ TPFF58 Black Galaxy | ☐ TPFF61 Seafoam Green | ☐ TPFF64 Daffodil Yellow |

Color Anodize Package

Standard:

☐ ZAN01 Black

Standard

Anodized Package:

☐ ZAN12 Gold
 ☐ ZAN15 Green
 ☐ ZAN13 Purple
 ☐ ZAN14 Orange
 ☐ ZAN10 Blue
 ☐ ZAN11 Red

No charge

Includes caster mount, caster fork and caster caps. If no anodize colour is selected, default will be black.

Frame Measurements

Please refer to Reference Diagrams at end of form for more information.

Seat Width 'A' | Measure inside width from side panel to side panel. Seat width grows 2" if selected width is between 8" and 12". 13" width will grow a maximum of 1", 14" will not grow any wider.

TPSW1: ☐ 8" ☐ 9" ☐ 10" ☐ 11" ☐ 12" ☐ 13" ☐ 14"

Seat Depth 'B' | The seat depth can be grown within the selected range.

Seat Depth Range: ☐ 9"-12" ☐ 11"-15"

TPSD1: ☐ 9" ☐ 10" ☐ 11" ☐ 12" ☐ 13" ☐ 14" ☐ 15"

Front Seat Height 'C' | Measure from the floor to the top edge of the front of the seat pan. Adjustable in 0.5" increments.

TPFSH1: ☐ 12" ☐ 12.5" ☐ 13" ☐ 13.5" ☐ 14" ☐ 14.5" ☐ 15" ☐ 15.5" ☐ 16" ☐ 16.5"

Rear Seat Height 'D' | Measure from the floor to the top edge of the seat pan at the intersection of the front of the back support pad and seat pan. Adjustable in 0.5" increments.

TPRSH1: ☐ 12" ☐ 12.5" ☐ 13" ☐ 13.5" ☐ 14" ☐ 14.5" ☐ 15" ☐ 15.5" ☐ 16" ☐ 16.5"

Seat Depth 'B'	20"			22"			24"			
	9"	12.5"-15"	12"-13.5"	0.5"-1.5"						
	10"	12"-15"	12"-13.5"	0.5"-2.25"						
	11"	12"-15"	12"-13.5"	0.5"-3"	13.5"-16.5"	13"-15"	0.5"-2.25"	14"-16.5"	13.5"-16"	0.5"-2.25"
	12"	12"-15"	12"-13.5"	0.5"-3.5"	13.5"-16"	13"-15"	0.5"-3.25"	14"-16.5"	13.5"-16"	0.5"-3.25"
	13"				13"-16"	13"-15"	0.5"-3.5"	13.5"-16.5"	13.5"-16"	0.5"-3.5"
	14"				13"-16"	13"-15"	0.5"-3.5"	13.5"-16.5"	13.5"-16"	0.5"-3.5"
	15"				13"-16"	13"-15"	1.5"-3.5"	13.5"-16.5"	13.5"-15.5"	1.5"-3.5"
		Front Height	Rear Height	COG	Front Height	Rear Height	COG	Front Height	Rear Height	COG
		Measurements								

Footrest

Seat to Footrest 'E': ☐ 3" ☐ 3.5" ☐ 4" ☐ 4.5" ☐ 5" ☐ 5.5" ☐ 6" ☐ 6.5" ☐ 7" ☐ 7.5" ☐ 8" ☐ 8.5"
☐ 9" ☐ 9.5" ☐ 10" ☐ 10.5" ☐ 11" ☐ 11.5" ☐ 12" ☐ 12.5" ☐ 13" ☐ 13.5" ☐ 14" ☐ 14.5"

Measure from the top of the footplate to the front of the seat pan. For angle adjustable footrest, measurement must be at least 2" shorter than front seat height. For flip-back footrest, measurement must be at least 3" shorter than front seat height.

- ☐ TPFT3 **Angle Adjustable Footrest** Standard
- ☐ TPFT24 **Flip-Back Depth and Angle Adjustable Footrest** Not available on 8", 9" or 14" seat widths. K0040 \$ 233.00

Back Supports

Comfort Back Supports

- ☐ SBRS19 **Comfort Unite 1 Back** Replaces Comfort Acta-Back. Available for 12"-14" seat widths. Inc 2-point quick release Compass hardware. Select shell height and cover. E2613 \$ 712
☐ 12" ☐ 14"
☐ Comfort-Tek Cover ☐ Diamond Spacer Cover
- ☐ SBRS22 **Comfort Embrace Back** Replaces Comfort Acta-Embrace Back. Available for 10"-14" seat widths. Inc 2-point quick release Compass hardware. Select shell height. E2613 \$ 700
☐ 10" ☐ 12" ☐ 14"

Comfort Back Options

- ☒ TPCHD1 **Cane Mount Backrest Mounting Hardware** WC-20 certified. Standard
- ☐ TPLAT1 **BodiLink 3"L x 4"D Lateral Trunk Support** Swing away, GlideWear cover, non-padded.
☐ Left E0956, E1028 \$ 423.00
☐ Right E0956, E1028 \$ 423.00
- ☐ TPHDR1 **Comfort BodiLink Head Support w/GlideWear® Cover** Uses piano taper midline hardware. Available on Acta-Back, Acta-Embrace and TiLite rigid back supports. The Comfort Head Support is WC-20 certified on the Comfort Unite and Comfort Embrace, WC-19 certified when used with the TiLite backrest. E0955, E1028 \$ 504.00
☐ 6" Width ☐ 8" Width

Backrest

Back Type

- | | | | |
|--------------------------------|--|-------|-----------|
| ☐ TPSB2 | Contoured Non-Folding Adjustable Height Rigid Back w/Pad 11"-13" height range. | E2613 | No charge |
| ☐ TPSB3 | Aluminum Canes for Aftermarket Back Supports, No Pad Not available with height adjustable fenders. | | No charge |

Push Handles | Height adjustable range is 33"-40" measured from the floor with 24" wheels.

- | | | |
|--------------------------------|---|-----------|
| ☐ TPSB6 | Bolt-On Height Adjustable Removable Push Handle, Single | \$ 196.00 |
| ☐ TPSB7 | Bolt-On Height Adjustable Removable Push Handles, Dual | \$ 239.00 |

Rear Wheel Positioning

Rear Wheel Spacing | Measure from outside of side panel to inside of rear tire at top of wheel. Not all rear wheel spacing is available with all accessories.

- | | | |
|----------------------------------|--|-----------|
| ☐ TPRWS1 | Fixed Rear Wheel Spacing | Standard |
| | ☐ 0.75" ☐ 1" ☐ 1.25" ☐ 1.5" ☐ 1.75" ☐ 2" | |
| ☐ TPCBR12 | Adjustable Rear Wheel Spacing Camber Plugs | Standard |
| ☐ TPCBR9 | Additional Spacers to Increase Rear Wheel Spacing | No charge |

Camber | Black aluminum camber tube standard.

- | | | |
|---------------|---|--|
| Camber Angle: | ☐ 2° (Minimum 1.5" rear wheel spacing) | ☐ 4° (Minimum 1.25" rear wheel spacing) |
| | ☐ 6° (Minimum 1" rear wheel spacing) | ☐ 8° |
| | ☐ 12° | |

Front Wheels

- | | | |
|---------------------------------|---|-----------|
| ☐ TPFW14 | 4"x0.75" Plastic Wheel w/Poly Tire | Standard |
| ☐ TPFW3 | 4"x0.75" Light-Up Micro Caster | \$ 88.00 |
| ☐ TPFW26 | 4"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire | No charge |
| ☐ TPFW27 | 4"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire Black hub. | \$ 217.00 |

Rear Wheels

20" rear wheels only available with seat depth range 9"-12". 22" and 24" rear wheels only available with 11"-15" seat depth range.

- | | | |
|---------------------------------|---|-----------|
| ☐ TPRW14 | TiLite SHADOW Rear Wheel Select size. | No charge |
| | ☐ 20" (451) ☐ 22" (501) ☐ 24" (540) | |
| ☐ TPRW7 | Spinergy SPOX Rear Wheel Select size and spoke color. Orange, pink, purple and green spokes only available in 20" size. | \$ 1,195 |
| | ☐ 20" (451) ☐ 22" (501) ☐ 24" (540) | |
| | ☐ Black ☐ Red ☐ Blue ☐ White ☐ Orange ☐ Pink ☐ Purple ☐ Green | |
| ☐ TPRW8 | Spinergy LX Rear Wheel Select size and spoke color. | \$ 1,295 |
| | ☐ 22" (501) ☐ 24" (540) | |
| | ☐ Black ☐ Red ☐ Blue ☐ White ☐ Orange ☐ Pink ☐ Purple ☐ Green | |
| ☐ TPRW9 | Spinergy Design Lab LX Rear Wheel Black spokes only. Select size and anodizing color. If same color is selected, anodized hubs and rims match colour anodized package. | \$ 1,995 |
| | ☐ 22" (501) ☐ 24" (540) | |
| | ☐ Gold ☐ Green ☐ Blue ☐ Red ☐ Purple ☐ Orange | |
| ☐ TPRW16 | Spinergy Carbon Blade Rear Wheel Select size. 2x2 Twill Weave Carbon Fiber. | \$ 1,795 |
| | ☐ 24" (540) | |
| ☐ TPRW17 | X-Core Mag Rear Wheel Select size. Not with short tabs, quad release axles, spoke guards or omit handrims options. Requires aluminum silver anodized handrims. | \$ 225 |
| | ☐ 20" (451) ☐ 22" (501) ☐ 24" (540) | |

Axles

- | | | |
|---------------------------------|---|-----------|
| ☐ TPAXL1 | Stainless Quick Release Axle | Standard |
| ☐ TPAXL3 | Quad Release Axle Not available with Mag rear wheels. | \$ 169.00 |

Tires

☐	TPRTR1	Treaded Tire 20"(451), 22"(501), 24"(540)		No charge
☐	TPRTR2	Airless Inserts	E2213	\$ 98.00
☐	TPRTR4	MBL SpeedLite Tire 20"(451), 22"(501), 24"(540)		No charge
☐	TPRTR30	Primo Sentinel Tire 20"(451), 22"(501), 24"(540)		\$ 136.00
☐	TPRTR8	Schwalbe Marathon Plus Evolution Tire 24"(540)		\$ 169

Handrims

☐	TPRIM1	Aluminum Handrim, Silver Anodized		Standard
☐	TPRIM21	Aluminum Handrim, Black Anodized	Not available with 20" rear wheels.	No charge
☐	TPRIM20	Q-Grip Handrim	Exclusive high-friction, anti-slip coating. Not available with 20" rear wheels.	\$ 358.00
☐	TPRIM5	Plastic Coated Handrim, Black	Not available with Mag rear wheels.	\$ 174.00
☐	TPRIM11	The Surge LT Handrim	Small profile. All-In-One oval with Gription™ strip. Not available with 20" rear wheels.	\$ 466.00
☐	TPRIM17	The Natural-Fit LT Handrim, Standard Grip	Reduced weight, small profile. Not available with 20" rear wheels, short tabs or Spinergy Carbon Blade rear wheels.	\$ 440.00
☐	TPRIM18	The Natural-Fit LT Handrim, Super Grip	Reduced weight, small profile. Not available with 20" rear wheels, short tabs or Spinergy Carbon Blade rear wheels.	\$ 514.00
☐	TPRIM19	The Natural-Fit LT Handrim, Without Thumb Grip	Reduced weight, small profile. Not available with 20" rear wheels.	\$ 423.00
☐	TPRIM22	Omit Handrims		No charge

Handrim Tab Length

☐	TPRIM7	Long Tabs		Standard
☐	TPRIM8	Short Tabs	Not available with Mag rear wheels or The Natural-Fit handrims with thumb grip.	\$ 28.00

Wheel Locks

☐	TPWLK30	Compact Push to Lock		No charge
		☐ Standard Location 'F' ☐ Forward Location 'G'		
☐	TPWLK31	Compact Pull to Lock		No charge
		☐ Standard Location 'F' ☐ Forward Location 'G'		

Seat Options

☒	TPSP1	Solid Seat Pan	Required.	E2231	\$ 163.00
-------------------------------------	-------	----------------	-----------	-------	-----------

Seat Cushions | Refer to reference charts at end of form for more information.

☐	TPCSH2	Comfort Curve Cushion	Stretch-Air cover.	E2601	\$ 122.00
☐	TPCSH23	Comfort Embrace Comfort-Tek Zero Elevation Cushion		E2607	\$ 428.00
☐	TPCSH11	ROHO Low Profile Single Compartment Cushion		E2622	\$ 534.00

Side Guards

☐	TPSDG30	Built In Side Guard with Height Adjustable Fender	Not available with 20" rear wheels.	E0956, E1034	\$ 487.00
☐	TPSDG31	Built In Side Guard without Fender			Standard

Armrests

☐	TPARM2	Height Adjustable Desk Arm	10" armpad. Works with fendered side guard attachment option.	E0973	\$ 363.00
☐	TPARM31	Height Adjustable Armrest	Attaches to built in side guards.		Standard

Positioning

Upper Body Positioning

☐	TPUBP1	Bodypoint PivotFit Dynamic Shoulder Harness	Cam buckle mounts, rear pull.	E0960	\$ 228.00
		☐ Extra Small ☐ Small			
☐	TPUBP2	Bodypoint Monoflex Chest Support		E0960	\$ 217.00
		☐ Extra Small ☐ Small			

Pelvic Positioning

☐	TPBLT1	Velcro Adjustable Belt, Black Nylon Webbing	E0978	\$ 61.00
☐	TPBLT8	Q'Straint WC-19 Wheelchair Anchored Pelvic Belt <i>Recommended for use with transit tie down.</i>	E0978	\$ 217.00
☐	TPBLT30	Bodypoint 1" Two-Point Padded Hip Belt <i>Swivel buckle.</i>	E0978	\$ 196.00
☐	TPBLT16	Bodypoint 1" Evoflex Belt <i>Small push-button buckle.</i>	E0978	\$ 201.00

Lower Body Positioning

☐	TPCLF15	Bodypoint Ankle Huggers	E0951	\$ 244.00
		☐ Extra Small ☐ Small		

Anti-Tips

☐	TPTIP1	Rear Anti-Tips, Aluminum, Black	E0971	\$ 169.00
☐	TPTIP7	Single Anti-Tip	E0971	\$ 136.00
		☐ Left Side ☐ Right Side		

Accessories


☐	TPGRD30	Clear Polycarbonate Spoke Guards w/TiLite Pilot Sticker Pack <i>Not available with Mag rear wheels.</i>	K0065	\$ 147.00
☐	TPBCK3	Permobil Medical Necessities Backpack		\$ 103.00
☐	TPCH1	Crutch Holder	E2207	\$ 244.00

General

Permobil manufactures TiLite manual wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. HCPCS Codes are for reference purposes only. Actual coding and coverage varies by payer. Proper HCPCS coding is not a guarantee of payment. Please verify the appropriate payer's coverage policy. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

Notes

TiFit Reference Diagrams



Seating Information

Comfort Acta-Embrace Cushion Compatibility

Seat Width 'A'	10"	10", 11", 12", 13", 14", 15"
	11"	10", 11", 12", 13", 14", 15"
	12"	10", 11", 12", 13", 14", 15"
	13"	10", 11", 12", 13", 14", 15"
	14"	10", 11", 12", 13", 14", 15"

Seat Depth 'B'

Roho Configuration Options | *Made to order options have a 2-3 week leadtime.*

Seat Width 'A'	8"	Made to Order	Made to Order	Made to Order	Made to Order	Made to Order	Made to Order
	9"	Made to Order	Made to Order	Made to Order	Made to Order	Made to Order	Made to Order
	10"	Made to Order	Made to Order	Made to Order	Made to Order	Made to Order	Made to Order
	11"	Made to Order	Stocked	Stocked	Stocked	Made to Order	Made to Order
	12"	Made to Order	Stocked	Stocked	Stocked	Made to Order	Made to Order
	13"	Made to Order	Made to Order	Made to Order	Stocked	Stocked	Stocked
	14"	Made to Order	Made to Order	Made to Order	Stocked	Stocked	Stocked
		10"	11"	12"	13"	14"	15"
Seat Depth 'B'							