



TiLite T Series - AADL W395



Dealer Name*:

Account #*:

City, Province, Post Code:

Phone:

Email:

ATP:

Rep Name:

P.O. Number:

Quote Number:

Client Reference*:

Client Weight:

Special considerations that need to be addressed here (e.g., diagnosis):

*Required field

Frame Style

- ☐ ATFS1 **Aero T Aluminum Frame** Full adjustability. 265lb weight limit. \$ 4,036.00

Frame Options

- ☐ THDR1 **Heavy Duty Frame** Not available with TR titanium Frame. 350lb weight limit. Includes tension adjustable by straps back upholstery. Requires folding aluminum adjustable height backrest, billet aluminum front wheels, Shadow or Spinergy SPOX rear wheels. LX rear wheels only available with 0° or 2° camber. Not available with Slipstream fork, velcro adjustable upholstery or Glide forks. No charge
- ☐ TQST1 **Q'Straint Adaptable Frame** Not available with TR titanium Frame. 350lb weight limit. Includes tension adjustable by straps back upholstery. Requires folding aluminum adjustable height backrest, billet aluminum front wheels, Shadow or Spinergy SPOX rear wheels. LX rear wheels only available with 0° or 2° camber. Not available with Slipstream fork, velcro adjustable upholstery, open loop footrest, V-front taper or Glide forks. No charge
- ☐ TPWR1 **Power Adaptable Frame** Requires LiteSpeed billet aluminum front wheels and 0° camber. Not available with single-sided forks or titanium camber tubes. Not required for SmartDrive. No charge
- ☐ TRFR1 **Reinforced Frame** No charge

SmartDrive Power Assist System

- ☐ TMAXSD1 **Permobil SmartDrive MX2+** Not available with rear seat heights less than 13.5" or 20" rear wheels. Maximum user weight capacity of 331 lbs. \$ 7,461.00
- ☐ TMAXSD9 **Battery Pack** Enclosed in SmartDrive system. \$ 528.00
- ☐ TMAXSD18 **Switch Control Buttons** No charge
- ☐ TMAXSD8 **SwitchControl Button with Mono Jack Option** No charge
- ☐ TMAXSD19 **Buddy Button** For use with SwitchControl button with mono jack option. No charge
- ☐ TMAXSD4 **Carrying Bag** \$ 110.00
- ☐ TMAXSD10 **Universal Padded Handle** \$ 31.00

Frame Finishes

- Frame Color:** ☐ TFF8 High Gloss Black ☐ TFF14 Ocean Blue Metallic
- Anodize Package:** ☒ Black Standard

TiFit: Tailored Individualized Frame

Measurements shown below are our most popular selections. Not seeing what you want? Contact customer service or view our Creative Mobility Division (CMD) catalog [here](#) to view additional options. Please refer to Reference Diagrams at end of form for more information.

Rear Seat Width 'A' | Measure from outside of seat tube at back post to the outside of the opposite seat tube at back post.

TSW1: ☐ 12" ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19"

TSW2: ☐ 20"

\$ 541.00

Front Seat Width (Seat Taper) 'T' | Measure from inside of front frame tube to inside of opposite front frame tube. Required with V-front end, not with standard front end.

- ☐ TFSW1 **No Taper** Standard
- ☐ TFSW2 **1" Taper** 3" narrower than Rear Seat Width 'A'. Includes tension adjustable by straps seat upholstery. Bolt on adjustable sling available through CMD #20. \$ 306.00
- ☐ TFSW4 **2" Taper** 4" narrower than Rear Seat Width 'A'. Includes tension adjustable by straps seat upholstery. Bolt on adjustable sling available through CMD #20. \$ 306.00

Seat Depth 'B' | Measure from front of back post to the front edge of seat sling. Frame bend will start 0.5" in front of seat sling unless Custom Frame Depth is specified or the combination of specifications affect the rolling dynamics of the chair.

TSD1: ☐ 12" ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19"

TSD2: ☐ 20"

\$ 590.00

TSD3 Custom Frame Depth: ☐ Specify:

No charge

Measure from front of seat sling to beginning of frame bend. Ranges from 0.5"-2" and can be given as specific as 0.25" increments. A custom frame depth outside the specified range may be available through our Creative Mobility Division. A custom frame depth of 20" will have an additional charge.

Ergonomic Seat 'U' | Measure from front of back post to beginning of upward bend of seat tube.

- ☐ TERG2 **Aero T and TRA Ergonomic Seat** Not available with V-front or seat taper, seat back that locks when folded, desk armrests, depth adjustable backrests, fold down rigid sideguards, fendered adjustable sideguards or Curve cushion. Bolt on adjustable sling available through CMD #20. Solid seat available through CMD #16. Minimum seat depth 12" for 6"-7". Minimum seat depth 13" for 8". \$ 568.00
- ☐ 6" ☐ 7" ☐ 8"

Front Seat Height 'C' | Measure from floor to top of seat tube at front edge of seat sling.

TFSH1: ☐ 15" ☐ 15.5" ☐ 16" ☐ 16.5" ☐ 17" ☐ 17.5" ☐ 18" ☐ 18.5" ☐ 19" ☐ 19.5" ☐ 20" ☐ 20.5" ☐ 21"

Rear Seat Height 'D' | Measure from floor to top of seat tube at back post.

TRSH1: ☐ 12.5" ☐ 13" ☐ 13.5" ☐ 14" ☐ 14.5" ☐ 15" ☐ 15.5" ☐ 16" ☐ 16.5" ☐ 17" ☐ 17.5" ☐ 18" ☐ 18.5" ☐ 19" ☐ 19.5" ☐ 20"

Rear Wheel	26"	15.5"-20"
	25"	15"-20"
	24"	14.5"-20" (15" on Aero T/TRA)
	22"	13.5"-20"
	20"	12.5"-20"

Rear Seat Height Range

Front Angle or Overall Frame Length | Select one. If Overall Frame Length is selected, front angle will be bent to achieve requested length.

- ☐ TFA1 **Frame Angle 'L'** Measure from frame front to floor. 85° not available with 6" front wheels. 90° not available with 5" or 6" front wheels. Standard
- ☐ 70° ☐ 75° ☐ 80° ☐ 85° ☐ 90°
- ☐ TFA2 **Overall Frame Length 'M'** Measure from front of backrest to front edge of footrest. No charge
- ☐ Specify:

Footrest

Standard Footrests

Seat to Footrest 'E': ☐ 13.5" ☐ 14" ☐ 14.5" ☐ 15" ☐ 15.5" ☐ 16" ☐ 16.5" ☐ 17" ☐ 17.5" ☐ 18" ☐ 18.5" ☐ 19"

Measure from front edge of seat sling to top rear of footrest. Measurement must be at least 2" shorter than Front Seat Height 'C'.

- ☐ TFTR3 **Angle Adjustable Footrest** Standard
- ☐ TFTR1 **Titanium Footrest w/Flat ABS Cover** \$ 219.00
- ☐ TFTR2 **Titanium Open Loop Footrest** \$ 279.00
- ☐ TFTR12 **Angle Adjustable Flip-Back Footrest** Not available with 10" or 11" rear seat width or heavy duty frame. \$ 257.00

High Mount Footrests

Seat to Footrest 'E': ☐ 7"-9" ☐ 9"-11" ☐ 11"-13" ☐ Set at

Measure from front edge of seat sling to top rear of footrest.

- ☐ TFTR4 **Angle Adjustable High Mount Footrest** \$ 160.00
- ☐ CMD040 **Custom Two-Piece Flip-Up Footplates** Please visit our CMD catalog and specify footrest type and rigidizer bar placement in Notes section of this order form. Some limitations may apply. CMD

Footrest Width | Select standard front end 'I' if no seat taper 'T' has been selected. Select V-front end 'V' if a seat taper 'T' has been selected.

- ☐ TFE1 **No V-Front Taper 'I'** Find your rear seat width 'A' in the chart below and select one of the options to the right. Measure from inside of front frame tube to inside of opposite front frame tube. Standard

Rear Seat Width 'A'	12"	☐ 10"	☐ 8"	
	13"	☐ 11"	☐ 9"	
	14"	☐ 12"	☐ 10"	☐ 8"
	15"	☐ 13"	☐ 11"	☐ 9"
	16"	☐ 14"	☐ 12"	☐ 10"
	17"	☐ 15"	☐ 13"	☐ 11"
	18"	☐ 16"	☐ 14"	☐ 12"
	19"	☐ 17"	☐ 15"	☐ 13"
	20"	☐ 18"	☐ 16"	☐ 14"
		No Taper	2" Taper	4" Taper
Footrest Width 'I'				

- ☐ TFE2 **V-Front Taper 'V'** Measure from inside of front frame tube to inside of opposite front frame tube 2.5" above footrest. Titanium footrests have 1" of adjustment up and down. Not available with ergonomic seat frames. No charge

☐ 9" ☐ 10" ☐ 11" ☐ 12" ☐ 13" ☐ 14" ☐ 15"

- ☐ CMD008 **Custom Abducted Front Frame Width** Adapted footrest and/or front seat width. Please specify desired front seat width and desired footrest width to the nearest 0.5" in the Notes section. Some limitations may apply. CMD

Back Supports

ROHO Back Supports

- ☐ SBRS14 **Agility Active Aluminum, Min Contour with Foam** Includes fixed hardware. Select shell height. \$ 979.00
- ☐ SBRS15 **Agility Active Aluminum, Mid Contour with Foam** Includes fixed hardware. Select shell height. \$ 979.00
- ☐ 8" ☐ 10"
- ☐ SBRS7 **Agility Carbon Min Contour Back with Foam/Air Insert** Includes fixed hardware. Select shell height. \$ 1,253.00
- ☐ ACB1 **Agility Carbon Min Upgrade** Required with Agility Carbon Backs. Not available with bolt-on push handles. \$ 234.00
- ☐ SBRS9 **Agility Carbon Mid Contour Back with Foam/Air Insert** Includes fixed hardware. Select shell height. \$ 1,278.00
- ☐ ACB2 **Agility Carbon Mid Upgrade** Required with Agility Carbon Backs. Not available with bolt-on push handles. \$ 221.00
- ☐ 10" ☐ 13" ☐ 16"

Comfort Back Supports

- ☐ SBRS22 **Comfort Embrace Back** Replaces Comfort Acta-Embrace Back. Includes 2-point quick release Compass hardware. Select shell height. \$ 1,099
- ☐ 10" ☐ 12" ☐ 14" ☐ 16"
- ☐ Comfort-Tek Cover ☐ Diamond Spacer Cover
- ☐ Privacy Shield \$ 70
- ☐ SBRS23 **Comfort Relief Back** Replaces Comfort Acta-Relief Back. Includes 2-point quick release Compass hardware. Select shell height. \$ 1,329
- ☐ 10" ☐ 12" ☐ 14" ☐ 16"
- ☐ Comfort-Tek Cover ☐ Diamond Spacer Cover
- ☐ Privacy Shield \$ 70

Back Upholstery

- ☐ TBUP2 **Velcro Adjustable** Not available with fold-down push handles. Standard
- ☐ TBUP3 **Tension Adjustable By Straps, Padded Nylon** Required for fold-down push handles. \$ 454.00
- ☐ TBUP5 **Tension Adjustable By Straps, AIR Breathable** Required for fold-down push handles. \$ 437.00
- ☐ TUPC3 **Ultrasuede Backrest Stripe, Black** Not available with AIR breathable upholstery. \$ 116.00
- ☐ TBUP0 **Omit Back Upholstery** No charge

Backrest

Back Type

- ☐ TSB1 **Folding Aluminum Adjustable Height Back** *Required for heavy duty frame.* Standard
- ☐ TSB2 **Folding Titanium Adjustable Height Back, Rounded Rigidizer Bar** *Titanium back will match frame finish.* \$ 323.00

Back Options

- ☐ TSB11 **Backrest Locks When Folded** *Available on all backs. Not available with ergonomic seat.* Standard
- ☐ TSB12 **4" Deep Backrest Rigidizer Bar** *Not a rounded style rigidizer bar. Not available with non-folding, fixed height titanium back type.* No charge
- ☐ TSB3 **TiShaft Back Release Bar, Aluminum** *Not available with non-folding or fixed height titanium backs.* Standard
- ☐ TSB5 **TiShaft Back Release Bar, 2x2 Twill Weave Carbon Fiber** *Not available with non-folding or fixed height titanium backs.* \$ 94.00

Push Handles

- ☐ TSB6 **Integrated Push Handles** No charge
- ☐ TSB7 **Integrated Fold-Down Push Handles** *Requires tension adjustable by straps back upholstery.* \$ 126.00
- ☐ TSB8 **Bolt-On Push Handles** *Includes 4" deep backrest rigidizer bar.* \$ 219.00
- ☐ TSB9 **Bolt-On Height-Adjustable Push Handles** *Includes 4" deep backrest rigidizer bar.* \$ 449.00

Folding Back Height 'F' | TSB1/TSB2. Set at 0.5" increments.

Push Handles	None / Bolt-On	☐ 8.5"-11.5"	☐ 12"-16.5"	☐ 16"-20.5"	Set at:
	Integrated	☐ 10.5"-12.5"	☐ 13.5"-16.5"	☐ 17.5"-20.5"	Set at:
	Fold-Down	☐ 13"-15"	☐ 14.5"-18.5"	☐ 18.5"-21"	Set at:
		Short	Medium	Tall	Specify
Folding Back Height					

Fixed Back Height 'F'

- Select height: ☐ 5" ☐ 5.5" ☐ 6" ☐ 6.5" ☐ 7" ☐ 7.5" ☐ 8" ☐ 8.5" ☐ 9" ☐ 9.5" ☐ 10" ☐ 10.5"
☐ 11" ☐ 11.5" ☐ 12" ☐ 12.5" ☐ 13" ☐ 13.5" ☐ 14" ☐ 14.5" ☐ 15" ☐ 15.5" ☐ 16" ☐ 16.5"
☐ 17" ☐ 17.5" ☐ 18" ☐ 18.5" ☐ 19" ☐ 19.5" ☐ 20"
- Required for fixed side guards. Not available with swing-away tubular or swing-away flip-back tubular armrests, integrated push handles, integrated fold-down push handles or bolt-on height adjustable push handles.*

Back Angle 'G' | TSBA1 Measure angle from front of back post to floor.

- Folding Back: ☐ 80° ☐ 83° ☐ 86° ☐ 89° ☐ 92° ☐ 95° ☐ 98° ☐ 101°
- Fixed Back: ☐ 80° ☐ 81° ☐ 82° ☐ 83° ☐ 84° ☐ 85° ☐ 86° ☐ 87° ☐ 88° ☐ 89° ☐ 90° ☐ 91° ☐ 92°
☐ 93° ☐ 94° ☐ 95° ☐ 96° ☐ 97° ☐ 98° ☐ 99° ☐ 100° ☐ 101°

Rear Wheel Positioning

Center of Gravity 'H' | Measure from front of back post at the seat tube to the center of the rear axle. If no center of gravity is selected the chair will ship with the center of gravity in the most rearward position. Ergo seat selection may limit available center of gravity.

- TCOG1: ☐ 1" ☐ 1.25" ☐ 1.5" ☐ 1.75" ☐ 2" ☐ 2.25" ☐ 2.5" ☐ 2.75" ☐ 3" ☐ 3.25" ☐ 3.5" ☐ 3.75"
☐ 4" ☐ 4.25" ☐ 4.5" ☐ 4.75" ☐ 5"

Rear Wheel Spacing 'J' | Measure from outside of back post to inside of rear tire. TRA 0.75" rear wheel spacing not with 0° or 2° camber.

- ☐ TRWS1 **Fixed Rear Wheel Spacing** Standard
- ☐ TCBR12 **Adjustable Rear Wheel Spacing Camber Plugs** No charge
- ☐ TCBR9 **Additional Spacers to Increase Rear Wheel Spacing** No charge
- ☐ TCBR11 **Zero Play Aluminum Camber Plugs** *Not available with power adaptable frame or more than 6° camber.* No charge

Camber 'K'

- TCBR1: ☐ 0° ☐ 2° ☐ 4° ☐ 6°
8° camber may be available through our Creative Mobility Division
- ☐ TCBR15 **Aluminum Camber Tube** Standard
- ☐ TCBR3 **Titanium Camber Tube** *Not available with power adaptable frame.* \$ 236.00

Seat Upholstery

☐ TSUP1	Tension Adjustable Bolt-On Upholstery, Black <i>Not available with ergonomic seat or tapered seat.</i>	Standard
☐ TSUP2	Tension Adjustable By Straps, Black <i>Standard with ergonomic seat and seat taper of 1"-3".</i>	\$ 203.00
☐ TSUP8	Aluminum Solid Seat <i>Ergonomic seat requires CMD016. Not with 10" or 11" rear seat width.</i>	\$ 312.00
☐ TSUP0	Omit Seat Upholstery	No charge
☐ CMD018	Custom iSleeve <i>Please provide sleeve width and depth (inches) in Notes section at the end of this order form. Add 1" to device depth and 0.25" to width. Some limitations may apply.</i>	CMD

Seat Cushions

☐ TCSH2	Comfort Curve Cushion <i>Stretch-Air cover.</i>	\$ 150.00
☐ TCSH23	Comfort Acta-Embrace™ Comfort-Tek Zero Elevation Cushion	\$ 522.00
☐ TCSH11	ROHO Single Compartment Cushion <i>Fits seat widths 14"-19".</i>	\$ 711.00
	☐ Low Profile ☐ Mid Profile ☐ High Profile	
☐ TCSH14	ROHO Quadtro Select Cushion <i>Fits seat widths 14"-19".</i>	\$ 711.00
	☐ Low Profile ☐ Mid Profile ☐ High Profile	
☐ TCSH16	ROHO Hybrid Select Cushion <i>Fits seat widths 14"-18", and 20"</i>	\$ 756.00

Front Wheels

☐ TFW24	3"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire	\$ 197.00
	☐ Silver Hub ☐ Black Hub	
☐ TFW14	4"x0.75" Plastic Wheel w/Poly Tire <i>Not available with Slipstream fork or power adaptable frame.</i>	No charge
☐ TFW3	4"x0.75" Light-Up Micro Caster <i>Not available with Slipstream fork or power adaptable frame.</i>	\$ 88.00
☐ TFW26	4"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire <i>Not available with power assist, power add-on system or heavy duty frame.</i>	\$ 104.00
☐ TFW27	4"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire	\$ 208.00
	☐ Silver Hub ☐ Black Hub	
☐ TFW6	5"x1" Plastic Wheel w/Poly Tire <i>Not available with 90° front angle, Slipstream fork or power adaptable frame.</i>	No charge
☐ TFW28	5"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire <i>Not available with 90° front angle or power adaptable frame.</i>	\$ 104.00
☐ TFW30	5"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire <i>Not with 90° front angle.</i>	\$ 274.00
	☐ Silver Hub ☐ Black Hub	
☐ TFW9	6"x1" Plastic Wheel w/Poly Tire <i>Not with 85° or 90° front angle, Slipstream fork or power adaptable frame.</i>	No charge
☐ TFW32	6"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire <i>Not with 85° or 90° front angle, Slipstream fork or power adaptable frame.</i>	\$ 110.00
☐ TFW33	6"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire <i>Not with 85° or 90° front angle or Slipstream fork.</i>	\$ 274.00
	☐ Silver Hub ☐ Black Hub	

Precision Lock Front Forks

☐ TFK1	TiLite Standard Fork	Standard
☐ TFK2	TiLite Slipstream Single-Sided Fork <i>Not available with power adaptable frame.</i>	\$ 203.00
☐ TGL1	Glide Suspension Fork <i>Black anodize only. Select rider weight.</i>	\$ 596.00
	☐ 125lbs ☐ 125-180lbs ☐ 180-265lbs	
☐ TTFS1	Titanium Fork Stem	\$ 66.00

Rear Wheels

See rear seat height ranges chart for rear seat height 'D' limitations.

☐ TRW14	TiLite SHADOW Rear Wheel Select size. 25" and 26" wheels require 15" minimum seat depth.	Standard
	☐ 20" (451) ☐ 22" (501) ☐ 24" (540) ☐ 25" (559) ☐ 26" (590)	
☐ TRW7	Spinergy SPOX Rear Wheel Select size and spoke color. 25" and 26" wheels require 15" minimum seat depth. Orange, pink, purple and green spokes only available in 20" size.	\$ 1,427
	☐ 20" (451) ☐ 22" (501) ☐ 24" (540) ☐ 25" (559) ☐ 26" (590)	
	☐ Black ☐ Red ☐ Blue ☐ White ☐ Orange ☐ Pink ☐ Purple ☐ Green	
☐ TRW8	Spinergy LX Rear Wheel Select size and spoke color. 25" and 26" wheels require 15" minimum seat depth.	\$ 1,547
	☐ 22" (501) ☐ 24" (540) ☐ 25" (559) ☐ 26" (590)	
	☐ Black ☐ Red ☐ Blue ☐ White ☐ Orange ☐ Pink ☐ Purple ☐ Green	
☐ TRW9	Spinergy Design Lab LX Rear Wheel Black spokes only. Select size and anodizing color. If same color is selected, anodized hubs and rims match colour anodized package. 25" and 26" wheels require 15" minimum seat depth.	\$ 2,382
	☐ 22" (501) ☐ 24" (540) ☐ 25" (559) ☐ 26" (590)	
	☐ Gold ☐ Green ☐ Blue ☐ Red ☐ Purple ☐ Orange	
☐ TRW16	Spinergy Carbon Blade Rear Wheel 2x2 twill weave carbon fiber. Not available with heavy duty frame. 25" wheels require 15" minimum seat depth.	\$ 2,003
	☐ 24" (540) ☐ 25" (559)	
☐ TRW17	X-Core Mag Rear Wheel Select size. Not available with heavy duty frame, short tabs, quad release axles, spoke guards or omit handrims options. Requires aluminum silver anodized handrims.	\$ 269
	☐ 20" (451) ☐ 22" (501) ☐ 24" (540)	

Axles

☐ TAXL1	Stainless Quick Release Axle	Standard
☐ TAXL2	Titanium Quick Release Axle	\$ 318.00
☐ TAXL3	Quad Release Not available with Mag or 26" (590) Spinergy Spox or LX rear wheels.	\$ 187.00

Pneumatic Rear Tires

☐ TRTR1	Treaded Tire 22"(501), 24"(540), 26"(590)	No charge
☐ TRTR2	Airless Inserts Not available with Spinergy Carbon Blade rear wheels.	\$ 99.00
☐ TRTR31	MBL SpeedLite Tire 20"(451), 22"(501), 24"(540), 25"(559), 26"(590)	No charge
☐ TRTR30	Primo Sentinel Tire 24"(540), 25"(559), 26"(590)	\$ 175.00
☐ TRTR8	Schwalbe Marathon Plus Evolution Tire 24"(540), 25"(559), 26"(590)	\$ 202
☐ TRTR32	Primo C-764 Knobby Tire 24"(540)	\$ 106.00

Solid Rear Tires

☐ TRTR22	Primo Orion Tire 20"(451), 22"(501), 24"(540), 25"(559) Not available with Spinergy Carbon Blade rear wheels.	\$ 121.00
☐ TRTR23	Primo Active Tire 24"(540), 25"(559) Not available with Spinergy Carbon Blade rear wheels.	\$ 148.00

Handrims

☐ TRIM1	Aluminum Handrim, Silver Anodized	Standard
☐ TRIM21	Aluminum Handrim, Black Anodized	No charge
☐ TRIM12	Titanium Handrim, Satin Finish Not available with 20"(451) or 22"(501) rear wheels or 24"(540) knobby tires.	\$ 536.00
☐ TRIM20	Q-Grip Handrim Exclusive high-friction, anti-slip coating. Not available with 20"(451) rear wheels.	\$ 410.00
☐ TRIM5	Plastic Coated Handrim, Black	\$ 175.00
☐ TRIM10	The Surge Handrim Large profile. All-In-One oval with Gription strip.	\$ 476.00
☐ TRIM11	The Surge LT Handrim Small profile. All-In-One oval with Gription strip.	\$ 514.00
☐ TRIM14	The Natural-Fit Handrim, Standard Grip Original, large profile. Not available with short tabs or Spinergy Carbon Blade rear wheels.	\$ 454.00
☐ TRIM15	The Natural-Fit Handrim, Super Grip Original, large profile. Not available with short tabs or Spinergy Carbon Blade rear wheels.	\$ 514.00
☐ TRIM16	The Natural-Fit Handrim, Without Thumb Grip Original, large profile.	\$ 415.00
☐ TRIM17	The Natural-Fit LT Handrim, Standard Grip Reduced weight, small profile. Not available with short tabs or Spinergy Carbon Blade rear wheels.	\$ 546.00
☐ TRIM18	The Natural-Fit LT Handrim, Super Grip Reduced weight, small profile. Not available with short tabs or Spinergy Carbon Blade rear wheels.	\$ 585.00
☐ TRIM19	The Natural-Fit LT Handrim, Without Thumb Grip Reduced weight, small profile.	\$ 503.00
☐ TRIM22	Omit Handrims	No charge

Handrim Tab Length

☐	TRIM7	Long Tabs	Standard
☐	TRIM8	Short Tabs <i>Not available with Mag rear wheels, The Natural-Fit handrims with thumb grip or Natural-Fit handrims on Shadow wheels.</i>	\$ 39.00

Wheel Locks

☐	TWLK7	Compact Composite Push to Lock <i>Available with 0.75"-1.5" rear wheel spacing only.</i>	No charge
☐	TWLK8	Compact Composite Pull to Lock <i>Available with 0.75"-1.5" rear wheel spacing only.</i>	No charge
☐	TWLK10	Extension Handles for Compact Composite Push/Pull to Lock	\$ 94.00
☐	TWLK11	MBL Composite Scissor Lock	\$ 225.00
☐	TWLK12	MBL Performance Aluminum Scissor Lock	\$ 245
☐	TWLK3	Uni-Lock	\$ 181.00
		☐ Push ☐ Pull	
☐	TWLK4	Extension Handles for Uni-Lock	\$ 94.00
☐	TWLK15	Grade Aids Push to Lock <i>Not available with solid rear tires.</i>	\$ 252.00

Side Guards

☐	CMD007	Custom Offset Profiled Side Guards (Removable or Fold-Down) <i>Please specify desired offset to nearest 0.5" and side guards type (rigid or fold-down) in the Notes section of this order form. Side guard material must be selected. Some limitations may apply.</i>	CMD
☐	TSDG13	Side Guards Custom Cut to Profile of Rear Wheels <i>Available for aluminum and carbon non-fendered side guards.</i>	\$ 345.00

Aluminum

☐	TSDG2	Aluminum Side Guards, Rigid Removable	\$ 318.00
		☐ Small ☐ Standard	
☐	TSDG3	Aluminum Side Guards, Quad Removable w/Hole	\$ 345.00
		☐ Small ☐ Standard	
☐	TSDG4	Aluminum Fold Down Side Guards, Rigid <i>Not available with ergonomic seats.</i>	\$ 410.00
		☐ Small ☐ Standard	
☐	TSDG6	Aluminum Fendered Side Guards, Adjustable <i>Available with 24"-25" rear wheels. Requires rear seat height 14.5"-18.5"</i>	\$ 437.00

2x2 Twill Weave Carbon Fiber

☐	TSDG7	Carbon Side Guards, Rigid Removable	\$ 372.00
		☐ Small ☐ Standard	
☐	TSDG8	Carbon Side Guards, Quad Removable w/Hole	\$ 410.00
		☐ Small ☐ Standard	
☐	TSDG9	Carbon Fold Down Side Guards, Rigid <i>Not available with ergonomic seats.</i>	\$ 558.00
		☐ Small ☐ Standard	
☐	TSDG12	Carbon Fendered Side Guards, Removable	\$ 650.00

Armrests

☐	TARM1	Swing Away Tubular Armrests <i>Select range. Not available with non-folding fixed height titanium back or height adjustable push handles.</i>	\$ 284.00
		☐ 6.5"-9.5" ☐ 9.5"-12.5"	
☐	TARM4	Swing Away Flip Back Tubular Armrests, 8.5"-11.5" <i>Not available with non-folding, fixed height titanium back.</i>	\$ 345.00
☐	TARM2	Removable Height Adjustable Desk Arm with Rigid Side Guard <i>Select range and pad size. Requires 14" minimum seat depth. Not available with ergonomic seat.</i>	\$ 394.00
		☐ 9"-12" ☐ 12"-15"	
		☐ Desk Length (10") ☐ Full Length (14")	
☐	TARM3	Removable Height Adjustable Desk Arm with Rigid Side Guard and Transfer Loop <i>Select range and pad size. Requires 15" minimum seat depth. Not available with ergonomic seat.</i>	\$ 410.00
		☐ 9"-12" ☐ 12"-15"	
		☐ Desk Length (10") ☐ Full Length (14")	
☐	TARM5	Armrests Covered in Black Ultrasuede <i>Not available with desk arms.</i>	\$ 126.00

Positioning Belts

☐	TBLT1	Velcro Adjustable Belt, Black Nylon Webbing	\$ 61.00
☐	TBLT2	Bodypoint 2" Belt, Push-Button Buckle (Auto Style)	\$ 104.00
☐	TBLT4	Bodypoint 1" Padded Belt, Push-Button Buckle (Auto Style) <i>Not available with 19" or 20" rear seat widths.</i>	\$ 187.00
☐	TBLT5	Bodypoint 1.5" Padded Belt, Push-Button Buckle (Auto Style) <i>Not available with 19" or 20" rear seat widths.</i>	\$ 187.00
☐	TBLT10	Bodypoint Optional Band Clamp Mounting Hardware	\$ 22.00
☐	TBLT16	Bodypoint EvoFlex Pelvic Stabilizer, Push Button	\$ 263.00

Calf Straps

☐	TCLF2	Velcro Adjustable	\$ 50.00
		☐ Small (9"-12" footrest widths) ☐ Medium (11"-15" footrest widths) ☐ Large (14"-17" footrest widths) ☐ Xtra Large (16"-19" footrest widths)	
☐	TCLF3	Bodypoint Padded Velcro Adjustable	\$ 88.00
		☐ Medium (11"-15" footrest widths) ☐ Large (16"-19" footrest widths)	
☐	TCLF1	Buckle Adjustable, Black Nylon Webbing	\$ 34.00

Residual Limb Supports

☐	TRSL1	Residual Limb Support, 8"W x 10"L, Individual	\$ 563.00
		☐ Left ☐ Right	
☐	TRSL2	Residual Limb Support, 8"W x 10"L, Pair	\$ 1,126.00

Anti-Tips

☐	TTIP1	Rear Anti-Tips, Aluminum Black	\$ 214.00
☐	TTIP2	User Friendly Anti-Tips, Flip-Up	\$ 394.00

Accessories

☐	TGRD1	Clear Lexan Spoke Guards <i>Not available with Mag rear wheels.</i>	\$ 148.00
☐	TBCK3	Permobil Medical Necessities Backpack	\$ 104.00
☐	TPCH3	Permobil Medical Necessities Seat Pouch	\$ 66.00
☐	TPCH4	Bodypoint Mobility Bag	\$ 143.00
☐	TUEST1	Upper Extremity Support Tray <i>Requires armrests.</i>	\$ 318.00
		☐ 12"-14" seat widths ☐ 13"-15" seat widths	
☐	TLUG1	Luggage Carriers	\$ 160.00
☐	TMPCT2	Black Neoprene Impact Guards with Ultrasuede Stripe	\$ 116.00
☐	TCH1	Crutch Holder	\$ 160.00








General

Permobil manufactures TiLite manual wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. HCPCS Codes are for reference purposes only. Actual coding and coverage varies by payer. Proper HCPCS coding is not a guarantee of payment. Please verify the appropriate payer's coverage policy. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.



Notes

TiFit Reference Diagrams

Note: The rear frame for a TR will be square when the centre of gravity 'H' measurement is less than 3.5". Aero T and TRA will always be a rounded rear frame.



Aero T and TRA will always have a rounded rear frame.



TR rear frame will be square when the center of gravity measurement is less than 3.5".

