permobil

Permobil Mattresses, Therapeutic Surfaces & Beds Guide

Retail Price List







About our offering

Permobil has mattresses, therapeutic surfaces, and beds to meet an even wider variety of client needs including different widths, lengths, and features!

Our newest offering is a line of powered active air surfaces for complex needs. Clients benefit from the Alternating Pressure and Pulsation Therapy systems with Low Air Loss, but with an added transfer edge for improving safety and reducing entrapment risk. We also have a new surface called the Atlas Hybrid that can be used with/without an electrical pump. Last but not least, we've expanded our Permobil Modular Systems with a 42" wide Adapt 2 and now a paediatric size for a crib. These customizable mattresses provide the optimal experience for each client.

Any of these mattresses or surfaces can be placed on our popular Tröst bed for long-term care clients or Hälsa bed for home care. We even have new personal beds for home with the Resten Hi-Lo and Select models. You can trust in our bed solutions to provide safety and comfort! For more information, please contact your local sales representative or our Customer Experience team.

ROHO® mattress section

We have expanded our line of ROHO compatible mattresses to give you even more options. These clinically proven, non-powered, and adjustable pressure relieving mattress sections can be inserted into our: Adapt and Adapt 2 Series, Acute, and leveling pads. These mattress sections are flame and fluid resistant and can

be used in a number of configurations whatever best suits your clients' needs!



To Place Your Order

To place your order, please contact our Customer Experience team at sales.ca@permobil.com or 800-265-9830.


To learn more, visit our website at www.permobil.com.

| Di la |
|--|
| |
| The state of the s |
| |
| |
| Table of Contents |
| Mattresses & therapeutic surfaces matrix |
| Key terms4 |
| MATTRESSES 5 Minds Mobility Adapt6 |
| 5 Minds Mobility Allay6 |
| 5 Minds Mobility Align & Align Bariatric |
| Permobil Optima 2 & Symmetry8 |
| Mattress replacement covers |
| Domus Edge11 |
| APEX Domus 4 |
| Serene Edge12 |
| APEX Serene Air 12 APEX Pro-Care Optima 13 |
| APEX Pro-Care Turn |
| THERAPEUTIC SURFACES - Non-Powered (Reactive) 5 Minds Mobility Atlas Hybrid14 |
| 5 Minds Mobility Atlas Auto Comfort15 |
| Therapeutic surfaces replacement covers |
| THERAPEUTIC SURFACES - Modular Non-Powered Adapt 2 Series, including NEW! Paediatric size |
| Acute |
| MATTRESS OVERLAYS ROHO Dry Floatation, SOFFLEX 2 & Prodigy overlays23 |
| Leveling pads |
| ACCESSORIES |
| Bolsters 11-12 Overlay base 24 |
| Gap fillers, including 8"H size for air surfaces24 Overbed tables |
| Halo Wing |
| BEDS |
| Tröst & Tröst Plus beds, packages & accessories |
| Resten Hi-Lo, NEW! Queen size |
| Resten Luxury & Hybrid mattresses |
| NOA Elite Riser 30 NOA Light 31 |
| NOA Bariatric Elite32 |
| NOA Elite accessories 32 Specifications & warranties 33 |

Mattresses & therapeutic surfaces matrix

| Features | , | | | | _ | | | | | | | | | | | |
|-------------------------|---|---|--|----------------------|-------------------------|---------------------------------------|----------------------------|---|-----------------|----------------------|---|---|----------------------------------|-----------------------------|-------------------|--------------|
| | Non- powered pressure redistrib. | Shear Reduction Technology (SRT) | Insert ROHO mattress sections | Gel swirl foam | Firm side borders | Raised bolster custom option | Special foot section | Alternating Pressure Therapy (APT) | Low Air Loss | Pulsation therapy | Continuous Lateral Rotation Technology | Recovery 5 or Recov- ery 5 Ultra cover | Fluid resis- tant cover | Plus/ bariatric sizes | Weight limit (lb) | Warranty |
| 5 Minds N | Nobility | Mattresse | S | | | | | ' | | | <u>'</u> | <u>'</u> | | _ | <u>'</u> | |
| Adapt | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | | | | | √ | | ✓ | 350- 500 | 24 months |
| Allay | ✓ | ✓ | | √ | ✓ | ✓ | | | | | | √ | | ✓ | 350- 500 | 24 months |
| Align & Align Bari. | √ | ✓ | | | √ | ✓ | ✓ | | | | | ✓ | | ✓ | 350- 800 | 24 months |
| Atlas Hybrid | ✓ | | | ✓ | ✓ | | ✓ | with pump | | | | | √ stretch | ✓ | 500- 700 | 36 months |
| Atlas Auto Comfort | √ | | | √ | ✓ | | √ | | | | | | √ stretch | ✓ | 500- 700 | 36 months |
| Permobil | Signatu | re Mattres | ses | | | <u>'</u> | | | | | | | | | <u>'</u> | |
| Sömn & Koppla | ✓ | | | | ✓ | Sömn only | | | | | | | √ Startex | ✓ | 350- 500 | 24 months |
| Optima 2 | √ | | | | | | | | | | | | | √ | 350- 500 | 24 months |
| Symmetry | √ | | | | √ | | ✓ | | | | | | √ stretch | ✓ | 350- 500 | 24 months |
| Resten Luxe & Hybrid | | | | ✓ | | | | | | | | | | ✓ | 350- 600 | 24 months |
| APEX & 5 | Minds | Mobility Po | owered | Surfa | ces | | | | | | | | | _ | | |
| Domus Edge | | | | | √ | | | √ | √ 8 L/min | | | | √ quilted | ✓ | 450 | 36 months |
| Domus 4 | | | | | | | | ✓ | √ 8 L/min | | | | √ quilted | | 450 | 36 months |
| Domus Auto | | | | | | | | ✓ | √ 8 L/min | | | | √ quilted | | 450 | 36 months |
| Serene Edge | | | | | √ | | | √ | √ 120 L/min | ✓ | | | quilted | ✓ | 400- 650 | 36 months |
| Serene Air | | | | | | | | ✓ | 120 L/min | √ | | | √ quilted | √ | 400- 920 | 36 months |
| Pro-Care Optima | | | | | | | √ | ✓ | 13 L/min | | | | √ stretch | | 550 | 36 months |
| Pro-Care Turn | | | | | | | ✓ | √ | 13 L/min | | √ | | ✓ | | 395 | 36 months |
| | Modula | r Systems | | | | | | | 13 [/111111 | | | | stretch | | | |
| Adapt 2 Series | √ | √ | √ | √ | √ | | ✓ | | | | | ✓ | | ✓ | 350 | 24 months |
| Acute | √ | | √ | √ | | | ✓ | | | | | ✓ | | | 350 | 24 months |
| Permobil | Overlay | s | | | | | | | | | | | | | | |
| ROHO Prodigy | ✓ | | | | | | | | | | | | | | 300 | 12 months |
| ROHO SOFFLEX 2 | √ | | | | | | | | | | | | | | 330 | 18 months |
| ROHO Dry Floatation | √ | | | | | | | | | | | | ✓ | ✓ | unlimited | 24 months |
| Leveling pads | ✓ | | ✓ | | | | | | | | | ✓ | | | 250- 1,000 | 24 months |

Key terms



Firm borders & bolsters

Firm perimeter borders on our many mattresses and surfaces ensure a solid transfer edge for improving safety and reducing entrapment risk. Domus Edge (shown above) and Serene Edge air surfaces provide extra support for clients. New add-on bolsters with 6" wide transfer edges

expands the 36" wide Domus 4 or Auto to 48" wide, and the Serene Air up to 60" wide.





High resiliency and **Hybrid foam**

High resiliency foam is a versatile foam used in many 5 Minds Mobility mattresses. It is considered to be a highly tensile surface, and can be cut to any size.

Our new Hybrid foam inserts for the Adapt 2 are even more supportive with 1.5" of high resiliency foam topped with 1.5" of Serene™ foam. Serene is similar to memory foam, but is even better for temperature regulation and breathability.



Shear Reduction Technology (SRT)

5 Minds Mobility's patented SRT is an innovative solution designed to address shearing forces present during client movement in bed, or when raising the head or foot section of the bed.

Low friction material separates two layers of foam - effectively moving the shearing to the core of the mattress and away from the client.



Gel swirl foam

Gel swirl foam provides enhanced comfort while reducing motion transfer across the mattress surface for an undisturbed sleep. The open cell structure promotes airflow allowing for excellent cooling. Gel swirl foam particles also:

- promote a cool, comfortable sleep
- ensure better heat conductivity than traditional visco elastic foam
- are non-toxic and hypo-allergenic



Recovery 5 covers

All modular foam bases and inserts, as well as the 5 Minds Mobility foam mattresses come standard with a Recovery 5 anti-microbial cover which provides incontinent protection and aids in keeping skin cool and dry.

It is removable, machine-washable, fluid-proof and fire-retardant. This material helps wick away moisture and increases air flow and breathability. The vinyl bottom under the mattress includes handles.

UPGRADE to the R5 Ultra! The new Recovery 5 Ultra is an additional zippered enclosure cover for the full mattress surface (and inserts) with 4-way stretch, is machine washable, and provides up to 6X more incontinent protection.

Note: For all modular surfaces, zippered top covers are required and sold separately.

ZIPPERED ENCLOSURE: Adapt 2 mattresses have zippered enclosure cover options in either Recovery 5 or R5 Ultra material. This removable top cover with its 360° zipper must be ordered separately.

For the Acute mattress, you can add a coverlet which has stretchy corners and fits over the entire mattress and inserts (shown above).

Mattresses

Adapt mattress

The Adapt modifiable foam surface with SRT addresses shearing forces present during client movement in bed or when raising the head or foot section of the bed, at the same time ensuring proper pressure redistribution. This mattress comes complete with three removable Gel swirl foam inserts, but allows for the addition of up to two ROHO® air mattress sections (sold separately). It comes with firm borders, a raised sleep surface, and full enclosure cover with a Recovery 5 top. Upgrade the mattress with a raised bolster or R5 Ultra cover (see below). **Note:** Custom sizes available, please call for details.



| Part Number | Description | Retail Price |
|--------------|---|--------------|
| AD3680 | ADAPT w/ SRT, 36" x 80" x 7" w/ R5 TOP COVER; 350 LB WT CAP. | \$1,497.00 |
| AD3684 | ADAPT w/ SRT, 36" x 84" x 7" w/ R5 TOP COVER; 350 LB WT CAP. | \$1,612.00 |
| AD3980 | ADAPT w/ SRT, 39" x 80" x 7" w/ R5 TOP COVER; 350 LB WT CAP. | \$1,612.00 |
| AD4280 | ADAPT w/ SRT, 42" x 80" x 7" w/ R5 TOP COVER; 500 LB WT CAP. | \$1,728.00 |
| AD4284 | ADAPT w/ SRT, 42" x 84" x 7" w/ R5 TOP COVER; 500 LB WT CAP. | \$1,843.00 |
| AD4880 | ADAPT w/ SRT, 48" x 80" x 7" w/ R5 TOP COVER; 500 LB WT CAP. | \$1,959.00 |
| AD5480 | ADAPT w/ SRT, 54" x 80" x 7" w/ R5 TOP COVER; 500 LB WT CAP. | \$2,074.00 |
| AD6080 | ADAPT w/ SRT, 60" x 80" x 7" w/ R5 TOP COVER; 500 LB WT CAP. | \$2,190.00 |
| Other | | |
| MATT1610 | 19" X 30" ROHO NEOPRENE MATTRESS SECTION (FOR 36" WIDE MATTRESSES) | \$1,171.00 |
| MATTSECT3400 | 19" X 34" ROHO NEOPRENE MATTRESS SECTION (FOR ALL OTHER MATTRESSES) | \$1,171.00 |
| RB-UPGRADE | CUSTOM* RAISED BOLSTER UPGRADE FOR ANY MATTRESS SIZE; MINIMUM QUANTITY ORDER OF FOUR OF THE SAME SIZE | \$210.00 |
| R5U-UPGRADE | ENHANCE THE COVER WITH R5 ULTRA FOR MAXIMUM INCONTINENCE PROTECTION | \$105.00 |

^{*}These systems are custom and made to order. They are non-refundable and non-exchangeable.

Allay mattress

The Allay Gel swirl foam replacement surface with SRT is our premium option for the client wanting optimal pressure redistribution and comfort. The cooling benefit of the gel-infused Visco foam top combined with a supportive core make this the best option for users <u>not</u> requiring ROHO air mattress sections. It is reversible, has firm borders, a raised sleep surface, and full enclosure cover with a Recovery 5 top. Upgrade the mattress with a raised bolster or R5 Ultra cover (see below). **Note**: Custom sizes available, please call for details.




| Part Number | Description | Retail Price |
|-------------|---|--------------|
| AL3680 | ALLAY w/ SRT, 36" x 80" x 7" w/ R5 TOP COVER; 350 LB WT CAP. | \$1,431.00 |
| AL3980 | ALLAY w/ SRT, 39" x 80" x 7" w/ R5 TOP COVER; 350 LB WT CAP. | \$1,562.00 |
| AL4280 | ALLAY w/ SRT, 42" x 80" x 7" w/ R5 TOP COVER; 500 LB WT CAP. | \$1,691.00 |
| AL4880 | ALLAY w/ SRT, 48" x 80" x 7" w/ R5 TOP COVER; 500 LB WT CAP. | \$1,952.00 |
| AL5480 | ALLAY w/ SRT, 54" x 80" x 7" w/ R5 TOP COVER; 500 LB WT CAP. | \$2,083.00 |
| AL6080 | ALLAY w/ SRT, 60" x 80" x 7" w/ R5 TOP COVER; 500 LB WT CAP. | \$2,213.00 |
| Other | | |
| RB-UPGRADE | CUSTOM* RAISED BOLSTER UPGRADE FOR ANY MATTRESS SIZE; MINIMUM QUANTITY ORDER OF FOUR OF THE SAME SIZE | \$210.00 |
| R5U-UPGRADE | ENHANCE THE COVER WITH R5 ULTRA FOR MAXIMUM INCONTINENCE PROTECTION | \$105.00 |
| NOU-UFGRADE | ENDANCE THE COVER WITH ROULTRA FOR IMAXIMOM INCOMPINENCE PROTECTION | \$100 |

^{*}These systems are custom and made to order. They are non-refundable and non-exchangeable.

Align mattress & Align Bariatric

The Align surface with SRT comes with firm borders, a raised sleep surface, and full enclosure cover with a Recovery 5 top. It is a unique multi-zone surface that provides differing immersion and envelopment levels that match the contours of the body. There is one piece of HR foam at the torso, with 4" square cuts in the foam under the sacrum, and 2" square cuts under the heels to release the surface tension (see pics below). The foam density and weight capacity are the differences for our Align & Align Bariatric options. Upgrade the mattress with a raised bolster or R5 Ultra cover (see below). **Note:** Custom sizes available, please call for details.



| Part Number | Description | Retail Price |
|-------------|---|--------------|
| A3680 | ALIGN w/ SRT - HR FOAM, 36" x 80" x 7" w/ R5 TOP COVER; 350 LB WT CAP. | \$1,190.00 |
| A3684 | ALIGN w/ SRT - HR FOAM, 36" x 84" x 7" w/ R5 TOP COVER; 350 LB WT CAP. | \$1,381.00 |
| A3980 | ALIGN w/ SRT - HR FOAM, 39" x 80" x 7" w/ R5 TOP COVER; 350 LB WT CAP. | \$1,323.00 |
| A4280 | ALIGN w/ SRT - HR FOAM, 42" x 80" x 7" w/ R5 TOP COVER; 500 LB WT CAP. | \$1,456.00 |
| A4284 | ALIGN w/ SRT - HR FOAM, 42" x 84" x 7" w/ R5 TOP COVER; 500 LB WT CAP. | \$1,612.00 |
| A4880 | ALIGN w/ SRT - HR FOAM, 48" x 80" x 7" w/ R5 TOP COVER; 500 LB WT CAP. | \$1,721.00 |
| A5480 | ALIGN w/ SRT - HR FOAM, 54" x 80" x 7" w/ R5 TOP COVER; 500 LB WT CAP. | \$1,855.00 |
| A6080 | ALIGN w/ SRT - HR FOAM, 60" x 80" x 7" w/ R5 TOP COVER; 500 LB WT CAP. | \$1,987.00 |
| Other | | |
| RB-UPGRADE | CUSTOM* RAISED BOLSTER UPGRADE FOR ANY MATTRESS SIZE; MINIMUM QUANTITY ORDER OF FOUR OF THE SAME SIZE | \$210.00 |
| R5U-UPGRADE | ENHANCE THE COVER WITH R5 ULTRA FOR MAXIMUM INCONTINENCE PROTECTION | \$105.00 |
| BARIATRIC | | |
| A3680B | BARIATRIC HIGH DENSITY FOAM w/ SRT, 36" x 80" x 7" w/ R5 TOP COVER; 600 LB WT CAP. | \$1,266.00 |
| A3980B | BARIATRIC HIGH DENSITY FOAM w/ SRT, 39" x 80" x 7" w/ R5 TOP COVER; 600 LB WT CAP. | \$1,381.00 |
| A4280B | BARIATRIC HIGH DENSITY FOAM w/ SRT, 42" x 80" x 7" w/ R5 TOP COVER; 600 LB WT CAP. | \$1,381.00 |
| A4880B | BARIATRIC HIGH DENSITY FOAM W/ SRT, 42 X 80" X 7" W/ R5 TOP COVER, 800 LB WT CAP. | |
| | | \$1,612.00 |
| A5480B | BARIATRIC HIGH DENSITY FOAM w/ SRT, 54" x 80" x 7" w/ R5 TOP COVER; 800 LB WT CAP. | \$1,728.00 |
| A6080B | BARIATRIC HIGH DENSITY FOAM w/ SRT, 60" x 80" x 7" w/ R5 TOP COVER; 800 LB WT CAP. | \$1,843.00 |
| Other | | |
| RB-UPGRADE | CUSTOM* RAISED BOLSTER UPGRADE FOR ANY MATTRESS SIZE; MINIMUM QUANTITY ORDER OF FOUR OF THE SAME SIZE | \$210.00 |
| R5U-UPGRADE | ENHANCE THE COVER WITH R5 ULTRA FOR MAXIMUM INCONTINENCE PROTECTION | \$105.00 |

^{*}These systems are custom and made to order. They are non-refundable and non-exchangeable.



One piece of HR foam at the torso



4" square cuts under the sacrum



2" square cuts under the heels

Permobil Signature mattresses

Permobil's newest Signature mattresses - the Sömn and Koppla - are customizable foam mattresses are specifically engineered to provide improved performance using pressure-reducing materials. We also offer the updated Optima 2 and Symmetry.

Sömn and Koppla premium mattresses are ideal for long-term care with their 3" wide transfer edges on both sides that help reduce the risk of falls. Each is comprised of a layered foam surface: **Sömn** is firmer with a 4" base layer and 2" top layer with Serene™ foam at the heel, while **Koppla** has a 2" thick foam base, plus a 2" middle and a 2" top layer of foam with decreasing density. Sömn and Koppla both come with fluid-resistant Startex full enclosure covers.



| Part Number | Description | Retail Price |
|-------------|---|--------------|
| SOMN3680 | SÖMN, 36" x 80" x 6" w/ STARTEX COVER; 350 LB WT CAP. | \$923.00 |
| SOMN3980 | SÖMN, 39" x 80" x 6" w/ STARTEX COVER; 350 LB WT CAP. | \$1,097.00 |
| SOMN4280 | SÖMN, 42" x 80" x 6" w/ STARTEX COVER; 350 LB WT CAP. | \$1,154.00 |
| SOMN4880 | SÖMN, 48" x 80" x 6" w/ STARTEX COVER; 500 LB WT CAP. | \$1,212.00 |
| KOPPLA3680 | KOPPLA, 36" x 80" x 6" w/ STARTEX COVER; 350 LB WT CAP. | \$1,097.00 |
| KOPPLA3980 | KOPPLA, 39" x 80" x 6" w/ STARTEX COVER; 350 LB WT CAP. | \$1,270.00 |
| KOPPLA4280 | KOPPLA, 42" x 80" x 6" w/ STARTEX COVER; 350 LB WT CAP. | \$1,328.00 |
| KOPPLA4880 | KOPPLA, 48" x 80" x 6" w/ STARTEX COVER; 500 LB WT CAP. | \$1,385.00 |

Optima 2 mattresses focus on comfort, safety, durability, and budget. They come with a waterproof nylon cover and a foam core constructed of 100% medical grade environmentally friendly non-allergenic foam. Mattress is reversible for longer product life.

| Part Number | Description | Retail Price |
|-------------|--|--------------|
| OPT2-3680 | OPTIMA 2, 36" x 80" x 6" w/ NYLON TOP COVER; 350 LB WT CAP. | \$472.00 |
| OPT2-3980 | OPTIMA 2, 39" x 80" x 6" w/ NYLON TOP COVER; 350 LB WT CAP. | \$629.00 |
| OPT2-4280 | OPTIMA 2 XT, 42" x 80" x 6" w/ NYLON TOP COVER; 350 LB WT CAP. | \$703.00 |
| OPT2-4880 | OPTIMA 2 XT, 48" x 80" x 6" w/ NYLON TOP COVER; 500 LB WT CAP. | \$818.00 |

Symmetry has three layers of foam with a convoluted high resiliency foam on top to improve client mobility, including a softer foam at the heel. It also includes firm side borders for secure transfers in/out of bed, and for keeping the client in a central position. Symmetry has a multiway stretch, fluid-resistant cover which provides incontinent protection.



| Part Number | Description | Retail Price |
|-------------|---|--------------|
| SYM-3680 | SYMMETRY, 36" x 80" x 6" w/ FLUID-RESISTANT TOP COVER & FIRM BORDERS; 350 LB WT CAP. | \$671.00 |
| SYM-3980 | SYMMETRY, 39" x 80" x 6" w/ FLUID-RESISTANT TOP COVER & FIRM BORDERS; 350 LB WT CAP. | \$808.00 |
| SYM-4280 | SYMMETRY XT, 42" x 80" x 6" w/ FLUID-RESISTANT TOP COVER & FIRM BORDERS; 350 LB WT CAP. | \$881.00 |
| SYM-4880 | SYMMETRY XT, 48" x 80" x 6" w/ FLUID-RESISTANT TOP COVER & FIRM BORDERS; 500 LB WT CAP. | \$986.00 |

Permobil Signature key terms

Foam borders (on the Symmetry, Sömn, and Koppla) - Constructed of 2.5 lb density and IFD 65. Provides firm edge support and helps to keep the patient in a central position and from falling off the mattress.

5 Minds Mobility replacement covers

Replacement enclosure covers with handles and vinyl bottom for 7" high: **Adapt, Allay, and Align** mattresses, including the new larger 84" length cover sizes. Custom sizes are also available. If ordering a custom raised bolster mattress, specialized raised bolster covers are available.



Raised bolster Recovery 5 cover shown above.

| Part Number | Description | Retail Price | |
|------------------|---|--------------|--|
| COV-REPL3680 | FULL ENCLOSURE COVER, ADAPT/ALLAY/ALIGN 36x80 (RECOVERY 5 FLUID-RESISTANCE, MULTI-WAY STRETCH) | \$508.00 | |
| COV-REPL3684 | FULL ENCLOSURE COVER, ADAPT/ALLAY/ALIGN 36x84 (RECOVERY 5 FLUID-RESISTANCE, MULTI-WAY STRETCH) | | |
| COV-REPL3980 | FULL ENCLOSURE COVER, ADAPT/ALLAY/ALIGN 39x80 (RECOVERY 5 FLUID-RESISTANCE, MULTI-WAY STRETCH) | \$573.00 | |
| COV-REPL4280 | FULL ENCLOSURE COVER, ADAPT/ALLAY/ALIGN 42x80 (RECOVERY 5 FLUID-RESISTANCE, MULTI-WAY STRETCH) | \$573.00 | |
| COV-REPL4284 | FULL ENCLOSURE COVER, ADAPT/ALLAY/ALIGN 42x84 (RECOVERY 5 FLUID-RESISTANCE, MULTI-WAY STRETCH) | \$630.00 | |
| COV-REPL4880 | FULL ENCLOSURE COVER, ADAPT/ALLAY/ALIGN 48x80 (RECOVERY 5 FLUID-RESISTANCE, MULTI-WAY STRETCH) | \$573.00 | |
| COV-REPL5480 | FULL ENCLOSURE COVER, ADAPT/ALLAY/ALIGN 54x80 (RECOVERY 5 FLUID-RESISTANCE, MULTI-WAY STRETCH) | \$625.00 | |
| COV-REPL6080 | FULL ENCLOSURE COVER, ADAPT/ALLAY/ALIGN 60x80 (RECOVERY 5 FLUID-RESISTANCE, MULTI-WAY STRETCH) | \$625.00 | |
| Recovery 5 Ultra | | | |
| COV-REPL3680-R5U | FULL ENCLOSURE COVER, ADAPT/ALLAY/ALIGN 36x80 (RECOVERY 5 ULTRA SUPERIOR FLUID-RESISTANCE, MULTI-WAY STRETCH) | \$573.00 | |
| COV-REPL3684-R5U | FULL ENCLOSURE COVER, ADAPT/ALLAY/ALIGN 36x84 (RECOVERY 5 ULTRA SUPERIOR FLUID-RESISTANCE, MULTI-WAY STRETCH) | \$630.00 | |
| COV-REPL3980-R5U | FULL ENCLOSURE COVER, ADAPT/ALLAY/ALIGN 39x80 (RECOVERY 5 ULTRA SUPERIOR FLUID-RESISTANCE, MULTI-WAY STRETCH) | \$688.00 | |
| COV-REPL4280-R5U | FULL ENCLOSURE COVER, ADAPT/ALLAY/ALIGN 42x80 (RECOVERY 5 ULTRA SUPERIOR FLUID-RESISTANCE, MULTI-WAY STRETCH) | \$688.00 | |
| COV-REPL4284-R5U | FULL ENCLOSURE COVER, ADAPT/ALLAY/ALIGN 42x84 (RECOVERY 5 ULTRA SUPERIOR FLUID-RESISTANCE, MULTI-WAY STRETCH) | \$688.00 | |
| COV-REPL4880-R5U | FULL ENCLOSURE COVER, ADAPT/ALLAY/ALIGN 48x80 (RECOVERY 5 ULTRA SUPERIOR FLUID-RESISTANCE, MULTI-WAY STRETCH) | \$688.00 | |
| COV-REPL5480-R5U | FULL ENCLOSURE COVER, ADAPT/ALLAY/ALIGN 54x80 (RECOVERY 5 ULTRA SUPERIOR FLUID-RESISTANCE, MULTI-WAY STRETCH) | \$746.00 | |
| COV-REPL6080-R5U | FULL ENCLOSURE COVER, ADAPT/ALLAY/ALIGN 60x80 (RECOVERY 5 ULTRA SUPERIOR FLUID-RESISTANCE, MULTI-WAY STRETCH) | \$746.00 | |

Permobil Signature replacement covers

Replacement enclosure covers with handles and vinyl bottom for 6" high: Sömn, Koppla, Optima 2, and Symmetry mattresses.

| Part Number | Description | Retail Price |
|-------------------|---|--------------|
| COV-REPL3680-STX | FULL ENCLOSURE COVER, SOMN/KOPPLA 36x80 (STARTEX, FLUID-RESISTANT, MULTI-WAY STRETCH) | \$524.00 |
| COV-REPL3980-STX | FULL ENCLOSURE COVER, SOMN/KOPPLA 39x80 (STARTEX, FLUID-RESISTANT, MULTI-WAY STRETCH) | \$556.00 |
| COV-REPL4280-STX | FULL ENCLOSURE COVER, SOMN/KOPPLA 42x80 (STARTEX, FLUID-RESISTANT, MULTI-WAY STRETCH) | \$577.00 |
| COV-REPL4880-STX | FULL ENCLOSURE COVER, SOMN/KOPPLA 48x80 (STARTEX, FLUID-RESISTANT, MULTI-WAY STRETCH) | \$598.00 |
| COV-REPL3680-OPT2 | FULL ENCLOSURE COVER, OPTIMA 2 36x80 (NYLON, WATERPROOF) | \$188.00 |
| COV-REPL3980-OPT2 | FULL ENCLOSURE COVER, OPTIMA 2 39x80 (NYLON, WATERPROOF) | \$220.00 |
| COV-REPL4280-OPT2 | FULL ENCLOSURE COVER, OPTIMA 2 42x80 (NYLON, WATERPROOF) | \$251.00 |
| COV-REPL4880-OPT2 | FULL ENCLOSURE COVER, OPTIMA 2 48x80 (NYLON, WATERPROOF) | \$283.00 |
| COV-REPL3680-SYM | FULL ENCLOSURE COVER, SYMMETRY 36x80 (FLUID-RESISTANT, MULTI-WAY STRETCH) | \$262.00 |
| COV-REPL3980-SYM | FULL ENCLOSURE COVER, SYMMETRY 39x80 (FLUID-RESISTANT, MULTI-WAY STRETCH) | \$293.00 |
| COV-REPL4280-SYM | FULL ENCLOSURE COVER, SYMMETRY 42x80 (FLUID-RESISTANT, MULTI-WAY STRETCH) | \$325.00 |
| COV-REPL4880-SYM | FULL ENCLOSURE COVER, SYMMETRY 48x80 (FLUID-RESISTANT, MULTI-WAY STRETCH) | \$356.00 |

Therapeutic Surfaces

- Powered (Active)
- Non-Powered (Reactive)
- Modular Non-Powered
- Mattress Overlays

Domus Edge

Low Air Loss Alternating Pressure Therapy system with transfer edge





Domus pump: Client weight references are incorporated into the control panel so caregivers can easily find the correct level of comfort and firmness needed for optimum comfort and healing.

New bolsters transform Domus 4 or Auto into a 48" wide edge surface!



At a glance

- Suitable for clients with medium to high risk for pressure injuries
- · Low Air Loss for microclimate control
- 3" wide transfer edges
- Lower 3" of air cells for the whole surface remain static to prevent bottoming out
- Fluid-resistant, four-way stretch, quilted cover
- Zippered cover allows for one-step removal
- . Maximum client weight of up to 450 lb
- 2 widths: 35.5" / 42" x 79"L x 8"H
- Gap filler extends 36" model to 84"L, see pg. 24

Domus Edge Alternating Pressure Therapy system takes powered therapeutic surfaces to the next level with the added safety feature of a 3" wide transfer edge. This helps enable safe and easy transfers on and off the surface, and also reduces falls and entrapment. Domus Edge has 10-minute operating cycle that allows the bladders to inflate/deflate to offload and transfer body pressure. Functions include Alternating Pressure Therapy and Constant Low Pressure (static), and the system also features max inflate and seat inflation for the Fowler position. The lower 3" of air cells remain static, even during power outages, to prevent clients from bottoming out. **New bolsters** with full enclosure R5 cover and 6" wide firm transfer edges on each side combine with Domus 4 or Auto to make a Domus Edge surface 48" wide*.

| Part Number | Description | Retail Price |
|-------------|--|--------------|
| DOM4EDG36 | DOMUS EDGE 36" SYSTEM: AIR SURFACE, PUMP, TOP & BOTTOM COVERS; 450 LB WT CAP. | \$3,495.00 |
| DOM4EDG42 | DOMUS EDGE 42" SYSTEM: AIR SURFACE, PUMP, TOP & BOTTOM COVERS; 450 LB WT CAP. | \$4,395.00 |
| DOMUS-PUMP | REPLACEMENT PUMP; USED ON ALL DOMUS SURFACES | \$1,295.00 |
| BOLSTER48 | FOAM BOLSTER WITH 6" FIRM TRANSFER EDGES ON EACH SIDE; FITS 36" DOMUS 4 or AUTO* | \$849.00 |

^{*}Domus system sold separately

Domus 4

Low Air Loss Alternating Pressure Therapy system



At a glance

- Suitable for clients with medium to high risk for pressure injuries
- · Low Air Loss for microclimate control
- Cell-on-cell bladder to avoid bottoming out
- Silver* nylon, fluid-resistant, quilted cover
- · Zippered cover allows for one-step removal
- Visual and audible alarm enhances clients' safety & security
- Maximum client weight of up to 450 lb
- 35.5"W x 79"L x 8"H; gap filler extends system to 84"L, see pg. 24

Domus 4 Alternating Pressure Therapy system has a 10-minute operating cycle that allows the bladders to inflate/deflate to offload and transfer body pressure. Functions include Alternating Pressure Therapy and Constant Low Pressure (static), and the system also features max inflate and seat inflation for the Fowler position. To protect your clients, the Domus 4 has a number of added safety features like a cell-on-cell bladder to avoid bottoming out, a highly visible CPR release valve, and an auto-lock on the settings after five (5) minutes.

| Part Number | Description | Retail Price |
|-------------|--|--------------|
| DOM436 | DOMUS 4 SYSTEM: AIR SURFACE, PUMP, TOP & BOTTOM COVERS; 450 LB WT CAP. | \$1,995.00 |

Domus Auto

Alternating Pressure Redistribution System

Domus Auto Alternating Pressure Therapy system has all the same features of the Domus 4 system (see above) with the added convenience of an automatic adjustment system that allows for fast and simple setup, along with finetuning to optimize client comfort.

| Part Number | Description | Retail Price |
|-------------|---|--------------|
| DOMAUTO36 | DOMUS AUTO SYSTEM: AIR SURFACE, PUMP, TOP & BOTTOM COVERS; 450 LB WT CAP. | \$2,495.00 |


Serene Edge

True Low Air Loss APT & Pulsation system with transfer edge





incorporated into the control panel so caregivers can easily find the correct level of comfort and



At a glance

- Suitable for vulnerable clients with high to very high risk for pressure injuries
- True Low Air Loss for microclimate control
- High performance blower offers >120 L/min of air loss for effective microclimate control
- 3" wide transfer edges plus 2" foam layer built into the surface base
- Fluid-resistant, four-way stretch, quilted cover
- · Zippered cover allows for one-step removal
- · Maximum client weight of up to 650 lb
- 2 widths: 36" / 42" x 79"L x 8"H
- Gap filler extends 36" model to 84"L, see pg. 24



Serene Edge is a dual mode system with Alternating Pressure Therapy (APT) or Pulsation options with the added safety feature of a 3" wide transfer edge. This helps enable safe and easy transfers on and off the surface, and also reduces falls and entrapment. APT modes assist in wound healing and prevention and there are four operating cycles (5/10/15/20 min.) based on therapy requirements. The Pulsation mode increases and decreases the surface's air volume every 15 seconds. Serene Edge also has a Constant Low Pressure (static) option, as well as a max inflate



| Part Number | Description | Retail Price |
|-------------|--|--------------|
| SERAEDG36 | SERENE EDGE 36" SYSTEM: AIR SURFACE, PUMP, TOP & BOTTOM COVERS; 400 LB WT CAP. | \$6,995.00 |
| SERAEDG42 | SERENE EDGE 42" SYSTEM: AIR SURFACE, PUMP, TOP & BOTTOM COVERS; 650 LB WT CAP. | \$7,995.00 |
| SERENE-PUMP | REPLACEMENT PUMP; USED ON ALL SERENE SURFACES | \$3,295.00 |
| BOLSTER48 | FOAM BOLSTER WITH 6" FIRM TRANSFER EDGES ON EACH SIDE; FITS 36" SERENE AIR* | \$849.00 |
| BOLSTER54 | FOAM BOLSTER WITH 6" FIRM TRANSFER EDGES ON EACH SIDE; FITS 42" SERENE AIR* | \$949.00 |
| BOLSTER60 | FOAM BOLSTER WITH 6" FIRM TRANSFER EDGES ON EACH SIDE; FITS 48" SERENE AIR* | \$1,049.00 |

and seat inflation for the Fowler position. A 2" foam layer is also built into the base of the surface to prevent clients from bottoming out. New bolsters (shown at right) with full enclosure R5 cover and 6" wide firm transfer edges on

each side combine with Serene Air systems to make Serene Edge surfaces in 48", 54", and 60" widths*.

Serene Air

True Low Air Loss Alternating Pressure Therapy & Pulsation system





At a glance

- Suitable for vulnerable clients with high to very high risk for pressure injuries
- True Low Air Loss for microclimate control
- High performance blower offers >120 L/min of air loss for effective microclimate control
- Silver* nylon, fluid-resistant, quilted cover
- · Zippered cover allows for one-step removal
- Maximum client weight of up to 920 lb (on 48" mattress)
- 3 widths: 36" / 42" / 48" x 79"L x 8"H
- Gap filler extends 36" model to 84"L, see pg. 24

Serene Air is a dual mode system with Alternating Pressure Therapy (APT) or Pulsation options. Serene Air also has a Constant Low Pressure (static) option, as well as a max inflate and seat inflation for the Fowler position. It includes a zippered full-length pocket with a 1" high foam base. Initial inflation takes less than three (3) minutes, and the surface has a quilted top cover with high moisture vapor transmission rate (MVTR). Serene Air has a number of added safety and convenience features: auto-lock on the settings after five (5) minutes, accessible and highly visible CPR release valve at the control pump, plus a handy carrying bag for compact storage and transport.

| Part Number | Description | Retail Price |
|-------------|--|--------------|
| SERA36 | SERENE AIR 36" SYSTEM: AIR SURFACE, PUMP, T&B COVERS, CARRYING BAG; 400 LB WT CAP. | \$5,995.00 |
| SERA42 | SERENE AIR 42" SYSTEM: AIR SURFACE, PUMP, T&B COVERS, CARRYING BAG; 650 LB WT CAP. | \$6,995.00 |
| SERA48 | SERENE AIR 48" SYSTEM: AIR SURFACE, PUMP, T&B COVERS, CARRYING BAG; 920 LB WT CAP. | \$7,995.00 |

^{*}Serene Air system sold separately

Pro-Care Optima

Alternating Pressure Therapy





Pro-Care Optima pump: One-touch set up, auto-sensing technology, and an intuitive graphic display. System has four operating cycles (10/15/20/25 min.) to meet clients' needs.

At a glance

- Suitable for clients with high to very high risk for pressure injuries
- Auto-sensing & auto-inflation technology, with intuitive one-touch set up
- · Comes standard with built-in battery backup
- Multi-zone air pressure for optimal comfort, including zero-pressure heel protection zone
- · Fluid-resistant and four-way stretch cover
- Maximum client weight of up to 550 lb
- 36"W x 79"L x 8"H

Pro-Care Optima Alternating Pressure Therapy system provides stability and comfort with three distinct zones. It fits your clients' contours offering offloading/redirecting pressure with distinct zones: torso, lower legs, and a static head pillow. The surface features a six cell-in-cell bladder for the sacrum area to provide excellent support when clients are sitting up. This prevents bottoming out, and the rectangular cells offer comfortable pressure relief. With a simple control knob located on the side of the five cylinders at the foot section, you can activate the zero pressure heel protection that allows the heels to be offloaded by deflating up to two bladders in the lower leg zone of the surface.

The auto-adjust feature provides optimal pressure based on the client's changing weight distribution. It also auto-inflates and adjusts pressure in the Fowler position when the head section is raised over 30 degrees to prevent bottoming out. Pro-Care Optima has the battery backup built-in which gives peace of mind and added protection. It can provide non-stop therapy for six hours during transport or in case of power failure.

| Part Number | Description | Retail Price |
|-------------|---|--------------|
| PCOPT | PRO-CARE OPTIMA 36": AIR SURFACE, PUMP, TOP & BOTTOM COVERS; 550 LB WT CAP. | \$9,995.00 |

Pro-Care Turn

Continuous Lateral Rotation Technology (CLRT)



At a glance

- Suitable for clients with high to very high risk for pressure injuries
- 30° of bilateral rotation therapy to enhance pressure relief
- 3-in-1 therapy available with effortless client repositioning
- Fluid-resistant and four-way stretch cover
- Zippered cover allows for one-step removal
- Visual and audible alarm enhances clients' safety & security
- Maximum client weight of up to 395 lb
- 35.5"W x 79"L x 5"H, with 10" H side bolsters

Bilateral turning can significantly reduce pressure on a client's body while increasing capillary circulation. For immobile or unconscious clients, caregivers must turn them up to 10x/day to achieve this. With the Pro-Care Turn surface, it saves caregivers' time, resources, and health while maintaining quality care for their clients. The 30° degrees of turning not only complies with clinical guidelines, but also provides optimal comfort.

There are three therapies available: CLRT, Alternating Pressure (APT), and Continuous Low Pressure modes, and they can be tailored to individual pressure care plans. The Turning dwell time can be selected between 10, 20, 30, and 60 minutes. APT mode that assists in wound healing and prevention is available in four operating cycles (10/15/20/25 min.) based on therapy requirements.

The high side bolsters reduce the risk of bed entrapment and accidental falling out of bed, and they support the client's arms while in turning therapy. The system also features seat inflation for the Fowler position, and single air-cell deflation for the vulnerable heel section allowing them to float with zero pressure. To protect your clients, the Pro-Care Turn has a number of added safety features like a rotation pause function to allow caregivers to stop a turn right away, transport mode that maintains static pressure relief for six hours during transport or a power failure, plus an accessible and highly visible CPR release valve.

| Part Number | Description | Retail Price |
|-------------|---|--------------|
| PCTURN | PRO-CARE TURN: AIR SURFACE, PUMP, TOP & BOTTOM COVERS; 395 LB WT CAP. | \$9,995.00 |

Atlas Hybrid



At a glance

- Suitable for vulnerable clients at risk for pressure injuries
- Interconnected chambers and valves for pressure redistribution
- Converts to a powered air surface with APT
- Sloped foot section helps reduce pressure from heels
- 4" wide firm borders for safe transfers
- · Fluid-resistant, four-way stretch cover
- Maximum client weight of up to 500-700 lb depending on size of surface
- Atlas Hybrid is available in three widths: 35" / 39" / 42" and 80"L x 7"H; custom sizes are available

The Atlas Hybrid is a non-powered, self-adjusting mattress that utilizes equalization technology to provide maximum pressure redistribution for the prevention and treatment of pressure injuries. Starting at the head, the convoluted high resiliency foam improves client upper body mobility. Next is a cool climate visco gel swirl foam top layer at the torso and middle section for immersion and heat dissipation, followed by a 5-degree sloped foot section made with Serene™ foam to help transfer excess pressure from the heels. Beneath the top layer, interconnected chambers provide pressure redistribution, and there are 4" wide firm borders to ensure a solid transfer edge.

The Atlas Hybrid is also convertible to a powered air surface that uses Alternating Pressure Therapy (APT) with the addition of the digital pump that is sold separately. When the pump is not being used, a valve closure kit is recommended and can be used to adjust the softness/firmness. A manual air pump can also be used to set the firmness of the chambers. Both the closure kit and manual pump come included with the Atlas Hybrid.

| Part Number | Description | Retail Price |
|-------------|--|--------------|
| ATL3580 | ATLAS 35" x 80" HYBRID SYSTEM: FOAM SURFACE, TOP & BOTTOM COVERS; 500 LB WT CAP. | \$1,945.00 |
| ATL3980 | ATLAS 39" x 80" HYBRID SYSTEM: FOAM SURFACE, TOP & BOTTOM COVERS; 600 LB WT CAP. | \$2,095.00 |
| ATL4280 | ATLAS 42" x 80" HYBRID SYSTEM: FOAM SURFACE, TOP & BOTTOM COVERS; 700 LB WT CAP. | \$2,295.00 |
| ATLPUMP | ATLAS DIGITAL PUMP FOR SIZES 35" - 42" | \$379.00 |

What's going on underneath?

- The 4" high interconnected chambers offer great pressure redistribution under the 2" high visco gel swirl foam and 1" standard foam top layers
- These chambers measure 30" wide by 52" long and fit within the cutout starting under the head area, continuing under the full torso, and stopping at the foot section
- When attached to the pump, you can see the alternating therapy provided by the foam/air chambers
- A valve closure kit can be used to keep the system closed when not attached to the digital pump
- The entire mattress is encased in a fire barrier sock fabric to meet fire retardant standards, and then a Startex cover to protect the foam from fluids



Atlas Auto Comfort



At a glance

- Suitable for vulnerable clients at risk for pressure injuries
- 2" high Gel swirl overlay for client comfort
- Interconnected chambers in a closed circuit for self-adjusting pressure redistribution
- Sloped foot section helps reduce pressure from heels
- 4" wide firm borders for safe transfers
- Fluid-resistant, four-way stretch cover
- Maximum client weight of up to 500-700 lb depending on size of surface
- Atlas Auto Comfort is available in three widths: 35" / 39" / 42" and 80"L x 7"H; custom sizes are available

New to Permobil, from the reactive surface family, Atlas Auto Comfort uses equalization technology to provide maximum pressure redistribution to prevent pressure injuries or to support the healing process. Beneath the 2" top layer of gel swirl foam there are interconnected chambers in a closed circuit provide self-adjusting pressure redistribution. Clients also benefit from 4" wide firm borders to ensure a solid transfer edge. Atlas Auto Comfort includes a 5-degree sloped foot section to help transfer excess pressure from the heels.

| Part Number | Description | Retail Price |
|-------------|--|--------------|
| ATLAC3580 | ATLAS AC 35" x 80" SYSTEM: FOAM SURFACE, TOP & BOTTOM COVERS; 500 LB WT CAP. | \$1,895.00 |
| ATLAC3584 | ATLAS AC 35" x 84" SYSTEM: FOAM SURFACE, TOP & BOTTOM COVERS; 500 LB WT CAP. | \$1,995.00 |
| ATLAC3980 | ATLAS AC 39" x 80" SYSTEM: FOAM SURFACE, TOP & BOTTOM COVERS; 600 LB WT CAP. | \$2,045.00 |
| ATLAC3984 | ATLAS AC 39" x 84" SYSTEM: FOAM SURFACE, TOP & BOTTOM COVERS; 600 LB WT CAP. | \$2,145.00 |
| ATLAC4280 | ATLAS AC 42" x 80" SYSTEM: FOAM SURFACE, TOP & BOTTOM COVERS; 700 LB WT CAP. | \$2,245.00 |
| ATLAC4284 | ATLAS AC 42" x 84" SYSTEM: FOAM SURFACE, TOP & BOTTOM COVERS; 700 LB WT CAP. | \$2,345.00 |

Therapeutic surfaces replacement covers

| Part Number | Description | Retail Price |
|---------------|--|--------------|
| Part Number | Description | Retail Price |
| COV-DOM4EDG36 | TOP COVER, DOMUS EDGE 36x80 (FLUID-RESISTANT, MULTI-WAY STRETCH, QUILTED, ZIPPERED) | \$479.00 |
| COV-DOM4EDG42 | TOP COVER, DOMUS EDGE 42x80 (FLUID-RESISTANT, MULTI-WAY STRETCH, QUILTED, ZIPPERED) | \$529.00 |
| COV-DOM36 | TOP COVER, DOMUS 4/AUTO 36x80 (SILVER* NYLON, FLUID-RESISTANT, QUILTED, ZIPPERED) | \$379.00 |
| COV-SERA3680 | TOP COVER, SERENE AIR 36x80 (SILVER* NYLON, FLUID-RESISTANT, QUILTED, ZIPPERED) | \$375.00 |
| COV-SERA4280 | TOP COVER, SERENE AIR 42x80 (SILVER* NYLON, FLUID-RESISTANT, QUILTED, ZIPPERED) | \$425.00 |
| COV-SERA4880 | TOP COVER, SERENE AIR 48x80 (SILVER* NYLON, FLUID-RESISTANT, QUILTED, ZIPPERED) | \$475.00 |
| COV-SERAEDG36 | TOP COVER, SERENE EDGE 36x80 (FLUID-RESISTANT, MULTI-WAY STRETCH, QUILTED, ZIPPERED) | \$479.00 |
| COV-SERAEDG42 | TOP COVER, SERENE EDGE 42x80 (FLUID-RESISTANT, MULTI-WAY STRETCH, QUILTED, ZIPPERED) | \$529.00 |
| COV-PCOPT | TOP COVER, PRO-CARE OPTIMA (FLUID-RESISTANT, MULTI-WAY STRETCH, ZIPPERED) | \$425.00 |
| COV-PCTURN | TOP COVER, PRO-CARE TURN (FLUID-RESISTANT, MULTI-WAY STRETCH, ZIPPERED) | \$725.00 |
| COV-ATL3580 | FULL ENCLOSURE COVER, ATLAS HYBRID 35x80 (FLUID-RESISTANT, MULTI-WAY STRETCH) | \$375.00 |
| COV-ATL3980 | FULL ENCLOSURE COVER, ATLAS HYBRID 39x80 (FLUID-RESISTANT, MULTI-WAY STRETCH) | \$395.00 |
| COV-ATL4280 | FULL ENCLOSURE COVER, ATLAS HYBRID 42x80 (FLUID-RESISTANT, MULTI-WAY STRETCH) | \$475.00 |
| COV-ATLAC3580 | FULL ENCLOSURE COVER, ATLAS AUTO COMFORT 35x80 (FLUID-RESISTANT, MULTI-WAY STRETCH) | \$375.00 |
| COV-ATLAC3584 | FULL ENCLOSURE COVER, ATLAS AUTO COMFORT 35x84 (FLUID-RESISTANT, MULTI-WAY STRETCH) | \$375.00 |
| COV-ATLAC3980 | FULL ENCLOSURE COVER, ATLAS AUTO COMFORT 39x80 (FLUID-RESISTANT, MULTI-WAY STRETCH) | \$395.00 |
| COV-ATLAC3984 | FULL ENCLOSURE COVER, ATLAS AUTO COMFORT 39x84 (FLUID-RESISTANT, MULTI-WAY STRETCH) | \$395.00 |
| COV-ATLAC4280 | FULL ENCLOSURE COVER, ATLAS AUTO COMFORT 42x80 (FLUID-RESISTANT, MULTI-WAY STRETCH) | \$475.00 |
| COV-ATLAC4284 | FULL ENCLOSURE COVER, ATLAS AUTO COMFORT 42x84 (FLUID-RESISTANT, MULTI-WAY STRETCH) | \$475.00 |
| | | |


Note: Adapt 2 and Acute modular systems have their covers listed on their product pages.

Adapt 2 Series modular system

36" x 80"


The Adapt 2 mattress with SRT is an innovative way to reduce the effects of pressure, friction and shear within the body of the mattress. The 3" wide firm perimeter borders and 4" permanent gel foam foot border hold the inserts in place without the need for snaps. The surface inserts can be configured based on the individual client's needs. The removable top cover with its 360° zipper must be ordered separately. Weight capacity: 350 lb

Note: Adapt 2 base and inserts include a fire barrier sock fabric to meet full-scale fire retardant standards (not pictured).



| Part Number | Description | Retail Price |
|--------------|--|--------------|
| AD2-SYS1-36 | ADAPT 2 SERIES SYSTEM 1 PKG. INCLUDES 1 ROHO MATTRESS SECTION (19" x 30") | \$3,292.00 |
| AD2-SYS2-36 | ADAPT 2 SERIES SYSTEM 2 PKG. INCLUDES 2 ROHO MATTRESS SECTIONS (19"x 30") | \$4,260.00 |
| AD2-SYS3-36 | ADAPT 2 SERIES SYSTEM 3 PKG. INCLUDES 3 ROHO MATTRESS SECTIONS (19" x 30") | \$5,278.00 |
| AD2-SYS4-36 | ADAPT 2 SERIES SYSTEM 4 PKG. INCLUDES 4 ROHO MATTRESS SECTIONS (19" x 30") | \$6,246.00 |
| AD2-SYS5-36 | ADAPT 2 SERIES SYSTEM 5 PKG. INCLUDES 2 ROHO MATTRESS SECTIONS (19" x 30") | \$4,260.00 |
| UPGRADE New! | SUBSTITUTE THE HIGH RESILIENCY FOAM INSERTS TO HYBRID FOAM | \$209.00 |











| System 1 AD2-SYS1-36 | System 2 AD2-SYS2-36 | System 3 AD2-SYS3-36 | System 4 AD2-SYS4-36 | System 5 AD2-SYS5-36 |
|---------------------------------------|--|---|--|--|
| Foam base covered with Recovery 5 | Foam base covered with Recovery 5 | Foam base covered with Recovery 5 | Foam base covered with Recovery 5 | Foam base covered with Recovery 5 |
| 1 x ROHO mattress section (19" x 30") | 2 x ROHO mattress sections (19" x 30") | 3 x ROHO mattress sections (19" x 30") | 4 x ROHO mattress sections (19" x 30") | 2 x ROHO mattress sections (19" x 30") |
| 1 x AD2-3030FS | 1 x AD2-3030FS | 1 x AD2-1930FS | | 1 x AD2-3030FS |
| 1 x AD2-1930FS | 1 x AD2-830FS | | | 1 x AD2-830FS |
| 1 x AD2-830FS | | | | |

Top cover options


| Part Number | Description | Retail Price |
|----------------|--|--------------|
| AD2-TC3680 | TOP COVER, ADAPT 2 36x80 (RECOVERY 5 FLUID-RESISTANCE, MULTI-WAY STRETCH, ZIPPERED) | \$403.00 |
| AD2-TC3680-R5U | TOP COVER, ADAPT 2 36x80 (RECOVERY 5 ULTRA SUPERIOR FLUID-RESISTANCE, MULTI-WAY STRETCH, ZIPPERED) | \$520.00 |

Adapt 2 Series modular system


39" x 80"

The Adapt 2 mattress with SRT is an innovative way to reduce the effects of pressure, friction and shear within the body of the mattress. The 3" wide firm perimeter borders and 4" permanent gel foam foot border hold the inserts in place without the need for snaps. The surface inserts can be configured based on the individual client's needs. The removable top cover with its 360° zipper must be ordered separately. Weight capacity: 350 lb

Note: Adapt 2 base and inserts include a fire barrier sock fabric to meet full-scale fire retardant standards (not pictured).



| Part Number | Description | Retail Price |
|--------------|--|--------------|
| AD2-SYS1-39 | ADAPT 2 SERIES SYSTEM 1 PKG. INCLUDES 1 ROHO MATTRESS SECTION (19" x 34") | \$3,562.00 |
| AD2-SYS2-39 | ADAPT 2 SERIES SYSTEM 2 PKG. INCLUDES 2 ROHO MATTRESS SECTIONS (19"x 34") | \$4,478.00 |
| AD2-SYS3-39 | ADAPT 2 SERIES SYSTEM 3 PKG. INCLUDES 3 ROHO MATTRESS SECTIONS (19" x 34") | \$5,451.00 |
| AD2-SYS4-39 | ADAPT 2 SERIES SYSTEM 4 PKG. INCLUDES 4 ROHO MATTRESS SECTIONS (19" x 34") | \$6,368.00 |
| AD2-SYS5-39 | ADAPT 2 SERIES SYSTEM 5 PKG. INCLUDES 2 ROHO MATTRESS SECTIONS (19" x 34") | \$4,478.00 |
| UPGRADE New! | SUBSTITUTE THE HIGH RESILIENCY FOAM INSERTS TO HYBRID FOAM | \$209.00 |











| System 1 AD2-SYS1-39 | System 2 AD2-SYS2-39 | System 3 AD2-SYS3-39 | System 4 AD2-SYS4-39 | System 5 AD2-SYS5-39 |
|--|--|---|--|---|
| Foam base covered with Recovery 5 | Foam base covered with Recovery 5 | Foam base covered with Recovery 5 | Foam base covered with Recovery 5 | Foam base covered with Recovery 5 |
| 1 x ROHO mattress section (19" x 34") | 2 x ROHO mattress sections (19" x 34") | 3 x ROHO mattress sections (19" x 34") | 4 x ROHO mattress sections (19" x 34") | 2 x ROHO mattress sections (19" x 34") |
| 1 x AD2-3033FS | 1 x AD2-3033FS | 1 x AD2-1933FS | | 1 x AD2-3033FS |
| 1 x AD2-1933FS | 1 x AD2-833FS | | | 1 x AD2-833FS |
| 1 x AD2-833FS | | | | |

Top cover options

| Part Number | Description | Retail Price |
|----------------|--|--------------|
| AD2-TC3980 | TOP COVER, ADAPT 2 39x80 (RECOVERY 5 FLUID-RESISTANCE, MULTI-WAY STRETCH, ZIPPERED) | \$480.00 |
| AD2-TC3980-R5U | TOP COVER, ADAPT 2 39x80 (RECOVERY 5 ULTRA SUPERIOR FLUID-RESISTANCE, MULTI-WAY STRETCH, ZIPPERED) | \$596.00 |

Adapt 2 Series modular system

42" x 80"

The Adapt 2 mattress with SRT is an innovative way to reduce the effects of pressure, friction and shear within the body of the mattress. The 4.5" wide firm perimeter borders and 4" permanent gel foam foot border hold the inserts in place without the need for snaps. The surface inserts can be configured based on the individual client's needs, and are the same size as the 39" insert options thanks to the addition of 1.5" of foam at each side border. The removable top cover with its 360° zipper must be ordered separately. Weight capacity: 350 lb

Note: Adapt 2 base and inserts include a fire barrier sock fabric to meet full-scale fire retardant standards (not pictured).




| Part Number | Description | Retail Price |
|--------------|--|--------------|
| AD2-SYS1-42 | ADAPT 2 SERIES SYSTEM 1 PKG. INCLUDES 1 ROHO MATTRESS SECTION (19" x 34") | \$3,793.00 |
| AD2-SYS2-42 | ADAPT 2 SERIES SYSTEM 2 PKG. INCLUDES 2 ROHO MATTRESS SECTIONS (19"x 34") | \$4,709.00 |
| AD2-SYS3-42 | ADAPT 2 SERIES SYSTEM 3 PKG. INCLUDES 3 ROHO MATTRESS SECTIONS (19" x 34") | \$5,682.00 |
| AD2-SYS4-42 | ADAPT 2 SERIES SYSTEM 4 PKG. INCLUDES 4 ROHO MATTRESS SECTIONS (19" x 34") | \$6,599.00 |
| AD2-SYS5-42 | ADAPT 2 SERIES SYSTEM 5 PKG. INCLUDES 2 ROHO MATTRESS SECTIONS (19" x 34") | \$4,709.00 |
| UPGRADE New! | SUBSTITUTE THE HIGH RESILIENCY FOAM INSERTS TO HYBRID FOAM | \$209.00 |











| System 1 AD2-SYS1-42 | System 2 AD2-SYS2-42 | System 3 AD2-SYS3-42 | System 4 AD2-SYS4-42 | System 5 AD2-SYS5-42 |
|--|---|---|--|---|
| Foam base covered with Recovery 5 | Foam base covered with Recovery 5 | Foam base covered with Recovery 5 | Foam base covered with Recovery 5 | Foam base covered with Recovery 5 |
| 1 x ROHO mattress section (19" x 34") | 2 x ROHO mattress sections (19" x 34") | 3 x ROHO mattress sections (19" x 34") | 4 x ROHO mattress sections (19" x 34") | 2 x ROHO mattress sections (19" x 34") |
| 1 x AD2-3033FS | 1 x AD2-3033FS | 1 x AD2-1933FS | | 1 x AD2-3033FS |
| 1 x AD2-1933FS | 1 x AD2-833FS | | | 1 x AD2-833FS |
| 1 x AD2-833FS | | | | |

Top cover options

| Part Number | Description | Retail Price |
|----------------|--|--------------|
| AD2-TC4280 | TOP COVER, ADAPT 2 42x80 (RECOVERY 5 FLUID-RESISTANCE, MULTI-WAY STRETCH, ZIPPERED) | \$573.00 |
| AD2-TC4280-R5U | TOP COVER, ADAPT 2 42x80 (RECOVERY 5 ULTRA SUPERIOR FLUID-RESISTANCE, MULTI-WAY STRETCH, ZIPPERED) | \$688.00 |

Adapt 2 Series mattress bases & components

Size: 36" / 39" / 42" x 80", can be customized to any size (call for a quote). Insert options: foam, hybrid foam, or ROHO® (foam & hybrid foam inserts are flippable for longevity and ease of use); each insert is encased in Recovery 5 material.

Top/Enclosure covers: **SOLD SEPARATELY.** The overall R5 top zippered cover provides extra incontinent protection. Or UPGRADE to the R5 Ultra top cover with 4-way stretch, is machine washable, and provides up to 6X more incontinent protection.

| Part Number | Description | Retail Price |
|-------------------------------|---|--------------|
| AD2-BASE3680 | 36" x 80" x 6" BASE COVERED WITH RECOVERY 5 | \$1,150.00 |
| Hybrid Foam Inserts | 30" Wide | |
| AD2-830HS | 8" x 30" HYBRID FOAM SECTION | \$272.00 |
| AD2-1930HS | 19" x 30" HYBRID FOAM SECTION | \$416.00 |
| AD2-3030HS | 30" x 30" HYBRID FOAM SECTION | \$520.00 |
| Foam Inserts | 30" Wide | |
| AD2-830FS | 8" x 30" FOAM SECTION | \$194.00 |
| AD2-1930FS | 19" x 30" FOAM SECTION | \$307.00 |
| AD2-3030FS | 30" x 30" FOAM SECTION | \$370.00 |
| ROHO Inserts | | |
| MATTSECT1610 | 19" x 30" ROHO NEOPRENE MATTRESS SECTION | \$1,171.00 |
| Covers | | |
| AD2-TC3680 | TOP COVER, ADAPT 2 36x80 (RECOVERY 5 FLUID-RESISTANCE, MULTI-WAY STRETCH, ZIPPERED) | \$403.00 |
| AD2-TC3680-R5U | TOP COVER, ADAPT 2 36x80 (RECOVERY 5 ULTRA SUPERIOR FLUID-RESISTANCE, MULTI-WAY STRETCH, ZIPPERED) | \$520.00 |
| AD2-BC3680 | BASE COVER, ADAPT 2 36x80 (RECOVERY 5 FLUID-RESISTANCE, MULTI-WAY STRETCH) | \$461.00 |
| COV-MATTSECT1610 | FULL ENCLOSURE COVER, ROHO MATTSECT 19x30 (RECOVERY 5 FLUID-RESISTANCE, MULTI-WAY STRETCH) | \$207.00 |
| Note: the same inserts fit bo | th the 39" and 42" width AD2 thanks to the 1.5" wider borders on the 42" | |
| AD2-BASE3980 | 39" x 80" x 6" BASE COVERED WITH RECOVERY 5 | \$1,271.00 |
| AD2-BASE4280 | 42" x 80" x 6" BASE COVERED WITH RECOVERY 5 | \$1,856.00 |
| Hybrid Foam Inserts | 33" Wide | |
| AD2-833HS | 8" x 33" HYBRID FOAM SECTION | \$324.00 |
| AD2-1933HS | 19" x 33" HYBRID FOAM SECTION | \$520.00 |
| AD2-3033HS | 30" x 33" HYBRID FOAM SECTION | \$596.00 |
| Foam Inserts | 33" Wide | |
| AD2-833FS | 8" x 33" FOAM SECTION | \$220.00 |
| AD2-1933FS | 19" x 33" FOAM SECTION | \$359.00 |
| AD2-3033FS | 30" x 33" FOAM SECTION | \$439.00 |
| ROHO Inserts | | |
| MATTSECT3400 | 19" x 34" ROHO NEOPRENE MATTRESS SECTION | \$1,171.00 |
| Covers | | |
| AD2-TC3980 | TOP COVER, ADAPT 2 39x80 (RECOVERY 5 FLUID-RESISTANCE, MULTI-WAY STRETCH, ZIPPERED) | \$480.00 |
| AD2-TC3980-R5U | TOP COVER, ADAPT 2 39x80 (RECOVERY 5 ULTRA SUPERIOR FLUID-RESISTANCE, MULTI-WAY STRETCH, ZIPPERED) | \$596.00 |
| AD2-BC3980 | BASE COVER, ADAPT 2 39x80 (RECOVERY 5 FLUID-RESISTANCE, MULTI-WAY STRETCH) | \$461.00 |
| AD2-TC4280 | TOP COVER, ADAPT 2 42x80 (RECOVERY 5 FLUID-RESISTANCE, MULTI-WAY STRETCH, ZIPPERED) | \$573.00 |
| AD2-TC4280-R5U | TOP COVER, ADAPT 2 42x80 (RECOVERY 5 ULTRA SUPERIOR FLUID-RESISTANCE, MULTI-WAY STRETCH, ZIPPERED) | \$688.00 |
| AD2-BC4280 | BASE COVER, ADAPT 2 42x80 (RECOVERY 5 FLUID-RESISTANCE, MULTI-WAY STRETCH) | \$577.00 |
| COV-MATTSECT3400 | FULL ENCLOSURE COVER, ROHO MATTSECT 19x34 (RECOVERY 5 FLUID-RESISTANCE, MULTI-WAY STRETCH) | \$207.00 |
| COV-MATTSECT3400-R5U | FULL ENCLOSURE COVER, ROHO MATTSECT 19x34 (RECOVERY 5 ULTRA SUPERIOR FLUID-RESISTANCE, MULTI-WAY STRETCH) | \$288.00 |

Adapt 2 Paediatric modular system

30" x 57" New option!

The Adapt 2 Paediatric mattress with optional ROHO® Isolette or Serene foam inserts is ideal for vulnerable children with high to very high risk for pressure injuries. The surface inserts can be configured based on the individual client's needs and can hold up to three (3) ROHO Isolette. This modified mattress fits exactly in a standard crib. The 1.5" wide firm perimeter borders and 4" permanent gel foam foot border hold the inserts in place without the need for snaps. The mattress base is protected by a Recovery 5 anti-microbial incontinent cover, as well as a Recovery 5 protective top cover with zipper. Weight capacity: 200 lb



| Description | Retail Price |
|---|--|
| ADAPT 2 SERIES PAEDIATRIC BASE WITH ZIPPERED TOP COVER | \$1,360.00 |
| 4" SERENE FOAM INSERT FOR PAEDIATRIC BASE | \$373.00 |
| 8" SERENE FOAM INSERT FOR PAEDIATRIC BASE | \$394.00 |
| 15" SERENE FOAM INSERT FOR PAEDIATRIC BASE | \$457.00 |
| ROHO NEONATAL/ISOLETTE (15" x 27") | \$1,885.00 |
| TOP COVER, ADAPT 2 PEDIATRIC (RECOVERY 5 FLUID-RESISTANCE, MULTI-WAY STRETCH, ZIPPERED) | \$289.00 |
| | ADAPT 2 SERIES PAEDIATRIC BASE WITH ZIPPERED TOP COVER 4" SERENE FOAM INSERT FOR PAEDIATRIC BASE 8" SERENE FOAM INSERT FOR PAEDIATRIC BASE 15" SERENE FOAM INSERT FOR PAEDIATRIC BASE ROHO NEONATAL/ISOLETTE (15" x 27") |

Possible configurations











Acute modular system

36" x 80"


The Acute mattress is our foam surface for the client at high risk of experiencing skin breakdown and may require the addition of ROHO® air components. There are multiple configuration system options available to best meet their needs.

Overall top coverlet must be ordered separately.


Weight capacity: 350 lb

Size: 36" x 80", can be customized to any size (call for a quote)

Insert options: foam or ROHO



| Part Number | Description | Retail Price |
|-------------|--|--------------|
| AC-SYS1-PKG | ACUTE SYSTEM 1 PKG. INCLUDES 1 ROHO MATTRESS SECTION (19" x 34") | \$2,310.00 |
| AC-SYS2-PKG | ACUTE SYSTEM 2 PKG. INCLUDES 2 ROHO MATTRESS SECTIONS (19" x 34") | \$3,639.00 |
| AC-SYS3-PKG | ACUTE SYSTEM 3 PKG. INCLUDES 2 ROHO MATTRESS SECTIONS (19" x 34"), MOVABLE DIVIDER | \$3,812.00 |
| AC-SYS4-PKG | ACUTE SYSTEM 4 PKG. INCLUDES 3 ROHO MATTRESS SECTIONS (19" x 34") | \$4,921.00 |
| AC-SYS5-PKG | ACUTE SYSTEM 5 PKG. INCLUDES 2 ROHO MATTRESS SECTIONS (19" x 34"), IMMOVABLE DIVIDER | \$3,650.00 |
| AC-SYS6-PKG | ACUTE SYSTEM 6 PKG. INCLUDES 2 ROHO MATTRESS SECTIONS (19" x 34"), T-SHAPE OPENING | \$3,968.00 |
| AC-SYS7-PKG | ACUTE SYSTEM 7 PKG. INCLUDES 4 ROHO MATTRESS SECTIONS (19" x 34") | \$4,026.00 |



| Acute system 1 | Acute system 2 | Acute system 3 | Acute system 4 | Acute system 5 | Acute system 6 | Acute system 7 |
|---------------------------------------|--|---|--|--|--|--|
| Foam base covered with Recovery 5 | Foam base covered with Recovery 5 | Foam base covered with Recovery 5, movable divider included | Foam base covered with Recovery 5 | Foam base covered with Recovery 5, immovable divider | Foam base covered with Recovery 5, T-shape opening | Foam base covered with Recovery 5 |
| 1 x ROHO mattress section (19" x 34") | 2 x ROHO mattress sections (19" x 34") | 2 x ROHO mattress sections (19" x 34") | 3 x ROHO mattress sections (19" x 34") | 2 x ROHO mattress sections (19" x 34") | 2 x ROHO mattress sections (19" x 34") | 4 x ROHO mattress sections (19" x 34") |

Top cover options

| Part Number | Description | Retail Price |
|------------------|---|--------------|
| AC-COVERLET9 | FITTED ELASTIC TOP COVER, 36x80x9 (RECOVERY 5 FLUID-RESISTANCE, MULTI-WAY STRETCH) | \$300.00 |
| AC-COVERLET9-R5U | FITTED ELASTIC TOP COVER, 36x80x9 (RECOVERY 5 ULTRA SUPERIOR FLUID-RESISTANCE, MULTI-WAY STRETCH) | \$369.00 |

Acute mattress bases & components

| Part Number | Description | Retail Price |
|----------------------|---|--------------|
| AC-SYS1 | ACUTE SYSTEM 1 - SECTIONAL BASE COVERED WITH RECOVERY 5 w/ 19" OPENING | \$1,039.00 |
| AC-SYS2 | ACUTE SYSTEM 2 - SECTIONAL BASE COVERED WITH RECOVERY 5 w/ 38" OPENING | \$1,098.00 |
| AC-SYS3 | C-SYS3 ACUTE SYSTEM 3 - SECTIONAL BASE COVERED WITH RECOVERY 5 w/ TWO 19" OPENINGS, MOVABLE DIVIDER | |
| AC-SYS4 | ACUTE SYSTEM 4 - SECTIONAL BASE COVERED WITH RECOVERY 5 w/ 57" OPENING | \$1,098.00 |
| AC-SYS5 | ACUTE SYSTEM 5 - SECTIONAL BASE COVERED WITH RECOVERY 5 w/ TWO 19" OPENINGS | \$1,098.00 |
| AC-SYS6 | ACUTE SYSTEM 6 - SECTIONAL BASE COVERED WITH RECOVERY 5 w/ TWO 19" OPENINGS (T-SHAPE) | \$1,418.00 |
| AC-SYS7 | ACUTE SYSTEM 7 - SECTIONAL BASE COVERED WITH RECOVERY 5 w/ 76" OPENING | \$972.00 |
| Inserts | | |
| AC-FI | ACUTE - FOAM INSERT 19" x 34" | \$265.00 |
| MATTSECT3400 | ROHO NEOPRENE MATTRESS SECTION 19" x 34" | \$1,171.00 |
| Covers | | |
| AC-COVERLET9 | FITTED ELASTIC TOP COVER, 36x80x9 (RECOVERY 5 FLUID-RESISTANCE, MULTI-WAY STRETCH) | \$300.00 |
| AC-COVERLET9-R5U | FITTED ELASTIC TOP COVER, 36x80x9 (RECOVERY 5 ULTRA SUPERIOR FLUID-RESISTANCE, MULTI-WAY STRETCH) | \$369.00 |
| AC-BC3680 | BASE COVER, ACUTE SYSTEM (RECOVERY 5 FLUID-RESISTANCE, MULTI-WAY STRETCH) *SPECIFY SYSTEM NUMBER | \$369.00 |
| COV-MATTSECT3400 | FULL ENCLOSURE COVER, ROHO MATTSECT 19x34 (RECOVERY 5 FLUID-RESISTANCE, MULTI-WAY STRETCH) | \$207.00 |
| COV-MATTSECT3400-R5U | FULL ENCLOSURE COVER, ROHO MATTSECT 19x34 (RECOVERY 5 ULTRA SUPERIOR FLUID-RESISTANCE, MULTI-WAY STRETCH) | \$288.00 |

ROHO mattress overlays



From left to right: Prodigy®, SOFFLEX 2®, and ROHO® Dry Floatation®

At a glance


- Provides skin/soft tissue protection
- Adjustable to client's shape and weight, and to changes in position over time
- Accommodates other positioning devices beneath the overlay
- Can be reused with multiple individuals
- Provides immersion and envelopment
- Sections can be easily cleaned and disinfected
- · Two-piece enclosure cover
- Maximum client weight of up to 300 lb for Prodigy, 330 lb for SOFFLEX 2, and no weight limit for ROHO Dry Floatation

These three mattress overlays are non-powered, air-filled, adjustable, zoned, reactive, full support surfaces that are intended to conform to an individual's shape to provide skin/soft tissue and deep tissue injury protection and an environment to facilitate wound healing. Flexible, interconnected air cells can be adjusted to match the individual's body contours and changing needs.

The first two are made of durable polyvinyl and comprise of three sections (27" long) to make an overlay that is 81" long, while the last one is made of neoprene and has four sections (19" long) for a 76" long overlay. As you compare the Prodigy, SOFFLEX 2, and ROHO DRY FLOATATION mattress overlays, you also see an increasing number of air cells.

With 720 air cells, the ROHO DRY FLOATATION neoprene mattress overlay allows for even greater conformity and better pressure distribution. Each cell is specially treated to minimize friction and resulting shear, and bariatric kit add-on sections are available to widen the overlay. Drainage holes in the base of the mattress overlay also allows for moisture to drain away from the body. The performance, efficiency and quality of this overlay system are unequaled in the industry.

Covers: The Prodigy & SOFFLEX 2 both use a two-piece zippered enclosure cover. It allows for the removal of the top portion for ease of cleaning if needed without removing the entire product from the bed mattress.






| Part Number | Description | Retail Price |
|----------------|---|--------------|
| 3800-1 | PRODIGY SINGLE SECTION (27" x 35"); NEED 3 FOR FULL-SIZE OVERLAY; 300 LB WT CAP. | \$198.00 |
| SFLX2EA | SOFFLEX 2 SINGLE SECTION (27" x 35"); NEED 3 FOR FULL-SIZE OVERLAY; 330 LB WT CAP. | \$732.00 |
| MATTSECT3400 | ROHO NEOPRENE MATTRESS SINGLE SECTION (19" x 34"); NEED 4 FOR FULL-SIZE OVERLAY; NO WT CAP. | \$1,171.00 |
| MATTSECTBARI | NEOPRENE SECTION W/ SNAPS FOR BARI ADD-ON (19" x 34"); XXX LB WT CAP. | \$1,226.00 |
| BARI310EA | BARIATRIC MATTRESS KIT, SINGLE SECTION (3 x 10 CELLS) | \$480.00 |
| BARI310SET | BARIATRIC MATTRESS KIT SYSTEM, FOUR SECTIONS (3 x 10 CELLS) | \$1,839.00 |
| BARI410EA | BARIATRIC MATTRESS KIT, SINGLE SECTION (4 x 10 CELLS) | \$480.00 |
| BARI410SET | BARIATRIC MATTRESS KIT SYSTEM, FOUR SECTIONS (4 x 10 CELLS) | \$1,839.00 |
| Covers | | |
| COV-MR2PCSTD-P | 2-PIECE ENCLOSURE COVER (36" x 81.75") FOR PRODIGY & SOFFLEX 2 | \$343.00 |
| 6162M | 2-PIECE ENCLOSURE INCONTINENT COVER (34" x 76") FOR ROHO NEOPRENE MATTRESS OVERLAY | \$343.00 |


Leveling pads

By using a combination of leveling pads and ROHO® air mattress sections to make an overlay, a high degree of pressure redistribution can be maintained. The pads are made with a high resiliency and stable foam base, and offer the utmost in comfort and support. Simply mix and match the four sections to best suit your needs! Weight capacity: 250 lb/piece


Below are three additional options using ROHO mattress sections.







| Option 1 | Option 2 | Option 3 |
|----------------------|---------------------|----------------------------------|
| fits 4 ROHO sections | fits 1 ROHO section | RP3480 (foam with precut centre) |
| 4 x MATTSECT3400 | 1 x MATTSECT3400 | fits 1 ROHO section |
| | 2 x RP3134 | 1 x MATTSECT3400 |



fits up to 4 ROHO sections

| Part Number | Description | Retail Price |
|----------------------|---|--------------|
| RP-2134 | FOAM LEVELING PAD 21" x 34" | \$368.00 |
| RP-3134 | FOAM LEVELING PAD 31" x 34" | \$357.00 |
| RP-3480 | FOAM LEVELING PAD 34" x 80" WITH PRECUT CENTRE | \$614.00 |
| MATTSECT3400 | ROHO NEOPRENE MATTRESS SECTION 19" x 34" | \$1,171.00 |
| Covers | | |
| AC-COVERLET3 | FITTED ELASTIC TOP COVER, 36x80x3 (RECOVERY 5 FLUID-RESISTANCE, MULTI-WAY STRETCH) | \$300.00 |
| AC-COVERLET3-R5U | FITTED ELASTIC TOP COVER, 36x80x3 (RECOVERY 5 ULTRA SUPERIOR FLUID-RESISTANCE, MULTI-WAY STRETCH) | \$369.00 |
| COV-MATTSECT3400 | FULL ENCLOSURE COVER, ROHO MATTSECT 19x34 (RECOVERY 5 FLUID-RESISTANCE, MULTI-WAY STRETCH) | \$207.00 |
| COV-MATTSECT3400-R5U | FULL ENCLOSURE COVER, ROHO MATTSECT 19x34 (RECOVERY 5 ULTRA SUPERIOR FLUID-RESISTANCE, MULTI-WAY STRETCH) | \$288.00 |

Foam accessories

| Part Number | Description | Retail Price |
|-------------|---|--------------|
| UL3680 | 3" OVERLAY BASE WITH NYLON COVER; 36" x 80" | \$346.00 |
| GAPF36 | $4"$ GAP FILLER TO EXTEND $36"$ WIDE FOAM MATTRESS FROM $80"$ LONG TO $84"$; $36"$ WIDE \times $6"$ HIGH; COMES IN FULL ENCLOSURE R5 COVER | \$550.00 |
| GAPF42 | 4" GAP FILLER TO EXTEND 42" WIDE FOAM MATTRESS FROM 80" LONG TO 84"; 42" WIDE x 6" HIGH; COMES IN FULL ENCLOSURE R5 COVER | \$607.00 |
| GAPF368 | 4" GAP FILLER TO EXTEND 36" WIDE AIR SURFACES FROM 80" LONG TO 84"; 36" WIDE x 8" HIGH; COMES IN FULL ENCLOSURE R5 COVER | \$550.00 |

Beds

Tröst & Tröst Plus

Adjustable long-term care bed





At a glance

- Standard and Plus sizes available
- Trendelenburg, Reverse Trendelenburg, and Comfort chair positioning
- 12-button hand pendant
- Low position: 6.75"; High position: 30"
- Flip down directional castor locks at head end, plus swivel locking castor
- · Centre step safety lock at foot end
- Mattress retainers and wall bumpers included
- Standard size: 36" x 80"; Plus size: 39" x 80" or 42" x 80" and length extension to 84"
- Weight capacity of up to 600 lb

The Tröst Bed is an electric, adjustable bed ideal for the long-term care setting and community. It offers a low height of 6.75" and a high height of 30" to meet all your caregiving and client needs. Tröst Bed has a 12-button hand pendant and optional footboard attendant controls, so customized comfort is an easy adjustment. The Trendelenburg (and Reverse Trendelenburg) feature comes standard and allows clients to lay flat on their back on a 15–30 degree incline with their feet elevated above the head. Available in a natural or wine colour head/footboard (sold separately, except when bought as a package), that are easy to use with quick release hardware.

Optional features:

- Plus model expandable to 39" and 42"
- Length extension from 80" to 84" (added at foot end only)
- · Embedded footboard attendant controls
- · Half side rails, interchangeable between head and foot sections
- Tall half side rails

- Assist bar (shown), interchangeable between head and foot sections
- Trapeze bar
- Battery backup
- Transport trolley
- Add an overbed table, see pg. 34



The **Tröst Bed standard** package includes: 36" head and footboard, and a set of half side rails. The **Tröst Plus Bariatric Bed** package options includes: 36", 39" or 42" head and footboard, and a set of half side rails. Please indicate your choice of natural or wine colour head/footboards.

| Part Number | Description | Retail Price |
|--------------------------|---|--------------|
| TROST-STD-HSR-PKG | TRÖST STANDARD BED w/ HALF SIDE RAILS, 36" WIDE HEAD/FOOTBOARD | \$3,575.00 |
| TROST-PLUS-HSR-PKG | TRÖST PLUS BARIATRIC BED w/ HALF SIDE RAILS, 36" WIDE HEAD/FOOTBOARD | \$3,975.00 |
| TROST-PLUS-HSR-PKG-39 | TRÖST PLUS BARIATRIC BED w/ HALF SIDE RAILS, 39" WIDE HEAD/FOOTBOARD | \$4,025.00 |
| TROST-PLUS-HSR-PKG-42 | TRÖST PLUS BARIATRIC BED w/ HALF SIDE RAILS, 42" WIDE HEAD/FOOTBOARD | \$4,075.00 |
| Components & Accessories | ; | |
| TROST-STD | TRÖST STANDARD BED, 36" x 80" (EXTENDABLE TO 84"), HEAD/FOOTBOARD SOLD SEP. | \$2,975.00 |
| TROST-PLUS | TRÖST PLUS BARIATRIC BED (EXTENDABLE TO 84" & EXPANDABLE TO 39" OR 42"), HEAD/FOOTBOARD SOLD SEP. | \$3,375.00 |
| TROST-HFB-NATURE | TRÖST BED 36" WIDE HEAD/FOOTBOARD - NATURE | \$275.00 |
| TROST-HFB-WINE | TRÖST BED 36" WIDE HEAD/FOOTBOARD - WINE | \$275.00 |
| TROST-HFB-NATURE39 | TRÖST BED 39" WIDE HEAD/FOOTBOARD - NATURE | \$325.00 |
| TROST-HFB-WINE39 | TRÖST BED 39" WIDE HEAD/FOOTBOARD - WINE | \$325.00 |
| TROST-HFB-NATURE42 | TRÖST BED 42" WIDE HEAD/FOOTBOARD - NATURE | \$385.00 |
| TROST-HFB-WINE42 | TRÖST BED 42" WIDE HEAD/FOOTBOARD - WINE | \$385.00 |
| TROST-SC | EMBEDDED STAFF CONTROL | \$365.00 |
| TROST-HSR | TRÖST BED HALF SIDE RAILS (HEAD OR FOOT, PAIR) | \$495.00 |
| TROST-THSR2 | TRÖST BED TALL HALF SIDE RAILS (HEAD OR FOOT, PAIR) | \$595.00 |
| TROST-AB | ASSIST BAR (HEAD OR FOOT, SINGLE) | \$155.00 |
| TROST-TT | TRANSPORT TROLLEY | \$675.00 |
| TROST-TB | TRÖST BED TRAPEZE BAR; 280 LB WT CAP. | \$495.00 |
| TROST-BB | BATTERY BACKUP | \$275.00 |
| TROST-EXT | TRÖST BED LENGTH EXTENSION (80"-84") | \$335.00 |

Hälsa & Hälsa Plus

Adjustable home care bed







At a glance

- · Standard and Plus sizes available
- Trendelenburg & Reverse Trendelenburg
- Dual posts on each side increases stability at all heights
- Low position: 9.875" (12.125" with large castor)
- High position: 25.625" (27.875" with large castor)
- Disassembles into four section for easy transport or storage
- Standard size: 36" x 80"; Plus size: four width options - 36", 39", 42", and 48"
- Weight capacity of up to 375 lb for standard and 500 lb for Plus/bariatric model

The Hälsa Bed is an electric, adjustable home care bed available in standard and Plus sizes and includes a 10-button hand pendant. The Trendelenburg (and Reverse Trendelenburg) feature comes standard and allows clients to lay flat on their back on a 15–30 degree incline with their feet elevated above the head. Additionally, it offers the capability of moving both the head and foot sections simultaneously, an upgraded large castor option, and easy disassembly for transportation and storage.

Hälsa & Hälsa Plus bed standard features:

- Natural or wine colour head/footboard
- · Hooks for incontinence bag
- 3 locking castors & 1 directional locking castor
- · Side and foot mattress retainers
- · Receptacle for trapeze bar
- · 2 receptacles for IV holder stems
- Virtually noise-free with low electrical consumption DC Linak motors

Optional features:

- Assist bar (shown) or half side rails, interchangeable between head and foot sections
- Full side rails
- · Trapeze bar
- Battery backup
- · Upgraded large castor option
- Add an overbed table, see pg. 34



The Hälsa Bed Half Side Rail (HSR) package includes: 36" head and footboard, and a set of half side rails. The Hälsa Bed Full Side Rail (FSR) package includes: 36" head and footboard, and a set of full size rails. The Hälsa Plus Bariatric Bed HSR package includes: 36", 39", 42" or 48" head and footboard, and a set of half side rails. Please indicate your choice of natural or wine colour head/footboards. Wine colour option only available in the 36" wide head/footboard.

| Part Number | Description | Retail Price |
|--------------------------|---|--------------|
| HALSA-STD-HSR-PKG | HÄLSA STD BED w/ HALF SIDE RAILS, 36" WIDE HEAD/FOOTBOARD; 375 LB WT CAP. | \$3,544.00 |
| HALSA-STD-FSR-PKG | HÄLSA STD BED w/ FULL SIDE RAILS, 36" WIDE HEAD/FOOTBOARD; 375 LB WT CAP. | \$3,544.00 |
| HALSA-PLUS-HSR-PKG | HÄLSA PLUS BED w/ HALF SIDE RAILS, 36" WIDE HEAD/FOOTBOARD; 500 LB WT CAP. | \$4,725.00 |
| HALSA-PLUS-HSR-PKG39 | HÄLSA PLUS BED w/ HALF SIDE RAILS, 39" WIDE HEAD/FOOTBOARD; 500 LB WT CAP. | \$4,778.00 |
| HALSA-PLUS-HSR-PKG42 | HÄLSA PLUS BED w/ HALF SIDE RAILS, 42" WIDE HEAD/FOOTBOARD; 500 LB WT CAP. | \$4,830.00 |
| HALSA-PLUS-HSR-PKG48 | HÄLSA PLUS BED w/ HALF SIDE RAILS, 48" WIDE HEAD/FOOTBOARD; 500 LB WT CAP. | \$4,883.00 |
| Components & Accessories | 3 | |
| HALSA-STD | HÄLSA STANDARD BED, 36" x 80", HEAD/FOOTBOARD INCLUDED | \$3,334.00 |
| HALSA-PLUS | HÄLSA PLUS BED (EXPANDABLE TO 39", 42" OR 48"), 36" HEAD/FOOTBOARD INCLUDED | \$4,515.00 |
| HALSA-HFBC-NATURE | HÄLSA BED 36" WIDE HEAD/FOOTBOARD - NATURE | \$352.00 |
| HALSA-HFBC-WINE | HÄLSA BED 36" WIDE HEAD/FOOTBOARD - WINE | \$348.00 |
| HALSA-HFBC-NATURE39 | HÄLSA BED 39" WIDE HEAD/FOOTBOARD - NATURE | \$405.00 |
| HALSA-HFBC-NATURE42 | HÄLSA BED 42" WIDE HEAD/FOOTBOARD - NATURE | \$457.00 |
| HALSA-HFBC-NATURE48 | HÄLSA BED 48" WIDE HEAD/FOOTBOARD - NATURE | \$510.00 |
| HALSA-BAT2 | BATTERY BACK UP | \$294.00 |
| HALSA-HSR | HÄLSA BED HALF SIDE RAILS (HEAD OR FOOT, PAIR) | \$520.00 |
| HALSA-FSR | HÄLSA BED FULL SIDE RAILS (PAIR) | \$573.00 |
| HALSA-AB | ASSIST BAR (HEAD OR FOOT, SINGLE) | \$163.00 |
| HALSA-TBB | HÄLSA BED TRAPEZE BAR BARIATRIC; 300 LB WT CAP. | \$520.00 |

Resten Hi-Lo

Adjustable elevating home comfort bed







At a glance

- Available in Twin XL and Queen with a height range of 16" - 30"
- Easy to adjust your position using a mobile app or the wireless remote control
- Programmable massage settings
- · Built-in USB charging ports, two per side
- Reduce back pressure using the zero-gravity function
- Head & foot articulation to raise your head and feet for a healthy night's sleep
- Four quiet motors to optimize your positioning (head, foot, lumbar, and hi-low)
- Sizes: 39" x 80" and 60" x 80"
- Weight capacity of 410-800 lb

Improve the way you sleep with an adjustable bed frame that can accommodate almost any sleeping position. Resten Hi-Lo elevating adjustable beds are remote-controlled and come with advanced features to help you personalize your sleep experience. Key features include lumbar support, programmable massage, LED underbed lighting, USB charging ports, endless head & foot positioning, mattress retainer to keep it in place as the bed adjusts, and individualized controls. Bed comes with legs and upgraded castors to suit needs.

Resten Hi-Lo bed standard features:

- Motorized lumbar support
- · Bed legs (shown right) and upgraded castors
- LED underbed lighting

Optional feature:

Fixed side rails

- · Memory settings for specific positions
- · Built-in USB charging ports, two per side

Add an overbed table, see pg. 34

Resten Luxury and Hybrid mattress:

- Multi-layer foam mattress, perfect for all types of sleepers
- Available in Twin XL and Queen to match the Resten Hi-Lo bed
- Luxury is 11" high, Hybrid is 12" high
- Hybrid has choice of plush or firm
- Removable zippered cover
- Compressed and rolled for ease of delivery
- Weight capacity is 350 lb for Twin XL and 600 lb for Queen



The **Resten Hi-Lo** bed is available in Twin XL and Queen sizes and includes castor upgrade. Side rails are sold separately. For couples wanting to adjust their own settings, put two Twin XL beds together with the free connector to make a King-sized bed!

| Part Number | Description | Retail Price |
|-----------------|--|--------------|
| RESTHL-TXL | RESTEN HI-LOW BED - TWIN XL w/ UPGRADED CASTORS, 39" x 80"; 410 LB WT CAP. | \$3,195.00 |
| RESTHL-QN New! | RESTEN HI-LOW BED - QUEEN w/ UPGRADED CASTORS, 60" x 80"; 800 LB WT CAP. | \$3,795.00 |
| RESTHL-SR | RESTEN HI-LOW - SIDE RAILS | \$279.00 |
| RESTLUX-TXL | LUXURY MATTRESS - TWIN XL; 350 LB WT CAP. | \$1,295.00 |
| RESTLUX-QN New! | LUXURY MATTRESS - QUEEN; 600 LB WT CAP. | \$1,695.00 |
| RESTHYP-TXL | HYBRID MATTRESS - TWIN XL (PLUSH); 350 LB WT CAP. | \$1,595.00 |
| RESTHYF-TXL | HYBRID MATTRESS - TWIN XL (FIRM); 350 LB WT CAP. | \$1,595.00 |
| RESTHYP-QN New! | HYBRID MATTRESS - QUEEN (PLUSH); 600 LB WT CAP. | \$1,995.00 |
| RESTHYF-QN New! | HYBRID MATTRESS - QUEEN (FIRM); 600 LB WT CAP. | \$1,995.00 |

Resten Select

Adjustable home comfort bed



At a glance

- Available in Twin XL, Double, and Queen
- Adjust your position using the wireless remote control with built-in storage pocket
- Programmable massage settings
- Built-in USB charging ports, two per side
- · Zero-gravity and anti-snore functions
- Head & foot articulation to raise your head and feet for a healthy night's sleep
- Sizes: 39" x 80", 54" x 75", and 60" x 80"
- Weight capacity of up to 600 lb

Adjustable bedframe that can accommodate almost any position to help with better sleep and decrease back pain. Key features include programmable massage, LED underbed lighting, USB charging ports, endless head & foot positioning, mattress retainer to keep it in place as the bed adjusts, and wireless controller. The fixed height legs are 12" high, but adjustable 12" legs can be added.

Resten Select bed standard features:

- LED underbed lighting
- · Memory settings for specific positions
- · Built-in USB charging ports, two per side


Optional feature:


- · Adjustable height legs
- · Add an overbed table, see pg. 34

Resten Luxury and Hybrid mattress:

- Multi-layer foam mattress, perfect for all types of sleepers
- Available in Twin XL, Double, and Queen to match the Resten Select bed
- Luxury is 11" high, Hybrid is 12" high
- Hybrid has choice of plush or firm







Resten Select bed comes in your choice of Twin XL, Double, or Queen. For couples wanting to adjust their own settings, put two Twin XL together to make a split King-sized bed! Connector comes standard for the two Twin XL

| Part Number | Description | Retail Price |
|------------------|--|--------------|
| RESTSEL-TXL | RESTEN SELECT BED - TWIN XL, 39" x 80"; 600 LB WT CAP. | \$1,995.00 |
| RESTSEL-DB New! | RESTEN SELECT BED - DOUBLE, 54" x 75"; 600 LB WT CAP. | \$2,195.00 |
| RESTSEL-QN | RESTEN SELECT BED - QUEEN, 60" x 80"; 600 LB WT CAP. | \$2,345.00 |
| RESTSEL-ADJLEG | RESTEN SELECT - 12" ADJUSTABLE HEIGHT LEGS | \$135.00 |
| RESTLUX-TXL | LUXURY MATTRESS - TWIN XL; 350 LB WT CAP. | \$1,295.00 |
| RESTLUX-DBL New! | LUXURY MATTRESS - DOUBLE; 500 LB WT CAP. | \$1,495.00 |
| RESTLUX-QN | LUXURY MATTRESS - QUEEN; 600 LB WT CAP. | \$1,695.00 |
| RESTHYP-TXL | HYBRID MATTRESS - TWIN XL (PLUSH); 350 LB WT CAP. | \$1,595.00 |
| RESTHYF-TXL | HYBRID MATTRESS - TWIN XL (FIRM); 350 LB WT CAP. | \$1,595.00 |
| RESTHYP-DBL New! | HYBRID MATTRESS - DOUBLE (PLUSH); 500 LB WT CAP. | \$1,795.00 |
| RESTHYF-DBL New! | HYBRID MATTRESS - DOUBLE (FIRM); 500 LB WT CAP. | \$1,795.00 |
| RESTHYP-QN | HYBRID MATTRESS - QUEEN (PLUSH); 600 LB WT CAP. | \$1,995.00 |
| RESTHYF-QN | HYBRID MATTRESS - QUEEN (FIRM); 600 LB WT CAP. | \$1,995.00 |

NOA Elite Riser





At a glance

- Low position of 7.5" and high position of 28"
- · Lighter than most high/low beds
- Can be easily moved at any height and in any position
- New safety lock option
- Four colour options for head/foot boards
- Standard size: 35" x 80"; as well as 42" width options
- Available with a shorter (76") and longer (84") sleep deck for increased comfort;
 Note: Extension options listed with Elite accessories on page 33
- Weight capacity of 600 lb

| Part Number | Description | Retail Price |
|--------------------|---|--------------|
| 4000 | NOA ELITE RISER BED, 35" x 76" | \$3,188.00 |
| 4001 | NOA ELITE RISER BED, 35" x 80" | \$3,188.00 |
| 4002 | NOA ELITE RISER BED, 35" x 84" | \$3,423.00 |
| 4005 | NOA ELITE TWIN BED, 42" x 80" | \$3,964.00 |
| 4006 | NOA ELITE TWIN BED, 42" x 84" | \$4,066.00 |
| 4050SEP | SAFETY LOCK KIT, FOR ELITE RISER SERIES | \$186.00 |
| Board Ontions - nl | ease select from these colour choices | |

| Board Options - please select from these colour choices XXX = LOK - LIGHT OAK, MAH - MAHOGANY, WAL - WALNUT, WCH - WILD CHERRY | | |
|--|--|----------|
| 4010XXX | NOA HEAD AND FOOT BOARDS MELAMINE - 36" | \$187.00 |
| 4012XXX | NOA HEAD AND FOOT BOARDS MELAMINE - 42" | \$337.00 |
| 4015XXX | NOA FOOTBOARD CONTROL KIT - 36" | \$581.00 |
| 4018XXX | NOA FOOTBOARD CONTROL KIT - 42" | \$722.00 |
| 4016 | HDPE FOOTBOARD CONTROL KIT - 36" (SAND COLOUR ONLY) | \$735.00 |
| 4022 | HDPE FOOTBOARD CONTROL KIT - 36" WITH COMFORT CHAIR, UNDERBED LIGHT (SAND COLOUR ONLY) | \$803.00 |
| 4025 | HDPE FOOTBOARD CONTROL KIT - 42" WITH COMFORT CHAIR, UNDERBED LIGHT (SAND COLOUR ONLY) | \$795.00 |
| 4021XXX | FOOTBOARD CONTROL KIT - 36" WITH COMFORT CHAIR, UNDERBED LIGHT | \$578.00 |
| 4024XXX | FOOTBOARD CONTROL KIT - 42" WITH COMFORT CHAIR, UNDERBED LIGHT | \$657.00 |
| | | |

Soft Touch Side Rails are designed to provide extra safety for users, as well as protection against bruises and abrasions. Meets HBSW/FDA guidelines on bedrail entrapment zones.



| Part Number | Description | Retail Price |
|-------------|--|--------------|
| 4026 | SOFT TOUCH SIDE RAIL SET FOR HEAD END | \$410.00 |
| 4028 | SOFT TOUCH SIDE RAIL SET FOR FOOT END | \$652.00 |
| 4034 | NOA ASSIST RAIL (FIXED) - AMBI COMPLETE | \$189.00 |
| 4035 | NOA ASSIST RAIL (FIXED) - AMBI COMPLETE TO FIT EXPANSION KIT | \$215.00 |
| 4073 | NOA ASSIST RAIL (FOLDING) - AMBI COMPLETE | \$292.00 |
| 4091 | NOA ASSIST RAIL (FOLDING) - AMBI COMPLETE TO FIT EXPANSION KIT | \$215.00 |

See additional NOA accessories on page 32

NOA Light




At a glance

- Low position of 8.5" and high position of 25"
- Fully adjustable
- Four motors provide end users with comfort positioning at the head and foot area, as well as height adjustment
- Four colour options for head/foot boards
- Standard size: 35" x 80" with expansion kit available to widen bed to 42"



• Weight capacity of 450 lb

| Part Number | Description | Retail Price |
|-------------|--|--------------|
| 4008 | NOA LIGHT BED, 35" x 80" | \$2,687.00 |
| 4075 | NOA LIGHT BED WITH BATTERY BACKUP, 35" x 80" | \$2,830.00 |
| 4051 | SAFETY LOCK KIT, FOR NOA LIGHT BED | \$186.00 |
| 4054 | HANDSET - 10 BUTTON | \$182.00 |
| 4708 | BATTERY BACKUP KIT, INCLUDES CONTROL BOX | \$590.00 |
| 4041 | EXPANSION KIT, FOR NOA LIGHT | \$443.00 |
| 4049 | 4" DECK EXTENSION FOR 35" NOA LIGHT | \$259.00 |
| | ease select from these colour choices AK, MAH - MAHOGANY, WAL - WALNUT, WCH - WILD CHERRY | |
| 4010XXX | NOA HEAD AND FOOT BOARDS MELAMINE - 36" | \$187.00 |
| 4012XXX | NOA HEAD AND FOOT BOARDS MELAMINE - 42" | \$337.00 |
| 4013 | HDPE HEAD AND FOOT BOARDS - 36" (SAND COLOUR ONLY) | \$510.00 |
| 4014TF | HDPE HEAD AND FOOT BOARDS - 42" (SAND COLOUR ONLY) | \$556.00 |
| 4020 | NOA FOOTBOARD CONTROL KIT, FOR NOA LIGHT BED | \$239.00 |
| | | (9.119) |





| Part Number | Description | Retail Price |
|-------------|--|--------------|
| 4026 | SOFT TOUCH SIDE RAIL SET FOR HEAD END | \$410.00 |
| 4034 | NOA ASSIST RAIL (FIXED) - AMBI COMPLETE | \$189.00 |
| 4035 | NOA ASSIST RAIL (FIXED) - AMBI COMPLETE TO FIT EXPANSION KIT | \$215.00 |
| 4073 | NOA ASSIST RAIL (FOLDING) - AMBI COMPLETE | \$292.00 |
| 4091 | NOA ASSIST RAIL (FOLDING) - AMBI COMPLETE TO FIT EXPANSION KIT | \$215.00 |

NOA Bariatric Elite



At a glance

- Low position of 9.5" and high position of 30"
- Designed for heavy work load and very sturdy throughout
- Four motors providing the high/low, head, and knee functions found on better resident room beds
- Standard size: 42" x 80"
- Weight capacity of 800 lb

| Part Number | Description | Retail Price |
|-------------|---|--------------|
| 4007 | NOA ELITE BARIATRIC BED PACKAGE, 42" x 80", INCLUDES HEAD/FOOTBOARDS, SET OF SOFT TOUCH SIDE RAILS, AND TRANSPORT TROLLEY | \$7,255.00 |
| 4019 | NOA FOOTBOARD CONTROL KIT, FOR BARIATRIC BED - 42" | \$785.00 |
| 4748 | HANDSET - 10 BUTTON BACKLIT NON-LOCKING | \$176.00 |



Soft Touch Side Rails are designed to provide extra safety for users, as well as protection against bruises and abrasions. Meets HBSW/FDA guidelines on bedrail entrapment zones.

| Part Number | Description | Retail Price |
|-------------|--|--------------|
| 4026 | SOFT TOUCH SIDE RAIL SET FOR HEAD END | \$410.00 |
| 4028 | SOFT TOUCH SIDE RAIL SET FOR FOOT END | \$652.00 |
| 4034 | NOA ASSIST RAIL (FIXED) - AMBI COMPLETE | \$189.00 |
| 4035 | NOA ASSIST RAIL (FIXED) - AMBI COMPLETE TO FIT EXPANSION KIT | \$215.00 |
| 4073 | NOA ASSIST RAIL (FOLDING) - AMBI COMPLETE | \$292.00 |
| 4091 | NOA ASSIST RAIL (FOLDING) - AMBI COMPLETE TO FIT EXPANSION KIT | \$215.00 |

NOA Elite accessories

| Part Number | Description | Retail Price |
|-------------|---|--------------|
| 4040 | NOA FULL OVERBED TRAPEZE SYSTEM | \$1,511.00 |
| 4042 | NOA EXPANSION KIT | \$443.00 |
| 4043 | NOA TRANSPORT TROLLEY | \$634.00 |
| 4045 | 4" DECK EXTENSION FOR 35" ELITE SERIES | \$259.00 |
| 4047 | 4" DECK EXTENSION FOR 39" ELITE SERIES | \$259.00 |
| 4048 | 4" DECK EXTENSION FOR 42" ELITE SERIES | \$259.00 |
| 4052 | UNDERBED LIGHT KIT (INCLUDES LIGHT & HANDSET) | \$325.00 |
| 4055 | HANDSET - 10 BUTTON LOCKING, TIMOTION | \$201.00 |
| 4057 | IN-RAIL CONTROL - 6 BUTTON, TIMOTION | \$182.00 |
| 4059 | BATTERY BACKUP KIT - TIMOTION | \$315.00 |



Overbed tables

Two new options for overbed tables! Each includes a rubberized border to reduce spillage, a powder-coated metal stand with adjustable height, and two locking castors.

Tilting table (far right) has a high-density laminate top that measures 16" x 24" and the pivot stops at 15°/30°/45°/90°. It comes fully assembled and the base is foldable for storage under the bed. Table has a max height of 37" and a minimum height of 25".

Composite table has an injection-molded top that measures 16" x 30" and includes a built-in cupholder. Automatic one-touch height adjustments to a max height of 42" and a minimum height of 30". Comes flat packed with tools and hardware for assembly.



| Part Number | Description | Retail Price |
|-------------|--|--------------|
| UNITYTILTOV | TILTING OVERBED TABLE WITH HIGH-DENSITY LAMINATE TOP | \$230.00 |
| UNITYCOMPOV | COMPOSITE OVERBED TABLE WITH CUPHOLDER | \$230.00 |

Halo Wing

The Halo Wing helps eliminate the four most dangerous entrapment zones creating a safe bed system environment. It can be used as a side rail replacement, for transferring or movement in bed, and preventing bed rollouts and falls. It exceeds the compressibility of mattress and 2 1/3" gap space along with the 60" exit point for Zones 2, 3, and 4.

Can be used on NOA and Tröst beds, and comes with a lifetime warranty.



| Part Number | Description | Retail Price |
|-------------|--------------------------------|--------------|
| 77701 | HALO SAFETY WING, SINGLE MOUNT | \$804.00 |
| 77702 | HALO SAFETY WING, DOUBLE MOUNT | \$1,000.00 |

Specifications & warranties

Weight capacity: 350 lb for standard sizes; 600 lb for bariatric; unless otherwise specified

Permobil mattresses & modular systems: 2 year limited warranty; unless otherwise specified

5 Minds Mobility mattresses/surfaces: 2 year limited warranty on foam mattresses; 3 year limited warranty on powered surfaces

ROHO® mattress sections & overlays: Prodigy (1 year), SOFFLEX 2 (18 months), mattress sections/overlay (2 years)

APEX surfaces: 3 year limited warranty, 1 year warranty on covers

Tröst and Hälsa bed: Bed frame (10 years), mechanical (4 years), electrical (2 years), wood products (1 year)

NOA beds: Welds (15 years), electrical & mechanical (4 years), side rails (2 years), wood products (1 year)

Resten beds: Bed frame (10 years), electrical (2 years)