

New product

# Dreamline STX2 Cushion

Learn more about the  
Dreamline STX2 cushion



The Dreamline STX2 cushion is a highly supportive, deep contoured, Gel Fluid cushion with abundant positioning options, for the less active user with mild to severe asymmetrical positioning requirements.

## Key Features

- **Adjustability.** With an array of positioning tools for on-the-spot adjustability and adaptability for various postures.
- **Stability.** A firm positioning base with a deep well gives high lateral stability accentuated with its deep contouring.
- **Pressure Protection.** Large dual Cluster Gel Fluid packs nestled in 3cm of Enduro foam provides effective protection from skin breakdown, all secured within the FLO2™ cover for extra protection during transfers.



## Unparalleled Positioning

The Dreamline STX2 cushion has a deep contour profile and incorporates advanced positioning options for symmetrical to severely asymmetrical postures. Pressure redistribution is effective through its patented Cluster Gel Fluid packs, which are neatly hidden in a unique secure zip-up compartment. A soft Comfort Wrap together with the FLO2™ fabric facilitates immersion and long-term comfort.



## Technical Specifications

<b>Weight Limit</b>	136kg (300lbs) - Up to 20"			
<b>Weight of Product</b>	3.8kg* (8.38lbs)			
<b>Width</b>	12" - 20" (30.5cm - 50.8cm)			
<b>Depth</b>	12" - 22" (30.5cm - 55.9cm)			
<b>Pelvic Well</b>	1.8" (4.5cm)			
<b>Positioning Base Height</b>	12"/14"	16"/18"	20"	22"
	4.5cm	5.5cm	6.0cm	6.5cm
<b>Wrap Thickness</b>	2.5cm	3.0cm	3.0cm	3.0cm
<b>Gel Packs</b>	2 x Gel Fluid Packs			
<b>Base Foam</b>	EVA Closed Cell Foam			
<b>Comfort Wrap</b>	Enduro Foam with FLO2™ Cover			
<b>Covers</b>	FLO2™ Incontinence and Outer Cover			
	3D Spacer Side Panel			
	Neoprene Front & Side Panel Non-Slip Rubberised Urethane Base Panel			
<b>Adjustable Positioning</b>	Pelvic Obliquity / Pelvic Rotation			
	Abduction/Adduction			
	Windswept			
	Hip Flexion			
<b>Warranty</b>	18 Months			

\*Weight Based on 18" x 18" (45.8cm x 45.8cm)

Product Code DX2C

