



ROHO® CONTOUR SELECT® Cushion

Designed for positioning and posture, the CONTOUR SELECT cushion keeps you centered in the middle of the cushion, with your pelvis back in the wheelchair. The ISOFLO® Memory Control locks air into four separate compartments, providing stability and positioning by minimizing side-to-side and front-to-back motion.



Pressure injury prevention and treatment



In-place individualization with ISOFLO®
Memory Control



Four independent compartments for
customized fit



Stability



Positioning

Dimensions

Width: 14.5 in. to 24.25 in. /
36.5 cm to 62.5 cm

Height: Approx. 2.25 in. to
4.25 in. / 5.5 cm to 10.5 cm

Range due to varying heights
of cells

Weight: Approx. 3.0 lbs. /
1.4 kg

Weight limit

No weight limit.
However, the cushion
must be properly sized to
the individual.

Limited warranty

Cushion: 36 months
Cover: 6 months



Our evidence shows the CONTOUR SELECT is ideal for people who:

- Currently have, have a history of or are at risk for skin/soft tissue breakdown and/or any category/stage pressure injury (including deep tissue injury)
- Need added positioning and stability
- Need leg alignment and pelvis stabilization
- Need additional relief under ischial tuberosities (ITs)