

# Take-Apart

The Take-Apart frame allows attendants to take the Fuze apart in seconds without tools. This feature reduces the weight of the heaviest piece of the wheelchair to 29 lbs or less. Perfect for transporting it in the trunk of a car or SUV. Available on Fuze™ T50, Fuze T50n and Fuze JR.



# Seating and Positioning

Permobil provides a wide range of seating and positioning products to address the unique needs of individuals using PDG wheelchairs.

Comfort and ROHO® offer a variety of cushions, back supports and accessories that are designed for pressure redistribution, stability, and comfort that can be tailored to meet the individual needs of users.



permobil

permobil.com | Ph 888.858.4422 | Fax 877.323.39093

MKTG004160B

# PDG Mobility™ Comparison matrix

permobil



												
	Fuze™ T50	Fuze T50n (NEW)	Fuze T50 Junior (peds)	Fuze M50* (NEW)	Stellar GL	Stellar GLT	Stellar Impact	Stellar HD	Bentley™	Bentley LT	Bentley LT-R (NEW)	Eclipse
Dynamic tilt range	Up to 50°	Up to 50°	Up to 50°	Up to 50°	-5° to 45°	Up to 45°	Up to 45°	Up to 45° (25° @ 600 lbs)	Up to 20°	Up to 20°	Up to 20°	None
Seat width	13" - 24"	16" to 24"	12" - 18"	15" - 24"	14" - 32"	12" - 24"	14" - 22"	22" - 32"	14" - 32"	12" - 24"	12" - 24"	20" - 48"
Seat depth	15" - 22"	15" - 22"	13" - 18"	15" - 22"	16" - 22"	16" - 22"	17" - 22"	18" - 22"	16" - 22"	16"- 22"	16" - 22"	16" - 24"
Seat to floor height	13" - 20"	16" - 20"	13" - 20"	13" - 19"	13" - 20"	13" - 20"	14" - 21"	16" - 18"	13" - 20"	13" - 20"	12" - 19"	14" - 20"
Overall width	+7.5" - 8.5"	+4.0"	+7.5" - 8.5"	+8.5"	+9.5 - 10.5"	+9.5 - 10.5"	+9.5 - 10.5"	+11"	+9.5" - 10.5"	+9.5" - 10.5"	+10.5"	+7.5"
Transport weight	44 lbs	44 lbs	42 lbs	52 lbs	51 lbs	34 lbs	64 lbs	62 lbs	51 lbs	33 lbs	34 lbs	44 lbs
User weight capacity	400 lbs	400 lbs	250 lbs	400 lbs	450 lbs	300 lbs	250 lbs	600 lbs	450 lbs	300 lbs	300 lbs	1000 lbs
Frame material	Aluminum	Aluminum	Aluminum	Aluminum	Steel	Aluminum	Steel	Steel	Steel	Aluminum	Aluminum	Steel
Power tilt option	Yes, up to 350 lbs	Yes, up to 350 lbs	No	No	Yes, up to 250 lbs	Yes	No	No	Yes, up to 250 lbs	Yes	Yes	No
Take-Apart	Yes	Yes	Yes	No	No	No	No	No	No	No	No	No
Hand or foot propulsion	Hand	No	Hand	Hand	Hand or Foot	Hand or Foot	Hand or Foot	Hand or Foot	Hand or Foot	Hand or Foot	Hand or Foot	Hand or Foot
Self tilt	No	No	No	No	Yes	Yes	Yes	Yes	Yes	Yes	Yes	No
Transit tiedowns	Yes, up to 300 lbs	Yes, up to 300 lbs	Yes	Yes, up to 300 lbs	Yes, up to 300 lbs	Yes	Yes	No	Yes, up to 300 lbs	Yes, up to 300 lbs	Yes, up to 300 lbs	Yes, up to 600 lbs

\*This product is exclusively available to our customers in Canada.