

Mototronik

Connect and ride



permobil

Set your own limits

Mototronik is a stylish and intuitive, frontmounted power assist device enabling you to travel over longer distances and different terrains, whilst conserving energy. **Simply connect and ride.**

- Cover longer distances at a higher speed of up to 20km/h*
- Front casters are lifted to make it easier to ride over challenging surfaces
- Easy on-screen configuration to adjust your preferred settings
- Convenient basket for storage
- Cruise control which can be engaged to safer downhill travel
- Range of up to 50km*

*dependent on configuration





Specifications

Max. speed	12 km/h
Range (battery 11.6 Ah)	up to 50 km*
Range (battery 5.8 Ah)	up to 25 km*
Max. user weight	150 kg
Obstacle climbing	50 mm
Length	400 mm
Width	540 mm
Height	900 mm
Weight (12" motor unit)	10 kg
Weight (14" motor unit)	13 kg

**With a charged battery, moderate speed, flat course, and user weighing 75 kg.*