

# Kombinationsavtal med Panthera rullstolar

	Panthera S3	Panthera U3	Panthera S3 Swing	Panthera U3 Light	Panthera X	Panthera X3	Panthera Bambino 3	Panthera Micro 3
Alber E-fix E35/E36	●	●	●					
Alber E-motion M15/M25	●	●	●					
Alber E-Pilot P15	●	●		●				
Alber Smoov O10*	●	●		●	●			
Alber Twion T24	●	●	●					
Anatomic Sitt Höftbälten	●	●	●				●	●
Anatomic Sitt Vristselar	●	●	●					
Anpassarna Robo Inn Attachments	●	●			●			
Bodypoint, bälten och bröstselar/bröstvåstar	●	●	●				●	●
Decon Wheel Adventus**	●	●	●	●	●			
FreeWheel Panthera**	●	●		●	●		●	
Gear Wheel	●	●	●					
Permobil Mototronik*****	●	●	●	●	●	●		
Permobil SmartDrive	●	●	●	●	●	●		
V-Trak ryggsystem***	●	●	●	●	●		●	
Rehasense HK Ltd PAWS City	●	●		●				
Rehasense HK Ltd PAWS Cruiser	●	●		●				
Rehasense HK Ltd PAWS Tourer	●	●		●				
Stricker Handbikes	●	●		●				
Sunrise Medical Jay J3 rygghöjd LT, MT, UT****	●	●		●			●	
Sunrise Medical Jay J3 Carbon rygg	●	●		●			●	
Sunrise Medical Attitude påhängscyklar	●	●		●				
Sunrise Medical Wheeldrive	●	●		●				

\*För Panthera S3 Swing individuell bedömning av funktionskontroll    \*\*Med förstärkt fotbåge på Panthera-rullstolar    \*\*\*Panthera S3 – för fällbara fotplattor funkar enbart bredd 39, 42 och 45, då med asymmetriska fotplattor

\*\*\*\*V-Trak på Panthera X endast i kombination med U3 Light rygg för fasta ryggsystem    \*\*\*\*\*UT ej godkänt på Light-rullstolar, se mont. anvisning