

Joystick Module

With Bluetooth



permobil

Features and benefits

The NEW Permobil Joystick Module (PJSM) unlocks the power of Bluetooth® wireless technology by allowing you to connect and control up to 4 supported devices, two of which can be iDevices®. The PJSM also allows you to take advantage of your iOS SwitchControl & VoiceOver functionality. All connected devices are listed in the universal Bluetooth mode and connected-device switching is possible without unpairing.

The PJSM also offers on screen mouse control of two applicable devices as well as built in Infra-Red [IR] Control with device command memory. Extend the functionality of the PJSM by connecting switches to the Profile/Mode Jack and programming them for click/navigation functions. This unprecedented amount of control is now standard on the Permobil M & F series Corpus/Corpus VS models and the M300 Corpus HD!

PLUS all the GREAT features from the original PJSM

- Compact, ergonomic design
- Durable metal housing
- Built-in light sensor to automatically adjust screen brightness
- Paddle switches for easy access to user controls
- Enhanced programming and customization options
- Speaker for enhanced horn volume



Programmable buttons

All buttons on the new PJSM can be reassigned for more customization. The soft keys are preprogrammed for lights and key functions with unique presets for standing/anterior tilt. A short press activates function 1 while a long press will activate

function 2. These key functions can adapt to mode. For example, when in Bluetooth mode, the keys are available to operate mouse clicks and iDevice® commands. Lastly, the Mode/Profile jack is preprogrammed to run IRM (long press).

