

Customer and Order Information			
Purchase Order *		Tag	v.0910
Dealer/Provider Name *			
Purchaser Name and Contact Information *			
Shipping Address *			
User Information (if available)		Weight _____ lbs	Height _____ ft/in
			*required
Order Notes/Instructions			

Base Chair & Frame Options			
☒ 2501000	Stellar GLT 45° Manual Tilt Wheelchair		W321 \$ 2,685
☒ 2501009	Lightweight Aluminum Frame Upgrade Required Option for the Stellar GLT.		\$ 300
Frame Type All frames offer 300 lbs Weight Capacity.			
☐ 2502010	Fixed Width Frame Select a Fixed Seat Width below.		
☐ 2502025	Adjustable Width Frame Select an Adjustable Seat Width below. \$ 200		
User Weight Range Required for setup of tilt uplift force.			
☐ 2502431	75 lbs to 100 lbs	☐ 2502432 100 lbs to 115 lbs	☐ 2502433 115 lbs to 150 lbs
☐ 2502434	150 lbs to 200 lbs	☐ 2502435 200 lbs to 250 lbs	☐ 2502436 250 lbs to 300 lbs
☐ 2502300	Transit Tie Down System (TTDS) TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details.		
PDG Premium Colours			
☐ 2505001	Black Shimmer	☐ 2505015 Sky Blue	
☐	Other Standard PDG Premium Colour Specify Colour _____		pdgmobility.com/colours \$ 125

Seat Width & Depth		Overall chair width is +9.5" for 12/16" RW and +10.5" for 20/22/24" RW	
Fixed Seat Widths Order with Fixed Frame.			
☐ 2503015	15" Seat Width	☐ 2503016 16" Seat Width	☐ 2503017 17" Seat Width
☐ 2503018	18" Seat Width	☐ 2503019 19" Seat Width	
☐ 2503020	20" Seat Width	☐ 2503021 21" Seat Width	☐ 2503022 22" Seat Width W111 \$ 250
☐ 2503023	23" Seat Width	☐ 2503024 24" Seat Width	W111 \$ 500
Adjustable Seat Widths Order with Adjustable Frame.			
☐ 2503215	15"-20" Range - Set at 15"	☐ 2503216 15"-20" Range - Set at 16"	☐ 2503217 15"-20" Range - Set at 17"
☐ 2503218	15"-20" Range - Set at 18"	☐ 2503219 15"-20" Range - Set at 19"	
☐ 2503220	15"-20" Range - Set at 20"		W111 \$ 250
Seat Depths			
☐ 2504016	16" Seat Depth	☐ 2504017 17" Seat Depth	☐ 2504018 18" Seat Depth
☐ 2504019	19" Seat Depth		
☐ 2504020	20" Seat Depth	☐ 2504021 21" Seat Depth	☐ 2504022 22" Seat Depth W111 \$ 250
☒ 2508001	Solid Aluminum Seat Pan		

Armrests			
Single Post Height Adjustable Armrests			
☐ 2506111	8"-11" Range	☐ 2506114 10"-14" Range	☐ 2506116 12"-16" Range
Regular Armpads			
☐ 2507100	Desk Length - 2" x 10"	☐ 2507101 Full Length - 2" x 14"	
☐ 2507310	Clothing Guards		\$ 125

Seat Height

➤ Don't forget

Select desired Seat-to-Floor Height

STF Height is adjustable within the range of the chosen Subframe Plate (Low/High).

6" Casters	Low							High						
	RW	13"	14"	15"	16"	17"	18"	RW	15"	16"	17"	18"	19"	20"
	12"	☐	☐	☐	☐	☐	☐	12"	☐	☐	☐	☐	☐	☐
	16"	☐	☐	☐	☐	☐	☐	16"	☐	☐	☐	☐	☐	☐
	20"	☐	☐	☐	☐	☐	☐	20"	☐	☐	☐	☐	☐	☐
	22"	☐	☐	☐	☐	☐	☐	22"	☐	☐	☐	☐	☐	☐
	24"	☐	☐	☐	☐	☐	☐	24"	☐	☐	☐	☐	☐	☐

8" Casters	Low							High						
	RW	13"	14"	15"	16"	17"	18"	RW	15"	16"	17"	18"	19"	20"
	16"	☐	☐	☐	☐	☐	☐	16"	☐	☐	☐	☐	☐	☐
24"	☐	☐	☐	☐	☐	☐	☐	24"	☐	☐	☐	☐	☐	☐

Casters & Rear Wheels

Urethane Casters

☐ 2511061 6" x 1.25" Caster size

☐ 2511081 8" x 1.25" Caster size

2514000 Composite Wheels with Urethane Tires

☐ 2512012 12" Wheel size

☐ 2512016 16" Wheel size

☐ 2512020 20" Wheel size

☐ 2512022 22" Wheel size

☐ 2512024 24" Wheel size

☒ 2513000 Quick Release Rear Axles

Handdrims

Handdrims are only available with 20", 22", and 24" Rear Wheels.

☐ 2515010 Standard Aluminum

☐ 2515000 No Handdrims

☐ 2515020 Plastic Coated

\$

175

Wheel Locks

☒ 2516010 Push to Lock Wheel Locks

☐ 2516100 Wheel Lock Extensions

Anti-Tippers

☒ 2517010 Rear Anti-Tippers

Required Safety Option.

Back Post

☒ 2518010 2519010 Angle/Height Adjustable Stroller Handle

Button on Outside.

Back Post Height

Measured from the Seat Pan to the Handle.

☐ 2518425 Short (22")

☐ 2518431 Tall (25")

☒ 2519090 90° Seat to Back Angle

☒ 2519220 Fold Down Back

Tilt/Recline Options

☒ 2523011 Attendant Tilt only

Trigger on Right Back Post.

Front Rigging

Plug In Front Rigging

☐ 2525410 Easy-Out Fixed 70°

Effortless insertion and removal. No Hanger Bracket remains on frame.

\$

95

☐ 2525110 Standard Plug In 70°

Legacy Design. Not compatible with Easy-Out.

Swing Away Front Rigging

☐ 2525140 Swing Away 70°

☐ 2525200 Swing Away Elevating Legrests

Includes Calf Pads.

W106

\$

375

☐ 2525300 1" Offset Hangers for Swing Away

\$

95

Permobil | 318 East Kent Ave South | Vancouver BC V5X 4N6 Canada | Phone 888 858 4422 | Email sales.pdg@permobil.com

2 / 3

Footplates & Options					
Footplates					
☐	2526010	Composite			
☐	2526020	Aluminum Angle Adjustable			W107 \$ 165
☒	2527040	Length Adjustable Heel Loops			
Front Rigging Extension Tubes					
☐	2528010	Short	14"-18" for Easy-Out Front Rigging 13"-17" for Swing Away and Standard Plug In Front Rigging		
☐	2528020	Long	18"-24" for Easy-Out Front Rigging 15"-19" for Swing Away and Standard Plug In Front Rigging		\$ 75

Please send your completed order form by email to sales.pdg@permobil.com.